



Recombinant Dictyostelium discoideum Serine/threonine-protein phosphatase 4 catalytic subunit (ppp4c)

Product Code	CSB-BP897511DKK
Abbreviation	ppp4c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y0B7
Product Type	Recombinant Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Purity	>85% (SDS-PAGE)
Sequence	MSSDLDRQIE QLKRCEIIKE SEVRALCSKA REILLEEGNV QRVDSPVTIC GDIHGQFYDL KELFKVGGDC PQTNYLFMGD FVDRGFYSVE TFLLLLALKV RYPDRITLIR GNHESRQITQ VYGFYEECVR KYGSVTVWKY CTEIFDYLSL SALVDGKIFC VHGGGLSPSIN TLDQIRAIDR KQEVPHGPM CDLMWSDPED IPGWNGSPRG AGFLFGEDVV QKFNDNMLE FICRAHQLVM EGFKYMNET LVTVWSAPNY CYRCGNVAAI LQLDENLKKK FAIFEAAPQE SRGAPAKKPA PEYFL
Source	Baculovirus
Target Names	ppp4c
Protein Names	Recommended name: Serine/threonine-protein phosphatase 4 catalytic subunit Short name= PP4C EC= 3.1.3.16
Expression Region	1-305
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.