



Recombinant Human Protein FAM50B (FAM50B)

Product Code	CSB-EP897521HU-B
Abbreviation	FAM50B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y247
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AQYKGTMRE AGRAMHLLKK RERQREQMEV LKQRIAEETI LKSQVDKRFS AHYDAVEAEL KSSTVGLVTL NDMKARQEAL VRERERQLAK RQHLEEQRLQ QERQREQEQR RERKRKISCL SFALDDLDDQ ADAAEARRAG NLGKNPDVDT SFLPDRDREE EENRLREELR QEWEAQREKV KDEEMEVTFS YWDGSGHRRT VRVRKGNTVQ QFLKKALQGL RKDFLELRSA GVEQLMFIKE DLILPHYHTF YDFIARARG KSGPLFSFDV HDDVRLLSDA TMEKDESHAG KVVLRSWYEK NKHIFPASRW EAYDPEKKWD KYTIR
Source	E.coli
Target Names	FAM50B
Protein Names	Recommended name: Protein FAM50B Alternative name(s): Protein XAP-5-like
Expression Region	2-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.