



Recombinant Human Afadin- and alpha-actinin-binding protein (SSX2IP)

Product Code	CSB-BP897526HU
Abbreviation	SSX2IP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y2D8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGDWMTVTDP GLSSESKTIS QYTSETKMSP SSLYSQQVLC SSIPLSKNVH SFFSAFCTED NIEQSISYLD QELTTFGFPS LYEESKGGKET KRELNIVAVL NCMNELLVLQ RKNLLAQENV ETQNLKLGSD MDHLQSCYSK LKEQLETSRR EMIGLQERDR QLQCKNRNLH QLLKNEKDEV QKLQNIASR ATQYNHDMKR KEREYNKLKE RLHQLVMNKK DKKIAMDILN YVGRADGKRG SWRTGKTEAR NEDEMYKILL NDYEYRQKQI LMENAELKKV LQQMKKEMIS LLSPQKKKPR ERVDDSTGTV ISDVEEDAGE LSRESMWDLS CETVREQLTN SIRKQWRILK SHVEKLDNQV SKVHLEGFND EDVISRQDHE QETEKLELEI QQCKEMIKTQ QQLLQQQLAT AYDDDTTSLR RDCYLLEEKE RLKEEWSLFK EQKKNFERER RSFTEAAIRL GLERKAFEEE RASWLKQQFL NMTTFDHQNS ENVKLFSAFS GSSDWDNLIV HSRQPQKKPH SVSNGSPVCM SKLTKSLPAS PSTSDFCQTR SCISEHSSIN VLNITAEIEK PNQVGGECTN QKWSVASRPG SQEGCYSGCS LSYTNSHVEK DDLP
Source	Baculovirus
Target Names	SSX2IP
Protein Names	Recommended name: Afadin- and alpha-actinin-binding protein Short name= ADIP Alternative name(s): Afadin DIL domain-interacting protein SSX2-interacting protein
Expression Region	1-614
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a protein that binds the cancer-testis antigen Synovial Sarcoma X breakpoint 2 protein. The encoded protein may regulate the activity of Synovial Sarcoma X breakpoint 2 protein in malignant cells. Alternate splicing results in multiple transcript variants. A pseudogene of this gene is found on



chromosome 3.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.