



# Recombinant Human Zinc finger and BTB domain-containing protein 32 (ZBTB32)

<b>Product Code</b>	CSB-BP897538HU
<b>Abbreviation</b>	ZBTB32
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9Y2Y4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSLPPIRLPS PYGSDRLVQL AARLRPALCD TLITVGSQEF PAHSLVLAGV SQQLGRRGQW ALGEGISPST FAQLLNFBVYG ESVELQPGEL RPLQEAARAL GVQSLEEACW RARGDRAKPP DPGLKKHQEE PEKPSRNPER ELGDPGEKQK PEQVSRTGGR EQEMLHKHSP PRGRPEMAGA TQEAQQEQTR SKEKRLQAPV GQRGADGKHG VLTWLRENPG GSEESLRKLP GPLPPAGSLQ TSVTPRPSWA EAPWLVGGQP ALWSILLMPP RYGIPFYHST PTTGAWQEVW REQRIPLSLN APKGLWSQNN LASSSPTPGS LPQGPAQLSP GEMEESDQGH TGALATCAGH EDKAGCPPRP HPPPAPPARS RPYACSVCGK RFSLKHQMET HYRVHTGEKP FSCSLCPQRS RDFSAMTKHL RTHGAAPYRC SLCGAGCPSL ASMQAHMRGH SPSQLPPGWT IRSTFLYSSS RPSRPSTSPC CPSSSTT
<b>Source</b>	Baculovirus
<b>Target Names</b>	ZBTB32
<b>Protein Names</b>	Recommended name: Zinc finger and BTB domain-containing protein 32 Alternative name(s): FANCC-interacting protein Fanconi anemia zinc finger protein Testis zinc finger protein Zinc finger protein 538
<b>Expression Region</b>	1-487
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.