



Recombinant Human Zinc finger and BTB domain-containing protein 32 (ZBTB32)

Product Code	CSB-MP897538HU
Abbreviation	ZBTB32
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y2Y4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSLPPIRLPS PYGSDRLVQL AARLRPALCD TLITVGSQEF PAHSLVLAGV SQQLGRRGQW ALGEGISPST FAQLLNFBVYG ESVELQPGEL RPLQEAARAL GVQSLEEACW RARGDRAKKP DPGLKKHQEE PEKPSRNPER ELGDPGEKQK PEQVSRTGGR EQEMLHKHSP PRGRPEMAGA TQEAQQEQTR SKEKRLQAPV GQRGADGKHG VLTWLRENPG GSEESLRKLP GPLPPAGSLQ TSVTPRPSWA EAPWLVGGQP ALWSILLMPP RYGIPFYHST PTTGAWQEVW REQRIPLSLN APKGLWSQNN LASSSPTPGS LPQGPAQLSP GEMEESDQGH TGALATCAGH EDKAGCPPRP HPPPAPPARS RPYACSVCGK RFSLKHQMET HYRVHTGEKP FSCSLCPQRS RDFSAMTKHL RTHGAAPYRC SLCGAGCPSL ASMQAHMRGH SPSQLPPGWT IRSTFLYSSS RPSRPSTSPC CPSSSTT
Source	Mammalian cell
Target Names	ZBTB32
Protein Names	Recommended name: Zinc finger and BTB domain-containing protein 32 Alternative name(s): FANCC-interacting protein Fanconi anemia zinc finger protein Testis zinc finger protein Zinc finger protein 538
Expression Region	1-487
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.