



Recombinant Human DNA-directed RNA polymerase III subunit RPC8 (POLR3H)

Product Code	CSB-MP897557HU
Abbreviation	POLR3H
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y535
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MFVLVEMVDT VRIPPWQFER KLNDSIAEEL NKKLANKVVY NVGLCICLFD ITKLEDAYVF PGDGASHTKV HFRCVVFHPF LDEILIGKIK GCSPEGVHVS LGFFDDILIP PESLQQPAKF DEAEQVWVWE YETEEGAHDL YMDTGEEIRF RVVDESFDVT SPTGPSSADA TTSSEELPKK EAPYTLVGS I SEPGLGLLSW WTSN
Source	Mammalian cell
Target Names	POLR3H
Protein Names	Recommended name: DNA-directed RNA polymerase III subunit RPC8 Short name= RNA polymerase III subunit C8 Alternative name(s): DNA-directed RNA polymerase III subunit H RNA polymerase III subunit 22.9 kDa subunit Short name= RPC22
Expression Region	1-204
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.