



# Recombinant Human Host cell factor 2 (HCFC2)

<b>Product Code</b>	CSB-EP897575HU
<b>Abbreviation</b>	HCFC2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9Y5Z7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAAPSLLNWR RVSSFTGPVP RARHGHRVA IRELMIFGG GNEGIADDELH VYNTATNQWF LPAVRGDIPP GCAAHGFVCD GTRILVFGGM VEYGRYSNEL YELQASRWLW KKVKPHPPPS GLPPCPRLGH SFSLYGNKCY LFGGLANESE DSNNNVPRYL NDFYELELQH GSGVVGWSIP VTKGVVPSPR ESHTAVIYCK KDSGSPKMYV FGGMCGARLD DLWQLDLETM SWSKPETKGT VPLPRSLHTA SVIGNKMYIF GGWVPHKGEN TETSPHDCEW RCTSSFSYLN LDTTEWTTLV SDSQEDKKNS RPRPRAGHCA VAIGTRLYFW SGRDGYKKAL NSQVCCKDLW YLDTEKPPAP SQVQLIKATT NSFHVKWDEV STVEGYLLQL STDLPYQAAS SDSSAAPNMQ GVRMDPHRQG SNNIVPNSIN DTINSTKTEQ PATKETSMTKN KPDFKALTDN NAILYPSLAS NASNHNSHV DMLRKNEGPH TSANVGVLS CLDVRTVIPE TSVSSTVSST QTMVTQQTIK TESSSTNGAV VKDETSLTTF STKSEVDETY ALPATKISR ETHATATPFS KETPSNPVAT VKAGERQWCD VGIFKNNTAL VSQFYLLPKG KQSISKVGNA DVPDYSLLKK QDLVPGTGYR FRVAANGCG IGPFSKISEF KTCIPGFPGA PSAVRISKNV EGIHLSWEPP TSPSGNILEY SAYLAIRTAQ IQDNPSQLVF MRIYCGLKTS CIVTAGQLAN AHIDYTSRPA IVFRISAKNE KGYGPATQVR WLQGNKKAP LN
<b>Source</b>	E.coli
<b>Target Names</b>	HCFC2
<b>Protein Names</b>	Recommended name: Host cell factor 2 Short name= HCF-2 Alternative name(s): C2 factor
<b>Expression Region</b>	1-792
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes one of two proteins which interact with VP16, a herpes simplex virus protein that initiates virus infection. Both the encoded protein and the original Herpes host cell factor interact with VP16 through a beta-propeller



domain. The original Herpes host cell factor, however, is effective at initiating viral infection while the encoded protein is not. Transcripts of varying length due to alternative polyadenylation signals have been described.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.