



Recombinant Human Transcription factor ETV7 (ETV7)

Product Code	CSB-EP897576HU-B
Abbreviation	ETV7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y603
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MQEGELAISP ISPVAAMPPL GTHVQARCEA QINLLGEGGI CKLPGRRLRIQ PALWSREDVL HWLRWAEQEY SLPCTAEHGF EMNGRALCIL TKDDFRHRAP SSGDVLYELL QYIKTQRRAL VCGPFFGGIF RLKTPTQHSP VPPEEVTGPS QMDTRRGHLL QPPDPGLTSN FGHLD DPGLA RWTPGKEESL NLCHCAELGC RTQGVCSFPA MPQAPIDGRI ADCRLLWDYV YQLLLDTRYE PYIKWEDKDA KIFRVVDPNG LARLWGNHKN RVNMTYEKMS RALRHYYKLN IIKKEPGQKL LFRFLKTPGK MVQDKHSHLE PLESQEQDRI EFKDKRPEIS P
Source	E.coli
Target Names	ETV7
Protein Names	Recommended name: Transcription factor ETV7 Alternative name(s): ETS translocation variant 7 ETS-related protein Tel2 Tel-related Ets factor Transcription factor Tel-2
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.