



Recombinant Human Testisin (PRSS21)

Product Code	CSB-BP897592HU
Abbreviation	PRSS21
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y6M0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	IVGGEDAEL GRWPWQGSLR LWDSHVCGVS LLSHRWALTA AHCFETYS DL SDPSGWMVQF GQLTSMPSFW SLQAYYTRYF VSNIYLSPRY LGNSPYDIAL VKLSAPVTYT KHIQPICLQA STFEFENRTD CWVTGWGYIK EDEALPSPH T LQEVQVAIIN NSMCNHLFLK YSFRKDIFGD MVCAGNAQGG KDACFGDSSG G PLACNKNGLW YQIGVVSXWGV GCGRPNRPGV YTNISHHFEW IQKLMAQS
Source	Baculovirus
Target Names	PRSS21
Protein Names	Recommended name: Testisin EC= 3.4.21.- Alternative name(s): Eosinophil serine protease 1 Short name= ESP-1 Serine protease 21
Expression Region	42-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a cell-surface anchored serine protease, which is a member of the trypsin family of serine proteases. It is predicted to be active on peptide linkages involving the carboxyl group of lysine or arginine. The protein localizes to the cytoplasm and the plasma membrane of premeiotic testicular germ cells and it may be involved in progression of testicular tumors of germ cell origin. Alternative splicing of this gene results in three transcript variants encoding three different isoforms.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.