




Ttr Antibody

Product Code	CSB-PA025270DA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07309
Immunogen	Recombinant Mouse Transthyretin protein (23-147AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transthyretin (Prealbumin), Ttr
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Ttr

Image



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA025270DA01MO at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.