



Sry Antibody, FITC conjugated

Product Code	CSB-PA022699LC01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05738
Immunogen	Recombinant Mouse Sex-determining region Y protein (1-143AA)
Raised In	Rabbit
Species Reactivity	Mouse
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sex-determining region Y protein (Testis-determining factor), Sry, Tdf Tdy
Immunogen Species	Mus musculus (Mouse)
Research Area	Epigenetics and Nuclear Signaling
Target Names	Sry
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.