



Hbb-b2 Antibody

Product Code	CSB-PA360812EA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02089
Immunogen	Recombinant Mouse Hemoglobin subunit beta-2 protein (2-147AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hemoglobin subunit beta-2 (Beta-2-globin) (Hemoglobin beta-2 chain) (Hemoglobin beta-minor chain), Hbb-b2
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Hbb-b2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.