



ADIPOQ Antibody, Biotin conjugated

Product Code	CSB-PA618019YD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15848
Immunogen	Recombinant Human Adiponectin protein (19-244AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adiponectin (30 kDa adipocyte complement-related protein) (Adipocyte complement-related 30 kDa protein) (ACRP30) (Adipocyte, C1q and collagen domain-containing protein) (Adipose most abundant gene transcript 1 protein) (apM-1) (Gelatin-binding protein), ADIPOQ, ACDC ACRP30 APM1 GBP28
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ADIPOQ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.