




ADIPOQ Antibody

Product Code	CSB-PA618019YA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15848
Immunogen	Recombinant Human Adiponectin protein (19-244AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adiponectin (30 kDa adipocyte complement-related protein) (Adipocyte complement-related 30 kDa protein) (ACRP30) (Adipocyte, C1q and collagen domain-containing protein) (Adipose most abundant gene transcript 1 protein) (apM-1) (Gelatin-binding protein), ADIPOQ, ACDC ACRP30 APM1 GBP28
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ADIPOQ

Image



Western blot

All lanes: ADIPOQ antibody at 2µg/ml

Lane 1: NIH/3T3 whole cell lysate


Lane 2: HL60 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 27 kDa

Observed band size: 27 kDa



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA618019YA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.