



Recombinant Human RNA-binding protein 34 (RBM34)

Product Code	CSB-EP019422HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P42696
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MALEGMSKRK RKRSVQEGEN PDDGVRGSPP EDYRLGQVAS SLFRGEHHSR GGTGRLASLF SSLEPQIQPV YVPVPKQTIK KTKRNEEEES TSQIERPLSQ EPAKKVKAKK KHTNAEKKLA DRESALASAD LEEEIHQKQG QKRKNSQPGV KVADRKILDD TEDTVVSQRK KIQINQEEER LKNERTVFVG NLPVTCNKKK LKSFFKEYGQ IESVRFRSLI PAEGTLSKKL AAIKRKIHPD QKNINAYVVF KEESAATQAL KRNGAQIADG FRIRVDLASE TSSRDKRSVF VGNLPYKVEE SAIEKHFLDC GSIMAVRIVR DKMTGIGKGF GYVLFENTDS VHLALKLNNS ELMGRKLRVM RSVNKEKFKQ QNSNPRLKNV SKPKQGLNFT SKTAEGHPKS LFIGEKAVLL KTKKKGQKKS GRPKKQRKQK
Source	E.coli
Target Names	RBM34
Protein Names	Recommended name: RNA-binding protein 34 Alternative name(s): RNA-binding motif protein 34
Expression Region	1-430
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.