



Recombinant Mouse Histone H1.1 (Hist1h1a)

Product Code	CSB-MP010376MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P43275
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MSETAPVAQA ASTATEKPAA AKKTKKPAKA AAPRKKPAGP SVSELIVQAV SSSKERSGVS LAALKKSLAA AGYDVEKNNS RIKLGLKSLV NKGTLVQTKG TGAAGSFKLN KKAESKAITT KVSVKAKASG AAKKPKKTAG AAKKTVKTP KKPKKPAVSK KTSKSPKKPK VVKAKKVAKS PAKAKAVKPK ASKAKVTKPK TPAKPKKAAP KKK
Source	Mammalian cell
Target Names	Hist1h1a
Protein Names	Recommended name: Histone H1.1 Alternative name(s): H1 VAR.3 Short name= H1a
Expression Region	1-213
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	Histones are basic nuclear proteins responsible for nucleosome structure of the chromosomal fiber in eukaryotes. Two molecules of each of the four core histones (H2A, H2B, H3, and H4) form an octamer, around which approximately 146 bp of DNA is wrapped in repeating units, called nucleosomes. The linker histone, H1, interacts with linker DNA between nucleosomes and functions in the compaction of chromatin into higher order structures. This gene is intronless and encodes a member of the histone H1 family. Transcripts from this gene lack polyA tails but instead contain a palindromic termination element. This gene is found in the large histone gene cluster on chromosome 6.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.