



Recombinant Human DNA repair protein RAD52 homolog (RAD52)

Product Code	CSB-YP019269HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P43351
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSGTEEAILG GRDSHPAAGG GSVLCFGQCQ YTAEEYQAIQ KALRQRLGPE YISSRMAGGG QKVCIYIEGHR VINLANEMFG YNGWAHSITQ QNVDFVDLNN GKFYVGVCAF VRVQLKDGSY HEDVGYGVSE GLKSKALSLE KARKEAVTDG LKRALRSFGN ALGNCILDKD YLRSLNKLPR QLPLEVDLTK AKRQDLEPSV EEARYNSCRP NMALGHPQLQ QVTSPSRPSH AVIPADQDCS SRSLSSSAVE SEATHQRKLR QKQLQQQFRE RMEKQQVRVS TPSAEKSEAA PPAPPVTHST PVTVSEPLLE KDFLAGVTQE LIKTLEDNSE KWAVTPDAGD GVVKPSSRAD PAQTSIDLAL NNQMVNTQNRTP HSHVCHQKPQ AKSGSWDLQT YSADQRTTGN WESHRKSQDM KKRKYDPS
Source	Yeast
Target Names	RAD52
Protein Names	Recommended name: DNA repair protein RAD52 homolog
Expression Region	1-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This protein shares similarity with <i>Saccharomyces cerevisiae</i> Rad52, a protein important for DNA double-strand break repair and homologous recombination. This gene product was shown to bind single-stranded DNA ends, and mediate the DNA-DNA interaction necessary for the annealing of complementary DNA strands. It was also found to interact with DNA recombination protein RAD51, which suggested its role in RAD51 related DNA recombination and repair.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.