



Recombinant Chicken Transcription factor GATA-4 (GATA4)

Product Code	CSB-EP009277CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P43691
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	>85% (SDS-PAGE)
Sequence	PVYVPTTRVP SMLPSLPYLP SSGSSQQASP VSSHSIWTQP GAESAAYNPG SSHPPVSPRF SFSTSTPIPS TSSRDAAAAY SSSLSLSANG REQYSRGFGS SYSSPYPAYV SPEMATTWTS SPFDSPMLHN LQSRGTPAAA RHANIEFFDD YSEGRECVNC GAMSTPLWRR DGTGHYLCNA CGLYHKMNGI NRPLFKPQRR LSASRRVGLS CANCHTTTTT LWRRNAEGEP VCNACGLYMK LHGVRPLAM RKEGIQTRKR KPKNLNKTKT PAGPSSSESL TPTTSSTSSS SSATTEEMR PIKTEPGLSS HYGHPSPISQ AFSVSAMSGH GSSIHPAISA LKLSPQAYQS AISQSPQASS KQDSWNSLVL AENHGDIITA
Source	E.coli
Target Names	GATA4
Protein Names	Recommended name: Transcription factor GATA-4 Alternative name(s): GATA-binding factor 4
Expression Region	1-380
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This gene encodes a member of the GATA family of zinc-finger transcription factors. Members of this family recognize the GATA motif which is present in the promoters of many genes. This protein is thought to regulate genes involved in embryogenesis and in myocardial differentiation and function. Mutations in this gene have been associated with cardiac septal defects.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.