



# Recombinant Chicken Transcription factor GATA-4 (GATA4)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-YP009277CH   |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.  |
| <b>Uniprot No.</b>       | P43691   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Gallus gallus (Chicken)  |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | PVYVPTTRVP SMLPSLPYLP SSGSSQQASP VSSHSIWTQP GAESAAYNPG<br>SSHPPVSPRF SFSTSTPIPS TSSRDAAAAY SSSLSLANG REQYSRGFGS<br>SYSSPYPAYV SPEMATTWTS SPFDSPMLHN LQSRGTPAAA RHANIEFFDD<br>YSEGRECVNC GAMSTPLWRR DGTGHYLCNA CGLYHKMNGI<br>NRPLFKPQRR LSASRRVGLS CANCHTTTTT LWRRNAEGEP VCNACGLYMK<br>LHGVRPLAM RKEGIQTRKR KPKNLNKTKT PAGPSSSESL TPTTSSTSSS<br>SSATTEEMR PIKTEPGLSS HYGHPSPISQ AFSVSAMSGH GSSIHPAISA<br>LKLSQAYQS AISQSPQASS KQDSWNSLVL AENHGDIITA |
| <b>Source</b>            | Yeast  |
| <b>Target Names</b>      | GATA4  |
| <b>Protein Names</b>     | Recommended name: Transcription factor GATA-4 Alternative name(s): GATA-binding factor 4   |
| <b>Expression Region</b> | 1-380  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | Full length protein  |
| <b>Target Details</b>    | This gene encodes a member of the GATA family of zinc-finger transcription factors. Members of this family recognize the GATA motif which is present in the promoters of many genes. This protein is thought to regulate genes involved in embryogenesis and in myocardial differentiation and function. Mutations in this gene have been associated with cardiac septal defects.  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life   |



of lyophilized form is 12 months at -20°C/-80°C.