



Recombinant Mouse Nuclear receptor subfamily 5 group A member 2 (Nr5a2)

Product Code	CSB-MP016066MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P45448
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MSASLDTGDF QEFLKHGLTA IASAPGSETR HSPKREEQLR EKRAGLPDRH RRPIPARSRL VMLPKVETEA PGLVRSHGEQ GQMPENMQVS QFKMVNYSYD EDLEELCPVC GDKVSGYHYG LLTCECKGF FKRTVQNQKR YTCIENQNCQ IDKTQRKRCP YCRFKKCIDV GMKLEAVRAD RMRGGRNKFG PMYKRDRALK QQKKALIRAN GLKLEAMSQV IQAMPSDLTS AIQNIHSASK GLPLSHVALP PTDYDRSPFV TSPISMTMPP HSSLHGYQPY GHFPSRAIKS EYPDPYSSSP ESMMGYSYMD GYQTNSPASI PHLILELLKC EPDEPQVQAK IMAYLQQEQS NRNRQEKLSA FGLLCKMADQ TLFSIVEWAR SSIFFRELKV DDQMKLLQNC WSELLILDHI YRQVAHGKEG TIFLVTGEHV DYSTIISHTE VAFNNLLSLA QELVVRLRSL QFDQREFVCL KFLVLFSSDV KNLENLQLVE GVQEQVNAAL LDYTVCNYPQ QTEKFGQLLL RLPEIRAIK QAEDYLYYKH VNGDVPYNNL LIEMLHAKRA
Source	Mammalian cell
Target Names	Nr5a2
Protein Names	Recommended name: Nuclear receptor subfamily 5 group A member 2 Alternative name(s): Liver receptor homolog 1 Short name= LRH-1
Expression Region	1-560
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.