



# Recombinant Mouse Nuclear receptor subfamily 5 group A member 2 (Nr5a2)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-YP016066MO   |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.  |
| <b>Uniprot No.</b>       | P45448   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Mus musculus (Mouse)   |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MSASLDTGDF QEFLKHGLTA IASAPGSETR HSPKREEQLR EKRAGLPDRH<br>RRPIPARSRL VMLPKVETEA PGLVRSHGEQ GQMPENMQVS QFKMVNYSYD<br>EDLEELCPVC GDKVSGYHYG LLTCECKGF FKRTVQNQKR YTCIENQNCQ<br>IDKTQRKRCP YCRFKKCIDV GMKLEAVRAD RMRGGRNKFG PMYKRDRALK<br>QQKKALIRAN GLKLEAMSQV IQAMPSDLTS AIQNIHSASK GLPLSHVALP<br>PTDYDRSPFV TSPISMTMPP HSSLHGYQPY GHFPSRAIKS EYPDPYSSSP<br>ESMMGYSYMD GYQTNSPASI PHLILELLKC EPDEPQVQAK IMAYLQQEQS<br>NRNRQEKLSA FGLLCKMADQ TLFSIVEWAR SSIFFRELKV DDQMKLLQNC<br>WSELLILDHI YRQVAHGKEG TIFLVTGEHV DYSTIISHTE VAFNNLLSLA<br>QELVVRLRSL QFDQREFVCL KFLVLFSSDV KNLENLQLVE GVQEQVNAAL<br>LDYTVCNYPQ QTEKFGQLLL RLPEIRAIK QAEDYLYYKH VNGDVPYNNL<br>LIEMLHAKRA |
| <b>Source</b>            | Yeast  |
| <b>Target Names</b>      | Nr5a2  |
| <b>Protein Names</b>     | Recommended name: Nuclear receptor subfamily 5 group A member 2<br>Alternative name(s): Liver receptor homolog 1 Short name= LRH-1   |
| <b>Expression Region</b> | 1-560  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | Full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |