



Recombinant Human Nuclear receptor subfamily 4 group A member 1 (NR4A1)

| | |
|--------------------------|--|
| Product Code | CSB-MP016062HU |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P22736 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MPCIQAQYGT PAPSPGPRDH LASDPLTPEF IKPTMDLASP EAAPAAPTAL PSFSTFMDGY TGEFDTFLYQ LPGTVQPCSS ASSSASSTSS SSATSPASAS FKFEDFQVYG CYPGPLSGPV DEALSSSSGSD YYGSPCSAPS PSTPSFQPPQ LSPWDGSFGH FSPSQTYEGL RAWTEQLPKA SGPPQPPAFF SFSPPTGPSP SLAQSPKLF PSQATHQLGE GESYSMPTAF PGLAPTSPHL EGSGILDTPV TSTKARSGAP GGSEGRCAVC GDNASCQHYG VRTCEGCKGF FKRTVQKNAK YICLANKDCP VDKRRRNRCQ FCRFQKCLAV GMVKEVVRTD SLKGRRGRLP SKPKQPPDAS PANLLTSLVR AHLDSGPSTA KLDYSKFQEL VLPHFGKEDA GDVQQFYDLL SGSLEVIRKW AEKIPGFAEL SPADQDLLLE SAFLELFILR LAYRSKPGEG KLIFCSGLVL HRLQCARGFG DWIDSILAFS RSLHLLVDV PAFACLSALV LITDRHGLQE PRRVEELQNR IASCLKEHVA AVAGEPQPAS CLSRLGKLP ELRTLCTQGL QRIFYLKLED LVPPPIIDK IFMDTLPF |
| Source | Mammalian cell |
| Target Names | NR4A1 |
| Protein Names | Recommended name: Nuclear receptor subfamily 4 group A member 1 Alternative name(s): Early response protein NAK1 Nuclear hormone receptor NUR/77 Short name= Nur77 Orphan nuclear receptor HMR Orphan nuclear receptor TR3 ST-59 |
| Expression Region | 1-598 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Target Details | This gene encodes a member of the steroid-thyroid hormone-retinoid receptor superfamily. Expression is induced by phytohemagglutinin in human lymphocytes and by serum stimulation of arrested fibroblasts. The encoded protein acts as a nuclear transcription factor. Translocation of the protein from the nucleus to mitochondria induces apoptosis. Multiple alternatively spliced variants, encoding the same protein, have been identified. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the |



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.