



# Recombinant Bovine Cytochrome b-c1 complex subunit 2, mitochondrial (UQCRC2)

<b>Product Code</b>	CSB-EP025668BO-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P23004
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SLKVAP KVKATEAPAG VPPHPQDLEF TRLPNGLVIA SLENYAPASR IGLFIKAGSR YENSNNLGTS HLLRLASSLT TKGASSFKIT RGIEAVGGKL SVTSTRENMA YTVECLRDDV DILMEFLNV TTAPEFRRWE VAALQPQLRI DKAVALQNPQ AHVIENLHAA AYRNALANSL YCPDYRIGKV TPVELHDYVQ NHFTSARMAL IGLGVSHPV LKQVAEQFLNI RGGLGLSGAK AKYHGGEIRE QNGDSL VHAA LVAESAAIGS AEANAFSVLQ HVLGAGPHVK RGSNATSSLY QAVAKGVHQP FDVSAFNASY SDSGLFGFYT ISQAASAGDV IKAAYNQVKT IAQGNLSNP D VQAAKNLKA GYLMSVESSE GFLDEVGSQA LAAGSYTPPS TVLQQIDAVA DADVINAACK FVSGRKSMAA SGNLGH TPF I DEL
<b>Source</b>	E.coli
<b>Target Names</b>	UQCRC2
<b>Protein Names</b>	Recommended name: Cytochrome b-c1 complex subunit 2, mitochondrial Alternative name(s): Complex III subunit 2 Core protein II Ubiquinol-cytochrome-c reductase complex core protein 2
<b>Expression Region</b>	15-453
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.