



Recombinant Human Splicing factor, proline- and glutamine-rich (SFPQ)

Product Code	CSB-EP021136HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P23246
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSRDRFRSRG GGGGGFHRRG GGGGRGGLHD FRSPPPGMGL NQNRGPMGPG PGQSGPKPPI PPPPHQQQQ QPPQPPPPQ QPPPHQPPPH PQPHQQQQPP PPPQDSSKPV VAQGPGPAPG VGSAPPASSS APPATPPTSG APPGSGPGPT PTPPPAVTSA PPGAPPPTPP SSGVPTTPPQ AGGPPPPAA VPGPGPGPKQ GPGPGGPKGG KMPGGPKPGG GPGLSTPGGH PKPPHRGGGE PRGGRQHHP YHQHHQGPP PGGPGRSEE KISDSEGFKA NLSLLRRPGE KTYTQRCRLF VGNLPADITE DEFKRLFAKY GEPGEVFINK GKGFGFIKLE SRALAEIAKA ELDDTPMRGR QLRVRFATHA AALSVRNLSP YVSNELLEEA FSQFGPIERA VVIVDDRGRS TGKGIVEFAS KPAARKAFER CSEGVLLTT TPRPVIVEPL EQLDDEDGLP EKLAQKNPMY QKERETPPRF AQHGTFEY SQRWKSLEDEM EKQREQVEK NMKDAKDKLE SEMEDAYHEH QANLLRQDLM RRQEELRRME ELHNQEMQKR KEMQLRQEEE RRRREEEMMI RQREMEEQMR RQREESYSRM GYMDPRERDM RMGGGGAMNM GDPYSGGGQK FPPLGGGGGI GYEANPGVPP ATMSGSMMS DMRTFRFGQG GAGPVGGQGP RGMGPGTPAG YGRGREEYEG PNKKPRF
Source	E.coli
Target Names	SFPQ
Protein Names	Recommended name: Splicing factor, proline- and glutamine-rich Alternative name(s): 100 kDa DNA-pairing protein Short name= hPOMp100 DNA-binding p52/p100 complex, 100 kDa subunit Polypyrimidine tract-binding protein-associated-splic
Expression Region	1-707
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.