



Recombinant Human COUP transcription factor 2 (NR2F2)

Product Code	CSB-MP016057HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P24468
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAMVVSTWRD PQDEVPGSQG SQASQAPPVP GPPPGAPHTP QTPGQGGPAS TPAQTAAGGQ GPGGGPGSDK QQQQHHIECV VCGDKSSGKH YGQFTCEGCK SFFKRSVRRN LSYTCRANRN CPIDQHHRNQ CQYCRLKKCL KVGMRREAVQ RGRMPPTQPT HGQFALTNGD PLNCHSYLSG YISLLLRAEP YPTSRFGSQC MQPNNIMGIE NICELAARML FSAVEWARNI PFFPDLQITD QVALLRLTWS ELFVLNAAQC SMPLHVAPLL AAAGLHASPMD SADRVAFMD HIRIFQEVE KLKALHVDSA EYSCLKAIVL FTSDACGLSD VAHVESLQEK SQCALEEYVR SQYPNQPTRF GKLLLRLPSL RTVSSSVIEQ LFFVRLVGKT PIETLIRDML LSGSSFNWPY MAIQ
Source	Mammalian cell
Target Names	NR2F2
Protein Names	Recommended name: COUP transcription factor 2 Short name= COUP-TF2 Alternative name(s): Apolipoprotein A-I regulatory protein 1 Short name= ARP-1 COUP transcription factor II Short name= COUP-TF II Nuclear receptor sub
Expression Region	1-414
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.