



Recombinant Mouse Protein C-ets-1 (Ets1)

Product Code	CSB-EP007852MO-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P27577
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MKAAVDLKPT LTIKTEKVD LELFSPDME CADVPLLTPS SKEMMSQALK ATFSGFTKEQ QRLGIPKDPR QWTETHVRDW VMWAVNEFSL KGVDVFQKFCM SGAALCALGK ECFLELAPDF VGDILWEHLE ILQKEDVKPY QVNGANPTYF ESCYTSDYFI SYGIEHAQCV PPSEFSEPSF ITESYQTLHP ISSEELLSLK YENDYPSVIL QDPLQTDTLQ TDYFAIKQEV LTPDNMCLGR ASRGKLGGQD SFESVESYDS CDRLTQSWSS QSSFNSLQRV PSYDSFDYED YPAALPNHKP KGTFKDYVRD RADLNKDKPV IPAAALAGYT GSGPIQLWQF LLELLTDKSC QSFISWTGDG WEFKLSDPDE VARRWGKRKN KPKMNYEKLS RGLRYYYDKN IHKTAGKRY VYRFVCDLQS LLGYTPEELH AMLDVKPDAD
Source	E.coli
Target Names	Ets1
Protein Names	Recommended name: Protein C-ets-1 Alternative name(s): p54
Expression Region	1-440
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.