



Recombinant Mouse Polyadenylate-binding protein 1 (Pabpc1)

Product Code	CSB-EP017352MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P29341
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MNPSAPSYPM ASLYVGD LHP DVTEAMLYEK FSPAGPILSI RVC RDMITRR SLGYAYVNFQ QPADAERALD TMNFDVIK GK PVRIMWSQRD PSLRKSGVGN IFIKNLDSI DNKALYDTFS AFGNILSCKV VCDENGSKGY GFVHFETQEA AERAIEKMNG MLLNDRKVFV GRFKSRKERE AELGARAKEF TNVYIKNFGE DMDDERL KEL FGKFGPALS V KVM TDES GKS KGFGFVS FER HEDAQKAVDE MNGKELNGKQ IYVGRAQKKV ERQTELKRKF EQMKQDRITR YQGVNLYVKN LDDGIDDERL RKEFSPFGTI TSAKVMMEGG RSKGFGFVCF SSPEEATKAV TEMNGRIVAT KPLYVALAQR KEERQAHLTN QYMQRMASVR AVPNPVINPY QPAPPSGYFM AAIPQTQNR AYYPPSQIAQ LRPSRWTAQ GARPHPFQNM PGAIRPAAPR PPFSTMRPAS SQVPRVMSTQ RVANTSTQTM GPRPAAAAA ATPAVRTVPQ YKYAAGVRNP QQHLNAQPQV TMQQPAVHVQ GQEPLTASML ASAPPQEQQK MLGERLFPLI QAMHPSLAGK ITGM LLEIDN SELLHMLES P ESLRSKVDEA VAVLQAHQAK EAAQKAVNSA TGVPTV
Source	E.coli
Target Names	Pabpc1
Protein Names	Recommended name: Polyadenylate-binding protein 1 Short name= PABP-1 Short name= Poly(A)-binding protein 1
Expression Region	1-636
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.