



Recombinant Rat Neutrophil gelatinase-associated lipocalin (Lcn2)

Product Code	CSB-EP012816RA-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P30152
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	QDSTQNLIPA PPLISVPLQP GFWTERFQGR WFVVGLAGNA VQKERQSRFT MYSTIYELQE DNSYNVTSIL VRGQGCRYWI RTFVPSSRPG QFTLGNIHSY PQIQSYDVQV ADTDYDQFAM VFFQKTSENK QYFKVTLYGR TKGLSDELKE RFVSAKSLG LKDNNIVFSV PTDQCIDN
Source	E.coli
Target Names	Lcn2
Protein Names	Recommended name: Neutrophil gelatinase-associated lipocalin Short name= NGAL Alternative name(s): Alpha-2-microglobulin-related protein Alpha-2U globulin-related protein Lipocalin-2 p25
Expression Region	21-198
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.