



CUSABIO catalogue of trial size antibodies

Product Name	Code	Species Reactivity	Applications
<u>2.5 Antibody</u>	CSB-PA366021HA01EEB	Enterobacteria phage T7	ELISA
<u>6*His Monoclonal Antibody</u>	CSB-MA000011M0m	All	ELISA, WB
<u>A1BG Antibody</u>	CSB-PA001001LA01HU	Human	ELISA, IHC
<u>A2M Antibody</u>	CSB-PA14579A0Rb	Human	ELISA, IHC, IF
<u>A2ML1 Antibody</u>	CSB-PA001005LA01HU	Human	ELISA, WB, IHC
<u>A4GNT Antibody</u>	CSB-PA001009LA01HU	Human	ELISA, IHC
<u>aac Antibody</u>	CSB-PA333489LA01ACY	Actinoplanes utahensis	ELISA
<u>AACS Antibody</u>	CSB-PA001012LA01HU	Human	ELISA, WB, IHC, IP
<u>AADAC Antibody</u>	CSB-PA001014LA01HU	Human	ELISA, IHC
<u>AADACL2 Antibody</u>	CSB-PA001016LA01HU	Human	ELISA, WB
<u>AADAT Antibody</u>	CSB-PA843276LA01HU	Human	ELISA, IHC
<u>AAGAB Antibody</u>	CSB-PA747468LA01HU	Human	ELISA, IHC, IF
<u>AAK1 Antibody</u>	CSB-PA643571LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>AANAT Antibody</u>	CSB-PA623090LA01HU	Human	ELISA, WB, IF
<u>AAO Antibody</u>	CSB-PA333022LA01DZK	Cucurbita maxima	ELISA, WB
<u>AAP1 Antibody</u>	CSB-PA327747LA01SVG	Saccharomyces cerevisiae	ELISA
<u>AAR2 Antibody</u>	CSB-PA896704LA01HU	Human	ELISA, IHC
<u>AARS Antibody</u>	CSB-PA001023LA01HU	Human	ELISA, IHC, IF
<u>AARSD1 Antibody</u>	CSB-PA863119LA01HU	Human	ELISA, WB, IHC
<u>AASDHPPT Antibody</u>	CSB-PA865121LA01HU	Human	ELISA, WB, IHC
<u>AASS Antibody</u>	CSB-PA883373LA01HU	Human	ELISA, WB, IHC
<u>AATF Antibody</u>	CSB-PA05264A0Rb	Human	ELISA, IHC
<u>ABCA12 Antibody</u>	CSB-PA774804LA01HU	Human	ELISA, IF
<u>ABCA2 Antibody</u>	CSB-PA887172LA01HU	Human	ELISA, IHC
<u>ABCA5 Antibody</u>	CSB-PA845163LA01HU	Human	ELISA, IHC
<u>ABCA9 Antibody</u>	CSB-PA812864LA01HU	Human	ELISA, IF
<u>ABCB1 Antibody</u>	CSB-PA11737A0Rb	Human	ELISA, IHC, IF
<u>Abcb4 Antibody</u>	CSB-PA001050LA01MO	Mouse	ELISA
<u>ABCB4 Antibody</u>	CSB-PA001050LA01HU	Human	ELISA, IHC, IF
<u>ABCB5 Antibody</u>	CSB-PA643574LA01HU	Human	ELISA, IHC, IF
<u>ABCB7 Antibody</u>	CSB-PA001053LA01HU	Human	ELISA, IHC
<u>ABCB9 Antibody</u>	CSB-PA889064LA01HU	Human, Mouse	ELISA, WB, IF
<u>ABCC2 Antibody</u>	CSB-PA852918LA01HU	Human	ELISA, IHC, IF
<u>ABCC3 Antibody</u>	CSB-PA001062LA01HU	Human	ELISA, IF
<u>ABCC4 Antibody</u>	CSB-PA001063LA01HU	Human	ELISA, IF
<u>ABCC5 Antibody</u>	CSB-PA001064LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>ABCC6 Antibody</u>	CSB-PA001065LA01HU	Human	ELISA, IF
<u>ABCC8 Antibody</u>	CSB-PA22439A0Rb	Human	ELISA, IHC, IF
<u>ABCC9 Antibody</u>	CSB-PA001067LA01HU	Human	ELISA, IHC, IF
<u>ABCD1 Antibody</u>	CSB-PA001068LA01HU	Human	ELISA, IHC
<u>ABCD2 Antibody</u>	CSB-PA866202LA01HU	Human	ELISA, IHC
<u>ABCD4 Antibody</u>	CSB-PA001075LA01HU	Human	ELISA, IHC
<u>ABCF3 Antibody</u>	CSB-PA889115LA01HU	Human	ELISA, IHC
<u>ABCG1 Antibody</u>	CSB-PA001080LA01HU	Human	ELISA, IHC
<u>ABCG2 Antibody</u>	CSB-PA891568LA01HU	Human	ELISA, IF
<u>ABCG5 Antibody</u>	CSB-PA863956LA01HU	Human	ELISA, WB, IHC, IF
<u>ABCG8 Antibody</u>	CSB-PA875651LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ABHD10 Antibody</u>	CSB-PA865140LA01HU	Human	ELISA, IHC
<u>Abhd11 Antibody</u>	CSB-PA811717LA01MO	Mouse	ELISA
<u>ABHD14A Antibody</u>	CSB-PA874812LA01HU	Human	ELISA, IHC, IF
<u>ABHD14B Antibody</u>	CSB-PA842687LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ABHD18 Antibody</u>	CSB-PA610370LA01HU	Human	ELISA, IHC
<u>ABHD2 Antibody</u>	CSB-PA001094LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>ABHD3 Antibody</u>	CSB-PA845142EA01HU	Human	ELISA, IHC, IF
<u>ABHD5 Antibody</u>	CSB-PA819884HA01HU	Human	ELISA, IF
<u>ABHD6 Antibody</u>	CSB-PA861144LA01HU	Human	ELISA, IHC, IF
<u>ABI2 Antibody</u>	CSB-PA001102LA01HU	Human	ELISA, IHC, IF
<u>ABI3 Antibody</u>	CSB-PA873732LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ABI3BP Antibody</u>	CSB-PA772035LA01HU	Human	ELISA, IHC
<u>ABI5 Antibody</u>	CSB-PA879765LA01DOA	Arabidopsis thaliana	ELISA
<u>ABL1 Antibody</u>	CSB-PA04355A0Rb	Human	ELISA, IHC, IP
<u>ABL2 Antibody</u>	CSB-PA001106LA01HU	Human	ELISA, WB
<u>ABLIM1 Antibody</u>	CSB-PA001107LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ABLIM2 Antibody</u>	CSB-PA753612LA01HU	Human	ELISA, IHC
<u>ABLIM3 Antibody</u>	CSB-PA001109LA01HU	Human	ELISA, IHC, IF
<u>ABO Antibody</u>	CSB-PA001110LA01HU	Human, Mouse	ELISA, WB
<u>ABR Antibody</u>	CSB-PA622521LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ABT1 Antibody</u>	CSB-PA001114LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ABTB1 Antibody</u>	CSB-PA839277LA01HU	Human	ELISA
<u>AC1 Antibody</u>	CSB-PA324434LA01ADX	African cassava mosaic virus	ELISA
<u>ACAA1 Antibody</u>	CSB-PA001117LA01HU	Human	ELISA, IHC
<u>ACAA2 Antibody</u>	CSB-PA001118LA01HU	Human, Mouse, Rat	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF, IP
<u>ACAD11 Antibody</u>	CSB-PA754597LA01HU	Human	ELISA, IHC, IF
<u>ACAD8 Antibody</u>	CSB-PA892457LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ACAD9 Antibody</u>	CSB-PA887999EA01HU	Human, Mouse	ELISA, WB, IHC
<u>ACADL Antibody</u>	CSB-PA001125LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ACADM Antibody</u>	CSB-PA001126LA01HU	Human	ELISA, IHC
<u>ACADSB Antibody</u>	CSB-PA001128LA01HU	Human, Rat	ELISA, WB, IHC
<u>ACADVL Antibody</u>	CSB-PA001129LA01HU	Human	ELISA, WB, IHC, IF
<u>ACAP1 Antibody</u>	CSB-PA622985EA01HU	Human, Mouse	ELISA, WB, IHC
<u>ACAT1 Antibody</u>	CSB-PA001134LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>ACBD3 Antibody</u>	CSB-PA863972LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>ACBD4 Antibody</u>	CSB-PA818758LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ACBD5 Antibody</u>	CSB-PA725893HA01HU	Human	ELISA, WB
<u>ACBD7 Antibody</u>	CSB-PA854056HA01HU	Human	ELISA
<u>ACD Antibody</u>	CSB-PA839289LA01HU	Human, Mouse	ELISA, WB
<u>Ace Antibody</u>	CSB-PA001149LA01RA	Rat	ELISA, IHC
<u>ACE2 Antibody</u>	CSB-PA866317LA01HU	Human	ELISA, IHC, IF
<u>ACHE Antibody</u>	CSB-PA001154EA01BO	Bovine	ELISA
<u>ACHE Antibody</u>	CSB-PA5619A0Rb	Human	ELISA, IHC, IF
<u>ache Antibody</u>	CSB-PA863598LA01DIL	Zebrafish	ELISA, WB
<u>ACKR1 Antibody</u>	CSB-PA624105LA01HU	Human, Mouse	ELISA, WB, IF
<u>ACLY Antibody</u>	CSB-PA001158LA01HU	Human	ELISA, WB, IHC, IF
<u>ACMSD Antibody</u>	CSB-PA851554LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>ACO1 Antibody</u>	CSB-PA001161LA01HU	Human	ELISA, WB, IHC, IF
<u>ACO2 Antibody</u>	CSB-PA859943LA01HU	Human, Rat, Mouse	ELISA, WB, IHC
<u>ACOD1 Antibody</u>	CSB-PA011825LA01HU	Human	ELISA, IHC
<u>ACOT1 Antibody</u>	CSB-PA768206LA01HU	Human	ELISA, IHC, IF
<u>ACOT11 Antibody</u>	CSB-PA855502LA01HU	Human	ELISA
<u>ACOT6 Antibody</u>	CSB-PA668210LA01HU	Human	ELISA, IHC, IF
<u>ACOT8 Antibody</u>	CSB-PA001170LA01HU	Human	ELISA, IHC
<u>ACOT9 Antibody</u>	CSB-PA897539LA01HU	Human	ELISA, IHC
<u>ACOX1 Antibody</u>	CSB-PA617998HA01HU	Human	ELISA, IHC, IF
<u>ACOX3 Antibody</u>	CSB-PA001174LA01HU	Human	ELISA
<u>ACP6 Antibody</u>	CSB-PA889072LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>acpP Antibody</u>	CSB-PA015636DA01SKX	Staphylococcus aureus	ELISA
<u>ACPP Antibody</u>	CSB-PA11309A0Rb	Human	ELISA, IHC
<u>acpS Antibody</u>	CSB-PA740005HA01SKV	Staphylococcus aureus	ELISA



Product Name	Code	Species Reactivity	Applications
<u>ACRV1 Antibody</u>	CSB-PA001186HA01HU	Human	ELISA, IHC
<u>ACS8 Antibody</u>	CSB-PA871226LA01DOA	Arabidopsis thaliana	ELISA, WB
<u>ACSBG1 Antibody</u>	CSB-PA853424LA01HU	Human	ELISA, IHC
<u>ACSBG2 Antibody</u>	CSB-PA685911LA01HU	Human	ELISA, WB, IHC
<u>ACSF2 Antibody</u>	CSB-PA853397LA01HU	Human	ELISA, IHC
<u>ACSF3 Antibody</u>	CSB-PA674579LA01HU	Human	ELISA, IHC
<u>ACSL1 Antibody</u>	CSB-PA001191LA01HU	Human	ELISA, IF
<u>ACSL4 Antibody</u>	CSB-PA060488EA01HU	Human	ELISA, WB, IHC
<u>ACSL5 Antibody</u>	CSB-PA891734HA01HU	Human, Rat	ELISA, WB, IHC
<u>ACSL5 Antibody</u>	CSB-PA891734LA01HU	Human	ELISA, IHC
<u>Acss2 Antibody</u>	CSB-PA882544LA01MO	Mouse, Rat	ELISA, WB
<u>ACTA1 Antibody</u>	CSB-PA001205LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>ACTA1 Antibody</u>	CSB-PA001205NJ01HU	Human	ELISA, IHC, IF
<u>ACTA1 Antibody</u>	CSB-PA001205OJ01HU	Human	ELISA
<u>Acta1 Antibody</u>	CSB-PA00750E0Rb	Mouse, Rat, Zebrafish	ELISA, WB
<u>acta1b Antibody</u>	CSB-PA20477A0Rb	Zebrafish	ELISA, WB
<u>ACTA2 Antibody</u>	CSB-PA12659A0Rb	Human	ELISA, IHC
<u>ACTB Antibody</u>	CSB-PA001207NJ01HU	Human	ELISA, IHC, IF
<u>ACTB Antibody</u>	CSB-PA01629A0Rb	Human	ELISA, WB, IHC, IF
<u>Actb Antibody</u>	CSB-PA07675A0Rb	Mouse, Human	ELISA, WB
<u>actba Antibody</u>	CSB-PA800252LA01DIL	Zebrafish	ELISA, WB
<u>actbb Antibody</u>	CSB-PA768804LA01DIL	Zebrafish	ELISA
<u>ACTG1 Antibody</u>	CSB-PA02619A0Rb	Human	ELISA, WB, IHC, IF
<u>ACTG2 Antibody</u>	CSB-PA12799A0Rb	Human	ELISA, WB, IHC
<u>ACTL7A Antibody</u>	CSB-PA896747LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ACTL7B Antibody</u>	CSB-PA897577LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ACTL8 Antibody</u>	CSB-PA887982LA01HU	Human	ELISA, WB, IHC
<u>ACTL9 Antibody</u>	CSB-PA840564LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ACTN1 Antibody</u>	CSB-PA001241LA01HU	Human	ELISA, WB, IF
<u>ACTN3 Antibody</u>	CSB-PA001243LA01HU	Human	ELISA, IF
<u>ACTN4 Antibody</u>	CSB-PA00814A0Rb	Human, Rat	ELISA, WB, IHC, IF
<u>ACTN4 Antibody</u>	CSB-PA00819A0Rb	Human, Mouse	ELISA, WB, IHC
<u>ACTR10 Antibody</u>	CSB-PA889156LA01HU	Human	ELISA, WB
<u>ACTR1A Antibody</u>	CSB-PA001246EA01HU	Human	ELISA, WB, IHC
<u>ACTR3 Antibody</u>	CSB-PA878955LA01HU	Human	ELISA, IHC, IF
<u>ACTR5 Antibody</u>	CSB-PA884496LA01HU	Human	ELISA, IHC, IF
<u>ACTR8 Antibody</u>	CSB-PA872523LA01HU	Human	ELISA, WB
<u>ACTRT1 Antibody</u>	CSB-PA001255LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>ACTRT2 Antibody</u>	CSB-PA001256LA01HU	Human	ELISA, WB
<u>ACTRT3 Antibody</u>	CSB-PA880135LA01HU	Human	ELISA, WB, IHC
<u>ACVR1B Antibody</u>	CSB-PA001258LA01HU	Human	ELISA, IF
<u>ACVRL1 Antibody</u>	CSB-PA00839A0Rb	Human	ELISA, WB, IF
<u>ACY3 Antibody</u>	CSB-PA842682LA01HU	Human	ELISA, WB, IF
<u>ACYP2 Antibody</u>	CSB-PA001266LA01HU	Human	ELISA, IHC, IF
<u>ADA Antibody</u>	CSB-PA001268EA01HU	Human	ELISA, WB, IHC, IF
<u>ADAL Antibody</u>	CSB-PA724932LA01HU	Human, Mouse	ELISA, WB
<u>ADAM10 Antibody</u>	CSB-PA001273LA01HU	Human	ELISA, IHC, IF
<u>ADAM12 Antibody</u>	CSB-PA001275LA01HU	Human	ELISA, IHC, IF
<u>ADAM15 Antibody</u>	CSB-PA01915A0Rb	Human	ELISA, IHC
<u>ADAM23 Antibody</u>	CSB-PA001285LA01HU	Human	ELISA, WB, IHC, IF
<u>ADAM29 Antibody</u>	CSB-PA891722LA01HU	Human	ELISA, IHC, IF
<u>ADAM7 Antibody</u>	CSB-PA863964LA01HU	Human	ELISA, IHC
<u>ADAM8 Antibody</u>	CSB-PA10259A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>ADAMDEC1 Antibody</u>	CSB-PA001297LA01HU	Human	ELISA, IHC
<u>ADAMTS12 Antibody</u>	CSB-PA001300LA01HU	Human	ELISA, IHC
<u>ADAMTS13 Antibody</u>	CSB-PA001301LA01HU	Human	ELISA, IHC
<u>ADAMTS15 Antibody</u>	CSB-PA837445LA01HU	Human	ELISA, IHC
<u>ADAMTS4 Antibody</u>	CSB-PA001311LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ADAMTS5 Antibody</u>	CSB-PA890943LA01HU	Human	ELISA, IHC
<u>ADAMTS7 Antibody</u>	CSB-PA891951LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ADAMTS9 Antibody</u>	CSB-PA868397LA01HU	Human	ELISA, IF
<u>ADAP2 Antibody</u>	CSB-PA889071LA01HU	Human	ELISA
<u>ADAR Antibody</u>	CSB-PA001324LA01HU	Human	ELISA, WB, IHC
<u>ADAT1 Antibody</u>	CSB-PA863123LA01HU	Human	ELISA, WB, IHC
<u>Adck2 Antibody</u>	CSB-PA761222LA01MO	Mouse, Human	ELISA, WB, IHC
<u>ADCY1 Antibody</u>	CSB-PA600142LA01HU	Human	ELISA, IHC, IF
<u>ADCY10 Antibody</u>	CSB-PA839393LA01HU	Human	ELISA, IHC
<u>ADCY5 Antibody</u>	CSB-PA001341LA01HU	Human	ELISA, IF
<u>ADCY6 Antibody</u>	CSB-PA001342LA01HU	Human	ELISA, IHC, IF
<u>ADCY8 Antibody</u>	CSB-PA001344LA01HU	Human	ELISA, IHC, IF
<u>add Antibody</u>	CSB-PA001268LA01ENV	Escherichia coli	ELISA
<u>ADGRA2 Antibody</u>	CSB-PA857001LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>ADGRA3 Antibody</u>	CSB-PA808998LA01HU	Human	ELISA, IF
<u>ADGRB2 Antibody</u>	CSB-PA523985LA01HU	Human, Rat	ELISA, WB, IF
<u>ADGRE2 Antibody</u>	CSB-PA887050LA01HU	Human	ELISA, IHC, IF
<u>ADGRG1 Antibody</u>	CSB-PA812970LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>ADGRG3 Antibody</u>	CSB-PA769807LA01HU	Human	ELISA
<u>ADGRL4 Antibody</u>	CSB-PA872541LA01HU	Human	ELISA, WB, IHC, IF
<u>ADH1A Antibody</u>	CSB-PA001353LA01HU	Human	ELISA, IHC, IF
<u>ADH1B Antibody</u>	CSB-PA001354HA01HU	Human	ELISA, WB, IHC, IF
<u>Adh5 Antibody</u>	CSB-PA001357LA01MO	Mouse	ELISA, WB
<u>ADH5 Antibody</u>	CSB-PA001357LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>ADH7 Antibody</u>	CSB-PA001362LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>Adh7 Antibody</u>	CSB-PA717544LA01MO	Mouse	ELISA, WB
<u>Adipoq Antibody</u>	CSB-PA07946A0Rb	Mouse	ELISA
<u>ADIPOQ Antibody</u>	CSB-PA07956A0Rb	Human, Mouse	ELISA, WB, IHC
<u>ADIPOQ Antibody</u>	CSB-PA07957A0Rb	Human	ELISA, IHC
<u>ADIPOQ Antibody</u>	CSB-PA618019YA01HU	Human, Mouse	ELISA, WB, IHC
<u>ADIPOQ Antibody</u>	CSB-PA661101LA01BO	Bovine	ELISA
<u>ADIPOQ Antibody</u>	CSB-PA661101YA01BO	Bovine	ELISA, WB
<u>ADIPOR1 Antibody</u>	CSB-PA001367LA01HU	Human	ELISA, WB, IHC, IF
<u>ADK Antibody</u>	CSB-PA346232LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>ADM Antibody</u>	CSB-PA001370HA01HU	Human	ELISA, IHC
<u>ADNP Antibody</u>	CSB-PA887964HA01HU	Human	ELISA, WB, IHC
<u>ADO Antibody</u>	CSB-PA839416LA01HU	Human	ELISA, IHC
<u>ADORA1 Antibody</u>	CSB-PA001375EA01HU	Human	ELISA, IF
<u>ADORA2A Antibody</u>	CSB-PA001376LA01HU	Human	ELISA, IHC, IF
<u>ADORA2B Antibody</u>	CSB-PA001377LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ADPRHL1 Antibody</u>	CSB-PA839860LA01HU	Human	ELISA
<u>ADRA1A Antibody</u>	CSB-PA001385PA01HU	Human	ELISA, WB, IF
<u>ADRA2A Antibody</u>	CSB-PA001388LA01HU	Human	ELISA, IF
<u>ADRB1 Antibody</u>	CSB-PA001391LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>Adrb2 Antibody</u>	CSB-PA001392LA01RA	Rat	ELISA
<u>ADRB3 Antibody</u>	CSB-PA001393LA01HU	Human	ELISA, IF
<u>ADRM1 Antibody</u>	CSB-PA001396HA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>AEBP1 Antibody</u>	CSB-PA814205LA01HU	Human	ELISA, IHC, IF
<u>AEN Antibody</u>	CSB-PA819883LA01HU	Human	ELISA, IF
<u>AFF2 Antibody</u>	CSB-PA001413LA01HU	Human	ELISA, WB
<u>AFF4 Antibody</u>	CSB-PA890687LA01HU	Human, Mouse	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>AFG1L Antibody</u>	CSB-PA819895LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>AFG3L2 Antibody</u>	CSB-PA04005A0Rb	Human	ELISA, WB, IHC
<u>Afp Antibody</u>	CSB-PA001421LA01MO	Mouse, Rat, Human	ELISA, WB
<u>AFP Antibody</u>	CSB-PA09987A0Rb	Human	ELISA, IHC
<u>AFP Antibody</u>	CSB-PA09989A0Rb	Human	ELISA, WB
<u>afp Antibody</u>	CSB-PA325933LA01APN	Aspergillus giganteus	ELISA
<u>AFP2 Antibody</u>	CSB-PA339081HA01RJP	Raphanus sativus	ELISA, WB
<u>AFTPH Antibody</u>	CSB-PA001422LA01HU	Human	ELISA, IHC
<u>AGA Antibody</u>	CSB-PA001423LA01HU	Human, Mouse	ELISA, WB
<u>AGA2 Antibody</u>	CSB-PA335766HA01SVG	Saccharomyces cerevisiae	ELISA
<u>AGBL2 Antibody</u>	CSB-PA732978LA01HU	Human, Mouse	ELISA, WB, IHC
<u>AGER Antibody</u>	CSB-PA14209A0Rb	Human	ELISA, IHC, IF
<u>AGK Antibody</u>	CSB-PA706634LA01HU	Human, Mouse	ELISA, WB, IHC
<u>AGL Antibody</u>	CSB-PA001446LA01HU	Human	ELISA, IF
<u>aglB Antibody</u>	CSB-PA381525LA01AVE	Aspergillus niger	ELISA
<u>aglB Antibody</u>	CSB-PA381525YA01AVE	Aspergillus niger	ELISA
<u>AGO2 Antibody</u>	CSB-PA891731LA01HU	Human, Mouse	ELISA, WB, IHC, ChIP
<u>AGPAT3 Antibody</u>	CSB-PA865125LA01HU	Human	ELISA, WB, IHC
<u>AGPAT5 Antibody</u>	CSB-PA001453LA01HU	Human, Mouse	ELISA, WB, IHC
<u>AGPS Antibody</u>	CSB-PA001457LA01HU	Human	ELISA, WB, IHC, IF
<u>AGR2 Antibody</u>	CSB-PA13057A0Rb	Human	ELISA, IHC
<u>AGR3 Antibody</u>	CSB-PA848404LA01HU	Human	ELISA
<u>AGRN Antibody</u>	CSB-PA001461LA01HU	Human	ELISA, WB, IHC, IF
<u>AGRP Antibody</u>	CSB-PA05699A0Rb	Human	ELISA, IHC
<u>AGT Antibody</u>	CSB-PA10079A0Rb	Human	ELISA, WB, IHC, IF
<u>AGT Antibody</u>	CSB-PA001463EA01HU	Human, Rat	ELISA, WB, IHC
<u>AGTR1 Antibody</u>	CSB-PA001465LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>AGTRAP Antibody</u>	CSB-PA744194LA01HU	Human	ELISA, IF
<u>Agxt Antibody</u>	CSB-PA001468LA01MO	Mouse	ELISA
<u>AGXT Antibody</u>	CSB-PA001468LA01HU	Human	ELISA, WB, IHC
<u>AGXT2 Antibody</u>	CSB-PA866322LA01HU	Human	ELISA, WB, IHC, IF
<u>AHCY Antibody</u>	CSB-PA001474HA01HU	Human	ELISA, IHC, IF
<u>AHCY Antibody</u>	CSB-PA001474LA01HU	Human	ELISA, IHC
<u>AHI1 Antibody</u>	CSB-PA836642LA01HU	Human	ELISA, WB, IHC
<u>AHNAK Antibody</u>	CSB-PA600257LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>AHNAK2 Antibody</u>	CSB-PA812872LA01HU	Human	ELISA, IHC, IF
<u>AHR Antibody</u>	CSB-PA001481LA01HU	Human	ELISA, IHC, IF
<u>AHRR Antibody</u>	CSB-PA001482LA01HU	Human	ELISA, IF
<u>AHSA1 Antibody</u>	CSB-PA05455A0Rb	Human	ELISA, WB, IHC
<u>AHSA2P Antibody</u>	CSB-PA741128LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>AHSG Antibody</u>	CSB-PA06369A0Rb	Human	ELISA, IHC, IF
<u>AICDA Antibody</u>	CSB-PA001487LA01HU	Human, Mouse	ELISA, WB
<u>AIDA Antibody</u>	CSB-PA839294LA01HU	Human	ELISA, IF
<u>Aif1 Antibody</u>	CSB-PA001490HA01MO	Mouse, Human	ELISA, WB
<u>Aif1 Antibody</u>	CSB-PA001490HA01RA	Rat	ELISA
<u>AIF1 Antibody</u>	CSB-PA00667A0Rb	Human	ELISA, IHC, IF
<u>AIFM1 Antibody</u>	CSB-PA001492HA01HU	Human	ELISA, WB, IHC, IF
<u>AIFM2 Antibody</u>	CSB-PA874794LA01HU	Human	ELISA, WB, IHC, IF
<u>aiiM Antibody</u>	CSB-PA512020HA01MTN	Microbacterium testaceum	ELISA
<u>AIM2 Antibody</u>	CSB-PA001499LA01HU	Human	ELISA, WB, IHC, IF
<u>AIMP1 Antibody</u>	CSB-PA618638EA01HU	Human	ELISA, WB, IHC, IF
<u>AIP Antibody</u>	CSB-PA516850LA01HU	Human	ELISA, IHC
<u>AIRE Antibody</u>	CSB-PA001502HA01HU	Human	ELISA, IHC, IF
<u>AJAP1 Antibody</u>	CSB-PA891947LA01HU	Human, Mouse	ELISA, WB, IHC
<u>AJUBA Antibody</u>	CSB-PA836229NA01HU	Human	ELISA, IHC
<u>AK Antibody</u>	CSB-PA511311HA01ETF	Penaeus monodon	ELISA, WB
<u>AK1 Antibody</u>	CSB-PA001508HA01HU	Human	ELISA, WB, IHC
<u>AK3 Antibody</u>	CSB-PA883438LA01HU	Human	ELISA, IHC
<u>AK5 Antibody</u>	CSB-PA897315LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>AKAP17A Antibody</u>	CSB-PA04299A0Rb	Human	ELISA, IF
<u>AKAP7 Antibody</u>	CSB-PA001529LA01HU	Human	ELISA, IHC,IF
<u>AKAP9 Antibody</u>	CSB-PA857882LA01HU	Human	ELISA, WB, IHC, IF
<u>AKIP1 Antibody</u>	CSB-PA868282LA01HU	Human	ELISA, IHC
<u>AKR1A1 Antibody</u>	CSB-PA001538LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>AKR1B1 Antibody</u>	CSB-PA001539LA01HU	Human	ELISA, WB, IHC, IF
<u>AKR1C1 Antibody</u>	CSB-PA001542LA01HU	Human	ELISA, IF
<u>AKR1C2 Antibody</u>	CSB-PA001543HA01HU	Human, Rat	ELISA, WB
<u>AKR1D1 Antibody</u>	CSB-PA001548LA01HU	Human	ELISA, IHC
<u>AKT1 Antibody</u>	CSB-PA15905A0Rb	Human	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF, IP
<u>AKT1S1 Antibody</u>	CSB-PA836185LA01HU	Human, Rat	ELISA, WB, IHC
<u>AKT2 Antibody</u>	CSB-PA15919A0Rb	Human, Rat	ELISA, WB, IHC, IF
<u>AKT3 Antibody</u>	CSB-PA15929A0Rb	Human, Mouse	ELISA, IP
<u>AKTIP Antibody</u>	CSB-PA862043LA01HU	Human	ELISA, IHC, IF
<u>alaA Antibody</u>	CSB-PA359250LA01ENV	Escherichia coli	ELISA, WB
<u>alaC Antibody</u>	CSB-PA302867HA01ENV	Escherichia coli	ELISA
<u>ALAS1 Antibody</u>	CSB-PA001559LA01HU	Human	ELISA, IHC
<u>ALAS2 Antibody</u>	CSB-PA001560LA01HU	Human	ELISA, WB, IHC
<u>Alb Antibody</u>	CSB-PA00010E1Rb	Rat	ELISA
<u>ALB Antibody</u>	CSB-PA00020E1Rb	Pigeon	ELISA
<u>ALB Antibody</u>	CSB-PA00030E1Rb	Chicken	ELISA
<u>ALB Antibody</u>	CSB-PA00040E1Rb	Horse	ELISA, WB
<u>ALB Antibody</u>	CSB-PA00050E1Rb	Sheep	ELISA
<u>ALB Antibody</u>	CSB-PA00060E1Rb	Human	ELISA
<u>ALB Antibody</u>	CSB-PA00070E0Gp	Rabbit	ELISA
<u>Alb Antibody</u>	CSB-PA00080E1Rb	Mouse	ELISA
<u>ALB Antibody</u>	CSB-PA001561JA01HU	Human	ELISA
<u>ALB Antibody</u>	CSB-PA001561LA01DO	Dog	ELISA
<u>ALB Antibody</u>	CSB-PA001561LA01MOV	Macaca fascicularis	ELISA
<u>ALB Antibody</u>	CSB-PA00180E1Rb	Goat	ELISA
<u>ALB Antibody</u>	CSB-PA00200E1Rb	Guinea Pig	ELISA
<u>ALB Antibody</u>	CSB-PA00210E1Rb	Pig	ELISA
<u>ALB Antibody</u>	CSB-PA00260E1Rb	Canine	ELISA
<u>ALB Antibody</u>	CSB-PA00950E1Rb	Bovine	ELISA
<u>ALB Antibody</u>	CSB-PA001561LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>Alb Antibody</u>	CSB-PA001561LA01RA	Rat	ELISA, WB
<u>albA Antibody</u>	CSB-PA840708HA01FHW	Pyrococcus furiosus	ELISA, WB
<u>ALCAM Antibody</u>	CSB-PA10209A0Rb	Human	ELISA, IHC, IF
<u>ald1 Antibody</u>	CSB-PA350554HA01SKX	Staphylococcus aureus	ELISA, WB
<u>ALDH16A1 Antibody</u>	CSB-PA811634LA01HU	Human	ELISA, IF
<u>Aldh1a1 Antibody</u>	CSB-PA001565LA01MO	Mouse	ELISA, WB
<u>ALDH1A2 Antibody</u>	CSB-PA001566EA01HU	Human, Mouse	ELISA, WB, IHC
<u>ALDH1A3 Antibody</u>	CSB-PA001567LA01HU	Human, Mouse	ELISA, WB, IF
<u>ALDH1L2 Antibody</u>	CSB-PA659640LA01HU	Human	ELISA, IHC
<u>ALDH3A1 Antibody</u>	CSB-PA001572EA01HU	Human	ELISA, WB, IHC, IF
<u>ALDH3B1 Antibody</u>	CSB-PA001574LA01HU	Human, Mouse, Rat	ELISA, WB, IF
<u>ALDH3B2 Antibody</u>	CSB-PA001575LA01HU	Human	ELISA, IHC
<u>ALDH5A1 Antibody</u>	CSB-PA001577LA01HU	Human, Mouse, Rat	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>ALDH6A1 Antibody</u>	CSB-PA001578LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ALDH7A1 Antibody</u>	CSB-PA001579LA01HU	Human, Rat, Mouse	ELISA, WB, IHC, IF, IP
<u>ALDH9A1 Antibody</u>	CSB-PA001582LA01HU	Human	ELISA, IHC, IF
<u>ALDOA Antibody</u>	CSB-PA00015A0Rb	Human, Rat	ELISA, WB, IHC
<u>ALDOC Antibody</u>	CSB-PA05659A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>ALG10 Antibody</u>	CSB-PA704720LA01HU	Human	ELISA, IHC
<u>ALG11 Antibody</u>	CSB-PA647192LA01HU	Human	ELISA, IHC, IF
<u>ALG12 Antibody</u>	CSB-PA887130LA01HU	Human	ELISA, IHC, IF
<u>ALG13 Antibody</u>	CSB-PA885677LA01HU	Human	ELISA, IF
<u>ALG14 Antibody</u>	CSB-PA850291LA01HU	Human	ELISA, IHC
<u>ALG2 Antibody</u>	CSB-PA863983LA01HU	Human	ELISA, IHC, IF
<u>ALG5 Antibody</u>	CSB-PA896924LA01HU	Human	ELISA, IHC
<u>ALG9 Antibody</u>	CSB-PA867161LA01HU	Human	ELISA, WB, IHC
<u>algL Antibody</u>	CSB-PA524422LA01DPX	Azotobacter vinelandii	ELISA, WB
<u>ALK Antibody</u>	CSB-PA21109A0Rb	Human	ELISA, IHC, IF
<u>ALKBH3 Antibody</u>	CSB-PA847209HA01HU	Human	ELISA, IHC, IF
<u>ALKBH4 Antibody</u>	CSB-PA885759LA01HU	Human, Mouse	ELISA, WB, IF
<u>Allergen Ara h 1, clone P17 Antibody</u>	CSB-PA331779LA01ANE	Peanut	ELISA
<u>Allergen Bla g 4 Antibody</u>	CSB-PA19857A0Rb	Blattella germanica	ELISA
<u>Allergen Fel d 4 Antibody</u>	CSB-PA19827A0Rb	Cat	ELISA
<u>ALMS1 Antibody</u>	CSB-PA844721LA01HU	Human	ELISA, IF
<u>ALMS1L Antibody</u>	CSB-PA853452LA01HU	Human	ELISA
<u>ALOX12 Antibody</u>	CSB-PA001617LA01HU	Human	ELISA, WB, IHC, IF
<u>ALOX12B Antibody</u>	CSB-PA001618LA01HU	Human	ELISA, IHC
<u>ALOX15 Antibody</u>	CSB-PA001621LA01HU	Human	ELISA, IHC, IF
<u>ALOX15B Antibody</u>	CSB-PA001622LA01HU	Human	ELISA, WB, IF
<u>ALOX5 Antibody</u>	CSB-PA001624LA01HU	Human	ELISA, IHC, IF
<u>Alpha-amanitin proprotein Antibody</u>	CSB-PA306453JO01AIF	Amanita phalloides	ELISA
<u>Alpha-amanitin proprotein Antibody</u>	CSB-PA306453JP01AIF	Amanita phalloides	ELISA
<u>Alpha-insect toxin LqhaIT Antibody</u>	CSB-PA322958LA01LDS	Leiurus quinquestriatus hebraeus	ELISA
<u>ALPI Antibody</u>	CSB-PA11989A0Rb	Human	ELISA, WB
<u>ALPL Antibody</u>	CSB-PA11999A0Rb	Human	ELISA, IHC
<u>Alpl Antibody</u>	CSB-PA17039A0Rb	Mouse	ELISA
<u>ALS2CR12 Antibody</u>	CSB-PA853485LA01HU	Human	ELISA, WB, IF



Product Name	Code	Species Reactivity	Applications
<u>ALTA1 Antibody</u>	CSB-PA303588LA01AZV	Alternaria alternata	ELISA
<u>ALX1 Antibody</u>	CSB-PA624031LA01HU	Human	ELISA, WB, IHC
<u>ALX3 Antibody</u>	CSB-PA001646LA01HU	Human	ELISA, IF
<u>ALX4 Antibody</u>	CSB-PA884437LA01HU	Human, Rat	ELISA, WB
<u>ALYREF Antibody</u>	CSB-PA774811LA01HU	Human	ELISA, IHC, IF
<u>AMA-1 Antibody</u>	CSB-PA338457HA01PLT	Plasmodium falciparum	ELISA
<u>AMACR Antibody</u>	CSB-PA001652LA01HU	Human, Mouse	ELISA, WB, IF
<u>AMBN Antibody</u>	CSB-PA882061LA01HU	Human, Rat	ELISA, WB
<u>AMBP Antibody</u>	CSB-PA001654YA01HU	Human	ELISA, WB, IHC
<u>AMBP Antibody</u>	CSB-PA10009A0Rb	Human	ELISA, IF
<u>AMD1 Antibody</u>	CSB-PA001568LA01HU	Human	ELISA, WB, IF
<u>AMDHD1 Antibody</u>	CSB-PA850402LA01HU	Human	ELISA
<u>AMFR Antibody</u>	CSB-PA333164LA01HU	Human	ELISA, WB, IHC
<u>Amh Antibody</u>	CSB-PA001666YA01MO	Mouse, Human	ELISA, WB, IHC
<u>AMH Antibody</u>	CSB-PA12639A0Rb	Human	ELISA, WB, IHC, IF
<u>amiF Antibody</u>	CSB-PA397462LA01BPF	Bradyrhizobium sp.	ELISA
<u>AMIGO1 Antibody</u>	CSB-PA773048LA01HU	Human	ELISA, IHC
<u>AMIGO2 Antibody</u>	CSB-PA769765LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>AMMECR1 Antibody</u>	CSB-PA22559A0Rb	Human	ELISA, IHC, IF
<u>AMN Antibody</u>	CSB-PA880128LA01HU	Human	ELISA, IHC, IF
<u>AMOTL2 Antibody</u>	CSB-PA809099NA01HU	Human	ELISA, WB
<u>ampC Antibody</u>	CSB-PA326492HA01EZX	Pseudomonas aeruginosa	ELISA
<u>AMPD1 Antibody</u>	CSB-PA001680LA01HU	Human	ELISA, WB
<u>AMPD2 Antibody</u>	CSB-PA001681LA01HU	Human	ELISA, WB, IHC, IF
<u>AMT Antibody</u>	CSB-PA001684LA01HU	Human, Mouse	ELISA, WB
<u>AMTN Antibody</u>	CSB-PA001685LA01HU	Human	ELISA, IHC
<u>ANAPC5 Antibody</u>	CSB-PA05104A0Rb	Human	ELISA, IHC
<u>ANG Antibody</u>	CSB-PA05729A0Rb	Human	ELISA, IHC
<u>ANG1 Antibody</u>	CSB-PA14679A0Rb	Bovine	ELISA
<u>ANGPT1 Antibody</u>	CSB-PA05739A0Rb	Human	ELISA, WB, IHC
<u>ANGPT2 Antibody</u>	CSB-PA05749A0Rb	Human, Mouse, Rat	ELISA, WB, IHC
<u>Angpt2 Antibody</u>	CSB-PA16959A0Rb	Mouse	ELISA, WB, IHC
<u>ANGPTL1 Antibody</u>	CSB-PA09689A0Rb	Human	ELISA, IHC
<u>ANGPTL2 Antibody</u>	CSB-PA891539LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ANGPTL2 Antibody</u>	CSB-PA891539YA01HU	Human	ELISA, IHC
<u>ANGPTL3 Antibody</u>	CSB-PA05759A0Rb	Human	ELISA, IHC
<u>ANGPTL4 Antibody</u>	CSB-PA866314EA01HU	Human	ELISA, WB
<u>ANGPTL7 Antibody</u>	CSB-PA001715LA01HU	Human	ELISA, WB, IHC
<u>ANGPTL8 Antibody</u>	CSB-PA757793LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>Angptl8 Antibody</u>	CSB-PA844436LA01MO	Mouse	ELISA
<u>Anionic trypsin Antibody</u>	CSB-PA356954LA01DO	Dog	ELISA
<u>ANK3 Antibody</u>	CSB-PA619642LA01HU	Human	ELISA, IHC, IF
<u>ANKDD1A Antibody</u>	CSB-PA669849LA01HU	Human	ELISA, IHC
<u>ANKMY1 Antibody</u>	CSB-PA885806LA01HU	Human	ELISA, IHC
<u>ANKMY2 Antibody</u>	CSB-PA815551LA01HU	Human	ELISA, IHC
<u>ANKRA2 Antibody</u>	CSB-PA872524LA01HU	Human	ELISA, IHC
<u>ANKRD10 Antibody</u>	CSB-PA878915LA01HU	Human, Mouse	ELISA, WB
<u>ANKRD11 Antibody</u>	CSB-PA757769LA01HU	Human	ELISA
<u>ANKRD2 Antibody</u>	CSB-PA875635LA01HU	Human	ELISA, IF
<u>ANKRD23 Antibody</u>	CSB-PA803108LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ANKRD33 Antibody</u>	CSB-PA768761LA01HU	Human	ELISA, IHC
<u>ANKRD42 Antibody</u>	CSB-PA847661LA01HU	Human	ELISA, IHC
<u>ANKRD49 Antibody</u>	CSB-PA837858LA01HU	Human	ELISA, WB, IHC
<u>ANKRD55 Antibody</u>	CSB-PA660753LA01HU	Human	ELISA, IHC, IF
<u>ANKRD7 Antibody</u>	CSB-PA842127LA01HU	Human	ELISA
<u>ANKS6 Antibody</u>	CSB-PA718074LA01HU	Human	ELISA, IHC, IF
<u>Annexin A5 Antibody</u>	CSB-PA001846HA01EQC	Cynops pyrrhogaster	ELISA
<u>ANO1 Antibody</u>	CSB-PA804805LA01HU	Human	ELISA, IHC, IF
<u>ANO2 Antibody</u>	CSB-PA885689LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>ANO4 Antibody</u>	CSB-PA653477LA01HU	Human	ELISA, IHC, IF
<u>ANO6 Antibody</u>	CSB-PA679761LA01HU	Human	ELISA, IF
<u>ANP32D Antibody</u>	CSB-PA001825LA01HU	Human	ELISA, IHC, IF
<u>ANPEP Antibody</u>	CSB-PA001827LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>Anxa1 Antibody</u>	CSB-PA001836HA01MO	Mouse, Human	ELISA, WB
<u>Anxa1 Antibody</u>	CSB-PA001836YA01MO	Human, Mouse	ELISA, WB, IP
<u>ANXA1 Antibody</u>	CSB-PA13457A0Rb	Human	ELISA, WB, IF, IP
<u>ANXA10 Antibody</u>	CSB-PA001837LA01HU	Human	ELISA, IHC, IF
<u>ANXA13 Antibody</u>	CSB-PA001839HA01DO	Dog	ELISA, WB
<u>ANXA2 Antibody</u>	CSB-PA001840HA01HU	Human, Mouse	ELISA, WB, IHC, IF, IP
<u>ANXA2 Antibody</u>	CSB-PA01064A0Rb	Human	ELISA, WB
<u>ANXA3 Antibody</u>	CSB-PA001844YA01HU	Human	ELISA, IHC
<u>ANXA3 Antibody</u>	CSB-PA17149A0Rb	Human	ELISA
<u>ANXA4 Antibody</u>	CSB-PA001845LA01DO	Dog	ELISA
<u>ANXA5 Antibody</u>	CSB-PA001846YA01HU	Human	ELISA, WB, IHC
<u>ANXA5 Antibody</u>	CSB-PA06384A0Rb	Human, Mouse	ELISA, WB, IHC
<u>ANXA6 Antibody</u>	CSB-PA00794A0Rb	Human, Rat	ELISA, WB
<u>ANXA9 Antibody</u>	CSB-PA001852LA01HU	Human, Mouse	ELISA, WB, IHC
<u>AOAH Antibody</u>	CSB-PA001853LA01HU	Human, Mouse, Rat	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>AOC3 Antibody</u>	CSB-PA624122LA01HU	Human	ELISA, IHC, IF
<u>AOX1 Antibody</u>	CSB-PA001858HA01HU	Human	ELISA, IHC, IF
<u>AOX1 Antibody</u>	CSB-PA24692A0Rb	Chlamydomonas reinhardtii	ELISA, WB
<u>AOX1A Antibody</u>	CSB-PA655106LA01DOA	Arabidopsis thaliana	ELISA, WB
<u>APIG1 Antibody</u>	CSB-PA001861LA01HU	Human	ELISA, IHC, IF
<u>APIS1 Antibody</u>	CSB-PA001866LA01HU	Human	ELISA
<u>APIS3 Antibody</u>	CSB-PA001868EA01HU	Human	ELISA, WB
<u>AP24 Antibody</u>	CSB-PA320333LA01NHE	Nicotiana tabacum	ELISA, WB
<u>AP2M1 Antibody</u>	CSB-PA00625A0Rb	Human	ELISA, WB, IHC
<u>AP2M1 Antibody</u>	CSB-PA00627A0Rb	Human, Mouse, Rat	ELISA, WB, IHC
<u>AP2S1 Antibody</u>	CSB-PA02604A0Rb	Human	ELISA, IHC, IF
<u>AP3M2 Antibody</u>	CSB-PA001878LA01HU	Human	ELISA, IHC
<u>AP3S2 Antibody</u>	CSB-PA001880LA01HU	Human	ELISA, WB
<u>AP4B1 Antibody</u>	CSB-PA897585LA01HU	Human	ELISA, WB, IHC
<u>APAF1 Antibody</u>	CSB-PA001885LA01HU	Human	ELISA, IHC, IF
<u>apaLIM Antibody</u>	CSB-PA523812LA01AVW	Acetobacter pasteurianus	ELISA
<u>APC2 Antibody</u>	CSB-PA001895LA01HU	Human	ELISA, IF
<u>APCDD1 Antibody</u>	CSB-PA814246LA01HU	Human	ELISA, WB, IHC
<u>APCS Antibody</u>	CSB-PA001898LA01HU	Human	ELISA, IHC
<u>APEH Antibody</u>	CSB-PA001899LA01HU	Human	ELISA, WB, IHC, IF
<u>APEX1 Antibody</u>	CSB-PA001900HA01HU	Human, Mouse	ELISA, WB, IHC, ChIP
<u>APH1A Antibody</u>	CSB-PA853390LA01HU	Human	ELISA, IF
<u>APH1B Antibody</u>	CSB-PA819900NA01HU	Human	ELISA, WB
<u>API5 Antibody</u>	CSB-PA05559A0Rb	Human, Mouse	ELISA, WB, IF
<u>APLN Antibody</u>	CSB-PA001909NJ01HU	Human	ELISA, IHC
<u>APLN Antibody</u>	CSB-PA001909OJ01HU	Human	ELISA, IHC
<u>APLNR Antibody</u>	CSB-PA001910EA01HU	Human	ELISA, IHC, IF
<u>APMAP Antibody</u>	CSB-PA884629LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>APOA1 Antibody</u>	CSB-PA001913HA01HU	Human	ELISA
<u>APOA1 Antibody</u>	CSB-PA001913LA01BO	Bovine, Rat, Mouse	ELISA, WB
<u>Apoa4 Antibody</u>	CSB-PA001916DA01MO	Mouse	ELISA
<u>APOA4 Antibody</u>	CSB-PA001916HA01HU	Human	ELISA, WB
<u>APOA5 Antibody</u>	CSB-PA001917LA01HU	Human	ELISA, IHC
<u>APOB Antibody</u>	CSB-PA001918LA01HU	Human	ELISA
<u>APOBEC1 Antibody</u>	CSB-PA001919LA01HU	Human	ELISA, IF
<u>APOBEC3A Antibody</u>	CSB-PA001921LA01HU	Human	ELISA
<u>APOBEC3D Antibody</u>	CSB-PA846585LA01HU	Human	ELISA, IHC, IF
<u>APOBEC3F Antibody</u>	CSB-PA816878LA01HU	Human, Mouse	ELISA, WB, IF



Product Name	Code	Species Reactivity	Applications
<u>APOBEC3G Antibody</u>	CSB-PA001926LA01HU	Human	ELISA, IHC, IF
<u>APOBEC4 Antibody</u>	CSB-PA823894LA01HU	Human	ELISA, IHC, IF
<u>APOC2 Antibody</u>	CSB-PA001932OA01HU	Human	ELISA, IHC
<u>APOC3 Antibody</u>	CSB-PA12197A0Rb	Human, Rat	ELISA, WB, IHC
<u>APOD Antibody</u>	CSB-PA001935EA01HU	Human	ELISA, IHC
<u>ApoE Antibody</u>	CSB-PA001936LA01MO	Mouse	ELISA, WB
<u>APOE Antibody</u>	CSB-PA07599A0Rb	Human	ELISA, IHC, IF
<u>APOH Antibody</u>	CSB-PA11079A0Rb	Human, Rat	ELISA, WB, IHC, IF
<u>APOLD1 Antibody</u>	CSB-PA853458LA01HU	Human	ELISA, IF
<u>APOOL Antibody</u>	CSB-PA001949LA01HU	Human	ELISA, IHC, IF
<u>APOPT1 Antibody</u>	CSB-PA850312LA01HU	Human, Mouse	ELISA, WB
<u>APP Antibody</u>	CSB-PA001950NJ01HU	Human	ELISA, WB
<u>APP Antibody</u>	CSB-PA001950OJ01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>APP Antibody</u>	CSB-PA07174A0Rb	Human	ELISA, IHC
<u>APP Antibody</u>	CSB-PA07179A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>APPL1 Antibody</u>	CSB-PA890917LA01HU	Human	ELISA, IHC
<u>APPL2 Antibody</u>	CSB-PA818777LA01HU	Human	ELISA, IHC, IF
<u>APT X Antibody</u>	CSB-PA001955LA01HU	Human	ELISA, IHC
<u>APX1 Antibody</u>	CSB-PA13479A0Rb	Arabidopsis thaliana	ELISA, WB
<u>APX2 Antibody</u>	CSB-PA14419A0Rb	Arabidopsis thaliana	ELISA, WB
<u>AQP1 Antibody</u>	CSB-PA16085A0Rb	Human, Mouse, Rat	ELISA, WB, IF, IP
<u>AQP1 Antibody</u>	CSB-PA16087A0Rb	Human	ELISA, IHC, IF
<u>AQP10 Antibody</u>	CSB-PA853482LA01HU	Human	ELISA, IHC, IF
<u>AQP2 Antibody</u>	CSB-PA16069A0Rb	Human	ELISA, IF
<u>AQP3 Antibody</u>	CSB-PA13909A0Rb	Human	ELISA, IHC, IF
<u>AQP4 Antibody</u>	CSB-PA11137A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>AQP6 Antibody</u>	CSB-PA16139A0Rb	Human	ELISA, IHC
<u>AR Antibody</u>	CSB-PA05979A0Rb	Human	ELISA, IHC
<u>Arachin Ahy-3 Antibody</u>	CSB-PA737391LA01ANE	Peanut	ELISA
<u>ARAF Antibody</u>	CSB-PA01679A0Rb	Human	ELISA, IHC, IF
<u>ARAP1 Antibody</u>	CSB-PA856998LA01HU	Human	ELISA, IHC, IF
<u>ARAP3 Antibody</u>	CSB-PA819904LA01HU	Human	ELISA, IHC
<u>AREG Antibody</u>	CSB-PA11339A0Rb	Human	ELISA, IHC, IF
<u>ARF1 Antibody</u>	CSB-PA00929A0Rb	Human	ELISA, IHC, IF
<u>ARF4 Antibody</u>	CSB-PA00915A0Rb	Human	ELISA, WB, IHC
<u>ARF5 Antibody</u>	CSB-PA00905A0Rb	Human	ELISA
<u>ARF6 Antibody</u>	CSB-PA001996LA01HU	Human	ELISA, WB, IHC
<u>ARFGEF2 Antibody</u>	CSB-PA896929LA01HU	Human, Rat, Mouse	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>ARFIP2 Antibody</u>	CSB-PA002003LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ARFRP1 Antibody</u>	CSB-PA621696LA01HU	Human	ELISA, IHC
<u>Arg1 Antibody</u>	CSB-PA733775HA01MO	Mouse	ELISA, IP
<u>ARG1 Antibody</u>	CSB-PA002005DA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>Arg2 Antibody</u>	CSB-PA002006LA01MO	Mouse, Human	ELISA, WB
<u>argF Antibody</u>	CSB-PA12907A0Rb	Escherichia coli	ELISA
<u>Arginine esterase Antibody</u>	CSB-PA362734LA01DO	Dog	ELISA, WB
<u>ARHGAP1 Antibody</u>	CSB-PA691829LA01HU	Human, Rat	ELISA, WB, IHC
<u>ARHGAP10 Antibody</u>	CSB-PA002012LA01HU	Human	ELISA, IHC, IF
<u>ARHGAP11A Antibody</u>	CSB-PA750800LA01HU	Human	ELISA
<u>ARHGAP11B Antibody</u>	CSB-PA664668LA01HU	Human	ELISA
<u>ARHGAP15 Antibody</u>	CSB-PA681513LA01HU	Human	ELISA, IHC, IF
<u>ARHGAP18 Antibody</u>	CSB-PA850791LA01HU	Human, Rat, Mouse	ELISA, WB, IHC, IF
<u>ARHGAP19 Antibody</u>	CSB-PA623954LA01HU	Human	ELISA, IHC
<u>ARHGAP25 Antibody</u>	CSB-PA002026LA01HU	Human	ELISA, IHC
<u>ARHGAP27 Antibody</u>	CSB-PA741070LA01HU	Human	ELISA, IHC
<u>ARHGAP29 Antibody</u>	CSB-PA681475LA01HU	Human	ELISA
<u>ARHGAP32 Antibody</u>	CSB-PA415667LA01HU	Human	ELISA, IF
<u>ARHGAP42 Antibody</u>	CSB-PA413088LA01HU	Human	ELISA, WB, IHC
<u>ARHGDIA Antibody</u>	CSB-PA00085A0Rb	Human	ELISA, WB, IHC
<u>ARHGDIB Antibody</u>	CSB-PA002038LA01HU	Human	ELISA, WB, IF
<u>ARHGDIG Antibody</u>	CSB-PA858727LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ARHGEF11 Antibody</u>	CSB-PA002043LA01HU	Human	ELISA, IHC
<u>ARHGEF12 Antibody</u>	CSB-PA868370LA01HU	Human	ELISA, IHC, IF
<u>ARHGEF18 Antibody</u>	CSB-PA019681HA01HU	Human	ELISA, IHC, IF
<u>ARHGEF18 Antibody</u>	CSB-PA05275A0Rb	Human	ELISA, WB
<u>ARHGEF39 Antibody</u>	CSB-PA854040LA01HU	Human	ELISA
<u>ARHGEF5 Antibody</u>	CSB-PA618628LA01HU	Human	ELISA, IHC
<u>ARHGEF6 Antibody</u>	CSB-PA618974LA01HU	Human	ELISA, IHC
<u>ARHGEF7 Antibody</u>	CSB-PA01445A0Rb	Human	ELISA, IF, IP
<u>ARID1A Antibody</u>	CSB-PA002058LA01HU	Human	ELISA, IHC
<u>ARID3A Antibody</u>	CSB-PA858728LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ARIH2 Antibody</u>	CSB-PA002069LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ARL1 Antibody</u>	CSB-PA002070LA01HU	Human	ELISA, WB, IHC, IF
<u>ARL11 Antibody</u>	CSB-PA850249LA01HU	Human	ELISA, IHC, IF
<u>ARL14 Antibody</u>	CSB-PA836664LA01HU	Human	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>ARL15 Antibody</u>	CSB-PA882136LA01HU	Human	ELISA, IHC
<u>ARL16 Antibody</u>	CSB-PA612678LA01HU	Human	ELISA, IHC
<u>ARL2 Antibody</u>	CSB-PA01265A0Rb	Human	ELISA, WB
<u>ARL2BP Antibody</u>	CSB-PA897097HA01HU	Human	ELISA, WB, IHC
<u>ARL4A Antibody</u>	CSB-PA002085HA01HU	Human	ELISA
<u>ARL4C Antibody</u>	CSB-PA002086LA01HU	Human	ELISA, IHC, IF
<u>ARL4D Antibody</u>	CSB-PA002087LA01HU	Human	ELISA, IHC
<u>ARL6IP1 Antibody</u>	CSB-PA614878LA01HU	Human	ELISA, IHC, IF
<u>ARL6IP5 Antibody</u>	CSB-PA002096LA01HU	Human	ELISA, IHC, IF
<u>ARL8A Antibody</u>	CSB-PA842642LA01HU	Human	ELISA
<u>ARL8B Antibody</u>	CSB-PA868330LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ARMC10 Antibody</u>	CSB-PA854026LA01HU	Human	ELISA, IHC, IF
<u>ARMC7 Antibody</u>	CSB-PA872505LA01HU	Human	ELISA, IHC
<u>ARMC8 Antibody</u>	CSB-PA816877LA01HU	Human	ELISA
<u>ARMC9 Antibody</u>	CSB-PA773583LA01HU	Human	ELISA, WB, IHC
<u>ARMCX1 Antibody</u>	CSB-PA878960LA01HU	Human	ELISA, IHC
<u>ARMCX2 Antibody</u>	CSB-PA755485LA01HU	Human	ELISA, WB, IHC
<u>ARMCX6 Antibody</u>	CSB-PA765977LA01HU	Human	ELISA, WB, IHC
<u>ARMT1 Antibody</u>	CSB-PA867178LA01HU	Human	ELISA, WB, IHC
<u>ARNTL Antibody</u>	CSB-PA002123LA01HU	Human	ELISA, IF
<u>aroA Antibody</u>	CSB-PA870830HA01ALAE	Agrobacterium sp.	ELISA
<u>aroD Antibody</u>	CSB-PA326162HA01SWW	Salmonella typhi	ELISA
<u>aroD Antibody</u>	CSB-PA326162YA01SWW	Salmonella typhi	ELISA
<u>ARPC1A Antibody</u>	CSB-PA856421LA01HU	Human	ELISA, WB, IF
<u>ARPC1B Antibody</u>	CSB-PA002126LA01HU	Human	ELISA, IHC
<u>ARPC2 Antibody</u>	CSB-PA03365A0Rb	Human	ELISA, IHC
<u>ARPC3 Antibody</u>	CSB-PA03324A0Rb	Human	ELISA, WB, IHC, IF
<u>ARPC4 Antibody</u>	CSB-PA002130LA01HU	Human	ELISA, IHC, IF
<u>ARPC5L Antibody</u>	CSB-PA861078LA01HU	Human	ELISA, IF
<u>ARPIN Antibody</u>	CSB-PA773801LA01HU	Human	ELISA, WB, IHC
<u>ARPP19 Antibody</u>	CSB-PA348885LA01HU	Human	ELISA, IHC, IF
<u>ARRB1 Antibody</u>	CSB-PA002134LA01HU	Human	ELISA, WB, IHC, IF
<u>ARRB2 Antibody</u>	CSB-PA002135LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>ARRDC1 Antibody</u>	CSB-PA850807LA01HU	Human	ELISA, IHC, IF
<u>ARRDC1-AS1 Antibody</u>	CSB-PA887963LA01HU	Human	ELISA, IHC
<u>ARRDC3 Antibody</u>	CSB-PA839293LA01HU	Human, Mouse	ELISA, WB
<u>ARSB Antibody</u>	CSB-PA002142LA01HU	Human, Mouse	ELISA, WB
<u>ARSE Antibody</u>	CSB-PA002146LA01HU	Human	ELISA, IHC, IF
<u>ARSG Antibody</u>	CSB-PA002150LA01HU	Human	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>ARSH Antibody</u>	CSB-PA688880LA01HU	Human	ELISA, WB, IF
<u>ARSI Antibody</u>	CSB-PA685925LA01HU	Human	ELISA, IHC
<u>ART1 Antibody</u>	CSB-PA002155LA01HU	Human	ELISA, IHC
<u>ART3 Antibody</u>	CSB-PA622665LA01HU	Human	ELISA, IHC, IF
<u>ART5 Antibody</u>	CSB-PA836246LA01HU	Human	ELISA, IHC
<u>ARTN Antibody</u>	CSB-PA002160LA01HU	Human	ELISA, IHC
<u>ARV1 Antibody</u>	CSB-PA867134EA01HU	Human	ELISA, IF
<u>ARX Antibody</u>	CSB-PA850418LA01HU	Human, Rat, Mouse	ELISA, WB, IHC
<u>AS3MT Antibody</u>	CSB-PA888021LA01HU	Human	ELISA, IHC, IF
<u>ASAH1 Antibody</u>	CSB-PA619774LA01HU	Human	ELISA, WB, IHC
<u>ASAH2 Antibody</u>	CSB-PA002170EA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>ASAP3 Antibody</u>	CSB-PA837444LA01HU	Human	ELISA, WB, IHC
<u>ASB10 Antibody</u>	CSB-PA845168LA01HU	Human	ELISA, IHC, IF
<u>ASB11 Antibody</u>	CSB-PA855501LA01HU	Human	ELISA, IHC
<u>ASB13 Antibody</u>	CSB-PA837876LA01HU	Human	ELISA, IHC
<u>ASB4 Antibody</u>	CSB-PA897287LA01HU	Human, Mouse, Rat	ELISA, WB
<u>ASB6 Antibody</u>	CSB-PA868341LA01HU	Human	ELISA, IHC
<u>ASB7 Antibody</u>	CSB-PA863988LA01HU	Human	ELISA, WB, IHC
<u>ASB8 Antibody</u>	CSB-PA863994LA01HU	Human	ELISA, WB, IHC
<u>ASB9 Antibody</u>	CSB-PA853405LA01HU	Human	ELISA
<u>ASCC1 Antibody</u>	CSB-PA002196LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ASCL1 Antibody</u>	CSB-PA002199LA01HU	Human, Rat	ELISA, WB, IHC
<u>ascl1a Antibody</u>	CSB-PA852654LA01DIL	Zebrafish	ELISA
<u>ASF1B Antibody</u>	CSB-PA878898LA01HU	Human	ELISA, WB, IHC
<u>ASGR1 Antibody</u>	CSB-PA17469A0Rb	Human	ELISA, IHC, IF
<u>ASGR2 Antibody</u>	CSB-PA002208LA01HU	Human	ELISA, IHC, IF
<u>ASH1L Antibody</u>	CSB-PA865114LA01HU	Human	ELISA, IHC, IF
<u>ASH2L Antibody</u>	CSB-PA883363LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ASIC1 Antibody</u>	CSB-PA001142LA01HU	Human, Mouse	ELISA, WB, IF
<u>ASIC2 Antibody</u>	CSB-PA613683LA01HU	Human	ELISA, IHC, IF
<u>ASIP Antibody</u>	CSB-PA002212LA01HU	Human	ELISA, IHC, IF
<u>ASL Antibody</u>	CSB-PA002213LA01HU	Human	ELISA, IHC, IF
<u>ASMT Antibody</u>	CSB-PA002216LA01HU	Human	ELISA, IF
<u>ASNS Antibody</u>	CSB-PA002219LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>ASNSD1 Antibody</u>	CSB-PA885745LA01HU	Human	ELISA, IHC
<u>aspf3 Antibody</u>	CSB-PA525004HA01NGS	Neosartorya fumigata	ELISA, WB
<u>ASPH Antibody</u>	CSB-PA002226HA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ASPHD1 Antibody</u>	CSB-PA716402LA01HU	Human	ELISA, IHC
<u>ASPHD2 Antibody</u>	CSB-PA760787LA01HU	Human	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>ASPM Antibody</u>	CSB-PA814245LA01HU	Human	ELISA, IHC, IF
<u>ASPN Antibody</u>	CSB-PA002230LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ASS1 Antibody</u>	CSB-PA002234LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ASTL Antibody</u>	CSB-PA760708LA01HU	Human	ELISA, IHC
<u>ASTN2 Antibody</u>	CSB-PA002252LA01HU	Human	ELISA, WB, IHC, IF
<u>ASZ1 Antibody</u>	CSB-PA845159LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>At1g32060 Antibody</u>	CSB-PA326237OA01DOA	Arabidopsis thaliana	ELISA, WB
<u>At2g24580 Antibody</u>	CSB-PA18349A0RB	Arabidopsis thaliana	ELISA, WB
<u>ATAD1 Antibody</u>	CSB-PA850855LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>ATAD3C Antibody</u>	CSB-PA716349LA01HU	Human	ELISA, IHC
<u>ATAT1 Antibody</u>	CSB-PA725542LA01HU	Human, Rat, Mouse	ELISA, WB, IHC
<u>ATE1 Antibody</u>	CSB-PA002268LA01HU	Human	ELISA, IHC, IF
<u>ATF3 Antibody</u>	CSB-PA002271EA01HU	Human	ELISA, WB, IHC
<u>ATF5 Antibody</u>	CSB-PA002276LA01HU	Human	ELISA, IHC
<u>ATG1 Antibody</u>	CSB-PA738902LA01CZI	Candida glabrata	ELISA
<u>ATG10 Antibody</u>	CSB-PA861988LA01HU	Human	ELISA, IF
<u>ATG101 Antibody</u>	CSB-PA883590LA01HU	Human	ELISA, IF
<u>ATG13 Antibody</u>	CSB-PA012203LA01HU	Human	ELISA, WB, IHC
<u>ATG14 Antibody</u>	CSB-PA741052LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>ATG16L1 Antibody</u>	CSB-PA734569LA01HU	Human	ELISA, IHC, IF
<u>ATG16L2 Antibody</u>	CSB-PA843310HA01HU	Human	ELISA, IHC, IF
<u>ATG2B Antibody</u>	CSB-PA822186LA01HU	Human	ELISA, IHC, IF
<u>ATG3 Antibody</u>	CSB-PA002288HA01HU	Human	ELISA, WB, IHC, IF
<u>ATG4A Antibody</u>	CSB-PA845174EA01HU	Human, Mouse	ELISA, WB
<u>ATG4B Antibody</u>	CSB-PA002290LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ATG4C Antibody</u>	CSB-PA836199LA01HU	Human	ELISA, IHC, IF
<u>ATG4D Antibody</u>	CSB-PA774797LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ATG5 Antibody</u>	CSB-PA15499A0Rb	Human	ELISA, IHC, IF
<u>ATG7 Antibody</u>	CSB-PA002294LA01HU	Human	ELISA, WB, IHC, IF
<u>ATG8 Antibody</u>	CSB-PA739788HA01CZI	Candida glabrata	ELISA
<u>ATIC Antibody</u>	CSB-PA002299LA01HU	Human	ELISA, IHC, IF
<u>ATL1 Antibody</u>	CSB-PA855500LA01HU	Human	ELISA, IHC
<u>ATL2 Antibody</u>	CSB-PA836776LA01HU	Human	ELISA, WB, IHC, IF
<u>ATL3 Antibody</u>	CSB-PA728503LA01HU	Human	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>Atoh1 Antibody</u>	CSB-PA002306LA01MO	Mouse	ELISA, WB
<u>ATOH1 Antibody</u>	CSB-PA846107LA01HU	Human	ELISA, WB, IHC
<u>ATOX1 Antibody</u>	CSB-PA002309LA01HU	Human	ELISA, IHC
<u>ATP12A Antibody</u>	CSB-PA002316LA01HU	Human	ELISA, IHC, IF
<u>ATP13A2 Antibody</u>	CSB-PA23379A0Rb	Human	ELISA, IF
<u>ATP1A1 Antibody</u>	CSB-PA18129A0Rb	Human	ELISA, IHC, IF
<u>ATP1A3 Antibody</u>	CSB-PA18139A0Rb	Human	ELISA
<u>ATP1B1 Antibody</u>	CSB-PA002326LA01HU	Human, Rat, Mouse	ELISA, WB, IHC, IF
<u>ATP1B3 Antibody</u>	CSB-PA002328LA01HU	Human	ELISA, IHC, IF
<u>ATP2A1 Antibody</u>	CSB-PA002332LA01HU	Human	ELISA, IHC, IF
<u>ATP2A2 Antibody</u>	CSB-PA002333LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ATP2B3 Antibody</u>	CSB-PA613698LA01HU	Human	ELISA, WB, IF
<u>ATP2B4 Antibody</u>	CSB-PA23389A0Rb	Human	ELISA, IHC, IF
<u>ATP2C1 Antibody</u>	CSB-PA002339LA01HU	Human	ELISA, IHC, IF
<u>ATP2C2 Antibody</u>	CSB-PA002340LA01HU	Human	ELISA, IF
<u>ATP4B Antibody</u>	CSB-PA002343EA01HU	Human	ELISA
<u>ATP5B Antibody</u>	CSB-PA002350EA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>ATP5F1A Antibody</u>	CSB-PA002344LA01HU	Human	ELISA, WB
<u>ATP5F1A Antibody</u>	CSB-PA002344HA01HU	Human, Mouse	ELISA, WB, IHC, IP
<u>Atp5f1b Antibody</u>	CSB-PA002350LA01MO	Mouse	ELISA
<u>ATP5F1C Antibody</u>	CSB-PA002353LA01HU	Human	ELISA, IHC
<u>ATP5F1D Antibody</u>	CSB-PA002355HA01HU	Human	ELISA, IHC, IF
<u>ATP5F1D Antibody</u>	CSB-PA00635A0Rb	Human	ELISA, WB, IHC
<u>ATP5L Antibody</u>	CSB-PA002374LA01HU	Human	ELISA, IHC
<u>ATP5ME Antibody</u>	CSB-PA05499A0Rb	Human	ELISA
<u>ATP6AP1L Antibody</u>	CSB-PA681474NA01HU	Human	ELISA, IHC, IF
<u>ATP6AP2 Antibody</u>	CSB-PA002384LA01HU	Human, Mouse	ELISA, WB
<u>ATP6V0A1 Antibody</u>	CSB-PA846292LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>ATP6V0A2 Antibody</u>	CSB-PA897489LA01HU	Human	ELISA, IF
<u>ATP6V0D2 Antibody</u>	CSB-PA818735LA01HU	Human	ELISA, WB, IHC, IF
<u>ATP6V1A Antibody</u>	CSB-PA002396LA01HU	Human	ELISA, IHC, IF
<u>ATP6V1C1 Antibody</u>	CSB-PA002399LA01HU	Human	ELISA, IHC, IF
<u>ATP6V1C2 Antibody</u>	CSB-PA822777LA01HU	Human, Mouse, Rat	ELISA, WB, IF
<u>ATP6V1D Antibody</u>	CSB-PA897507LA01HU	Human	ELISA, WB
<u>ATP6V1E1 Antibody</u>	CSB-PA02295A0Rb	Human	ELISA, WB, IF
<u>ATP6V1E2 Antibody</u>	CSB-PA836175LA01HU	Human	ELISA, IF



Product Name	Code	Species Reactivity	Applications
<u>ATP6V1F Antibody</u>	CSB-PA01574A0Rb	Human	ELISA, IHC, IF
<u>ATP6V1G1 Antibody</u>	CSB-PA00135A0Rb	Human	ELISA, WB, IHC, IF
<u>ATP6V1G2 Antibody</u>	CSB-PA002409LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>ATP7A Antibody</u>	CSB-PA002414LA01HU	Human	ELISA, IHC, IF
<u>ATP8A1 Antibody</u>	CSB-PA896697LA01HU	Human	ELISA, IHC, IF
<u>ATP8A2 Antibody</u>	CSB-PA865134LA01HU	Human	ELISA, IHC, IF
<u>ATP8B1 Antibody</u>	CSB-PA002418LA01HU	Human	ELISA, IF
<u>ATP9A Antibody</u>	CSB-PA002422LA01HU	Human	ELISA, WB, IHC, IF
<u>ATPAF1 Antibody</u>	CSB-PA002424LA01HU	Human	ELISA, IHC
<u>ATPAF2 Antibody</u>	CSB-PA002425LA01HU	Human	ELISA
<u>atpB Antibody</u>	CSB-PA13999A0Rb	Lemna minor	ELISA
<u>atpB Antibody</u>	CSB-PA20469A0Rb	Arabidopsis thaliana	ELISA
<u>atpE Antibody</u>	CSB-PA428118LA01AUQ	strain OhILAs	ELISA
<u>ATR Antibody</u>	CSB-PA622666LA01HU	Human	ELISA, WB, IHC
<u>ATRX Antibody</u>	CSB-PA002437LA01HU	Human	ELISA, IHC, IF
<u>ATXN1L Antibody</u>	CSB-PA002440LA01HU	Human	ELISA, IHC
<u>ATXN2 Antibody</u>	CSB-PA859520LA01HU	Human	ELISA, WB, IF
<u>ATXN3 Antibody</u>	CSB-PA002443LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ATXN7 Antibody</u>	CSB-PA002445LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>ATXN7L1 Antibody</u>	CSB-PA892147LA01HU	Human	ELISA, IHC
<u>ATXN7L3 Antibody</u>	CSB-PA613538LA01HU	Human	ELISA, IHC
<u>AURKA Antibody</u>	CSB-PA002454HA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>AURKB Antibody</u>	CSB-PA002458EA01HU	Human	ELISA, IHC, IF
<u>AURKC Antibody</u>	CSB-PA002460LA01HU	Human	ELISA, WB, IHC, IF
<u>AVIL Antibody</u>	CSB-PA002464LA01HU	Human	ELISA, WB
<u>AVL9 Antibody</u>	CSB-PA843316LA01HU	Human	ELISA, IHC, IF
<u>AVP Antibody</u>	CSB-PA11069A0Rb	Human, Rat	ELISA, WB, IHC
<u>AVPI1 Antibody</u>	CSB-PA732925LA01HU	Human	ELISA, IHC
<u>Avpr1a Antibody</u>	CSB-PA15624A0Rb	Rat, Mouse	ELISA, WB, IHC
<u>AVPR1A Antibody</u>	CSB-PA002468LA01HU	Human	ELISA, IHC, IF
<u>Avpr1b Antibody</u>	CSB-PA15639A0Rb	Rat	ELISA, WB
<u>AVPR2 Antibody</u>	CSB-PA002470LA01PI	Pig	ELISA
<u>AVPR2 Antibody</u>	CSB-PA002470LA01HU	Human	ELISA, IHC, IF
<u>AXIN1 Antibody</u>	CSB-PA517675LA01HU	Human	ELISA, IHC
<u>AXIN2 Antibody</u>	CSB-PA897472LA01HU	Human	ELISA, IHC, IF
<u>AXIN2 Antibody</u>	CSB-PA897472OA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>AXL Antibody</u>	CSB-PA326981LA01HU	Human	ELISA, IHC, IF
<u>AZGP1 Antibody</u>	CSB-PA11039A0Rb	Human	ELISA, WB, IHC
<u>AZI2 Antibody</u>	CSB-PA863992LA01HU	Human	ELISA, IHC, IF
<u>AZIN1 Antibody</u>	CSB-PA002484LA01HU	Human, Mouse	ELISA, WB, IHC
<u>B1 Antibody</u>	CSB-PA882501LA01HKA	Human herpesvirus 6B	ELISA, WB
<u>B2M Antibody</u>	CSB-PA002486YA01HU	Human, Mouse	ELISA, WB, IHC
<u>B2M Antibody</u>	CSB-PA05786A0Rb	Human, Mouse	ELISA, WB
<u>B2m Antibody</u>	CSB-PA002486LA01MO	Mouse	ELISA, WB, IHC
<u>B3GALT2 Antibody</u>	CSB-PA002491LA01HU	Human	ELISA, IHC, IF
<u>B3GALT4 Antibody</u>	CSB-PA002492LA01HU	Human, Mouse	ELISA, WB, IHC
<u>B3GNT2 Antibody</u>	CSB-PA002500LA01HU	Human	ELISA, WB, IHC
<u>B3GNT5 Antibody</u>	CSB-PA002503EA01HU	Human	ELISA, IHC
<u>B3GNT6 Antibody</u>	CSB-PA744392LA01HU	Human	ELISA, WB, IHC
<u>B4GALNT2 Antibody</u>	CSB-PA854161LA01HU	Human	ELISA, IHC
<u>B4GALT2 Antibody</u>	CSB-PA002514LA01HU	Human	ELISA, WB
<u>B4GALT3 Antibody</u>	CSB-PA00704A0Rb	Human	ELISA, WB, IHC
<u>B4GALT5 Antibody</u>	CSB-PA527356LA01HU	Human, Mouse	ELISA, WB, IHC
<u>B4GALT6 Antibody</u>	CSB-PA871391LA01HU	Human	ELISA, WB
<u>B4GALT7 Antibody</u>	CSB-PA890666LA01HU	Human, Mouse	ELISA, WB
<u>B4GAT1 Antibody</u>	CSB-PA002499EA01HU	Human	ELISA, IHC
<u>BAAT Antibody</u>	CSB-PA613484LA01HU	Human	ELISA, IHC, IF
<u>BABAM2 Antibody</u>	CSB-PA873684LA01HU	Human	ELISA, IHC, IF
<u>BACE1 Antibody</u>	CSB-PA002524LA01HU	Human	ELISA, IHC, IF
<u>BACE2 Antibody</u>	CSB-PA896746LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>BACH2 Antibody</u>	CSB-PA883639LA01HU	Human	ELISA, IHC, IF
<u>Bag3 Antibody</u>	CSB-PA864372HA01MO	Mouse	ELISA, IP
<u>BAG3 Antibody</u>	CSB-PA002531LA01HU	Human	ELISA, IHC, IF
<u>BAG4 Antibody</u>	CSB-PA002532LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>BAG5 Antibody</u>	CSB-PA890743HA01HU	Human, Mouse	ELISA, WB, IHC
<u>BAGE3 Antibody</u>	CSB-PA768236LA01HU	Human	ELISA, IHC
<u>BAIAP2 Antibody</u>	CSB-PA891988LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>BAIAP2L1 Antibody</u>	CSB-PA887045LA01HU	Human	ELISA, IHC
<u>BAIAP3 Antibody</u>	CSB-PA002547LA01HU	Human	ELISA, IHC
<u>BALF5 Antibody</u>	CSB-PA355978LA01EFA	Epstein-Barr virus	ELISA
<u>BANF1 Antibody</u>	CSB-PA002550LA01HU	Human	ELISA, IHC, IF
<u>bar Antibody</u>	CSB-PA321441HA01SOD	Streptomyces hygroscopicus	ELISA
<u>BARD1 Antibody</u>	CSB-PA860351LA01HU	Human	ELISA, WB
<u>BARHL1 Antibody</u>	CSB-PA880147LA01HU	Human	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>BARHL2 Antibody</u>	CSB-PA878919LA01HU	Human	ELISA, IHC
<u>BARX1 Antibody</u>	CSB-PA888023LA01HU	Human	ELISA, WB, IHC
<u>BAS1 Antibody</u>	CSB-PA822134HA01DOA	Arabidopsis thaliana	ELISA, WB
<u>BATF Antibody</u>	CSB-PA619074LA01HU	Human, Mouse	ELISA, WB, IHC
<u>BATF2 Antibody</u>	CSB-PA822685LA01HU	Human	ELISA, IHC
<u>Batf3 Antibody</u>	CSB-PA002572LA01RA	Rat	ELISA
<u>BATF3 Antibody</u>	CSB-PA882082LA01HU	Human	ELISA, IF
<u>BAZ2A Antibody</u>	CSB-PA883428HA01HU	Human	ELISA, IF
<u>BAZ2B Antibody</u>	CSB-PA883427LA01HU	Human	ELISA, IHC, IF
<u>BBC3 Antibody</u>	CSB-PA861160NA01HU	Human, Rat	ELISA, WB
<u>BBE1 Antibody</u>	CSB-PA341191LA01EOT	Eschscholzia californica	ELISA
<u>BBOX1 Antibody</u>	CSB-PA530450LA01HU	Human	ELISA, IHC, IF
<u>BBS1 Antibody</u>	CSB-PA836736LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>BBS10 Antibody</u>	CSB-PA823443LA01HU	Human	ELISA, IHC
<u>BBS5 Antibody</u>	CSB-PA847625LA01HU	Human	ELISA, IHC
<u>BBS7 Antibody</u>	CSB-PA814219LA01HU	Human	ELISA, IF
<u>BBS9 Antibody</u>	CSB-PA659641LA01HU	Human	ELISA, IHC, IF
<u>BBX Antibody</u>	CSB-PA845172LA01HU	Human	ELISA, IHC, IF
<u>BCAM Antibody</u>	CSB-PA341920LA01HU	Human	ELISA, IHC, IF
<u>BCAP31 Antibody</u>	CSB-PA03425A0Rb	Human, Mouse, Rat	ELISA, WB, IHC, IF, IP
<u>BCAR1 Antibody</u>	CSB-PA002593LA01HU	Human	ELISA, WB, IHC, IF
<u>BCAT1 Antibody</u>	CSB-PA002600EA01HU	Human	ELISA, IHC, IF
<u>BCAT2 Antibody</u>	CSB-PA002601LA01HU	Human, Mouse	ELISA, WB, IHC
<u>BCDIN3D Antibody</u>	CSB-PA767206LA01HU	Human	ELISA
<u>BCHE Antibody</u>	CSB-PA002604LA01DO	Dog	ELISA, IHC
<u>BCHE Antibody</u>	CSB-PA002604LA01HU	Human	ELISA, WB, IF
<u>BCKDHA Antibody</u>	CSB-PA002605LA01HU	Human	ELISA, WB, IHC
<u>BCL10 Antibody</u>	CSB-PA002608HA01HU	Human	ELISA, IHC, IF
<u>BCL11B Antibody</u>	CSB-PA880159HA01HU	Human	ELISA, IHC
<u>Bcl2 Antibody</u>	CSB-PA002611LA01MO	Mouse	ELISA, IHC
<u>BCL2 Antibody</u>	CSB-PA12959A0Rb	Human	ELISA, IHC, IF
<u>BCL2A1 Antibody</u>	CSB-PA614974LA01HU	Human	ELISA, IHC, IF
<u>BCL2L1 Antibody</u>	CSB-PA002613HA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>BCL2L10 Antibody</u>	CSB-PA884627LA01HU	Human	ELISA, IHC, IF
<u>BCL2L11 Antibody</u>	CSB-PA002615LA01HU	Human	ELISA, IHC, IF
<u>BCL2L2 Antibody</u>	CSB-PA821707HA01HU	Human	ELISA, IHC
<u>BCL3 Antibody</u>	CSB-PA002623LA01HU	Human, Mouse	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>BCL9 Antibody</u>	CSB-PA002631LA01HU	Human	ELISA, WB, IHC, IF
<u>BCLAF1 Antibody</u>	CSB-PA889149LA01HU	Human	ELISA, IF
<u>BCO1 Antibody</u>	CSB-PA864015LA01HU	Human	ELISA, IHC
<u>BCOR Antibody</u>	CSB-PA761519LA01HU	Human	ELISA
<u>BDH1 Antibody</u>	CSB-PA002648LA01HU	Human, Mouse	ELISA, WB, IHC
<u>BDH2 Antibody</u>	CSB-PA874814LA01HU	Human	ELISA, IHC
<u>BDKRB2 Antibody</u>	CSB-PA002652LA01HU	Human	ELISA, WB, IHC
<u>BDNF Antibody</u>	CSB-PA002655LA01HU	Human	ELISA, IHC
<u>BDNF Antibody</u>	CSB-PA05775A0Rb	Human, Mouse	ELISA, WB, IHC
<u>Bdnf Antibody</u>	CSB-PA16639A0Rb	Rat, Human	ELISA, WB
<u>BECN1 Antibody</u>	CSB-PA617917LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>benC Antibody</u>	CSB-PA27287A0Rb	Streptomyces sp	ELISA, WB
<u>BEND7 Antibody</u>	CSB-PA850823LA01HU	Human	ELISA, IHC
<u>bepA Antibody</u>	CSB-PA355426LA01ENV	Escherichia coli	ELISA, WB
<u>BET1 Antibody</u>	CSB-PA002670LA01HU	Human, Rat	ELISA, WB
<u>Beta-Actin Monoclonal Antibody</u>	CSB-MA000091M0m		ELISA
<u>Beta-lactamase inhibitory protein Antibody</u>	CSB-PA327444LA01FOA	Streptomyces clavuligerus	ELISA, WB
<u>Beta-lytic metalloendopeptidase Antibody</u>	CSB-PA355563LA01LLB	Lysobacter enzymogenes	ELISA, WB
<u>Beta-mammal toxin C_{ss}4 Antibody</u>	CSB-PA350214LA01DRC	Centruroides suffusus	ELISA
<u>BETV4 Antibody</u>	CSB-PA19907A0Rb	Betula pendula	ELISA, WB
<u>BETVIA Antibody</u>	CSB-PA322213HA01BSS	Betula pendula	ELISA, WB
<u>BETVII Antibody</u>	CSB-PA335151HA01BSS	Betula pendula	ELISA
<u>BEX1 Antibody</u>	CSB-PA862062LA01HU	Human	ELISA, IHC
<u>BEX2 Antibody</u>	CSB-PA883628LA01HU	Human	ELISA, IHC, IF
<u>BEX5 Antibody</u>	CSB-PA705160LA01HU	Human	ELISA, IHC
<u>BFSP2 Antibody</u>	CSB-PA615697LA01HU	Human	ELISA, WB, IHC, IF
<u>BG2 Antibody</u>	CSB-PA327203LA01DOA	Arabidopsis thaliana	ELISA, WB
<u>Bglap Antibody</u>	CSB-PA002682DA01RA	Rat	ELISA
<u>BGLAP Antibody</u>	CSB-PA11834A0Rb	Human	ELISA, IHC
<u>BGN Antibody</u>	CSB-PA002683EA01HU	Human	ELISA, IHC, IF
<u>BHLHA15 Antibody</u>	CSB-PA755911LA01HU	Human	ELISA, IHC
<u>BHLHE40 Antibody</u>	CSB-PA002691LA01HU	Human	ELISA, IHC, IF
<u>BHLHE41 Antibody</u>	CSB-PA883658LA01HU	Human	ELISA, WB
<u>BHMT2 Antibody</u>	CSB-PA11829A0Rb	Human, Mouse	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>BHRF1 Antibody</u>	CSB-PA314488LA01EFC	Epstein-Barr virus	ELISA, WB
<u>BICC1 Antibody</u>	CSB-PA872503LA01HU	Human	ELISA, IHC, IF
<u>BICD2 Antibody</u>	CSB-PA844723LA01HU	Human	ELISA, WB, IHC, IF
<u>BID Antibody</u>	CSB-PA002698HA01HU	Human	ELISA, IHC
<u>BIN2 Antibody</u>	CSB-PA002701LA01HU	Human	ELISA, IHC, IF
<u>BIN3 Antibody</u>	CSB-PA002702LA01HU	Human	ELISA, IF
<u>BIRC2 Antibody</u>	CSB-PA618777EA01HU	Human, Mouse, Rat	ELISA, WB
<u>BIRC3 Antibody</u>	CSB-PA623808LA01HU	Human, Mouse	ELISA, WB, IHC
<u>Birc5 Antibody</u>	CSB-PA002706HA01MO	Mouse, Human	ELISA, WB
<u>BIRC5 Antibody</u>	CSB-PA07249A0Rb	Human	ELISA, IF
<u>bla Antibody</u>	CSB-PA352353YA01ENL	Escherichia coli	ELISA
<u>bla Antibody</u>	CSB-PA321251LA01EZW	Pseudomonas aeruginosa	ELISA
<u>bla Antibody</u>	CSB-PA333368LA01ENL	Escherichia coli	ELISA, WB
<u>blaNDM-1 Antibody</u>	CSB-PA07617A0Rb	Klebsiella pneumoniae	ELISA
<u>blaNDM-1 Antibody</u>	CSB-PA510513HA01KBG	Klebsiella pneumoniae	ELISA
<u>BLLF1 Antibody</u>	CSB-PA360986LA01EFA	Epstein-Barr virus	ELISA
<u>BLLF3 Antibody</u>	CSB-PA319032LA01EFC	Epstein-Barr virus	ELISA, WB
<u>BLNK Antibody</u>	CSB-PA840984LA01HU	Human	ELISA, WB, IHC, IF
<u>BLOC1S1 Antibody</u>	CSB-PA00495A0Rb	Human	ELISA, WB, IHC, IF
<u>BLOC1S5 Antibody</u>	CSB-PA837439LA01HU	Human	ELISA, IHC
<u>BLVRB Antibody</u>	CSB-PA002722HA01HU	Human	ELISA, WB, IHC, IF
<u>BLZF1 Antibody</u>	CSB-PA872478LA01HU	Human, Mouse	ELISA, WB, IHC
<u>BMI1 Antibody</u>	CSB-PA03345A0Rb	Human	ELISA, WB, IHC
<u>BMP10 Antibody</u>	CSB-PA002734LA01HU	Human	ELISA, IHC
<u>BMP2 Antibody</u>	CSB-PA09419A0Rb	Human	ELISA, WB, IHC, IF
<u>bmp2b Antibody</u>	CSB-PA20849A0Rb	Zebrafish	ELISA, WB
<u>BMP3 Antibody</u>	CSB-PA002739LA01HU	Human	ELISA, WB, IHC
<u>BMP4 Antibody</u>	CSB-PA05799A0Rb	Human, Mouse, Zebrafish	ELISA, WB, IHC
<u>BMP6 Antibody</u>	CSB-PA10709A0Rb	Human	ELISA, IHC
<u>BMP7 Antibody</u>	CSB-PA05805A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>BMPER Antibody</u>	CSB-PA839829LA01HU	Human	ELISA, IHC
<u>BMPR1A Antibody</u>	CSB-PA002748LA01HU	Human	ELISA, WB, IHC, IF
<u>BMPR1B Antibody</u>	CSB-PA002751LA01HU	Human, Rat	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>BOK Antibody</u>	CSB-PA892353LA01HU	Human	ELISA, IHC
<u>BOLA1 Antibody</u>	CSB-PA002774LA01HU	Human	ELISA
<u>BOLA2 Antibody</u>	CSB-PA867143LA01HU	Human	ELISA, IHC, IF
<u>BOLA3 Antibody</u>	CSB-PA690486LA01HU	Human	ELISA, IHC
<u>BORA Antibody</u>	CSB-PA750873LA01HU	Human	ELISA, WB, IHC
<u>BORCS5 Antibody</u>	CSB-PA822164LA01HU	Human	ELISA, IHC
<u>BORCS7 Antibody</u>	CSB-PA842639LA01HU	Human	ELISA, IHC, IF
<u>botA Antibody</u>	CSB-PA320799LA01CLQ	Clostridium botulinum	ELISA
<u>botB Antibody</u>	CSB-PA318083LA01CLQ	Clostridium botulinum	ELISA
<u>botD Antibody</u>	CSB-PA323098LA01CLQ	Clostridium botulinum	ELISA
<u>botF Antibody</u>	CSB-PA341192HA01CLQ	Clostridium botulinum	ELISA, WB
<u>Botulinum neurotoxin type C1 Antibody</u>	CSB-PA325596LA01CLQ	Clostridium botulinum	ELISA
<u>Botulinum neurotoxin type E Antibody</u>	CSB-PA327013LA01DUK	Clostridium butyricum	ELISA
<u>BPESC1 Antibody</u>	CSB-PA884420LA01HU	Human	ELISA, IHC
<u>BPI Antibody</u>	CSB-PA002783LA01HU	Human	ELISA, WB, IF
<u>BPIFB3 Antibody</u>	CSB-PA003693LA01HU	Human, Mouse	ELISA, WB
<u>BRAF Antibody</u>	CSB-PA002791LA01HU	Human	ELISA, IHC, IF
<u>BRAP Antibody</u>	CSB-PA773793LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>BRAT1 Antibody</u>	CSB-PA757659LA01HU	Human	ELISA, WB, IHC
<u>BRCA1 Antibody</u>	CSB-PA09659A0Rb	Human	ELISA
<u>brca2 Antibody</u>	CSB-PA24309A0Rb	Zebrafish	ELISA, WB
<u>BRCC3 Antibody</u>	CSB-PA002798LA01HU	Human	ELISA, WB, IHC
<u>BRD1 Antibody</u>	CSB-PA002799LA01HU	Human	ELISA, IHC
<u>BRD3 Antibody</u>	CSB-PA624013HA01HU	Human	ELISA, IHC, IF
<u>BRD4 Antibody</u>	CSB-PA002802LA01HU	Human	ELISA, IF
<u>BRD7 Antibody</u>	CSB-PA865097LA01HU	Human	ELISA, IHC, IF
<u>BRD8 Antibody</u>	CSB-PA880933LA01HU	Human	ELISA, IHC
<u>BRINP1 Antibody</u>	CSB-PA006515LA01HU	Human	ELISA, IHC, IF
<u>BRINP2 Antibody</u>	CSB-PA866338LA01HU	Human	ELISA, IHC
<u>BRINP3 Antibody</u>	CSB-PA755233LA01HU	Human	ELISA, WB, IHC
<u>BRIP1 Antibody</u>	CSB-PA866300LA01HU	Human	ELISA, IHC, IF
<u>BROX Antibody</u>	CSB-PA730136LA01HU	Human	ELISA, WB, IHC, IF
<u>BRSK1 Antibody</u>	CSB-PA837438LA01HU	Human, Mouse	ELISA, WB, IHC
<u>BSG Antibody</u>	CSB-PA11759A0Rb	Human	ELISA, WB, IHC
<u>BSX Antibody</u>	CSB-PA664240LA01HU	Human	ELISA, IHC
<u>BTAF1 Antibody</u>	CSB-PA002839LA01HU	Human	ELISA, IF
<u>BTBD2 Antibody</u>	CSB-PA05095A0Rb	Human	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>BTD Antibody</u>	CSB-PA002854LA01HU	Human	ELISA, IHC
<u>BTF3 Antibody</u>	CSB-PA01714A0Rb	Human	ELISA, IHC, IF
<u>btfP Antibody</u>	CSB-PA346537LA01BDP	Bacteroides fragilis	ELISA, WB
<u>BTG4 Antibody</u>	CSB-PA878918LA01HU	Human, Mouse	ELISA, WB, IHC
<u>BTK Antibody</u>	CSB-PA002867LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>BTN2A1 Antibody</u>	CSB-PA755482LA01HU	Human	ELISA, IHC
<u>BTN2A2 Antibody</u>	CSB-PA848825LA01HU	Human, Mouse	ELISA, WB, IHC
<u>BTN3A1 Antibody</u>	CSB-PA002873LA01HU	Human, Mouse	ELISA, WB, IHC
<u>BTN3A2 Antibody</u>	CSB-PA002874LA01HU	Human	ELISA
<u>BTN3A3 Antibody</u>	CSB-PA002875LA01HU	Human	ELISA
<u>BTNL2 Antibody</u>	CSB-PA883443LA01HU	Human, Mouse	ELISA, WB, IHC
<u>BTNL3 Antibody</u>	CSB-PA764943LA01HU	Human	ELISA, IHC
<u>BTNL9 Antibody</u>	CSB-PA744270LA01HU	Human	ELISA
<u>BTRC Antibody</u>	CSB-PA897464LA01HU	Human	ELISA, WB, IHC, IF
<u>BUB1 Antibody</u>	CSB-PA002882LA01HU	Human	ELISA, IHC
<u>BUB1B Antibody</u>	CSB-PA002883LA01HU	Human	ELISA, IHC, IF, IP
<u>BUB3 Antibody</u>	CSB-PA002884LA01HU	Human, Mouse	ELISA, WB, IHC
<u>BUD13 Antibody</u>	CSB-PA863104LA01HU	Human	ELISA, IHC, IF
<u>BUD31 Antibody</u>	CSB-PA01749A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>BZLF2 Antibody</u>	CSB-PA365876LA01EFA	Epstein-Barr virus	ELISA, WB
<u>C Antibody</u>	CSB-PA323345LA01RAI	Rabies virus	ELISA, WB
<u>C10orf111 Antibody</u>	CSB-PA843145LA01HU	Human	ELISA, IHC, IF
<u>C10orf53 Antibody</u>	CSB-PA818719LA01HU	Human	ELISA, IHC
<u>C10orf62 Antibody</u>	CSB-PA719409LA01HU	Human	ELISA, IHC
<u>C10orf82 Antibody</u>	CSB-PA837861LA01HU	Human, Mouse	ELISA, WB, IHC
<u>C10orf91 Antibody</u>	CSB-PA722562LA01HU	Human	ELISA, IHC
<u>C11orf1 Antibody</u>	CSB-PA880979LA01HU	Human	ELISA, IHC
<u>C11orf45 Antibody</u>	CSB-PA819448LA01HU	Human	ELISA, IHC
<u>C11orf53 Antibody</u>	CSB-PA809008LA01HU	Human, Mouse	ELISA, WB, IHC
<u>C11orf68 Antibody</u>	CSB-PA862008LA01HU	Human, Mouse	ELISA, WB, IHC
<u>C11orf84 Antibody</u>	CSB-PA887125LA01HU	Human	ELISA
<u>C12orf43 Antibody</u>	CSB-PA850269LA01HU	Human	ELISA, IHC
<u>C12orf49 Antibody</u>	CSB-PA867163LA01HU	Human	ELISA, IHC
<u>C12orf50 Antibody</u>	CSB-PA822743LA01HU	Human	ELISA
<u>C12orf57 Antibody</u>	CSB-PA859105LA01HU	Human	ELISA
<u>C12orf60 Antibody</u>	CSB-PA719466LA01HU	Human	ELISA, IHC, IF
<u>C12orf65 Antibody</u>	CSB-PA880965LA01HU	Human	ELISA, IHC
<u>C14orf177 Antibody</u>	CSB-PA703583LA01HU	Human	ELISA, IHC
<u>C15orf32 Antibody</u>	CSB-PA654885LA01HU	Human	ELISA



Product Name	Code	Species Reactivity	Applications
<u>C15orf40 Antibody</u>	CSB-PA837852LA01HU	Human	ELISA, IHC
<u>C15orf41 Antibody</u>	CSB-PA897474LA01HU	Human	ELISA, WB
<u>C15orf53 Antibody</u>	CSB-PA850848LA01HU	Human	ELISA
<u>C16orf62 Antibody</u>	CSB-PA773584HA01HU	Human	ELISA, WB, IF
<u>C16orf78 Antibody</u>	CSB-PA848813LA01HU	Human	ELISA, IHC
<u>C16orf87 Antibody</u>	CSB-PA764800LA01HU	Human	ELISA, WB, IP
<u>C17orf47 Antibody</u>	CSB-PA850873LA01HU	Human	ELISA, IHC
<u>C17orf58 Antibody</u>	CSB-PA643573LA01HU	Human	ELISA, IHC, IF
<u>C17orf64 Antibody</u>	CSB-PA771458LA01HU	Human, Mouse	ELISA, WB, IHC
<u>C17orf75 Antibody</u>	CSB-PA881009LA01HU	Human	ELISA
<u>C17orf78 Antibody</u>	CSB-PA847630LA01HU	Human	ELISA, IHC, IF
<u>C19orf12 Antibody</u>	CSB-PA003434LA01HU	Human	ELISA, IHC, IF
<u>C19orf44 Antibody</u>	CSB-PA884480LA01HU	Human	ELISA, IHC
<u>C19orf47 Antibody</u>	CSB-PA003463LA01HU	Human, Mouse	ELISA, WB, IHC
<u>C19orf53 Antibody</u>	CSB-PA891978LA01HU	Human	ELISA
<u>C1D Antibody</u>	CSB-PA03474A0Rb	Human	ELISA, WB, IHC, IF
<u>C1GALT1 Antibody</u>	CSB-PA882094LA01HU	Human	ELISA, WB, IHC, IF
<u>C1orf105 Antibody</u>	CSB-PA003497LA01HU	Human	ELISA, IHC
<u>C1orf112 Antibody</u>	CSB-PA878881LA01HU	Human	ELISA, IHC
<u>C1orf123 Antibody</u>	CSB-PA873673LA01HU	Human, Mouse	ELISA, WB, IHC
<u>C1orf127 Antibody</u>	CSB-PA854077LA01HU	Human	ELISA, IHC
<u>C1orf158 Antibody</u>	CSB-PA839778LA01HU	Human	ELISA
<u>C1orf189 Antibody</u>	CSB-PA730130LA01HU	Human	ELISA, IHC, IF
<u>C1orf210 Antibody</u>	CSB-PA808993LA01HU	Human	ELISA
<u>C1orf53 Antibody</u>	CSB-PA719622LA01HU	Human	ELISA, IHC
<u>C1orf54 Antibody</u>	CSB-PA840992LA01HU	Human	ELISA, IHC
<u>C1orf64 Antibody</u>	CSB-PA843340LA01HU	Human	ELISA
<u>C1orf87 Antibody</u>	CSB-PA847603LA01HU	Human	ELISA
<u>C1QA Antibody</u>	CSB-PA003637LA01HU	Human	ELISA, IHC, IF
<u>C1QC Antibody</u>	CSB-PA003641LA01HU	Human	ELISA, IHC
<u>C1QTNF1 Antibody</u>	CSB-PA871599EA01HU	Human	ELISA, IHC
<u>C1QTNF3 Antibody</u>	CSB-PA003648LA01HU	Human	ELISA, IHC
<u>C1QTNF6 Antibody</u>	CSB-PA874831EA01HU	Human, Mouse	ELISA, WB, IHC
<u>C1QTNF9B Antibody</u>	CSB-PA450719LA01HU	Human	ELISA, IHC, IF
<u>C1R Antibody</u>	CSB-PA003655LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>C1S Antibody</u>	CSB-PA003657LA01HU	Human	ELISA, IHC
<u>C20orf196 Antibody</u>	CSB-PA812897LA01HU	Human	ELISA
<u>C20orf197 Antibody</u>	CSB-PA822688DA01HU	Human	ELISA
<u>C20orf27 Antibody</u>	CSB-PA875631LA01HU	Human	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>C21orf62-AS1 Antibody</u>	CSB-PA613779LA01HU	Human	ELISA
<u>C22orf23 Antibody</u>	CSB-PA871624LA01HU	Human	ELISA
<u>C2CD3 Antibody</u>	CSB-PA673367LA01HU	Human	ELISA, IHC, IF
<u>C2CD6 Antibody</u>	CSB-PA681515LA01HU	Human	ELISA, IHC
<u>C2orf15 Antibody</u>	CSB-PA823878LA01HU	Human	ELISA, WB, IHC
<u>C2orf49 Antibody</u>	CSB-PA861147LA01HU	Human	ELISA, IHC
<u>C2orf50 Antibody</u>	CSB-PA839365LA01HU	Human	ELISA
<u>C2orf69 Antibody</u>	CSB-PA822731LA01HU	Human	ELISA, IHC
<u>C2orf76 Antibody</u>	CSB-PA659452LA01HU	Human	ELISA, IHC
<u>C2orf80 Antibody</u>	CSB-PA603916LA01HU	Human	ELISA
<u>C2orf83 Antibody</u>	CSB-PA709610LA01HU	Human	ELISA, IHC
<u>C3 Antibody</u>	CSB-PA10599A0Rb	Human	ELISA, IHC, IF
<u>C3 Antibody</u>	CSB-PA365507LA01MO	Mouse, Human	ELISA, WB
<u>C3 Antibody</u>	CSB-PA360475HA01RA	Rat	ELISA
<u>C3orf14 Antibody</u>	CSB-PA884512LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>C3orf18 Antibody</u>	CSB-PA892447LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>C3orf20 Antibody</u>	CSB-PA818767LA01HU	Human	ELISA, IF
<u>C3orf36 Antibody</u>	CSB-PA662245LA01HU	Human	ELISA, IHC
<u>C3orf49 Antibody</u>	CSB-PA856904LA01HU	Human	ELISA, WB
<u>C4A Antibody</u>	CSB-PA003949LA01HU	Human	ELISA, WB, IHC
<u>C4B Antibody</u>	CSB-PA313364LA01HU	Human	ELISA, IHC
<u>C4BPA Antibody</u>	CSB-PA003951LA01HU	Human	ELISA, IHC
<u>C4orf3 Antibody</u>	CSB-PA823892HA01HU	Human	ELISA, IHC, IF
<u>C4orf3 Antibody</u>	CSB-PA823892LA01HU	Human	ELISA, WB
<u>C4orf33 Antibody</u>	CSB-PA822682LA01HU	Human, Mouse	ELISA, WB
<u>C4orf46 Antibody</u>	CSB-PA687370LA01HU	Human	ELISA, IHC
<u>C5 Antibody</u>	CSB-PA003995LA01HU	Human	ELISA, IHC
<u>C5 Antibody</u>	CSB-PA003995YA01PI	Pig	ELISA, IHC
<u>C5AR1 Antibody</u>	CSB-PA003996LA01HU	Human	ELISA, WB, IF
<u>C5AR2 Antibody</u>	CSB-PA868390PA01HU	Human	ELISA, IHC
<u>C5orf22 Antibody</u>	CSB-PA675468LA01HU	Human, Mouse	ELISA, WB, IHC
<u>C5orf24 Antibody</u>	CSB-PA773601LA01HU	Human, Mouse	ELISA, WB, IHC
<u>C5orf51 Antibody</u>	CSB-PA004030LA01HU	Human	ELISA, IHC, IF
<u>C6 Antibody</u>	CSB-PA004036LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>C6orf136 Antibody</u>	CSB-PA713102LA01HU	Human	ELISA
<u>C6orf141 Antibody</u>	CSB-PA725558LA01HU	Human	ELISA, IHC
<u>C6orf203 Antibody</u>	CSB-PA889170LA01HU	Human	ELISA, IHC
<u>C6orf223 Antibody</u>	CSB-PA836654LA01HU	Human	ELISA, IHC, IF
<u>C6orf89 Antibody</u>	CSB-PA764939LA01HU	Human	ELISA, IHC
<u>C7orf31 Antibody</u>	CSB-PA839824LA01HU	Human	ELISA, WB, IF
<u>C7orf33 Antibody</u>	CSB-PA840977LA01HU	Human	ELISA



Product Name	Code	Species Reactivity	Applications
<u>C7orf50 Antibody</u>	CSB-PA866265LA01HU	Human	ELISA, WB
<u>C7orf61 Antibody</u>	CSB-PA814238LA01HU	Human	ELISA, IHC
<u>C8G Antibody</u>	CSB-PA004196LA01HU	Human	ELISA, WB
<u>C8orf31 Antibody</u>	CSB-PA843300LA01HU	Human	ELISA
<u>C8orf34 Antibody</u>	CSB-PA675467LA01HU	Human, Mouse	ELISA, WB, IHC
<u>C8orf48 Antibody</u>	CSB-PA853456LA01HU	Human, Mouse	ELISA, WB, IHC
<u>C8orf59 Antibody</u>	CSB-PA836635LA01HU	Human	ELISA
<u>C8orf74 Antibody</u>	CSB-PA764703LA01HU	Human	ELISA, IHC
<u>C9 Antibody</u>	CSB-PA004259LA01HU	Human	ELISA, IHC, IF
<u>C9orf16 Antibody</u>	CSB-PA05339A0Rb	Human	ELISA, IHC
<u>C9orf163 Antibody</u>	CSB-PA822739LA01HU	Human	ELISA, IHC
<u>C9orf43 Antibody</u>	CSB-PA819445LA01HU	Human	ELISA, IHC
<u>C9orf57 Antibody</u>	CSB-PA730146LA01HU	Human	ELISA, IHC
<u>C9orf64 Antibody</u>	CSB-PA713128LA01HU	Human	ELISA, WB, IHC
<u>C9orf72 Antibody</u>	CSB-PA842710LA01HU	Human	ELISA, IHC, IF
<u>CA1 Antibody</u>	CSB-PA004364HA01HU	Human, Mouse	ELISA, WB, IHC, IF, IP
<u>Ca1 Antibody</u>	CSB-PA004364HA01RA	Rat	ELISA, WB
<u>Ca1 Antibody</u>	CSB-PA004364YA01RA	Rat	ELISA, WB
<u>CA10 Antibody</u>	CSB-PA878879EA01HU	Human, Mouse	ELISA, WB
<u>CA11 Antibody</u>	CSB-PA004366LA01HU	Human	ELISA, IHC
<u>CA12 Antibody</u>	CSB-PA004367LA01HU	Human, Mouse	ELISA, WB
<u>CA12 Antibody</u>	CSB-PA004367YA01HU	Human, Mouse	ELISA, WB, IHC
<u>CA13 Antibody</u>	CSB-PA847613LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CA2 Antibody</u>	CSB-PA004370HA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>CA3 Antibody</u>	CSB-PA004371LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CA4 Antibody</u>	CSB-PA004372LA01HU	Human	ELISA, IHC, IF
<u>CA5B Antibody</u>	CSB-PA896871EA01HU	Human, Mouse	ELISA, WB, IHC
<u>CA6 Antibody</u>	CSB-PA004377LA01HU	Human	ELISA, IHC
<u>CAB39 Antibody</u>	CSB-PA640939LA01HU	Human	ELISA, IF
<u>CAB39L Antibody</u>	CSB-PA872527LA01HU	Human	ELISA, WB, IF
<u>CABIN1 Antibody</u>	CSB-PA23399A0Rb	Human, Rat	ELISA, WB, IHC
<u>CABP2 Antibody</u>	CSB-PA873608LA01HU	Human	ELISA, IF
<u>CABP4 Antibody</u>	CSB-PA004390LA01HU	Human	ELISA, IHC, IF
<u>CABP5 Antibody</u>	CSB-PA865092LA01HU	Human, Mouse	ELISA, WB
<u>CACNA1A Antibody</u>	CSB-PA004397LA01HU	Human	ELISA, IHC, IF
<u>CACNA1B Antibody</u>	CSB-PA004398LA01HU	Human	ELISA, IHC, IF
<u>CACNA1C Antibody</u>	CSB-PA615724LA01HU	Human	ELISA, IHC, IF
<u>CACNA1D Antibody</u>	CSB-PA23409A0Rb	Human	ELISA, IF
<u>CACNA1H Antibody</u>	CSB-PA004404EA01HU	Human	ELISA, IHC, IF
<u>CACNA2D1 Antibody</u>	CSB-PA004407LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>CACNA2D4 Antibody</u>	CSB-PA801802LA01HU	Human	ELISA, IHC, IF
<u>CACNB2 Antibody</u>	CSB-PA804482LA01HU	Human	ELISA, IHC, IF
<u>CACNB4 Antibody</u>	CSB-PA004414LA01HU	Human	ELISA, IHC, IF
<u>CACNG3 Antibody</u>	CSB-PA004417LA01HU	Human	ELISA, WB, IHC, IF
<u>CACNG6 Antibody</u>	CSB-PA883625OA01HU	Human	ELISA, IF
<u>CACNG8 Antibody</u>	CSB-PA823907OA01HU	Human	ELISA, IHC
<u>CACYBP Antibody</u>	CSB-PA872536LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>CADM3 Antibody</u>	CSB-PA004427LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>cagA Antibody</u>	CSB-PA14367A0Rb	Helicobacter pylori	ELISA
<u>cagA Antibody</u>	CSB-PA301307LA01HUS	Helicobacter pylori	ELISA
<u>cagA Antibody</u>	CSB-PA345531YA01HUV	Helicobacter pylori	ELISA
<u>CALB1 Antibody</u>	CSB-PA004432LA01HU	Human	ELISA, IHC, IF
<u>CALCA Antibody</u>	CSB-PA12134A0Rb	Human	ELISA
<u>CALCA Antibody</u>	CSB-PA004434HA01DO	Dog	ELISA, WB
<u>CALHM3 Antibody</u>	CSB-PA768231LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CALM1 Antibody</u>	CSB-PA004445HA01HU	Human	ELISA, IF
<u>CALM1 Antibody</u>	CSB-PA08284A0Rb	Human	ELISA, IHC
<u>CALM1 Antibody</u>	CSB-PA08285A0Rb	Human, Mouse	ELISA, WB
<u>CALM1 Antibody</u>	CSB-PA08287A0Rb	Human	ELISA
<u>CALML4 Antibody</u>	CSB-PA836218EA01HU	Human	ELISA, IHC
<u>CALML5 Antibody</u>	CSB-PA873713HA01HU	Human	ELISA, WB, IHC
<u>CALN1 Antibody</u>	CSB-PA871605LA01HU	Human, Mouse	ELISA, WB
<u>CALR Antibody</u>	CSB-PA07864A0Rb	Human	ELISA, WB, IHC
<u>CALR Antibody</u>	CSB-PA07867A0Rb	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>CALR Antibody</u>	CSB-PA004458YA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CALR Antibody</u>	CSB-PA004458YA01PI	Pig	ELISA
<u>CAMK1 Antibody</u>	CSB-PA614509LA01HU	Human, Mouse, Rat	ELISA, WB, IF
<u>Camk2n1 Antibody</u>	CSB-PA004469HA01MO	Mouse	ELISA
<u>CAMK2N1 Antibody</u>	CSB-PA004469NA01HU	Human	ELISA, IF
<u>CAMK2N2 Antibody</u>	CSB-PA004470NA01HU	Human	ELISA, IF
<u>CAMKK2 Antibody</u>	CSB-PA822299LA01HU	Human	ELISA, IHC
<u>CAMLG Antibody</u>	CSB-PA004475LA01HU	Human	ELISA, WB, IF, IP
<u>CAMTA1 Antibody</u>	CSB-PA897322LA01HU	Human	ELISA, IF
<u>CAND1 Antibody</u>	CSB-PA774815LA01HU	Human	ELISA, WB, IHC, IF
<u>CANT1 Antibody</u>	CSB-PA13099A0Rb	Human	ELISA



Product Name	Code	Species Reactivity	Applications
<u>CANX Antibody</u>	CSB-PA00899A0Rb	Human	ELISA, WB, IHC, IF
<u>CAP1 Antibody</u>	CSB-PA004486LA01HU	Human	ELISA, WB, IHC, IF
<u>CAP2 Antibody</u>	CSB-PA004487LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CAPN1 Antibody</u>	CSB-PA004490LA01HU	Human	ELISA, WB, IHC, IF
<u>CAPN1 Antibody</u>	CSB-PA004490EA01HU	Human	ELISA, WB, IHC
<u>CAPN10 Antibody</u>	CSB-PA864027LA01HU	Human	ELISA, IF
<u>CAPN14 Antibody</u>	CSB-PA004495HA01HU	Human	ELISA, IHC, IF
<u>CAPN2 Antibody</u>	CSB-PA004496LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>CAPN2 Antibody</u>	CSB-PA004496LA01PI	Pig	ELISA
<u>CAPN6 Antibody</u>	CSB-PA897594LA01HU	Human	ELISA, WB, IHC
<u>CAPNS1 Antibody</u>	CSB-PA03179A0Rb	Human	ELISA, IF
<u>CAPNS2 Antibody</u>	CSB-PA822248LA01HU	Human, Mouse	ELISA, WB, IF
<u>Caprin1 Antibody</u>	CSB-PA711733HA01RA	Rat	ELISA
<u>CAPRIN1 Antibody</u>	CSB-PA614529LA01HU	Human	ELISA, IHC, IF
<u>CAPRIN2 Antibody</u>	CSB-PA740238OA01HU	Human	ELISA, IHC, IF
<u>CAPS Antibody</u>	CSB-PA622689LA01HU	Human	ELISA, IHC
<u>CAPS2 Antibody</u>	CSB-PA861167LA01HU	Human	ELISA, WB, IHC
<u>Capsid protein VP1 Antibody</u>	CSB-PA362102LA01HPM	Human parvovirus B19	ELISA
<u>CAPZA1 Antibody</u>	CSB-PA004510LA01HU	Human	ELISA, WB, IHC, IF
<u>CAPZA2 Antibody</u>	CSB-PA004512LA01HU	Human	ELISA, IHC
<u>CARD14 Antibody</u>	CSB-PA883622LA01HU	Human	ELISA, IHC, IF
<u>CARD16 Antibody</u>	CSB-PA682727LA01HU	Human	ELISA, IF
<u>CARD18 Antibody</u>	CSB-PA004521LA01HU	Human	ELISA, IHC
<u>CARD19 Antibody</u>	CSB-PA842711LA01HU	Human	ELISA, IHC
<u>CARD6 Antibody</u>	CSB-PA863138LA01HU	Human	ELISA, IHC, IF
<u>CARHSP1 Antibody</u>	CSB-PA896497LA01HU	Human	ELISA
<u>CARM1 Antibody</u>	CSB-PA771462LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CARMIL1 Antibody</u>	CSB-PA713356LA01HU	Human	ELISA, IF
<u>CARMIL2 Antibody</u>	CSB-PA739364LA01HU	Human	ELISA, WB, IP
<u>CARS Antibody</u>	CSB-PA004529LA01HU	Human, Mouse	ELISA, WB, IHC
<u>cas3 Antibody</u>	CSB-PA24029A0Rb	Thermus thermophilus	ELISA, WB
<u>CASC1 Antibody</u>	CSB-PA754371LA01HU	Human, Rat	ELISA, WB, IHC
<u>Casein Antibody</u>	CSB-PA00880E0Rb	Human	ELISA
<u>CASK Antibody</u>	CSB-PA004539LA01HU	Human	ELISA, IHC
<u>CASP1 Antibody</u>	CSB-PA05824A0Rb	Human	ELISA, IHC
<u>CASP1 Antibody</u>	CSB-PA05829A0Rb	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>CASP10 Antibody</u>	CSB-PA838815LA01HU	Human	ELISA, IF
<u>Casp12 Antibody</u>	CSB-PA004545LA01MO	Mouse, Human	ELISA, WB
<u>CASP12 Antibody</u>	CSB-PA747641LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CASP2 Antibody</u>	CSB-PA004547LA01HU	Human	ELISA, IHC, IF
<u>CASP3 Antibody</u>	CSB-PA05689A0Rb	Human	ELISA, IHC, IF
<u>CASP4 Antibody</u>	CSB-PA004549LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CASP5 Antibody</u>	CSB-PA004550LA01HU	Human	ELISA
<u>CASP6 Antibody</u>	CSB-PA004551EA01HU	Human, Mouse	ELISA, WB, IHC
<u>CASP7 Antibody</u>	CSB-PA15479A0Rb	Human	ELISA, WB, IHC, IF, IP
<u>CASP8 Antibody</u>	CSB-PA18109A0Rb	Human	ELISA, WB, IHC, IF
<u>CASP9 Antibody</u>	CSB-PA004555LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CASQ1 Antibody</u>	CSB-PA004556LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CASQ2 Antibody</u>	CSB-PA004557LA01HU	Human	ELISA, WB, IHC
<u>CASR Antibody</u>	CSB-PA004558EA01HU	Human	ELISA, IF
<u>CAST Antibody</u>	CSB-PA004561YA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CASTOR2 Antibody</u>	CSB-PA28239A0Rb	Human	ELISA, WB
<u>CASTOR3 Antibody</u>	CSB-PA822746LA01HU	Human	ELISA, IHC
<u>CASZ1 Antibody</u>	CSB-PA773036LA01HU	Human	ELISA, WB, IHC
<u>CAT1 Antibody</u>	CSB-PA822149LA01DOA	Arabidopsis thaliana	ELISA, WB
<u>Cationic trypsin Antibody</u>	CSB-PA361971LA01DO	Dog	ELISA
<u>CATSPER2 Antibody</u>	CSB-PA842731LA01HU	Human	ELISA, WB, IHC, IF
<u>CATSPER4 Antibody</u>	CSB-PA745763LA01HU	Human	ELISA, IHC, IF
<u>CATSPERE Antibody</u>	CSB-PA719393LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CAV1 Antibody</u>	CSB-PA004571OA01HU	Human	ELISA, IF
<u>CAV2 Antibody</u>	CSB-PA004572LA01HU	Human	ELISA, WB, IHC, IF
<u>CAVIN1 Antibody</u>	CSB-PA019070LA01HU	Human	ELISA, IF
<u>CBFA2T3 Antibody</u>	CSB-PA004576LA01HU	Human	ELISA, WB, IHC
<u>CBFB Antibody</u>	CSB-PA623835LA01HU	Human	ELISA
<u>CBL Antibody</u>	CSB-PA004578LA01HU	Human	ELISA, WB, IHC, IF
<u>CBLIF Antibody</u>	CSB-PA009417LA01HU	Human	ELISA, IHC
<u>CBLIF Antibody</u>	CSB-PA009417YA01HU	Human, Mouse	ELISA, WB, IHC
<u>CBLL2 Antibody</u>	CSB-PA836684LA01HU	Human, Mouse	ELISA, WB, IF



Product Name	Code	Species Reactivity	Applications
<u>CBR1 Antibody</u>	CSB-PA004586LA01HU	Human, Mouse	ELISA, WB, IF
<u>Cbs Antibody</u>	CSB-PA849665LA01MO	Mouse	ELISA, IF
<u>CBS Antibody</u>	CSB-PA12719A0Rb	Human, Mouse	ELISA, WB, IHC, IF, IP
<u>CBWD2 Antibody</u>	CSB-PA814202LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CBX2 Antibody</u>	CSB-PA613521LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CBX3 Antibody</u>	CSB-PA613383LA01HU	Human	ELISA, IF
<u>CBX4 Antibody</u>	CSB-PA004600EA01HU	Human	ELISA, IHC, IF
<u>CBX5 Antibody</u>	CSB-PA004601HA01HU	Human	ELISA, WB, IHC
<u>CBY1 Antibody</u>	CSB-PA896889LA01HU	Human	ELISA, WB
<u>CCDC103 Antibody</u>	CSB-PA812877LA01HU	Human, Mouse	ELISA, WB
<u>CCDC106 Antibody</u>	CSB-PA863135LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CCDC112 Antibody</u>	CSB-PA818774LA01HU	Human	ELISA, WB, IHC, IF
<u>CCDC113 Antibody</u>	CSB-PA880934LA01HU	Human	ELISA, IHC
<u>CCDC116 Antibody</u>	CSB-PA812904LA01HU	Human	ELISA, WB, IHC
<u>CCDC127 Antibody</u>	CSB-PA842643LA01HU	Human, Mouse	ELISA, WB
<u>CCDC134 Antibody</u>	CSB-PA863989EA01HU	Human	ELISA, IHC
<u>CCDC14 Antibody</u>	CSB-PA672101EA01HU	Human	ELISA, WB, IHC
<u>CCDC150 Antibody</u>	CSB-PA836720LA01HU	Human	ELISA, IHC
<u>CCDC155 Antibody</u>	CSB-PA818716LA01HU	Human	ELISA
<u>CCDC170 Antibody</u>	CSB-PA815585LA01HU	Human	ELISA, WB
<u>CCDC173 Antibody</u>	CSB-PA613149LA01HU	Human	ELISA
<u>CCDC181 Antibody</u>	CSB-PA725903LA01HU	Human	ELISA, IHC
<u>CCDC184 Antibody</u>	CSB-PA681476LA01HU	Human	ELISA, IHC
<u>CCDC24 Antibody</u>	CSB-PA836665LA01HU	Human	ELISA, WB, IHC
<u>CCDC26 Antibody</u>	CSB-PA823440LA01HU	Human, Rat	ELISA, WB
<u>CCDC28B Antibody</u>	CSB-PA866287LA01HU	Human	ELISA, IHC
<u>CCDC32 Antibody</u>	CSB-PA883607LA01HU	Human	ELISA, WB
<u>CCDC36 Antibody</u>	CSB-PA816910LA01HU	Human	ELISA, IHC
<u>CCDC42 Antibody</u>	CSB-PA842714LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CCDC57 Antibody</u>	CSB-PA652685LA01HU	Human	ELISA, WB, IHC
<u>CCDC59 Antibody</u>	CSB-PA882165LA01HU	Human	ELISA, IHC
<u>CCDC6 Antibody</u>	CSB-PA623068LA01HU	Human	ELISA, WB, IHC, IF
<u>CCDC66 Antibody</u>	CSB-PA004728LA01HU	Human	ELISA, IHC
<u>CCDC68 Antibody</u>	CSB-PA880954LA01HU	Human	ELISA, IHC
<u>CCDC8 Antibody</u>	CSB-PA880938LA01HU	Human, Mouse	ELISA, WB, IF
<u>CCDC81 Antibody</u>	CSB-PA744397LA01HU	Human	ELISA, IHC
<u>CCDC82 Antibody</u>	CSB-PA847633LA01HU	Human	ELISA
<u>CCDC84 Antibody</u>	CSB-PA771438LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>CCDC87 Antibody</u>	CSB-PA878896LA01HU	Human	ELISA
<u>CCDC89 Antibody</u>	CSB-PA818737LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CCDC91 Antibody</u>	CSB-PA801816LA01HU	Human	ELISA, IHC
<u>CCDC94 Antibody</u>	CSB-PA866295LA01HU	Human	ELISA, IHC
<u>CCDC97 Antibody</u>	CSB-PA836209LA01HU	Human	ELISA, IHC
<u>CCIN Antibody</u>	CSB-PA004770LA01HU	Human	ELISA, IHC, IF
<u>CCKAR Antibody</u>	CSB-PA004772LA01HU	Human	ELISA, WB, IHC, IF
<u>CCKBR Antibody</u>	CSB-PA12619A0Rb	Human	ELISA, IHC, IF
<u>Ccl12 Antibody</u>	CSB-PA09139A0Rb	Mouse	ELISA
<u>CCL13 Antibody</u>	CSB-PA08939A0Rb	Human	ELISA, IHC
<u>CCL14 Antibody</u>	CSB-PA08947A0Rb	Human	ELISA, IHC
<u>CCL15 Antibody</u>	CSB-PA614984LA01HU	Human	ELISA, IHC
<u>CCL16 Antibody</u>	CSB-PA08967A0Rb	Human	ELISA
<u>CCL17 Antibody</u>	CSB-PA856406LA01HU	Human	ELISA, IHC
<u>CCL2 Antibody</u>	CSB-PA05864A0Rb	Human	ELISA, IHC, IF
<u>CCL2 Antibody</u>	CSB-PA05865A0Rb	Human	ELISA, WB, IHC, IF
<u>Ccl2 Antibody</u>	CSB-PA09047A0Rb	Mouse	ELISA, IF
<u>CCL21 Antibody</u>	CSB-PA004785DA01HU	Human	ELISA, IHC
<u>CCL22 Antibody</u>	CSB-PA05895A0Rb	Human	ELISA, WB
<u>CCL24 Antibody</u>	CSB-PA05909A0Rb	Human	ELISA, IHC
<u>Ccl25 Antibody</u>	CSB-PA09189A0Rb	Mouse	ELISA
<u>CCL27 Antibody</u>	CSB-PA896726LA01HU	Human	ELISA, IHC
<u>CCL3 Antibody</u>	CSB-PA05939A0Rb	Human, Mouse	ELISA, WB, IHC
<u>Ccl4 Antibody</u>	CSB-PA09069A0Rb	Mouse	ELISA, IHC
<u>CCL4 Antibody</u>	CSB-PA05949A0Rb	Human	ELISA, IHC
<u>CCL5 Antibody</u>	CSB-PA05955A0Rb	Human	ELISA, WB
<u>CCL5 Antibody</u>	CSB-PA05959A0Rb	Human	ELISA, WB
<u>CCL7 Antibody</u>	CSB-PA05965A0Rb	Human	ELISA, WB
<u>CCL8 Antibody</u>	CSB-PA07629A0Rb	Human	ELISA, IHC
<u>Ccl9 Antibody</u>	CSB-PA09119A0Rb	Mouse	ELISA
<u>CCNA1 Antibody</u>	CSB-PA004804LA01HU	Human	ELISA, IHC
<u>CCNB2 Antibody</u>	CSB-PA004808LA01HU	Human	ELISA, WB, IF
<u>CCNC Antibody</u>	CSB-PA004810LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CCND1 Antibody</u>	CSB-PA004811LA01HU	Human	ELISA
<u>CCND2 Antibody</u>	CSB-PA004812LA01HU	Human, Mouse	ELISA, WB, IF
<u>CCNE1 Antibody</u>	CSB-PA004817LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CCNH Antibody</u>	CSB-PA004822LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>CCNJL Antibody</u>	CSB-PA004826LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CCNO Antibody</u>	CSB-PA004831LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>CCNT1 Antibody</u>	CSB-PA004832LA01HU	Human	ELISA, IHC
<u>CCNY Antibody</u>	CSB-PA847689LA01HU	Human	ELISA, WB, IHC, IF
<u>CCR1 Antibody</u>	CSB-PA004839LA01HU	Human	ELISA, IHC, IF
<u>CCR10 Antibody</u>	CSB-PA004840LA01HU	Human	ELISA, WB, IHC, IF
<u>CCR2 Antibody</u>	CSB-PA004841LA01HU	Human	ELISA, IF
<u>CCR3 Antibody</u>	CSB-PA004842LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CCR4 Antibody</u>	CSB-PA004843LA01HU	Human	ELISA, IHC, IF
<u>Ccr5 Antibody</u>	CSB-PA004844LA01MO	Mouse	ELISA
<u>CCR5 Antibody</u>	CSB-PA004844LA01HU	Human, Mouse	ELISA, WB, IF
<u>CCR6 Antibody</u>	CSB-PA004845EA01HU	Human	ELISA, WB
<u>CCR7 Antibody</u>	CSB-PA004846LA01HU	Human	ELISA, IHC, IF
<u>CCRL2 Antibody</u>	CSB-PA004852LA01HU	Human	ELISA, WB, IHC
<u>CCS Antibody</u>	CSB-PA004854LA01HU	Human	ELISA, IHC
<u>CCT5 Antibody</u>	CSB-PA004862EA01HU	Human	ELISA, WB, IF
<u>CCT6A Antibody</u>	CSB-PA004865EA01HU	Human	ELISA, WB
<u>CCT6B Antibody</u>	CSB-PA846076LA01HU	Human	ELISA, IHC
<u>CCT7 Antibody</u>	CSB-PA859112LA01HU	Human	ELISA, WB, IHC, IF
<u>CCT8 Antibody</u>	CSB-PA004873LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CCZ1 Antibody</u>	CSB-PA004156LA01HU	Human, Rat	ELISA, WB
<u>CD109 Antibody</u>	CSB-PA740988LA01HU	Human	ELISA, IHC, IF
<u>CD14 Antibody</u>	CSB-PA07085A0Rb	Human	ELISA, IHC, IF
<u>CD14 Monoclonal Antibody</u>	CSB-MA004879A0m	Human, Mouse, Rabbit	ELISA, WB, IHC, IF, FC
<u>CD14 Monoclonal Antibody</u>	CSB-MA004879A1m	Human, Mouse, Rat, Rabbit	ELISA, WB, IHC, IF
<u>CD146 Monoclonal Antibody</u>	CSB-MA013563A0m	Human	ELISA, WB, IF, FC
<u>CD151 Antibody</u>	CSB-PA004880LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CD160 Antibody</u>	CSB-PA004881EA01HU	Human	ELISA, IHC
<u>Cd163 Antibody</u>	CSB-PA650452LA01MO	Mouse, Human	ELISA, WB
<u>CD163 Antibody</u>	CSB-PA801238LA01HU	Human	ELISA, IF
<u>CD164 Antibody</u>	CSB-PA22579A0Rb	Human	ELISA, IHC, IF
<u>CD177 Antibody</u>	CSB-PA850816LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CD19 Monoclonal Antibody</u>	CSB-MA004888A0m	Human	ELISA, WB, IHC, IF, FC
<u>CD19 Monoclonal Antibody</u>	CSB-MA004888A1m	Human	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			FC
<u>CD1C Antibody</u>	CSB-PA004891LA01HU	Human, Mouse	ELISA, WB, IF
<u>CD1E Antibody</u>	CSB-PA004893LA01HU	Human	ELISA, IHC
<u>CD2 Antibody</u>	CSB-PA004894LA01HU	Human	ELISA, IHC, IF
<u>CD200 Antibody</u>	CSB-PA004895EA01HU	Human	ELISA, IF
<u>Cd207 Antibody</u>	CSB-PA837695LA01MO	Mouse	ELISA
<u>CD207 Antibody</u>	CSB-PA892321LA01HU	Human	ELISA, IHC, IF
<u>CD209 Antibody</u>	CSB-PA873604LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CD24 Antibody</u>	CSB-PA10468A0Rb	Human	ELISA, IF
<u>CD247 Antibody</u>	CSB-PA10499A0Rb	Human	ELISA, WB, IHC, IF, IP
<u>CD27 Antibody</u>	CSB-PA08235A0Rb	Human, Rat	ELISA, WB, IHC, IF
<u>CD2AP Antibody</u>	CSB-PA004914LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>CD300LB Antibody</u>	CSB-PA004919LA01HU	Human	ELISA
<u>CD300LF Antibody</u>	CSB-PA819470LA01HU	Human	ELISA, IHC
<u>CD300LG Antibody</u>	CSB-PA751016LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>CD302 Antibody</u>	CSB-PA818249LA01HU	Human	ELISA
<u>CD31 Monoclonal Antibody</u>	CSB-MA017767A0m	Human	ELISA, WB, IHC, IF, FC
<u>CD33 Antibody</u>	CSB-PA004925LA01HU	Human	ELISA, WB, IHC
<u>CD34 Antibody</u>	CSB-PA10509A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>CD36 Antibody</u>	CSB-PA09479A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>CD37 Antibody</u>	CSB-PA004928LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CD3E Antibody</u>	CSB-PA10484A0Rb	Human	ELISA, WB, IHC, IF
<u>CD3E Antibody</u>	CSB-PA004931HA01HU	Human	ELISA, IHC, IF
<u>CD3EAP Antibody</u>	CSB-PA004932LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>CD3G Antibody</u>	CSB-PA004933LA01HU	Human, Mouse, Rat	ELISA, WB, IF
<u>CD40 Antibody</u>	CSB-PA004936LA01HU	Human, Mouse	ELISA, WB, IF
<u>CD40LG Antibody</u>	CSB-PA06005A0Rb	Human	ELISA
<u>CD44 Monoclonal Antibody</u>	CSB-MA004938A0m	Human	ELISA, WB, IHC, IF, FC
<u>CD46 Antibody</u>	CSB-PA004939EA01HU	Human, Mouse	ELISA, WB, IHC
<u>Cd46 Antibody</u>	CSB-PA004939LA01MO	Mouse	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>CD47 Antibody</u>	CSB-PA004940LA01HU	Human	ELISA, IHC, IF
<u>CD4L-2 Antibody</u>	CSB-PA29721A0Rb	Grass carp	ELISA
<u>CD4-like protein Antibody</u>	CSB-PA29711A0Rb	Grass carp	ELISA
<u>CD5 Antibody</u>	CSB-PA004942LA01HU	Human	ELISA, IF
<u>CD55 Antibody</u>	CSB-PA004945LA01HU	Human	ELISA, IHC, IF
<u>CD59 Antibody</u>	CSB-PA004947EA01HU	Human	ELISA, IHC, IF
<u>CD59 Antibody</u>	CSB-PA004947YA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>Cd63 Antibody</u>	CSB-PA004950EA01MO	Mouse	ELISA
<u>Cd68 Antibody</u>	CSB-PA004951LA01MO	Mouse	ELISA
<u>CD7 Antibody</u>	CSB-PA004953YA01HU	Human	ELISA, IHC
<u>CD70 Antibody</u>	CSB-PA08009A0Rb	Human	ELISA, IF
<u>CD72 Antibody</u>	CSB-PA004955LA01HU	Human	ELISA, WB, IF
<u>Cd74 Antibody</u>	CSB-PA004956LA01MO	Mouse	ELISA, WB
<u>Cd74 Antibody</u>	CSB-PA004956LA01RA	Rat, Mouse	ELISA, WB
<u>CD79A Antibody</u>	CSB-PA004957LA01HU	Human	ELISA, IHC, IF
<u>CD79B Antibody</u>	CSB-PA004958LA01HU	Human	ELISA, IF
<u>Cd81 Antibody</u>	CSB-PA004960DA01MO	Mouse	ELISA
<u>CD81 Antibody</u>	CSB-PA004960NA01HU	Human	ELISA, IHC
<u>Cd82 Antibody</u>	CSB-PA004961LA01MO	Mouse, Human	ELISA, WB
<u>CD82 Antibody</u>	CSB-PA13269A0Rb	Human	ELISA, IHC, IF
<u>CD83 Antibody</u>	CSB-PA004962LA01HU	Human	ELISA, IF
<u>CD86 Antibody</u>	CSB-PA004965EA01HU	Human	ELISA
<u>CD8A Antibody</u>	CSB-PA004966LA01HU	Human	ELISA, IHC
<u>CD9 Antibody</u>	CSB-PA10559A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>CD96 Antibody</u>	CSB-PA004971LA01HU	Human	ELISA, WB, IHC, IF
<u>CD99L2 Antibody</u>	CSB-PA840569LA01HU	Human	ELISA, WB, IHC, IF
<u>CDA Antibody</u>	CSB-PA00845A0Rb	Human	ELISA, WB, IHC, IF
<u>CDADC1 Antibody</u>	CSB-PA874826LA01HU	Human	ELISA, IHC
<u>CDAN1 Antibody</u>	CSB-PA810285LA01HU	Human	ELISA, IHC, IF
<u>CDC123 Antibody</u>	CSB-PA004984LA01HU	Human	ELISA, IHC
<u>CDC23 Antibody</u>	CSB-PA891530LA01HU	Human	ELISA, IF
<u>CDC25A Antibody</u>	CSB-PA12709A0Rb	Human	ELISA, IF
<u>CDC25B Antibody</u>	CSB-PA004995LA01HU	Human	ELISA, WB, IHC, IF
<u>CDC25C Antibody</u>	CSB-PA004996LA01HU	Human	ELISA, IHC
<u>CDC26 Antibody</u>	CSB-PA822822LA01HU	Human	ELISA, IHC
<u>CDC27 Antibody</u>	CSB-PA004999LA01HU	Human	ELISA, IF



Product Name	Code	Species Reactivity	Applications
<u>CDC34 Antibody</u>	CSB-PA005004LA01HU	Human	ELISA, IHC
<u>CDC37 Antibody</u>	CSB-PA613684LA01HU	Human	ELISA, WB, IHC, IF
<u>CDC42 Antibody</u>	CSB-PA005008DA01HU	Human	ELISA, WB
<u>CDC42BPB Antibody</u>	CSB-PA22709A0Rb	Human	ELISA, IHC, IF
<u>CDC42EP2 Antibody</u>	CSB-PA005013LA01HU	Human	ELISA, IHC
<u>CDC42EP3 Antibody</u>	CSB-PA890918LA01HU	Human	ELISA, IHC
<u>CDC6 Antibody</u>	CSB-PA857873LA01HU	Human	ELISA, WB, IHC
<u>CDC73 Antibody</u>	CSB-PA750783LA01HU	Human	ELISA, IF
<u>CDCA7L Antibody</u>	CSB-PA850301LA01HU	Human	ELISA, IHC
<u>CDH1 Antibody</u>	CSB-PA06239A0Rb	Human	ELISA, IHC, IF
<u>CDH1 Antibody</u>	CSB-PA005034MA01HU	Human	ELISA, WB, IHC, IF, IP
<u>CDH-1 Antibody</u>	CSB-PA312474YA01EUK	Phanerochaete chrysosporium	ELISA
<u>CDH1 Monoclonal Antibody</u>	CSB-MA005034A0m	Human	ELISA, WB, IF, FC
<u>CDH11 Antibody</u>	CSB-PA005036EA01HU	Human	ELISA, IHC, IF
<u>CDH12 Antibody</u>	CSB-PA13109A0Rb	Human	ELISA, WB, IHC, IF
<u>CDH15 Antibody</u>	CSB-PA005040LA01HU	Human, Rat	ELISA, WB, IF
<u>CDH18 Antibody</u>	CSB-PA615709LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CDH2 Antibody</u>	CSB-PA005045LA01HU	Human, Mouse	ELISA, WB, IF
<u>CDH26 Antibody</u>	CSB-PA818250LA01HU	Human	ELISA
<u>CDH3 Antibody</u>	CSB-PA10715A0Rb	Human	ELISA, IHC, IF
<u>CDH5 Antibody</u>	CSB-PA07519A0Rb	Human, Mouse	ELISA, WB, IF
<u>CDH7 Antibody</u>	CSB-PA805303LA01HU	Human	ELISA, WB, IHC, IF
<u>CDH8 Antibody</u>	CSB-PA005057LA01HU	Human	ELISA, IF
<u>CDH9 Antibody</u>	CSB-PA892342LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CDK1 Antibody</u>	CSB-PA361848LA01HU	Human	ELISA, IHC
<u>CDK11A Antibody</u>	CSB-PA05204A0Rb	Human, Rat	ELISA, WB, IHC, IF
<u>CDK12 Antibody</u>	CSB-PA882147LA01HU	Human	ELISA, WB
<u>CDK13 Antibody</u>	CSB-PA613482LA01HU	Human	ELISA, IHC
<u>CDK17 Antibody</u>	CSB-PA28369A0Rb	Human	ELISA, WB, IHC
<u>CDK2 Antibody</u>	CSB-PA005061LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>CDK2AP1 Antibody</u>	CSB-PA005062DA01HU	Human	ELISA, IHC
<u>CDK2AP2 Antibody</u>	CSB-PA005063LA01HU	Human	ELISA, IHC
<u>CDK3 Antibody</u>	CSB-PA005064LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>CDK4 Antibody</u>	CSB-PA11659A0Rb	Human	ELISA, IHC
<u>CDK5 Antibody</u>	CSB-PA005067LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CDK5RAP2 Antibody</u>	CSB-PA847226LA01HU	Human	ELISA, IHC, IF
<u>CDK5RAP3 Antibody</u>	CSB-PA842691LA01HU	Human, Rat	ELISA, WB, IF
<u>CDK7 Antibody</u>	CSB-PA005075LA01HU	Human	ELISA, IHC, IF
<u>CDKAL1 Antibody</u>	CSB-PA716574LA01HU	Human	ELISA, IHC, IF
<u>CDKB1-2 Antibody</u>	CSB-PA651993LA01DOA	Arabidopsis thaliana	ELISA
<u>CDKL5 Antibody</u>	CSB-PA005085LA01HU	Human, Rat	ELISA, WB
<u>CDKN1A Antibody</u>	CSB-PA005086DA01HU	Human	ELISA, IHC, IF
<u>CDKN1C Antibody</u>	CSB-PA005088LA01HU	Human	ELISA, IF
<u>CDKN2A Antibody</u>	CSB-PA005089LA01HU	Human	ELISA, IHC
<u>CDKN2AIPNL Antibody</u>	CSB-PA005091LA01HU	Human	ELISA, IHC
<u>CDKN2B Antibody</u>	CSB-PA005092LA01HU	Human	ELISA, IHC, IF
<u>CDKN2C Antibody</u>	CSB-PA005094LA01HU	Human	ELISA, IHC, IF
<u>CDKN3 Antibody</u>	CSB-PA005096LA01HU	Human	ELISA, IHC, IF
<u>CDON Antibody</u>	CSB-PA679696LA01HU	Human	ELISA, WB, IF
<u>CDPF1 Antibody</u>	CSB-PA764676LA01HU	Human	ELISA
<u>CDR2 Antibody</u>	CSB-PA005103LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CDR2L Antibody</u>	CSB-PA005104LA01HU	Human	ELISA, WB, IHC
<u>CDS1 Antibody</u>	CSB-PA842170LA01HU	Human	ELISA, IHC, IF
<u>CDS2 Antibody</u>	CSB-PA005123LA01HU	Human	ELISA, WB, IHC
<u>CDSN Antibody</u>	CSB-PA005124LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CDT1 Antibody</u>	CSB-PA005125LA01HU	Human	ELISA, WB
<u>CDX1 Antibody</u>	CSB-PA005127LA01HU	Human, Rat	ELISA, WB, IF
<u>CDX2 Antibody</u>	CSB-PA11699A0Rb	Human	ELISA
<u>CDX4 Antibody</u>	CSB-PA005129LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CDYL Antibody</u>	CSB-PA896685LA01HU	Human	ELISA, WB, IHC
<u>CEACAM1 Antibody</u>	CSB-PA10749A0Rb	Human	ELISA, WB, IF
<u>CEACAM16 Antibody</u>	CSB-PA005158LA01HU	Human	ELISA, IHC
<u>CEACAM4 Antibody</u>	CSB-PA005164LA01HU	Human	ELISA
<u>CEACAM5 Antibody</u>	CSB-PA005165LA01HU	Human	ELISA, WB, IHC, IF
<u>CEACAM6 Antibody</u>	CSB-PA005166LA01HU	Human	ELISA, IHC, IF
<u>CEACAM7 Antibody</u>	CSB-PA005167LA01HU	Human	ELISA, IHC
<u>CEBPA Antibody</u>	CSB-PA005180HA01HU	Human	ELISA, WB, IHC
<u>CEBPB Antibody</u>	CSB-PA005181LA01HU	Human	ELISA, IHC
<u>CEBPE Antibody</u>	CSB-PA618016LA01HU	Human, Rat	ELISA, WB
<u>CEL Antibody</u>	CSB-PA325708LA01HU	Human	ELISA, IHC
<u>Cela2a Antibody</u>	CSB-PA005196LA01MO	Mouse	ELISA, WB
<u>CELA3A Antibody</u>	CSB-PA007588LA01HU	Human, Mouse	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>CELF2 Antibody</u>	CSB-PA006217LA01HU	Human	ELISA, WB, IHC, IF
<u>CELF4 Antibody</u>	CSB-PA883643LA01HU	Human	ELISA, IHC
<u>CELF6 Antibody</u>	CSB-PA822233LA01HU	Human	ELISA, WB
<u>CELSR1 Antibody</u>	CSB-PA868360NA01HU	Human	ELISA, IHC
<u>CELSR3 Antibody</u>	CSB-PA852817LA01HU	Human	ELISA, IHC
<u>CENPA Antibody</u>	CSB-PA005205LA01HU	Human	ELISA, IF
<u>CENPB Antibody</u>	CSB-PA002171LA01HU	Human	ELISA, IHC
<u>CENPJ Antibody</u>	CSB-PA005213LA01HU	Human	ELISA, WB, IF
<u>CENPX Antibody</u>	CSB-PA022871LA01HU	Human	ELISA, IHC
<u>CEP104 Antibody</u>	CSB-PA012199LA01HU	Human	ELISA, IHC
<u>CEP126 Antibody</u>	CSB-PA868394LA01HU	Human	ELISA, IHC
<u>CEP135 Antibody</u>	CSB-PA721063LA01HU	Human	ELISA, IF
<u>CEP164 Antibody</u>	CSB-PA892491LA01HU	Human	ELISA, IHC
<u>CEP19 Antibody</u>	CSB-PA850384LA01HU	Human	ELISA, IHC
<u>CEP192 Antibody</u>	CSB-PA851560LA01HU	Human	ELISA, IHC, IF
<u>CEP250 Antibody</u>	CSB-PA871581LA01HU	Human	ELISA, IF
<u>CEP290 Antibody</u>	CSB-PA005241LA01HU	Human	ELISA, IHC, IF
<u>CEP63 Antibody</u>	CSB-PA839376LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CEP68 Antibody</u>	CSB-PA765741LA01HU	Human	ELISA
<u>CEP70 Antibody</u>	CSB-PA836778LA01HU	Human	ELISA, WB, IHC, IF
<u>CEP76 Antibody</u>	CSB-PA855035LA01HU	Human	ELISA, WB, IHC
<u>CEP83 Antibody</u>	CSB-PA896730LA01HU	Human	ELISA
<u>CEP89 Antibody</u>	CSB-PA853505LA01HU	Human	ELISA
<u>CEP97 Antibody</u>	CSB-PA808994LA01HU	Human	ELISA, IF
<u>CERK Antibody</u>	CSB-PA005256LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CERKL Antibody</u>	CSB-PA669888LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CERS1 Antibody</u>	CSB-PA012761LA01HU	Human, Rat	ELISA, WB, IHC
<u>CERS3 Antibody</u>	CSB-PA811597LA01HU	Human, Mouse	ELISA, WB
<u>CERS5 Antibody</u>	CSB-PA847636LA01HU	Human	ELISA, WB, IHC, IF
<u>CERS6 Antibody</u>	CSB-PA751145LA01HU	Human	ELISA, WB, IHC
<u>CES1 Antibody</u>	CSB-PA005258LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>Ces1c Antibody</u>	CSB-PA338557YA01MO	Mouse, Human, Rat	ELISA, WB, IHC
<u>CES3 Antibody</u>	CSB-PA005260LA01HU	Human	ELISA, IHC, IF
<u>CESA4 Antibody</u>	CSB-PA771138PA01DOA	Arabidopsis thaliana	ELISA, WB
<u>CETN2 Antibody</u>	CSB-PA005265HA01HU	Human	ELISA
<u>cfaE Antibody</u>	CSB-PA328563LA01ENL	Escherichia coli	ELISA



Product Name	Code	Species Reactivity	Applications
<u>CFAP100 Antibody</u>	CSB-PA669847LA01HU	Human	ELISA, IHC, IF
<u>CFAP157 Antibody</u>	CSB-PA686386LA01HU	Human	ELISA, IHC
<u>CFAP410 Antibody</u>	CSB-PA003758LA01HU	Human	ELISA, IHC
<u>CFAP52 Antibody</u>	CSB-PA839781LA01HU	Human	ELISA, WB, IHC
<u>CFAP61 Antibody</u>	CSB-PA850915LA01HU	Human	ELISA, IHC
<u>CFAP65 Antibody</u>	CSB-PA757963LA01HU	Human	ELISA, WB, IHC, IF
<u>CFB Antibody</u>	CSB-PA11049A0Rb	Human	ELISA, WB, IHC
<u>Cfb Antibody</u>	CSB-PA005268LA01MO	Mouse, Human	ELISA, WB
<u>CFHR1 Antibody</u>	CSB-PA13809A0Rb	Human	ELISA, IHC
<u>CFHR2 Antibody</u>	CSB-PA13817A0Rb	Human	ELISA, WB, IHC
<u>CFHR5 Antibody</u>	CSB-PA883624LA01HU	Human	ELISA, WB, IHC
<u>CFI Antibody</u>	CSB-PA005279LA01HU	Human	ELISA, WB, IHC
<u>CFL1 Antibody</u>	CSB-PA04709A0Rb	Human	ELISA, IHC, IF
<u>CFLAR Antibody</u>	CSB-PA005282LA01HU	Human	ELISA, WB, IHC, IF
<u>CGGBP1 Antibody</u>	CSB-PA890676LA01HU	Human	ELISA, IHC
<u>CGN Antibody</u>	CSB-PA005302LA01HU	Human	ELISA, IHC
<u>CH25H Antibody</u>	CSB-PA005307LA01HU	Human	ELISA, IHC
<u>CHAC1 Antibody</u>	CSB-PA883606LA01HU	Human	ELISA, IHC, IF
<u>CHAD Antibody</u>	CSB-PA005310LA01HU	Human	ELISA, IHC
<u>CHAF1A Antibody</u>	CSB-PA618758LA01HU	Human	ELISA, IHC
<u>CHAT Antibody</u>	CSB-PA005314LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CHCHD4 Antibody</u>	CSB-PA850802HA01HU	Human	ELISA, IHC
<u>chd Antibody</u>	CSB-PA005377LA01DIL	Zebrafish	ELISA
<u>CHD1L Antibody</u>	CSB-PA803136LA01HU	Human	ELISA, WB, IF
<u>CHD4 Antibody</u>	CSB-PA01097A0Rb	Human	ELISA, IHC, IF
<u>CHD5 Antibody</u>	CSB-PA848409LA01HU	Human	ELISA, IHC, IF
<u>CHDH Antibody</u>	CSB-PA839862LA01HU	Human	ELISA, IHC
<u>CHEK1 Antibody</u>	CSB-PA005338LA01HU	Human	ELISA
<u>CHEK2 Antibody</u>	CSB-PA005339LA01HU	Human	ELISA, IHC
<u>CHFR Antibody</u>	CSB-PA839315LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CHGA Antibody</u>	CSB-PA005344YA01HU	Human	ELISA, IHC
<u>CHGB Antibody</u>	CSB-PA005345LA01HU	Human, Mouse	ELISA, WB
<u>CHI3L2 Antibody</u>	CSB-PA620893LA01HU	Human	ELISA, IHC
<u>CHID1 Antibody</u>	CSB-PA883614LA01HU	Human	ELISA, WB, IHC
<u>CHKA Antibody</u>	CSB-PA005353LA01HU	Human	ELISA, WB, IHC, IF
<u>CHKB Antibody</u>	CSB-PA005354LA01HU	Human, Rat, Mouse	ELISA, WB, IHC, IF
<u>CHLP Antibody</u>	CSB-PA874935LA01DOA	Arabidopsis thaliana	ELISA



Product Name	Code	Species Reactivity	Applications
<u>CHM Antibody</u>	CSB-PA005356LA01HU	Human	ELISA, IHC, IF
<u>CHML Antibody</u>	CSB-PA005357LA01HU	Human	ELISA, IF
<u>CHMP1A Antibody</u>	CSB-PA881032LA01HU	Human	ELISA, WB, IHC, IF
<u>CHMP2A Antibody</u>	CSB-PA04524A0Rb	Human	ELISA, WB, IHC
<u>CHMP4C Antibody</u>	CSB-PA853396LA01HU	Human	ELISA, IHC
<u>CHMP5 Antibody</u>	CSB-PA878944LA01HU	Human	ELISA, IHC, IF
<u>CHMP6 Antibody</u>	CSB-PA846622LA01HU	Human, Rat	ELISA, WB, IHC
<u>CHMP7 Antibody</u>	CSB-PA855478LA01HU	Human	ELISA, IHC
<u>CHN2 Antibody</u>	CSB-PA005369LA01HU	Human	ELISA, IHC
<u>CHN48 Antibody</u>	CSB-PA362364LA01NHE	Nicotiana tabacum	ELISA, WB
<u>choD Antibody</u>	CSB-PA18185A0Rb	Mycobacterium tuberculosis	ELISA
<u>choD Antibody</u>	CSB-PA704121EA01MVZ	Mycobacterium tuberculosis	ELISA
<u>CHODL Antibody</u>	CSB-PA005370LA01HU	Human	ELISA, WB, IHC, IF
<u>CHP1 Antibody</u>	CSB-PA860773LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>CHRD1 Antibody</u>	CSB-PA866283LA01HU	Human	ELISA, WB, IHC
<u>Chrm1 Antibody</u>	CSB-PA005381HA01RA	Rat	ELISA, WB
<u>CHRM1 Antibody</u>	CSB-PA005381PA01HU	Human	ELISA, IF
<u>CHRM2 Antibody</u>	CSB-PA005382PA01HU	Human	ELISA, IF
<u>CHRM3 Antibody</u>	CSB-PA005383LA01HU	Human	ELISA, IHC, IF
<u>CHRM5 Antibody</u>	CSB-PA005385LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>Chrna1 Antibody</u>	CSB-PA005386LA01MO	Mouse, Human	ELISA, WB
<u>CHRNA1 Antibody</u>	CSB-PA005386LA01TOT	Torpedo californica	ELISA, WB
<u>CHRNA1 Antibody</u>	CSB-PA11789A0Rb	Human, Rat	ELISA, WB, IHC
<u>CHRN1 Antibody</u>	CSB-PA15049A0Rb	Human	ELISA, WB, IHC, IF
<u>CHRN2 Antibody</u>	CSB-PA005396LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>CHRN3 Antibody</u>	CSB-PA005397LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CHRND Antibody</u>	CSB-PA005399LA01HU	Human	ELISA, IF
<u>CHS1 Antibody</u>	CSB-PA16277A0Rb	Baker's yeast	ELISA
<u>CHST1 Antibody</u>	CSB-PA005402LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>CHST10 Antibody</u>	CSB-PA005403LA01HU	Human	ELISA
<u>CHST11 Antibody</u>	CSB-PA882066LA01HU	Human	ELISA, IHC
<u>CHST14 Antibody</u>	CSB-PA839851LA01HU	Human	ELISA, WB, IHC
<u>CHST3 Antibody</u>	CSB-PA759175LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>CHST9 Antibody</u>	CSB-PA745332LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CHSY1 Antibody</u>	CSB-PA005416LA01HU	Human	ELISA, IHC
<u>CHTOP Antibody</u>	CSB-PA897109HA01HU	Human, Mouse	ELISA, WB, IHC
<u>CIAO2A Antibody</u>	CSB-PA872502LA01HU	Human, Mouse	ELISA, WB
<u>CIAO3 Antibody</u>	CSB-PA863991LA01HU	Human	ELISA, WB, IHC, IF
<u>CIAPIN1 Antibody</u>	CSB-PA738750HA01HU	Human	ELISA
<u>CIART Antibody</u>	CSB-PA850790LA01HU	Human	ELISA, IHC
<u>CIB2 Antibody</u>	CSB-PA005427LA01HU	Human	ELISA, IF
<u>CIDEA Antibody</u>	CSB-PA005431LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CIDEB Antibody</u>	CSB-PA12869A0Rb	Human	ELISA, IHC
<u>CINP Antibody</u>	CSB-PA861151LA01HU	Human	ELISA
<u>CIPC Antibody</u>	CSB-PA880155LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CISD2 Antibody</u>	CSB-PA839803LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CKAP2L Antibody</u>	CSB-PA809015LA01HU	Human	ELISA, IHC, IF
<u>CKAP4 Antibody</u>	CSB-PA005453LA01HU	Human	ELISA, WB, IF
<u>CKAP5 Antibody</u>	CSB-PA617899LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>CKB Antibody</u>	CSB-PA01495A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>CKM Antibody</u>	CSB-PA09809A0Rb	Human, Mouse	ELISA, WB, IHC
<u>CKMT1A Antibody</u>	CSB-PA09817A0Rb	Human	ELISA, IHC, IF
<u>CKS2 Antibody</u>	CSB-PA005468LA01HU	Human	ELISA, IHC
<u>CLASP2 Antibody</u>	CSB-PA005472LA01HU	Human	ELISA, IHC, IP
<u>CLC Antibody</u>	CSB-PA22679A0Rb	Human	ELISA, IHC, IF
<u>CLCA1 Antibody</u>	CSB-PA005474LA01HU	Human	ELISA, IHC
<u>CLCC1 Antibody</u>	CSB-PA857024LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>CLCF1 Antibody</u>	CSB-PA906951LA01HU	Human	ELISA, IHC
<u>CLCN1 Antibody</u>	CSB-PA005480LA01HU	Human	ELISA, IF
<u>CLCN2 Antibody</u>	CSB-PA005481LA01HU	Human	ELISA, IF
<u>CLCN5 Antibody</u>	CSB-PA005484LA01HU	Human	ELISA, IHC
<u>CLCNKB Antibody</u>	CSB-PA005488LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>CLDN1 Antibody</u>	CSB-PA005490LA01HU	Human	ELISA, IHC, IF
<u>CLDN10 Antibody</u>	CSB-PA005491LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CLDN12 Antibody</u>	CSB-PA005493LA01HU	Human	ELISA, IHC
<u>CLDN15 Antibody</u>	CSB-PA005495EA01HU	Human	ELISA, IHC, IF
<u>CLDN20 Antibody</u>	CSB-PA005501LA01HU	Human	ELISA, WB, IF



Product Name	Code	Species Reactivity	Applications
<u>CLDN3 Antibody</u>	CSB-PA13119A0Rb	Human	ELISA, WB, IHC
<u>CLDN5 Antibody</u>	CSB-PA005507LA01HU	Human	ELISA, IHC, IF
<u>CLDN6 Antibody</u>	CSB-PA005508LA01HU	Human	ELISA, IHC, IF
<u>CLDND1 Antibody</u>	CSB-PA873687LA01HU	Human	ELISA, IF
<u>CLEC10A Antibody</u>	CSB-PA811601LA01HU	Human	ELISA, IHC, IF
<u>CLEC14A Antibody</u>	CSB-PA801222LA01HU	Human	ELISA, IHC
<u>CLEC18A Antibody</u>	CSB-PA005521LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CLEC1A Antibody</u>	CSB-PA843319LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CLEC2D Antibody</u>	CSB-PA883401LA01HU	Human	ELISA
<u>CLEC3A Antibody</u>	CSB-PA005530LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CLEC4C Antibody</u>	CSB-PA855470LA01HU	Human	ELISA, IHC
<u>CLEC4D Antibody</u>	CSB-PA841219LA01HU	Human	ELISA, IHC
<u>CLEC4E Antibody</u>	CSB-PA891553LA01HU	Human	ELISA, IHC, IF
<u>CLEC4G Antibody</u>	CSB-PA744269LA01HU	Human	ELISA, IHC
<u>CLEC4M Antibody</u>	CSB-PA862003LA01HU	Human, Mouse	ELISA, WB, IF
<u>CLEC7A Antibody</u>	CSB-PA883623LA01HU	Human	ELISA, IHC, IF
<u>CLEC9A Antibody</u>	CSB-PA740922HA01HU	Human, Mouse	ELISA, WB, IF
<u>clfB Antibody</u>	CSB-PA753442HA01SKV	Staphylococcus aureus	ELISA
<u>CLGN Antibody</u>	CSB-PA005544LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CLIC2 Antibody</u>	CSB-PA005547LA01HU	Human	ELISA, WB, IHC, IF
<u>CLINT1 Antibody</u>	CSB-PA005552LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>CLIP2 Antibody</u>	CSB-PA871393LA01HU	Human, Mouse, Rat	ELISA, WB, IF
<u>CLN5 Antibody</u>	CSB-PA005567LA01HU	Human	ELISA, IHC, IF
<u>CLN6 Antibody</u>	CSB-PA882125LA01HU	Human	ELISA, IHC
<u>CLN8 Antibody</u>	CSB-PA866210LA01HU	Human	ELISA, IHC
<u>CLOCK Antibody</u>	CSB-PA005574EA01HU	Human	ELISA, WB, IHC, IF
<u>CLPB Antibody</u>	CSB-PA884428LA01HU	Human	ELISA, WB, IHC
<u>clpP Antibody</u>	CSB-PA401627LA01CWV	Clostridium botulinum	ELISA
<u>CLPP Antibody</u>	CSB-PA618094LA01HU	Human	ELISA, IHC
<u>clpS Antibody</u>	CSB-PA16649A0Rb	Escherichia coli	ELISA, WB
<u>CLPS Antibody</u>	CSB-PA005578LA01HU	Human	ELISA, IHC
<u>CLPTMIL Antibody</u>	CSB-PA850323LA01HU	Human	ELISA, WB, IHC, IF
<u>CLRN2 Antibody</u>	CSB-PA005584LA01HU	Human	ELISA, IHC, IF
<u>CLRN3 Antibody</u>	CSB-PA818764LA01HU	Human	ELISA, IF
<u>CLSTN2 Antibody</u>	CSB-PA863978LA01HU	Human	ELISA, IHC, IF
<u>CLSTN3 Antibody</u>	CSB-PA871555LA01HU	Human	ELISA, IF
<u>CLTB Antibody</u>	CSB-PA03147A0Rb	Human	ELISA, IHC
<u>CLTCL1 Antibody</u>	CSB-PA005594LA01HU	Human, Mouse	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>CLU Antibody</u>	CSB-PA06029A0Rb	Human	ELISA, IHC, IF
<u>CLUAP1 Antibody</u>	CSB-PA836181LA01HU	Human	ELISA, IHC
<u>CLVS2 Antibody</u>	CSB-PA713118HA01HU	Human, Mouse	ELISA, WB
<u>CLYBL Antibody</u>	CSB-PA850772LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CMA1 Antibody</u>	CSB-PA005599LA01HU	Human, Mouse	ELISA, WB
<u>Cma1 Antibody</u>	CSB-PA005599LA01MO	Mouse	ELISA
<u>CMBL Antibody</u>	CSB-PA846604LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CMC1 Antibody</u>	CSB-PA005603LA01HU	Human	ELISA, IHC
<u>CMC4 Antibody</u>	CSB-PA349766LA01HU	Human	ELISA, IHC
<u>CMIP Antibody</u>	CSB-PA815577LA01HU	Human	ELISA, IHC
<u>cmk Antibody</u>	CSB-PA363865HA01EVN	Escherichia coli	ELISA
<u>CMKLR1 Antibody</u>	CSB-PA005614LA01HU	Human	ELISA, IHC, IF
<u>CMTM2 Antibody</u>	CSB-PA844708LA01HU	Human	ELISA, IF
<u>CMTM5 Antibody</u>	CSB-PA842660LA01HU	Human	ELISA, IHC
<u>CMTM8 Antibody</u>	CSB-PA809032LA01HU	Human	ELISA, IF
<u>CMTR2 Antibody</u>	CSB-PA812902LA01HU	Human	ELISA, IHC, IF
<u>c-MYC Monoclonal Antibody</u>	CSB-MA000041M0m	All	ELISA, WB, IF, IP
<u>CNBP Antibody</u>	CSB-PA005637HA01HU	Human	ELISA, IHC
<u>CNDP2 Antibody</u>	CSB-PA856967HA01HU	Human	ELISA, IHC
<u>CNFN Antibody</u>	CSB-PA005641LA01HU	Human	ELISA, IHC, IF
<u>CNGA1 Antibody</u>	CSB-PA005642LA01HU	Human	ELISA, IF
<u>CNGA4 Antibody</u>	CSB-PA808540LA01HU	Human	ELISA, IHC
<u>CNGB1 Antibody</u>	CSB-PA623897LA01HU	Human	ELISA, WB, IHC, IF
<u>CNIH2 Antibody</u>	CSB-PA754273LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CNIH3 Antibody</u>	CSB-PA848394LA01HU	Human	ELISA, IHC, IF
<u>CNIH4 Antibody</u>	CSB-PA885782LA01HU	Human	ELISA, WB, IHC, IF
<u>CNKSR2 Antibody</u>	CSB-PA005653LA01HU	Human	ELISA, IHC, IF
<u>CNKSR3 Antibody</u>	CSB-PA757602LA01HU	Human, Mouse	ELISA, WB, IF
<u>CNMD Antibody</u>	CSB-PA012854LA01HU	Human	ELISA
<u>CNN2 Antibody</u>	CSB-PA860764LA01HU	Human, Mouse	ELISA, WB, IF
<u>CNN3 Antibody</u>	CSB-PA005657LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CNNM1 Antibody</u>	CSB-PA878873LA01HU	Human	ELISA, IF
<u>CNNM2 Antibody</u>	CSB-PA872518LA01HU	Human	ELISA, IHC, IF
<u>CNNM3 Antibody</u>	CSB-PA843333LA01HU	Human	ELISA, WB, IHC, IF
<u>CNNM4 Antibody</u>	CSB-PA764727LA01HU	Human	ELISA, IHC, IF
<u>CNOT1 Antibody</u>	CSB-PA005663LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>CNOT10 Antibody</u>	CSB-PA884495LA01HU	Human	ELISA, IHC
<u>CNOT2 Antibody</u>	CSB-PA865193LA01HU	Human	ELISA, WB
<u>CNOT3 Antibody</u>	CSB-PA005666LA01HU	Human	ELISA, IHC, IF
<u>CNOT9 Antibody</u>	CSB-PA842140LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CNP Antibody</u>	CSB-PA21939A0Rb	Human	ELISA, IHC
<u>CNPY2 Antibody</u>	CSB-PA05034A0Rb	Human	ELISA, WB, IHC
<u>CNPY4 Antibody</u>	CSB-PA843134LA01HU	Human	ELISA, IHC
<u>CNR2 Antibody</u>	CSB-PA005679LA01HU	Human	ELISA, IHC, IF
<u>CNTF Antibody</u>	CSB-PA06039A0Rb	Human	ELISA, IHC, IF
<u>CNTRF Antibody</u>	CSB-PA005684EA01HU	Human	ELISA, WB, IHC, IF
<u>CNTN1 Antibody</u>	CSB-PA623657LA01HU	Human	ELISA, IF
<u>CNTN6 Antibody</u>	CSB-PA891801LA01HU	Human	ELISA, WB, IHC
<u>CNTNAP1 Antibody</u>	CSB-PA005692LA01HU	Human	ELISA, IHC
<u>COA1 Antibody</u>	CSB-PA004170LA01HU	Human	ELISA, IHC, IF
<u>COA3 Antibody</u>	CSB-PA896698LA01HU	Human	ELISA, IHC
<u>COA6 Antibody</u>	CSB-PA683360LA01HU	Human	ELISA, IHC
<u>COASY Antibody</u>	CSB-PA618754LA01HU	Human	ELISA, WB, IF
<u>cobB Antibody</u>	CSB-PA536037LA01XAZ	Xanthomonas campestris pv. campestris	ELISA, WB
<u>COBLL1 Antibody</u>	CSB-PA706643LA01HU	Human	ELISA, IHC, IF
<u>Cobra venom factor Antibody</u>	CSB-PA856221LA01NAF	Naja kaouthia	ELISA, WB
<u>COG3 Antibody</u>	CSB-PA839348LA01HU	Human	ELISA, IHC
<u>COG6 Antibody</u>	CSB-PA23419A0Rb	Human	ELISA, IHC
<u>COG8 Antibody</u>	CSB-PA846673LA01HU	Human	ELISA, IHC
<u>COL11A1 Antibody</u>	CSB-PA005716LA01HU	Human	ELISA, IHC
<u>COL12A1 Antibody</u>	CSB-PA005719LA01CH	Chicken	ELISA
<u>COL12A1 Antibody</u>	CSB-PA859521LA01HU	Human	ELISA, IHC
<u>COL17A1 Antibody</u>	CSB-PA18099A0Rb	Human	ELISA, IHC
<u>Col1a1 Antibody</u>	CSB-PA005727YA01RA	Rat	ELISA, IHC
<u>COL1A1 Antibody</u>	CSB-PA005727OA01HU	Human	ELISA, IF
<u>COL20A1 Antibody</u>	CSB-PA873727LA01HU	Human	ELISA, IHC
<u>COL3A1 Antibody</u>	CSB-PA11459A0Rb	Human	ELISA
<u>COL4A1 Antibody</u>	CSB-PA11464A0Rb	Human	ELISA, IHC, IF
<u>COL4A1 Antibody</u>	CSB-PA11469A0Rb	Human	ELISA, IHC, IF
<u>COL4A2 Antibody</u>	CSB-PA11947A0Rb	Human	ELISA, IHC
<u>COL4A3 Antibody</u>	CSB-PA005743LA01HU	Human	ELISA, IHC
<u>COL4A3BP Antibody</u>	CSB-PA896912LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>COL4A5 Antibody</u>	CSB-PA005746LA01HU	Human	ELISA, IHC
<u>COL6A1 Antibody</u>	CSB-PA005751NA01HU	Human	ELISA, IF
<u>COL6A1 Antibody</u>	CSB-PA005751OA01HU	Human	ELISA, IF



Product Name	Code	Species Reactivity	Applications
<u>COL6A2 Antibody</u>	CSB-PA005752LA01HU	Human	ELISA, IHC
<u>COL6A6 Antibody</u>	CSB-PA409861LA01HU	Human	ELISA, IHC
<u>COL6A6 Antibody</u>	CSB-PA409861YA01HU	Human	ELISA, IHC
<u>COL8A1 Antibody</u>	CSB-PA005755LA01HU	Human, Mouse	ELISA, WB, IF
<u>COL9A2 Antibody</u>	CSB-PA617900OA01HU	Human, Rat	ELISA, WB
<u>COLEC10 Antibody</u>	CSB-PA896556LA01HU	Human	ELISA, IHC
<u>COLGALT2 Antibody</u>	CSB-PA818260LA01HU	Human, Mouse	ELISA, WB, IHC
<u>COMMD1 Antibody</u>	CSB-PA818712LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>COMMD10 Antibody</u>	CSB-PA809070LA01HU	Human	ELISA, IF
<u>COMMD2 Antibody</u>	CSB-PA774827LA01HU	Human	ELISA, IHC
<u>COMMD4 Antibody</u>	CSB-PA04965A0Rb	Human	ELISA, IHC
<u>COMMD6 Antibody</u>	CSB-PA772019LA01HU	Human	ELISA, IHC
<u>COMMD9 Antibody</u>	CSB-PA882164LA01HU	Human	ELISA, IHC
<u>Conglutin Antibody</u>	CSB-PA720709HA01AVE	Peanut	ELISA
<u>COPB1 Antibody</u>	CSB-PA005783LA01HU	Human, Mouse	ELISA, WB, IHC
<u>COPE Antibody</u>	CSB-PA03669A0Rb	Human	ELISA, WB, IHC
<u>COPG2 Antibody</u>	CSB-PA866200LA01HU	Human	ELISA, IHC, IF
<u>Copper-metallothionein Antibody</u>	CSB-PA345921DA01HDB	Roman snail	ELISA, WB
<u>Copper-metallothionein Antibody</u>	CSB-PA345921EA01HDB	Roman snail	ELISA, WB
<u>COPS7A Antibody</u>	CSB-PA880068HA01HU	Human, Mouse	ELISA, WB, IHC
<u>COPS7B Antibody</u>	CSB-PA872526LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>COPS8 Antibody</u>	CSB-PA860346LA01HU	Human	ELISA, IHC
<u>COPZ1 Antibody</u>	CSB-PA005797LA01HU	Human	ELISA, WB, IHC
<u>COQ3 Antibody</u>	CSB-PA882156LA01HU	Human	ELISA, IHC
<u>COQ4 Antibody</u>	CSB-PA896885LA01HU	Human	ELISA, IHC
<u>COQ6 Antibody</u>	CSB-PA897098LA01HU	Human, Mouse	ELISA, WB, IHC
<u>COQ8B Antibody</u>	CSB-PA839305LA01HU	Human	ELISA, WB, IHC, IF
<u>CORO1A Antibody</u>	CSB-PA005813LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CORO1B Antibody</u>	CSB-PA880087LA01HU	Human	ELISA, IHC, IF
<u>CORO1C Antibody</u>	CSB-PA891968LA01HU	Human, Mouse	ELISA, WB, IF
<u>CORO7 Antibody</u>	CSB-PA005819LA01HU	Human	ELISA
<u>COTL1 Antibody</u>	CSB-PA622751HA01HU	Human, Mouse	ELISA, WB, IHC
<u>COX11 Antibody</u>	CSB-PA896759LA01HU	Human	ELISA, IHC, IF
<u>COX15 Antibody</u>	CSB-PA759167LA01HU	Human, Rat	ELISA, WB, IHC
<u>COX17 Antibody</u>	CSB-PA619858LA01HU	Human	ELISA, WB, IHC
<u>COX19 Antibody</u>	CSB-PA677720LA01HU	Human	ELISA, IF



Product Name	Code	Species Reactivity	Applications
<u>COX4I1 Antibody</u>	CSB-PA00509A0Rb	Human, Rat, Zebrafish	ELISA, WB, IHC
<u>COX5A Antibody</u>	CSB-PA00414A0Rb	Human	ELISA, WB, IHC, IF
<u>COX5B Antibody</u>	CSB-PA00404A0Rb	Human	ELISA, WB, IF
<u>COX6A1 Antibody</u>	CSB-PA02315A0Rb	Human	ELISA, WB, IHC, IF
<u>COX6B1 Antibody</u>	CSB-PA00395A0Rb	Human	ELISA, IF
<u>CP Antibody</u>	CSB-PA14599A0Rb	Human	ELISA, IHC
<u>CPA2 Antibody</u>	CSB-PA005876LA01HU	Human	ELISA, WB, IHC
<u>CPA6 Antibody</u>	CSB-PA005880LA01HU	Human	ELISA
<u>CPB2 Antibody</u>	CSB-PA005884LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CPE Antibody</u>	CSB-PA005886LA01HU	Human	ELISA, WB
<u>CPEB1 Antibody</u>	CSB-PA861184LA01HU	Human	ELISA
<u>CPEB2 Antibody</u>	CSB-PA767205LA01HU	Human, Mouse, Rat	ELISA, WB, IF
<u>CPEB3 Antibody</u>	CSB-PA818773LA01HU	Human	ELISA, IF
<u>CPEB4 Antibody</u>	CSB-PA618177LA01HU	Human	ELISA, WB, IF
<u>CPED1 Antibody</u>	CSB-PA004182LA01HU	Human	ELISA, IHC
<u>cpg2 Antibody</u>	CSB-PA361895DA01PW	Pseudomonas sp.	ELISA, WB
<u>CPLX1 Antibody</u>	CSB-PA005893LA01HU	Human	ELISA, IHC, IF
<u>CPLX3 Antibody</u>	CSB-PA851972LA01HU	Human	ELISA
<u>CPLX4 Antibody</u>	CSB-PA767213LA01HU	Human	ELISA
<u>CPNE1 Antibody</u>	CSB-PA859524LA01HU	Human	ELISA, IHC
<u>CPNE2 Antibody</u>	CSB-PA839323LA01HU	Human	ELISA, IHC
<u>CPNE3 Antibody</u>	CSB-PA005902LA01HU	Human	ELISA, WB, IHC, IF
<u>CPNE4 Antibody</u>	CSB-PA850253LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>CPNE7 Antibody</u>	CSB-PA866203LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>CPNE8 Antibody</u>	CSB-PA773062LA01HU	Human	ELISA, IHC
<u>CPNE9 Antibody</u>	CSB-PA812898LA01HU	Human	ELISA, IHC
<u>CPO Antibody</u>	CSB-PA005909LA01HU	Human	ELISA, IF
<u>CPOX Antibody</u>	CSB-PA005910LA01HU	Human	ELISA, WB, IHC, IF
<u>CPPED1 Antibody</u>	CSB-PA871558LA01HU	Human	ELISA, IHC
<u>CPS1 Antibody</u>	CSB-PA005913LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>CPSF6 Antibody</u>	CSB-PA621966LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CPT1A Antibody</u>	CSB-PA005922LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CPT1B Antibody</u>	CSB-PA821674LA01HU	Human	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>CPT1C Antibody</u>	CSB-PA855048LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CPT2 Antibody</u>	CSB-PA005925LA01HU	Human	ELISA, IHC
<u>CPTP Antibody</u>	CSB-PA719416LA01HU	Human	ELISA, IHC, IF
<u>CPVL Antibody</u>	CSB-PA875664LA01HU	Human	ELISA, IHC
<u>CPXCR1 Antibody</u>	CSB-PA822678LA01HU	Human, Mouse	ELISA, WB, IHC
<u>Cpz Antibody</u>	CSB-PA005931LA01RA	Rat	ELISA, WB
<u>CR1 Antibody</u>	CSB-PA005932EA01HU	Human	ELISA, IHC
<u>Cr2 Antibody</u>	CSB-PA005934LA01MO	Mouse	ELISA
<u>CR2 Antibody</u>	CSB-PA005934LA01HU	Human	ELISA, IHC, IF
<u>CRABP1 Antibody</u>	CSB-PA005935DA01HU	Human, Mouse	ELISA, WB, IHC
<u>CRABP2 Antibody</u>	CSB-PA005936LA01HU	Human	ELISA, WB, IHC
<u>CRADD Antibody</u>	CSB-PA005938HA01HU	Human	ELISA, IF
<u>CRAT Antibody</u>	CSB-PA14169A0Rb	Human	ELISA, IHC
<u>CRB2 Antibody</u>	CSB-PA692212LA01HU	Human	ELISA, IHC
<u>CRCP Antibody</u>	CSB-PA005945LA01HU	Human	ELISA, IHC, IF
<u>CRD1 Antibody</u>	CSB-PA862988LA01DOA	Arabidopsis thaliana	ELISA, WB
<u>CREB1 Antibody</u>	CSB-PA005947HA01HU	Human, Mouse	ELISA, WB, IHC
<u>CREB3L1 Antibody</u>	CSB-PA856900LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>CREB3L4 Antibody</u>	CSB-PA844735LA01HU	Human	ELISA, IHC
<u>CREBBP Antibody</u>	CSB-PA838805LA01HU	Human	ELISA, IHC, IF
<u>CREM Antibody</u>	CSB-PA005961LA01HU	Human	ELISA, IHC, IF
<u>creN7 Antibody</u>	CSB-PA502491LA01FPB	Sulfolobus islandicus	ELISA
<u>crfb12 Antibody</u>	CSB-PA29691A0Rb	Zebrafish	ELISA
<u>CRH Antibody</u>	CSB-PA005963OA01HU	Human	ELISA, IHC
<u>CRHR1 Antibody</u>	CSB-PA005965EA01HU	Human	ELISA, WB, IHC, IF
<u>CRHR1 Antibody</u>	CSB-PA005965YA01HU	Human	ELISA, IHC, IF
<u>CRHR2 Antibody</u>	CSB-PA623797LA01HU	Human, Rat, Mouse	ELISA, WB, IHC, IF
<u>CRIP2 Antibody</u>	CSB-PA005970HA01HU	Human	ELISA
<u>CRIPT Antibody</u>	CSB-PA005973LA01HU	Human	ELISA, IHC, IF
<u>CRISPLD1 Antibody</u>	CSB-PA884457LA01HU	Human	ELISA, IHC
<u>CRISPLD2 Antibody</u>	CSB-PA884430LA01HU	Human	ELISA, WB
<u>CRK Antibody</u>	CSB-PA005979LA01HU	Human, Rat	ELISA, WB, IF
<u>CRKL Antibody</u>	CSB-PA005980LA01HU	Human	ELISA, WB, IF
<u>CRLF1 Antibody</u>	CSB-PA005982LA01HU	Human	ELISA, IHC
<u>CRMP1 Antibody</u>	CSB-PA623907HA01HU	Human, Mouse, Rat	ELISA, WB
<u>CRNN Antibody</u>	CSB-PA005988HA01HU	Human, Mouse, Rat	ELISA, WB
<u>CRP Antibody</u>	CSB-PA06079A0Rb	Human	ELISA, WB, IHC
<u>crp Antibody</u>	CSB-PA359785LA01ENV	Escherichia coli	ELISA



Product Name	Code	Species Reactivity	Applications
<u>CRTAC1 Antibody</u>	CSB-PA005997LA01HU	Human	ELISA, WB
<u>CRTC1 Antibody</u>	CSB-PA757785LA01HU	Human	ELISA, WB, IHC, IF
<u>CRTC2 Antibody</u>	CSB-PA706631LA01HU	Human	ELISA, IHC, IF
<u>CRTC3 Antibody</u>	CSB-PA747623LA01HU	Human	ELISA, IHC
<u>cry1Ab Antibody</u>	CSB-PA363358HA01BDC	Bacillus thuringiensis subsp. kurstaki	ELISA
<u>cry1Ab Antibody</u>	CSB-PA363358YA01BDC	Bacillus thuringiensis subsp. kurstaki	ELISA
<u>cry1Ac Antibody</u>	CSB-PA356448EA01BDC	Bacillus thuringiensis subsp. kurstaki	ELISA, WB
<u>cry1Ac Antibody</u>	CSB-PA356448YA01BDC	Bacillus thuringiensis subsp.	ELISA
<u>cry1Fb Antibody</u>	CSB-PA528987LA01BRS	Bacillus thuringiensis subsp. morrisoni	ELISA
<u>cry3Aa Antibody</u>	CSB-PA358337LA01BDG	Bacillus thuringiensis subsp. Tenebrionis	ELISA
<u>CRYAA Antibody</u>	CSB-PA006007LA01HU	Human	ELISA, IF
<u>CRYAB Antibody</u>	CSB-PA006008HA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>CRYBA1 Antibody</u>	CSB-PA006009LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CRYBA4 Antibody</u>	CSB-PA006011LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CRYGN Antibody</u>	CSB-PA841000LA01HU	Human	ELISA, IHC
<u>CRYM Antibody</u>	CSB-PA622791LA01HU	Human	ELISA, IHC, IF
<u>CRYZ Antibody</u>	CSB-PA600843LA01HU	Human	ELISA, IHC
<u>CRYZL1 Antibody</u>	CSB-PA006029LA01HU	Human	ELISA, WB
<u>CS Antibody</u>	CSB-PA006031LA01HU	Human, Zebrafish	ELISA, WB, IHC
<u>Csad Antibody</u>	CSB-PA006032LA01MO	Mouse	ELISA, WB
<u>CSAD Antibody</u>	CSB-PA006032LA01HU	Human	ELISA, IHC
<u>CSAG1 Antibody</u>	CSB-PA006033LA01HU	Human	ELISA, IHC
<u>CSDE1 Antibody</u>	CSB-PA006040LA01HU	Human	ELISA, IHC, IF
<u>CSE1L Antibody</u>	CSB-PA006042LA01HU	Human	ELISA, WB, IHC, IF
<u>CSF1 Antibody</u>	CSB-PA06859A0Rb	Human	ELISA, IHC, IF
<u>Csf1 Antibody</u>	CSB-PA07777A0Rb	Mouse	ELISA
<u>Csf1 Antibody</u>	CSB-PA07779A0Rb	Mouse	ELISA, IHC
<u>CSF1R Antibody</u>	CSB-PA006044LA01CA	Felis silvestris catus	ELISA
<u>CSF2 Antibody</u>	CSB-PA06475A0Rb	Human	ELISA, IHC
<u>CSF2RA Antibody</u>	CSB-PA006046LA01HU	Human	ELISA, IHC, IF
<u>CSF2RB Antibody</u>	CSB-PA006047LA01HU	Human	ELISA, IHC, IF
<u>CSF3 Antibody</u>	CSB-PA06459A0Rb	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>CSF3R Antibody</u>	CSB-PA860321LA01HU	Human	ELISA, IHC, IF
<u>CSGALNACT2 Antibody</u>	CSB-PA006052LA01HU	Human	ELISA, WB, IHC
<u>CSH1 Antibody</u>	CSB-PA11377A0Rb	Human	ELISA, WB
<u>CSH2 Antibody</u>	CSB-PA006054LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CSHL1 Antibody</u>	CSB-PA006055LA01HU	Human	ELISA, WB, IHC
<u>CSK Antibody</u>	CSB-PA006056LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CSMD2 Antibody</u>	CSB-PA770346LA01HU	Human	ELISA, IHC
<u>CSMD3 Antibody</u>	CSB-PA800102LA01HU	Human	ELISA, IHC, IF
<u>CSN2 Antibody</u>	CSB-PA006063LA01HU	Human	ELISA, IHC
<u>CSN3 Antibody</u>	CSB-PA006064LA01HU	Human	ELISA, IHC
<u>CSNK1A1 Antibody</u>	CSB-PA006065LA01HU	Human	ELISA, IHC
<u>CSNK1D Antibody</u>	CSB-PA006067LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF, IP
<u>CSNK1E Antibody</u>	CSB-PA01315A0Rb	Human	ELISA, WB, IHC
<u>CSNK1E Antibody</u>	CSB-PA01319A0Rb	Human	ELISA, WB, IHC, IF
<u>CSNK1G1 Antibody</u>	CSB-PA888030EA01HU	Human	ELISA, WB, IHC
<u>CSNK1G3 Antibody</u>	CSB-PA896549LA01HU	Human	ELISA, IF
<u>CSNK2A1 Antibody</u>	CSB-PA006072LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>CSNK2A2 Antibody</u>	CSB-PA526319HA01HU	Human	ELISA, IHC, IF
<u>CSNK2A3 Antibody</u>	CSB-PA086427LA01HU	Human	ELISA, WB, IHC, IP
<u>CSNK2B Antibody</u>	CSB-PA01655A0Rb	Human	ELISA, IHC
<u>csrA Antibody</u>	CSB-PA543019HA01ENT	Escherichia coli	ELISA, WB
<u>CSRNP3 Antibody</u>	CSB-PA855508NA01HU	Human	ELISA, WB
<u>CSRP2 Antibody</u>	CSB-PA06937A0Rb	Human	ELISA, IHC
<u>cssB Antibody</u>	CSB-PA347417HA01ENL	Escherichia coli	ELISA
<u>CST11 Antibody</u>	CSB-PA887950LA01HU	Human	ELISA, IHC
<u>CST3 Antibody</u>	CSB-PA06196A0Rb	Human	ELISA, WB, IHC
<u>CST3 Antibody</u>	CSB-PA06197A0Rb	Human	ELISA, WB, IHC
<u>CST5 Antibody</u>	CSB-PA006093LA01HU	Human	ELISA, IHC
<u>CST9L Antibody</u>	CSB-PA862015LA01HU	Human	ELISA, IHC
<u>CSTA Antibody</u>	CSB-PA006099DA01HU	Human	ELISA, WB, IHC, IF
<u>CSTB Antibody</u>	CSB-PA00224A0Rb	Human	ELISA, WB, IHC, IF, IP
<u>CTAGE1 Antibody</u>	CSB-PA836295LA01HU	Human	ELISA, IF
<u>CTAGE5 Antibody</u>	CSB-PA006131LA01HU	Human	ELISA, IHC
<u>Ctcf1 Antibody</u>	CSB-PA006139LA01MO	Mouse, Human	ELISA, WB
<u>CTDNEP1 Antibody</u>	CSB-PA05415A0Rb	Human	ELISA, IHC
<u>CTF1 Antibody</u>	CSB-PA621965LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>CTGF Antibody</u>	CSB-PA09429A0Rb	Human	ELISA, IHC, IF
<u>Ctgf Antibody</u>	CSB-PA006147LA01MO	Mouse, Rat	ELISA, WB
<u>CTH Antibody</u>	CSB-PA006160LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>CTHRC1 Antibody</u>	CSB-PA822189EA01HU	Human	ELISA, IHC
<u>CTHRC1 Antibody</u>	CSB-PA822189LA01HU	Human	ELISA, IHC
<u>CTLA4 Antibody</u>	CSB-PA006163LA01HU	Human	ELISA, IF
<u>CTNNA1 Antibody</u>	CSB-PA006164EA01HU	Human	ELISA, WB, IHC, IF
<u>CTNNA2 Antibody</u>	CSB-PA006165LA01HU	Human	ELISA, IF
<u>CTNNA3 Antibody</u>	CSB-PA883421LA01HU	Human	ELISA, IHC
<u>CTNNB1 Antibody</u>	CSB-PA00174A0Rb	Human	ELISA, WB, IHC, IF, IP
<u>CTNND1 Antibody</u>	CSB-PA006172LA01HU	Human	ELISA, IHC, IF
<u>CTNND2 Antibody</u>	CSB-PA892496LA01HU	Human	ELISA, IHC, IF
<u>CTPS1 Antibody</u>	CSB-PA006176LA01HU	Human	ELISA, WB, IHC
<u>CTPS2 Antibody</u>	CSB-PA865117EA01HU	Human	ELISA, WB, IHC, IP
<u>CTRB1 Antibody</u>	CSB-PA324640LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CTRL Antibody</u>	CSB-PA006182LA01HU	Human	ELISA, IHC, IF
<u>CTSA Antibody</u>	CSB-PA006184LA01HU	Human	ELISA, IHC, IF
<u>CTSB Antibody</u>	CSB-PA06974A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>CTSC Antibody</u>	CSB-PA006186LA01HU	Human	ELISA, IHC
<u>CTSD Antibody</u>	CSB-PA02129A0Rb	Human	ELISA, IHC, IP
<u>CTSE Antibody</u>	CSB-PA006188LA01HU	Human	ELISA, IHC
<u>CTSF Antibody</u>	CSB-PA15279A0Rb	Human	ELISA
<u>Ctsg Antibody</u>	CSB-PA006190LA01MO	Mouse	ELISA
<u>CTSH Antibody</u>	CSB-PA006191LA01HU	Human	ELISA, WB, IHC
<u>Ctsk Antibody</u>	CSB-PA006192LA01RA	Rat	ELISA, IHC
<u>CTSK Antibody</u>	CSB-PA15119A0Rb	Human	ELISA
<u>CTSL Antibody</u>	CSB-PA006193LA01DSU	Chlorocebus aethiops, Human	ELISA, WB
<u>CTSL Antibody</u>	CSB-PA11499A0Rb	Human	ELISA, IHC, IF
<u>Ctss Antibody</u>	CSB-PA006204LA01MO	Mouse	ELISA
<u>CTSS Antibody</u>	CSB-PA10729A0Rb	Human	ELISA, WB, IHC
<u>CTSZ Antibody</u>	CSB-PA006206LA01HU	Human	ELISA, IHC
<u>CUEDC2 Antibody</u>	CSB-PA884464HA01HU	Human	ELISA
<u>CUL2 Antibody</u>	CSB-PA623822LA01HU	Human	ELISA, WB, IHC
<u>CUL4A Antibody</u>	CSB-PA619782LA01HU	Human	ELISA, IF
<u>CUL7 Antibody</u>	CSB-PA617942LA01HU	Human	ELISA, IHC, IF
<u>CUP1-1 Antibody</u>	CSB-PA317846DA01SVG	Saccharomyces cerevisiae	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>CUTA Antibody</u>	CSB-PA006226LA01HU	Human	ELISA, IHC, IF
<u>CUTC Antibody</u>	CSB-PA873651LA01HU	Human	ELISA, IHC, IF
<u>CUX1 Antibody</u>	CSB-PA614451LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>CUX1 Antibody</u>	CSB-PA006228LA01HU	Human	ELISA, IF
<u>CWC22 Antibody</u>	CSB-PA875721LA01HU	Human	ELISA, IHC, IF
<u>CWF19L1 Antibody</u>	CSB-PA727965LA01HU	Human	ELISA, WB, IHC
<u>cwlM Antibody</u>	CSB-PA339780HA01BQT	Bacillus licheniformis	ELISA
<u>CX3CL1 Antibody</u>	CSB-PA06087A0Rb	Human	ELISA, IHC, IF
<u>CX3CR1 Antibody</u>	CSB-PA006236LA01HU	Human	ELISA, IHC, IF
<u>CXADR Antibody</u>	CSB-PA006237LA01HU	Human	ELISA, IHC
<u>CXCL1 Antibody</u>	CSB-PA06097A0Rb	Human	ELISA
<u>Cxcl1 Antibody</u>	CSB-PA09217A0Rb	Mouse	ELISA
<u>CXCL10 Antibody</u>	CSB-PA06107A0Rb	Human	ELISA
<u>CXCL10 Antibody</u>	CSB-PA822646LA01MOW	Macaca mulatta	ELISA
<u>Cxcl11 Antibody</u>	CSB-PA09289A0Rb	Mouse	ELISA, IHC
<u>CXCL11 Antibody</u>	CSB-PA06119A0Rb	Human	ELISA, IHC, IF
<u>CXCL12 Antibody</u>	CSB-PA343703LA01HU	Human	ELISA, IHC
<u>CXCL13 Antibody</u>	CSB-PA06139A0Rb	Human	ELISA, IHC
<u>Cxcl14 Antibody</u>	CSB-PA09319A0Rb	Mouse	ELISA
<u>Cxcl16 Antibody</u>	CSB-PA09337A0Rb	Mouse	ELISA
<u>CXCL2 Antibody</u>	CSB-PA08989A0Rb	Human	ELISA
<u>CXCL3 Antibody</u>	CSB-PA08999A0Rb	Human	ELISA, IHC
<u>Cxcl3 Antibody</u>	CSB-PA09239A0Rb	Mouse	ELISA
<u>CXCL5 Antibody</u>	CSB-PA06155A0Rb	Human	ELISA, WB
<u>Cxcl5 Antibody</u>	CSB-PA09257A0Rb	Mouse	ELISA, WB
<u>CXCL6 Antibody</u>	CSB-PA06169A0Rb	Human	ELISA, IHC
<u>CXCL8 Antibody</u>	CSB-PA011671YA01DO	Dog	ELISA
<u>CXCL8 Antibody</u>	CSB-PA08327A0Rb	Human	ELISA, WB, IHC
<u>CXCL8 Antibody</u>	CSB-PA15789A0Rb	Gallus gallus	ELISA
<u>CXCL8 Antibody</u>	CSB-PA08329A0Rb	Human	ELISA, IF
<u>Cxcl9 Antibody</u>	CSB-PA09269A0Rb	Mouse	ELISA
<u>CXCL9 Antibody</u>	CSB-PA06189A0Rb	Human	ELISA, IHC
<u>CXCR2 Antibody</u>	CSB-PA011673HA01HU	Human	ELISA, IHC, IF
<u>CXCR3 Antibody</u>	CSB-PA006253LA01HU	Human	ELISA, WB, IHC, IF
<u>CXCR4 Antibody</u>	CSB-PA006254YA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CXCR5 Antibody</u>	CSB-PA006255LA01HU	Human	ELISA, IHC, IF
<u>CXorf1 Antibody</u>	CSB-PA006258LA01HU	Human	ELISA, IHC
<u>CXorf21 Antibody</u>	CSB-PA867184LA01HU	Human	ELISA, IHC
<u>CXorf36 Antibody</u>	CSB-PA887998LA01HU	Human	ELISA, WB, IHC
<u>CXorf40A Antibody</u>	CSB-PA837446LA01HU	Human	ELISA



Product Name	Code	Species Reactivity	Applications
<u>CXorf40B Antibody</u>	CSB-PA842656LA01HU	Human	ELISA
<u>CXorf56 Antibody</u>	CSB-PA880981LA01HU	Human	ELISA, IHC
<u>CXorf58 Antibody</u>	CSB-PA856973LA01HU	Human	ELISA
<u>CXXC5 Antibody</u>	CSB-PA741975LA01HU	Human	ELISA, IHC
<u>CYB5A Antibody</u>	CSB-PA03134A0Rb	Human	ELISA, IF
<u>CYB5B Antibody</u>	CSB-PA006310LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CYB5R1 Antibody</u>	CSB-PA883402LA01HU	Human	ELISA, IHC
<u>CYB5R2 Antibody</u>	CSB-PA738178LA01HU	Human, Mouse	ELISA, WB, IHC
<u>CYB5R3 Antibody</u>	CSB-PA02655A0Rb	Human	ELISA, WB, IHC, IF
<u>CYBA Antibody</u>	CSB-PA006323LA01HU	Human	ELISA, IHC, IF
<u>CYBB Antibody</u>	CSB-PA006325LA01HU	Human	ELISA, IHC, IF
<u>CYBC1 Antibody</u>	CSB-PA859832LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CYBRD1 Antibody</u>	CSB-PA006326LA01HU	Human	ELISA, IF
<u>CYC1 Antibody</u>	CSB-PA01515A0Rb	Human	ELISA, WB, IHC, IF
<u>CYCB1-1 Antibody</u>	CSB-PA339075LA01DOA	Arabidopsis thaliana	ELISA, WB
<u>Cyclomaltodextrin glucanotransferase Antibody</u>	CSB-PA327084LA01EQJ	Bacillus macerans	ELISA, WB
<u>CYCS Antibody</u>	CSB-PA006328EA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>CYFIP1 Antibody</u>	CSB-PA748722LA01HU	Human	ELISA, IF
<u>CYGB Antibody</u>	CSB-PA006376LA01HU	Human	ELISA, IHC, IF
<u>CYLC2 Antibody</u>	CSB-PA618873LA01HU	Human	ELISA, WB
<u>CYLD Antibody</u>	CSB-PA865105LA01HU	Human	ELISA, IHC, IF
<u>CYP11A1 Antibody</u>	CSB-PA006389LA01HU	Human, Rat	ELISA, WB, IHC
<u>CYP11B1 Antibody</u>	CSB-PA006390LA01HU	Human, Mouse	ELISA, WB
<u>CYP11B2 Antibody</u>	CSB-PA006391LA01HU	Human	ELISA, IHC
<u>cyp125 Antibody</u>	CSB-PA351308HA01MVZ	Mycobacterium tuberculosis	ELISA
<u>CYP17A1 Antibody</u>	CSB-PA006392EA01HU	Human	ELISA, WB, IHC, IF
<u>CYP17A1 Antibody</u>	CSB-PA006392LA01HU	Human	ELISA, WB, IHC, IF
<u>cyp1a Antibody</u>	CSB-PA20839A0Rb	Zebrafish	ELISA
<u>cyp1a Antibody</u>	CSB-PA24439A0Rb	Zebrafish	ELISA, WB
<u>CYP1A1 Antibody</u>	CSB-PA006395LA01HU	Human	ELISA, IHC
<u>CYP1A2 Antibody</u>	CSB-PA006396LA01HU	Human	ELISA, IHC
<u>CYP1B1 Antibody</u>	CSB-PA618091LA01HU	Human	ELISA, WB
<u>CYP20A1 Antibody</u>	CSB-PA747625LA01HU	Human	ELISA, IHC, IF
<u>CYP21A2 Antibody</u>	CSB-PA006400LA01HU	Human	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>cyp26a1 Antibody</u>	CSB-PA006402LA01DIL	Zebrafish	ELISA
<u>cyp26b1 Antibody</u>	CSB-PA721935LA01DIL	Zebrafish	ELISA, WB
<u>CYP26C1 Antibody</u>	CSB-PA751026LA01HU	Human	ELISA, IF
<u>CYP27A1 Antibody</u>	CSB-PA880560LA01HU	Human	ELISA, IHC
<u>CYP2A6 Antibody</u>	CSB-PA006410LA01HU	Human, Rat	ELISA, WB
<u>CYP2B6 Antibody</u>	CSB-PA006414LA01HU	Human	ELISA, IHC, IF
<u>CYP2C19 Antibody</u>	CSB-PA006417LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CYP2E1 Antibody</u>	CSB-PA006425EA01HU	Human, Mouse	ELISA, WB, IHC
<u>CYP2U1 Antibody</u>	CSB-PA768766LA01HU	Human	ELISA, IF
<u>CYP39A1 Antibody</u>	CSB-PA882143LA01HU	Human	ELISA, WB, IHC
<u>Cyp3a13 Antibody</u>	CSB-PA731056LA01MO	Mouse	ELISA
<u>CYP3A4 Antibody</u>	CSB-PA006439LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>CYP3A7 Antibody</u>	CSB-PA006444LA01HU	Human	ELISA, IHC, IF
<u>CYP4F11 Antibody</u>	CSB-PA872539LA01HU	Human, Mouse	ELISA, WB, IF
<u>CYP4F2 Antibody</u>	CSB-PA006451LA01HU	Human	ELISA, IHC
<u>CYP4F22 Antibody</u>	CSB-PA740654LA01HU	Human	ELISA, WB
<u>CYP4X1 Antibody</u>	CSB-PA847606LA01HU	Human	ELISA, IHC, IF
<u>CYP51A1 Antibody</u>	CSB-PA614993LA01HU	Human	ELISA, WB, IHC
<u>CYP7A1 Antibody</u>	CSB-PA006460LA01HU	Human	ELISA, WB, IHC
<u>CYP7B1 Antibody</u>	CSB-PA006461LA01HU	Human	ELISA, IHC
<u>CYP8B1 Antibody</u>	CSB-PA006462LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>CYR61 Antibody</u>	CSB-PA10349A0Rb	Human	ELISA, IHC
<u>CYREN Antibody</u>	CSB-PA883613HA01HU	Human	ELISA, IHC
<u>CYS4 Antibody</u>	CSB-PA14487A0Rb	Saccharomyces cerevisiae	ELISA
<u>CYSLTR1 Antibody</u>	CSB-PA006465LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>cvsP Antibody</u>	CSB-PA325847HA01ENV	Escherichia coli	ELISA, WB
<u>CYTL1 Antibody</u>	CSB-PA873640LA01HU	Human	ELISA, IHC
<u>Cytotoxin 4 Antibody</u>	CSB-PA355645LA01NFG	Naja atra	ELISA
<u>d2hgdh Antibody</u>	CSB-PA374600LA01DIL	Zebrafish	ELISA
<u>D2HGDH Antibody</u>	CSB-PA843153LA01HU	Human, Mouse	ELISA, WB
<u>D8L Antibody</u>	CSB-PA322653LA01VAA	Vaccinia virus	ELISA, WB
<u>Da-3 Antibody</u>	CSB-PA22309A0Rb	Trifolium repens	ELISA
<u>DAAM2 Antibody</u>	CSB-PA774794LA01HU	Human	ELISA, IHC
<u>DAB2IP Antibody</u>	CSB-PA730138OA01HU	Human	ELISA, IF
<u>DACH1 Antibody</u>	CSB-PA006482LA01HU	Human	ELISA, IHC, IF
<u>DACT1 Antibody</u>	CSB-PA873691LA01HU	Human	ELISA, IHC, IF
<u>DAD1 Antibody</u>	CSB-PA03124A0Rb	Human	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>DAG1 Antibody</u>	CSB-PA613486NA01HU	Human	ELISA, IHC
<u>DALRD3 Antibody</u>	CSB-PA707703LA01HU	Human	ELISA
<u>DAND5 Antibody</u>	CSB-PA850835LA01HU	Human	ELISA, IHC
<u>DAO Antibody</u>	CSB-PA006494LA01HU	Human	ELISA, WB, IHC
<u>DAOA Antibody</u>	CSB-PA006495EA01HU	Human, Mouse	ELISA, WB
<u>DAP Antibody</u>	CSB-PA01034A0Rb	Human	ELISA, IHC
<u>DAPK2 Antibody</u>	CSB-PA006500EA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>DAPK3 Antibody</u>	CSB-PA006501LA01HU	Human	ELISA, WB, IF
<u>DAPL1 Antibody</u>	CSB-PA006502LA01HU	Human	ELISA, IHC
<u>DARS Antibody</u>	CSB-PA006505LA01HU	Human	ELISA
<u>DARS2 Antibody</u>	CSB-PA747490LA01HU	Human	ELISA, WB, IHC
<u>DAXX Antibody</u>	CSB-PA871395LA01HU	Human	ELISA, IHC
<u>DAZAP1 Antibody</u>	CSB-PA822205LA01HU	Human	ELISA, WB, IF
<u>DAZAP2 Antibody</u>	CSB-PA622986LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>DAZL Antibody</u>	CSB-PA006514LA01HU	Human	ELISA, IHC, IF
<u>DBF4 Antibody</u>	CSB-PA887010LA01HU	Human	ELISA, WB, IHC
<u>DBF4B Antibody</u>	CSB-PA822784LA01HU	Human	ELISA, WB, IHC
<u>DBI Antibody</u>	CSB-PA006519HA01HU	Human	ELISA, IHC
<u>DBNDD2 Antibody</u>	CSB-PA863102HA01HU	Human	ELISA, IHC
<u>DBR1 Antibody</u>	CSB-PA890914LA01HU	Human	ELISA, IHC
<u>DBT Antibody</u>	CSB-PA006530LA01HU	Human, Mouse	ELISA, WB
<u>DCAF1 Antibody</u>	CSB-PA896722LA01HU	Human	ELISA, IHC, IF
<u>DCAF4 Antibody</u>	CSB-PA855479LA01HU	Human	ELISA, WB
<u>DCAF8 Antibody</u>	CSB-PA725896HA01HU	Human, Mouse	ELISA, WB, IHC
<u>DCAKD Antibody</u>	CSB-PA855481LA01HU	Human	ELISA
<u>DCD Antibody</u>	CSB-PA006537DA01HU	Human	ELISA, IHC
<u>DCHS1 Antibody</u>	CSB-PA850320OA01HU	Human	ELISA, IHC, IF
<u>DCK Antibody</u>	CSB-PA006547LA01HU	Human	ELISA, IHC
<u>DCLK1 Antibody</u>	CSB-PA006548LA01HU	Human	ELISA, IHC, IF
<u>DCLRE1B Antibody</u>	CSB-PA880993LA01HU	Human	ELISA, IF
<u>DCN Antibody</u>	CSB-PA006554LA01BO	Bovine, Human	ELISA, WB
<u>Dcn Antibody</u>	CSB-PA006554LA01MO	Mouse	ELISA
<u>DCN Antibody</u>	CSB-PA10289A0Rb	Human, Mouse	ELISA, WB, IHC
<u>DCP2 Antibody</u>	CSB-PA810265LA01HU	Human	ELISA, IHC
<u>DCSTAMP Antibody</u>	CSB-PA023625LA01HU	Human	ELISA, IHC, IF
<u>DCT Antibody</u>	CSB-PA006562LA01HU	Human	ELISA, IHC
<u>DCTN2 Antibody</u>	CSB-PA614427LA01HU	Human	ELISA, IHC, IF
<u>DCTN3 Antibody</u>	CSB-PA006566LA01HU	Human	ELISA, IHC
<u>DCTN4 Antibody</u>	CSB-PA892445LA01HU	Human	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>DCTN5 Antibody</u>	CSB-PA871573LA01HU	Human	ELISA, IHC
<u>DCUN1D1 Antibody</u>	CSB-PA839329HA01HU	Human	ELISA, IHC
<u>DCUN1D2 Antibody</u>	CSB-PA754270LA01HU	Human	ELISA
<u>DCX Antibody</u>	CSB-PA006576DA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>DCX Antibody</u>	CSB-PA006576LA01HU	Human	ELISA, IHC, IF
<u>dcvD Antibody</u>	CSB-PA18017A0Rb	Escherichia coli	ELISA, WB
<u>DDA1 Antibody</u>	CSB-PA880116LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>DDAH2 Antibody</u>	CSB-PA006580LA01HU	Human, Mouse	ELISA, WB, IHC
<u>DDB1 Antibody</u>	CSB-PA618079LA01HU	Human	ELISA, IHC, IF
<u>DDB2 Antibody</u>	CSB-PA846067LA01HU	Human, Rat, Mouse	ELISA, WB
<u>DDC Antibody</u>	CSB-PA006583LA01HU	Human	ELISA, IHC, IF
<u>DDHD2 Antibody</u>	CSB-PA006586LA01HU	Human	ELISA, WB, IHC
<u>DDI2 Antibody</u>	CSB-PA732937LA01HU	Human	ELISA, WB, IHC
<u>DDIT3 Antibody</u>	CSB-PA006589LA01HU	Human	ELISA, IHC, IF
<u>DDIT4L Antibody</u>	CSB-PA853399LA01HU	Human	ELISA, IHC, IF
<u>DDO Antibody</u>	CSB-PA006593LA01HU	Human, Rat	ELISA, WB, IHC
<u>DDR2 Antibody</u>	CSB-PA620985LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>Ddt Antibody</u>	CSB-PA006598HA01MO	Mouse	ELISA
<u>DDX17 Antibody</u>	CSB-PA846106LA01HU	Human	ELISA, IF
<u>DDX18 Antibody</u>	CSB-PA882116LA01HU	Human	ELISA, IF
<u>DDX19A Antibody</u>	CSB-PA889116HA01HU	Human, Mouse	ELISA, WB, IHC
<u>DDX27 Antibody</u>	CSB-PA836220HA01HU	Human	ELISA
<u>DDX28 Antibody</u>	CSB-PA873656LA01HU	Human	ELISA, WB
<u>DDX31 Antibody</u>	CSB-PA872517LA01HU	Human	ELISA, IHC
<u>DDX39B Antibody</u>	CSB-PA00304A0Rb	Human	ELISA, WB, IHC, IF
<u>DDX3X Antibody</u>	CSB-PA006621LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF, IP
<u>DDX4 Antibody</u>	CSB-PA882074LA01HU	Human	ELISA, IHC, IF
<u>DDX42 Antibody</u>	CSB-PA773053LA01HU	Human	ELISA, IHC, IF
<u>DDX47 Antibody</u>	CSB-PA867123LA01HU	Human	ELISA, WB
<u>DDX50 Antibody</u>	CSB-PA861080LA01HU	Human	ELISA, WB, IHC
<u>DDX55 Antibody</u>	CSB-PA13569A0Rb	Human	ELISA, IHC
<u>DDX58 Antibody</u>	CSB-PA006638LA01HU	Human	ELISA
<u>DDX59 Antibody</u>	CSB-PA722563LA01HU	Human	ELISA, IHC
<u>DDX6 Antibody</u>	CSB-PA006640LA01HU	Human	ELISA, WB, IHC, IF
<u>DEAF1 Antibody</u>	CSB-PA006644LA01HU	Human	ELISA, WB
<u>def Antibody</u>	CSB-PA017707HA01EOD	Escherichia coli	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>DEF6 Antibody</u>	CSB-PA884467LA01HU	Human, Mouse	ELISA, WB, IF
<u>DEF8 Antibody</u>	CSB-PA765060LA01HU	Human	ELISA
<u>Defa5 Antibody</u>	CSB-PA006657EA01MO	Mouse	ELISA
<u>DEFB1 Antibody</u>	CSB-PA09854A0Rb	Human	ELISA, IHC, IF
<u>DEFB104A Antibody</u>	CSB-PA840974LA01HU	Human	ELISA, IHC
<u>DEFB126 Antibody</u>	CSB-PA006696EA01EXP	Pongo pygmaeus	ELISA
<u>DEFB126 Antibody</u>	CSB-PA866323LA01HU	Human	ELISA, IHC
<u>DEFB129 Antibody</u>	CSB-PA863952LA01HU	Human	ELISA, WB, IHC
<u>Defb14 Antibody</u>	CSB-PA749254HA01MO	Mouse, Rat	ELISA, WB
<u>Defensin D2 Antibody</u>	CSB-PA305599LA01FKI	Spinacia oleracea	ELISA
<u>Defensin-like protein 1 Antibody</u>	CSB-PA316614LA01DAE	Dahlia merckii	ELISA
<u>DEG15 Antibody</u>	CSB-PA819818LA01DOA	Arabidopsis thaliana	ELISA, WB
<u>DEGP12 Antibody</u>	CSB-PA873288LA01DOA	Arabidopsis thaliana	ELISA, WB
<u>Delta-ctenitoxin-Pn1a Antibody</u>	CSB-PA349440HA01EUV	Phoneutria nigriventer	ELISA
<u>DENND1A Antibody</u>	CSB-PA823472LA01HU	Human	ELISA, IHC
<u>DENND1C Antibody</u>	CSB-PA810272LA01HU	Human, Rat	ELISA, WB, IHC
<u>DENND2C Antibody</u>	CSB-PA724250LA01HU	Human	ELISA, WB
<u>DENND2D Antibody</u>	CSB-PA887987LA01HU	Human	ELISA, WB, IHC
<u>DENND5A Antibody</u>	CSB-PA757108LA01HU	Human	ELISA, IHC
<u>DENR Antibody</u>	CSB-PA006724HA01HU	Human	ELISA
<u>deoD Antibody</u>	CSB-PA18357A0Rb	Escherichia coli	ELISA, WB
<u>DEPDC1B Antibody</u>	CSB-PA851968LA01HU	Human	ELISA, WB, IHC, IF
<u>DEPP1 Antibody</u>	CSB-PA865135LA01HU	Human	ELISA, IHC, IF
<u>DES Antibody</u>	CSB-PA13149A0Rb	Human	ELISA, IHC, IF
<u>DET1 Antibody</u>	CSB-PA762572LA01HU	Human	ELISA, IHC
<u>DEUP1 Antibody</u>	CSB-PA004729LA01HU	Human	ELISA, IHC
<u>DEXI Antibody</u>	CSB-PA529744LA01HU	Human	ELISA, IHC
<u>DFFB Antibody</u>	CSB-PA006738LA01HU	Human	ELISA, IHC
<u>DGAT1 Antibody</u>	CSB-PA006817LA01HU	Human	ELISA, IHC
<u>DGAT2 Antibody</u>	CSB-PA850407LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>dgcQ Antibody</u>	CSB-PA08524A0Rb	Escherichia coli	ELISA, WB
<u>DGCR6L Antibody</u>	CSB-PA880133HA01HU	Human, Mouse	ELISA, WB
<u>DGKG Antibody</u>	CSB-PA006836LA01HU	Human, Mouse	ELISA, WB
<u>DGKZ Antibody</u>	CSB-PA618788LA01HU	Human	ELISA, IF
<u>DGUOK Antibody</u>	CSB-PA621974HA01HU	Human	ELISA, IHC
<u>dharma Antibody</u>	CSB-PA20959A0Rb	Zebrafish	ELISA
<u>DHCR7 Antibody</u>	CSB-PA006844LA01HU	Human	ELISA, WB, IHC
<u>DHFR Antibody</u>	CSB-PA006847HA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>DHH Antibody</u>	CSB-PA006851LA01HU	Human	ELISA, IF
<u>DHODH Antibody</u>	CSB-PA006852LA01HU	Human	ELISA, IHC
<u>DHPS Antibody</u>	CSB-PA006854LA01HU	Human	ELISA, IHC
<u>DHRS12 Antibody</u>	CSB-PA006857LA01HU	Human	ELISA, WB, IHC, IF
<u>DHRS13 Antibody</u>	CSB-PA761493LA01HU	Human	ELISA, WB, IHC
<u>DHRS3 Antibody</u>	CSB-PA006860LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>DHRS7 Antibody</u>	CSB-PA897480LA01HU	Human	ELISA, WB, IHC
<u>DHRS7B Antibody</u>	CSB-PA760785LA01HU	Human	ELISA
<u>DHRS9 Antibody</u>	CSB-PA006867DA01HU	Human	ELISA, IHC, IF
<u>DHRSX Antibody</u>	CSB-PA836670LA01HU	Human	ELISA, IHC, IF
<u>DHX15 Antibody</u>	CSB-PA006870LA01HU	Human	ELISA, WB
<u>DHX30 Antibody</u>	CSB-PA752098LA01HU	Human	ELISA, IHC
<u>DHX32 Antibody</u>	CSB-PA748724LA01HU	Human, Mouse	ELISA, WB, IHC
<u>DHX33 Antibody</u>	CSB-PA884479LA01HU	Human	ELISA, WB, IHC, IF
<u>DHX34 Antibody</u>	CSB-PA621761LA01HU	Human	ELISA, IHC, IF
<u>DHX35 Antibody</u>	CSB-PA887985LA01HU	Human, Mouse	ELISA, WB, IHC
<u>DHX36 Antibody</u>	CSB-PA862002LA01HU	Human	ELISA, IF
<u>DHX38 Antibody</u>	CSB-PA835687LA01HU	Human	ELISA, IHC
<u>DHX9 Antibody</u>	CSB-PA600078LA01HU	Human	ELISA, IHC
<u>DIABLO Antibody</u>	CSB-PA865113LA01HU	Human	ELISA, IHC, IF
<u>DIAPH3 Antibody</u>	CSB-PA885723LA01HU	Human	ELISA, IF
<u>DICER1 Antibody</u>	CSB-PA891575LA01HU	Human	ELISA, IHC
<u>DIMT1 Antibody</u>	CSB-PA892161LA01HU	Human, Mouse	ELISA, WB, IHC
<u>DIO1 Antibody</u>	CSB-PA343445LA01HU	Human	ELISA
<u>DIO2 Antibody</u>	CSB-PA821705LA01HU	Human	ELISA, IHC, IF
<u>Dio3 Antibody</u>	CSB-PA344617LA01RA	Rat	ELISA
<u>DIO3 Antibody</u>	CSB-PA344617LA01HU	Human	ELISA, IF
<u>DIP2A Antibody</u>	CSB-PA623930LA01HU	Human	ELISA, IHC
<u>DIRAS1 Antibody</u>	CSB-PA006905LA01HU	Human	ELISA, IHC, IF
<u>DIRC1 Antibody</u>	CSB-PA846572LA01HU	Human	ELISA, IHC
<u>DIS3 Antibody</u>	CSB-PA897532LA01HU	Human	ELISA, IHC
<u>DIS3L Antibody</u>	CSB-PA823476LA01HU	Human	ELISA, WB, IHC
<u>DIS3L2 Antibody</u>	CSB-PA818258LA01HU	Human	ELISA, IHC
<u>DISC1 Antibody</u>	CSB-PA889093LA01HU	Human	ELISA, WB, IF
<u>Disintegrin EC3B Antibody</u>	CSB-PA305523HA01EAH	Echis carinatus	ELISA
<u>DIXDC1 Antibody</u>	CSB-PA621878LA01HU	Human	ELISA, IHC
<u>DKC1 Antibody</u>	CSB-PA006919LA01HU	Human, Mouse, Rat	ELISA, WB, IF
<u>DKK1 Antibody</u>	CSB-PA06219A0Rb	Human	ELISA, IHC, IF
<u>DKK2 Antibody</u>	CSB-PA08279A0Rb	Human	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>DKK3 Antibody</u>	CSB-PA890660LA01HU	Human	ELISA
<u>DKKL1 Antibody</u>	CSB-PA892450LA01HU	Human	ELISA, IHC
<u>DLAT Antibody</u>	CSB-PA006926LA01HU	Human	ELISA, WB, IP
<u>DLD Antibody</u>	CSB-PA006928LA01HU	Human, Rat, Mouse	ELISA, WB, IHC
<u>DLEU1 Antibody</u>	CSB-PA006930LA01HU	Human	ELISA, IHC, IF
<u>DLG1 Antibody</u>	CSB-PA614264LA01HU	Human	ELISA, IHC, IF
<u>DLG2 Antibody</u>	CSB-PA621882LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>DLG4 Antibody</u>	CSB-PA006938LA01HU	Human	ELISA, IHC
<u>DLG5 Antibody</u>	CSB-PA823464LA01HU	Human	ELISA, IF
<u>DLGAP2 Antibody</u>	CSB-PA882169LA01HU	Human	ELISA, IF
<u>DLK1 Antibody</u>	CSB-PA006945LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>Dlk2 Antibody</u>	CSB-PA812967LA01MO	Mouse	ELISA, WB
<u>DLL3 Antibody</u>	CSB-PA882142LA01HU	Human	ELISA, IHC
<u>DLL4 Antibody</u>	CSB-PA006949LA01HU	Human	ELISA, WB, IF
<u>DLST Antibody</u>	CSB-PA861526LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>DLX1 Antibody</u>	CSB-PA006952LA01HU	Human	ELISA, IF
<u>dlx2a Antibody</u>	CSB-PA342691LA01DIL	Zebrafish	ELISA
<u>DLX5 Antibody</u>	CSB-PA006956LA01HU	Human	ELISA, IHC, IF
<u>DLX6 Antibody</u>	CSB-PA006957LA01HU	Human	ELISA, WB
<u>DMAPI Antibody</u>	CSB-PA873611LA01HU	Human	ELISA, WB, IHC, IF
<u>DMC1 Antibody</u>	CSB-PA617920LA01HU	Human	ELISA, WB, IHC, IF
<u>DMD Antibody</u>	CSB-PA006963LA01HU	Human	ELISA, IHC
<u>DME Antibody</u>	CSB-PA839701LA01DOA	Arabidopsis thaliana	ELISA
<u>DMGDH Antibody</u>	CSB-PA883415LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>DMKN Antibody</u>	CSB-PA006965LA01HU	Human, Mouse	ELISA, WB
<u>DMPK Antibody</u>	CSB-PA006968LA01HU	Human	ELISA, IHC, IF
<u>DMRT1 Antibody</u>	CSB-PA896914LA01HU	Human, Mouse	ELISA, WB, IHC
<u>DMRT3 Antibody</u>	CSB-PA885692LA01HU	Human	ELISA, IHC
<u>DMRTA1 Antibody</u>	CSB-PA730143LA01HU	Human	ELISA, IHC
<u>DMXL2 Antibody</u>	CSB-PA819469LA01HU	Human	ELISA, IHC, IF
<u>DNA2 Antibody</u>	CSB-PA006982LA01HU	Human	ELISA, IHC
<u>DNAAF3 Antibody</u>	CSB-PA850845LA01HU	Human	ELISA, IHC, IF
<u>DNAAF4 Antibody</u>	CSB-PA855504LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>DNAAF5 Antibody</u>	CSB-PA771471LA01HU	Human	ELISA
<u>DNAH1 Antibody</u>	CSB-PA878961LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>DNAH12 Antibody</u>	CSB-PA747792LA01HU	Human	ELISA, WB, IF
<u>DNAH3 Antibody</u>	CSB-PA823461LA01HU	Human	ELISA, IHC
<u>DNAH5 Antibody</u>	CSB-PA819476LA01HU	Human	ELISA, IHC, IF
<u>DNAH7 Antibody</u>	CSB-PA841220LA01HU	Human	ELISA, IF
<u>DNAJA1 Antibody</u>	CSB-PA006999LA01HU	Human	ELISA, WB, IHC, IF
<u>DNAJA3 Antibody</u>	CSB-PA836207LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>DNAJB1 Antibody</u>	CSB-PA03934A0Rb	Human, Mouse	ELISA, WB, IHC
<u>DNAJB11 Antibody</u>	CSB-PA890663LA01HU	Human	ELISA, IHC, IF
<u>DNAJB14 Antibody</u>	CSB-PA819454LA01HU	Human	ELISA, IHC, IF
<u>DNAJB2 Antibody</u>	CSB-PA007008HA01HU	Human	ELISA, WB, IHC
<u>DNAJB4 Antibody</u>	CSB-PA10929A0Rb	Human, Mouse	ELISA, WB, IHC
<u>DNAJB5 Antibody</u>	CSB-PA007011LA01HU	Human	ELISA, WB, IHC
<u>DNAJB6 Antibody</u>	CSB-PA03239A0Rb	Human, Mouse	ELISA, WB, IF
<u>DNAJB8 Antibody</u>	CSB-PA836779LA01HU	Human, Mouse	ELISA, WB, IHC
<u>DNAJC10 Antibody</u>	CSB-PA809006LA01HU	Human, Rat	ELISA, WB, IHC
<u>DNAJC15 Antibody</u>	CSB-PA22849A0Rb	Human	ELISA, IHC
<u>DNAJC17 Antibody</u>	CSB-PA873664LA01HU	Human	ELISA, WB, IHC
<u>DNAJC3 Antibody</u>	CSB-PA615670LA01HU	Human	ELISA, IHC, IF
<u>DNAJC30 Antibody</u>	CSB-PA846659HA01HU	Human	ELISA, WB, IHC
<u>DNAJC4 Antibody</u>	CSB-PA865087LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>DNAJC6 Antibody</u>	CSB-PA007040LA01HU	Human	ELISA, IHC, IF
<u>dnaK Antibody</u>	CSB-PA633459HA01EGW	Escherichia coli	ELISA, WB
<u>DNAL4 Antibody</u>	CSB-PA007046LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>DNALI1 Antibody</u>	CSB-PA007047LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>DNASE1 Antibody</u>	CSB-PA007049LA01HU	Human, Mouse	ELISA, WB, IHC
<u>DNASE1L3 Antibody</u>	CSB-PA621686EA01HU	Human	ELISA, IHC
<u>DNASE2 Antibody</u>	CSB-PA007053LA01HU	Human	ELISA, WB, IHC
<u>DNASE2B Antibody</u>	CSB-PA841227LA01HU	Human	ELISA, WB
<u>DND1 Antibody</u>	CSB-PA809023LA01HU	Human	ELISA, IF
<u>DNER Antibody</u>	CSB-PA818784NA01HU	Human	ELISA, IF
<u>DNM1 Antibody</u>	CSB-PA007062LA01HU	Human	ELISA, IHC, IF
<u>DNM1L Antibody</u>	CSB-PA007078EA01HU	Human, Mouse	ELISA, WB, IHC
<u>DNMT1 Antibody</u>	CSB-PA007082HA01HU	Human	ELISA, IHC, IF
<u>DNMT3A Antibody</u>	CSB-PA007083LA01HU	Human	ELISA, IHC, IF
<u>DNMT3B Antibody</u>	CSB-PA883358HA01HU	Human	ELISA, WB, IHC, IF
<u>DNMT3L Antibody</u>	CSB-PA891717LA01HU	Human	ELISA, WB, IF



Product Name	Code	Species Reactivity	Applications
<u>DNPEP Antibody</u>	CSB-PA007087LA01HU	Human	ELISA, IHC, IF
<u>DNTT Antibody</u>	CSB-PA007088LA01HU	Human	ELISA, IHC, IF
<u>DNTTIP1 Antibody</u>	CSB-PA007089LA01HU	Human	ELISA, IHC
<u>DOC2B Antibody</u>	CSB-PA618877LA01HU	Human	ELISA, IHC, IF
<u>DOCK1 Antibody</u>	CSB-PA619864LA01HU	Human	ELISA, IHC, IF
<u>DOCK10 Antibody</u>	CSB-PA842645LA01HU	Human	ELISA, IHC
<u>DOCK3 Antibody</u>	CSB-PA818270LA01HU	Human	ELISA, IHC, IF
<u>DOCK5 Antibody</u>	CSB-PA875689LA01HU	Human	ELISA
<u>DOCK6 Antibody</u>	CSB-PA846633LA01HU	Human	ELISA, IHC
<u>DOCK7 Antibody</u>	CSB-PA836266LA01HU	Human	ELISA, WB, IF
<u>DODA Antibody</u>	CSB-PA767105LA01PIAD	Portulaca grandiflora	ELISA
<u>DOHH Antibody</u>	CSB-PA880104LA01HU	Human	ELISA, IF
<u>DOK1 Antibody</u>	CSB-PA858723HA01HU	Human	ELISA, IHC, IF
<u>DOK2 Antibody</u>	CSB-PA007108LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>DOK3 Antibody</u>	CSB-PA007019LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>DOK5 Antibody</u>	CSB-PA007111LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>DOK6 Antibody</u>	CSB-PA007112LA01HU	Human	ELISA
<u>DOK7 Antibody</u>	CSB-PA623281LA01HU	Human	ELISA, WB, IF
<u>DOT1L Antibody</u>	CSB-PA848418LA01HU	Human	ELISA, WB, IHC
<u>DPEP1 Antibody</u>	CSB-PA007123LA01HU	Human	ELISA, WB, IF
<u>DPH1 Antibody</u>	CSB-PA871625LA01HU	Human	ELISA, IHC
<u>DPH2 Antibody</u>	CSB-PA859031LA01HU	Human, Mouse	ELISA, WB, IHC
<u>DPH3 Antibody</u>	CSB-PA822214LA01HU	Human	ELISA
<u>DPH5 Antibody</u>	CSB-PA862001LA01HU	Human	ELISA, IHC, IF
<u>DPH6 Antibody</u>	CSB-PA752102LA01HU	Human	ELISA, IHC, IF
<u>DPH7 Antibody</u>	CSB-PA866281LA01HU	Human	ELISA, IHC, IF
<u>DPM1 Antibody</u>	CSB-PA007134LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>DPP4 Antibody</u>	CSB-PA06229A0Rb	Human, Rat	ELISA, WB, IHC, IF
<u>DPP6 Antibody</u>	CSB-PA007140YA01HU	Human	ELISA
<u>DPP9 Antibody</u>	CSB-PA769770LA01HU	Human	ELISA, IHC
<u>DPPA3 Antibody</u>	CSB-PA744290LA01HU	Human	ELISA, IHC
<u>DPPA4 Antibody</u>	CSB-PA765975EA01HU	Human	ELISA, IHC
<u>DPPA5 Antibody</u>	CSB-PA007147LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>dps Antibody</u>	CSB-PA314354HA01MVW	Mycobacterium smegmatis	ELISA, WB
<u>DPY19L1 Antibody</u>	CSB-PA23079A0Rb	Human	ELISA, IF
<u>DPY19L2 Antibody</u>	CSB-PA757514LA01HU	Human	ELISA, IF



Product Name	Code	Species Reactivity	Applications
<u>DPY19L3 Antibody</u>	CSB-PA754545LA01HU	Human	ELISA, WB, IF
<u>DPY30 Antibody</u>	CSB-PA861193LA01HU	Human	ELISA, IHC
<u>DPYD Antibody</u>	CSB-PA622516LA01HU	Human	ELISA, IF
<u>DPYSL2 Antibody</u>	CSB-PA614975YA01HU	Human	ELISA, WB, IHC, IF
<u>DPYSL4 Antibody</u>	CSB-PA007172LA01HU	Human	ELISA, IHC, IF
<u>DRAM1 Antibody</u>	CSB-PA836675LA01HU	Human	ELISA, IHC
<u>DRAP1 Antibody</u>	CSB-PA05144A0Rb	Human	ELISA
<u>DRC7 Antibody</u>	CSB-PA809014LA01HU	Human	ELISA, IHC
<u>DRD1 Antibody</u>	CSB-PA007178LA01HU	Human	ELISA, WB, IHC, IF
<u>DRD5 Antibody</u>	CSB-PA007182EA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>DRG1 Antibody</u>	CSB-PA897524LA01HU	Human	ELISA, WB, IHC, IF
<u>DRG2 Antibody</u>	CSB-PA007186LA01HU	Human	ELISA
<u>DROSHA Antibody</u>	CSB-PA889096LA01HU	Human	ELISA, IHC
<u>DRP2 Antibody</u>	CSB-PA623806LA01HU	Human	ELISA, IF
<u>dsbA Antibody</u>	CSB-PA528735HA01FGP	Pseudomonas syringae pv. tomato	ELISA
<u>dsbA Antibody</u>	CSB-PA528735YA01FGP	Pseudomonas syringae pv. tomato	ELISA
<u>DSC3 Antibody</u>	CSB-PA613507LA01HU	Human	ELISA, IHC, IF
<u>DSCAM Antibody</u>	CSB-PA007192LA01HU	Human	ELISA, IHC
<u>DSCC1 Antibody</u>	CSB-PA887133LA01HU	Human	ELISA, WB
<u>DSCR4 Antibody</u>	CSB-PA007197LA01HU	Human	ELISA
<u>DSCR9 Antibody</u>	CSB-PA007200LA01HU	Human	ELISA, IHC
<u>dsdA Antibody</u>	CSB-PA472744HA01SWQ	Salmonella heidelberg	ELISA
<u>dsdA Antibody</u>	CSB-PA472744YA01SWQ	Salmonella heidelberg	ELISA
<u>DSE Antibody</u>	CSB-PA890925LA01HU	Human	ELISA, IHC, IF
<u>DSG1 Antibody</u>	CSB-PA861938LA01HU	Human	ELISA, IHC, IF
<u>DSG2 Antibody</u>	CSB-PA622752LA01HU	Human	ELISA, IHC, IF
<u>DSG3 Antibody</u>	CSB-PA007205LA01HU	Human	ELISA, IHC, IF
<u>Dsg3 Antibody</u>	CSB-PA007205LA01MO	Mouse	ELISA
<u>DSG3 Antibody</u>	CSB-PA007205YA01HU	Human, Mouse	ELISA, WB
<u>DSP Antibody</u>	CSB-PA007208LA01HU	Human	ELISA
<u>DST Antibody</u>	CSB-PA17019A0Rb	Human	ELISA, WB, IF
<u>DSTN Antibody</u>	CSB-PA007211LA01HU	Human	ELISA, IHC
<u>DTL Antibody</u>	CSB-PA889160LA01HU	Human	ELISA, IHC
<u>DTNBP1 Antibody</u>	CSB-PA822207LA01HU	Human	ELISA, IHC, IF
<u>DTWD1 Antibody</u>	CSB-PA007219EA01HU	Human	ELISA, IHC, IF
<u>DTX1 Antibody</u>	CSB-PA773055LA01HU	Human	ELISA, IF



Product Name	Code	Species Reactivity	Applications
<u>DTX3 Antibody</u>	CSB-PA839833LA01HU	Human	ELISA, IHC, IF
<u>DTX4 Antibody</u>	CSB-PA764784LA01HU	Human	ELISA, IF
<u>DUOXA1 Antibody</u>	CSB-PA04949A0Rb	Human	ELISA, IF
<u>DUS3L Antibody</u>	CSB-PA842675LA01HU	Human	ELISA, IHC
<u>DUS4L Antibody</u>	CSB-PA007237LA01HU	Human	ELISA, IHC
<u>DUSP1 Antibody</u>	CSB-PA007238LA01HU	Human	ELISA, WB
<u>DUSP11 Antibody</u>	CSB-PA007240LA01HU	Human	ELISA, WB
<u>DUSP13 Antibody</u>	CSB-PA734923HA01HU	Human	ELISA, WB, IHC, IF
<u>DUSP13 Antibody</u>	CSB-PA883436LA01HU	Human	ELISA, IHC
<u>DUSP15 Antibody</u>	CSB-PA861994LA01HU	Human	ELISA, IHC, IF
<u>DUSP19 Antibody</u>	CSB-PA855469LA01HU	Human	ELISA, IHC, IF
<u>DUSP21 Antibody</u>	CSB-PA875677LA01HU	Human	ELISA, IHC
<u>DUSP26 Antibody</u>	CSB-PA863128LA01HU	Human	ELISA, IHC, IF
<u>DUSP3 Antibody</u>	CSB-PA007255LA01HU	Human	ELISA, IHC
<u>DUSP4 Antibody</u>	CSB-PA619752LA01HU	Human, Mouse	ELISA, WB, IHC
<u>DUSP5 Antibody</u>	CSB-PA620976LA01HU	Human	ELISA, IF
<u>DUSP6 Antibody</u>	CSB-PA621973LA01HU	Human	ELISA, IHC, IF
<u>DUSP9 Antibody</u>	CSB-PA859528LA01HU	Human	ELISA, IHC, IF
<u>DUX4 Antibody</u>	CSB-PA890667LA01HU	Human	ELISA, IHC
<u>DVL1P1 Antibody</u>	CSB-PA007285OA01HU	Human	ELISA, WB
<u>DVL3 Antibody</u>	CSB-PA856442LA01HU	Human	ELISA, IHC, IF
<u>DXO Antibody</u>	CSB-PA007116LA01HU	Human	ELISA, IHC, IF
<u>DYDC1 Antibody</u>	CSB-PA819901LA01HU	Human	ELISA, WB
<u>dye Antibody</u>	CSB-PA772049LA01DIL	Zebrafish	ELISA, WB
<u>DYM Antibody</u>	CSB-PA007291LA01HU	Human, Mouse	ELISA, WB, IHC
<u>DYNC1H1 Antibody</u>	CSB-PA613491LA01HU	Human	ELISA, IHC, IF
<u>DYNC1LI1 Antibody</u>	CSB-PA897312LA01HU	Human	ELISA, IHC
<u>DYNC2H1 Antibody</u>	CSB-PA818763LA01HU	Human	ELISA, WB, IHC, IF
<u>DYNLL1 Antibody</u>	CSB-PA01284A0Rb	Human	ELISA, WB, IHC, IF
<u>DYNLL2 Antibody</u>	CSB-PA839322LA01HU	Human	ELISA
<u>DYNLRB1 Antibody</u>	CSB-PA865093HA01HU	Human	ELISA, IHC
<u>DYRK3 Antibody</u>	CSB-PA007312LA01HU	Human	ELISA, IHC
<u>DYRK4 Antibody</u>	CSB-PA882080LA01HU	Human	ELISA, IHC, IF
<u>DYSF Antibody</u>	CSB-PA007314LA01HU	Human	ELISA, IHC, IF
<u>E2 Antibody</u>	CSB-PA361949LA01HMN	Human type 18 papillomavirus	ELISA
<u>E2F1 Antibody</u>	CSB-PA007340LA01HU	Human, Mouse	ELISA, WB, IHC, IP
<u>E2F2 Antibody</u>	CSB-PA623909LA01HU	Human	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>E2F3 Antibody</u>	CSB-PA007342LA01HU	Human	ELISA, IF
<u>E2F4 Antibody</u>	CSB-PA00759A0Rb	Human	ELISA, WB, IHC
<u>E2F5 Antibody</u>	CSB-PA624021LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>E6 Antibody</u>	CSB-PA334186LA01HNY	Human papillomavirus type 52	ELISA
<u>E6 Antibody</u>	CSB-PA365854LA01HML	Human papillomavirus type 16	ELISA
<u>E7 Antibody</u>	CSB-PA361948LA01HMN	Human papillomavirus type 18	ELISA
<u>EAF2 Antibody</u>	CSB-PA856910LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>EAPP Antibody</u>	CSB-PA007353LA01HU	Human	ELISA, WB, IHC
<u>EARS2 Antibody</u>	CSB-PA705509LA01HU	Human	ELISA, IHC
<u>EBF2 Antibody</u>	CSB-PA884503LA01HU	Human	ELISA, IHC
<u>EBF3 Antibody</u>	CSB-PA880977LA01HU	Human	ELISA, WB
<u>EBI3 Antibody</u>	CSB-PA007360LA01HU	Human	ELISA, WB, IHC
<u>EBP Antibody</u>	CSB-PA624016EA01HU	Human	ELISA, WB, IHC, IF
<u>ECD Antibody</u>	CSB-PA007370LA01HU	Human	ELISA, WB
<u>ECEL1 Antibody</u>	CSB-PA007373LA01HU	Human	ELISA, IHC, IF
<u>ECH1 Antibody</u>	CSB-PA621643LA01HU	Human	ELISA, WB, IHC, IF
<u>ECHDC1 Antibody</u>	CSB-PA882103LA01HU	Human	ELISA
<u>ECI2 Antibody</u>	CSB-PA527983LA01HU	Human, Rat, Mouse	ELISA, WB, IHC
<u>ECM1 Antibody</u>	CSB-PA007383LA01HU	Human	ELISA, WB, IHC
<u>ecoRVM Antibody</u>	CSB-PA366147LA01ENL	Escherichia coli	ELISA
<u>ECSCR Antibody</u>	CSB-PA007385LA01HU	Human	ELISA, IHC, IF
<u>ECT2 Antibody</u>	CSB-PA888001LA01HU	Human	ELISA, WB, IHC, IP
<u>EDAR Antibody</u>	CSB-PA892357LA01HU	Human	ELISA, WB, IHC, IF
<u>EDC4 Antibody</u>	CSB-PA744034LA01HU	Human	ELISA, IHC, IF
<u>EDEM1 Antibody</u>	CSB-PA838793LA01HU	Human	ELISA, IHC
<u>EDEM2 Antibody</u>	CSB-PA861146LA01HU	Human	ELISA, IHC
<u>EDNRB Antibody</u>	CSB-PA007404LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>EEA1 Antibody</u>	CSB-PA007405LA01HU	Human	ELISA, WB, IHC, IF
<u>EEF1A1 Antibody</u>	CSB-PA007409LA01HU	Human	ELISA, IHC, IF
<u>EEF1AKMT1 Antibody</u>	CSB-PA819896HA01HU	Human	ELISA, WB
<u>EEF1D Antibody</u>	CSB-PA007431LA01HU	Human	ELISA, WB
<u>EEF1E1 Antibody</u>	CSB-PA007432LA01HU	Human	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>EEF2 Antibody</u>	CSB-PA007434LA01HU	Human	ELISA, IHC, IF
<u>EEF2K Antibody</u>	CSB-PA007435LA01HU	Human	ELISA, IHC, IF
<u>EEF2KMT Antibody</u>	CSB-PA853420LA01HU	Human	ELISA, WB, IHC
<u>EFCAB1 Antibody</u>	CSB-PA862053LA01HU	Human	ELISA
<u>EFCAB12 Antibody</u>	CSB-PA740681LA01HU	Human	ELISA, IHC
<u>EFCAB14 Antibody</u>	CSB-PA04499A0Rb	Human	ELISA, WB, IHC
<u>EFCAB7 Antibody</u>	CSB-PA007447LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>EFEMP1 Antibody</u>	CSB-PA007450LA01HU	Human	ELISA, IHC, IF
<u>EFEMP2 Antibody</u>	CSB-PA007451LA01HU	Human, Mouse	ELISA, WB, IHC
<u>EFHB Antibody</u>	CSB-PA822723LA01HU	Human	ELISA, WB
<u>EFHD1 Antibody</u>	CSB-PA880107LA01HU	Human	ELISA, WB, IHC, IF
<u>EFHD2 Antibody</u>	CSB-PA856907LA01HU	Human	ELISA, IHC, IF
<u>Efna1 Antibody</u>	CSB-PA007460LA01RA	Rat	ELISA
<u>Efna5 Antibody</u>	CSB-PA007464LA01RA	Rat	ELISA
<u>Efnb3 Antibody</u>	CSB-PA007467LA01MO	Mouse, Human	ELISA, WB
<u>EFR3B Antibody</u>	CSB-PA897089LA01HU	Human	ELISA, IHC
<u>EGF Antibody</u>	CSB-PA06244A0Rb	Human	ELISA
<u>egf Antibody</u>	CSB-PA24449A0Rb	Zebrafish	ELISA, WB
<u>EGFL6 Antibody</u>	CSB-PA811602LA01HU	Human	ELISA, WB, IHC
<u>EGFL7 Antibody</u>	CSB-PA890692LA01HU	Human	ELISA, IHC
<u>EGFLAM Antibody</u>	CSB-PA723649LA01HU	Human, Mouse	ELISA, WB
<u>EGFR Antibody</u>	CSB-PA10279A0Rb	Human	ELISA, IHC, IF
<u>EGFR Antibody</u>	CSB-PA005571LA01HU	Human	ELISA, WB, IF
<u>egfra Antibody</u>	CSB-PA20899A0Rb	Zebrafish	ELISA, WB
<u>EGLN2 Antibody</u>	CSB-PA007482EA01HU	Human	ELISA, IHC, IF
<u>Egln3 Antibody</u>	CSB-PA14797A0Rb	Mouse, Rat	ELISA, WB
<u>EGR1 Antibody</u>	CSB-PA007484LA01HU	Human	ELISA
<u>EGR2 Antibody</u>	CSB-PA007485LA01HU	Human, Rat	ELISA, WB, IF
<u>EGR3 Antibody</u>	CSB-PA007486LA01HU	Human	ELISA, WB
<u>egsA Antibody</u>	CSB-PA455307HA01CFM	Geobacillus stearothermophilus	ELISA
<u>EHBP1 Antibody</u>	CSB-PA836723LA01HU	Human	ELISA
<u>EHD1 Antibody</u>	CSB-PA884470LA01HU	Human	ELISA, WB, IHC
<u>EHD2 Antibody</u>	CSB-PA873710LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>EHD3 Antibody</u>	CSB-PA878940LA01HU	Human	ELISA, IHC, IF
<u>EHHADH Antibody</u>	CSB-PA600104LA01HU	Human, Mouse	ELISA, WB, IHC
<u>EIF1 Antibody</u>	CSB-PA007503LA01HU	Human	ELISA, IF
<u>EIF1AD Antibody</u>	CSB-PA850842LA01HU	Human	ELISA



Product Name	Code	Species Reactivity	Applications
<u>EIF2AK1 Antibody</u>	CSB-PA007510LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>EIF2AK2 Antibody</u>	CSB-PA007511LA01HU	Human	ELISA, IHC, IF
<u>EIF2AK3 Antibody</u>	CSB-PA007512LA01HU	Human	ELISA, IHC, IF
<u>EIF2S1 Antibody</u>	CSB-PA007523LA01HU	Human	ELISA, IHC, IF
<u>EIF3A Antibody</u>	CSB-PA623903LA01HU	Human	ELISA, WB, IF
<u>EIF3C Antibody</u>	CSB-PA857867LA01HU	Human	ELISA, WB, IHC, IF
<u>EIF3G Antibody</u>	CSB-PA03099A0Rb	Human	ELISA, IHC, IF
<u>EIF3I Antibody</u>	CSB-PA03105A0Rb	Human, Rat	ELISA, WB, IF
<u>EIF3J Antibody</u>	CSB-PA007539LA01HU	Human	ELISA, IHC
<u>EIF3M Antibody</u>	CSB-PA745333LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>EIF4A2 Antibody</u>	CSB-PA00384A0Rb	Human	ELISA, WB, IHC, IF
<u>EIF4A2 Antibody</u>	CSB-PA00385A0Rb	Human	ELISA, WB
<u>EIF4A3 Antibody</u>	CSB-PA007553LA01HU	Human	ELISA, WB, IHC, IF
<u>EIF4B Antibody</u>	CSB-PA007554LA01HU	Human	ELISA, IHC, IF
<u>EIF4E Antibody</u>	CSB-PA007556LA01HU	Human	ELISA, IHC, IF
<u>EIF4EBP1 Antibody</u>	CSB-PA623815LA01HU	Human	ELISA, IHC, IF
<u>EIF4EBP2 Antibody</u>	CSB-PA007563DA01HU	Human	ELISA, IHC
<u>EIF4EBP2 Antibody</u>	CSB-PA007563YA01HU	Human, Mouse	ELISA, WB, IHC
<u>EIF4EBP3 Antibody</u>	CSB-PA007565HA01HU	Human	ELISA, IHC
<u>EIF4ENIF1 Antibody</u>	CSB-PA873633LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>EIF4G2 Antibody</u>	CSB-PA007569LA01HU	Human	ELISA, WB, IHC, IF
<u>EIF4G3 Antibody</u>	CSB-PA007570OA01HU	Human	ELISA, IHC, IF
<u>EIF5 Antibody</u>	CSB-PA007572HA01HU	Human	ELISA, IHC, IF
<u>EIF5A2 Antibody</u>	CSB-PA872418LA01HU	Human	ELISA, IHC, IF
<u>EIF6 Antibody</u>	CSB-PA007582LA01HU	Human	ELISA, IHC, IF
<u>EIPR1 Antibody</u>	CSB-PA681508LA01HU	Human	ELISA, IHC
<u>ELAC1 Antibody</u>	CSB-PA880988LA01HU	Human	ELISA, IHC
<u>ELANE Antibody</u>	CSB-PA007587LA01HU	Human	ELISA, IHC, IF
<u>ELAVL1 Antibody</u>	CSB-PA618993HA01HU	Human	ELISA, WB
<u>ELAVL2 Antibody</u>	CSB-PA007593LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ELAVL4 Antibody</u>	CSB-PA007595LA01HU	Human	ELISA, WB, IHC, IF
<u>ELF3 Antibody</u>	CSB-PA007598LA01HU	Human	ELISA, IHC, IF
<u>ELK1 Antibody</u>	CSB-PA007603LA01HU	Human, Mouse	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>ELK3 Antibody</u>	CSB-PA007606LA01HU	Human	ELISA, IHC
<u>ELMO1 Antibody</u>	CSB-PA856404HA01HU	Human	ELISA, WB
<u>ELMO2 Antibody</u>	CSB-PA836237LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ELMO3 Antibody</u>	CSB-PA846590LA01HU	Human	ELISA, IHC, IF
<u>ELMOD2 Antibody</u>	CSB-PA815588LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ELMOD3 Antibody</u>	CSB-PA842670LA01HU	Human	ELISA, WB, IHC, IF
<u>ELOC Antibody</u>	CSB-PA023281EA01HU	Human	ELISA
<u>ELOVL1 Antibody</u>	CSB-PA887138LA01HU	Human	ELISA, IHC, IF
<u>ELOVL2 Antibody</u>	CSB-PA889139LA01HU	Human	ELISA, IHC, IF
<u>ELOVL4 Antibody</u>	CSB-PA880925LA01HU	Human	ELISA, WB, IHC, IF
<u>ELOVL5 Antibody</u>	CSB-PA882144LA01HU	Human	ELISA, IF
<u>ELP1 Antibody</u>	CSB-PA011571LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>ELP2 Antibody</u>	CSB-PA740224LA01HU	Human	ELISA, IHC, IF
<u>ELP3 Antibody</u>	CSB-PA884499LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>ELP4 Antibody</u>	CSB-PA853408LA01HU	Human	ELISA, IHC, IF
<u>ELP6 Antibody</u>	CSB-PA608122LA01HU	Human	ELISA, IHC, IF
<u>eltB Antibody</u>	CSB-PA315643LA01ENL	Escherichia coli	ELISA
<u>EMB Antibody</u>	CSB-PA754239LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>EMC4 Antibody</u>	CSB-PA023887OA01HU	Human	ELISA, IF
<u>EMC6 Antibody</u>	CSB-PA023898LA01HU	Human	ELISA, IHC, IF
<u>EMC8 Antibody</u>	CSB-PA005835LA01HU	Human	ELISA, WB, IHC, IF
<u>EMC9 Antibody</u>	CSB-PA897482LA01HU	Human	ELISA, IF
<u>EMD Antibody</u>	CSB-PA007633LA01HU	Human, Mouse	ELISA, WB, IHC, IF, IP
<u>EME1 Antibody</u>	CSB-PA856897LA01HU	Human	ELISA, IHC
<u>EMILIN1 Antibody</u>	CSB-PA896753LA01HU	Human	ELISA, IHC
<u>EML2 Antibody</u>	CSB-PA007643LA01HU	Human	ELISA, IHC, IF
<u>EML4 Antibody</u>	CSB-PA862066LA01HU	Human	ELISA, WB
<u>EMP1 Antibody</u>	CSB-PA007648LA01HU	Human	ELISA, IHC, IF
<u>EMP2 Antibody</u>	CSB-PA007649LA01HU	Human	ELISA, IHC, IF
<u>EMP3 Antibody</u>	CSB-PA007650LA01HU	Human, Mouse	ELISA, WB, IF
<u>EMX2 Antibody</u>	CSB-PA007657LA01HU	Human, Rat, Mouse	ELISA, WB, IHC
<u>EN1 Antibody</u>	CSB-PA007659LA01CH	Chicken	ELISA
<u>EN2 Antibody</u>	CSB-PA007660LA01HU	Human	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>ENC1 Antibody</u>	CSB-PA007663LA01HU	Human, Mouse	ELISA, WB
<u>ENDOG Antibody</u>	CSB-PA007665LA01HU	Human, Mouse	ELISA, WB
<u>ENG Antibody</u>	CSB-PA007667LA01HU	Human	ELISA, IF
<u>ENHO Antibody</u>	CSB-PA744265OA01HU	Human	ELISA, IHC
<u>ENKD1 Antibody</u>	CSB-PA867120LA01HU	Human, Mouse	ELISA, WB, IHC
<u>eno Antibody</u>	CSB-PA351974LA01SKX	Staphylococcus aureus	ELISA, WB
<u>Eno1 Antibody</u>	CSB-PA007670HA01MO	Human, Mouse	ELISA, WB
<u>ENO1 Antibody</u>	CSB-PA02395A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>ENO2 Antibody</u>	CSB-PA007672YA01HU	Human	ELISA
<u>ENO2 Antibody</u>	CSB-PA09764A0Rb	Human, Mouse	ELISA, WB, IHC
<u>ENO3 Antibody</u>	CSB-PA13169A0Rb	Human	ELISA, IHC, IF
<u>ENOSF1 Antibody</u>	CSB-PA748723LA01HU	Human	ELISA, IHC
<u>ENOX1 Antibody</u>	CSB-PA007676LA01HU	Human	ELISA, WB, IHC
<u>ENPEP Antibody</u>	CSB-PA14759A0Rb	Human	ELISA, IHC
<u>ENPP1 Antibody</u>	CSB-PA007679LA01HU	Human	ELISA, IF
<u>ENPP6 Antibody</u>	CSB-PA740915LA01HU	Human	ELISA, IHC, IF
<u>ENPP7 Antibody</u>	CSB-PA751012LA01HU	Human, Mouse	ELISA, WB, IF
<u>ENSA Antibody</u>	CSB-PA01704A0Rb	Human	ELISA, IHC, IF
<u>entC1 Antibody</u>	CSB-PA365590HA01FKZ	Staphylococcus aureus	ELISA
<u>entC2 Antibody</u>	CSB-PA339410HA01FKZ	Staphylococcus aureus	ELISA
<u>entC3 Antibody</u>	CSB-PA357923HA01FKZ	Staphylococcus aureus	ELISA
<u>entD Antibody</u>	CSB-PA324966LA01FKZ	Staphylococcus aureus	ELISA
<u>entE Antibody</u>	CSB-PA320170HA01FKZ	Staphylococcus aureus	ELISA
<u>entG Antibody</u>	CSB-PA357924LA01SKY	Staphylococcus aureus	ELISA
<u>entG Antibody</u>	CSB-PA362954HA01SKX	Staphylococcus aureus	ELISA
<u>entH Antibody</u>	CSB-PA357926LA01FKZ	Staphylococcus aureus	ELISA
<u>ENTPD1 Antibody</u>	CSB-PA007690LA01HU	Human	ELISA, WB, IHC, IF
<u>ENTPD3 Antibody</u>	CSB-PA007692LA01HU	Human	ELISA, IHC, IF
<u>ENTPD7 Antibody</u>	CSB-PA007696LA01HU	Human, Mouse	ELISA, WB, IF
<u>ENTPD8 Antibody</u>	CSB-PA007697LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ENTR1 Antibody</u>	CSB-PA020897LA01HU	Human	ELISA, IHC, IF
<u>ENY2 Antibody</u>	CSB-PA889067LA01HU	Human	ELISA
<u>EOMES Antibody</u>	CSB-PA007701LA01HU	Human, Rat	ELISA, WB, IHC
<u>EP300 Antibody</u>	CSB-PA007702HA01HU	Human	ELISA, IHC, IF
<u>EPAS1 Antibody</u>	CSB-PA860780LA01HU	Human	ELISA, IHC, IF
<u>EPB2 Antibody</u>	CSB-PA326194LA01HWQ	Hordeum vulgare	ELISA
<u>EPCAM Antibody</u>	CSB-PA10799A0Rb	Human, Rat	ELISA, WB, IHC, IF
<u>EPHA1 Antibody</u>	CSB-PA007720LA01HU	Human, Mouse	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>EPHA2 Antibody</u>	CSB-PA007722LA01HU	Human	ELISA, IHC, IF
<u>EPHA4 Antibody</u>	CSB-PA007724LA01HU	Human	ELISA, IHC, IF
<u>EPHA5 Antibody</u>	CSB-PA007725LA01HU	Human	ELISA, WB, IHC, IF
<u>EPHA6 Antibody</u>	CSB-PA890673LA01HU	Human	ELISA, IHC
<u>EPHA7 Antibody</u>	CSB-PA618004LA01HU	Human	ELISA, IHC, IF
<u>EPHA8 Antibody</u>	CSB-PA007728LA01HU	Human	ELISA, IF
<u>EPHB2 Antibody</u>	CSB-PA007730LA01HU	Human	ELISA, IHC, IF
<u>EPHB3 Antibody</u>	CSB-PA007731LA01HU	Human, Mouse, Rat	ELISA, WB, IF
<u>EPHB6 Antibody</u>	CSB-PA007733LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>EPHX1 Antibody</u>	CSB-PA007734LA01HU	Human	ELISA, WB, IHC
<u>EPHX2 Antibody</u>	CSB-PA007735LA01HU	Human	ELISA, WB, IHC
<u>EPM2A Antibody</u>	CSB-PA007738LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>EPM2AIP1 Antibody</u>	CSB-PA759172LA01HU	Human	ELISA, IHC
<u>EPN2 Antibody</u>	CSB-PA007741LA01HU	Human	ELISA, IF
<u>EPN3 Antibody</u>	CSB-PA867130LA01HU	Human	ELISA, IHC, IF
<u>EPO Antibody</u>	CSB-PA007743LA01BO	Bovine, Human	ELISA, WB
<u>EPO Antibody</u>	CSB-PA007743LA01MOV	Macaca fascicularis	ELISA
<u>EPO Antibody</u>	CSB-PA06337A0Rb	Human	ELISA, IHC, IF
<u>EPO Antibody</u>	CSB-PA06339A0Rb	Human	ELISA, IHC
<u>EPOR Antibody</u>	CSB-PA007744LA01HU	Human	ELISA, WB, IHC, IF
<u>EPPK1 Antibody</u>	CSB-PA007745LA01HU	Human	ELISA
<u>EPRS Antibody</u>	CSB-PA007747LA01HU	Human	ELISA, WB, IHC, IF
<u>EPS8 Antibody</u>	CSB-PA615539LA01HU	Human	ELISA, IHC
<u>EPS8L2 Antibody</u>	CSB-PA887991LA01HU	Human	ELISA, IHC
<u>EPS8L3 Antibody</u>	CSB-PA823469LA01HU	Human	ELISA, IHC
<u>EPSTI1 Antibody</u>	CSB-PA839347LA01HU	Human	ELISA, WB, IHC
<u>EPX Antibody</u>	CSB-PA007756LA01HU	Human	ELISA, IHC
<u>ERAL1 Antibody</u>	CSB-PA007759LA01HU	Human	ELISA, IHC, IF
<u>ERAP1 Antibody</u>	CSB-PA007760LA01HU	Human	ELISA, WB, IHC, IF
<u>ERAP2 Antibody</u>	CSB-PA007761LA01HU	Human	ELISA, IHC
<u>ERAS Antibody</u>	CSB-PA13979A0Rb	Human	ELISA, WB
<u>ERBB2 Antibody</u>	CSB-PA09759A0Rb	Human	ELISA, IHC
<u>ERBB4 Antibody</u>	CSB-PA007766LA01HU	Human	ELISA, IHC, IF
<u>ERC1 Antibody</u>	CSB-PA007767LA01HU	Human	ELISA, WB, IF
<u>ERCC1 Antibody</u>	CSB-PA007769LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>ERCC4 Antibody</u>	CSB-PA849795HA01HU	Human	ELISA, IHC, IF
<u>ERCC6L2 Antibody</u>	CSB-PA004261LA01HU	Human	ELISA, IHC
<u>ERG Antibody</u>	CSB-PA007781HA01HU	Human, Rat	ELISA, IHC, IF, IP
<u>ERGIC1 Antibody</u>	CSB-PA842630LA01HU	Human	ELISA, WB, IHC, IF
<u>ERGIC2 Antibody</u>	CSB-PA857020EA01HU	Human, Rat, Mouse	ELISA, WB, IHC, IF
<u>ERGIC3 Antibody</u>	CSB-PA896688LA01HU	Human	ELISA, WB, IHC, IF
<u>ERH Antibody</u>	CSB-PA02745A0Rb	Human	ELISA, IHC
<u>ERI1 Antibody</u>	CSB-PA22449A0Rb	Human, Rat	ELISA, WB, IHC, IF
<u>ERLIN1 Antibody</u>	CSB-PA007790LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ERLIN2 Antibody</u>	CSB-PA007791LA01HU	Human	ELISA, IHC, IF
<u>ERN1 Antibody</u>	CSB-PA007795LA01HU	Human	ELISA, IHC
<u>ERO1B Antibody</u>	CSB-PA774834LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ERP27 Antibody</u>	CSB-PA842657LA01HU	Human	ELISA, IHC
<u>ERP29 Antibody</u>	CSB-PA05344A0Rb	Human	ELISA, WB, IHC
<u>ERP29 Antibody</u>	CSB-PA05345A0Rb	Human	ELISA, WB, IHC
<u>ERVK-6 Antibody</u>	CSB-PA724336LA01HU	Human	ELISA, IHC, IF
<u>ERVW-1 Antibody</u>	CSB-PA891578LA01HU	Human	ELISA, IHC, IF
<u>ESAM Antibody</u>	CSB-PA850258LA01HU	Human	ELISA, IHC
<u>ESD Antibody</u>	CSB-PA007823LA01HU	Human	ELISA, IF
<u>ESF1 Antibody</u>	CSB-PA867152LA01HU	Human	ELISA
<u>ESM1 Antibody</u>	CSB-PA007825LA01HU	Human	ELISA, IHC
<u>ESM1 Antibody</u>	CSB-PA007825HA01HU	Human	ELISA, IHC
<u>ESPL1 Antibody</u>	CSB-PA621785LA01HU	Human	ELISA, IF
<u>ESR1 Antibody</u>	CSB-PA007830YA01HU	Human	ELISA, IHC, IF
<u>ESR1 Antibody</u>	CSB-PA11399A0Rb	Human, Rat	ELISA, WB, IHC, ChIP
<u>ESR2 Antibody</u>	CSB-PA007831LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ESRG Antibody</u>	CSB-PA632275LA01HU	Human	ELISA, IHC
<u>ESRP2 Antibody</u>	CSB-PA007833LA01HU	Human	ELISA, WB, IHC, IF
<u>ESRRB Antibody</u>	CSB-PA007836LA01HU	Human	ELISA, WB
<u>ESRRG Antibody</u>	CSB-PA007837LA01HU	Human	ELISA
<u>esxA Antibody</u>	CSB-PA763968HA01SKW	Staphylococcus aureus	ELISA, WB
<u>esxB Antibody</u>	CSB-PA363658YA01MVZ	Mycobacterium tuberculosis	ELISA
<u>ESYT2 Antibody</u>	CSB-PA008279LA01HU	Human	ELISA, IHC, IF
<u>eta Antibody</u>	CSB-PA357644LA01FKZ	Staphylococcus aureus	ELISA, WB
<u>E-Tag Monoclonal Antibody</u>	CSB-MA000151M0m	All	ELISA, WB, IP



Product Name	Code	Species Reactivity	Applications
<u>ETFDH Antibody</u>	CSB-PA007846LA01PI	Pig	ELISA
<u>ETFDH Antibody</u>	CSB-PA619056LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ETFRF1 Antibody</u>	CSB-PA753708LA01HU	Human	ELISA, IHC
<u>ETHE1 Antibody</u>	CSB-PA007847HA01HU	Human, Mouse	ELISA, WB, IHC
<u>ETNK1 Antibody</u>	CSB-PA007850LA01HU	Human	ELISA, WB, IHC, IF
<u>ETS1 Antibody</u>	CSB-PA007852EA01HU	Human	ELISA
<u>ETS2 Antibody</u>	CSB-PA007853LA01HU	Human, Mouse	ELISA, WB
<u>ETV1 Antibody</u>	CSB-PA007854LA01HU	Human, Mouse, Rat	ELISA, WB
<u>ETV2 Antibody</u>	CSB-PA007855LA01HU	Human	ELISA, IF
<u>ETV5 Antibody</u>	CSB-PA007859LA01HU	Human	ELISA, WB, IHC
<u>ETV7 Antibody</u>	CSB-PA897576LA01HU	Human, Rat	ELISA, WB, IHC
<u>eutC Antibody</u>	CSB-PA18239A0Rb	Escherichia coli	ELISA, WB
<u>EVA1C Antibody</u>	CSB-PA003777LA01HU	Human	ELISA, WB, IHC
<u>EVC Antibody</u>	CSB-PA007862LA01HU	Human	ELISA, WB, IHC
<u>EVC2 Antibody</u>	CSB-PA769777LA01HU	Human	ELISA, IHC, IF
<u>EVI2A Antibody</u>	CSB-PA007865LA01HU	Human	ELISA, IHC, IF
<u>EXD2 Antibody</u>	CSB-PA873662LA01HU	Human, Mouse	ELISA, WB, IHC
<u>EXO1 Antibody</u>	CSB-PA007876LA01HU	Human	ELISA, WB, IHC
<u>EXOC3 Antibody</u>	CSB-PA527663LA01HU	Human	ELISA, WB
<u>EXOC7 Antibody</u>	CSB-PA007886LA01HU	Human	ELISA, WB, IHC, IF
<u>EXOC8 Antibody</u>	CSB-PA007887LA01HU	Human	ELISA, IF
<u>EXOSC1 Antibody</u>	CSB-PA897103LA01HU	Human	ELISA, IHC
<u>EXOSC10 Antibody</u>	CSB-PA007889LA01HU	Human	ELISA, IHC
<u>EXOSC3 Antibody</u>	CSB-PA889083LA01HU	Human	ELISA, IHC
<u>EXOSC5 Antibody</u>	CSB-PA878857LA01HU	Human	ELISA, WB, IHC
<u>EXOSC8 Antibody</u>	CSB-PA853387LA01HU	Human, Mouse	ELISA, WB
<u>EXPH5 Antibody</u>	CSB-PA843342LA01HU	Human	ELISA, IHC
<u>EXT1 Antibody</u>	CSB-PA620959LA01HU	Human	ELISA, WB
<u>EXT2 Antibody</u>	CSB-PA849974LA01HU	Human	ELISA, IHC
<u>EXTL3 Antibody</u>	CSB-PA03075A0Rb	Human	ELISA, WB, IHC, IF
<u>EYA1 Antibody</u>	CSB-PA857862LA01HU	Human	ELISA, IHC, IF
<u>EYA2 Antibody</u>	CSB-PA007907LA01HU	Human	ELISA, IHC
<u>EYA4 Antibody</u>	CSB-PA007909LA01HU	Human	ELISA, IHC, IF
<u>EZH1 Antibody</u>	CSB-PA852913DA01HU	Human	ELISA, IHC
<u>EZR Antibody</u>	CSB-PA02159A0Rb	Human	ELISA, IHC, IF
<u>F Antibody</u>	CSB-PA318261LA01CCQ	Canine distemper virus	ELISA
<u>F Antibody</u>	CSB-PA356041LA01HTO	Human respiRatory syncytial virus A	ELISA
<u>F11 Antibody</u>	CSB-PA007916LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>F11R Antibody</u>	CSB-PA897579LA01HU	Human	ELISA, IHC, IF
<u>F12 Antibody</u>	CSB-PA007918LA01PI	Pig	ELISA
<u>F12 Antibody</u>	CSB-PA007918LA01HU	Human	ELISA, IHC
<u>F13A1 Antibody</u>	CSB-PA007919LA01HU	Human	ELISA, IHC, IF
<u>F13B Antibody</u>	CSB-PA007922LA01HU	Human	ELISA, WB
<u>F2 Antibody</u>	CSB-PA007923LA01HU	Human	ELISA, WB, IF
<u>F2R Antibody</u>	CSB-PA007924LA01HU	Human	ELISA, IHC, IF
<u>F2RL1 Antibody</u>	CSB-PA007925LA01HU	Human	ELISA, IHC, IF
<u>F3 Antibody</u>	CSB-PA007928LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>F5 Antibody</u>	CSB-PA007929HA01HU	Human	ELISA, IHC, IF
<u>F7 Antibody</u>	CSB-PA007930HA01HU	Human	ELISA, IHC
<u>F8 Antibody</u>	CSB-PA007932LA01HU	Human	ELISA, IHC, IF
<u>F9 Antibody</u>	CSB-PA007936LA01HU	Human	ELISA, IHC, IF
<u>FAAH Antibody</u>	CSB-PA007938LA01HU	Human	ELISA, IHC, IF
<u>FAAP20 Antibody</u>	CSB-PA750774HA01HU	Human	ELISA, IHC
<u>fabD Antibody</u>	CSB-PA686101LA01FLB	Staphylococcus aureus	ELISA
<u>fabF Antibody</u>	CSB-PA359466LA01ENV	Escherichia coli	ELISA
<u>fabF Antibody</u>	CSB-PA692016LA01FLB	Staphylococcus aureus	ELISA
<u>fabI Antibody</u>	CSB-PA750140HA01SKV	Staphylococcus aureus, Mouse	ELISA, WB
<u>FABP1 Antibody</u>	CSB-PA007940YA01HU	Human, Mouse	ELISA, WB, IHC
<u>FABP3 Antibody</u>	CSB-PA007943YA01HU	Human	ELISA
<u>FABP3 Antibody</u>	CSB-PA07764A0Rb	Human, Mouse	ELISA, WB, IHC
<u>FABP3 Antibody</u>	CSB-PA07767A0Rb	Human	ELISA
<u>FABP4 Antibody</u>	CSB-PA007945LA01HU	Human	ELISA, WB, IHC
<u>FABP7 Antibody</u>	CSB-PA007956LA01HU	Human	ELISA, IHC, IF
<u>FABP9 Antibody</u>	CSB-PA604514HA01HU	Human	ELISA, IF
<u>FADS1 Antibody</u>	CSB-PA007959LA01HU	Human	ELISA, WB, IHC, IF
<u>FAF1 Antibody</u>	CSB-PA892359LA01HU	Human, Mouse	ELISA, WB, IHC
<u>FAH Antibody</u>	CSB-PA007965LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>FAHD2A Antibody</u>	CSB-PA822218LA01HU	Human	ELISA, WB, IHC
<u>FAIM Antibody</u>	CSB-PA889122LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>FAM102B Antibody</u>	CSB-PA735976LA01HU	Human	ELISA, IHC, IF
<u>FAM105A Antibody</u>	CSB-PA873659LA01HU	Human	ELISA, WB, IHC
<u>FAM107A Antibody</u>	CSB-PA007987LA01HU	Human	ELISA, IHC
<u>FAM107B Antibody</u>	CSB-PA884429LA01HU	Human, Mouse	ELISA, WB
<u>FAM109B Antibody</u>	CSB-PA753688LA01HU	Human	ELISA
<u>FAM110B Antibody</u>	CSB-PA008005LA01HU	Human, Mouse	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>FAM110C Antibody</u>	CSB-PA637287LA01HU	Human	ELISA, IHC
<u>FAM111A Antibody</u>	CSB-PA853484LA01HU	Human	ELISA, IHC, IF
<u>FAM114A2 Antibody</u>	CSB-PA878875LA01HU	Human, Mouse	ELISA, WB, IHC
<u>FAM117B Antibody</u>	CSB-PA754170LA01HU	Human	ELISA, IHC
<u>FAM118A Antibody</u>	CSB-PA885747LA01HU	Human	ELISA, WB, IHC, IF
<u>FAM118B Antibody</u>	CSB-PA858411LA01HU	Human	ELISA, WB, IHC
<u>FAM122B Antibody</u>	CSB-PA768758LA01HU	Human, Mouse	ELISA, WB, IHC
<u>FAM124A Antibody</u>	CSB-PA801237LA01HU	Human	ELISA, IHC
<u>FAM124B Antibody</u>	CSB-PA880982LA01HU	Human	ELISA, WB, IHC
<u>FAM126B Antibody</u>	CSB-PA814227LA01HU	Human	ELISA, WB, IHC
<u>FAM131A Antibody</u>	CSB-PA744268LA01HU	Human	ELISA, IHC
<u>FAM131B Antibody</u>	CSB-PA774828LA01HU	Human	ELISA
<u>FAM134B Antibody</u>	CSB-PA875683LA01HU	Human	ELISA, IHC
<u>FAM135B Antibody</u>	CSB-PA677719LA01HU	Human	ELISA, IHC
<u>FAM136A Antibody</u>	CSB-PA836190LA01HU	Human	ELISA, WB, IHC
<u>FAM13A Antibody</u>	CSB-PA008075LA01HU	Human	ELISA, IHC
<u>FAM149B1 Antibody</u>	CSB-PA839295LA01HU	Human, Mouse	ELISA, WB, IHC
<u>FAM151B Antibody</u>	CSB-PA740923LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>FAM153C Antibody</u>	CSB-PA676570LA01HU	Human	ELISA, IHC
<u>FAM156A Antibody</u>	CSB-PA847690LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>FAM160A2 Antibody</u>	CSB-PA822711LA01HU	Human	ELISA, IHC
<u>FAM160B1 Antibody</u>	CSB-PA726129LA01HU	Human	ELISA, IHC
<u>FAM161A Antibody</u>	CSB-PA661553LA01HU	Human, Mouse	ELISA, WB, IHC
<u>FAM163A Antibody</u>	CSB-PA853423LA01HU	Human	ELISA, IHC
<u>FAM166A Antibody</u>	CSB-PA753727LA01HU	Human	ELISA, IHC
<u>FAM167A Antibody</u>	CSB-PA836244LA01HU	Human	ELISA, IHC
<u>FAM168B Antibody</u>	CSB-PA008125LA01HU	Human	ELISA, IHC, IF
<u>FAM171B Antibody</u>	CSB-PA757601LA01HU	Human	ELISA, IHC, IF
<u>FAM174A Antibody</u>	CSB-PA851535LA01HU	Human	ELISA, IHC
<u>FAM182A Antibody</u>	CSB-PA729200LA01HU	Human	ELISA, IF
<u>FAM187B Antibody</u>	CSB-PA623168LA01HU	Human	ELISA, IHC
<u>FAM19A1 Antibody</u>	CSB-PA770353LA01HU	Human	ELISA
<u>FAM19A3 Antibody</u>	CSB-PA773794LA01HU	Human	ELISA
<u>FAM19A4 Antibody</u>	CSB-PA836253LA01HU	Human	ELISA
<u>FAM20C Antibody</u>	CSB-PA816901LA01HU	Human	ELISA, WB, IHC
<u>FAM222B Antibody</u>	CSB-PA851963LA01HU	Human, Mouse	ELISA, WB, IHC
<u>FAM234A Antibody</u>	CSB-PA887949LA01HU	Human	ELISA, IF
<u>FAM234B Antibody</u>	CSB-PA012262LA01HU	Human	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>FAM32A Antibody</u>	CSB-PA897549LA01HU	Human	ELISA, IHC
<u>FAM3A Antibody</u>	CSB-PA008226LA01HU	Human	ELISA, IHC
<u>FAM3B Antibody</u>	CSB-PA008227LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>FAM3C Antibody</u>	CSB-PA856401LA01HU	Human	ELISA, WB
<u>FAM43A Antibody</u>	CSB-PA854028LA01HU	Human, Mouse	ELISA, WB, IHC
<u>FAM45A Antibody</u>	CSB-PA840565LA01HU	Human	ELISA
<u>FAM49B Antibody</u>	CSB-PA873658LA01HU	Human	ELISA
<u>FAM50B Antibody</u>	CSB-PA897521LA01HU	Human	ELISA, IHC, IF
<u>FAM53B Antibody</u>	CSB-PA617903LA01HU	Human	ELISA, WB, IHC
<u>FAM60A Antibody</u>	CSB-PA882059LA01HU	Human	ELISA, IHC
<u>FAM69A Antibody</u>	CSB-PA732928LA01HU	Human	ELISA, WB
<u>FAM71A Antibody</u>	CSB-PA810307LA01HU	Human	ELISA, IHC, IF
<u>FAM71C Antibody</u>	CSB-PA843337LA01HU	Human	ELISA, WB
<u>FAM71D Antibody</u>	CSB-PA818742LA01HU	Human	ELISA, WB
<u>FAM71E2 Antibody</u>	CSB-PA847639LA01HU	Human	ELISA
<u>FAM71F1 Antibody</u>	CSB-PA853446LA01HU	Human	ELISA, IHC
<u>FAM71F2 Antibody</u>	CSB-PA750765LA01HU	Human	ELISA, IHC
<u>FAM76B Antibody</u>	CSB-PA702405LA01HU	Human	ELISA
<u>FAM78A Antibody</u>	CSB-PA711454LA01HU	Human	ELISA, WB, IHC
<u>FAM83D Antibody</u>	CSB-PA887979LA01HU	Human	ELISA, WB, IF
<u>FAM83G Antibody</u>	CSB-PA008345LA01HU	Human	ELISA, IHC, IF
<u>FAM84A Antibody</u>	CSB-PA846652LA01HU	Human	ELISA, IHC
<u>FAM89B Antibody</u>	CSB-PA854045LA01HU	Human	ELISA, IHC, IF
<u>FAM92A Antibody</u>	CSB-PA008396LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>FAN1 Antibody</u>	CSB-PA22859A0Rb	Human	ELISA, WB, IHC
<u>FANCA Antibody</u>	CSB-PA008413LA01HU	Human	ELISA, IF
<u>FANCB Antibody</u>	CSB-PA818751LA01HU	Human	ELISA, IF
<u>FANCC Antibody</u>	CSB-PA008415LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>FANCD2 Antibody</u>	CSB-PA866311LA01HU	Human	ELISA, IHC
<u>FANCF Antibody</u>	CSB-PA873614LA01HU	Human	ELISA, WB
<u>FANCI Antibody</u>	CSB-PA889120LA01HU	Human	ELISA, IF
<u>FANCL Antibody</u>	CSB-PA882120LA01HU	Human	ELISA, IHC, IF
<u>FANK1 Antibody</u>	CSB-PA819458HA01HU	Human	ELISA, IHC, IF
<u>FAR1 Antibody</u>	CSB-PA837860LA01HU	Human	ELISA, IHC
<u>FAR2 Antibody</u>	CSB-PA822239LA01HU	Human	ELISA, IHC, IF
<u>FARP1 Antibody</u>	CSB-PA897492LA01HU	Human	ELISA, WB, IHC, IF
<u>FARP2 Antibody</u>	CSB-PA008428LA01HU	Human	ELISA, IHC, IF
<u>FARSA Antibody</u>	CSB-PA897084LA01HU	Human	ELISA, IHC, IF
<u>FARSB Antibody</u>	CSB-PA882099LA01HU	Human, Mouse	ELISA, WB, IF



Product Name	Code	Species Reactivity	Applications
<u>FAS Antibody</u>	CSB-PA06354A0Rb	Human	ELISA, IHC, IF
<u>fas Antibody</u>	CSB-PA24199A0Rb	Zebrafish	ELISA, WB
<u>FASLG Antibody</u>	CSB-PA008434LA01HU	Human	ELISA, IHC, IF
<u>FASN Antibody</u>	CSB-PA13179A0Rb	Human	ELISA, IHC, IF
<u>FASTKD2 Antibody</u>	CSB-PA882149LA01HU	Human, Rat, Mouse	ELISA, WB, IHC
<u>FASTKD5 Antibody</u>	CSB-PA759173LA01HU	Human	ELISA, IHC
<u>FAT3 Antibody</u>	CSB-PA819472LA01HU	Human	ELISA, IHC, IF
<u>FAT4 Antibody</u>	CSB-PA740927LA01HU	Human	ELISA, IF
<u>FATE1 Antibody</u>	CSB-PA839274LA01HU	Human	ELISA, IHC, IF
<u>FAU Antibody</u>	CSB-PA01825A0Rb	Human	ELISA, IHC, IF
<u>FAU Antibody</u>	CSB-PA008447LA01HU	Human	ELISA, IHC, IF
<u>FBLN2 Antibody</u>	CSB-PA008453LA01HU	Human	ELISA, WB, IHC
<u>FBLN7 Antibody</u>	CSB-PA687502LA01HU	Human	ELISA, IHC
<u>FBN1 Antibody</u>	CSB-PA008456LA01HU	Human	ELISA, IHC
<u>FBN2 Antibody</u>	CSB-PA008457LA01HU	Human	ELISA, IF
<u>FBP2 Antibody</u>	CSB-PA008460LA01HU	Human, Rat	ELISA, WB, IF
<u>fbpA Antibody</u>	CSB-PA18279A0Rb	Mycobacterium tuberculosis	ELISA, WB
<u>fbpB Antibody</u>	CSB-PA324138LA01MKZ	Mycobacterium kansasii	ELISA, WB
<u>fbpC Antibody</u>	CSB-PA363614LA01MVZ	Mycobacterium tuberculosis	ELISA, WB
<u>FBXL17 Antibody</u>	CSB-PA866216LA01HU	Human	ELISA, IHC
<u>FBXL18 Antibody</u>	CSB-PA008469LA01HU	Human	ELISA, IHC, IF
<u>FBXL19 Antibody</u>	CSB-PA754242LA01HU	Human	ELISA, IF
<u>FBXL2 Antibody</u>	CSB-PA892134LA01HU	Human	ELISA, IF
<u>FBXL20 Antibody</u>	CSB-PA822228LA01HU	Human	ELISA, IHC
<u>FBXL21 Antibody</u>	CSB-PA891953LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>FBXL22 Antibody</u>	CSB-PA757543LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>FBXO10 Antibody</u>	CSB-PA890915LA01HU	Human	ELISA, WB, IHC, IF
<u>FBXO15 Antibody</u>	CSB-PA843324LA01HU	Human	ELISA, IHC
<u>FBXO24 Antibody</u>	CSB-PA008491LA01HU	Human, Mouse	ELISA, WB, IHC
<u>FBXO25 Antibody</u>	CSB-PA848401LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>FBXO28 Antibody</u>	CSB-PA889119LA01HU	Human	ELISA, IHC
<u>FBXO3 Antibody</u>	CSB-PA892132LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>FBXO30 Antibody</u>	CSB-PA844709LA01HU	Human, Mouse	ELISA, WB, IHC
<u>FBXO40 Antibody</u>	CSB-PA883385LA01HU	Human	ELISA, IHC
<u>FBXO44 Antibody</u>	CSB-PA872497LA01HU	Human	ELISA



Product Name	Code	Species Reactivity	Applications
<u>FBXO9 Antibody</u>	CSB-PA891719LA01HU	Human	ELISA, IHC, IF
<u>FBXW10 Antibody</u>	CSB-PA726395LA01HU	Human	ELISA, IHC
<u>Fbxw10 Antibody</u>	CSB-PA725549LA01MO	Mouse	ELISA
<u>FBXW12 Antibody</u>	CSB-PA751062LA01HU	Human	ELISA, IHC
<u>FBXW2 Antibody</u>	CSB-PA891730LA01HU	Human, Mouse	ELISA, WB, IHC
<u>FBXW5 Antibody</u>	CSB-PA846577LA01HU	Human	ELISA, WB, IHC, IF
<u>FBXW7 Antibody</u>	CSB-PA822163LA01HU	Human	ELISA, IHC
<u>FCER1A Antibody</u>	CSB-PA008532LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>FCER1G Antibody</u>	CSB-PA008533LA01HU	Human	ELISA, IHC, IF
<u>FCER2 Antibody</u>	CSB-PA02219A0Rb	Human	ELISA, IHC, IF
<u>FCGR2A Antibody</u>	CSB-PA04459A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>FCGR2B Antibody</u>	CSB-PA10839A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>Fcgr3 Antibody</u>	CSB-PA357412LA01MO	Mouse	ELISA
<u>FCGR3A Antibody</u>	CSB-PA008543LA01HU	Human	ELISA, IHC, IF
<u>FCGR3B Antibody</u>	CSB-PA008544LA01HU	Human	ELISA, IF
<u>FCGRT Antibody</u>	CSB-PA008545EA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>FCGRT Antibody</u>	CSB-PA851443LA01MOV	Macaca fascicularis	ELISA
<u>FCHO1 Antibody</u>	CSB-PA008546LA01HU	Human	ELISA, IHC, IF
<u>FCHO2 Antibody</u>	CSB-PA606615LA01HU	Human, Mouse	ELISA, WB, IF
<u>FCHSD2 Antibody</u>	CSB-PA008549LA01HU	Human	ELISA, WB, IHC
<u>Fcmmr Antibody</u>	CSB-PA007972LA01MO	Mouse, Rat	ELISA, WB
<u>FCRL3 Antibody</u>	CSB-PA853475LA01HU	Human	ELISA, IHC, IF
<u>FCRL4 Antibody</u>	CSB-PA836277LA01HU	Human	ELISA, WB
<u>FCRL5 Antibody</u>	CSB-PA853494LA01HU	Human	ELISA, IF
<u>FCRL6 Antibody</u>	CSB-PA718779LA01HU	Human	ELISA
<u>FDFT1 Antibody</u>	CSB-PA008562LA01HU	Human	ELISA, WB, IHC, IF
<u>FDPS Antibody</u>	CSB-PA008563LA01HU	Human	ELISA, WB, IHC, IP
<u>FDXR Antibody</u>	CSB-PA008575LA01HU	Human, Mouse	ELISA, WB, IHC
<u>FEM1B Antibody</u>	CSB-PA891718LA01HU	Human, Mouse	ELISA, WB
<u>FEN1 Antibody</u>	CSB-PA008585LA01DOA	Arabidopsis thaliana	ELISA, WB
<u>FERD3L Antibody</u>	CSB-PA850424LA01HU	Human	ELISA, IHC, IF
<u>FERMT1 Antibody</u>	CSB-PA874786LA01HU	Human	ELISA, IHC, IF
<u>FERMT2 Antibody</u>	CSB-PA853382LA01HU	Human	ELISA, IHC, IF
<u>FERMT3 Antibody</u>	CSB-PA771441LA01HU	Human	ELISA, IHC, IF
<u>Ferritin-2, chloroplastic</u>	CSB-PA836022LA01GGV	Glycine max	ELISA



Product Name	Code	Species Reactivity	Applications
<u>Antibody</u>			
<u>FEV Antibody</u>	CSB-PA858717LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>FEZ1 Antibody</u>	CSB-PA03389A0Rb	Human	ELISA, IHC, IF
<u>FEZF1 Antibody</u>	CSB-PA008602LA01HU	Human	ELISA, WB, IHC
<u>FFAR1 Antibody</u>	CSB-PA008604LA01HU	Human	ELISA, IHC, IF
<u>FFAR2 Antibody</u>	CSB-PA008605LA01HU	Human	ELISA, IHC, IF
<u>FFAR3 Antibody</u>	CSB-PA008606LA01HU	Human	ELISA, IHC, IF
<u>FGA Antibody</u>	CSB-PA008607LA01HU	Human	ELISA, IHC
<u>FGB Antibody</u>	CSB-PA008608LA01HU	Human	ELISA, IHC
<u>FGD5 Antibody</u>	CSB-PA757939LA01HU	Human	ELISA, IHC, IF
<u>FGD6 Antibody</u>	CSB-PA761649EA01HU	Human	ELISA, WB, IHC
<u>Fgf1 Antibody</u>	CSB-PA008615HA01MO	Mouse, Human	ELISA, WB, IHC
<u>FGF1 Antibody</u>	CSB-PA06379A0Rb	Human	ELISA
<u>FGF13 Antibody</u>	CSB-PA008619LA01HU	Human, Mouse	ELISA, WB, IHC
<u>Fgf15 Antibody</u>	CSB-PA522052LA01MO	Mouse, Rat	ELISA, WB
<u>FGF16 Antibody</u>	CSB-PA008621LA01HU	Human	ELISA, IHC
<u>FGF17 Antibody</u>	CSB-PA008622LA01HU	Human, Mouse	ELISA, WB, IHC
<u>Fgf2 Antibody</u>	CSB-PA06259A0Rb	Mouse, Human	ELISA, WB
<u>FGF2 Antibody</u>	CSB-PA12847A0Rb	Human	ELISA, IHC
<u>FGF2 Antibody</u>	CSB-PA12849A0Rb	Human, Mouse	ELISA, WB, IHC
<u>Fgf2 Antibody</u>	CSB-PA319495LA01RA	Rat	ELISA, WB
<u>Fgf21 Antibody</u>	CSB-PA008627LA01MO	Mouse	ELISA
<u>FGF21 Antibody</u>	CSB-PA06404A0Rb	Human	ELISA, WB, IHC
<u>FGF22 Antibody</u>	CSB-PA008628LA01HU	Human	ELISA, IHC
<u>Fgf23 Antibody</u>	CSB-PA008629LA01RA	Rat	ELISA, WB
<u>FGF23 Antibody</u>	CSB-PA13919A0Rb	Human	ELISA, IF
<u>fgf3 Antibody</u>	CSB-PA008630LA01DIL	Zebrafish	ELISA
<u>FGF4 Antibody</u>	CSB-PA008631OA01HU	Human	ELISA, IHC
<u>FGF7 Antibody</u>	CSB-PA008634LA01HU	Human	ELISA
<u>FGF8 Antibody</u>	CSB-PA008637LA01HU	Human	ELISA, IHC
<u>FGFBP2 Antibody</u>	CSB-PA887162LA01HU	Human	ELISA, WB, IHC
<u>FGFBP3 Antibody</u>	CSB-PA819447OA01HU	Human	ELISA, IHC
<u>FGFR1 Antibody</u>	CSB-PA008642LA01HU	Human	ELISA, WB, IF
<u>FGFR1OP2 Antibody</u>	CSB-PA008644LA01HU	Human, Mouse	ELISA, WB
<u>FGFR2 Antibody</u>	CSB-PA008645LA01HU	Human	ELISA, IHC
<u>FGFR3 Antibody</u>	CSB-PA008646LA01HU	Human	ELISA, IHC, IF
<u>FGFR4 Antibody</u>	CSB-PA008648LA01HU	Human	ELISA, IHC, IF
<u>FGFRL1 Antibody</u>	CSB-PA008650LA01HU	Human	ELISA, IHC, IF
<u>FGG Antibody</u>	CSB-PA008651LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>FGGY Antibody</u>	CSB-PA853393LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
FGL1 Antibody	CSB-PA008653LA01HU	Human	ELISA, WB, IHC
FGL2 Antibody	CSB-PA619871LA01HU	Human	ELISA, IHC
FGR Antibody	CSB-PA008655LA01HU	Human	ELISA, IHC, IF
FH Antibody	CSB-PA008659LA01HU	Human, Mouse	ELISA, WB, IHC, IF
FHLT Antibody	CSB-PA11547A0Rb	Human	ELISA, IHC, IF
Fhit Antibody	CSB-PA530819HA01MO	Mouse	ELISA
FHL1 Antibody	CSB-PA619784LA01HU	Human, Mouse, Rat	ELISA, WB, IF
FHL2 Antibody	CSB-PA617906LA01HU	Human	ELISA, IHC, IF
FHL5 Antibody	CSB-PA725899LA01HU	Human	ELISA, IF
FHOD1 Antibody	CSB-PA896540LA01HU	Human	ELISA, IHC
FIBIN Antibody	CSB-PA851526LA01HU	Human	ELISA
FIBP Antibody	CSB-PA008671LA01HU	Human	ELISA, WB, IHC, IF
FICD Antibody	CSB-PA880112LA01HU	Human, Rat	ELISA, WB, IHC, IF
FIGN Antibody	CSB-PA008676LA01HU	Human, Mouse	ELISA, WB, IHC, IF
FILIP1 Antibody	CSB-PA800226LA01HU	Human	ELISA, IHC, IF
FILIP1L Antibody	CSB-PA679703LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
fim2 Antibody	CSB-PA356626LA01BUA	Bordetella pertussis	ELISA
fim3 Antibody	CSB-PA322971LA01BUA	Bordetella pertussis	ELISA
FIS1 Antibody	CSB-PA008684LA01HU	Human	ELISA, IHC, IF
FITC Antibody	CSB-PA0051110Rb		ELISA
FIZ1 Antibody	CSB-PA839415LA01HU	Human	ELISA, WB, IHC
FKBP10 Antibody	CSB-PA822180LA01HU	Human	ELISA, IHC
FKBP11 Antibody	CSB-PA878925LA01HU	Human, Mouse	ELISA, WB, IHC, IF
FKBP14 Antibody	CSB-PA865159LA01HU	Human, Mouse	ELISA, WB
FKBP1A Antibody	CSB-PA02325A0Rb	Human	ELISA, WB, IHC, IF
FKBP3 Antibody	CSB-PA008699HA01HU	Human	ELISA, WB, IHC, IF
FKBP4 Antibody	CSB-PA008700EA01HU	Human	ELISA, WB, IHC, IF
FKBP6 Antibody	CSB-PA008702LA01HU	Human	ELISA, WB, IHC, IF
FKBP8 Antibody	CSB-PA617913LA01HU	Human, Mouse	ELISA, WB, IHC
FKBP9P1 Antibody	CSB-PA758809LA01HU	Human	ELISA
fkpA Antibody	CSB-PA355178HA01EGX	Escherichia coli O6:H1	ELISA
fkpA Antibody	CSB-PA517042HA01AXW	Aeromonas hydrophila	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>FKRP Antibody</u>	CSB-PA862050LA01HU	Human	ELISA, IHC
<u>FKTN Antibody</u>	CSB-PA008709LA01HU	Human	ELISA, WB
<u>fla Antibody</u>	CSB-PA319964LA01BUD	Borrelia burgdorferi	ELISA
<u>FLAD1 Antibody</u>	CSB-PA854116LA01HU	Human	ELISA, WB, IHC
<u>Flag Tag Monoclonal Antibody</u>	CSB-MA000021M0m	All	ELISA, WB, IP
<u>FLCN Antibody</u>	CSB-PA008711HA01HU	Human	ELISA, IHC
<u>FLG Antibody</u>	CSB-PA008712HA01HU	Human	ELISA, IHC
<u>FLII Antibody</u>	CSB-PA008714LA01HU	Human	ELISA, WB, IHC, IF
<u>fliC Antibody</u>	CSB-PA09734A0Rb	Escherichia coli	ELISA
<u>Flii Antibody</u>	CSB-PA884922LA01MO	Mouse	ELISA
<u>FLII Antibody</u>	CSB-PA615661LA01HU	Human	ELISA, IHC, IF
<u>FLNA Antibody</u>	CSB-PA008724LA01HU	Human	ELISA, IHC, IF
<u>FLOT1 Antibody</u>	CSB-PA008727LA01HU	Human, Rat, Mouse	ELISA, WB, IHC, IF, IP
<u>FLRT1 Antibody</u>	CSB-PA865196EA01HU	Human, Mouse	ELISA, WB, IHC
<u>FLRT3 Antibody</u>	CSB-PA882162LA01HU	Human, Mouse	ELISA, WB, IF
<u>FLT1 Antibody</u>	CSB-PA07269A0Rb	Human	ELISA, IHC, IF
<u>FLT3 Antibody</u>	CSB-PA008733LA01HU	Human	ELISA, WB
<u>FLT3LG Antibody</u>	CSB-PA008734LA01HU	Human	ELISA, IHC, IF
<u>FLT4 Antibody</u>	CSB-PA07539A0Rb	Human	ELISA, IHC, IF
<u>FMN1 Antibody</u>	CSB-PA715012LA01HU	Human	ELISA, WB, IHC, IF
<u>FMNL1 Antibody</u>	CSB-PA008743LA01HU	Human	ELISA, IHC
<u>FMNL3 Antibody</u>	CSB-PA815554LA01HU	Human	ELISA, IHC, IF
<u>FMO3 Antibody</u>	CSB-PA008749LA01HU	Human	ELISA, IHC
<u>FMO4 Antibody</u>	CSB-PA008750LA01HU	Human	ELISA, IHC, IF
<u>FMOD Antibody</u>	CSB-PA008755LA01HU	Human, Rat	ELISA, WB, IHC
<u>FN1 Antibody</u>	CSB-PA13187A0Rb	Human	ELISA
<u>FNDC11 Antibody</u>	CSB-PA883611LA01HU	Human	ELISA, IHC
<u>FNDC3B Antibody</u>	CSB-PA709599LA01HU	Human	ELISA, IHC, IF
<u>FNDC4 Antibody</u>	CSB-PA872504LA01HU	Human	ELISA
<u>FNDC7 Antibody</u>	CSB-PA736197LA01HU	Human	ELISA
<u>FNIP1 Antibody</u>	CSB-PA819483LA01HU	Human	ELISA, IHC, IF
<u>FNTB Antibody</u>	CSB-PA008780LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>FOLH1 Antibody</u>	CSB-PA11729A0Rb	Human	ELISA, WB, IHC
<u>FOLR1 Antibody</u>	CSB-PA008784LA01HU	Human	ELISA, IHC, IF
<u>FOLR2 Antibody</u>	CSB-PA008786LA01HU	Human	ELISA, IHC, IF
<u>FOSB Antibody</u>	CSB-PA008791LA01HU	Human	ELISA, IHC
<u>FOSL2 Antibody</u>	CSB-PA008793HA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>FOXA1 Antibody</u>	CSB-PA008794LA01HU	Human	ELISA, IHC
<u>FOXA3 Antibody</u>	CSB-PA008796LA01HU	Human	ELISA, IHC
<u>FOXC1 Antibody</u>	CSB-PA618640LA01HU	Human	ELISA, WB, IF
<u>FOXC2 Antibody</u>	CSB-PA857880LA01HU	Human	ELISA, IHC
<u>FOXD3 Antibody</u>	CSB-PA890729LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>FOXE1 Antibody</u>	CSB-PA008811LA01HU	Human, Mouse	ELISA, WB, IF
<u>FOXE3 Antibody</u>	CSB-PA618775LA01HU	Human, Rat, Mouse	ELISA, WB, IF
<u>FOXF2 Antibody</u>	CSB-PA623661LA01HU	Human	ELISA
<u>FOXG1 Antibody</u>	CSB-PA008815HA01HU	Human	ELISA, IHC, IF
<u>FOXH1 Antibody</u>	CSB-PA008817LA01HU	Human	ELISA, WB
<u>FOXI1 Antibody</u>	CSB-PA622520LA01HU	Human	ELISA, WB, IF
<u>FOXK1 Antibody</u>	CSB-PA008824LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>FOXL2 Antibody</u>	CSB-PA008827LA01HU	Human	ELISA, WB
<u>FOXM1 Antibody</u>	CSB-PA008828HA01HU	Human	ELISA, IHC
<u>FOXN1 Antibody</u>	CSB-PA008829LA01HU	Human, Mouse	ELISA, WB, IHC
<u>FOXN4 Antibody</u>	CSB-PA850403LA01HU	Human	ELISA, IF
<u>FOXO1 Antibody</u>	CSB-PA16019A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>FOXO4 Antibody</u>	CSB-PA008839OA01HU	Human	ELISA, IHC, IF
<u>FOXO6 Antibody</u>	CSB-PA008840LA01HU	Human	ELISA, IF
<u>FOXO6 Antibody</u>	CSB-PA008840OA01HU	Human	ELISA, IHC, IF
<u>FOXP1 Antibody</u>	CSB-PA008841EA01HU	Human	ELISA, WB, IHC, IF, ChIP
<u>FOXP1 Antibody</u>	CSB-PA008841LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>FOXP3 Antibody</u>	CSB-PA15229A0Rb	Human, Mouse	ELISA, WB, IHC
<u>FOXQ1 Antibody</u>	CSB-PA874860LA01HU	Human	ELISA, WB, IHC, IF
<u>FOXR2 Antibody</u>	CSB-PA757660LA01HU	Human	ELISA, WB, IHC
<u>FOXRED2 Antibody</u>	CSB-PA816892HA01HU	Human	ELISA, IHC
<u>FPGS Antibody</u>	CSB-PA23719A0Rb	Human	ELISA, WB
<u>FPR1 Antibody</u>	CSB-PA07129A0Rb	Human	ELISA, WB, IHC, IF
<u>FPR2 Antibody</u>	CSB-PA008855LA01HU	Human	ELISA, WB, IF
<u>FPR3 Antibody</u>	CSB-PA008856LA01HU	Human	ELISA, IHC, IF
<u>fre Antibody</u>	CSB-PA18537A0Rb	Escherichia coli	ELISA, WB
<u>FREM1 Antibody</u>	CSB-PA686013LA01HU	Human	ELISA, IHC
<u>FREM2 Antibody</u>	CSB-PA722558LA01HU	Human	ELISA, IHC, IF
<u>FRMD5 Antibody</u>	CSB-PA773602LA01HU	Human, Mouse	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>FRMD7 Antibody</u>	CSB-PA744425LA01HU	Human	ELISA, WB, IHC
<u>FRS2 Antibody</u>	CSB-PA009004LA01HU	Human	ELISA, IHC
<u>FRS3 Antibody</u>	CSB-PA009005LA01HU	Human	ELISA, IHC, IF
<u>FRYL Antibody</u>	CSB-PA009007LA01HU	Human	ELISA, IHC, IF
<u>FSCB Antibody</u>	CSB-PA708110LA01HU	Human, Mouse	ELISA, WB, IHC
<u>FSCN1 Antibody</u>	CSB-PA619084LA01HU	Human	ELISA, WB, IHC, IF
<u>FSCN2 Antibody</u>	CSB-PA009011LA01HU	Human	ELISA
<u>FSD1 Antibody</u>	CSB-PA861136LA01HU	Human	ELISA, IHC
<u>FSHB Antibody</u>	CSB-PA12335A0Rb	Human	ELISA
<u>FSIP1 Antibody</u>	CSB-PA843306LA01HU	Human	ELISA, IHC
<u>FSIP2 Antibody</u>	CSB-PA688567LA01HU	Human	ELISA, WB, IF
<u>FSTL1 Antibody</u>	CSB-PA613266LA01HU	Human	ELISA, IHC, IF
<u>FTH1 Antibody</u>	CSB-PA02084A0Rb	Human	ELISA
<u>FTH1 Antibody</u>	CSB-PA02089A0Rb	Human	ELISA
<u>FTL Antibody</u>	CSB-PA11407A0Rb	Human	ELISA, IHC
<u>FTL Antibody</u>	CSB-PA11409A0Rb	Human	ELISA, IHC, IF
<u>FTL3 Antibody</u>	CSB-PA21589A0Rb	Soybean	ELISA, WB
<u>FTSJ1 Antibody</u>	CSB-PA009055LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ftsZ Antibody</u>	CSB-PA343512LA01WAP	Wolbachia pipientis	ELISA, WB
<u>ftsZ Antibody</u>	CSB-PA359270HA01EGX	Escherichia coli O6:H1	ELISA, WB
<u>FUCA1 Antibody</u>	CSB-PA009062LA01HU	Human	ELISA, IHC
<u>Fuca1 Antibody</u>	CSB-PA858759LA01MO	Mouse, Rat	ELISA, WB
<u>FUCA2 Antibody</u>	CSB-PA12529A0Rb	Human	ELISA, WB, IHC
<u>FUK Antibody</u>	CSB-PA847604LA01HU	Human	ELISA, IHC, IF
<u>FUNDC1 Antibody</u>	CSB-PA816886LA01HU	Human, Mouse	ELISA, WB, IHC
<u>FUOM Antibody</u>	CSB-PA002922LA01HU	Human	ELISA, IHC, IF
<u>FURIN Antibody</u>	CSB-PA009068LA01HU	Human	ELISA, IHC
<u>FUS Antibody</u>	CSB-PA02704A0Rb	Human	ELISA, WB, IHC, IF
<u>FUT3 Antibody</u>	CSB-PA009077LA01HU	Human, Rat	ELISA, WB, IHC
<u>FUT8 Antibody</u>	CSB-PA887159LA01HU	Human, Mouse	ELISA, WB, IHC
<u>FUZ Antibody</u>	CSB-PA883593LA01HU	Human	ELISA, IHC, IF
<u>Fxn Antibody</u>	CSB-PA009086LA01MO	Mouse, Human	ELISA, WB
<u>FXR1 Antibody</u>	CSB-PA009087LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>FXYD1 Antibody</u>	CSB-PA009089LA01HU	Human	ELISA
<u>FXYD6 Antibody</u>	CSB-PA863943LA01HU	Human	ELISA, IHC
<u>FXYD7 Antibody</u>	CSB-PA009097LA01HU	Human	ELISA, IHC, IF
<u>FYCO1 Antibody</u>	CSB-PA866262LA01HU	Human, Mouse	ELISA, WB, IHC
<u>FYN Antibody</u>	CSB-PA009101LA01HU	Human	ELISA, IHC, IF
<u>FYTTD1 Antibody</u>	CSB-PA842744LA01HU	Human	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>FZD1 Antibody</u>	CSB-PA891570LA01HU	Human	ELISA, IHC, IF
<u>FZD10 Antibody</u>	CSB-PA891782LA01HU	Human	ELISA, WB, IHC, IF
<u>FZD4 Antibody</u>	CSB-PA890755LA01HU	Human	ELISA, WB, IHC, IF
<u>FZD7 Antibody</u>	CSB-PA009110LA01HU	Human	ELISA, WB, IHC
<u>FZD8 Antibody</u>	CSB-PA880970LA01HU	Human, Mouse	ELISA, WB, IHC
<u>FZD9 Antibody</u>	CSB-PA009112LA01HU	Human, Mouse	ELISA, WB, IHC
<u>FZR1 Antibody</u>	CSB-PA892473LA01HU	Human	ELISA, WB, IHC
<u>G Antibody</u>	CSB-PA14899A0Rb	Rabies virus	ELISA, WB
<u>G Antibody</u>	CSB-PA365937LA01HPO	Human respiRatory syncytial virus A	ELISA
<u>G Antibody</u>	CSB-PA389005LA01RIE	Rabies virus	ELISA
<u>G0S2 Antibody</u>	CSB-PA009114EA01HU	Human, Mouse	ELISA, WB, IHC
<u>G2E3 Antibody</u>	CSB-PA741972LA01HU	Human	ELISA, IF
<u>G6PC Antibody</u>	CSB-PA009118LA01HU	Human	ELISA, IHC
<u>G6PC2 Antibody</u>	CSB-PA873624LA01HU	Human	ELISA, IHC, IF
<u>G6PD Antibody</u>	CSB-PA009121LA01HU	Human	ELISA, IHC, IF
<u>GAA Antibody</u>	CSB-PA009125LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>GAB1 Antibody</u>	CSB-PA614416LA01HU	Human	ELISA, IHC, IF
<u>GAB2 Antibody</u>	CSB-PA892175LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>GABARAP Antibody</u>	CSB-PA05484A0Rb	Human	ELISA, IF
<u>GABARAPL1 Antibody</u>	CSB-PA861986LA01HU	Human	ELISA, IHC, IF
<u>GABARAPL2 Antibody</u>	CSB-PA03484A0Rb	Human, Rat	ELISA, WB, IHC
<u>GABARAPL2 Antibody</u>	CSB-PA03485A0Rb	Human	ELISA, WB
<u>GABPA Antibody</u>	CSB-PA22769A0Rb	Human	ELISA, IHC, IF
<u>GABPB2 Antibody</u>	CSB-PA844704LA01HU	Human	ELISA, WB, IHC, IF
<u>GABRA2 Antibody</u>	CSB-PA009141LA01HU	Human	ELISA, IHC, IF
<u>GABRA4 Antibody</u>	CSB-PA009143LA01HU	Human	ELISA, WB, IHC, IF
<u>GABRA6 Antibody</u>	CSB-PA623078LA01HU	Human	ELISA, IHC
<u>GABRB2 Antibody</u>	CSB-PA009147LA01HU	Human	ELISA, IHC, IF
<u>GABRB3 Antibody</u>	CSB-PA009148LA01HU	Human	ELISA, WB, IHC, IF
<u>GABRG2 Antibody</u>	CSB-PA009152LA01HU	Human	ELISA, IF
<u>GABRQ Antibody</u>	CSB-PA891564LA01HU	Human	ELISA, IHC, IF
<u>Gad2 Antibody</u>	CSB-PA009160HA01MO	Mouse	ELISA
<u>GAD2 Antibody</u>	CSB-PA11159A0Rb	Human, Mouse	ELISA, WB, IHC
<u>gadA Antibody</u>	CSB-PA302526HA01ENV	Escherichia coli	ELISA



Product Name	Code	Species Reactivity	Applications
<u>GADD45A Antibody</u>	CSB-PA009161HA01HU	Human	ELISA
<u>GADD45B Antibody</u>	CSB-PA009163LA01HU	Human	ELISA, IHC, IF
<u>GADD45G Antibody</u>	CSB-PA009164HA01HU	Human	ELISA, IHC
<u>gag Antibody</u>	CSB-PA15019A0Rb	Human immunodeficiency virus type 1 group M subtype B	ELISA
<u>gag Antibody</u>	CSB-PA356017DA01MHK	Moloney murine leukemia virus	ELISA
<u>gag Antibody</u>	CSB-PA356017EA01MHK	Moloney murine leukemia virus	ELISA
<u>GAGE5 Antibody</u>	CSB-PA009186EA01HU	Human	ELISA, IHC
<u>gag-pol Antibody</u>	CSB-PA15387A0Rb	Human immunodeficiency virus type 1 group M subtype B	ELISA
<u>GAK Antibody</u>	CSB-PA009190LA01HU	Human	ELISA, IHC, IF
<u>GALE Antibody</u>	CSB-PA613499HA01HU	Human	ELISA, WB, IHC, IF
<u>GALK1 Antibody</u>	CSB-PA009198LA01HU	Human	ELISA, WB, IF
<u>GALK2 Antibody</u>	CSB-PA009199LA01HU	Human	ELISA, IHC
<u>GALM Antibody</u>	CSB-PA850268HA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>GALNS Antibody</u>	CSB-PA009202EA01HU	Human	ELISA
<u>GALNT11 Antibody</u>	CSB-PA847687LA01HU	Human	ELISA, IHC, IF
<u>GALNT12 Antibody</u>	CSB-PA814224LA01HU	Human	ELISA, IHC
<u>GALNT14 Antibody</u>	CSB-PA836212LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>GALNT15 Antibody</u>	CSB-PA822698LA01HU	Human, Mouse	ELISA, WB, IHC
<u>GALNT16 Antibody</u>	CSB-PA850796LA01HU	Human, Mouse	ELISA, WB
<u>GALNT17 Antibody</u>	CSB-PA743592LA01HU	Human	ELISA, IHC
<u>GALNT7 Antibody</u>	CSB-PA773016LA01HU	Human	ELISA, WB, IHC
<u>GALNT9 Antibody</u>	CSB-PA862072LA01HU	Human	ELISA, WB
<u>GALR1 Antibody</u>	CSB-PA009223LA01HU	Human	ELISA, WB
<u>GALR2 Antibody</u>	CSB-PA009224LA01HU	Human	ELISA, IF
<u>GALT Antibody</u>	CSB-PA009226LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF, IP
<u>Gamma-glutamyl hydrolase Antibody</u>	CSB-PA009389LA01GGV	Glycine max	ELISA, WB
<u>Gamt Antibody</u>	CSB-PA009227LA01RA	Rat	ELISA, WB
<u>GANC Antibody</u>	CSB-PA837451LA01HU	Human	ELISA, WB, IHC
<u>GAPDH Antibody</u>	CSB-PA00025A0Rb	Human	ELISA, WB, IHC, IF
<u>GAPDH Antibody</u>	CSB-PA009232HA01CH	Chicken	ELISA
<u>GAPDH Antibody</u>	CSB-PA00029A0Rb	Human, Mouse	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>GAPDH Antibody</u>	CSB-PA009232HA01DXU	Cricetulus griseus	ELISA, WB
<u>GAPDH Monoclonal Antibody</u>	CSB-MA000071M0m		ELISA, WB, IHC, IP, IF, FC
<u>Gapdh1 Antibody</u>	CSB-PA362145LA01DLU	Drosophila melanogaster	ELISA
<u>GAR1 Antibody</u>	CSB-PA878917LA01HU	Human	ELISA, IHC
<u>GAREM1 Antibody</u>	CSB-PA863993LA01HU	Human	ELISA, IHC, IF
<u>GARS Antibody</u>	CSB-PA009262EA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF, IP
<u>GAS1 Antibody</u>	CSB-PA009264LA01HU	Human	ELISA, WB, IHC, IF
<u>GAS2 Antibody</u>	CSB-PA009265LA01HU	Human, Rat	ELISA, WB, IHC
<u>GAS6 Antibody</u>	CSB-PA619872LA01HU	Human	ELISA, WB
<u>GAS7 Antibody</u>	CSB-PA009271LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>GATA1 Antibody</u>	CSB-PA009274LA01HU	Human	ELISA, WB, IHC
<u>gata1-a Antibody</u>	CSB-PA326431LA01XBE	Xenopus laevis	ELISA
<u>GATA5 Antibody</u>	CSB-PA883616HA01HU	Human	ELISA, IHC, IF
<u>GATA6 Antibody</u>	CSB-PA846113LA01HU	Human	ELISA, IHC
<u>GATAD2A Antibody</u>	CSB-PA803156LA01HU	Human, Rat	ELISA, WB, IHC
<u>GATAD2B Antibody</u>	CSB-PA823905LA01HU	Human	ELISA, IHC
<u>GATB Antibody</u>	CSB-PA017792LA01HU	Human	ELISA, WB, IHC, IF
<u>GATD3A Antibody</u>	CSB-PA01814A0Rb	Human	ELISA, WB, IHC
<u>gB Antibody</u>	CSB-PA362479LA01HJX	Human herpesvirus 2	ELISA
<u>Gba Antibody</u>	CSB-PA009289HA01MO	Mouse	ELISA, WB
<u>GBA Antibody</u>	CSB-PA009289LA01HU	Human	ELISA, WB, IHC
<u>GBA2 Antibody</u>	CSB-PA881025LA01HU	Human	ELISA, IHC
<u>GBP2 Antibody</u>	CSB-PA009298LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>GBP4 Antibody</u>	CSB-PA853481LA01HU	Human	ELISA, WB, IHC, IF
<u>GBP5 Antibody</u>	CSB-PA847203LA01HU	Human, Mouse	ELISA, WB, IHC
<u>GBP6 Antibody</u>	CSB-PA751154LA01HU	Human	ELISA, IHC
<u>GBX2 Antibody</u>	CSB-PA009305LA01HU	Human	ELISA, IHC
<u>GC Antibody</u>	CSB-PA07569A0Rb	Human	ELISA, WB, IHC
<u>GCC1 Antibody</u>	CSB-PA846598LA01HU	Human	ELISA, WB, IHC
<u>GCDH Antibody</u>	CSB-PA849798LA01HU	Human, Mouse	ELISA, WB, IHC, IP
<u>GCGR Antibody</u>	CSB-PA009316LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>GCHFR Antibody</u>	CSB-PA009318LA01HU	Human	ELISA, IHC, IF
<u>GCKR Antibody</u>	CSB-PA622767LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>GCLC Antibody</u>	CSB-PA009321LA01HU	Human	ELISA, IHC
<u>GCLM Antibody</u>	CSB-PA009322LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>GCM2 Antibody</u>	CSB-PA009324LA01HU	Human, Rat	ELISA, WB, IHC
<u>GCNA Antibody</u>	CSB-PA836285LA01HU	Human	ELISA, WB, IHC
<u>GCNT2 Antibody</u>	CSB-PA009328LA01HU	Human	ELISA
<u>GCNT3 Antibody</u>	CSB-PA009329LA01HU	Human	ELISA, WB, IHC
<u>GCSAM Antibody</u>	CSB-PA854052LA01HU	Human	ELISA, IHC
<u>GCSAML Antibody</u>	CSB-PA689332LA01HU	Human	ELISA, IHC
<u>GCY1 Antibody</u>	CSB-PA321197HA01SVG	Saccharomyces cerevisiae	ELISA
<u>GCY1 Antibody</u>	CSB-PA321197YA01SVG	Saccharomyces cerevisiae	ELISA
<u>gD Antibody</u>	CSB-PA18929A0Rb	Human herpesvirus 1	ELISA
<u>GDAP1L1 Antibody</u>	CSB-PA009339LA01HU	Human	ELISA, WB, IHC
<u>GDAP2 Antibody</u>	CSB-PA868350LA01HU	Human	ELISA, WB
<u>GDF10 Antibody</u>	CSB-PA009343LA01HU	Human, Rat	ELISA, WB, IHC
<u>GDF11 Antibody</u>	CSB-PA009344LA01PI	Pig, Mouse	ELISA, WB
<u>GDF15 Antibody</u>	CSB-PA10029A0Rb	Human	ELISA, WB, IHC
<u>GDF2 Antibody</u>	CSB-PA892325LA01HU	Human	ELISA, IHC
<u>GDF3 Antibody</u>	CSB-PA868295NA01HU	Human	ELISA, IF
<u>Gdf5 Antibody</u>	CSB-PA009349LA01MO	Mouse	ELISA, IHC
<u>GDF7 Antibody</u>	CSB-PA768769LA01HU	Human	ELISA, WB, IHC
<u>GDF9 Antibody</u>	CSB-PA009352EA01HU	Human	ELISA, IHC
<u>GDI1 Antibody</u>	CSB-PA01229A0Rb	Human	ELISA, WB, IHC, IF
<u>GDI2 Antibody</u>	CSB-PA009354LA01HU	Human, Mouse	ELISA, WB, IHC
<u>GDI2 Antibody</u>	CSB-PA00945A0Rb	Human	ELISA, WB
<u>GDPD2 Antibody</u>	CSB-PA864028LA01HU	Human	ELISA, WB, IHC, IF
<u>GDV1 Antibody</u>	CSB-PA28187A0Rb	Plasmodium relictum	ELISA, WB
<u>geIE Antibody</u>	CSB-PA769310LA01ELW	Enterococcus faecalis	ELISA, WB
<u>GEMIN2 Antibody</u>	CSB-PA021329LA01HU	Human	ELISA, WB, IF, IP
<u>GEMIN4 Antibody</u>	CSB-PA009363LA01HU	Human	ELISA, IHC, IF
<u>GEMIN6 Antibody</u>	CSB-PA823902LA01HU	Human	ELISA, IHC
<u>GEMIN7 Antibody</u>	CSB-PA862039EA01HU	Human	ELISA, WB
<u>GEMIN8 Antibody</u>	CSB-PA873675LA01HU	Human, Mouse	ELISA, WB, IHC
<u>Genome polyprotein Antibody</u>	CSB-PA18549A0Rb	Dengue virus type 1	ELISA
<u>Genome polyprotein Antibody</u>	CSB-PA333180LA01HFD	Hepatitis C virus genotype 1a	ELISA
<u>Genome polyprotein Antibody</u>	CSB-PA362073LA01HQD	Human rhinovirus A serotype 89	ELISA
<u>Genome polyprotein Antibody</u>	CSB-PA714772LA01HYH	Human enterovirus 71	ELISA
<u>GFAP Antibody</u>	CSB-PA009369LA01HU	Human	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>Gfap Antibody</u>	CSB-PA009369LA01RA	Rat, Mouse	ELISA, WB
<u>gfap Antibody</u>	CSB-PA009369LA01DIL	Zebrafish, Human	ELISA, WB
<u>GFAP Monoclonal Antibody</u>	CSB-MA009369A0m	Human, Mouse, Rat	ELISA, WB, IHC, IF, FC
<u>Gfer Antibody</u>	CSB-PA733907HA01RA	Rat	ELISA
<u>GFOD2 Antibody</u>	CSB-PA665476LA01HU	Human	ELISA, WB, IHC
<u>GFP Monoclonal Antibody</u>	CSB-MA000051M0m	All	ELISA, WB, IP
<u>GFRA1 Antibody</u>	CSB-PA009379LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>GFRA2 Antibody</u>	CSB-PA009380LA01HU	Human, Mouse	ELISA, WB, IF
<u>GFRA3 Antibody</u>	CSB-PA009381LA01HU	Human	ELISA, IHC, IF
<u>GFRA4 Antibody</u>	CSB-PA009382LA01HU	Human	ELISA, IHC, IF
<u>GFRA5 Antibody</u>	CSB-PA009383LA01HU	Human	ELISA, IHC, IF
<u>GFRA6 Antibody</u>	CSB-PA009384LA01HU	Human	ELISA, IHC, IF
<u>GFRA7 Antibody</u>	CSB-PA009385LA01HU	Human	ELISA, IHC, IF
<u>GFRA8 Antibody</u>	CSB-PA009386LA01HU	Human	ELISA, IHC, IF
<u>GFRA9 Antibody</u>	CSB-PA009387LA01HU	Human	ELISA, IHC, IF
<u>GFRA10 Antibody</u>	CSB-PA009388LA01HU	Human	ELISA, IHC, IF
<u>GFRA11 Antibody</u>	CSB-PA009389LA01HU	Human	ELISA, IHC, IF
<u>GFRA12 Antibody</u>	CSB-PA009390LA01HU	Human	ELISA, IHC, IF
<u>GFRA13 Antibody</u>	CSB-PA009391LA01HU	Human	ELISA, IHC, IF
<u>GFRA14 Antibody</u>	CSB-PA009392LA01HU	Human	ELISA, IHC, IF
<u>GFRA15 Antibody</u>	CSB-PA009393HA01HU	Human	ELISA, IHC
<u>GG Antibody</u>	CSB-PA361844HA01HWY	Human herpesvirus 1	ELISA
<u>GGA1 Antibody</u>	CSB-PA890911LA01HU	Human	ELISA, IHC, IF
<u>GGA2 Antibody</u>	CSB-PA890912LA01HU	Human	ELISA, IHC, IF
<u>GGA3 Antibody</u>	CSB-PA882151LA01HU	Human	ELISA, IHC
<u>GGCT Antibody</u>	CSB-PA009387LA01HU	Human	ELISA, IHC
<u>GGCX Antibody</u>	CSB-PA009388LA01HU	Human	ELISA, IHC
<u>GGH Antibody</u>	CSB-PA009389LA01HU	Human	ELISA, WB
<u>GGPS1 Antibody</u>	CSB-PA009393HA01HU	Human	ELISA, IHC
<u>GGT1 Antibody</u>	CSB-PA10329A0Rb	Human	ELISA, IHC, IF
<u>GGTLC1 Antibody</u>	CSB-PA009404LA01HU	Human	ELISA, IHC
<u>GGTLC2 Antibody</u>	CSB-PA009405LA01HU	Human, Mouse	ELISA, WB, IHC
<u>gH Antibody</u>	CSB-PA319015LA01HWV	Human cytomegalovirus	ELISA
<u>GH1 Antibody</u>	CSB-PA009407LA01HU	Human	ELISA, IHC
<u>GH2 Antibody</u>	CSB-PA009408LA01HU	Human, Mouse	ELISA, WB, IHC
<u>GHITM Antibody</u>	CSB-PA872487LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>GHR Antibody</u>	CSB-PA009411LA01HU	Human	ELISA, IHC
<u>ghrA Antibody</u>	CSB-PA482379HA01ENM	Escherichia coli	ELISA, WB
<u>GHRH Antibody</u>	CSB-PA009412EA01SH	Sheep	ELISA
<u>GHRHR Antibody</u>	CSB-PA009413LA01HU	Human	ELISA, IHC, IF
<u>GHSR Antibody</u>	CSB-PA842164LA01HU	Human	ELISA, IHC
<u>GID8 Antibody</u>	CSB-PA882124LA01HU	Human	ELISA, IHC
<u>GIMAP4 Antibody</u>	CSB-PA885732LA01HU	Human	ELISA, IF
<u>GIMAP7 Antibody</u>	CSB-PA818820LA01HU	Human	ELISA, WB
<u>GINS1 Antibody</u>	CSB-PA617926LA01HU	Human	ELISA, IHC
<u>GINS2 Antibody</u>	CSB-PA896481LA01HU	Human	ELISA, IHC
<u>GINS4 Antibody</u>	CSB-PA871563LA01HU	Human	ELISA, IHC
<u>GIP Antibody</u>	CSB-PA009434PA01HU	Human	ELISA, IHC
<u>GIPC1 Antibody</u>	CSB-PA009435LA01HU	Human	ELISA, IHC, IF
<u>GIT1 Antibody</u>	CSB-PA009439LA01HU	Human, Mouse	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>GIT2 Antibody</u>	CSB-PA613489LA01HU	Human	ELISA, IHC
<u>GJA4 Antibody</u>	CSB-PA009447LA01HU	Human	ELISA, IHC, IF
<u>GJA5 Antibody</u>	CSB-PA009448LA01HU	Human	ELISA, IHC, IF
<u>GJA8 Antibody</u>	CSB-PA009449LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>GJB1 Antibody</u>	CSB-PA009451LA01HU	Human	ELISA, IHC, IF
<u>GJB2 Antibody</u>	CSB-PA009452LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>GJB3 Antibody</u>	CSB-PA009453LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>GJB4 Antibody</u>	CSB-PA865136LA01HU	Human	ELISA
<u>GJB5 Antibody</u>	CSB-PA009455DA01HU	Human	ELISA, IHC, IF
<u>GJB6 Antibody</u>	CSB-PA009456LA01HU	Human	ELISA, IHC, IF
<u>GJB7 Antibody</u>	CSB-PA757637LA01HU	Human	ELISA, IF
<u>GJC1 Antibody</u>	CSB-PA009458EA01HU	Human	ELISA, WB, IF
<u>GJD2 Antibody</u>	CSB-PA891536LA01HU	Human	ELISA, WB, IHC, IF
<u>GK Antibody</u>	CSB-PA300097LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>GK2 Antibody</u>	CSB-PA009466LA01HU	Human	ELISA, WB, IHC
<u>GK3P Antibody</u>	CSB-PA614527LA01HU	Human	ELISA, WB, IHC
<u>GK5 Antibody</u>	CSB-PA744414LA01HU	Human	ELISA, WB
<u>GKAP1 Antibody</u>	CSB-PA719617LA01HU	Human	ELISA, IHC
<u>GKN1 Antibody</u>	CSB-PA10589A0Rb	Human	ELISA
<u>Gkn2 Antibody</u>	CSB-PA866591LA01MO	Mouse	ELISA
<u>gL Antibody</u>	CSB-PA515269LA01HWW	Human cytomegalovirus	ELISA
<u>GLA Antibody</u>	CSB-PA009474LA01HU	Human	ELISA, IHC
<u>glaB1 Antibody</u>	CSB-PA406363HA01BOP	Bacteroides vulgatus	ELISA
<u>GLB1 Antibody</u>	CSB-PA009476LA01HU	Human	ELISA, WB, IHC, IF
<u>GLB1L2 Antibody</u>	CSB-PA811610LA01HU	Human	ELISA, IHC
<u>glcB Antibody</u>	CSB-PA327661LA01ENV	Escherichia coli	ELISA
<u>GLCE Antibody</u>	CSB-PA009492LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>GLDC Antibody</u>	CSB-PA328400LA01HU	Human	ELISA, IHC
<u>GLE1 Antibody</u>	CSB-PA700651EA01HU	Human	ELISA, WB, IHC, IF, IP
<u>GLI1 Antibody</u>	CSB-PA12989A0Rb	Human	ELISA, IHC, IF
<u>GLI2 Antibody</u>	CSB-PA009500OA01HU	Human	ELISA, WB, IHC
<u>GLI4 Antibody</u>	CSB-PA009502LA01HU	Human	ELISA, IHC
<u>GLIPR1 Antibody</u>	CSB-PA009503LA01HU	Human	ELISA, IHC, IF
<u>GLIPR2 Antibody</u>	CSB-PA880974LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>GLIS2 Antibody</u>	CSB-PA861185LA01HU	Human	ELISA, IHC, IF
<u>GLIS3 Antibody</u>	CSB-PA839863LA01HU	Human	ELISA, IHC, IF
<u>GLO1 Antibody</u>	CSB-PA862860LA01DOA	Arabidopsis thaliana	ELISA, WB
<u>GLP1R Antibody</u>	CSB-PA009514LA01HU	Human	ELISA, WB, IHC, IF
<u>GLP1R Antibody</u>	CSB-PA009514YA01HU	Human, Mouse	ELISA, WB, IF
<u>Glp1r Antibody</u>	CSB-PA009514HA01MO	Mouse	ELISA, WB
<u>glpD Antibody</u>	CSB-PA320174HA01ENV	Escherichia coli	ELISA
<u>glpK Antibody</u>	CSB-PA18407A0Rb	Escherichia coli	ELISA, WB
<u>GLRA1 Antibody</u>	CSB-PA009516LA01HU	Human	ELISA, WB, IHC, IF
<u>GLRB Antibody</u>	CSB-PA009520LA01HU	Human	ELISA, IF
<u>GLRX Antibody</u>	CSB-PA009521EA01HU	Human	ELISA, IHC, IF
<u>GLRX5 Antibody</u>	CSB-PA769767LA01HU	Human	ELISA, IHC
<u>GLS Antibody</u>	CSB-PA009528DA01HU	Human, Mouse	ELISA, WB, IHC
<u>GLS Antibody</u>	CSB-PA009528LA01HU	Human	ELISA
<u>GLS2 Antibody</u>	CSB-PA887058LA01HU	Human	ELISA, IF
<u>glsA1 Antibody</u>	CSB-PA302868EA01ENV	Escherichia coli	ELISA, WB
<u>GLT1D1 Antibody</u>	CSB-PA842719LA01HU	Human	ELISA, IHC
<u>GLT8D1 Antibody</u>	CSB-PA715010LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>GLT8D2 Antibody</u>	CSB-PA009536LA01HU	Human, Mouse	ELISA, WB
<u>GLTP Antibody</u>	CSB-PA009539LA01BO	Bovine	ELISA
<u>GLU1 Antibody</u>	CSB-PA24649A0Rb	Arabidopsis thaliana	ELISA, WB
<u>gluD Antibody</u>	CSB-PA335284HA01DUN	Clostridium difficile	ELISA
<u>GLUL Antibody</u>	CSB-PA009553LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>GM2A Antibody</u>	CSB-PA02375A0Rb	Human, Rat	ELISA, WB
<u>GM2A Antibody</u>	CSB-PA02379A0Rb	Human, Mouse, Rat	ELISA, WB, IHC, IP
<u>GMDS Antibody</u>	CSB-PA009569LA01HU	Human	ELISA, IHC, IF
<u>GMEB1 Antibody</u>	CSB-PA896751LA01HU	Human	ELISA
<u>GMFG Antibody</u>	CSB-PA009574LA01HU	Human	ELISA
<u>GMIP Antibody</u>	CSB-PA878952LA01HU	Human	ELISA, IHC, IF
<u>GMNN Antibody</u>	CSB-PA009577LA01HU	Human	ELISA, IHC, IF
<u>GMPPB Antibody</u>	CSB-PA896740LA01HU	Human	ELISA, IHC
<u>GMPR Antibody</u>	CSB-PA339756HA01HU	Human	ELISA, IHC
<u>GMPR2 Antibody</u>	CSB-PA873740LA01HU	Human	ELISA
<u>GNA11 Antibody</u>	CSB-PA009583LA01HU	Human	ELISA, IHC, IF
<u>GNA12 Antibody</u>	CSB-PA009584LA01HU	Human	ELISA, WB, IHC, IF
<u>GNA14 Antibody</u>	CSB-PA009586LA01HU	Human, Mouse	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>GNAI3 Antibody</u>	CSB-PA009591LA01HU	Human	ELISA
<u>GNAL Antibody</u>	CSB-PA009592LA01HU	Human	ELISA, WB, IHC
<u>Gnao1 Antibody</u>	CSB-PA009593LA01MO	Mouse, Rat	ELISA, WB
<u>GNAQ Antibody</u>	CSB-PA009594LA01HU	Human	ELISA, IHC, IF
<u>GNAS Antibody</u>	CSB-PA352389LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>GNAS Antibody</u>	CSB-PA307619LA01HU	Human	ELISA, WB, IHC, IF
<u>GNAT1 Antibody</u>	CSB-PA009598LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>GNAT2 Antibody</u>	CSB-PA009599LA01HU	Human	ELISA, IF
<u>GNAT3 Antibody</u>	CSB-PA009600LA01HU	Human	ELISA, IHC, IF
<u>GNAZ Antibody</u>	CSB-PA009601LA01HU	Human	ELISA, IHC, IF
<u>GNB1 Antibody</u>	CSB-PA009602LA01HU	Human, Mouse, Rat	ELISA, WB, IF
<u>GNB4 Antibody</u>	CSB-PA862056LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>GNB5 Antibody</u>	CSB-PA009608LA01HU	Human	ELISA, IHC, IF
<u>GNG10 Antibody</u>	CSB-PA009610LA01HU	Human	ELISA, IHC, IF
<u>GNG11 Antibody</u>	CSB-PA009611LA01HU	Human	ELISA, IHC, IF
<u>GNG12 Antibody</u>	CSB-PA883362OA01HU	Human	ELISA, IHC, IF
<u>GNG2 Antibody</u>	CSB-PA009614OA01HU	Human	ELISA, IHC, IF
<u>GNG4 Antibody</u>	CSB-PA009616LA01HU	Human	ELISA, IHC, IF
<u>GNG5 Antibody</u>	CSB-PA009617LA01HU	Human	ELISA, IHC, IF
<u>GNG7 Antibody</u>	CSB-PA009619LA01HU	Human	ELISA, IF
<u>GNGT2 Antibody</u>	CSB-PA009622LA01HU	Human	ELISA, IF
<u>GNL1 Antibody</u>	CSB-PA009623LA01HU	Human	ELISA, WB, IF
<u>GNL2 Antibody</u>	CSB-PA618806LA01HU	Human	ELISA
<u>GNL3L Antibody</u>	CSB-PA865150LA01HU	Human	ELISA
<u>GNMT Antibody</u>	CSB-PA009628LA01HU	Human	ELISA, IHC, IF
<u>GNPAT Antibody</u>	CSB-PA009629LA01HU	Human	ELISA, IHC
<u>GNPDA2 Antibody</u>	CSB-PA855058LA01HU	Human	ELISA, IF
<u>GNPNAT1 Antibody</u>	CSB-PA856926LA01HU	Human, Mouse	ELISA, WB, IHC
<u>GNRH2 Antibody</u>	CSB-PA009636LA01HU	Human	ELISA, IHC
<u>GOLGA3 Antibody</u>	CSB-PA600096LA01HU	Human	ELISA, IHC, IF
<u>GOLGB1 Antibody</u>	CSB-PA622788LA01HU	Human	ELISA, IHC, IF
<u>GOLM1 Antibody</u>	CSB-PA10644A0Rb	Human	ELISA, IHC
<u>GOLPH3 Antibody</u>	CSB-PA872494LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>GOLPH3L Antibody</u>	CSB-PA009669LA01HU	Human	ELISA, IHC
<u>GOSR2 Antibody</u>	CSB-PA009678LA01HU	Human	ELISA, WB, IHC
<u>GOT1 Antibody</u>	CSB-PA009679HA01HU	Human, Mouse	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>GOT2 Antibody</u>	CSB-PA03065A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>GP Antibody</u>	CSB-PA310843LA01ZAA	Zaire ebolavirus	ELISA
<u>GP Antibody</u>	CSB-PA314099LA01ZAT	Zaire ebolavirus	ELISA
<u>GP Antibody</u>	CSB-PA350772LA01ZAT	Zaire ebolavirus	ELISA
<u>GP Antibody</u>	CSB-PA308957LA01ZAA	Zaire ebolavirus	ELISA, WB
<u>Gp2 Antibody</u>	CSB-PA880511LA01MO	Mouse	ELISA
<u>GP6 Antibody</u>	CSB-PA872550LA01HU	Human	ELISA, IHC, IF
<u>GP9 Antibody</u>	CSB-PA009690LA01HU	Human	ELISA, IHC, IF
<u>GPA33 Antibody</u>	CSB-PA857459EA01HU	Human	ELISA, WB, IHC, IF
<u>GPAM Antibody</u>	CSB-PA875724LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>GPBAR1 Antibody</u>	CSB-PA819471LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>GPC Antibody</u>	CSB-PA357830LA01LKW	Lymphocytic choriomeningitis virus	ELISA
<u>GPC1 Antibody</u>	CSB-PA009703LA01HU	Human	ELISA, IHC, IF
<u>GPC3 Antibody</u>	CSB-PA12749A0Rb	Human	ELISA, IF
<u>Gpc4 Antibody</u>	CSB-PA009706EA01MO	Mouse	ELISA
<u>GPC4 Antibody</u>	CSB-PA009706LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>GPC6 Antibody</u>	CSB-PA009708LA01HU	Human	ELISA, IHC
<u>GPD1 Antibody</u>	CSB-PA009709HA01HU	Human	ELISA, IHC, IF
<u>gpd-1 Antibody</u>	CSB-PA361431LA01CXY	Caenorhabditis elegans	ELISA
<u>GPD2 Antibody</u>	CSB-PA009711LA01HU	Human	ELISA, IHC
<u>GPER1 Antibody</u>	CSB-PA859933LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>GPHN Antibody</u>	CSB-PA009716LA01HU	Human	ELISA, IHC
<u>GPI Antibody</u>	CSB-PA00367A0Rb	Human, Mouse	ELISA, WB, IHC, IF, IP
<u>Gpi Antibody</u>	CSB-PA009717LA01MO	Mouse	ELISA
<u>GPIHBP1 Antibody</u>	CSB-PA811603LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>GPLD1 Antibody</u>	CSB-PA009721LA01HU	Human	ELISA, WB, IHC
<u>GPM6A Antibody</u>	CSB-PA009722LA01HU	Human, Mouse	ELISA, WB, IF
<u>GPM6B Antibody</u>	CSB-PA614417LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>GPN3 Antibody</u>	CSB-PA009726LA01HU	Human	ELISA, IHC
<u>GPR1 Antibody</u>	CSB-PA009728LA01HU	Human	ELISA, IHC, IF
<u>GPR101 Antibody</u>	CSB-PA846682LA01HU	Human	ELISA, WB, IF
<u>GPR107 Antibody</u>	CSB-PA736202LA01HU	Human	ELISA, IHC, IF
<u>GPR12 Antibody</u>	CSB-PA009742LA01HU	Human	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>GPR141 Antibody</u>	CSB-PA773798LA01HU	Human	ELISA
<u>GPR143 Antibody</u>	CSB-PA009758LA01HU	Human, Rat	ELISA, WB
<u>GPR146 Antibody</u>	CSB-PA850271LA01HU	Human	ELISA, IHC, IF
<u>GPR156 Antibody</u>	CSB-PA844046LA01HU	Human	ELISA, IHC
<u>GPR158 Antibody</u>	CSB-PA735975LA01HU	Human	ELISA, IHC, IF
<u>GPR161 Antibody</u>	CSB-PA04114A0Rb	Human	ELISA, WB, IF
<u>GPR17 Antibody</u>	CSB-PA619758LA01HU	Human, Mouse	ELISA, WB, IF
<u>GPR173 Antibody</u>	CSB-PA882095LA01HU	Human	ELISA, IHC, IF
<u>GPR174 Antibody</u>	CSB-PA887146LA01HU	Human	ELISA, WB, IHC, IF
<u>GPR176 Antibody</u>	CSB-PA618889LA01HU	Human	ELISA, WB, IF
<u>GPR179 Antibody</u>	CSB-PA757664LA01HU	Human	ELISA, IHC, IF
<u>GPR18 Antibody</u>	CSB-PA614525LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>GPR180 Antibody</u>	CSB-PA769784LA01HU	Human	ELISA, IF
<u>GPR183 Antibody</u>	CSB-PA009792LA01HU	Human	ELISA, IHC, IF
<u>GPR21 Antibody</u>	CSB-PA860348LA01HU	Human	ELISA, IHC, IF
<u>GPR22 Antibody</u>	CSB-PA859519LA01HU	Human, Mouse	ELISA, WB, IF
<u>GPR27 Antibody</u>	CSB-PA885717LA01HU	Human	ELISA, IHC, IF
<u>GPR3 Antibody</u>	CSB-PA009800LA01HU	Human	ELISA, IHC, IF
<u>GPR32 Antibody</u>	CSB-PA009803LA01HU	Human	ELISA, WB, IHC, IF
<u>GPR34 Antibody</u>	CSB-PA891571LA01HU	Human, Mouse	ELISA, WB, IF
<u>GPR37 Antibody</u>	CSB-PA009810NA01HU	Human	ELISA, WB, IHC, IF
<u>GPR37 Antibody</u>	CSB-PA009810OA01HU	Human	ELISA, IHC
<u>GPR39 Antibody</u>	CSB-PA009812LA01HU	Human	ELISA, IHC, IF
<u>GPR4 Antibody</u>	CSB-PA009813LA01HU	Human	ELISA, IHC, IF
<u>GPR52 Antibody</u>	CSB-PA896700LA01HU	Human	ELISA, IF
<u>GPR55 Antibody</u>	CSB-PA009820OA01HU	Human	ELISA, IHC
<u>GPR75 Antibody</u>	CSB-PA009829LA01HU	Human	ELISA, WB, IHC, IF
<u>GPR82 Antibody</u>	CSB-PA856999LA01HU	Human	ELISA, IHC, IF
<u>GPR83 Antibody</u>	CSB-PA865180LA01HU	Human	ELISA, WB, IF
<u>GPR84 Antibody</u>	CSB-PA865107LA01HU	Human	ELISA, IF
<u>GPR85 Antibody</u>	CSB-PA009837LA01HU	Human	ELISA, IHC
<u>GPR98 Antibody</u>	CSB-PA23089A0Rb	Human	ELISA, IF
<u>GPRC5B Antibody</u>	CSB-PA868367LA01HU	Human	ELISA, IF
<u>GPRC6A Antibody</u>	CSB-PA722571LA01HU	Human	ELISA, WB, IHC, IF
<u>GPS2 Antibody</u>	CSB-PA009859LA01HU	Human	ELISA, IF
<u>GPSM1 Antibody</u>	CSB-PA801266LA01HU	Human	ELISA, IF



Product Name	Code	Species Reactivity	Applications
<u>GP3M3 Antibody</u>	CSB-PA896517LA01HU	Human	ELISA, IHC
<u>GPT Antibody</u>	CSB-PA009864LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF, IP
<u>GPT2 Antibody</u>	CSB-PA823460LA01HU	Human	ELISA, WB, IHC
<u>GPX2 Antibody</u>	CSB-PA009867LA01HU	Human	ELISA, IF
<u>GPX4 Antibody</u>	CSB-PA03049A0Rb	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>GPX5 Antibody</u>	CSB-PA009870YA01PI	Pig	ELISA
<u>GPX5 Antibody</u>	CSB-PA009870LA01PI	Pig	ELISA
<u>GPX7 Antibody</u>	CSB-PA009872LA01HU	Human	ELISA, IHC
<u>GPX8 Antibody</u>	CSB-PA848417LA01HU	Human	ELISA, IHC, IF
<u>GRA1 Antibody</u>	CSB-PA319538YA01TOV	Toxoplasma gondii	ELISA
<u>GRAMD1B Antibody</u>	CSB-PA662055LA01HU	Human	ELISA, IHC, IF
<u>GRAMD4 Antibody</u>	CSB-PA764225LA01HU	Human, Mouse	ELISA, WB, IHC
<u>GRAP Antibody</u>	CSB-PA618789LA01HU	Human	ELISA, IHC
<u>GRAP2 Antibody</u>	CSB-PA009884LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>GRB10 Antibody</u>	CSB-PA009887LA01HU	Human	ELISA, IHC, IF
<u>GRB14 Antibody</u>	CSB-PA009888LA01HU	Human	ELISA, IHC, IF
<u>GRB2 Antibody</u>	CSB-PA009889EA01HU	Human	ELISA
<u>GREM1 Antibody</u>	CSB-PA009892LA01HU	Human	ELISA, IHC
<u>GRHL1 Antibody</u>	CSB-PA868368LA01HU	Human	ELISA, WB, IHC
<u>GRHL2 Antibody</u>	CSB-PA750350LA01HU	Human	ELISA, IHC, IF
<u>GRHL3 Antibody</u>	CSB-PA844730LA01HU	Human	ELISA, IF
<u>GRHPR Antibody</u>	CSB-PA009897HA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>GRIA3 Antibody</u>	CSB-PA009900LA01HU	Human	ELISA, IHC, IF
<u>GRIA4 Antibody</u>	CSB-PA009901LA01HU	Human	ELISA, IHC, IF
<u>GRID2IP Antibody</u>	CSB-PA009904LA01HU	Human	ELISA, IHC, IF
<u>GRIK1 Antibody</u>	CSB-PA009906LA01HU	Human	ELISA, IHC, IF
<u>GRIK2 Antibody</u>	CSB-PA618751LA01HU	Human	ELISA, IHC, IF
<u>GRIK4 Antibody</u>	CSB-PA619054LA01HU	Human	ELISA, WB
<u>GRIN1 Antibody</u>	CSB-PA009911LA01HU	Human	ELISA, IF
<u>GRIN2A Antibody</u>	CSB-PA14129A0Rb	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>GRIN2C Antibody</u>	CSB-PA009914LA01HU	Human	ELISA, IHC, IF
<u>GRK2 Antibody</u>	CSB-PA00969A0Rb	Human	ELISA, IHC
<u>Grk2 Antibody</u>	CSB-PA858153LA01MO	Mouse	ELISA
<u>GRK3 Antibody</u>	CSB-PA001395LA01HU	Human	ELISA, WB, IHC, IF
<u>GRK4 Antibody</u>	CSB-PA009925LA01HU	Human	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>GRK5 Antibody</u>	CSB-PA009926LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>GRK7 Antibody</u>	CSB-PA009929LA01HU	Human	ELISA, IF
<u>GRM1 Antibody</u>	CSB-PA009931LA01HU	Human	ELISA, IHC, IF
<u>GRM2 Antibody</u>	CSB-PA621775LA01HU	Human, Rat	ELISA, WB, IF
<u>GRM3 Antibody</u>	CSB-PA622790LA01HU	Human	ELISA, IF
<u>GRM6 Antibody</u>	CSB-PA009936LA01HU	Human	ELISA, IF
<u>groL Antibody</u>	CSB-PA394210HA01YAL	Yersinia pestis	ELISA
<u>groL Antibody</u>	CSB-PA330104HA01DRZ	Chlamydia pneumoniae	ELISA
<u>GRP Antibody</u>	CSB-PA009940EA01HU	Human	ELISA, IHC, IF
<u>GRP Antibody</u>	CSB-PA009940HA01HU	Human	ELISA, IHC
<u>GRPR Antibody</u>	CSB-PA009943LA01HU	Human	ELISA, IHC, IF
<u>grxD Antibody</u>	CSB-PA359726LA01ENV	Escherichia coli	ELISA, WB
<u>grxD Antibody</u>	CSB-PA359726LA11ENV	Escherichia coli	ELISA, WB
<u>GSC Antibody</u>	CSB-PA009950LA01HU	Human	ELISA, IHC
<u>GSDMA Antibody</u>	CSB-PA850415LA01HU	Human	ELISA, WB, IHC, IF
<u>GSDMB Antibody</u>	CSB-PA823445LA01HU	Human	ELISA, WB, IHC
<u>GSDMC Antibody</u>	CSB-PA883633LA01HU	Human	ELISA, IHC
<u>GSDMD Antibody</u>	CSB-PA009956LA01HU	Human	ELISA, WB
<u>GSDME Antibody</u>	CSB-PA006766LA01HU	Human	ELISA, IHC, IF
<u>GSG1L Antibody</u>	CSB-PA740924LA01HU	Human	ELISA, IF
<u>gshA Antibody</u>	CSB-PA363915HA01ENV	Escherichia coli	ELISA
<u>GSK3A Antibody</u>	CSB-PA05369A0Rb	Human	ELISA, WB, IF
<u>GSK3B Antibody</u>	CSB-PA009963LA01HU	Human, Mouse	ELISA, WB, IF
<u>GSN Antibody</u>	CSB-PA009965HA01HU	Human	ELISA, WB, IHC, IF
<u>GSPT1 Antibody</u>	CSB-PA009966LA01HU	Human	ELISA, IHC, IF
<u>GSR Antibody</u>	CSB-PA009968LA01HU	Human	ELISA, IHC, IF
<u>GSS Antibody</u>	CSB-PA009969LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>GST Monoclonal Antibody</u>	CSB-MA000031M0m	All	ELISA, WB
<u>GSTA1 Antibody</u>	CSB-PA009970LA01HU	Human	ELISA, IHC, IF
<u>GSTA1 Antibody</u>	CSB-PA009970HA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>Gsta1 Antibody</u>	CSB-PA009970HA01MO	Mouse	ELISA, WB
<u>GSTA4 Antibody</u>	CSB-PA009973HA01HU	Human	ELISA, IHC, IF
<u>GSTK1 Antibody</u>	CSB-PA009978HA01HU	Human	ELISA, WB
<u>GSTM1 Antibody</u>	CSB-PA009979LA01HU	Human	ELISA, IHC, IF
<u>GSTM3 Antibody</u>	CSB-PA009982LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>Gstp1 Antibody</u>	CSB-PA009989HA01MO	Mouse, Human	ELISA, WB
<u>Gstp1 Antibody</u>	CSB-PA009989HA01RA	Rat, Human	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>GSTP1 Antibody</u>	CSB-PA10334A0Rb	Human	ELISA, WB
<u>GSTP1 Antibody</u>	CSB-PA10337A0Rb	Human	ELISA
<u>GSTT1 Antibody</u>	CSB-PA009991LA01HU	Human	ELISA, IHC, IF
<u>GSX1 Antibody</u>	CSB-PA009995LA01HU	Human	ELISA, IHC
<u>GSX2 Antibody</u>	CSB-PA009996LA01HU	Human	ELISA, IHC
<u>GT-1 Antibody</u>	CSB-PA861892LA01DOA	Arabidopsis thaliana	ELISA
<u>GTDC1 Antibody</u>	CSB-PA676737LA01HU	Human	ELISA, WB
<u>GTF2B Antibody</u>	CSB-PA010001HA01HU	Human	ELISA, WB, IHC, IF
<u>GTF2E2 Antibody</u>	CSB-PA010003LA01HU	Human, Mouse	ELISA, WB
<u>GTF2F1 Antibody</u>	CSB-PA010004HA01HU	Human	ELISA, IHC
<u>GTF2F2 Antibody</u>	CSB-PA010005LA01HU	Human, Mouse	ELISA, WB
<u>GTF2H2 Antibody</u>	CSB-PA618808LA01HU	Human	ELISA
<u>GTF2H5 Antibody</u>	CSB-PA010014EA01HU	Human	ELISA, IHC, IF
<u>GTF2IRD1 Antibody</u>	CSB-PA888303LA01HU	Human	ELISA, IHC
<u>GTF3C5 Antibody</u>	CSB-PA896534LA01HU	Human	ELISA, IHC
<u>GTF3C6 Antibody</u>	CSB-PA853368LA01HU	Human	ELISA, IHC
<u>GTPBP1 Antibody</u>	CSB-PA010031LA01HU	Human	ELISA, IHC, IF
<u>GTPBP4 Antibody</u>	CSB-PA887173LA01HU	Human	ELISA, WB, IF
<u>GTPBP6 Antibody</u>	CSB-PA04439A0Rb	Human	ELISA, WB, IHC
<u>GTSF1 Antibody</u>	CSB-PA848827LA01HU	Human	ELISA
<u>Guanine nucleotide-binding protein subunit beta-like protein Antibody</u>	CSB-PA659907LA01GGV	Glycine max	ELISA, WB
<u>GUCA1A Antibody</u>	CSB-PA010044LA01HU	Human	ELISA
<u>Guca2b Antibody</u>	CSB-PA010048LA01MO	Mouse	ELISA
<u>GUCA2B Antibody</u>	CSB-PA613693LA01HU	Human	ELISA, IHC, IF
<u>GUCD1 Antibody</u>	CSB-PA839384LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>GUCY1A1 Antibody</u>	CSB-PA010050LA01HU	Human	ELISA, IHC, IF
<u>GUCY2C Antibody</u>	CSB-PA010053LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>GUCY2D Antibody</u>	CSB-PA010054LA01HU	Human	ELISA, WB, IF
<u>GULP1 Antibody</u>	CSB-PA887007LA01HU	Human, Rat	ELISA, WB, IF
<u>GV51_0914 Antibody</u>	CSB-PA333053LA01GHL	Gardnerella vaginalis 5-1	ELISA
<u>GYPA Antibody</u>	CSB-PA010074DA01HU	Human	ELISA
<u>GYPE Antibody</u>	CSB-PA010077LA01HU	Human	ELISA
<u>Gzma Antibody</u>	CSB-PA010081LA01MO	Mouse, Human	ELISA, WB
<u>GZMB Antibody</u>	CSB-PA010082LA01HU	Human	ELISA, IHC, IF
<u>GZMK Antibody</u>	CSB-PA010084LA01HU	Human	ELISA, IHC, IF
<u>GZMM Antibody</u>	CSB-PA010085LA01HU	Human	ELISA, IHC, IF
<u>H2AFB1 Antibody</u>	CSB-PA010091LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>H2AFJ Antibody</u>	CSB-PA010094DA01HU	Human, Mouse	ELISA, WB, IHC
<u>H2AFV Antibody</u>	CSB-PA03709A0Rb	Human	ELISA
<u>H2AFX Antibody</u>	CSB-PA15429A0Rb	Human	ELISA, IHC
<u>H2AFY2 Antibody</u>	CSB-PA885785LA01HU	Human	ELISA
<u>H2-D1 Antibody</u>	CSB-PA365684LA01MO	Mouse, Human	ELISA, WB
<u>H2-K1 Antibody</u>	CSB-PA355780EA01MO	Mouse	ELISA
<u>H3F3A Antibody</u>	CSB-PA010109LA01HU	Human	ELISA, WB, IHC
<u>H3F3A Antibody</u>	CSB-PA01435A0Rb	Human	ELISA, WB, IHC
<u>H671_3g9701 Antibody</u>	CSB-PA006187LA01DXU	Cricetulus griseus	ELISA
<u>HA Antibody</u>	CSB-PA18939A0Rb	Influenza A virus(H7N9)	ELISA
<u>HA Antibody</u>	CSB-PA356317LA01IDG	Influenza A virus	ELISA
<u>HA Antibody</u>	CSB-PA607773LA01IER	Influenza A virus	ELISA, WB
<u>HA-33 Antibody</u>	CSB-PA338424LA01CLQ	Clostridium botulinum	ELISA, WB
<u>HAAO Antibody</u>	CSB-PA010112LA01HU	Human	ELISA, IHC, IF
<u>HABP2 Antibody</u>	CSB-PA12159A0Rb	Human	ELISA, IHC
<u>HACD2 Antibody</u>	CSB-PA761567LA01HU	Human	ELISA, WB, IHC, IF
<u>HACD3 Antibody</u>	CSB-PA878945LA01HU	Human	ELISA, IHC
<u>HACE1 Antibody</u>	CSB-PA818264LA01HU	Human	ELISA, IHC
<u>HADHB Antibody</u>	CSB-PA03035A0Rb	Human	ELISA, IHC
<u>HAGH Antibody</u>	CSB-PA11195A0Rb	Human	ELISA, IHC
<u>HAL Antibody</u>	CSB-PA010123LA01HU	Human	ELISA, IHC, IF
<u>HAO1 Antibody</u>	CSB-PA891938LA01HU	Human	ELISA, IHC, IF
<u>HAO2 Antibody</u>	CSB-PA882145LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>HAP1 Antibody</u>	CSB-PA010129LA01HU	Human	ELISA, IF
<u>HARS2 Antibody</u>	CSB-PA010138HA01HU	Human, Mouse	ELISA, WB, IHC
<u>HAS1 Antibody</u>	CSB-PA010139LA01HU	Human	ELISA, IF
<u>Has2 Antibody</u>	CSB-PA010140LA01MO	Mouse	ELISA
<u>HAS2 Antibody</u>	CSB-PA010140LA01HU	Human	ELISA, IF
<u>HAT1 Antibody</u>	CSB-PA010143EA01HU	Human, Rat	ELISA, WB, IF
<u>HA-Tag</u>	CSB-MA000141M0m	Human,Mouse	ELISA, WB, IF, FC, IP
<u>HAVCR1 Antibody</u>	CSB-PA09995A0Rb	Human, Mouse	ELISA, WB, IHC
<u>Havcr1 Antibody</u>	CSB-PA17109A0Rb	Mouse	ELISA, IHC
<u>HAVCR2 Antibody</u>	CSB-PA010145LA01HU	Human	ELISA, IHC, IF
<u>HAX1 Antibody</u>	CSB-PA03515A0Rb	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>Hba Antibody</u>	CSB-PA010147LA01MO	Mouse, Human	ELISA, WB
<u>HBA1 Antibody</u>	CSB-PA12175A0Rb	Human	ELISA
<u>HBB Antibody</u>	CSB-PA010150LA01HU	Human, Rat	ELISA, WB, IHC
<u>Hbb-b2 Antibody</u>	CSB-PA360812EA01MO	Mouse	ELISA
<u>HBD Antibody</u>	CSB-PA010152LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>HBE1 Antibody</u>	CSB-PA010153LA01HU	Human	ELISA, IHC
<u>HBEGF Antibody</u>	CSB-PA857429LA01HU	Human	ELISA, IHC, IF
<u>HBG1 Antibody</u>	CSB-PA010155HA01HU	Human	ELISA, IHC
<u>HBG2 Antibody</u>	CSB-PA010156LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>HBQ1 Antibody</u>	CSB-PA010159LA01HU	Human	ELISA
<u>HCFC1 Antibody</u>	CSB-PA010166LA01HU	Human	ELISA, IHC, IF
<u>HCFC1R1 Antibody</u>	CSB-PA865161LA01HU	Human	ELISA, IHC
<u>HCN1 Antibody</u>	CSB-PA010216LA01HU	Human	ELISA, IHC, IF
<u>HCP5 Antibody</u>	CSB-PA754020LA01HU	Human	ELISA, IHC
<u>HCRT1 Antibody</u>	CSB-PA010231HA01HU	Human	ELISA, IF
<u>HDAC1 Antibody</u>	CSB-PA17869A0Rb	Human	ELISA, IHC, IF, ChIP
<u>HDAC10 Antibody</u>	CSB-PA010236LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>HDAC3 Antibody</u>	CSB-PA010239LA01HU	Human, Mouse	ELISA, WB, IHC, IF, ChIP
<u>HDAC4 Antibody</u>	CSB-PA17884A0Rb	Human	ELISA, IHC
<u>HDAC5 Antibody</u>	CSB-PA010241HA01HU	Human	ELISA, IHC, IF
<u>HDAC7 Antibody</u>	CSB-PA17899A0Rb	Human	ELISA
<u>HDAC8 Antibody</u>	CSB-PA010244LA01HU	Human	ELISA, IHC, IF
<u>HDAC8 Antibody</u>	CSB-PA010244EA01HU	Human, Rat	ELISA, WB
<u>HDAC9 Antibody</u>	CSB-PA17907A0Rb	Human	ELISA, IF
<u>HDC Antibody</u>	CSB-PA13409A0Rb	Human	ELISA, WB, IHC
<u>HDDC3 Antibody</u>	CSB-PA836666HA01HU	Human, Mouse	ELISA, WB, IHC
<u>HDGF Antibody</u>	CSB-PA01004A0Rb	Human	ELISA
<u>HDHD2 Antibody</u>	CSB-PA880936LA01HU	Human	ELISA, IHC, IF
<u>HDHD3 Antibody</u>	CSB-PA883592LA01HU	Human, Mouse	ELISA, WB, IHC
<u>HDT2 Antibody</u>	CSB-PA684969YA01DOA	Arabidopsis thaliana	ELISA
<u>HEATR9 Antibody</u>	CSB-PA003369LA01HU	Human	ELISA, WB, IHC
<u>HEBP2 Antibody</u>	CSB-PA896539HA01HU	Human	ELISA, IHC
<u>HECTD2 Antibody</u>	CSB-PA732977LA01HU	Human	ELISA, IHC, IF
<u>HECTD3 Antibody</u>	CSB-PA725567LA01HU	Human	ELISA, IF
<u>HECW1 Antibody</u>	CSB-PA748484HA01HU	Human	ELISA, IF
<u>HECW2 Antibody</u>	CSB-PA889188LA01HU	Human	ELISA, WB, IHC, IF
<u>HELLS Antibody</u>	CSB-PA868307HA01HU	Human	ELISA, IHC
<u>HEMK1 Antibody</u>	CSB-PA010285LA01HU	Human	ELISA, IHC
<u>Hemoglobin Antibody</u>	CSB-PA00240E0Rb	Human	ELISA, WB
<u>Hemoglobin Antibody</u>	CSB-PA00480E0Rb	Guinea Pig	ELISA, WB
<u>Hemoglobin Antibody</u>	CSB-PA00490E0Rb	Mouse	ELISA, WB
<u>Hemoglobin Antibody</u>	CSB-PA00500E0Gp	Rabbit	ELISA



Product Name	Code	Species Reactivity	Applications
<u>Hemoglobin Antibody</u>	CSB-PA00790E0Gp	Pigeon	ELISA
<u>HEPACAM Antibody</u>	CSB-PA010287LA01HU	Human	ELISA, WB, IF
<u>HEPACAM2 Antibody</u>	CSB-PA010288LA01HU	Human	ELISA, IHC
<u>her6 Antibody</u>	CSB-PA20989A0Rb	Zebrafish	ELISA
<u>HERC3 Antibody</u>	CSB-PA618972LA01HU	Human	ELISA
<u>HERC5 Antibody</u>	CSB-PA883435LA01HU	Human	ELISA, IHC, IF
<u>HERC6 Antibody</u>	CSB-PA812876LA01HU	Human	ELISA, IHC
<u>HES2 Antibody</u>	CSB-PA897496LA01HU	Human	ELISA, IF
<u>HES5 Antibody</u>	CSB-PA729212LA01HU	Human	ELISA, IHC
<u>HES7 Antibody</u>	CSB-PA887160LA01HU	Human	ELISA, IHC
<u>HESX1 Antibody</u>	CSB-PA871390LA01HU	Human	ELISA
<u>HEXA Antibody</u>	CSB-PA12259A0Rb	Human, Rat	ELISA, WB, IHC, IF
<u>HEXB Antibody</u>	CSB-PA010316LA01HU	Human	ELISA, IHC
<u>HEXIM2 Antibody</u>	CSB-PA010319LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>HEY1 Antibody</u>	CSB-PA896909LA01HU	Human	ELISA, IHC
<u>hfb1 Antibody</u>	CSB-PA347143HA01HYE	Hypocrea jecorina	ELISA, WB
<u>HFE Antibody</u>	CSB-PA653744LA01HU	Human	ELISA
<u>HFM1 Antibody</u>	CSB-PA010325LA01HU	Human	ELISA, IHC, IF
<u>HHEX Antibody</u>	CSB-PA010335OA01HU	Human	ELISA, IHC
<u>HHIPL1 Antibody</u>	CSB-PA822236LA01HU	Human	ELISA
<u>HHLA2 Antibody</u>	CSB-PA891555LA01HU	Human	ELISA, IF
<u>HHLA3 Antibody</u>	CSB-PA895005LA01HU	Human	ELISA
<u>HIC1 Antibody</u>	CSB-PA623921LA01HU	Human	ELISA, IHC
<u>HIF1A Antibody</u>	CSB-PA12807A0Rb	Human	ELISA, WB, IHC, IF
<u>HIF1AN Antibody</u>	CSB-PA865160LA01HU	Human, Mouse	ELISA, WB, IHC, IF, IP
<u>HIF3A Antibody</u>	CSB-PA896696LA01HU	Human, Mouse	ELISA, WB, IHC
<u>HIGD1B Antibody</u>	CSB-PA889181LA01HU	Human	ELISA, IHC, IF
<u>HIKESHI Antibody</u>	CSB-PA706632HA01HU	Human	ELISA, IHC
<u>HINFP Antibody</u>	CSB-PA010361LA01HU	Human, Rat	ELISA, WB
<u>HINT1 Antibody</u>	CSB-PA05639A0Rb	Human	ELISA, IHC, IF
<u>HIP1R Antibody</u>	CSB-PA010366LA01HU	Human, Rat	ELISA, WB, IHC
<u>HIPK1 Antibody</u>	CSB-PA803159LA01HU	Human	ELISA, IHC, IF
<u>HIPK2 Antibody</u>	CSB-PA867138LA01HU	Human	ELISA, IHC
<u>HIPK3 Antibody</u>	CSB-PA872492LA01HU	Human	ELISA, IHC, IF
<u>HIST1H1C Antibody</u>	CSB-PA010378HA01HU	Human	ELISA, IHC, IF
<u>HIST1H2AB Antibody</u>	CSB-PA010385HA01HU	Human	ELISA, IHC
<u>HIST1H2AC Antibody</u>	CSB-PA846293LA01HU	Human	ELISA, IHC, IF
<u>HIST1H2AG Antibody</u>	CSB-PA010389LA01HU	Human	ELISA, WB, IHC, IF, ChIP



Product Name	Code	Species Reactivity	Applications
<u>HIST1H2BB Antibody</u>	CSB-PA010402HA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>HIST1H2BC Antibody</u>	CSB-PA01694A0Rb	Human	ELISA, IHC, IF
<u>HIST1H2BD Antibody</u>	CSB-PA010404LA01HU	Human	ELISA, IHC, IF
<u>HIST1H2BK Antibody</u>	CSB-PA010411LA01HU	Human	ELISA, IHC, IF
<u>HIST1H2BL Antibody</u>	CSB-PA860357EA01HU	Human	ELISA, IHC
<u>Hist1h2bm Antibody</u>	CSB-PA010413LA01MO	Mouse	ELISA
<u>HIST1H2BN Antibody</u>	CSB-PA857877LA01HU	Human	ELISA, IHC
<u>HIST1H4A Antibody</u>	CSB-PA010429LA01HU	Human, Mouse	ELISA, WB, IHC
<u>HIST2H2AC Antibody</u>	CSB-PA010445LA01HU	Human	ELISA, IHC
<u>HIST2H2BE Antibody</u>	CSB-PA620981LA01HU	Human	ELISA, IHC
<u>hist2h2l Antibody</u>	CSB-PA740767LA01DIL	Zebrafish	ELISA, WB
<u>HIST2H3A Antibody</u>	CSB-PA744479HA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>HIST3H3 Antibody</u>	CSB-PA623092HA01HU	Human	ELISA, IHC
<u>Histidine-rich protein PFHRP-II Antibody</u>	CSB-PA356476HA01PLT	Plasmodium falciparum	ELISA
<u>Histone H2B Antibody</u>	CSB-PA355820LA01CSJ	Cairina moschata	ELISA
<u>HJURP Antibody</u>	CSB-PA010467LA01HU	Human	ELISA, WB
<u>HK1 Antibody</u>	CSB-PA010468HA01HU	Human, Rat	ELISA, WB
<u>HK1 Antibody</u>	CSB-PA010468LA01HU	Human	ELISA
<u>HK2 Antibody</u>	CSB-PA010469LA01HU	Human	ELISA, WB, IF
<u>HKDC1 Antibody</u>	CSB-PA647983LA01HU	Human	ELISA, WB, IHC
<u>HLA-B Antibody</u>	CSB-PA15089A0Rb	Human	ELISA, WB, IHC, IF
<u>HLA-B Antibody</u>	CSB-PA642773LA01HU	Human	ELISA, WB, IHC
<u>HLA-C Antibody</u>	CSB-PA15099A0Rb	Human	ELISA, WB, IHC, IF
<u>HLA-C Antibody</u>	CSB-PA320739LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>HLA-C Antibody</u>	CSB-PA356221LA01HU	Human	ELISA, IHC, IF
<u>HLA-DMA Antibody</u>	CSB-PA338892LA01HU	Human	ELISA, IHC
<u>HLA-DPA1 Antibody</u>	CSB-PA14819A0Rb	Human	ELISA, WB, IF
<u>HLA-DPB1 Antibody</u>	CSB-PA14809A0Rb	Human	ELISA, WB, IHC
<u>HLA-DQA1 Antibody</u>	CSB-PA14839A0Rb	Human	ELISA, IHC
<u>HLA-DQA2 Antibody</u>	CSB-PA14829A0Rb	Human	ELISA, IHC
<u>HLA-DQB1 Antibody</u>	CSB-PA14849A0Rb	Human, Mouse	ELISA, WB, IHC
<u>HLA-DRA Antibody</u>	CSB-PA17939A0Rb	Human	ELISA, IHC
<u>HLA-DRB1 Antibody</u>	CSB-PA17949A0Rb	Human	ELISA
<u>HLA-DRB1 Antibody</u>	CSB-PA19179A0Rb	Human	ELISA, WB, IF
<u>HLA-DRB3 Antibody</u>	CSB-PA302920LA01HU	Human, Mouse	ELISA, WB, IF
<u>HLA-DRB4 Antibody</u>	CSB-PA03909A0Rb	Human	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>HLA-E Antibody</u>	CSB-PA320269LA01HU	Human	ELISA
<u>HLA-F Antibody</u>	CSB-PA335552LA01HU	Human	ELISA, IF
<u>HLA-G Antibody</u>	CSB-PA02109A0Rb	Human	ELISA, WB, IHC, IF
<u>HLCS Antibody</u>	CSB-PA010518LA01HU	Human	ELISA, IHC
<u>hlgB Antibody</u>	CSB-PA357862HA01SKY	Staphylococcus aureus	ELISA
<u>hlgC Antibody</u>	CSB-PA739984LA01SKV	Staphylococcus aureus	ELISA
<u>HLTF Antibody</u>	CSB-PA010520LA01HU	Human	ELISA, WB, IHC, IP
<u>hly Antibody</u>	CSB-PA639324HA01FLF	Staphylococcus aureus	ELISA
<u>HM13 Antibody</u>	CSB-PA851544LA01HU	Human, Mouse	ELISA, WB, IF
<u>HMBS Antibody</u>	CSB-PA010524LA01HU	Human	ELISA, WB, IHC, IF
<u>HMGB1 Antibody</u>	CSB-PA01604A0Rb	Human, Mouse	ELISA, WB, IHC
<u>HMGB1 Antibody</u>	CSB-PA010553YA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>HMGB2 Antibody</u>	CSB-PA862013LA01HU	Human	ELISA, IHC, IF
<u>HMGB3 Antibody</u>	CSB-PA010561LA01HU	Human	ELISA, WB, IHC
<u>HMGB4 Antibody</u>	CSB-PA851975LA01HU	Human	ELISA, IHC
<u>HMGCLL1 Antibody</u>	CSB-PA851531LA01HU	Human	ELISA
<u>HMGCR Antibody</u>	CSB-PA010565LA01HU	Human	ELISA, WB, IHC, IF
<u>Hmgcs1 Antibody</u>	CSB-PA812948HA01MO	Mouse	ELISA
<u>Hmgcs1 Antibody</u>	CSB-PA812948YA01MO	Mouse	ELISA, WB
<u>HMGCS1 Antibody</u>	CSB-PA010566LA01HU	Human	ELISA, IHC, IF
<u>HMGCS2 Antibody</u>	CSB-PA010567LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>HMGN3 Antibody</u>	CSB-PA623009LA01HU	Human	ELISA, IHC
<u>HMGXB4 Antibody</u>	CSB-PA871514PA01HU	Human	ELISA, WB
<u>HMMR Antibody</u>	CSB-PA14189A0Rb	Human	ELISA, IHC, IF
<u>HMOX1 Antibody</u>	CSB-PA010583LA01HU	Human	ELISA, IHC, IF
<u>HMOX2 Antibody</u>	CSB-PA010584HA01HU	Human	ELISA, IF
<u>HMX2 Antibody</u>	CSB-PA010587LA01HU	Human	ELISA, IHC, IF
<u>HNF4A Antibody</u>	CSB-PA010594LA01HU	Human, Mouse	ELISA, WB
<u>HNRNPA1 Antibody</u>	CSB-PA010600HA01HU	Human	ELISA, WB, IHC, IF
<u>HNRNPA2B1 Antibody</u>	CSB-PA010602DA01HU	Human	ELISA, IHC, IF
<u>HNRNPAB Antibody</u>	CSB-PA010604LA01HU	Human	ELISA, WB, IHC, IF
<u>HNRNPD Antibody</u>	CSB-PA02544A0Rb	Human	ELISA, WB, IHC
<u>HNRNPF Antibody</u>	CSB-PA010608LA01HU	Human	ELISA, IHC
<u>HNRNPH1 Antibody</u>	CSB-PA010609HA01HU	Human	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>HNRNPH1 Antibody</u>	CSB-PA04264A0Rb	Human	ELISA, WB
<u>HNRNPL Antibody</u>	CSB-PA010612HA01HU	Human	ELISA, IHC, IF
<u>HNRNPL Antibody</u>	CSB-PA010612LA01HU	Human, Mouse	ELISA, WB, IHC
<u>HNRNPLL Antibody</u>	CSB-PA010630LA01HU	Human	ELISA, IHC, IF
<u>hns Antibody</u>	CSB-PA359764LA01ENV	Escherichia coli	ELISA, WB
<u>hns Antibody</u>	CSB-PA359764EA11ENV	Escherichia coli	ELISA, WB
<u>HOGA1 Antibody</u>	CSB-PA771464LA01HU	Human	ELISA, IHC, IF
<u>HOMER1 Antibody</u>	CSB-PA03007A0Rb	Human	ELISA, IHC
<u>HOMER2 Antibody</u>	CSB-PA885720LA01HU	Human	ELISA, IHC
<u>HOMER3 Antibody</u>	CSB-PA878880LA01HU	Human, Mouse	ELISA, WB, IHC
<u>HOOK1 Antibody</u>	CSB-PA891526LA01HU	Human	ELISA, WB, IHC
<u>HOOK3 Antibody</u>	CSB-PA769790LA01HU	Human	ELISA, WB, IHC, IP
<u>HOPX Antibody</u>	CSB-PA859411LA01HU	Human	ELISA, IHC
<u>HORMAD1 Antibody</u>	CSB-PA803140LA01HU	Human	ELISA
<u>HORMAD2 Antibody</u>	CSB-PA836683LA01HU	Human	ELISA, IHC
<u>HOXA11 Antibody</u>	CSB-PA010649LA01HU	Human	ELISA, IHC, IF
<u>HOXA2 Antibody</u>	CSB-PA010652LA01HU	Human, Mouse	ELISA, WB, IHC
<u>HOXA4 Antibody</u>	CSB-PA010654LA01HU	Human, Rat	ELISA, WB, IHC
<u>HOXA5 Antibody</u>	CSB-PA010655LA01HU	Human, Rat	ELISA, WB, IHC
<u>HOXA6 Antibody</u>	CSB-PA010656EA01HU	Human, Mouse	ELISA, WB, IHC
<u>HOXB13 Antibody</u>	CSB-PA842162LA01HU	Human, Mouse	ELISA, WB
<u>HOXB3 Antibody</u>	CSB-PA010663LA01HU	Human	ELISA, IHC
<u>HOXB4 Antibody</u>	CSB-PA010664LA01HU	Human	ELISA, IHC, IF
<u>HOXC11 Antibody</u>	CSB-PA010672LA01HU	Human	ELISA, IHC, IF
<u>HOXC6 Antibody</u>	CSB-PA010677LA01HU	Human	ELISA, IF
<u>HOXC9 Antibody</u>	CSB-PA010679LA01HU	Human	ELISA
<u>HOXD1 Antibody</u>	CSB-PA884426LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>HOXD12 Antibody</u>	CSB-PA010684LA01HU	Human	ELISA, IF
<u>HOXD13 Antibody</u>	CSB-PA010685LA01HU	Human	ELISA, IHC
<u>HOXD3 Antibody</u>	CSB-PA010686LA01HU	Human, Rat	ELISA, WB
<u>HOXD8 Antibody</u>	CSB-PA010688LA01HU	Human, Mouse	ELISA, WB
<u>HOXD9 Antibody</u>	CSB-PA010689LA01HU	Human, Mouse	ELISA, WB, IHC
<u>HP Antibody</u>	CSB-PA14609A0Rb	Human, Mouse	ELISA, WB, IHC
<u>Hp Antibody</u>	CSB-PA17069A0Rb	Rat	ELISA, WB, IHC
<u>HPCAL1 Antibody</u>	CSB-PA02269A0Rb	Human	ELISA, IHC, IF
<u>HPCAL4 Antibody</u>	CSB-PA891785LA01HU	Human, Mouse	ELISA, WB
<u>HPD Antibody</u>	CSB-PA010698LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>HPGD Antibody</u>	CSB-PA010702LA01HU	Human	ELISA, IHC, IF
<u>HPN Antibody</u>	CSB-PA010704LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>HPR Antibody</u>	CSB-PA24659A0Rb	Arabidopsis thaliana	ELISA, WB
<u>HPRT1 Antibody</u>	CSB-PA15827A0Rb	Human	ELISA, IHC
<u>HPS6 Antibody</u>	CSB-PA803158LA01HU	Human	ELISA, IHC, IF
<u>HPSE Antibody</u>	CSB-PA010716LA01HU	Human	ELISA, WB, IHC
<u>Hpx Antibody</u>	CSB-PA010723LA01RA	Rat, Human	ELISA, WB
<u>HPX Antibody</u>	CSB-PA14619A0Rb	Human	ELISA, IHC, IF
<u>Hras Antibody</u>	CSB-PA15539A0Rb	Rat, Mouse, Human	ELISA, WB, IHC
<u>HRAS Antibody</u>	CSB-PA010726LA01HU	Human	ELISA, IHC, IF
<u>HRG Antibody</u>	CSB-PA010736LA01HU	Human	ELISA, WB, IHC
<u>HRH1 Antibody</u>	CSB-PA010737LA01HU	Human	ELISA, WB, IHC, IF
<u>HRH3 Antibody</u>	CSB-PA897568LA01HU	Human	ELISA, WB, IHC, IF
<u>HRH4 Antibody</u>	CSB-PA887972LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>HS2ST1 Antibody</u>	CSB-PA010746LA01HU	Human, Rat	ELISA, WB, IHC
<u>HS3ST1 Antibody</u>	CSB-PA010747LA01HU	Human	ELISA, IHC
<u>HS6ST1 Antibody</u>	CSB-PA010756LA01HU	Human, Mouse	ELISA, WB, IF
<u>HSBP1 Antibody</u>	CSB-PA00454A0Rb	Human	ELISA, IHC
<u>HSCB Antibody</u>	CSB-PA816893LA01HU	Human, Rat	ELISA, WB, IHC
<u>HSD11B1 Antibody</u>	CSB-PA010763LA01HU	Human	ELISA, WB, IHC, IF
<u>HSD11B1L Antibody</u>	CSB-PA770354LA01HU	Human	ELISA, IHC
<u>HSD11B2 Antibody</u>	CSB-PA010765LA01HU	Human	ELISA, WB, IHC
<u>HSD17B1 Antibody</u>	CSB-PA010766LA01HU	Human	ELISA, IHC, IF
<u>HSD17B10 Antibody</u>	CSB-PA860776HA01HU	Human	ELISA, WB, IF, IP
<u>HSD17B11 Antibody</u>	CSB-PA843318HA01HU	Human, Mouse	ELISA, WB
<u>HSD17B12 Antibody</u>	CSB-PA687494LA01HU	Human	ELISA, WB, IHC, IF
<u>HSD17B14 Antibody</u>	CSB-PA859026HA01HU	Human, Mouse	ELISA, WB
<u>HSD17B7 Antibody</u>	CSB-PA010776LA01HU	Human	ELISA, IHC, IF
<u>HSD3B1 Antibody</u>	CSB-PA010781LA01HU	Human	ELISA, WB, IHC
<u>HSD3B7 Antibody</u>	CSB-PA861999LA01HU	Human, Mouse	ELISA, WB, IF
<u>HSDL2 Antibody</u>	CSB-PA761576LA01HU	Human	ELISA, WB, IHC
<u>HSF1 Antibody</u>	CSB-PA010791LA01HU	Human	ELISA, WB, IHC, IF
<u>HSF2 Antibody</u>	CSB-PA010792LA01HU	Human	ELISA, IHC, IF
<u>HSP90AB1 Antibody</u>	CSB-PA00109A0Rb	Human	ELISA, WB, IHC, IF
<u>HSP90B1 Antibody</u>	CSB-PA10887A0Rb	Human	ELISA, WB, IHC, IF
<u>HSPA12A Antibody</u>	CSB-PA010817LA01HU	Human	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>HSPA12B Antibody</u>	CSB-PA10989A0Rb	Human	ELISA, IHC
<u>HSPA13 Antibody</u>	CSB-PA10964A0Rb	Human	ELISA, IHC
<u>HSPA1A Antibody</u>	CSB-PA10899A0Rb	Human	ELISA, WB, IHC, IF
<u>HSPA1A Antibody</u>	CSB-PA357302LA01HU	Human	ELISA, WB, IHC
<u>Hspa1b Antibody</u>	CSB-PA010822HA01MO	Mouse, Human	ELISA, WB
<u>HSPA1B Antibody</u>	CSB-PA28047A0Rb	Human	ELISA, WB, IF
<u>HSPA1L Antibody</u>	CSB-PA010823HA01HU	Human	ELISA, WB, IHC, IF
<u>HSPA4 Antibody</u>	CSB-PA010825LA01HU	Human	ELISA, IF
<u>HSPA5 Antibody</u>	CSB-PA05027A0Rb	Human	ELISA, WB, IHC, IF
<u>Hspa5 Antibody</u>	CSB-PA010827YA01MO	Mouse, Rat	ELISA, WB, IF
<u>HSPA6 Antibody</u>	CSB-PA863082HA01HU	Human	ELISA, WB
<u>HSPA8 Antibody</u>	CSB-PA03249A0Rb	Human, Mouse, Rat	ELISA, WB, IHC
<u>HSPA9 Antibody</u>	CSB-PA010830HA01HU	Human	ELISA, IHC
<u>HSPB1 Antibody</u>	CSB-PA10869A0Rb	Human	ELISA, WB, IHC
<u>HSPB6 Antibody</u>	CSB-PA010837LA01HU	Human, Mouse	ELISA, WB
<u>HSPB7 Antibody</u>	CSB-PA10979A0Rb	Human	ELISA, IHC
<u>HSPBAP1 Antibody</u>	CSB-PA853412LA01HU	Human	ELISA, IHC
<u>HSPD1 Antibody</u>	CSB-PA010847LA01HU	Human, Mouse	ELISA, WB, IHC
<u>HSPE1 Antibody</u>	CSB-PA10857A0Rb	Human	ELISA, IHC
<u>htpG Antibody</u>	CSB-PA358897LA01ENV	Escherichia coli	ELISA, WB
<u>HTR1A Antibody</u>	CSB-PA010881LA01HU	Human	ELISA, IHC, IF
<u>HTR1E Antibody</u>	CSB-PA010885LA01HU	Human	ELISA, WB, IHC, IF
<u>HTR2A Antibody</u>	CSB-PA010887LA01HU	Human	ELISA, IHC, IF
<u>HTR3A Antibody</u>	CSB-PA010890LA01HU	Human, Mouse, Rat	ELISA, WB, IF
<u>HTR3D Antibody</u>	CSB-PA747864LA01HU	Human	ELISA, WB, IF
<u>HTR3E Antibody</u>	CSB-PA010894LA01HU	Human	ELISA, WB
<u>HTR4 Antibody</u>	CSB-PA010895LA01HU	Human	ELISA, IHC, IF
<u>HTR7 Antibody</u>	CSB-PA010899LA01HU	Human	ELISA, IHC
<u>HTRA1 Antibody</u>	CSB-PA835695LA01HU	Human	ELISA, IHC, IF
<u>HTRA2 Antibody</u>	CSB-PA010902EA01HU	Human	ELISA, WB, IHC, IP
<u>hup1 Antibody</u>	CSB-PA358373DA01FQI	Streptomyces lividans	ELISA
<u>hupB Antibody</u>	CSB-PA363134HA01SWW	Salmonella typhi	ELISA, WB
<u>HUS1 Antibody</u>	CSB-PA010909LA01HU	Human	ELISA, WB, IHC, IF
<u>HUS1B Antibody</u>	CSB-PA844090LA01HU	Human	ELISA, IHC
<u>HYDIN Antibody</u>	CSB-PA670104LA01HU	Human	ELISA
<u>HYI Antibody</u>	CSB-PA719395LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>HYPM Antibody</u>	CSB-PA006268LA01HU	Human	ELISA, IHC
<u>I79_001071 Antibody</u>	CSB-PA20599A0Rb	Cricetulus griseus	ELISA
<u>IAPP Antibody</u>	CSB-PA010931YA01HU	Human	ELISA
<u>IAPP Antibody</u>	CSB-PA12224A0Rb	Human, Mouse	ELISA, WB, IHC
<u>IBA57 Antibody</u>	CSB-PA722566LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ICA1 Antibody</u>	CSB-PA12109A0Rb	Human, Mouse	ELISA, WB, IHC
<u>ICA1L Antibody</u>	CSB-PA847691LA01HU	Human	ELISA, WB
<u>ICAM1 Antibody</u>	CSB-PA07149A0Rb	Human, Rat	ELISA, WB, IHC, IF
<u>ICE1 Antibody</u>	CSB-PA897466LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ICE2 Antibody</u>	CSB-PA717659LA01HU	Human	ELISA, WB
<u>ICK Antibody</u>	CSB-PA891800LA01HU	Human	ELISA, WB, IHC
<u>ID2 Antibody</u>	CSB-PA010967LA01HU	Human	ELISA, IHC, IF
<u>ID3 Antibody</u>	CSB-PA01959A0Rb	Human	ELISA, IF
<u>ID4 Antibody</u>	CSB-PA010970LA01HU	Human	ELISA, IF
<u>IDE Antibody</u>	CSB-PA010988LA01HU	Human	ELISA, IHC, IF
<u>Idh1 Antibody</u>	CSB-PA010989LA01MO	Mouse	ELISA, IHC, IF
<u>IDH1 Antibody</u>	CSB-PA010989EA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>IDH3A Antibody</u>	CSB-PA010991LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF, IP
<u>IDH3B Antibody</u>	CSB-PA05515A0Rb	Human, Mouse	ELISA, WB
<u>IDH3G Antibody</u>	CSB-PA010993LA01HU	Human	ELISA
<u>IDI2 Antibody</u>	CSB-PA871603LA01HU	Human	ELISA, IHC
<u>IDO1 Antibody</u>	CSB-PA010996HA01HU	Human	ELISA, WB, IHC, IF
<u>IDO2 Antibody</u>	CSB-PA010997LA01HU	Human	ELISA, IHC, IF
<u>IDUA Antibody</u>	CSB-PA011000LA01HU	Human	ELISA, IHC
<u>IER2 Antibody</u>	CSB-PA01729A0Rb	Human	ELISA
<u>IER3 Antibody</u>	CSB-PA011002LA01HU	Human	ELISA, IF
<u>IFI27L2 Antibody</u>	CSB-PA884454LA01HU	Human	ELISA, IHC, IF
<u>IFI44 Antibody</u>	CSB-PA011014LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>IFI44L Antibody</u>	CSB-PA011015LA01HU	Human, Mouse	ELISA, WB, IF
<u>IFIH1 Antibody</u>	CSB-PA880143LA01HU	Human	ELISA, IHC, IF
<u>IFIT1 Antibody</u>	CSB-PA011018LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>IFIT2 Antibody</u>	CSB-PA011021LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>IFIT3 Antibody</u>	CSB-PA011022HA01HU	Human	ELISA, WB, IHC, IF, IP



Product Name	Code	Species Reactivity	Applications
<u>IFIT5 Antibody</u>	CSB-PA011023LA01HU	Human	ELISA, WB, IHC, IF
<u>IFITM1 Antibody</u>	CSB-PA01204A0Rb	Human	ELISA, WB, IF
<u>IFITM2 Antibody</u>	CSB-PA011025LA01HU	Human	ELISA, WB, IHC, IF
<u>IFITM3 Antibody</u>	CSB-PA011026LA01HU	Human	ELISA, WB, IHC, IP
<u>IFITM5 Antibody</u>	CSB-PA011028HA01HU	Human	ELISA, IHC, IF
<u>IFNA1 Antibody</u>	CSB-PA365592HA01HU	Human	ELISA, IHC
<u>IFNA10 Antibody</u>	CSB-PA011033LA01HU	Human, Mouse	ELISA, WB, IHC
<u>IFNA14 Antibody</u>	CSB-PA011035EA01HU	Human, Mouse	ELISA, WB, IHC
<u>IFNA14 Antibody</u>	CSB-PA011035LA01HU	Human	ELISA
<u>IFNA16 Antibody</u>	CSB-PA011036LA01HU	Human	ELISA, IHC
<u>IFNA2 Antibody</u>	CSB-PA06504A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>IFNA21 Antibody</u>	CSB-PA011039EA01HU	Human, Mouse	ELISA, WB, IP
<u>IFNA4 Antibody</u>	CSB-PA011040LA01HU	Human, Mouse	ELISA, WB, IHC
<u>IFNA5 Antibody</u>	CSB-PA011041NA01HU	Human	ELISA, IHC
<u>IFNAG Antibody</u>	CSB-PA344226LA01BO	Bovine	ELISA, WB
<u>IFNAR2 Antibody</u>	CSB-PA11439A0Rb	Human	ELISA, IHC
<u>IFNB1 Antibody</u>	CSB-PA06515A0Rb	Human	ELISA, IHC
<u>Ifne Antibody</u>	CSB-PA800374LA01MO	Mouse	ELISA
<u>IFNG Antibody</u>	CSB-PA011050HA01HU	Human	ELISA, WB
<u>IFNG Antibody</u>	CSB-PA06524A0Rb	Human	ELISA
<u>IFNG Antibody</u>	CSB-PA06525A0Rb	Human	ELISA, WB
<u>IFNG Antibody</u>	CSB-PA14709A0Rb	Guinea Pig	ELISA, WB
<u>Ifng Antibody</u>	CSB-PA15574A0Rb	Rat	ELISA
<u>Ifng Antibody</u>	CSB-PA15587A0Rb	Mouse	ELISA, IHC
<u>IFNG Antibody</u>	CSB-PA15594A0Rb	Chicken	ELISA
<u>IFNG Antibody</u>	CSB-PA011050YA01BO	Bovine	ELISA, WB
<u>IFNGR2 Antibody</u>	CSB-PA011052LA01HU	Human	ELISA, IHC, IF
<u>IFNL1 Antibody</u>	CSB-PA13249A0Rb	Human	ELISA, IHC
<u>IFNL4 Antibody</u>	CSB-PA20629A0Rb	Human	ELISA, IF
<u>ifnlr1 Antibody</u>	CSB-PA29701A0Rb	Zebrafish	ELISA
<u>IFNT1 Antibody</u>	CSB-PA322795LA01BO	Bovine	ELISA
<u>IFNW1 Antibody</u>	CSB-PA011061LA01BO	Bovine	ELISA
<u>IFRD1 Antibody</u>	CSB-PA011069LA01HU	Human	ELISA, IHC
<u>IFRD2 Antibody</u>	CSB-PA011070LA01HU	Human	ELISA, WB, IHC
<u>IFT140 Antibody</u>	CSB-PA842756LA01HU	Human	ELISA, IHC
<u>IFT172 Antibody</u>	CSB-PA887019LA01HU	Human	ELISA, IHC, IF
<u>IFT22 Antibody</u>	CSB-PA863998LA01HU	Human	ELISA, IHC
<u>IFT27 Antibody</u>	CSB-PA874822LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>IFT46 Antibody</u>	CSB-PA873621LA01HU	Human, Mouse	ELISA, WB, IHC
<u>IFT52 Antibody</u>	CSB-PA896883LA01HU	Human	ELISA, WB
<u>IFT74 Antibody</u>	CSB-PA846657LA01HU	Human	ELISA, IHC
<u>IFT80 Antibody</u>	CSB-PA885802LA01HU	Human	ELISA, IHC
<u>IGF1 Antibody</u>	CSB-PA08399A0Rb	Human	ELISA, IHC
<u>IGF1R Antibody</u>	CSB-PA08395A0Rb	Human	ELISA, IF
<u>IGF2BP1 Antibody</u>	CSB-PA878937LA01HU	Human	ELISA, WB, IHC
<u>IGF2BP3 Antibody</u>	CSB-PA011092LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>IGF2R Antibody</u>	CSB-PA011093LA01HU	Human	ELISA, IHC, IF
<u>Igfbp1 Antibody</u>	CSB-PA011095LA01MO	Mouse, Human	ELISA, WB
<u>IGFBP1 Antibody</u>	CSB-PA08419A0Rb	Human	ELISA, WB, IHC
<u>IGFBP6 Antibody</u>	CSB-PA011100LA01HU	Human	ELISA, IHC
<u>IGFBP7 Antibody</u>	CSB-PA16259A0Rb	Human, Rat	ELISA, WB
<u>IgG Antibody</u>	CSB-PA00760E0Gp	Goat	ELISA
<u>IgG Antibody</u>	CSB-PA00770E0Gp	Goat	ELISA, WB
<u>IgG Antibody</u>	CSB-PA00780E0Gp	Human	ELISA
<u>IgG Antibody</u>	CSB-PA00790E0Rb	Human	ELISA, WB
<u>IgG2a Antibody</u>	CSB-PA00820E0Rb	Mouse	ELISA
<u>IGHA2 Antibody</u>	CSB-PA011112LA01HU	Human	ELISA, IHC, IF
<u>IGHD Antibody</u>	CSB-PA011113LA01HU	Human, Mouse	ELISA, WB, IHC
<u>IGHG1 Antibody</u>	CSB-PA011156LA01HU	Human	ELISA, IHC
<u>IGHM Antibody</u>	CSB-PA360788LA01HU	Human	ELISA, WB, IHC
<u>IGHMBP2 Antibody</u>	CSB-PA011172LA01HU	Human	ELISA, IHC
<u>IGKC Antibody</u>	CSB-PA011340LA01HU	Human	ELISA, IHC, IF
<u>IGLL5 Antibody</u>	CSB-PA494948LA01HU	Human	ELISA
<u>IGSF11 Antibody</u>	CSB-PA691547LA01HU	Human, Mouse	ELISA, WB
<u>IGSF6 Antibody</u>	CSB-PA011564LA01HU	Human	ELISA
<u>IGSF9 Antibody</u>	CSB-PA23039A0Rb	Human	ELISA, IHC
<u>IHH Antibody</u>	CSB-PA622779LA01HU	Human	ELISA, IHC, IF
<u>III Antibody</u>	CSB-PA303116LA01ECY	Enterobacteria phage M13	ELISA
<u>IKBIP Antibody</u>	CSB-PA751229LA01HU	Human	ELISA, WB
<u>IKBKB Antibody</u>	CSB-PA16309A0Rb	Human	ELISA, WB, IHC, IF
<u>IKZF1 Antibody</u>	CSB-PA621668LA01HU	Human	ELISA
<u>IKZF1 Antibody</u>	CSB-PA621668HA01HU	Human	ELISA, WB, IHC, IF
<u>IL10 Antibody</u>	CSB-PA011580LA01BO	Bovine	ELISA
<u>IL10 Antibody</u>	CSB-PA011580LA01MOW	Macaca mulatta	ELISA
<u>IL10RA Antibody</u>	CSB-PA621688EA01HU	Human	ELISA, WB, IF
<u>IL10RB Antibody</u>	CSB-PA600091LA01HU	Human, Rat	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>IL12A Antibody</u>	CSB-PA16819A0Rb	Human	ELISA, IHC
<u>IL12B Antibody</u>	CSB-PA011587LA01SH	Sheep	ELISA
<u>IL12B Antibody</u>	CSB-PA011587YA01SH	Sheep	ELISA, WB, IHC
<u>IL12B Antibody</u>	CSB-PA06615A0Rb	Human, Mouse	ELISA, WB, IHC
<u>IL12RB1 Antibody</u>	CSB-PA011588LA01HU	Human	ELISA, WB
<u>IL12RB2 Antibody</u>	CSB-PA859518LA01HU	Human	ELISA, IF
<u>IL13 Antibody</u>	CSB-PA011590LA01HU	Human	ELISA, IHC, IF
<u>IL13RA1 Antibody</u>	CSB-PA011591LA01HU	Human	ELISA, IHC, IF
<u>IL15 Antibody</u>	CSB-PA06635A0Rb	Human	ELISA, IHC
<u>IL16 Antibody</u>	CSB-PA06649A0Rb	Human	ELISA, IF
<u>IL17A Antibody</u>	CSB-PA06659A0Rb	Human	ELISA, IHC, IF
<u>IL17D Antibody</u>	CSB-PA848386LA01HU	Human	ELISA, IHC
<u>IL17F Antibody</u>	CSB-PA842733LA01HU	Human	ELISA, IHC
<u>IL17RA Antibody</u>	CSB-PA12929A0Rb	Human, Mouse	ELISA, WB, IF
<u>IL18 Antibody</u>	CSB-PA08377A0Rb	Human	ELISA
<u>IL18 Antibody</u>	CSB-PA614514YA01HU	Human	ELISA, IHC
<u>IL18 Antibody</u>	CSB-PA895937LA01DO	Dog	ELISA
<u>IL18BP Antibody</u>	CSB-PA16739A0Rb	Human	ELISA, IHC
<u>IL18bp Antibody</u>	CSB-PA16769A0Rb	Mouse	ELISA, IHC
<u>IL18R1 Antibody</u>	CSB-PA623807LA01HU	Human	ELISA, WB, IHC
<u>IL1A Antibody</u>	CSB-PA011613LA01MOW	Macaca mulatta, Human	ELISA, WB
<u>IL1A Antibody</u>	CSB-PA06565A0Rb	Human	ELISA
<u>IL1B Antibody</u>	CSB-PA011614HA01SH	Sheep	ELISA, WB
<u>IL1B Antibody</u>	CSB-PA011614YA01SH	Sheep	ELISA, WB
<u>IL1B Antibody</u>	CSB-PA06574A0Rb	Human	ELISA
<u>IL1b Antibody</u>	CSB-PA17347A0Rb	Mouse	ELISA
<u>IL1F10 Antibody</u>	CSB-PA837869HA01HU	Human	ELISA, IHC
<u>IL1R1 Antibody</u>	CSB-PA011621EA01HU	Human	ELISA, WB, IF
<u>IL1RAPL1 Antibody</u>	CSB-PA011624LA01HU	Human	ELISA, WB, IHC, IF
<u>IL1RL1 Antibody</u>	CSB-PA07209A0Rb	Human	ELISA
<u>IL1RL2 Antibody</u>	CSB-PA888016LA01HU	Human	ELISA, IHC
<u>IL1RN Antibody</u>	CSB-PA06677A0Rb	Human	ELISA, WB
<u>IL2 Antibody</u>	CSB-PA011629LA01HU	Human	ELISA, IHC
<u>IL22RA1 Antibody</u>	CSB-PA822717LA01HU	Human	ELISA, IHC, IF
<u>IL22RA2 Antibody</u>	CSB-PA850247LA01HU	Human	ELISA, IHC
<u>IL23A Antibody</u>	CSB-PA06719A0Rb	Human	ELISA, IHC, IF
<u>IL23R Antibody</u>	CSB-PA733162LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>IL24 Antibody</u>	CSB-PA613378LA01HU	Human	ELISA, IHC
<u>IL2RA Antibody</u>	CSB-PA011649LA01HU	Human, Rat	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>IL2rb Antibody</u>	CSB-PA011650LA01MO	Mouse	ELISA
<u>IL2RG Antibody</u>	CSB-PA011651EA01HU	Human	ELISA, WB, IHC
<u>IL3 Antibody</u>	CSB-PA06725A0Rb	Human	ELISA, IHC
<u>IL32 Antibody</u>	CSB-PA08309A0Rb	Human	ELISA, IHC, IF
<u>IL33 Antibody</u>	CSB-PA011656LA01HU	Human	ELISA, IHC, IF
<u>IL-33 Antibody</u>	CSB-PA011656DA01PI	Pig, Mouse, Human	ELISA, WB
<u>IL34 Antibody</u>	CSB-PA751146LA01HU	Human	ELISA, WB, IHC
<u>IL36RN Antibody</u>	CSB-PA866201HA01HU	Human, Mouse, Rat	ELISA, WB
<u>IL4 Antibody</u>	CSB-PA011659LA01DO	Dog	ELISA
<u>IL4 Antibody</u>	CSB-PA06737A0Rb	Human	ELISA, IHC
<u>IL4 Antibody</u>	CSB-PA011659YA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>IL4I1 Antibody</u>	CSB-PA850426LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>IL6 Antibody</u>	CSB-PA06757A0Rb	Human, Mouse	ELISA, WB, IHC
<u>IL7R Antibody</u>	CSB-PA011670EA01HU	Human	ELISA, IHC
<u>ILDR1 Antibody</u>	CSB-PA768201LA01HU	Human	ELISA, WB, IHC
<u>ILF2 Antibody</u>	CSB-PA00114A0Rb	Human	ELISA, WB, IHC
<u>ILK Antibody</u>	CSB-PA00999A0Rb	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>ILKAP Antibody</u>	CSB-PA011686HA01HU	Human, Mouse, Rat	ELISA, WB
<u>ILVBL Antibody</u>	CSB-PA011687LA01HU	Human	ELISA, WB, IHC, IF
<u>IMMT Antibody</u>	CSB-PA623101LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>IMPACT Antibody</u>	CSB-PA889191LA01HU	Human	ELISA, WB, IHC
<u>IMPDH1 Antibody</u>	CSB-PA22459A0Rb	Human	ELISA, IHC, IF
<u>IMPDH2 Antibody</u>	CSB-PA011701LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>IMPDH2 Antibody</u>	CSB-PA01195A0Rb	Human	ELISA, WB, IHC, IF
<u>IMPG2 Antibody</u>	CSB-PA011704LA01HU	Human	ELISA, IF
<u>IMUP Antibody</u>	CSB-PA003450DA01HU	Human	ELISA
<u>INCA1 Antibody</u>	CSB-PA605764LA01HU	Human	ELISA, IHC
<u>INF2 Antibody</u>	CSB-PA637740LA01HU	Human	ELISA, IHC
<u>ING1 Antibody</u>	CSB-PA891532LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ING3 Antibody</u>	CSB-PA865171LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ING4 Antibody</u>	CSB-PA891792LA01HU	Human	ELISA, IHC
<u>ING5 Antibody</u>	CSB-PA820185LA01HU	Human	ELISA, IF
<u>INHA Antibody</u>	CSB-PA11599A0Rb	Human	ELISA, IHC
<u>inhA Antibody</u>	CSB-PA363781HA01MVZ	Mycobacterium tuberculosis	ELISA



Product Name	Code	Species Reactivity	Applications
<u>INHBA Antibody</u>	CSB-PA07589A0Rb	Human	ELISA, WB, IHC
<u>INHBB Antibody</u>	CSB-PA13929A0Rb	Human	ELISA, WB, IHC, IF
<u>inIA Antibody</u>	CSB-PA758331HA01AAD	Listeria monocytogenes serotype 4b	ELISA, WB
<u>INMT Antibody</u>	CSB-PA011723LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>INO80 Antibody</u>	CSB-PA892466LA01HU	Human	ELISA, WB, IF
<u>INPP1 Antibody</u>	CSB-PA011729LA01HU	Human	ELISA, IHC
<u>INPP4B Antibody</u>	CSB-PA011732HA01HU	Human	ELISA, IHC
<u>INPP5A Antibody</u>	CSB-PA622780HA01HU	Human	ELISA, IHC, IF
<u>INPP5B Antibody</u>	CSB-PA011734LA01HU	Human	ELISA, IHC, IF
<u>INPP5E Antibody</u>	CSB-PA873641LA01HU	Human	ELISA, IHC, IF
<u>INPP5F Antibody</u>	CSB-PA896692NA01HU	Human	ELISA, IF
<u>INS Antibody</u>	CSB-PA12089A0Rb	Human	ELISA, IHC
<u>INSIG1 Antibody</u>	CSB-PA011745LA01HU	Human	ELISA, WB, IHC, IF
<u>INSL4 Antibody</u>	CSB-PA621784LA01HU	Human	ELISA, IHC
<u>INSR Antibody</u>	CSB-PA011753LA01HU	Human	ELISA, IHC, IF
<u>INSRR Antibody</u>	CSB-PA011755LA01HU	Human	ELISA, WB, IHC
<u>Intermediate capsid protein VP6 Antibody</u>	CSB-PA321618LA01PGK	Rotavirus A	ELISA
<u>Intermediate capsid protein VP6 Antibody</u>	CSB-PA541481LA01ROG	Rotavirus A	ELISA
<u>Intermediate capsid protein VP6 Antibody</u>	CSB-PA542165LA01ROH	Rotavirus A	ELISA
<u>INTS10 Antibody</u>	CSB-PA011758EA01HU	Human	ELISA, IHC, IF
<u>INTS13 Antibody</u>	CSB-PA868331LA01HU	Human	ELISA, WB, IP
<u>INTS14 Antibody</u>	CSB-PA847227LA01HU	Human, Mouse	ELISA, WB, IHC
<u>INTS3 Antibody</u>	CSB-PA011761LA01HU	Human	ELISA, IF
<u>INTS6 Antibody</u>	CSB-PA011764LA01HU	Human	ELISA, IF
<u>INTS9 Antibody</u>	CSB-PA011767LA01HU	Human	ELISA, WB, IHC, IF
<u>INTU Antibody</u>	CSB-PA891544LA01HU	Human	ELISA, IHC, IF
<u>Invasin Antibody</u>	CSB-PA321672LA01Y AQ	Yersinia enterocolitica	ELISA
<u>INVS Antibody</u>	CSB-PA011769LA01HU	Human	ELISA, WB, IHC, IF
<u>iolG Antibody</u>	CSB-PA329675HA01BRG	Bacillus subtilis	ELISA, WB
<u>IP6K1 Antibody</u>	CSB-PA011770LA01HU	Human	ELISA, WB
<u>ipaH9.8 Antibody</u>	CSB-PA513237HA01SZC	Shigella flexneri serotype X	ELISA
<u>IPO11 Antibody</u>	CSB-PA887057LA01HU	Human	ELISA



Product Name	Code	Species Reactivity	Applications
<u>IPO13 Antibody</u>	CSB-PA011776LA01HU	Human	ELISA, WB, IHC, IF
<u>IPP Antibody</u>	CSB-PA011782LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>IQCB1 Antibody</u>	CSB-PA621862LA01HU	Human, Rat, Mouse	ELISA, WB, IHC
<u>IQCC Antibody</u>	CSB-PA679698LA01HU	Human	ELISA, IHC
<u>IQCF2 Antibody</u>	CSB-PA815574LA01HU	Human	ELISA
<u>IQCG Antibody</u>	CSB-PA867116LA01HU	Human, Mouse	ELISA, WB, IF
<u>IQGAP2 Antibody</u>	CSB-PA615701LA01HU	Human	ELISA, IHC, IF
<u>IREB2 Antibody</u>	CSB-PA011813LA01HU	Human, Mouse	ELISA, WB, IHC
<u>IRF2BP1 Antibody</u>	CSB-PA815541LA01HU	Human	ELISA, IHC, IF
<u>IRF3 Antibody</u>	CSB-PA011818LA01HU	Human	ELISA, WB, IHC, IF
<u>Irf3 Antibody</u>	CSB-PA011818YA01MO	Mouse	ELISA
<u>Irf5 Antibody</u>	CSB-PA011820LA01MO	Mouse	ELISA
<u>IRF5 Antibody</u>	CSB-PA011820LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>IRF6 Antibody</u>	CSB-PA011821LA01HU	Human, Mouse	ELISA, WB, IF
<u>IRF7 Antibody</u>	CSB-PA011822LA01HU	Human	ELISA, IHC, IF
<u>Irf8 Antibody</u>	CSB-PA011823HA01MO	Mouse, Human	ELISA, WB
<u>IRF9 Antibody</u>	CSB-PA011824LA01HU	Human, Mouse	ELISA, WB, IF
<u>IRS2 Antibody</u>	CSB-PA011830LA01HU	Human	ELISA, IHC, IF
<u>IRX3 Antibody</u>	CSB-PA011835LA01HU	Human, Rat	ELISA, WB
<u>IRX4 Antibody</u>	CSB-PA011836LA01HU	Human	ELISA, WB, IF
<u>IRX6 Antibody</u>	CSB-PA011838LA01HU	Human, Rat	ELISA, WB, IHC
<u>ISCA1 Antibody</u>	CSB-PA887126LA01HU	Human	ELISA, IHC
<u>ISCA2 Antibody</u>	CSB-PA769773LA01HU	Human	ELISA, WB, IHC
<u>ISCU Antibody</u>	CSB-PA887955DA01HU	Human	ELISA
<u>ISG15 Antibody</u>	CSB-PA09744A0Rb	Human	ELISA, WB, IP
<u>isl1 Antibody</u>	CSB-PA011846LA01DIL	Zebrafish	ELISA
<u>ISL1 Antibody</u>	CSB-PA011846LA01HU	Human	ELISA, IHC
<u>ISOC1 Antibody</u>	CSB-PA839302LA01HU	Human	ELISA, WB, IHC
<u>ISPD Antibody</u>	CSB-PA388827LA01HU	Human	ELISA, IHC, IF
<u>IST1 Antibody</u>	CSB-PA012171LA01HU	Human	ELISA, IHC, IF
<u>ISX Antibody</u>	CSB-PA011854LA01HU	Human	ELISA, IHC
<u>ITCH Antibody</u>	CSB-PA839346HA01HU	Human	ELISA, WB, IF
<u>ITGA1 Antibody</u>	CSB-PA011861LA01HU	Human	ELISA, IHC, IF
<u>ITGA2 Antibody</u>	CSB-PA011864LA01HU	Human	ELISA, IHC, IF
<u>ITGA2B Antibody</u>	CSB-PA011865LA01HU	Human	ELISA, IHC, IF
<u>ITGA3 Antibody</u>	CSB-PA011866LA01HU	Human	ELISA, IHC, IF
<u>ITGA4 Antibody</u>	CSB-PA011867LA01HU	Human	ELISA, IHC, IF
<u>ITGA5 Antibody</u>	CSB-PA011868LA01HU	Human, Mouse	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>ITGA6 Antibody</u>	CSB-PA011869LA01HU	Human, Rat, Mouse	ELISA, WB, IF
<u>ITGA7 Antibody</u>	CSB-PA618797LA01HU	Human	ELISA, IF
<u>ITGA8 Antibody</u>	CSB-PA011871NA01HU	Human	ELISA, IHC
<u>ITGAL Antibody</u>	CSB-PA011875LA01HU	Human	ELISA, WB, IHC, IF
<u>Itgal Antibody</u>	CSB-PA011875LA01MO	Mouse	ELISA
<u>ITGAM Antibody</u>	CSB-PA011876LA01HU	Human	ELISA, IHC, IF
<u>ITGAV Antibody</u>	CSB-PA011877LA01HU	Human	ELISA, WB, IHC, IF
<u>Itgb1 Antibody</u>	CSB-PA011880LA01MO	Mouse, Human	ELISA, WB, IF
<u>Itgb1 Antibody</u>	CSB-PA011880LA01RA	Rat	ELISA
<u>ITGB1BP1 Antibody</u>	CSB-PA011881LA01HU	Human	ELISA, WB, IHC
<u>ITGB1BP2 Antibody</u>	CSB-PA892333LA01HU	Human	ELISA, IHC
<u>ITGB2 Antibody</u>	CSB-PA011884LA01HU	Human	ELISA, IHC, IF
<u>ITGB4 Antibody</u>	CSB-PA011887LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ITGB5 Antibody</u>	CSB-PA011889LA01HU	Human	ELISA, IHC
<u>ITGB6 Antibody</u>	CSB-PA011890LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>ITIH1 Antibody</u>	CSB-PA4779A0Rb	Human	ELISA, IHC
<u>ITIH3 Antibody</u>	CSB-PA22869A0Rb	Human	ELISA, IHC
<u>ITIH4 Antibody</u>	CSB-PA011897LA01HU	Human	ELISA, IHC
<u>ITIH4 Antibody</u>	CSB-PA011897YA01HU	Human	ELISA, IHC
<u>ITIH5 Antibody</u>	CSB-PA768213LA01HU	Human	ELISA, IHC
<u>ITM2A Antibody</u>	CSB-PA011903LA01HU	Human	ELISA, WB, IHC, IF
<u>ITM2B Antibody</u>	CSB-PA011904LA01HU	Human	ELISA, IHC, IF
<u>ITPR1 Antibody</u>	CSB-PA614538LA01HU	Human	ELISA, IHC, IF
<u>ITPRID2 Antibody</u>	CSB-PA022710LA01HU	Human	ELISA, IHC
<u>ITSN1 Antibody</u>	CSB-PA621885LA01HU	Human	ELISA, WB, IHC
<u>ITSN2 Antibody</u>	CSB-PA873709LA01HU	Human	ELISA, WB, IF
<u>IVD Antibody</u>	CSB-PA011921LA01HU	Human	ELISA, IHC
<u>IVL Antibody</u>	CSB-PA011922LA01HU	Human	ELISA, WB, IHC, IF
<u>IVNS1ABP Antibody</u>	CSB-PA896937LA01HU	Human	ELISA, IHC, IF
<u>IX Antibody</u>	CSB-PA301675LA01ECY	Enterobacteria phage M13	ELISA
<u>IYD Antibody</u>	CSB-PA754272LA01HU	Human	ELISA, IF
<u>IZUMO1 Antibody</u>	CSB-PA812903LA01HU	Human	ELISA, IHC, IF
<u>IZUMO4 Antibody</u>	CSB-PA636327LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>JAG1 Antibody</u>	CSB-PA01949A0Rb	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>JAGN1 Antibody</u>	CSB-PA818709LA01HU	Human	ELISA, WB, IHC, IF
<u>JAK2 Antibody</u>	CSB-PA011931LA01HU	Human	ELISA
<u>JAK3 Antibody</u>	CSB-PA011932LA01HU	Human	ELISA
<u>JAM2 Antibody</u>	CSB-PA011936LA01HU	Human	ELISA, WB
<u>JAM3 Antibody</u>	CSB-PA874828LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>JARID2 Antibody</u>	CSB-PA852915HA01HU	Human	ELISA, IHC, IF
<u>JMJD1C Antibody</u>	CSB-PA624030HA01HU	Human	ELISA, IHC
<u>JMY Antibody</u>	CSB-PA850838LA01HU	Human	ELISA, IHC, IF
<u>JPH1 Antibody</u>	CSB-PA872553LA01HU	Human	ELISA, IHC, IF
<u>JPT1 Antibody</u>	CSB-PA890732HA01HU	Human, Rat	ELISA, WB
<u>JRK Antibody</u>	CSB-PA525639NA01HU	Human, Mouse	ELISA, WB
<u>JUN Antibody</u>	CSB-PA16929A0Rb	Human, Mouse, Rat	ELISA, WB, IF
<u>JUND Antibody</u>	CSB-PA04469A0Rb	Human	ELISA, WB
<u>JUP Antibody</u>	CSB-PA12839A0Rb	Human	ELISA, WB, IHC
<u>KALRN Antibody</u>	CSB-PA011981LA01HU	Human	ELISA, IHC, IF
<u>KANK1 Antibody</u>	CSB-PA614803LA01HU	Human	ELISA, IHC
<u>KANK2 Antibody</u>	CSB-PA720638LA01HU	Human	ELISA, IHC
<u>KANSL1 Antibody</u>	CSB-PA770344LA01HU	Human	ELISA, IHC, IF
<u>KARS Antibody</u>	CSB-PA04689A0Rb	Human	ELISA, WB, IHC, IF
<u>KAT2A Antibody</u>	CSB-PA821706LA01HU	Human	ELISA, IHC
<u>KAT2B Antibody</u>	CSB-PA835706LA01HU	Human	ELISA, IHC
<u>KAT5 Antibody</u>	CSB-PA17849A0Rb	Human	ELISA, WB, IHC, IF
<u>KAT7 Antibody</u>	CSB-PA03865A0Rb	Human, Rat	ELISA, WB, IHC
<u>KAT8 Antibody</u>	CSB-PA17859A0Rb	Human, Mouse, Rat	ELISA, WB
<u>KATNAL2 Antibody</u>	CSB-PA809022LA01HU	Human	ELISA, WB, IHC
<u>KBTBD2 Antibody</u>	CSB-PA811624LA01HU	Human	ELISA, WB
<u>KBTBD4 Antibody</u>	CSB-PA882119LA01HU	Human, Mouse	ELISA, WB, IHC
<u>KBTBD6 Antibody</u>	CSB-PA774812LA01HU	Human	ELISA, WB, IHC, IF
<u>KCMF1 Antibody</u>	CSB-PA873718LA01HU	Human	ELISA, IHC
<u>KCNA1 Antibody</u>	CSB-PA012005HA01HU	Human	ELISA
<u>KCNA1 Antibody</u>	CSB-PA012005LA01HU	Human, Mouse, Rat	ELISA, WB, IF
<u>KCNA6 Antibody</u>	CSB-PA012011LA01HU	Human	ELISA, IHC
<u>KCNA7 Antibody</u>	CSB-PA847218LA01HU	Human	ELISA, IHC
<u>KCNAB1 Antibody</u>	CSB-PA012013LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>KCNB1 Antibody</u>	CSB-PA613519LA01HU	Human	ELISA, IF
<u>KCNC1 Antibody</u>	CSB-PA012018LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>KCNC3 Antibody</u>	CSB-PA621758LA01HU	Human	ELISA, IHC, IF
<u>KCND1 Antibody</u>	CSB-PA012022LA01HU	Human	ELISA, IHC, IF
<u>KCNE1 Antibody</u>	CSB-PA012025LA01HU	Human	ELISA, IHC, IF
<u>Kcne2 Antibody</u>	CSB-PA012027LA01RA	Rat	ELISA
<u>KCNE3 Antibody</u>	CSB-PA22569A0Rb	Human	ELISA, IHC, IF
<u>KCNG4 Antibody</u>	CSB-PA844728LA01HU	Human	ELISA, IHC, IF
<u>KCNH2 Antibody</u>	CSB-PA619634LA01HU	Human	ELISA, IHC, IF
<u>KCNH7 Antibody</u>	CSB-PA878877LA01HU	Human	ELISA, IHC, IF
<u>KCNIP3 Antibody</u>	CSB-PA012045LA01HU	Human	ELISA, WB, IHC, IF
<u>KCNIP4 Antibody</u>	CSB-PA750882LA01HU	Human	ELISA, IHC, IF
<u>KCNJ1 Antibody</u>	CSB-PA012047LA01HU	Human	ELISA, WB, IF
<u>KCNJ10 Antibody</u>	CSB-PA012048LA01HU	Human	ELISA, IHC, IF
<u>KCNJ11 Antibody</u>	CSB-PA622781LA01HU	Human, Mouse	ELISA, WB, IHC
<u>KCNJ3 Antibody</u>	CSB-PA012056LA01HU	Human	ELISA, IHC, IF
<u>KCNJ5 Antibody</u>	CSB-PA012058LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>KCNJ8 Antibody</u>	CSB-PA613604LA01HU	Human	ELISA, IHC, IF
<u>KCNK1 Antibody</u>	CSB-PA012062LA01HU	Human	ELISA, IHC, IF
<u>KCNK13 Antibody</u>	CSB-PA862058LA01HU	Human	ELISA, WB, IHC, IF
<u>KCNK15 Antibody</u>	CSB-PA867147LA01HU	Human	ELISA, IF
<u>KCNK16 Antibody</u>	CSB-PA850434LA01HU	Human	ELISA, IHC
<u>KCNK3 Antibody</u>	CSB-PA012071LA01HU	Human	ELISA, IHC, IF
<u>KCNK9 Antibody</u>	CSB-PA889068LA01HU	Human	ELISA, IF
<u>KCNMA1 Antibody</u>	CSB-PA614255LA01HU	Human	ELISA, IHC, IF
<u>KCNMB1 Antibody</u>	CSB-PA621958LA01HU	Human	ELISA, IF
<u>KCNMB3 Antibody</u>	CSB-PA878838LA01HU	Human	ELISA, IF
<u>KCNN3 Antibody</u>	CSB-PA880079LA01HU	Human	ELISA, IF
<u>KCNQ3 Antibody</u>	CSB-PA012091LA01HU	Human	ELISA, WB, IF
<u>KCNQ4 Antibody</u>	CSB-PA012092LA01HU	Human	ELISA, IF
<u>KCNRG Antibody</u>	CSB-PA847638LA01HU	Human	ELISA, IHC
<u>KCNS2 Antibody</u>	CSB-PA892470LA01HU	Human, Mouse	ELISA, WB, IF
<u>KCNT1 Antibody</u>	CSB-PA686387LA01HU	Human	ELISA, IHC, IF
<u>KCNV1 Antibody</u>	CSB-PA747493LA01HU	Human	ELISA, IF
<u>KCNV2 Antibody</u>	CSB-PA848410LA01HU	Human	ELISA, IHC, IF
<u>KCTD1 Antibody</u>	CSB-PA744477LA01HU	Human	ELISA
<u>KCTD16 Antibody</u>	CSB-PA727857LA01HU	Human	ELISA, IHC
<u>KCTD18 Antibody</u>	CSB-PA757654LA01HU	Human	ELISA
<u>KCTD19 Antibody</u>	CSB-PA624192LA01HU	Human	ELISA, IHC
<u>KCTD4 Antibody</u>	CSB-PA837857LA01HU	Human	ELISA, IHC
<u>KCTD6 Antibody</u>	CSB-PA822756LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>KCTD7 Antibody</u>	CSB-PA850395LA01HU	Human	ELISA, IF
<u>KCTD9 Antibody</u>	CSB-PA741970LA01HU	Human	ELISA, WB
<u>KDM1A Antibody</u>	CSB-PA525314HA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>KDM1B Antibody</u>	CSB-PA822748LA01HU	Human	ELISA, IHC
<u>KDM2B Antibody</u>	CSB-PA836777LA01HU	Human	ELISA, IF
<u>KDM4A Antibody</u>	CSB-PA012137LA01HU	Human	ELISA, IHC, IF
<u>KDM4C Antibody</u>	CSB-PA872489HA01HU	Human	ELISA, IHC
<u>KDM4D Antibody</u>	CSB-PA718355LA01HU	Human	ELISA, IHC, IF
<u>KDM5A Antibody</u>	CSB-PA012141LA01HU	Human	ELISA, IHC, IF
<u>KDM5C Antibody</u>	CSB-PA011939HA01HU	Human	ELISA, IF
<u>KDM6A Antibody</u>	CSB-PA025774LA01HU	Human	ELISA, IHC
<u>KDM6B Antibody</u>	CSB-PA517674LA01HU	Human	ELISA, IF
<u>KDM8 Antibody</u>	CSB-PA839789LA01HU	Human	ELISA, IHC
<u>KDR Antibody</u>	CSB-PA07279A0Rb	Human	ELISA, IHC
<u>KEAP1 Antibody</u>	CSB-PA012147LA01HU	Human, Mouse	ELISA, WB, IHC
<u>Keap1 Antibody</u>	CSB-PA012147LA01MO	Mouse, Human	ELISA, WB
<u>KEL Antibody</u>	CSB-PA012148LA01HU	Human	ELISA, WB, IHC
<u>kgp Antibody</u>	CSB-PA464342LA01EXZ	Porphyromonas gingivalis	ELISA, WB
<u>KHDC1 Antibody</u>	CSB-PA709298LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>KHDC4 Antibody</u>	CSB-PA770364LA01HU	Human, Mouse	ELISA, WB, IHC
<u>KHDRBS2 Antibody</u>	CSB-PA730139LA01HU	Human, Mouse	ELISA, WB, IHC
<u>KIAA0319 Antibody</u>	CSB-PA22689A0Rb	Human	ELISA, IHC, IF
<u>KIAA0319L Antibody</u>	CSB-PA816919LA01HU	Human	ELISA, IHC, IF
<u>KIAA0513 Antibody</u>	CSB-PA012195LA01HU	Human, Mouse	ELISA, WB, IF
<u>KIAA1191 Antibody</u>	CSB-PA846582LA01HU	Human, Mouse	ELISA, WB, IHC
<u>KIAA1257 Antibody</u>	CSB-PA891961LA01HU	Human	ELISA, IHC
<u>KIAA1324L Antibody</u>	CSB-PA012251LA01HU	Human	ELISA, IHC
<u>KIAA1841 Antibody</u>	CSB-PA750733LA01HU	Human	ELISA
<u>KIDINS220 Antibody</u>	CSB-PA012306LA01HU	Human	ELISA, WB, IHC, IF
<u>KIF11 Antibody</u>	CSB-PA012307LA01HU	Human, Mouse	ELISA, WB, IF
<u>KIF12 Antibody</u>	CSB-PA850294LA01HU	Human	ELISA, IHC
<u>KIF13A Antibody</u>	CSB-PA880944LA01HU	Human	ELISA, IHC
<u>KIF13B Antibody</u>	CSB-PA873625LA01HU	Human	ELISA, IHC, IF
<u>KIF14 Antibody</u>	CSB-PA613582LA01HU	Human	ELISA, IHC, IF
<u>KIF15 Antibody</u>	CSB-PA868310LA01HU	Human	ELISA, IF
<u>KIF17 Antibody</u>	CSB-PA865208LA01HU	Human	ELISA, IHC, IF
<u>KIF18A Antibody</u>	CSB-PA847753LA01HU	Human	ELISA, IF
<u>Kif1a Antibody</u>	CSB-PA012318LA01MO	Mouse	ELISA
<u>KIF1BP Antibody</u>	CSB-PA836204LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>Kif1c Antibody</u>	CSB-PA012320LA01RA	Rat	ELISA
<u>KIF20A Antibody</u>	CSB-PA012321LA01HU	Human	ELISA, IHC, IP
<u>KIF20B Antibody</u>	CSB-PA839399LA01HU	Human	ELISA, IHC, IF
<u>KIF21B Antibody</u>	CSB-PA012324LA01HU	Human	ELISA, IHC
<u>KIF22 Antibody</u>	CSB-PA619894LA01HU	Human	ELISA, IHC, IF
<u>KIF23 Antibody</u>	CSB-PA23569A0Rb	Human	ELISA, IHC
<u>KIF24 Antibody</u>	CSB-PA716354LA01HU	Human	ELISA, IHC
<u>KIF25 Antibody</u>	CSB-PA883440LA01HU	Human	ELISA, IHC
<u>KIF26B Antibody</u>	CSB-PA642548LA01HU	Human	ELISA, IHC
<u>KIF2A Antibody</u>	CSB-PA012332LA01HU	Human	ELISA, IHC
<u>KIF3C Antibody</u>	CSB-PA012337LA01HU	Human	ELISA, WB, IHC, IF
<u>KIF5B Antibody</u>	CSB-PA012342LA01HU	Human	ELISA, WB, IHC
<u>KIFAP3 Antibody</u>	CSB-PA012347LA01HU	Human	ELISA, IHC, IF
<u>KIFC1 Antibody</u>	CSB-PA871587LA01HU	Human	ELISA, IF
<u>KIR2DL1 Antibody</u>	CSB-PA012352LA01HU	Human	ELISA, IHC, IF
<u>KIR2DL2 Antibody</u>	CSB-PA012353LA01HU	Human	ELISA, IF
<u>KIR2DL3 Antibody</u>	CSB-PA012354LA01HU	Human	ELISA, IHC, IF
<u>KIR2DS1 Antibody</u>	CSB-PA623946LA01HU	Human	ELISA
<u>KIR2DS2 Antibody</u>	CSB-PA012360LA01HU	Human	ELISA, IHC, IF
<u>KIR2DS3 Antibody</u>	CSB-PA613530LA01HU	Human	ELISA, IF
<u>KIR2DS5 Antibody</u>	CSB-PA619902LA01HU	Human	ELISA, IF
<u>KIR3DL1 Antibody</u>	CSB-PA012364LA01HU	Human	ELISA, WB, IHC, IF
<u>KIR3DL2 Antibody</u>	CSB-PA012365LA01HU	Human	ELISA, IF
<u>KIR3DS1 Antibody</u>	CSB-PA614817LA01HU	Human	ELISA, IF
<u>KIRREL2 Antibody</u>	CSB-PA754407LA01HU	Human, Mouse	ELISA, WB
<u>KIRREL3 Antibody</u>	CSB-PA811640LA01HU	Human	ELISA, IF
<u>KISS1 Antibody</u>	CSB-PA621883LA01HU	Human	ELISA, IHC
<u>KIT Antibody</u>	CSB-PA07109A0Rb	Human	ELISA, IF
<u>KL Antibody</u>	CSB-PA12729A0Rb	Human	ELISA, IHC, IF
<u>KLB Antibody</u>	CSB-PA769815LA01HU	Human	ELISA, IHC, IF
<u>KLC4 Antibody</u>	CSB-PA012382LA01HU	Human	ELISA
<u>KLF1 Antibody</u>	CSB-PA619760LA01HU	Human	ELISA, IF
<u>KLF11 Antibody</u>	CSB-PA012385LA01HU	Human	ELISA, IHC
<u>KLF12 Antibody</u>	CSB-PA897495LA01HU	Human, Mouse, Rat	ELISA, WB
<u>KLF13 Antibody</u>	CSB-PA897476LA01HU	Human	ELISA, WB, IHC
<u>KLF2 Antibody</u>	CSB-PA896743LA01HU	Human	ELISA, IF
<u>KLF5 Antibody</u>	CSB-PA613429LA01HU	Human	ELISA, WB
<u>KLF7 Antibody</u>	CSB-PA012397LA01HU	Human, Mouse	ELISA, WB
<u>KLHDC8A Antibody</u>	CSB-PA816911LA01HU	Human, Mouse	ELISA, WB, IHC
<u>KLHL12 Antibody</u>	CSB-PA700650LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>KLHL2 Antibody</u>	CSB-PA012421LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>KLHL20 Antibody</u>	CSB-PA22879A0Rb	Human	ELISA, IHC, IF
<u>KLHL22 Antibody</u>	CSB-PA684469LA01HU	Human, Mouse	ELISA, WB
<u>KLHL23 Antibody</u>	CSB-PA822749LA01HU	Human, Mouse	ELISA, WB
<u>KLHL3 Antibody</u>	CSB-PA890684HA01HU	Human	ELISA, IHC
<u>KLHL34 Antibody</u>	CSB-PA850783LA01HU	Human	ELISA, IHC
<u>KLHL38 Antibody</u>	CSB-PA652063LA01HU	Human	ELISA, IHC
<u>KLHL41 Antibody</u>	CSB-PA011995LA01HU	Human	ELISA
<u>KLHL42 Antibody</u>	CSB-PA885803LA01HU	Human	ELISA, IHC, IF
<u>KLK1 Antibody</u>	CSB-PA012446LA01HU	Human	ELISA, IHC, IF
<u>KLK10 Antibody</u>	CSB-PA14039A0Rb	Human	ELISA
<u>KLK14 Antibody</u>	CSB-PA868375LA01HU	Human	ELISA, IHC
<u>Klk1b22 Antibody</u>	CSB-PA325391LA01MO	Mouse	ELISA
<u>KLK2 Antibody</u>	CSB-PA012453LA01HU	Human	ELISA, IHC
<u>KLK3 Antibody</u>	CSB-PA06779A0Rb	Human	ELISA, IHC
<u>KLK5 Antibody</u>	CSB-PA897100LA01HU	Human	ELISA, WB, IHC
<u>KLK6 Antibody</u>	CSB-PA842165LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>Klk7 Antibody</u>	CSB-PA012458LA01RA	Rat	ELISA
<u>KLK8 Antibody</u>	CSB-PA012459LA01HU	Human	ELISA, IHC, IF
<u>Klkb1 Antibody</u>	CSB-PA012461LA01RA	Rat, Human	ELISA, WB
<u>Klkb1 Antibody</u>	CSB-PA012461YA01RA	Rat	ELISA
<u>KLKB1 Antibody</u>	CSB-PA012461LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>KLRC2 Antibody</u>	CSB-PA012467LA01HU	Human	ELISA, IHC
<u>KLRC3 Antibody</u>	CSB-PA23099A0Rb	Human	ELISA, IHC, IF
<u>KLRD1 Antibody</u>	CSB-PA615672LA01HU	Human	ELISA, IHC, IF
<u>KLRG1 Antibody</u>	CSB-PA846610LA01HU	Human	ELISA, IHC
<u>KLRK1 Antibody</u>	CSB-PA012474LA01HU	Human	ELISA, WB, IHC, IF
<u>KMO Antibody</u>	CSB-PA012475LA01HU	Human	ELISA, IHC, IF
<u>KMT2B Antibody</u>	CSB-PA890939LA01HU	Human	ELISA, IF
<u>KMT2E Antibody</u>	CSB-PA818269LA01HU	Human, Mouse	ELISA, WB
<u>KMT5A Antibody</u>	CSB-PA889082LA01HU	Human	ELISA, IHC
<u>KMT5B Antibody</u>	CSB-PA674569LA01HU	Human	ELISA, WB
<u>Kng1 Antibody</u>	CSB-PA012479LA01RA	Rat	ELISA
<u>KNG1 Antibody</u>	CSB-PA14647A0Rb	Human	ELISA, WB, IHC, IF
<u>KP6 killer toxin Antibody</u>	CSB-PA325875LA01UBC	Ustilago maydis P6 virus	ELISA
<u>KPNA1 Antibody</u>	CSB-PA012483LA01HU	Human	ELISA, IHC
<u>KPNA2 Antibody</u>	CSB-PA012484LA01HU	Human	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>KPNA4 Antibody</u>	CSB-PA012486LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>KPNB1 Antibody</u>	CSB-PA622929LA01HU	Human	ELISA, WB, IHC
<u>KRAS Antibody</u>	CSB-PA13969A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>KRAS Antibody</u>	CSB-PA13967A0Rb	Human	ELISA, IHC, IF
<u>KRBOX4 Antibody</u>	CSB-PA711455LA01HU	Human, Mouse	ELISA, WB, IF
<u>KRCC1 Antibody</u>	CSB-PA868279LA01HU	Human	ELISA, IF
<u>KRIT1 Antibody</u>	CSB-PA012501LA01HU	Human	ELISA, IHC, IF
<u>KRT1 Antibody</u>	CSB-PA012503LA01HU	Human	ELISA, IHC, IF
<u>KRT10 Antibody</u>	CSB-PA012504LA01HU	Human	ELISA, WB, IHC, IF
<u>KRT12 Antibody</u>	CSB-PA859510LA01HU	Human	ELISA, IF
<u>KRT13 Antibody</u>	CSB-PA012511LA01HU	Human	ELISA, IHC
<u>KRT14 Antibody</u>	CSB-PA012512LA01HU	Human	ELISA, IHC, IF
<u>KRT16 Antibody</u>	CSB-PA012515LA01HU	Human	ELISA, WB, IHC, IF
<u>KRT17 Antibody</u>	CSB-PA012516LA01BO	Bovine	ELISA
<u>KRT18 Antibody</u>	CSB-PA10629A0Rb	Human, Zebrafish	ELISA, WB, IHC, IF
<u>Krt18 Antibody</u>	CSB-PA012519LA01MO	Mouse, Human	ELISA, WB, IHC
<u>KRT19 Antibody</u>	CSB-PA11259A0Rb	Human	ELISA
<u>KRT20 Antibody</u>	CSB-PA012536LA01HU	Human	ELISA, WB, IHC, IF
<u>KRT33B Antibody</u>	CSB-PA621778LA01HU	Human	ELISA, WB, IHC
<u>KRT35 Antibody</u>	CSB-PA856422LA01HU	Human	ELISA, WB, IHC
<u>KRT38 Antibody</u>	CSB-PA012554LA01HU	Human	ELISA, WB
<u>KRT4 Antibody</u>	CSB-PA012556LA01HU	Human	ELISA, IHC, IF
<u>KRT40 Antibody</u>	CSB-PA718197LA01HU	Human	ELISA, IHC
<u>KRT5 Antibody</u>	CSB-PA11689A0Rb	Human	ELISA, IF
<u>KRT6A Antibody</u>	CSB-PA012561EA01HU	Human	ELISA, WB, IHC
<u>KRT6A Antibody</u>	CSB-PA012561LA01HU	Human	ELISA, WB, IHC, IF
<u>KRT7 Antibody</u>	CSB-PA012564LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>KRT71 Antibody</u>	CSB-PA660951LA01HU	Human	ELISA, IHC
<u>KRT72 Antibody</u>	CSB-PA619911LA01HU	Human	ELISA, WB
<u>KRT73 Antibody</u>	CSB-PA774833LA01HU	Human	ELISA, WB
<u>KRT76 Antibody</u>	CSB-PA012570LA01HU	Human	ELISA, WB, IHC, IF
<u>KRT77 Antibody</u>	CSB-PA767211LA01HU	Human	ELISA, IHC
<u>KRT79 Antibody</u>	CSB-PA719894LA01HU	Human, Mouse	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>KRT8 Antibody</u>	CSB-PA11269A0Rb	Human	ELISA, WB, IHC, IF
<u>KRT80 Antibody</u>	CSB-PA760849LA01HU	Human, Mouse	ELISA, WB, IHC
<u>KRT81 Antibody</u>	CSB-PA012576LA01HU	Human, Mouse	ELISA, WB, IHC
<u>KRTAP11-1 Antibody</u>	CSB-PA810266LA01HU	Human	ELISA
<u>KRTAP13-1 Antibody</u>	CSB-PA815543LA01HU	Human	ELISA
<u>KRTAP13-3 Antibody</u>	CSB-PA666167LA01HU	Human	ELISA
<u>KRTAP13-4 Antibody</u>	CSB-PA667309LA01HU	Human	ELISA
<u>KRTAP19-5 Antibody</u>	CSB-PA659465LA01HU	Human	ELISA
<u>KRTAP23-1 Antibody</u>	CSB-PA012641LA01HU	Human	ELISA, IHC
<u>KRTAP26-1 Antibody</u>	CSB-PA764791LA01HU	Human	ELISA
<u>KRTAP3-1 Antibody</u>	CSB-PA863160LA01HU	Human	ELISA
<u>KRTAP3-2 Antibody</u>	CSB-PA871617LA01HU	Human	ELISA, IHC
<u>KRTAP3-3 Antibody</u>	CSB-PA874846LA01HU	Human	ELISA
<u>KRTAP8-1 Antibody</u>	CSB-PA814201LA01HU	Human	ELISA, IHC
<u>KSR1 Antibody</u>	CSB-PA012695LA01HU	Human	ELISA, WB, IF
<u>KSR2 Antibody</u>	CSB-PA012696LA01HU	Human	ELISA, IHC, IF
<u>KTI12 Antibody</u>	CSB-PA839314LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>KTNI Antibody</u>	CSB-PA773030LA01HU	Human	ELISA, WB, IHC, IP
<u>KY Antibody</u>	CSB-PA012702LA01HU	Human	ELISA, WB, IHC, IF
<u>KYAT1 Antibody</u>	CSB-PA613699LA01HU	Human	ELISA, IHC
<u>KYAT3 Antibody</u>	CSB-PA740996LA01HU	Human	ELISA, IHC
<u>KYNU Antibody</u>	CSB-PA621970LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>L Antibody</u>	CSB-PA806333LA01VBE	Vesicular stomatitis Indiana virus	ELISA, WB
<u>L1 Antibody</u>	CSB-PA08349A0Rb	Human papillomavirus type 16	ELISA
<u>L1 Antibody</u>	CSB-PA08359A0Rb	Human papillomavirus type 18	ELISA, WB
<u>L1CAM Antibody</u>	CSB-PA07439A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>L1TD1 Antibody</u>	CSB-PA729208LA01HU	Human	ELISA, IHC
<u>L2 Antibody</u>	CSB-PA356939LA01HMN	Human papillomavirus type 18	ELISA
<u>L2 Antibody</u>	CSB-PA365848LA01HML	Human papillomavirus type 16	ELISA
<u>L2HGDH Antibody</u>	CSB-PA012709LA01BO	Bovine	ELISA
<u>L2HGDH Antibody</u>	CSB-PA864008LA01HU	Human	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>L3HYPDH Antibody</u>	CSB-PA846613LA01HU	Human	ELISA, WB, IHC
<u>L3MBTL2 Antibody</u>	CSB-PA839279LA01HU	Human, Rat	ELISA, WB, IF
<u>L3MBTL4 Antibody</u>	CSB-PA836703LA01HU	Human, Rat	ELISA, WB, IHC
<u>LACTB2 Antibody</u>	CSB-PA706635LA01HU	Human	ELISA, IHC, IF
<u>Lactoferrin Antibody</u>	CSB-PA00860E0Rb	Sheep	ELISA
<u>lacZ Antibody</u>	CSB-PA009476LA01ENV	Escherichia coli	ELISA, WB
<u>LAG3 Antibody</u>	CSB-PA012719LA01HU	Human	ELISA, WB, IHC, IF
<u>LAIR1 Antibody</u>	CSB-PA012721LA01HU	Human	ELISA, IHC, IF
<u>LAIR2 Antibody</u>	CSB-PA012722LA01HU	Human	ELISA, IHC
<u>LAMA2 Antibody</u>	CSB-PA012726LA01HU	Human	ELISA, IHC, IF
<u>LAMA3 Antibody</u>	CSB-PA623096LA01HU	Human	ELISA, IHC
<u>LAMA4 Antibody</u>	CSB-PA621948LA01HU	Human	ELISA, WB, IF
<u>LAMA5 Antibody</u>	CSB-PA012729LA01HU	Human	ELISA, IHC
<u>LAMB1 Antibody</u>	CSB-PA012730LA01HU	Human	ELISA, IHC
<u>LAMB2 Antibody</u>	CSB-PA14279A0Rb	Human	ELISA
<u>LAMB3 Antibody</u>	CSB-PA615717LA01HU	Human	ELISA, IHC
<u>LAMB4 Antibody</u>	CSB-PA012734LA01HU	Human	ELISA, IHC, IF
<u>LAMC1 Antibody</u>	CSB-PA11477A0Rb	Human	ELISA, IHC, IF
<u>LAMC2 Antibody</u>	CSB-PA618803LA01HU	Human	ELISA, WB, IHC
<u>LAMP2 Antibody</u>	CSB-PA012740EA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>LAMP3 Antibody</u>	CSB-PA892177LA01HU	Human	ELISA, WB, IHC
<u>LAMTOR1 Antibody</u>	CSB-PA757083LA01HU	Human	ELISA, IHC, IF
<u>LANCL1 Antibody</u>	CSB-PA012742LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>LANCL2 Antibody</u>	CSB-PA865129LA01HU	Human	ELISA, WB, IHC, IF
<u>lapB Antibody</u>	CSB-PA359577HA01ENV	Escherichia coli	ELISA, WB
<u>LAPTM4A Antibody</u>	CSB-PA621858LA01HU	Human	ELISA, IHC, IF
<u>LAPTM4B Antibody</u>	CSB-PA774814LA01HU	Human	ELISA, WB, IHC, IF
<u>LAPTM5 Antibody</u>	CSB-PA622667LA01HU	Human	ELISA, IHC
<u>LARP1B Antibody</u>	CSB-PA734222HA01HU	Human	ELISA
<u>LARP6 Antibody</u>	CSB-PA887112LA01HU	Human	ELISA, IHC
<u>LARP7 Antibody</u>	CSB-PA672326LA01HU	Human	ELISA, IHC
<u>LARS Antibody</u>	CSB-PA873736LA01HU	Human	ELISA, WB, IHC, IF
<u>Larval cuticle protein A1A Antibody</u>	CSB-PA304533LA01TBF	Tenebrio molitor	ELISA
<u>LAS1L Antibody</u>	CSB-PA897285LA01HU	Human	ELISA, WB, IF
<u>lasA Antibody</u>	CSB-PA318588LA01EZC	Pseudomonas aeruginosa	ELISA



Product Name	Code	Species Reactivity	Applications
<u>LASP1 Antibody</u>	CSB-PA04629A0Rb	Human, Mouse	ELISA, WB, IHC, IF, IP
<u>LATS1 Antibody</u>	CSB-PA012769LA01HU	Human	ELISA, IHC, IF, IP
<u>lats2 Antibody</u>	CSB-PA24329A0Rb	Zebrafish	ELISA, WB
<u>LAX1 Antibody</u>	CSB-PA811614LA01HU	Human, Mouse, Rat	ELISA, WB, IF
<u>LAYN Antibody</u>	CSB-PA012772LA01HU	Human	ELISA, IF
<u>LBH Antibody</u>	CSB-PA700659LA01HU	Human	ELISA, IF
<u>LBP Antibody</u>	CSB-PA012775LA01HU	Human	ELISA, IHC
<u>LCA5 Antibody</u>	CSB-PA769789LA01HU	Human, Mouse	ELISA, WB
<u>LCAT Antibody</u>	CSB-PA11369A0Rb	Human	ELISA, WB
<u>LCE1B Antibody</u>	CSB-PA012785LA01HU	Human	ELISA, IHC, IF
<u>LCE2A Antibody</u>	CSB-PA719417LA01HU	Human	ELISA, IHC
<u>LCE3C Antibody</u>	CSB-PA725885LA01HU	Human	ELISA, IF
<u>LCHN Antibody</u>	CSB-PA012231LA01HU	Human	ELISA, WB, IHC
<u>LCK Antibody</u>	CSB-PA012806LA01HU	Human	ELISA, IHC, IF, IP
<u>LCN2 Antibody</u>	CSB-PA06834A0Rb	Human	ELISA, WB
<u>LCN2 Antibody</u>	CSB-PA06836A0Rb	Human	ELISA, WB
<u>LCN2 Antibody</u>	CSB-PA06837A0Rb	Human	ELISA
<u>Lcn3 Antibody</u>	CSB-PA730794HA01MO	Mouse	ELISA
<u>LCOR Antibody</u>	CSB-PA856960LA01HU	Human	ELISA, WB
<u>LCP1 Antibody</u>	CSB-PA012824HA01HU	Human	ELISA, WB, IHC, IP
<u>LCP1 Antibody</u>	CSB-PA012824YA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>LCT Antibody</u>	CSB-PA012827LA01HU	Human	ELISA, IHC
<u>LCY1 Antibody</u>	CSB-PA657129LA01DOA	Arabidopsis thaliana	ELISA
<u>LDAH Antibody</u>	CSB-PA875685LA01HU	Human	ELISA, WB, IHC
<u>LDB1 Antibody</u>	CSB-PA771433LA01HU	Human	ELISA, WB, IHC, IF
<u>ldcA Antibody</u>	CSB-PA303469HA01ENV	Escherichia coli	ELISA
<u>ldcC Antibody</u>	CSB-PA16349A0Rb	Escherichia coli	ELISA, WB
<u>ldh1 Antibody</u>	CSB-PA348716EA01LMS	Lactobacillus plantarum	ELISA
<u>ldh1 Antibody</u>	CSB-PA354501HA01SKX	Staphylococcus aureus	ELISA, WB
<u>LDHA Antibody</u>	CSB-PA00045A0Rb	Human	ELISA, WB, IHC, IF
<u>LDHC Antibody</u>	CSB-PA012844LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>LDHD Antibody</u>	CSB-PA801250LA01HU	Human	ELISA, WB, IHC, IF
<u>LDLR Antibody</u>	CSB-PA09949A0Rb	Human	ELISA, IF
<u>LDLRAD1 Antibody</u>	CSB-PA729206LA01HU	Human	ELISA, IF
<u>lecA Antibody</u>	CSB-PA313124HA01EZX	Pseudomonas aeruginosa	ELISA



Product Name	Code	Species Reactivity	Applications
<u>LECT2 Antibody</u>	CSB-PA012855LA01HU	Human	ELISA, IHC
<u>LEF1 Antibody</u>	CSB-PA012856LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>LEFTY1 Antibody</u>	CSB-PA012857LA01HU	Human, Mouse, Rat	ELISA, WB
<u>LEMD1 Antibody</u>	CSB-PA715035LA01HU	Human	ELISA, IHC, IF
<u>LEMD2 Antibody</u>	CSB-PA854097LA01HU	Human	ELISA, WB, IHC
<u>LEO1 Antibody</u>	CSB-PA845150NA01HU	Human	ELISA, WB, IHC
<u>LEP Antibody</u>	CSB-PA06799A0Rb	Human	ELISA, IHC, IF
<u>LEP Antibody</u>	CSB-PA012870LA01BO	Bovine	ELISA, WB
<u>lepA Antibody</u>	CSB-PA350834LA01ENV	Escherichia coli	ELISA
<u>lepB Antibody</u>	CSB-PA360447LA01ENV	Escherichia coli	ELISA
<u>lepB Antibody</u>	CSB-PA604637LA01MVZ	Mycobacterium tuberculosis	ELISA
<u>LEPR Antibody</u>	CSB-PA06809A0Rb	Human	ELISA, IF
<u>LETM1 Antibody</u>	CSB-PA012877LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>lexA Antibody</u>	CSB-PA363985HA01ENV	Escherichia coli	ELISA, WB
<u>LEXM Antibody</u>	CSB-PA670268LA01HU	Human, Rat	ELISA, WB, IHC
<u>LGALS1 Antibody</u>	CSB-PA012882HA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>LGALS12 Antibody</u>	CSB-PA846606LA01HU	Human, Rat	ELISA, WB, IHC
<u>LGALS14 Antibody</u>	CSB-PA848400LA01HU	Human	ELISA, WB
<u>LGALS3 Antibody</u>	CSB-PA012887YA01HU	Human	ELISA, WB, IHC
<u>Lgals3bp Antibody</u>	CSB-PA012888LA01RA	Rat	ELISA
<u>Lgals4 Antibody</u>	CSB-PA012889HA01RA	Rat	ELISA, IF
<u>LGALS4 Antibody</u>	CSB-PA012889LA01HU	Human	ELISA, IHC, IF
<u>LGALS7 Antibody</u>	CSB-PA09727A0Rb	Human	ELISA, IHC
<u>Lgals7 Antibody</u>	CSB-PA09877A0Rb	Mouse	ELISA
<u>LGALS9 Antibody</u>	CSB-PA02649A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>LGALS9C Antibody</u>	CSB-PA728553LA01HU	Human, Mouse	ELISA, WB
<u>LGB Antibody</u>	CSB-PA00280E0Rb	Bovine	ELISA, WB
<u>LGI2 Antibody</u>	CSB-PA822674LA01HU	Human	ELISA, IHC
<u>LGI4 Antibody</u>	CSB-PA818678LA01HU	Human	ELISA, IHC
<u>LGMN Antibody</u>	CSB-PA012903LA01HU	Human	ELISA, IHC
<u>LGR4 Antibody</u>	CSB-PA883618LA01HU	Human	ELISA, IHC, IF
<u>LGR5 Antibody</u>	CSB-PA012906LA01HU	Human	ELISA, IHC, IF
<u>LGR6 Antibody</u>	CSB-PA884515LA01HU	Human	ELISA, WB, IHC
<u>LHB Antibody</u>	CSB-PA14659A0Rb	Human	ELISA, IHC
<u>LHCGR Antibody</u>	CSB-PA012911LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>Lhegr Antibody</u>	CSB-PA012911LA01RA	Rat	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>LHFPL1 Antibody</u>	CSB-PA774822LA01HU	Human	ELISA, IHC, IF
<u>LHFPL5 Antibody</u>	CSB-PA844703LA01HU	Human	ELISA, IHC, IF
<u>LHFPL6 Antibody</u>	CSB-PA012912LA01HU	Human	ELISA
<u>LHPP Antibody</u>	CSB-PA867114HA01HU	Human	ELISA, IHC
<u>LHX2 Antibody</u>	CSB-PA012920LA01HU	Human, Rat	ELISA, WB, IF
<u>LHX3 Antibody</u>	CSB-PA871386LA01HU	Human	ELISA, WB
<u>LIAS Antibody</u>	CSB-PA012927LA01HU	Human	ELISA, WB, IHC
<u>LIFR Antibody</u>	CSB-PA012929LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>LIG1 Antibody</u>	CSB-PA012930YA01HU	Human	ELISA, IHC, IF
<u>LIG3 Antibody</u>	CSB-PA012931LA01HU	Human	ELISA, IHC, IF
<u>ligA Antibody</u>	CSB-PA322745HA01ENV	Escherichia coli	ELISA, WB
<u>ligA Antibody</u>	CSB-PA718299LA01DAAA	Desulfotalea psychrophila	ELISA, WB
<u>LILRA2 Antibody</u>	CSB-PA854017LA01HU	Human	ELISA, WB, IHC
<u>LILRA3 Antibody</u>	CSB-PA836677LA01HU	Human	ELISA, WB, IHC
<u>LILRA4 Antibody</u>	CSB-PA012936LA01HU	Human	ELISA, IHC, IF
<u>LILRB1 Antibody</u>	CSB-PA847747LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>LILRB2 Antibody</u>	CSB-PA013000LA01HU	Human	ELISA, WB, IF
<u>Lilrb3 Antibody</u>	CSB-PA012941YA01MO	Mouse	ELISA
<u>LILRB3 Antibody</u>	CSB-PA012941LA01HU	Human	ELISA, WB, IF
<u>LILRB5 Antibody</u>	CSB-PA012943LA01HU	Human	ELISA, IHC
<u>LIMD2 Antibody</u>	CSB-PA012950LA01HU	Human	ELISA, IHC, IF
<u>LIME1 Antibody</u>	CSB-PA884463LA01HU	Human	ELISA, IHC, IF
<u>LIMK1 Antibody</u>	CSB-PA012952LA01HU	Human	ELISA, IHC
<u>LIMS1 Antibody</u>	CSB-PA012955EA01HU	Human, Rat	ELISA, WB, IHC
<u>LIMS3 Antibody</u>	CSB-PA012957LA01HU	Human	ELISA, IHC, IF
<u>LIN28A Antibody</u>	CSB-PA13619A0Rb	Human	ELISA, IHC, IF
<u>LIN28B Antibody</u>	CSB-PA741049LA01HU	Human	ELISA, IHC, IF
<u>LIN37 Antibody</u>	CSB-PA836222LA01HU	Human	ELISA, IHC
<u>LIN7A Antibody</u>	CSB-PA012963LA01HU	Human, Mouse	ELISA, WB
<u>LIN9 Antibody</u>	CSB-PA716370LA01HU	Human	ELISA, IHC
<u>LINC00299 Antibody</u>	CSB-PA765077LA01HU	Human	ELISA, IF
<u>LINC00846 Antibody</u>	CSB-PA882112LA01HU	Human	ELISA
<u>LINC01553 Antibody</u>	CSB-PA002950LA01HU	Human, Mouse, Rat	ELISA, WB
<u>LIPA Antibody</u>	CSB-PA012972EA01HU	Human	ELISA, WB, IHC, IF
<u>LIPC Antibody</u>	CSB-PA012974LA01HU	Human	ELISA, IHC
<u>LIPE Antibody</u>	CSB-PA012975LA01HU	Human	ELISA, IHC, IF
<u>LIPG Antibody</u>	CSB-PA897572LA01HU	Human	ELISA, IHC, IF, IP
<u>LIPH Antibody</u>	CSB-PA848833LA01HU	Human	ELISA, IHC, IF
<u>LIX1 Antibody</u>	CSB-PA822702LA01HU	Human, Mouse	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>IktA Antibody</u>	CSB-PA314743EA01ESE	Mannheimia haemolytica	ELISA
<u>L-lactate dehydrogenase Antibody</u>	CSB-PA632935LA01PLS	Plasmodium falciparum	ELISA
<u>LLCFC1 Antibody</u>	CSB-PA839358LA01HU	Human	ELISA, WB, IHC
<u>LLGL1 Antibody</u>	CSB-PA622993LA01HU	Human	ELISA, IHC, IF
<u>LLGL2 Antibody</u>	CSB-PA764711LA01HU	Human	ELISA, IF
<u>LLPH Antibody</u>	CSB-PA861123LA01HU	Human	ELISA
<u>LMAN1 Antibody</u>	CSB-PA012991LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>LMAN1L Antibody</u>	CSB-PA862055LA01HU	Human	ELISA, IHC
<u>LMAN2L Antibody</u>	CSB-PA863945LA01HU	Human	ELISA
<u>LMBRD1 Antibody</u>	CSB-PA885730LA01HU	Human, Mouse	ELISA, WB, IHC
<u>LMCD1 Antibody</u>	CSB-PA873714LA01HU	Human	ELISA, IHC, IF
<u>LMF2 Antibody</u>	CSB-PA013001LA01HU	Human, Mouse	ELISA, WB, IF
<u>LMNA Antibody</u>	CSB-PA013003HA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>LMNB1 Antibody</u>	CSB-PA10169A0Rb	Human	ELISA, IHC
<u>LMNB2 Antibody</u>	CSB-PA013005LA01HU	Human	ELISA, IHC, IF
<u>LMNTD1 Antibody</u>	CSB-PA854083LA01HU	Human, Mouse	ELISA, WB
<u>LMO1 Antibody</u>	CSB-PA02055A0Rb	Human	ELISA, WB
<u>LMO2 Antibody</u>	CSB-PA013009HA01HU	Human, Mouse	ELISA, WB, IHC
<u>LMO3 Antibody</u>	CSB-PA844705LA01HU	Human	ELISA, WB, IHC
<u>LMO7 Antibody</u>	CSB-PA823897LA01HU	Human	ELISA, IHC, IF
<u>LMP2 Antibody</u>	CSB-PA321086HA01EFA	Epstein-Barr virus	ELISA, WB
<u>LMX1A Antibody</u>	CSB-PA848414LA01HU	Human	ELISA, WB, IHC
<u>LNX1 Antibody</u>	CSB-PA848393LA01HU	Human	ELISA, IHC
<u>LNX2 Antibody</u>	CSB-PA854035LA01HU	Human	ELISA, IHC
<u>lolA Antibody</u>	CSB-PA419575LA01EJF	Escherichia coli O9:H4	ELISA
<u>LOLPIB Antibody</u>	CSB-PA670318HA01LQK	Lolium perenne	ELISA
<u>LONP1 Antibody</u>	CSB-PA013032LA01HU	Human	ELISA, IHC, IF
<u>LONRF1 Antibody</u>	CSB-PA615067LA01HU	Human	ELISA, IHC
<u>Lox Antibody</u>	CSB-PA013038LA01RA	Rat, Human	ELISA, WB
<u>LOX Antibody</u>	CSB-PA16909A0Rb	Human	ELISA, WB, IHC
<u>LOX Antibody</u>	CSB-PA013038LA01HU	Human, Rat	ELISA, WB, IHC
<u>LOXHD1 Antibody</u>	CSB-PA814211LA01HU	Human	ELISA, IHC, IF
<u>LOXL1 Antibody</u>	CSB-PA013040LA01HU	Human	ELISA, IHC
<u>Loxl1 Antibody</u>	CSB-PA013040LA01MO	Mouse	ELISA
<u>Loxl2 Antibody</u>	CSB-PA013041LA01MO	Mouse, Human	ELISA, WB
<u>LOXL2 Antibody</u>	CSB-PA897283LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>LOXL3 Antibody</u>	CSB-PA013042LA01HU	Human	ELISA, IHC
<u>Loxl4 Antibody</u>	CSB-PA846069LA01MO	Mouse, Rat	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>LPA Antibody</u>	CSB-PA362446YA01HU	Human	ELISA
<u>LPAL2 Antibody</u>	CSB-PA618087LA01HU	Human	ELISA, IHC, IF
<u>LPAR1 Antibody</u>	CSB-PA013047LA01HU	Human	ELISA, IHC, IF
<u>LPAR2 Antibody</u>	CSB-PA013048OA01HU	Human	ELISA, WB, IHC
<u>LPAR3 Antibody</u>	CSB-PA013049OA01HU	Human	ELISA, IHC
<u>LPAR4 Antibody</u>	CSB-PA013050LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>LPAR5 Antibody</u>	CSB-PA013051OA01HU	Human	ELISA, IHC, IF
<u>LPAR6 Antibody</u>	CSB-PA017336OA01HU	Human	ELISA, WB, IHC, IF
<u>LPCAT2 Antibody</u>	CSB-PA013053LA01HU	Human	ELISA, WB, IHC
<u>LPCAT3 Antibody</u>	CSB-PA744028LA01HU	Human	ELISA, IHC
<u>LPGAT1 Antibody</u>	CSB-PA821685LA01HU	Human	ELISA, WB, IHC
<u>LPL Antibody</u>	CSB-PA013065LA01HU	Human	ELISA, IHC, IF
<u>lplA Antibody</u>	CSB-PA510777HA01ENU	Escherichia coli	ELISA, WB
<u>LPO Antibody</u>	CSB-PA013066LA01HU	Human	ELISA, WB, IHC
<u>LPP Antibody</u>	CSB-PA856611LA01HU	Human	ELISA, WB, IHC
<u>lptA Antibody</u>	CSB-PA360010HA01ENV	Escherichia coli	ELISA
<u>lpxC Antibody</u>	CSB-PA358910HA01ENV	Escherichia coli	ELISA
<u>lpxC Antibody</u>	CSB-PA342468LA01EZXX	Pseudomonas aeruginosa	ELISA
<u>LPXN Antibody</u>	CSB-PA013068LA01HU	Human	ELISA, WB, IHC, IF
<u>LRAT Antibody</u>	CSB-PA013069LA01HU	Human	ELISA, WB
<u>LRBA Antibody</u>	CSB-PA013070LA01HU	Human	ELISA, IHC, IF
<u>LRCH1 Antibody</u>	CSB-PA897471LA01HU	Human	ELISA, WB, IHC, IF
<u>LRG1 Antibody</u>	CSB-PA11059A0Rb	Human, Mouse	ELISA, WB, IHC
<u>LRIG2 Antibody</u>	CSB-PA013085LA01HU	Human	ELISA, WB, IF
<u>LRIG3 Antibody</u>	CSB-PA751017LA01HU	Human	ELISA, IF
<u>LRIT1 Antibody</u>	CSB-PA879074LA01HU	Human	ELISA, WB
<u>LRP1 Antibody</u>	CSB-PA013091NA01HU	Human	ELISA, IHC, IF
<u>LRP2 Antibody</u>	CSB-PA009078LA01HU	Human	ELISA, IHC, IF
<u>LRP2 Antibody</u>	CSB-PA013096LA01HU	Human	ELISA, IHC
<u>LRP8 Antibody</u>	CSB-PA623900LA01HU	Human	ELISA, WB, IF
<u>LRPAP1 Antibody</u>	CSB-PA013104HA01HU	Human	ELISA, WB
<u>LRPPRC Antibody</u>	CSB-PA013105LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>LRR1 Antibody</u>	CSB-PA836247LA01HU	Human	ELISA, IHC, IF
<u>LRRC15 Antibody</u>	CSB-PA819484LA01HU	Human	ELISA, IHC, IF
<u>LRRC18 Antibody</u>	CSB-PA839794LA01HU	Human	ELISA, IHC
<u>LRRC19 Antibody</u>	CSB-PA872509LA01HU	Human	ELISA, IHC, IF
<u>LRRC2 Antibody</u>	CSB-PA866321LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>LRRC20 Antibody</u>	CSB-PA823454LA01HU	Human, Mouse	ELISA, WB, IHC
<u>LRRC29 Antibody</u>	CSB-PA819894LA01HU	Human	ELISA
<u>LRRC31 Antibody</u>	CSB-PA744274LA01HU	Human, Mouse	ELISA, WB, IHC
<u>LRRC32 Antibody</u>	CSB-PA618886LA01HU	Human	ELISA, IHC
<u>LRRC39 Antibody</u>	CSB-PA822195LA01HU	Human	ELISA, IHC
<u>LRRC3B Antibody</u>	CSB-PA857000LA01HU	Human	ELISA, IHC, IF
<u>LRRC3C Antibody</u>	CSB-PA408988LA01HU	Human	ELISA, IHC
<u>LRRC40 Antibody</u>	CSB-PA881002LA01HU	Human	ELISA, WB
<u>LRRC42 Antibody</u>	CSB-PA896899LA01HU	Human	ELISA, WB, IHC
<u>LRRC45 Antibody</u>	CSB-PA850272LA01HU	Human	ELISA, IHC, IF
<u>LRRC4C Antibody</u>	CSB-PA888027LA01HU	Human	ELISA, IHC, IF
<u>LRRC57 Antibody</u>	CSB-PA843302LA01HU	Human, Mouse	ELISA, WB
<u>LRRC61 Antibody</u>	CSB-PA871582LA01HU	Human	ELISA, IHC
<u>LRRC8E Antibody</u>	CSB-PA743982LA01HU	Human	ELISA, IHC, IF
<u>LRRFIP1 Antibody</u>	CSB-PA013175HA01HU	Human	ELISA, IHC, IF
<u>LRRFIP2 Antibody</u>	CSB-PA896921LA01HU	Human	ELISA, WB, IHC
<u>LRRIQ3 Antibody</u>	CSB-PA013178LA01HU	Human	ELISA, IHC, IF
<u>LRRK1 Antibody</u>	CSB-PA655748LA01HU	Human	ELISA, WB
<u>LRRK2 Antibody</u>	CSB-PA722493LA01HU	Human	ELISA, IHC, IF
<u>LRRN1 Antibody</u>	CSB-PA747638LA01HU	Human	ELISA, IHC
<u>LRRN2 Antibody</u>	CSB-PA013183LA01HU	Human	ELISA, IHC, IF
<u>LRRN3 Antibody</u>	CSB-PA880968LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>LRRTM2 Antibody</u>	CSB-PA013188LA01HU	Human	ELISA
<u>LRRTM4 Antibody</u>	CSB-PA774813LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>LRTM1 Antibody</u>	CSB-PA884513LA01HU	Human	ELISA, WB, IHC, IF
<u>LRTM2 Antibody</u>	CSB-PA847660LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>LRTOMT Antibody</u>	CSB-PA02667A0Rb	Human	ELISA, IF
<u>LSAMP Antibody</u>	CSB-PA013197LA01HU	Human, Mouse, Rat	ELISA, WB, IF
<u>LSM1 Antibody</u>	CSB-PA013199LA01HU	Human	ELISA, IHC
<u>LSM10 Antibody</u>	CSB-PA836169LA01HU	Human	ELISA, IHC
<u>LSM11 Antibody</u>	CSB-PA013201LA01HU	Human	ELISA, WB, IF
<u>LSM3 Antibody</u>	CSB-PA013206LA01HU	Human	ELISA
<u>LSM4 Antibody</u>	CSB-PA013207HA01HU	Human	ELISA, WB, IHC, IF
<u>LSM6 Antibody</u>	CSB-PA013209LA01HU	Human	ELISA
<u>LSM7 Antibody</u>	CSB-PA890731LA01HU	Human	ELISA, WB, IHC
<u>LSM8 Antibody</u>	CSB-PA013212LA01HU	Human	ELISA, IHC
<u>LSMEM1 Antibody</u>	CSB-PA818730LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>LSMEM2 Antibody</u>	CSB-PA822677LA01HU	Human	ELISA, IHC, IF
<u>lspL Antibody</u>	CSB-PA526307HA01RKU	Rhizobium meliloti	ELISA
<u>LSR Antibody</u>	CSB-PA801252LA01HU	Human	ELISA, IHC, IF
<u>lta Antibody</u>	CSB-PA24269A0Rb	Zebrafish	ELISA, WB
<u>LTA4H Antibody</u>	CSB-PA013219LA01HU	Human	ELISA, WB, IHC
<u>LTB Antibody</u>	CSB-PA771253LA01EQV	Chimpanzee	ELISA
<u>Ltbp1 Antibody</u>	CSB-PA807018LA01MO	Mouse	ELISA
<u>LTBP4 Antibody</u>	CSB-PA02049A0Rb	Human	ELISA, IHC, IF
<u>LTBR Antibody</u>	CSB-PA08129A0Rb	Human	ELISA, IHC, IF
<u>LTF Antibody</u>	CSB-PA00870E0Rb	Human	ELISA, IHC, IF
<u>LTF Antibody</u>	CSB-PA013229JA01HU	Human	ELISA
<u>LTF Antibody</u>	CSB-PA013229LA01HU	Human	ELISA, IHC
<u>LTK Antibody</u>	CSB-PA013230LA01HU	Human	ELISA, IF
<u>LTN1 Antibody</u>	CSB-PA026646LA01HU	Human	ELISA, IHC
<u>LTV1 Antibody</u>	CSB-PA853422LA01HU	Human	ELISA, IHC
<u>Lum Antibody</u>	CSB-PA013234LA01RA	Rat	ELISA, WB
<u>luxS Antibody</u>	CSB-PA541717HA01ENV	Escherichia coli	ELISA, WB
<u>LUZP2 Antibody</u>	CSB-PA013236LA01HU	Human	ELISA, WB, IHC, IF
<u>LXN Antibody</u>	CSB-PA013240LA01HU	Human	ELISA, IHC, IF
<u>LY6D Antibody</u>	CSB-PA613492LA01HU	Human	ELISA
<u>LY6E Antibody</u>	CSB-PA619076LA01HU	Human	ELISA, IHC, IF
<u>Ly6e Antibody</u>	CSB-PA717533LA01MO	Mouse, Human	ELISA, WB
<u>LY6G6C Antibody</u>	CSB-PA013245LA01HU	Human	ELISA
<u>LY6H Antibody</u>	CSB-PA013249LA01HU	Human	ELISA, IHC, IF
<u>LY6K Antibody</u>	CSB-PA621061LA01HU	Human	ELISA, IHC
<u>Ly75 Antibody</u>	CSB-PA013251HA01MO	Mouse	ELISA
<u>LY96 Antibody</u>	CSB-PA013254LA01HU	Human	ELISA, IHC
<u>LYAR Antibody</u>	CSB-PA882130LA01HU	Human	ELISA, WB, IHC, IF
<u>LYG1 Antibody</u>	CSB-PA843138LA01HU	Human	ELISA, IHC
<u>LYG2 Antibody</u>	CSB-PA773017LA01HU	Human	ELISA, IHC
<u>LYN Antibody</u>	CSB-PA013259LA01HU	Human	ELISA, IF
<u>LYPD2 Antibody</u>	CSB-PA761498LA01HU	Human	ELISA, IHC
<u>LYPD3 Antibody</u>	CSB-PA013263LA01HU	Human, Mouse	ELISA, WB
<u>LYPD4 Antibody</u>	CSB-PA751011LA01HU	Human	ELISA
<u>Lypla1 Antibody</u>	CSB-PA013268DA01MO	Mouse	ELISA
<u>Lypla1 Antibody</u>	CSB-PA013268LA01MO	Mouse	ELISA, IF
<u>LYPLA2 Antibody</u>	CSB-PA05074A0Rb	Human	ELISA, IF
<u>LYPLAL1 Antibody</u>	CSB-PA736204LA01HU	Human	ELISA, IHC, IF
<u>LYRM1 Antibody</u>	CSB-PA013272LA01HU	Human	ELISA, IHC
<u>LYRM2 Antibody</u>	CSB-PA885727LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>LYRM4 Antibody</u>	CSB-PA875729LA01HU	Human	ELISA, IHC
<u>LYSMD1 Antibody</u>	CSB-PA836299LA01HU	Human	ELISA, IHC
<u>LYSMD3 Antibody</u>	CSB-PA768760LA01HU	Human, Mouse	ELISA, WB
<u>Lysozyme C-1 Antibody</u>	CSB-PA360572LA01BZD	Anas platyrhynchos	ELISA, WB
<u>LYVE1 Antibody</u>	CSB-PA897574EA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>LYZ Antibody</u>	CSB-PA02769A0Rb	Human	ELISA, IHC
<u>LYZL1 Antibody</u>	CSB-PA747630LA01HU	Human	ELISA, IHC
<u>LYZL2 Antibody</u>	CSB-PA773594LA01HU	Human	ELISA, IHC
<u>LZTFL1 Antibody</u>	CSB-PA885687LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>LZTS1 Antibody</u>	CSB-PA22699A0Rb	Human	ELISA, IHC, IF
<u>M Antibody</u>	CSB-PA322192LA01RAI	Rabies virus	ELISA, WB
<u>M Antibody</u>	CSB-PA356390LA01VBG	Vesicular stomatitis Indiana virus	ELISA, WB
<u>M6PR Antibody</u>	CSB-PA02625A0Rb	Human	ELISA, IHC, IF
<u>MAB21L1 Antibody</u>	CSB-PA623800LA01HU	Human	ELISA, WB, IHC
<u>MAB21L2 Antibody</u>	CSB-PA896901LA01HU	Human, Mouse	ELISA, WB, IF
<u>MACC1 Antibody</u>	CSB-PA761627LA01HU	Human	ELISA, WB, IHC, IF
<u>MACF1 Antibody</u>	CSB-PA890946LA01HU	Human	ELISA, IHC
<u>MACROD1 Antibody</u>	CSB-PA860244LA01HU	Human	ELISA, IHC
<u>MAD2L1 Antibody</u>	CSB-PA619756LA01HU	Human	ELISA, IHC
<u>MAD2L2 Antibody</u>	CSB-PA887064LA01HU	Human	ELISA
<u>MADCAM1 Antibody</u>	CSB-PA618776HA01HU	Human	ELISA, WB, IF
<u>MAEA Antibody</u>	CSB-PA013310LA01HU	Human	ELISA, WB, IHC, IF
<u>MAEL Antibody</u>	CSB-PA842695LA01HU	Human	ELISA, IHC, IF
<u>MAF Antibody</u>	CSB-PA013312LA01HU	Human	ELISA, IHC
<u>MAF1 Antibody</u>	CSB-PA872421HA01HU	Human, Mouse	ELISA, WB
<u>MAFB Antibody</u>	CSB-PA897570LA01HU	Human, Mouse	ELISA, WB
<u>MAG Antibody</u>	CSB-PA013322LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>MAGEA10 Antibody</u>	CSB-PA013324LA01HU	Human	ELISA, WB, IHC
<u>MAGEA2 Antibody</u>	CSB-PA013327LA01HU	Human	ELISA, WB
<u>MAGEA3 Antibody</u>	CSB-PA013329LA01HU	Human	ELISA, WB, IHC
<u>MAGEA4 Antibody</u>	CSB-PA013330LA01HU	Human	ELISA, IHC
<u>MAGEA5 Antibody</u>	CSB-PA013331LA01HU	Human	ELISA, IHC
<u>MAGEA8 Antibody</u>	CSB-PA013334LA01HU	Human	ELISA, IHC, IF
<u>MAGEA9 Antibody</u>	CSB-PA013335LA01HU	Human	ELISA, IHC
<u>MAGEB10 Antibody</u>	CSB-PA853459LA01HU	Human, Mouse	ELISA, WB
<u>MAGEB18 Antibody</u>	CSB-PA853461LA01HU	Human, Mouse	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>MAGEC3 Antibody</u>	CSB-PA819466LA01HU	Human, Rat	ELISA, WB
<u>MAGED2 Antibody</u>	CSB-PA891976LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>MAGI1 Antibody</u>	CSB-PA857014LA01HU	Human	ELISA, IF
<u>MAGI2 Antibody</u>	CSB-PA771436LA01HU	Human	ELISA, IHC, IF
<u>MAGOHB Antibody</u>	CSB-PA822175LA01HU	Human	ELISA, IHC
<u>MAIP1 Antibody</u>	CSB-PA855488LA01HU	Human, Mouse	ELISA, WB, IHC
<u>MAJIN Antibody</u>	CSB-PA659444LA01HU	Human	ELISA, IHC
<u>Major allergen Api g 1, isoallergen 2 Antibody</u>	CSB-PA310896LA01DNL	Apium graveolens	ELISA
<u>Major allergen Bla g 1.0101 Antibody</u>	CSB-PA19837A0Rb	Blattella germanica	ELISA, WB
<u>Major pepsin inhibitor 3 Antibody</u>	CSB-PA322527LA01DOS	Ascaris suum	ELISA
<u>Major pollen allergen Aln g 1 Antibody</u>	CSB-PA19897A0Rb	Alnus glutinosa	ELISA, WB
<u>Major pollen allergen Art v 1 Antibody</u>	CSB-PA774628HA01AOH	Artemisia vulgaris	ELISA, WB
<u>Major pollen allergen Car b 1 isoforms 1A and 1B Antibody</u>	CSB-PA330956LA01CEC	Carpinus betulus	ELISA, WB
<u>Major pollen allergen Dac g 4 Antibody</u>	CSB-PA307868LA01DAC	Dactylis glomeRata	ELISA
<u>MAK Antibody</u>	CSB-PA013369LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>MAL Antibody</u>	CSB-PA013372LA01HU	Human	ELISA, IHC, IF
<u>MALD1 Antibody</u>	CSB-PA331771HA01MPS	Apple	ELISA
<u>MALD3 Antibody</u>	CSB-PA862991HA01MPS	Apple	ELISA, WB
<u>MAMDC2 Antibody</u>	CSB-PA801798LA01HU	Human	ELISA
<u>MAML1 Antibody</u>	CSB-PA838784LA01HU	Human	ELISA, IHC
<u>MAML2 Antibody</u>	CSB-PA809028LA01HU	Human	ELISA, WB, IHC
<u>MAN2A1 Antibody</u>	CSB-PA613697LA01HU	Human	ELISA, WB, IHC, IF
<u>MAN2B2 Antibody</u>	CSB-PA896489LA01HU	Human	ELISA, WB, IHC
<u>MAN2C1 Antibody</u>	CSB-PA878885LA01HU	Human	ELISA, IHC, IF
<u>manA Antibody</u>	CSB-PA014754HA01ENV	Escherichia coli	ELISA, WB
<u>MANBAL Antibody</u>	CSB-PA889080LA01HU	Human	ELISA, IHC, IF
<u>MAOA Antibody</u>	CSB-PA013398LA01HU	Human	ELISA, WB, IHC
<u>MAOB Antibody</u>	CSB-PA013399LA01HU	Human, Mouse	ELISA, WB, IHC
<u>MAP10 Antibody</u>	CSB-PA879069LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>MAP1A Antibody</u>	CSB-PA013400LA01HU	Human	ELISA, IHC, IF
<u>MAP1LC3A Antibody</u>	CSB-PA880971DA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>MAP1LC3B Antibody</u>	CSB-PA887936LA01HU	Human	ELISA, IHC
<u>MAP1LC3B2 Antibody</u>	CSB-PA013404LA01HU	Human	ELISA, IF
<u>MAP1LC3C Antibody</u>	CSB-PA887153LA01HU	Human	ELISA, IHC, IF
<u>MAP1S Antibody</u>	CSB-PA13549A0Rb	Human	ELISA, IHC
<u>MAP2K1 Antibody</u>	CSB-PA013409LA01HU	Human, Mouse	ELISA, WB, IHC
<u>MAP2K2 Antibody</u>	CSB-PA04909A0Rb	Human	ELISA, WB, IHC, IF
<u>MAP2K3 Antibody</u>	CSB-PA013412LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>Map2k4 Antibody</u>	CSB-PA013413LA01MO	Mouse, Rat, Human	ELISA, WB
<u>MAP2K4 Antibody</u>	CSB-PA013413LA01HU	Human	ELISA, WB, IHC, IF
<u>MAP2K5 Antibody</u>	CSB-PA615667HA01HU	Human, Mouse, Rat	ELISA, WB, IF
<u>MAP2K6 Antibody</u>	CSB-PA013415LA01HU	Human	ELISA, IF
<u>MAP2K7 Antibody</u>	CSB-PA013416LA01HU	Human, Mouse	ELISA, WB, IF
<u>Map3k1 Antibody</u>	CSB-PA15289A0Rb	Mouse	ELISA
<u>MAP3K11 Antibody</u>	CSB-PA01469A0Rb	Human	ELISA, IHC, IF
<u>MAP3K12 Antibody</u>	CSB-PA013420LA01HU	Human	ELISA
<u>MAP3K13 Antibody</u>	CSB-PA013421LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>MAP3K19 Antibody</u>	CSB-PA707114LA01HU	Human	ELISA
<u>MAP3K2 Antibody</u>	CSB-PA013424LA01HU	Human	ELISA, IHC
<u>MAP3K21 Antibody</u>	CSB-PA719421LA01HU	Human	ELISA, WB
<u>MAP3K4 Antibody</u>	CSB-PA013426LA01HU	Human	ELISA, IHC
<u>MAP3K5 Antibody</u>	CSB-PA013427LA01HU	Human	ELISA, IHC, IF
<u>MAP3K6 Antibody</u>	CSB-PA013428LA01HU	Human	ELISA, IHC, IF
<u>MAP3K7 Antibody</u>	CSB-PA013429LA01HU	Human	ELISA, IHC
<u>MAP3K7CL Antibody</u>	CSB-PA349825HA01HU	Human	ELISA, IHC, IF
<u>MAP3K8 Antibody</u>	CSB-PA013433LA01HU	Human	ELISA, IHC, IF
<u>MAP3K9 Antibody</u>	CSB-PA013434LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>MAP4 Antibody</u>	CSB-PA013435LA01HU	Human	ELISA, IHC, IF
<u>MAP4K2 Antibody</u>	CSB-PA013437LA01HU	Human	ELISA, IHC, IF
<u>MAP4K4 Antibody</u>	CSB-PA013439LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>MAPK1 Antibody</u>	CSB-PA013448LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>MAPK10 Antibody</u>	CSB-PA013449LA01HU	Human	ELISA, IHC, IF
<u>MAPK11 Antibody</u>	CSB-PA013450LA01HU	Human	ELISA, IHC, IF
<u>MAPK12 Antibody</u>	CSB-PA013451LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>MAPK13 Antibody</u>	CSB-PA013452LA01HU	Human	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>MAPK14 Antibody</u>	CSB-PA013453LA01HU	Human	ELISA, WB, IHC
<u>MAPK3 Antibody</u>	CSB-PA013456LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>MAPK6 Antibody</u>	CSB-PA013458LA01HU	Human	ELISA, IHC
<u>MAPK8 Antibody</u>	CSB-PA16759A0Rb	Human	ELISA, WB, IHC, IF
<u>MAPK9 Antibody</u>	CSB-PA013471LA01HU	Human	ELISA, IF, IP
<u>MAPKAPK5 Antibody</u>	CSB-PA818241LA01HU	Human	ELISA, IF
<u>MAPRE1 Antibody</u>	CSB-PA620890LA01HU	Human	ELISA, IHC
<u>MAPRE2 Antibody</u>	CSB-PA618987LA01HU	Human, Mouse	ELISA, WB, IF
<u>Mapt Antibody</u>	CSB-PA013481YA01RA	Rat	ELISA
<u>MAPT Antibody</u>	CSB-PA11169A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>MARK2 Antibody</u>	CSB-PA762568LA01HU	Human	ELISA, WB, IHC
<u>MARK3 Antibody</u>	CSB-PA013498LA01HU	Human	ELISA, IHC, IF
<u>MARS2 Antibody</u>	CSB-PA836221LA01HU	Human	ELISA, WB, IHC
<u>MARVELD3 Antibody</u>	CSB-PA842633LA01HU	Human	ELISA, IF
<u>MAS1 Antibody</u>	CSB-PA013505LA01HU	Human	ELISA, WB, IHC, IF
<u>MAS1L Antibody</u>	CSB-PA013506LA01HU	Human	ELISA, WB, IHC, IF
<u>MASP1 Antibody</u>	CSB-PA013508LA01HU	Human	ELISA, WB
<u>MASP2 Antibody</u>	CSB-PA013510LA01HU	Human	ELISA, WB, IHC
<u>MAST1 Antibody</u>	CSB-PA897529LA01HU	Human	ELISA, IHC, IF
<u>MAST3 Antibody</u>	CSB-PA013513LA01HU	Human, Rat, Mouse	ELISA, WB, IHC, IF
<u>MASTL Antibody</u>	CSB-PA846629LA01HU	Human	ELISA, WB, IHC
<u>MAT2A Antibody</u>	CSB-PA013517LA01HU	Human	ELISA, IHC, IF
<u>MAT2B Antibody</u>	CSB-PA882157LA01HU	Human, Rat	ELISA, WB, IF
<u>MATN2 Antibody</u>	CSB-PA013521LA01HU	Human	ELISA, IHC
<u>MATN3 Antibody</u>	CSB-PA013522EA01HU	Human, Mouse	ELISA, WB
<u>MATN4 Antibody</u>	CSB-PA013523LA01HU	Human	ELISA
<u>MAX Antibody</u>	CSB-PA013527LA01HU	Human	ELISA, IHC, IF
<u>MAZ Antibody</u>	CSB-PA013528LA01HU	Human	ELISA, WB, IHC
<u>MB Antibody</u>	CSB-PA013529YA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>MB Antibody</u>	CSB-PA07917A0Rb	Human, Mouse, Rat	ELISA, WB, IHC
<u>MBD1 Antibody</u>	CSB-PA883444LA01HU	Human	ELISA, IHC, IF
<u>MBD3 Antibody</u>	CSB-PA013532LA01HU	Human, Mouse	ELISA, WB, IHC
<u>MBD5 Antibody</u>	CSB-PA868388LA01HU	Human	ELISA
<u>MBIP Antibody</u>	CSB-PA889101LA01HU	Human	ELISA, IHC
<u>MBNL1 Antibody</u>	CSB-PA889088LA01HU	Human	ELISA, IHC
<u>MBNL3 Antibody</u>	CSB-PA868322LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>MBP Antibody</u>	CSB-PA20287A0Rb	Macaca fascicularis	ELISA
<u>MBP Antibody</u>	CSB-PA013551LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>MBP Monoclonal Antibody</u>	CSB-MA000061M0m	All	ELISA, WB
<u>MC1R Antibody</u>	CSB-PA013558LA01HU	Human	ELISA, WB, IHC, IF
<u>MC3R Antibody</u>	CSB-PA013560LA01HU	Human	ELISA, IHC, IF
<u>MC4R Antibody</u>	CSB-PA013561LA01HU	Human	ELISA, IHC
<u>MC5R Antibody</u>	CSB-PA013562LA01HU	Human	ELISA, IHC, IF
<u>MCAM Antibody</u>	CSB-PA10439A0Rb	Human	ELISA, IHC
<u>MCAT Antibody</u>	CSB-PA815557LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>MCC Antibody</u>	CSB-PA013571LA01HU	Human	ELISA, IHC, IF
<u>MCCC1 Antibody</u>	CSB-PA853497EA01HU	Human	ELISA, IF
<u>MCF2 Antibody</u>	CSB-PA013579LA01HU	Human	ELISA, IHC, IF
<u>MCHR2 Antibody</u>	CSB-PA013585LA01HU	Human	ELISA, IHC, IF
<u>MCIDAS Antibody</u>	CSB-PA518689LA01HU	Human	ELISA, IHC
<u>Mcl1 Antibody</u>	CSB-PA013588HA01MO	Mouse	ELISA
<u>MCL1 Antibody</u>	CSB-PA03829A0Rb	Human	ELISA, WB, IHC, IF
<u>MCM2 Antibody</u>	CSB-PA013590LA01HU	Human	ELISA, WB, IHC, IF
<u>MCM3AP Antibody</u>	CSB-PA013592LA01HU	Human	ELISA, IHC, IF
<u>MCM4 Antibody</u>	CSB-PA013594LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>MCM8 Antibody</u>	CSB-PA890907LA01HU	Human	ELISA, IHC
<u>MCM9 Antibody</u>	CSB-PA882135LA01HU	Human, Mouse	ELISA, WB, IHC
<u>Mept4 Antibody</u>	CSB-PA325157LA01MO	Mouse	ELISA
<u>mcr1 Antibody</u>	CSB-PA745804LA01ENL	Escherichia coli	ELISA
<u>mcr1 Antibody</u>	CSB-PA745804FA01ENL	Escherichia coli	ELISA
<u>MCRS1 Antibody</u>	CSB-PA013611LA01HU	Human	ELISA, WB, IHC
<u>MCTP2 Antibody</u>	CSB-PA738576LA01HU	Human	ELISA, IHC, IF
<u>MCTS1 Antibody</u>	CSB-PA892143HA01HU	Human	ELISA, IHC, IF
<u>MCU Antibody</u>	CSB-PA847694LA01HU	Human	ELISA, WB, IHC
<u>MCUB Antibody</u>	CSB-PA873671LA01HU	Human	ELISA, IHC, IF
<u>MDC1 Antibody</u>	CSB-PA623929LA01HU	Human	ELISA, IHC, IF
<u>MDFIC Antibody</u>	CSB-PA889174OA01HU	Human	ELISA, IHC, IF
<u>MDH1 Antibody</u>	CSB-PA013621EA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>MDK Antibody</u>	CSB-PA10574A0Rb	Human	ELISA, IHC, IF
<u>MDM2 Antibody</u>	CSB-PA11628A0Rb	Human	ELISA, WB, IHC
<u>MDM4 Antibody</u>	CSB-PA013627LA01HU	Human, Mouse	ELISA, WB
<u>MDR1 Antibody</u>	CSB-PA001046HA01EKM	Entamoeba histolytica	ELISA



Product Name	Code	Species Reactivity	Applications
<u>MDR1 Antibody</u>	CSB-PA318386LA01PLT	Plasmodium falciparum	ELISA
<u>ME1 Antibody</u>	CSB-PA013632LA01HU	Human	ELISA, IHC, IF
<u>ME2 Antibody</u>	CSB-PA013633HA01HU	Human	ELISA, IHC
<u>ME3 Antibody</u>	CSB-PA613700LA01HU	Human, Mouse	ELISA, WB, IHC
<u>MEAF6 Antibody</u>	CSB-PA872530LA01HU	Human	ELISA, IHC
<u>MECOM Antibody</u>	CSB-PA621671LA01HU	Human	ELISA, IHC
<u>MECR Antibody</u>	CSB-PA887132LA01HU	Human, Rat	ELISA, WB, IHC
<u>MED1 Antibody</u>	CSB-PA621880HA01HU	Human	ELISA, IHC
<u>MED12 Antibody</u>	CSB-PA853101LA01HU	Human	ELISA, IHC
<u>MED15 Antibody</u>	CSB-PA842753LA01HU	Human	ELISA, WB, IF
<u>MED17 Antibody</u>	CSB-PA865147LA01HU	Human, Mouse	ELISA, WB, IHC
<u>MED18 Antibody</u>	CSB-PA871577LA01HU	Human	ELISA, WB, IHC
<u>MED20 Antibody</u>	CSB-PA862046LA01HU	Human	ELISA, IHC, IF
<u>MED22 Antibody</u>	CSB-PA620882LA01HU	Human	ELISA, IHC, IF
<u>MED23 Antibody</u>	CSB-PA891548LA01HU	Human	ELISA, IHC, IF
<u>MED24 Antibody</u>	CSB-PA013657LA01HU	Human	ELISA, IHC
<u>MED25 Antibody</u>	CSB-PA741144HA01HU	Human	ELISA, IHC
<u>MED26 Antibody</u>	CSB-PA013659LA01HU	Human	ELISA, IHC
<u>MED27 Antibody</u>	CSB-PA764716LA01HU	Human	ELISA, WB
<u>MED28 Antibody</u>	CSB-PA863955LA01HU	Human	ELISA, IHC
<u>MED30 Antibody</u>	CSB-PA822225LA01HU	Human	ELISA, IHC, IF
<u>MED31 Antibody</u>	CSB-PA013664LA01HU	Human	ELISA, IHC
<u>MED37A Antibody</u>	CSB-PA888748LA01DOA	Arabidopsis thaliana	ELISA, WB
<u>MED4 Antibody</u>	CSB-PA865098LA01HU	Human	ELISA, IF
<u>MED7 Antibody</u>	CSB-PA013667LA01HU	Human, Mouse	ELISA, WB
<u>MEF2A Antibody</u>	CSB-PA013670LA01HU	Human	ELISA, WB
<u>MEF2D Antibody</u>	CSB-PA613522LA01HU	Human, Mouse	ELISA, WB, IHC
<u>MEFV Antibody</u>	CSB-PA013675LA01HU	Human	ELISA, IHC
<u>MEGF10 Antibody</u>	CSB-PA822242LA01HU	Human	ELISA, IHC, IF
<u>MEGF8 Antibody</u>	CSB-PA773606LA01HU	Human	ELISA, IHC, IF
<u>MEGF9 Antibody</u>	CSB-PA872434LA01HU	Human	ELISA, WB, IHC, IF
<u>MEI1 Antibody</u>	CSB-PA713139LA01HU	Human	ELISA, IHC
<u>MEI4 Antibody</u>	CSB-PA427046LA01HU	Human, Mouse	ELISA, WB
<u>MEIKIN Antibody</u>	CSB-PA20649A0Rb	Human	ELISA, IHC
<u>MEIS1 Antibody</u>	CSB-PA013686HA01HU	Human, Mouse	ELISA, WB, IHC, IP
<u>MEIS2 Antibody</u>	CSB-PA013687LA01HU	Human, Mouse	ELISA, WB, IF
<u>MEIS3 Antibody</u>	CSB-PA013689LA01HU	Human	ELISA, IHC
<u>meIA Antibody</u>	CSB-PA523525HA01BRG	Bacillus subtilis	ELISA
<u>MELTF Antibody</u>	CSB-PA013754LA01HU	Human	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>MEN1 Antibody</u>	CSB-PA013694HA01HU	Human	ELISA, IHC, IF
<u>MENT Antibody</u>	CSB-PA871579LA01HU	Human, Rat	ELISA, WB, IHC
<u>MEOX2 Antibody</u>	CSB-PA013696LA01HU	Human	ELISA, WB
<u>MET Antibody</u>	CSB-PA013714LA01HU	Human	ELISA, IHC, IF
<u>met Antibody</u>	CSB-PA24469A0Rb	Zebrafish	ELISA, WB
<u>Met Antibody</u>	CSB-PA27899A0Rb	Locusta migratoria	ELISA, WB
<u>METAP1 Antibody</u>	CSB-PA013715LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>METAP2 Antibody</u>	CSB-PA013716LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>metC Antibody</u>	CSB-PA14479A0Rb	Escherichia coli	ELISA
<u>METRNL Antibody</u>	CSB-PA731035LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>Metrnl Antibody</u>	CSB-PA855340LA01MO	Mouse, Rat	ELISA, WB
<u>METTL14 Antibody</u>	CSB-PA884519LA01HU	Human	ELISA, WB, IHC
<u>METTL21C Antibody</u>	CSB-PA733169LA01HU	Human	ELISA, IHC, IF
<u>METTL25 Antibody</u>	CSB-PA839813LA01HU	Human	ELISA, IHC
<u>METTL2A Antibody</u>	CSB-PA856955LA01HU	Human, Mouse	ELISA, WB
<u>METTL7A Antibody</u>	CSB-PA862041LA01HU	Human	ELISA, IHC
<u>MEX3A Antibody</u>	CSB-PA013740LA01HU	Human	ELISA, WB, IHC
<u>MEX3B Antibody</u>	CSB-PA751152LA01HU	Human	ELISA, IHC, IF
<u>MFAP3 Antibody</u>	CSB-PA013747LA01HU	Human	ELISA
<u>MFAP3L Antibody</u>	CSB-PA013748LA01HU	Human	ELISA, IHC, IF
<u>MFAP4 Antibody</u>	CSB-PA013749LA01HU	Human	ELISA, IHC
<u>MFAP5 Antibody</u>	CSB-PA623799LA01HU	Human	ELISA, IHC
<u>MFF Antibody</u>	CSB-PA013751LA01HU	Human	ELISA, WB
<u>MFGE8 Antibody</u>	CSB-PA07974A0Rb	Human, Mouse	ELISA, WB, IHC
<u>MFN2 Antibody</u>	CSB-PA013756LA01HU	Human, Rat	ELISA, WB, IHC
<u>MFNG Antibody</u>	CSB-PA013757LA01HU	Human	ELISA, IHC
<u>MFSD1 Antibody</u>	CSB-PA863974LA01HU	Human	ELISA, IHC, IF
<u>MFSD11 Antibody</u>	CSB-PA013761LA01HU	Human	ELISA, IHC
<u>MFSD14B Antibody</u>	CSB-PA713104LA01HU	Human	ELISA, IHC, IF
<u>MFSD4A Antibody</u>	CSB-PA850797LA01HU	Human	ELISA, WB, IF
<u>MFSD5 Antibody</u>	CSB-PA747224LA01HU	Human	ELISA, IHC, IF
<u>MFSD8 Antibody</u>	CSB-PA844087LA01HU	Human	ELISA, IHC, IF
<u>MGAT3 Antibody</u>	CSB-PA22969A0Rb	Human, Mouse	ELISA, WB, IHC
<u>MGAT4B Antibody</u>	CSB-PA890781LA01HU	Human	ELISA, WB, IHC
<u>MGAT5 Antibody</u>	CSB-PA600992LA01HU	Human	ELISA, IF
<u>MglI Antibody</u>	CSB-PA013787LA01RA	Rat, Mouse	ELISA, WB
<u>MGME1 Antibody</u>	CSB-PA880084LA01HU	Human	ELISA, WB, IHC
<u>MGMT Antibody</u>	CSB-PA11859A0Rb	Human	ELISA, IHC
<u>MGP Antibody</u>	CSB-PA013789LA01HU	Human	ELISA, IHC
<u>mgsA Antibody</u>	CSB-PA507015HA01BQK	Bacillus cereus	ELISA



Product Name	Code	Species Reactivity	Applications
<u>MGST1 Antibody</u>	CSB-PA013791LA01HU	Human	ELISA, IHC
<u>MGST2 Antibody</u>	CSB-PA013792LA01HU	Human	ELISA, IHC, IF
<u>MGST3 Antibody</u>	CSB-PA013793LA01HU	Human	ELISA, WB, IF
<u>MHT1 Antibody</u>	CSB-PA14269A0Rb	Saccharomyces cerevisiae	ELISA, WB
<u>MIA Antibody</u>	CSB-PA613694LA01HU	Human	ELISA, IHC
<u>MIA2 Antibody</u>	CSB-PA822278LA01HU	Human	ELISA, IHC
<u>MIA3 Antibody</u>	CSB-PA711448LA01HU	Human	ELISA
<u>MIB1 Antibody</u>	CSB-PA13289A0Rb	Human	ELISA, WB, IHC, IF
<u>mib1 Antibody</u>	CSB-PA896152LA01DIL	Zebrafish	ELISA
<u>MIB2 Antibody</u>	CSB-PA839291LA01HU	Human	ELISA, IHC
<u>MICAL1 Antibody</u>	CSB-PA851555LA01HU	Human	ELISA, WB, IHC, IF
<u>MICAL2 Antibody</u>	CSB-PA013808LA01HU	Human	ELISA, IHC, IF
<u>MICU1 Antibody</u>	CSB-PA860652LA01HU	Human	ELISA, WB, IHC
<u>MICU2 Antibody</u>	CSB-PA818265LA01HU	Human	ELISA, IHC
<u>MID1IP1 Antibody</u>	CSB-PA873607LA01HU	Human	ELISA, IHC, IF
<u>MIEF1 Antibody</u>	CSB-PA873622LA01HU	Human, Mouse	ELISA, WB
<u>MIEN1 Antibody</u>	CSB-PA874795LA01HU	Human	ELISA, IHC
<u>MIER2 Antibody</u>	CSB-PA818692LA01HU	Human	ELISA, WB, IHC
<u>MIF Antibody</u>	CSB-PA06867A0Rb	Human, Rat	ELISA, WB, IHC, IF
<u>Mif Antibody</u>	CSB-PA013826LA01MO	Mouse	ELISA, WB
<u>MIF4GD Antibody</u>	CSB-PA013827LA01HU	Human	ELISA
<u>MILR1 Antibody</u>	CSB-PA801819LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>MINDY1 Antibody</u>	CSB-PA843272LA01HU	Human, Mouse	ELISA, WB
<u>Minor capsid protein VP2 Antibody</u>	CSB-PA355948LA01JAK	JC polyomavirus	ELISA
<u>Minpp1 Antibody</u>	CSB-PA013831LA01RA	Rat	ELISA
<u>Minpp1 Antibody</u>	CSB-PA013831YA01RA	Rat	ELISA
<u>MINPP1 Antibody</u>	CSB-PA891977EA01HU	Human	ELISA, WB, IHC
<u>MIOS Antibody</u>	CSB-PA865167LA01HU	Human	ELISA, IHC
<u>MIOX Antibody</u>	CSB-PA866220LA01HU	Human, Mouse	ELISA, WB, IP
<u>mip Antibody</u>	CSB-PA401660HA01LLJ	Legionella pneumophila	ELISA
<u>MIP Antibody</u>	CSB-PA013834LA01HU	Human	ELISA, IHC, IF
<u>MIS12 Antibody</u>	CSB-PA875639LA01HU	Human	ELISA, WB, IHC
<u>mitfa Antibody</u>	CSB-PA20909A0Rb	Zebrafish	ELISA
<u>MKI67 Antibody</u>	CSB-PA11617A0Rb	Human	ELISA, IHC
<u>MKL1 Antibody</u>	CSB-PA853375LA01HU	Human	ELISA, IF
<u>MKL2 Antibody</u>	CSB-PA891962LA01HU	Human, Mouse	ELISA, WB, IHC
<u>MKLN1 Antibody</u>	CSB-PA890926LA01HU	Human, Rat	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>MKRN1 Antibody</u>	CSB-PA883390LA01HU	Human	ELISA, IHC
<u>MKRN2 Antibody</u>	CSB-PA872420LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>MKX Antibody</u>	CSB-PA818257LA01HU	Human	ELISA, WB, IHC, IF
<u>MLANA Antibody</u>	CSB-PA621968LA01HU	Human	ELISA, WB, IHC
<u>MLC1 Antibody</u>	CSB-PA613581LA01HU	Human	ELISA, IHC, IF
<u>MLH1 Antibody</u>	CSB-PA17009A0Rb	Human	ELISA, WB, IHC, IF
<u>MLKL Antibody</u>	CSB-PA850851LA01HU	Human	ELISA, IHC
<u>MLLT10 Antibody</u>	CSB-PA014632LA01HU	Human	ELISA, IHC
<u>MLLT11 Antibody</u>	CSB-PA618752LA01HU	Human	ELISA, IF
<u>MLLT4 Antibody</u>	CSB-PA014636LA01HU	Human	ELISA, IHC, IF
<u>MLN Antibody</u>	CSB-PA11109A0Rb	Human	ELISA, IHC
<u>MLPH Antibody</u>	CSB-PA014640LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>MLST8 Antibody</u>	CSB-PA863130LA01HU	Human	ELISA, IHC, IF
<u>MLX Antibody</u>	CSB-PA014643LA01HU	Human	ELISA, IHC, IF
<u>MLXIP Antibody</u>	CSB-PA862054LA01HU	Human	ELISA, IHC
<u>MLYCD Antibody</u>	CSB-PA014646LA01HU	Human	ELISA, WB, IHC, IF
<u>MMADHC Antibody</u>	CSB-PA875665LA01HU	Human	ELISA, IHC, IF
<u>MMD2 Antibody</u>	CSB-PA814229LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>MMGT1 Antibody</u>	CSB-PA014655LA01HU	Human	ELISA, WB, IF
<u>MMP1 Antibody</u>	CSB-PA07009A0Rb	Human	ELISA, WB, IHC
<u>MMP10 Antibody</u>	CSB-PA07019A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>MMP11 Antibody</u>	CSB-PA12369A0Rb	Human	ELISA, IF
<u>MMP12 Antibody</u>	CSB-PA12379A0Rb	Human, Mouse	ELISA, WB, IHC
<u>MMP13 Antibody</u>	CSB-PA07029A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>Mmp13 Antibody</u>	CSB-PA17739A0Rb	Mouse	ELISA, IHC
<u>MMP14 Antibody</u>	CSB-PA12389A0Rb	Human	ELISA, IHC
<u>Mmp2 Antibody</u>	CSB-PA014666LA01RA	Rat	ELISA
<u>MMP2 Antibody</u>	CSB-PA06879A0Rb	Human	ELISA, IHC, IF
<u>MMP2 Antibody</u>	CSB-PA014666LA01HU	Human	ELISA, IHC
<u>MMP28 Antibody</u>	CSB-PA880950LA01HU	Human	ELISA, WB, IHC, IF
<u>MMP3 Antibody</u>	CSB-PA07449A0Rb	Human, Mouse	ELISA, WB, IHC
<u>Mmp3 Antibody</u>	CSB-PA17509A0Rb	Mouse	ELISA, WB, IHC
<u>Mmp7 Antibody</u>	CSB-PA014677EA01RA	Rat, Mouse, Human	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>MMP7 Antibody</u>	CSB-PA07455A0Rb	Human	ELISA, IHC
<u>MMP8 Antibody</u>	CSB-PA07469A0Rb	Human, Mouse	ELISA, WB, IF
<u>MNAT1 Antibody</u>	CSB-PA014686LA01HU	Human, Mouse	ELISA, WB, IHC
<u>MND1 Antibody</u>	CSB-PA880120LA01HU	Human	ELISA, WB, IHC
<u>MNDA Antibody</u>	CSB-PA014688LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>MNX1 Antibody</u>	CSB-PA014692LA01HU	Human	ELISA, WB
<u>MOB1A Antibody</u>	CSB-PA872519HA01HU	Human	ELISA, WB, IHC, IF, IP
<u>MOB1B Antibody</u>	CSB-PA762574HA01HU	Human	ELISA, IHC, IF, IP
<u>MOBP Antibody</u>	CSB-PA014701LA01HU	Human	ELISA, WB, IF
<u>MOCS2 Antibody</u>	CSB-PA014706LA01HU	Human	ELISA, WB
<u>MOCS2 Antibody</u>	CSB-PA529746LA01HU	Human	ELISA, IF
<u>MOGAT1 Antibody</u>	CSB-PA836275LA01HU	Human	ELISA, IHC, IF
<u>MOGAT2 Antibody</u>	CSB-PA663570LA01HU	Human	ELISA, IHC
<u>MON1A Antibody</u>	CSB-PA801242LA01HU	Human	ELISA, WB, IHC
<u>MORC1 Antibody</u>	CSB-PA803127LA01HU	Human	ELISA, IHC, IF
<u>MORC2 Antibody</u>	CSB-PA896555LA01HU	Human	ELISA, WB, IHC
<u>MORN3 Antibody</u>	CSB-PA747476LA01HU	Human, Mouse	ELISA, WB
<u>MORN5 Antibody</u>	CSB-PA736208LA01HU	Human, Mouse	ELISA, WB, IHC
<u>MOS Antibody</u>	CSB-PA014733LA01HU	Human	ELISA, WB, IHC
<u>MOV10 Antibody</u>	CSB-PA862068LA01HU	Human	ELISA, WB, IHC
<u>MOV10L1 Antibody</u>	CSB-PA871604LA01HU	Human	ELISA, WB, IHC, IF
<u>MP68 Antibody</u>	CSB-PA00675A0Rb	Human	ELISA, WB, IHC, IF
<u>MPC2 Antibody</u>	CSB-PA002816LA01HU	Human, Mouse	ELISA, WB, IHC
<u>MPDU1 Antibody</u>	CSB-PA014744LA01HU	Human	ELISA, IHC, IF
<u>MPEG1 Antibody</u>	CSB-PA642654LA01HU	Human	ELISA, IF
<u>MPI Antibody</u>	CSB-PA014754LA01HU	Human	ELISA, IF
<u>MPL Antibody</u>	CSB-PA014755LA01HU	Human	ELISA, IHC, IF
<u>MPN_083 Antibody</u>	CSB-PA303444LA01Mlw	Mycoplasma pneumoniae	ELISA, WB
<u>MPO Antibody</u>	CSB-PA06899A0Rb	Human	ELISA, IHC, IF
<u>MPP1 Antibody</u>	CSB-PA014758LA01HU	Human	ELISA, IHC
<u>MPP3 Antibody</u>	CSB-PA615679LA01HU	Human	ELISA, IHC, IF
<u>MPP5 Antibody</u>	CSB-PA847626LA01HU	Human	ELISA, IF
<u>MPPED1 Antibody</u>	CSB-PA014767LA01HU	Human	ELISA, IHC, IF
<u>MPPED2 Antibody</u>	CSB-PA618997LA01HU	Human, Mouse	ELISA, WB
<u>MPST Antibody</u>	CSB-PA014770LA01HU	Human	ELISA, IHC, IF
<u>mpt51 Antibody</u>	CSB-PA14937A0Rb	Mycobacterium tuberculosis	ELISA
<u>mpt63 Antibody</u>	CSB-PA14907A0Rb	Mycobacterium	ELISA



Product Name	Code	Species Reactivity	Applications
		tuberculosis	
<u>mpt63 Antibody</u>	CSB-PA358712YA01MVZ	Mycobacterium tuberculosis	ELISA, WB
<u>mpt64 Antibody</u>	CSB-PA14949A0Rb	Mycobacterium tuberculosis	ELISA, WB
<u>mpt64 Antibody</u>	CSB-PA358713YA01MVZ	Mycobacterium tuberculosis	ELISA
<u>mpt64 Antibody</u>	CSB-PA14947A0Rb	Mycobacterium tuberculosis	ELISA, WB
<u>MPV17 Antibody</u>	CSB-PA014771LA01HU	Human, Mouse	ELISA, WB, IF
<u>MPZ Antibody</u>	CSB-PA07844A0Rb	Human	ELISA, WB
<u>MPZ Antibody</u>	CSB-PA07845A0Rb	Human	ELISA
<u>MPZL1 Antibody</u>	CSB-PA014775LA01HU	Human	ELISA, IHC, IF
<u>MPZL2 Antibody</u>	CSB-PA014776LA01HU	Human	ELISA, IHC
<u>MPZL3 Antibody</u>	CSB-PA740916LA01HU	Human	ELISA, IF
<u>MR1 Antibody</u>	CSB-PA853268LA01HU	Human	ELISA, WB, IHC, IF
<u>MRAP Antibody</u>	CSB-PA851545LA01HU	Human	ELISA, IF
<u>MRAP2 Antibody</u>	CSB-PA853421LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>MRC1 Antibody</u>	CSB-PA014783LA01HU	Human	ELISA
<u>mrcA Antibody</u>	CSB-PA365805LA01ENV	Escherichia coli	ELISA
<u>mrcA Antibody</u>	CSB-PA363014LA01NGG	Neisseria meningitidis serogroup B	ELISA
<u>MRE11A Antibody</u>	CSB-PA014786EA01HU	Human	ELISA, IHC, IF
<u>MREG Antibody</u>	CSB-PA014788LA01HU	Human, Mouse	ELISA, WB
<u>MRFAP1 Antibody</u>	CSB-PA896920EA01HU	Human	ELISA
<u>MRFAP1L1 Antibody</u>	CSB-PA014790EA01HU	Human	ELISA, WB
<u>MRGPRX2 Antibody</u>	CSB-PA839361LA01HU	Human, Mouse	ELISA, WB, IF
<u>MRGPRX3 Antibody</u>	CSB-PA850382LA01HU	Human	ELISA, IF
<u>MRJP1 Antibody</u>	CSB-PA522725EA01DNK	Apis mellifera	ELISA
<u>mrkA Antibody</u>	CSB-PA318939LA01KBG	Klebsiella pneumoniae	ELISA
<u>MRM2 Antibody</u>	CSB-PA883420LA01HU	Human	ELISA, IHC
<u>MRM3 Antibody</u>	CSB-PA867196LA01HU	Human	ELISA, IHC
<u>MRNIP Antibody</u>	CSB-PA747338LA01HU	Human, Mouse	ELISA, WB
<u>MRO Antibody</u>	CSB-PA874842LA01HU	Human	ELISA, WB, IF
<u>MRPL1 Antibody</u>	CSB-PA874841HA01HU	Human	ELISA, WB, IHC
<u>MRPL11 Antibody</u>	CSB-PA896709LA01HU	Human	ELISA, WB, IF
<u>MRPL14 Antibody</u>	CSB-PA014819LA01HU	Human, Mouse	ELISA, WB, IHC
<u>MRPL16 Antibody</u>	CSB-PA882128LA01HU	Human	ELISA, IHC, IF
<u>MRPL18 Antibody</u>	CSB-PA861987LA01HU	Human	ELISA, WB, IHC
<u>MRPL19 Antibody</u>	CSB-PA014826EA01HU	Human, Mouse	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>MRPL2 Antibody</u>	CSB-PA735972LA01HU	Human	ELISA, IHC
<u>MRPL20 Antibody</u>	CSB-PA863153LA01HU	Human	ELISA, IHC
<u>MRPL21 Antibody</u>	CSB-PA773581LA01HU	Human	ELISA, IHC
<u>MRPL22 Antibody</u>	CSB-PA878907LA01HU	Human	ELISA, WB, IHC
<u>MRPL27 Antibody</u>	CSB-PA868378LA01HU	Human	ELISA, IHC
<u>MRPL3 Antibody</u>	CSB-PA014836LA01HU	Human	ELISA, WB, IHC
<u>MRPL30 Antibody</u>	CSB-PA819459LA01HU	Human	ELISA
<u>MRPL33 Antibody</u>	CSB-PA014839LA01HU	Human	ELISA, IHC
<u>MRPL36 Antibody</u>	CSB-PA885784LA01HU	Human	ELISA
<u>MRPL38 Antibody</u>	CSB-PA856920LA01HU	Human	ELISA, IHC
<u>MRPL41 Antibody</u>	CSB-PA811620LA01HU	Human	ELISA, IHC
<u>MRPL42 Antibody</u>	CSB-PA014852EA01HU	Human	ELISA
<u>MRPL44 Antibody</u>	CSB-PA864007LA01HU	Human	ELISA, WB, IHC
<u>MRPL46 Antibody</u>	CSB-PA880959LA01HU	Human	ELISA
<u>MRPL48 Antibody</u>	CSB-PA836217LA01HU	Human, Mouse	ELISA, WB, IHC
<u>MRPL49 Antibody</u>	CSB-PA613399LA01HU	Human, Mouse	ELISA, WB, IHC
<u>MRPL50 Antibody</u>	CSB-PA839804LA01HU	Human	ELISA, IHC
<u>MRPL51 Antibody</u>	CSB-PA687045LA01HU	Human	ELISA, WB, IHC
<u>MRPL53 Antibody</u>	CSB-PA853410LA01HU	Human	ELISA, IHC
<u>MRPL54 Antibody</u>	CSB-PA014871LA01HU	Human	ELISA, IHC
<u>MRPL55 Antibody</u>	CSB-PA014872EA01HU	Human	ELISA, WB, IHC
<u>MRPL58 Antibody</u>	CSB-PA613490HA01HU	Human	ELISA, WB, IHC
<u>MRPL9 Antibody</u>	CSB-PA04864A0Rb	Human	ELISA, WB
<u>MRPS10 Antibody</u>	CSB-PA014875LA01HU	Human	ELISA, IHC
<u>MRPS14 Antibody</u>	CSB-PA014882LA01HU	Human	ELISA
<u>MRPS18A Antibody</u>	CSB-PA882117LA01HU	Human	ELISA, WB
<u>MRPS18C Antibody</u>	CSB-PA014893LA01HU	Human	ELISA
<u>MRPS2 Antibody</u>	CSB-PA897102LA01HU	Human	ELISA, WB
<u>MRPS28 Antibody</u>	CSB-PA896495LA01HU	Human	ELISA, IHC
<u>MRPS31 Antibody</u>	CSB-PA846089LA01HU	Human	ELISA, WB, IHC
<u>MRPS36 Antibody</u>	CSB-PA014927LA01HU	Human	ELISA, IHC
<u>MRPS5 Antibody</u>	CSB-PA014934LA01HU	Human	ELISA, WB
<u>MRRF Antibody</u>	CSB-PA850282LA01HU	Human	ELISA, IHC
<u>Ms4a1 Antibody</u>	CSB-PA015007LA01MO	Mouse, Human	ELISA, WB
<u>MS4A2 Antibody</u>	CSB-PA015013LA01HU	Human	ELISA, IHC, IF
<u>MS4A3 Antibody</u>	CSB-PA856947LA01HU	Human	ELISA, IF
<u>MS4A4A Antibody</u>	CSB-PA822238LA01HU	Human	ELISA, IHC, IF
<u>MS4A6A Antibody</u>	CSB-PA875659LA01HU	Human	ELISA, IF
<u>MS4A6E Antibody</u>	CSB-PA856919LA01HU	Human	ELISA, IHC
<u>MS4A7 Antibody</u>	CSB-PA875636LA01HU	Human	ELISA, IHC, IF
<u>MS4A8 Antibody</u>	CSB-PA883629LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>MSA2 Antibody</u>	CSB-PA343481DA01EWP	Plasmodium falciparum	ELISA
<u>MSC Antibody</u>	CSB-PA015025LA01HU	Human	ELISA, IHC
<u>mscL Antibody</u>	CSB-PA363948LA01EOD	Escherichia coli O157:H7	ELISA
<u>MSH2 Antibody</u>	CSB-PA16889A0Rb	Human, Mouse	ELISA, WB, IHC, IP
<u>MSH4 Antibody</u>	CSB-PA015030LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>MSH6 Antibody</u>	CSB-PA015032HA01HU	Human	ELISA, IHC, IF
<u>MSI1 Antibody</u>	CSB-PA015033LA01HU	Human	ELISA, IHC, IF
<u>MSL2 Antibody</u>	CSB-PA875722LA01HU	Human	ELISA, WB
<u>Msln Antibody</u>	CSB-PA015044LA01MO	Mouse	ELISA
<u>MSLN Antibody</u>	CSB-PA10159A0Rb	Human	ELISA, IHC, IF, IP
<u>MSMB Antibody</u>	CSB-PA015046LA01HU	Human	ELISA, IHC
<u>MSN Antibody</u>	CSB-PA01189A0Rb	Human	ELISA, WB, IHC, IF
<u>MSR1 Antibody</u>	CSB-PA015050LA01HU	Human, Mouse	ELISA, WB, IHC
<u>MSRA Antibody</u>	CSB-PA883456LA01HU	Human	ELISA, IHC
<u>MSRB1 Antibody</u>	CSB-PA882163LA01HU	Human	ELISA, IHC, IF
<u>MSS51 Antibody</u>	CSB-PA681195LA01HU	Human	ELISA, WB
<u>MST1 Antibody</u>	CSB-PA015055LA01HU	Human	ELISA, IHC
<u>MST1R Antibody</u>	CSB-PA015056LA01HU	Human	ELISA, IF
<u>MSTN Antibody</u>	CSB-PA11869A0Rb	Human	ELISA, WB, IHC
<u>MSTO1 Antibody</u>	CSB-PA880106LA01HU	Human	ELISA, IHC
<u>MSX1 Antibody</u>	CSB-PA015068LA01HU	Human	ELISA, IHC
<u>MSX2 Antibody</u>	CSB-PA015069LA01HU	Human, Rat	ELISA, WB, IHC
<u>MT1A Antibody</u>	CSB-PA17444A0Rb	Human	ELISA, IHC
<u>MT1B Antibody</u>	CSB-PA17419A0Rb	Human	ELISA
<u>MT1E Antibody</u>	CSB-PA17364A0Rb	Human	ELISA, IHC, IF
<u>MT1F Antibody</u>	CSB-PA17409A0Rb	Human	ELISA
<u>MT1G Antibody</u>	CSB-PA17384A0Rb	Human	ELISA, IHC
<u>MT1H Antibody</u>	CSB-PA17399A0Rb	Human	ELISA
<u>MT1M Antibody</u>	CSB-PA17429A0Rb	Human	ELISA
<u>MT1X Antibody</u>	CSB-PA17374A0Rb	Human	ELISA, IHC
<u>MT2731 Antibody</u>	CSB-PA300323LA01MVZ	Mycobacterium tuberculosis	ELISA
<u>MT2A Antibody</u>	CSB-PA17354A0Rb	Human	ELISA, IHC
<u>MT3 Antibody</u>	CSB-PA04639A0Rb	Human	ELISA, IHC, IF
<u>MT4 Antibody</u>	CSB-PA17459A0Rb	Human	ELISA
<u>MTAP Antibody</u>	CSB-PA622639LA01HU	Human, Mouse	ELISA, WB, IF
<u>MTBP Antibody</u>	CSB-PA856921LA01HU	Human	ELISA, IHC
<u>MTCH2 Antibody</u>	CSB-PA897586LA01HU	Human	ELISA, IHC
<u>MT-CO1 Antibody</u>	CSB-PA015072LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>MT-CO2 Antibody</u>	CSB-PA015073LA01HU	Human	ELISA, WB, IHC
<u>MTERF1 Antibody</u>	CSB-PA811872LA01HU	Human	ELISA, WB
<u>MTERF4 Antibody</u>	CSB-PA800224LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>MTHFD2L Antibody</u>	CSB-PA872521LA01HU	Human, Mouse	ELISA, WB
<u>MTHFR Antibody</u>	CSB-PA015158LA01HU	Human	ELISA, WB, IHC, IF
<u>MTHFR Antibody</u>	CSB-PA015158YA01HU	Human	ELISA
<u>MTHFS Antibody</u>	CSB-PA015159LA01HU	Human	ELISA, WB, IHC, IF
<u>MTMR11 Antibody</u>	CSB-PA015176LA01HU	Human	ELISA
<u>MTMR12 Antibody</u>	CSB-PA880158LA01HU	Human	ELISA, IF
<u>MTMR2 Antibody</u>	CSB-PA613418LA01HU	Human	ELISA, WB, IHC, IF
<u>MTMR3 Antibody</u>	CSB-PA622675HA01HU	Human	ELISA, IHC, IF
<u>MTMR6 Antibody</u>	CSB-PA897519LA01HU	Human	ELISA
<u>MTMR7 Antibody</u>	CSB-PA896478LA01HU	Human, Mouse	ELISA, WB
<u>MTMR9 Antibody</u>	CSB-PA847211LA01HU	Human	ELISA, IF
<u>MT-ND1 Antibody</u>	CSB-PA015076LA01HU	Human, Mouse	ELISA, WB, IHC
<u>MT-ND5 Antibody</u>	CSB-PA015081LA01HU	Human	ELISA, WB
<u>MTNR1B Antibody</u>	CSB-PA015190LA01HU	Human	ELISA, IHC, IF
<u>MTOR Antibody</u>	CSB-PA008968HA01HU	Human	ELISA, IHC, IF
<u>MTR Antibody</u>	CSB-PA859939LA01HU	Human	ELISA, IHC, IP
<u>MTRF1 Antibody</u>	CSB-PA015197LA01HU	Human	ELISA, WB, IHC, IF
<u>MTRR Antibody</u>	CSB-PA890659EA01HU	Human	ELISA, WB, IHC
<u>MTSS1 Antibody</u>	CSB-PA015201LA01HU	Human	ELISA, WB, IHC
<u>MTTP Antibody</u>	CSB-PA015209LA01HU	Human	ELISA, IHC, IF
<u>mturn Antibody</u>	CSB-PA24389A0Rb	Zebrafish	ELISA, WB
<u>MUC1 Antibody</u>	CSB-PA11319A0Rb	Human	ELISA, IHC, IF
<u>MUC16 Antibody</u>	CSB-PA10657A0Rb	Human	ELISA, IHC
<u>MUC2 Antibody</u>	CSB-PA015222LA01HU	Human	ELISA
<u>MUC20 Antibody</u>	CSB-PA822693LA01HU	Human	ELISA, WB, IHC, IF
<u>Muc5ac Antibody</u>	CSB-PA17229A0Rb	Mouse	ELISA
<u>Mug1 Antibody</u>	CSB-PA328826LA01MO	Mouse	ELISA
<u>Mup11 Antibody</u>	CSB-PA361419LA01MO	Mouse	ELISA, WB
<u>Mup11 Antibody</u>	CSB-PA361419YA01MO	Mouse	ELISA
<u>Mup2 Antibody</u>	CSB-PA319316HA01MO	Mouse	ELISA
<u>Mup2 Antibody</u>	CSB-PA319316YA01MO	Mouse	ELISA
<u>Mup6 Antibody</u>	CSB-PA360881HA01MO	Mouse	ELISA, IF
<u>MUS81 Antibody</u>	CSB-PA836272LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>Mut Antibody</u>	CSB-PA015243YA01MO	Mouse	ELISA
<u>mutL Antibody</u>	CSB-PA543397HA01ENV	Escherichia coli	ELISA, WB
<u>mutT Antibody</u>	CSB-PA357355HA01ENV	Escherichia coli	ELISA
<u>MVD Antibody</u>	CSB-PA015246LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>MVP Antibody</u>	CSB-PA015248LA01HU	Human, Rat	ELISA, WB
<u>MX1 Antibody</u>	CSB-PA015249LA01HU	Human, Mouse	ELISA, WB, IHC, IF, IP
<u>MX2 Antibody</u>	CSB-PA015250LA01HU	Human	ELISA, IHC, IF
<u>MXD3 Antibody</u>	CSB-PA015252LA01HU	Human	ELISA, IHC
<u>MXI1 Antibody</u>	CSB-PA015254LA01HU	Human	ELISA, IF
<u>MYADM Antibody</u>	CSB-PA015258LA01HU	Human	ELISA, IHC, IF
<u>MYB75 Antibody</u>	CSB-PA863725LA01DOA	Arabidopsis thaliana	ELISA, WB
<u>MYBL2 Antibody</u>	CSB-PA015264LA01HU	Human	ELISA, IHC, IF
<u>MYBPC1 Antibody</u>	CSB-PA717389LA01HU	Human	ELISA, IHC, IF
<u>MYBPC2 Antibody</u>	CSB-PA618884LA01HU	Human	ELISA, IHC, IF
<u>Mybpc3 Antibody</u>	CSB-PA015267YA01RA	Rat	ELISA
<u>MYBPC3 Antibody</u>	CSB-PA613524HA01HU	Human, Mouse	ELISA, WB
<u>MYBPH Antibody</u>	CSB-PA618765LA01HU	Human	ELISA, IF
<u>Myc Antibody</u>	CSB-PA015270EA01MO	Mouse	ELISA
<u>MYC Antibody</u>	CSB-PA02589A0Rb	Human, Mouse	ELISA, WB, IHC
<u>MYCBP Antibody</u>	CSB-PA858710HA01HU	Human	ELISA, IHC
<u>mycn Antibody</u>	CSB-PA868688LA01DIL	Zebrafish	ELISA, WB
<u>MYD88 Antibody</u>	CSB-PA859945LA01HU	Human	ELISA, IF
<u>MYEF2 Antibody</u>	CSB-PA015283LA01HU	Human	ELISA, IHC, IF
<u>MYF5 Antibody</u>	CSB-PA015287LA01HU	Human, Mouse, Rat	ELISA, WB
<u>MYF6 Antibody</u>	CSB-PA015288LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>MYH11 Antibody</u>	CSB-PA015291LA01HU	Human	ELISA, IHC
<u>MYH14 Antibody</u>	CSB-PA773786LA01HU	Human	ELISA, IHC, IF
<u>MYH2 Antibody</u>	CSB-PA892458LA01HU	Human	ELISA, WB
<u>MYH3 Antibody</u>	CSB-PA015297LA01HU	Human	ELISA, IHC, IF
<u>MYH4 Antibody</u>	CSB-PA22719A0Rb	Human	ELISA, IHC, IF
<u>MYH6 Antibody</u>	CSB-PA015299HA01HU	Human	ELISA, IHC
<u>MYH7 Antibody</u>	CSB-PA015300LA01HU	Human	ELISA, IHC
<u>MYH7B Antibody</u>	CSB-PA015301LA01HU	Human	ELISA, IHC
<u>MYH9 Antibody</u>	CSB-PA01405A0Rb	Human	ELISA, IHC, IF
<u>MYL2 Antibody</u>	CSB-PA10059A0Rb	Human	ELISA, IHC, IF
<u>MYL6 Antibody</u>	CSB-PA05154A0Rb	Human	ELISA, IHC, IF
<u>MYL9 Antibody</u>	CSB-PA015318LA01HU	Human	ELISA, IHC, IF
<u>MYLK Antibody</u>	CSB-PA618995LA01HU	Human	ELISA, IHC, IF
<u>MYLPF Antibody</u>	CSB-PA856892HA01HU	Human, Mouse	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>MYNN Antibody</u>	CSB-PA015330LA01HU	Human	ELISA, IHC
<u>MYO10 Antibody</u>	CSB-PA884628LA01HU	Human	ELISA, IHC, IF
<u>MYO18A Antibody</u>	CSB-PA856414LA01HU	Human	ELISA, IHC
<u>MYO18B Antibody</u>	CSB-PA815545LA01HU	Human	ELISA, IF
<u>MYO19 Antibody</u>	CSB-PA836223LA01HU	Human	ELISA, WB
<u>MYO1E Antibody</u>	CSB-PA621526LA01HU	Human	ELISA, IHC
<u>MYO1G Antibody</u>	CSB-PA015344LA01HU	Human	ELISA, IHC, IF
<u>MYO3B Antibody</u>	CSB-PA848841LA01HU	Human	ELISA, IHC
<u>MYO5A Antibody</u>	CSB-PA857430LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>MYO7A Antibody</u>	CSB-PA622657LA01HU	Human	ELISA, IHC, IF
<u>MYOC Antibody</u>	CSB-PA859950LA01HU	Human	ELISA, WB, IHC
<u>MYOD1 Antibody</u>	CSB-PA015358LA01HU	Human, Mouse	ELISA, WB
<u>MYOG Antibody</u>	CSB-PA015360LA01HU	Human	ELISA
<u>MYORG Antibody</u>	CSB-PA764657LA01HU	Human	ELISA, IHC
<u>MYOZ2 Antibody</u>	CSB-PA873609LA01HU	Human	ELISA, WB, IHC, IF
<u>MYRF Antibody</u>	CSB-PA897527LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>MYRIP Antibody</u>	CSB-PA839876LA01HU	Human	ELISA, IHC, IF
<u>MYSM1 Antibody</u>	CSB-PA716576LA01HU	Human	ELISA, IHC
<u>MYT1 Antibody</u>	CSB-PA015387LA01HU	Human	ELISA, IHC
<u>MZF1 Antibody</u>	CSB-PA015389LA01HU	Human	ELISA, WB, IHC
<u>N Antibody</u>	CSB-PA321352LA01RAI	Rabies virus	ELISA
<u>N Antibody</u>	CSB-PA321352YA01RAI	Rabies virus	ELISA
<u>N Antibody</u>	CSB-PA318492LA01JAM	Junin mammarenavirus	ELISA, WB
<u>N Antibody</u>	CSB-PA334828LA01PXE	Puumala virus	ELISA
<u>N Antibody</u>	CSB-PA318401LA01LNP	Lassa virus	ELISA
<u>N Antibody</u>	CSB-PA339612HA01BKT	Bovine respiratory syncytial virus	ELISA, WB
<u>N4BP2L2 Antibody</u>	CSB-PA838813LA01HU	Human	ELISA, IHC
<u>NAA15 Antibody</u>	CSB-PA863143LA01HU	Human, Rat, Mouse	ELISA, WB
<u>NAA38 Antibody</u>	CSB-PA871557LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>NAA60 Antibody</u>	CSB-PA875691LA01HU	Human	ELISA, IHC
<u>NAAA Antibody</u>	CSB-PA22919A0Rb	Human	ELISA, IHC
<u>NAALADL2 Antibody</u>	CSB-PA704507LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>NABP1 Antibody</u>	CSB-PA850257HA01HU	Human	ELISA, WB
<u>NACA2 Antibody</u>	CSB-PA887940LA01HU	Human	ELISA, WB
<u>NADK Antibody</u>	CSB-PA015409LA01HU	Human	ELISA, WB
<u>NADK2 Antibody</u>	CSB-PA676825LA01HU	Human, Rat, Mouse	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>NAGK Antibody</u>	CSB-PA883457HA01HU	Human	ELISA, WB
<u>NAGPA Antibody</u>	CSB-PA015417LA01HU	Human	ELISA, WB, IHC
<u>NAIF1 Antibody</u>	CSB-PA715122LA01HU	Human, Mouse	ELISA, WB, IHC
<u>NAIP Antibody</u>	CSB-PA615662LA01HU	Human	ELISA, WB, IHC, IF
<u>NALCN Antibody</u>	CSB-PA815590LA01HU	Human	ELISA, IF
<u>NAMPT Antibody</u>	CSB-PA015422EA01HU	Human	ELISA, WB, IHC, IF
<u>NANOG Antibody</u>	CSB-PA888008LA01HU	Human	ELISA, WB, IHC
<u>NANOG Monoclonal Antibody</u>	CSB-MA888008A0m	Human, Mouse, Rat	ELISA, WB, ICC, IF, FC
<u>NANOGP8 Antibody</u>	CSB-PA754125LA01HU	Human, Mouse, Rat	ELISA, WB
<u>NANOS1 Antibody</u>	CSB-PA837878LA01HU	Human	ELISA, IF
<u>NANOS2 Antibody</u>	CSB-PA015437LA01HU	Human	ELISA, IHC, IF
<u>NANOS3 Antibody</u>	CSB-PA015438LA01HU	Human	ELISA, IHC, IF
<u>NANP Antibody</u>	CSB-PA855041LA01HU	Human	ELISA, IHC
<u>NANS Antibody</u>	CSB-PA885699LA01HU	Human	ELISA, WB, IHC, IF
<u>NAP1L5 Antibody</u>	CSB-PA015445LA01HU	Human	ELISA, IHC
<u>NAPA Antibody</u>	CSB-PA015447LA01HU	Human	ELISA, WB, IHC, IF
<u>NAPEPLD Antibody</u>	CSB-PA764247LA01HU	Human, Mouse	ELISA, WB
<u>NAPRT Antibody</u>	CSB-PA015451EA01HU	Human, Mouse	ELISA, WB
<u>NAPSA Antibody</u>	CSB-PA11529A0Rb	Human, Mouse	ELISA, WB, IHC
<u>NARF Antibody</u>	CSB-PA015454LA01HU	Human, Mouse	ELISA, WB, IF
<u>NARS Antibody</u>	CSB-PA015459LA01HU	Human	ELISA, IHC, IF
<u>NARS2 Antibody</u>	CSB-PA856950LA01HU	Human	ELISA, IHC
<u>NAT1 Antibody</u>	CSB-PA015464LA01HU	Human	ELISA, IHC, IF
<u>NAT10 Antibody</u>	CSB-PA863939LA01HU	Human	ELISA, WB, IHC, IP
<u>NAT2 Antibody</u>	CSB-PA015471LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>NAT8 Antibody</u>	CSB-PA887033LA01HU	Human	ELISA, IHC, IF
<u>NBEAL2 Antibody</u>	CSB-PA751159LA01HU	Human	ELISA, IHC, IF
<u>NBPF3 Antibody</u>	CSB-PA872422LA01HU	Human	ELISA, WB, IHC
<u>NCAM2 Antibody</u>	CSB-PA015512LA01HU	Human	ELISA, WB, IHC, IF
<u>NCDN Antibody</u>	CSB-PA015524LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>NCF1 Antibody</u>	CSB-PA015525LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>NCF2 Antibody</u>	CSB-PA015528LA01HU	Human	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>NCF4 Antibody</u>	CSB-PA617999LA01HU	Human	ELISA, IHC, IF
<u>NCK2 Antibody</u>	CSB-PA015531HA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>NCLN Antibody</u>	CSB-PA822169LA01HU	Human	ELISA, IHC
<u>NCMAP Antibody</u>	CSB-PA735965LA01HU	Human	ELISA, IHC, IF
<u>NCOA1 Antibody</u>	CSB-PA015539HA01HU	Human	ELISA, IHC
<u>NCOA5 Antibody</u>	CSB-PA015544LA01HU	Human, Mouse	ELISA, WB, IF
<u>NCOA6 Antibody</u>	CSB-PA015545LA01HU	Human	ELISA, IHC
<u>NCOR1 Antibody</u>	CSB-PA015547LA01HU	Human	ELISA, IF
<u>NCR1 Antibody</u>	CSB-PA015549LA01HU	Human	ELISA, IHC, IF
<u>NCR3 Antibody</u>	CSB-PA015551LA01HU	Human	ELISA, IHC, IF
<u>NCR3LG1 Antibody</u>	CSB-PA737873LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ND-75 Antibody</u>	CSB-PA821986LA01DLU	Drosophila melanogaster	ELISA
<u>NDC1 Antibody</u>	CSB-PA874808LA01HU	Human	ELISA, IHC
<u>NDC80 Antibody</u>	CSB-PA015601LA01HU	Human	ELISA, IHC
<u>NDFIP1 Antibody</u>	CSB-PA887119LA01HU	Human	ELISA, WB, IHC, IF
<u>NDN Antibody</u>	CSB-PA859515LA01HU	Human	ELISA, IHC, IF
<u>NDP Antibody</u>	CSB-PA015609LA01HU	Human	ELISA, IHC
<u>ndr1 Antibody</u>	CSB-PA24539A0Rb	Zebrafish	ELISA, WB
<u>NDRG1 Antibody</u>	CSB-PA835678LA01HU	Human	ELISA, WB, IHC, IF
<u>NDRG4 Antibody</u>	CSB-PA890932LA01HU	Human	ELISA, WB, IHC, IF
<u>NDST1 Antibody</u>	CSB-PA015614LA01HU	Human	ELISA, WB, IHC
<u>NDUFA1 Antibody</u>	CSB-PA01665A0Rb	Human	ELISA, IHC, IF
<u>NDUFA2 Antibody</u>	CSB-PA02989A0Rb	Human	ELISA, IHC
<u>NDUFA3 Antibody</u>	CSB-PA015624EA01HU	Human	ELISA, IHC, IF
<u>NDUFA4L2 Antibody</u>	CSB-PA878874LA01HU	Human	ELISA, IHC
<u>NDUFA5 Antibody</u>	CSB-PA614988LA01HU	Human	ELISA, IHC
<u>NDUFA6 Antibody</u>	CSB-PA015631EA01HU	Human	ELISA, IHC
<u>NDUFAF1 Antibody</u>	CSB-PA897101LA01HU	Human	ELISA, IHC
<u>NDUFAF2 Antibody</u>	CSB-PA836643LA01HU	Human	ELISA, WB, IHC
<u>NDUFAF3 Antibody</u>	CSB-PA880103LA01HU	Human	ELISA, IHC
<u>NDUFB11 Antibody</u>	CSB-PA868343LA01HU	Human	ELISA, IHC, IF
<u>NDUFB3 Antibody</u>	CSB-PA015645LA01HU	Human	ELISA, IHC
<u>NDUFB5 Antibody</u>	CSB-PA015652LA01HU	Human	ELISA, WB
<u>NDUFB6 Antibody</u>	CSB-PA015653LA01HU	Human, Mouse	ELISA, WB, IHC, IP
<u>NDUFB7 Antibody</u>	CSB-PA01594A0Rb	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>NDUFB8 Antibody</u>	CSB-PA015655LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>NDUFB9 Antibody</u>	CSB-PA897317LA01HU	Human, Rat	ELISA, WB
<u>NDUFS2 Antibody</u>	CSB-PA015661LA01HU	Human, Rat, Mouse	ELISA, WB, IHC
<u>NDUFS4 Antibody</u>	CSB-PA015663LA01HU	Human	ELISA, WB, IHC
<u>NDUFS5 Antibody</u>	CSB-PA02735A0Rb	Human	ELISA, WB
<u>NDUFS8 Antibody</u>	CSB-PA015667LA01HU	Human, Rat	ELISA, WB
<u>NDUFV2 Antibody</u>	CSB-PA05385A0Rb	Human	ELISA, WB
<u>NDUFV3 Antibody</u>	CSB-PA015671LA01HU	Human	ELISA
<u>NEB Antibody</u>	CSB-PA324989LA01HU	Human	ELISA, IHC, IF
<u>NEBL Antibody</u>	CSB-PA015673LA01HU	Human	ELISA, IF
<u>NECAB1 Antibody</u>	CSB-PA843298LA01HU	Human	ELISA, IF
<u>NECTIN2 Antibody</u>	CSB-PA835690LA01HU	Human	ELISA, WB
<u>NECTIN4 Antibody</u>	CSB-PA822274LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>NEDD4L Antibody</u>	CSB-PA850411LA01HU	Human	ELISA, WB, IF
<u>NEDD8 Antibody</u>	CSB-PA618018DA01HU	Human	ELISA, IHC
<u>NEDD9 Antibody</u>	CSB-PA622771LA01HU	Human	ELISA, IHC
<u>nef Antibody</u>	CSB-PA27307A0Rb	Human immunodeficiency virus 1	ELISA, WB
<u>NEFL Antibody</u>	CSB-PA015688LA01HU	Human	ELISA, IHC, IF, IP
<u>NEGR1 Antibody</u>	CSB-PA015692LA01HU	Human	ELISA, WB, IHC, IF
<u>NEIL1 Antibody</u>	CSB-PA846619LA01HU	Human	ELISA, IHC
<u>NEK1 Antibody</u>	CSB-PA836281LA01HU	Human	ELISA, IHC, IF
<u>NEK10 Antibody</u>	CSB-PA747816LA01HU	Human	ELISA, IHC
<u>NEK11 Antibody</u>	CSB-PA844052LA01HU	Human	ELISA, WB, IHC
<u>NEK4 Antibody</u>	CSB-PA015702LA01HU	Human	ELISA, IF
<u>NEK5 Antibody</u>	CSB-PA740713LA01HU	Human	ELISA, IHC
<u>NEK7 Antibody</u>	CSB-PA840576LA01HU	Human	ELISA, IF
<u>NEK9 Antibody</u>	CSB-PA837435LA01HU	Human, Rat	ELISA, WB, IHC
<u>NELFB Antibody</u>	CSB-PA03765A0Rb	Human, Mouse	ELISA, WB, IHC, ChIP
<u>NELL2 Antibody</u>	CSB-PA859098LA01HU	Human	ELISA, IHC
<u>NENF Antibody</u>	CSB-PA892354NA01HU	Human	ELISA, IHC
<u>NEO1 Antibody</u>	CSB-PA852917LA01HU	Human	ELISA, IHC, IF
<u>NES Antibody</u>	CSB-PA015713YA01HU	Human	ELISA, IHC, IF
<u>NES Monoclonal Antibody</u>	CSB-MA0157131A0m	Human	ELISA, WB, IHC, IF, FC
<u>NET1 Antibody</u>	CSB-PA768775LA01HU	Human	ELISA
<u>NETO1 Antibody</u>	CSB-PA840571LA01HU	Human	ELISA, WB, IF
<u>NETO2 Antibody</u>	CSB-PA818759LA01HU	Human	ELISA, WB, IHC
<u>NEU1 Antibody</u>	CSB-PA17249A0Rb	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>NEU3 Antibody</u>	CSB-PA17267A0Rb	Human	ELISA, WB
<u>NEU4 Antibody</u>	CSB-PA848831LA01HU	Human	ELISA, IHC, IF
<u>NEURL1B Antibody</u>	CSB-PA015722LA01HU	Human	ELISA, IHC, IF
<u>NEUROD1 Antibody</u>	CSB-PA615700LA01HU	Human, Mouse	ELISA, WB, IF
<u>NEUROD2 Antibody</u>	CSB-PA614900LA01HU	Human	ELISA
<u>NEUROD6 Antibody</u>	CSB-PA846677LA01HU	Human, Mouse	ELISA, WB
<u>NEUROG2 Antibody</u>	CSB-PA887960LA01HU	Human, Rat	ELISA, WB
<u>NF1 Antibody</u>	CSB-PA015734LA01HU	Human	ELISA, IHC, IF
<u>NF2 Antibody</u>	CSB-PA015741LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>NFAM1 Antibody</u>	CSB-PA850874LA01HU	Human	ELISA, IHC
<u>NFAT5 Antibody</u>	CSB-PA015744LA01HU	Human	ELISA, IF
<u>NFATC1 Antibody</u>	CSB-PA17827A0Rb	Human	ELISA, IHC, IF
<u>NFATC3 Antibody</u>	CSB-PA615540LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>NFE2L1 Antibody</u>	CSB-PA015751LA01HU	Human	ELISA, IHC, IF
<u>NFE2L3 Antibody</u>	CSB-PA896721LA01HU	Human	ELISA, WB, IHC
<u>NFIC Antibody</u>	CSB-PA015756LA01HU	Human	ELISA, IHC, IF
<u>NFIX Antibody</u>	CSB-PA617941LA01HU	Human	ELISA, IHC
<u>NFKB1 Antibody</u>	CSB-PA10354A0Rb	Human	ELISA, WB, IHC, IF, ChIP
<u>NFKB2 Antibody</u>	CSB-PA10449A0Rb	Human	ELISA, IHC, IF
<u>NFKBID Antibody</u>	CSB-PA844092LA01HU	Human	ELISA, WB
<u>NFKBIL1 Antibody</u>	CSB-PA887001LA01HU	Human	ELISA, IHC
<u>nfuA Antibody</u>	CSB-PA540344LA01ENT	Escherichia coli	ELISA, WB
<u>nfuA Antibody</u>	CSB-PA540344EA11ENT	Escherichia coli	ELISA, WB
<u>NGB Antibody</u>	CSB-PA015776LA01HU	Human	ELISA
<u>NGEF Antibody</u>	CSB-PA850809LA01HU	Human	ELISA, IF
<u>Ngf Antibody</u>	CSB-PA16859A0Rb	Rat	ELISA, IP
<u>NGF Antibody</u>	CSB-PA643662LA01PI	Pig	ELISA
<u>NGF Antibody</u>	CSB-PA643662LA01HU	Human	ELISA, IHC
<u>NGFR Antibody</u>	CSB-PA10479A0Rb	Human, Mouse	ELISA, WB, IHC
<u>NHEJ1 Antibody</u>	CSB-PA881005LA01HU	Human	ELISA, IF
<u>NHLRC2 Antibody</u>	CSB-PA850853LA01HU	Human	ELISA, WB, IHC
<u>NHP2 Antibody</u>	CSB-PA873677LA01HU	Human	ELISA, IHC, IF
<u>NICN1 Antibody</u>	CSB-PA866273LA01HU	Human	ELISA, WB, IHC
<u>Nid1 Antibody</u>	CSB-PA015802LA01MO	Mouse, Rat	ELISA, WB
<u>Nid1 Antibody</u>	CSB-PA015802YA01MO	Mouse	ELISA
<u>NID1 Antibody</u>	CSB-PA015802LA01HU	Human	ELISA, IHC
<u>NIF3L1 Antibody</u>	CSB-PA880926LA01HU	Human	ELISA, WB, IHC
<u>nikR Antibody</u>	CSB-PA363929HA01EOD	Escherichia coli	ELISA
<u>NIM1K Antibody</u>	CSB-PA810298LA01HU	Human	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>NIN Antibody</u>	CSB-PA843155LA01HU	Human	ELISA, IHC
<u>NINJ1 Antibody</u>	CSB-PA856441LA01HU	Human	ELISA, IHC, IF
<u>NINJ2 Antibody</u>	CSB-PA885775LA01HU	Human	ELISA, IHC, IF
<u>NIP7 Antibody</u>	CSB-PA897077HA01HU	Human	ELISA, IHC, IF
<u>NIPBL Antibody</u>	CSB-PA750379LA01HU	Human	ELISA, WB, IHC
<u>NIPSNAP3B Antibody</u>	CSB-PA880093LA01HU	Human	ELISA, IHC
<u>NIT1 Antibody</u>	CSB-PA768228LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>NIT2 Antibody</u>	CSB-PA882075LA01HU	Human, Mouse	ELISA, WB
<u>NKAIN4 Antibody</u>	CSB-PA818240LA01HU	Human	ELISA, IHC, IF
<u>NKAPL Antibody</u>	CSB-PA702797NA01HU	Human	ELISA, WB
<u>NKD1 Antibody</u>	CSB-PA853369LA01HU	Human, Rat	ELISA, WB, IF
<u>NKG7 Antibody</u>	CSB-PA613690LA01HU	Human	ELISA, IF
<u>NKIRAS1 Antibody</u>	CSB-PA885767LA01HU	Human	ELISA, WB, IHC
<u>NKIRAS2 Antibody</u>	CSB-PA873696LA01HU	Human	ELISA, IHC
<u>NKTR Antibody</u>	CSB-PA015838LA01HU	Human	ELISA
<u>nkx2.5 Antibody</u>	CSB-PA20879A0Rb	Zebrafish	ELISA
<u>Nkx2-2 Antibody</u>	CSB-PA015842LA01MO	Mouse	ELISA
<u>NKX2-3 Antibody</u>	CSB-PA848391LA01HU	Human	ELISA, WB
<u>NKX2-8 Antibody</u>	CSB-PA015847LA01HU	Human	ELISA, IHC
<u>NKX3-2 Antibody</u>	CSB-PA015849LA01HU	Human, Rat	ELISA, WB
<u>NLE1 Antibody</u>	CSB-PA889125LA01HU	Human	ELISA
<u>NLGN1 Antibody</u>	CSB-PA839786LA01HU	Human	ELISA, IHC, IF
<u>NLGN3 Antibody</u>	CSB-PA873703LA01HU	Human	ELISA, IF
<u>NLGN4X Antibody</u>	CSB-PA843131LA01HU	Human	ELISA, IF
<u>NLN Antibody</u>	CSB-PA863161LA01HU	Human	ELISA, WB, IHC, IF
<u>NLRC5 Antibody</u>	CSB-PA771457LA01HU	Human	ELISA, IHC, IF
<u>NLRP1 Antibody</u>	CSB-PA871630EA01HU	Human	ELISA, IF
<u>NLRP10 Antibody</u>	CSB-PA773044LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>Nlrp3 Antibody</u>	CSB-PA823181LA01MO	Mouse	ELISA, IF
<u>NLRP3 Antibody</u>	CSB-PA015871LA01HU	Human	ELISA, IHC, IF
<u>NLRP6 Antibody</u>	CSB-PA015874LA01HU	Human	ELISA, IHC, IF
<u>NLRP7 Antibody</u>	CSB-PA840998LA01HU	Human	ELISA, IHC
<u>NLRX1 Antibody</u>	CSB-PA015878LA01HU	Human	ELISA, IHC
<u>NMBR Antibody</u>	CSB-PA015881LA01HU	Human	ELISA, WB
<u>NME1 Antibody</u>	CSB-PA02775A0Rb	Human	ELISA, IHC, IF
<u>NME5 Antibody</u>	CSB-PA015890LA01HU	Human	ELISA, IF
<u>NME6 Antibody</u>	CSB-PA015891LA01HU	Human	ELISA, WB, IHC
<u>NME7 Antibody</u>	CSB-PA897289LA01HU	Human, Mouse	ELISA, WB, IHC
<u>NMI Antibody</u>	CSB-PA618768LA01HU	Human, Rat	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>NMNAT2 Antibody</u>	CSB-PA866332HA01HU	Human	ELISA, IHC
<u>NMNAT3 Antibody</u>	CSB-PA839418LA01HU	Human, Mouse	ELISA, WB, IHC
<u>NMRAL1 Antibody</u>	CSB-PA867194LA01HU	Human	ELISA, IF
<u>NMT1 Antibody</u>	CSB-PA015900LA01HU	Human	ELISA, IHC
<u>NMUR1 Antibody</u>	CSB-PA884509LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>NNMT Antibody</u>	CSB-PA015906LA01HU	Human	ELISA, IF
<u>NOA1 Antibody</u>	CSB-PA836715LA01HU	Human, Mouse	ELISA, WB, IHC
<u>NOB1 Antibody</u>	CSB-PA892472LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>NOD2 Antibody</u>	CSB-PA881021LA01HU	Human	ELISA, IHC, IF
<u>NOG Antibody</u>	CSB-PA015917LA01HU	Human	ELISA, IHC
<u>NOL10 Antibody</u>	CSB-PA015918LA01HU	Human	ELISA, WB, IHC
<u>NOL3 Antibody</u>	CSB-PA015921HA01HU	Human	ELISA, IHC, IF
<u>NOL4L Antibody</u>	CSB-PA856989LA01HU	Human	ELISA, WB, IHC, IF
<u>NOL6 Antibody</u>	CSB-PA872506LA01HU	Human	ELISA, IHC
<u>NOMO1 Antibody</u>	CSB-PA015931LA01HU	Human	ELISA, WB, IHC, IF
<u>NOP10 Antibody</u>	CSB-PA873610LA01HU	Human	ELISA, WB, IHC
<u>NOP53 Antibody</u>	CSB-PA885778LA01HU	Human, Rat	ELISA, WB, IHC
<u>NOP58 Antibody</u>	CSB-PA896703LA01HU	Human	ELISA, WB, IF
<u>NOS1 Antibody</u>	CSB-PA14099A0Rb	Human, Mouse	ELISA, WB, IF
<u>NOS1AP Antibody</u>	CSB-PA015942LA01HU	Human	ELISA, IHC, IF
<u>NOS2 Antibody</u>	CSB-PA015943LA01HU	Human	ELISA, IF
<u>Nos3 Antibody</u>	CSB-PA16579A0Rb	Rat	ELISA
<u>NOSTRIN Antibody</u>	CSB-PA818236LA01HU	Human, Mouse	ELISA, WB, IHC
<u>nosZ Antibody</u>	CSB-PA703536LA01PBZ	Paracoccus denitrificans	ELISA
<u>NOTCH1 Antibody</u>	CSB-PA015949LA01HU	Human	ELISA, IHC, IF
<u>notch1a Antibody</u>	CSB-PA015949LA01DIL	Zebrafish	ELISA
<u>NOTCH2 Antibody</u>	CSB-PA015950LA01HU	Human	ELISA, IHC, IF
<u>NOTCH3 Antibody</u>	CSB-PA892153LA01HU	Human	ELISA, IHC
<u>notchl Antibody</u>	CSB-PA24369A0Rb	Zebrafish	ELISA, WB
<u>NOTUM Antibody</u>	CSB-PA747435LA01HU	Human, Mouse	ELISA, WB
<u>NOVA1 Antibody</u>	CSB-PA015957LA01HU	Human	ELISA, WB, IHC
<u>NOX1 Antibody</u>	CSB-PA015959LA01HU	Human	ELISA, IHC, IF
<u>NOX4 Antibody</u>	CSB-PA015961LA01HU	Human	ELISA, IHC, IF
<u>Nox4 Antibody</u>	CSB-PA015961LA01RA	Rat, Human	ELISA, WB
<u>NOXA1 Antibody</u>	CSB-PA015963LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>NOXO1 Antibody</u>	CSB-PA015964LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>NP Antibody</u>	CSB-PA356056LA01IFZ	Influenza A virus	ELISA
<u>NP Antibody</u>	CSB-PA356318HA01IJU	Influenza B virus	ELISA
<u>NP Antibody</u>	CSB-PA361330LA01IJK	Influenza B virus	ELISA
<u>NP Antibody</u>	CSB-PA529597LA01IMP	Influenza A virus	ELISA
<u>NP Antibody</u>	CSB-PA841983HA01RCI	Reston ebolavirus	ELISA, WB
<u>NP Antibody</u>	CSB-PA325769LA01IIM	Influenza A virus	ELISA, WB
<u>NPAS2 Antibody</u>	CSB-PA860777LA01HU	Human	ELISA, WB, IHC, IF
<u>NPAS4 Antibody</u>	CSB-PA814204LA01HU	Human	ELISA, WB, IHC
<u>NPBWR2 Antibody</u>	CSB-PA015973LA01HU	Human	ELISA, IF
<u>NPC1 Antibody</u>	CSB-PA890689LA01HU	Human	ELISA, WB, IHC
<u>NPC2 Antibody</u>	CSB-PA015976LA01HU	Human	ELISA, IHC, IF
<u>NPC2 Antibody</u>	CSB-PA634688LA01DO	Dog	ELISA
<u>NPFFR2 Antibody</u>	CSB-PA015983NA01HU	Human	ELISA, IHC
<u>NPFFR2 Antibody</u>	CSB-PA015983OA01HU	Human	ELISA, IF
<u>NPHP1 Antibody</u>	CSB-PA015985DA01HU	Human	ELISA, IHC, IF
<u>NPHP4 Antibody</u>	CSB-PA015987LA01HU	Human	ELISA, IHC, IF
<u>NPHS1 Antibody</u>	CSB-PA015988LA01HU	Human	ELISA, WB, IHC, IF
<u>NPL Antibody</u>	CSB-PA874830LA01HU	Human, Mouse	ELISA, WB, IF
<u>NPM2 Antibody</u>	CSB-PA774790LA01HU	Human	ELISA, IHC
<u>NPPB Antibody</u>	CSB-PA016021YA01HU	Human	ELISA
<u>NPPB Antibody</u>	CSB-PA07924A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>NPPB Antibody</u>	CSB-PA07926A0Rb	Human	ELISA, WB, IHC
<u>NPPB Antibody</u>	CSB-PA07927A0Rb	Human	ELISA
<u>NPR1 Antibody</u>	CSB-PA016023LA01HU	Human	ELISA
<u>NPR2 Antibody</u>	CSB-PA016024LA01HU	Human	ELISA
<u>Npr3 Antibody</u>	CSB-PA016025LA01RA	Rat	ELISA
<u>NPR3 Antibody</u>	CSB-PA016025LA01HU	Human	ELISA, IHC
<u>NPRL3 Antibody</u>	CSB-PA619643LA01HU	Human	ELISA, WB, IHC
<u>NPSR1 Antibody</u>	CSB-PA757818LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>NPTX1 Antibody</u>	CSB-PA614901LA01HU	Human	ELISA, IHC, IF
<u>NPY2R Antibody</u>	CSB-PA016036LA01HU	Human	ELISA, WB, IF
<u>NPY4R Antibody</u>	CSB-PA018593LA01HU	Human	ELISA, IHC
<u>NPY6R Antibody</u>	CSB-PA857447LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>NQO1 Antibody</u>	CSB-PA016039LA01HU	Human	ELISA, IHC, IF
<u>Nqo2 Antibody</u>	CSB-PA721470LA01RA	Rat, Mouse	ELISA, WB
<u>NQO2 Antibody</u>	CSB-PA016040LA01HU	Human, Mouse	ELISA, WB, IHC, IF, IP
<u>NR0B1 Antibody</u>	CSB-PA016041LA01HU	Human	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>NR0B2 Antibody</u>	CSB-PA624025LA01HU	Human	ELISA
<u>NR1D1 Antibody</u>	CSB-PA016043LA01HU	Human	ELISA, IHC
<u>NR1D2 Antibody</u>	CSB-PA623948LA01HU	Human	ELISA, WB
<u>NR1H2 Antibody</u>	CSB-PA016045HA01HU	Human	ELISA, IHC, IF
<u>NR1H3 Antibody</u>	CSB-PA619753LA01HU	Human, Mouse	ELISA, WB, IHC
<u>NR1H4 Antibody</u>	CSB-PA839406LA01HU	Human	ELISA, WB, IHC
<u>NR1I3 Antibody</u>	CSB-PA13019A0Rb	Human	ELISA, IHC
<u>NR2C2AP Antibody</u>	CSB-PA016053LA01HU	Human	ELISA, WB, IHC
<u>NR2E3 Antibody</u>	CSB-PA865585LA01HU	Human	ELISA, IHC
<u>NR2F1 Antibody</u>	CSB-PA016056LA01HU	Human	ELISA, WB
<u>NR2F2 Antibody</u>	CSB-PA016057LA01HU	Human	ELISA, IHC
<u>NR2F6 Antibody</u>	CSB-PA016058LA01HU	Human	ELISA, IHC
<u>NR3C1 Antibody</u>	CSB-PA016059LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>NR5A2 Antibody</u>	CSB-PA016066LA01HU	Human, Mouse	ELISA, WB, IHC
<u>NRARP Antibody</u>	CSB-PA016069LA01HU	Human	ELISA, IHC
<u>NRAS Antibody</u>	CSB-PA016070LA01HU	Human	ELISA, IHC, IF
<u>NRBP2 Antibody</u>	CSB-PA878883LA01HU	Human	ELISA, IHC, IF, IP
<u>NRCAM Antibody</u>	CSB-PA016074LA01HU	Human	ELISA, IHC, IF
<u>nrdA Antibody</u>	CSB-PA355522HA01EVN	Escherichia coli	ELISA
<u>nrdB Antibody</u>	CSB-PA303206HA01ENV	Escherichia coli	ELISA
<u>NRDC Antibody</u>	CSB-PA016075LA01HU	Human	ELISA, WB, IF
<u>nrdD Antibody</u>	CSB-PA15989A0Rb	Escherichia coli	ELISA
<u>NREP Antibody</u>	CSB-PA620971LA01HU	Human	ELISA, IHC, IF
<u>NRG1 Antibody</u>	CSB-PA09589A0Rb	Human	ELISA, IHC, IF
<u>NRG2 Antibody</u>	CSB-PA016078LA01HU	Human	ELISA, WB, IHC, IF
<u>NRG4 Antibody</u>	CSB-PA10569A0Rb	Human	ELISA, IHC, IF
<u>NRIP2 Antibody</u>	CSB-PA880083LA01HU	Human, Rat	ELISA, WB, IHC
<u>NRIP3 Antibody</u>	CSB-PA868283LA01HU	Human	ELISA, IHC
<u>NRL Antibody</u>	CSB-PA016086LA01HU	Human	ELISA, IF
<u>NRN1 Antibody</u>	CSB-PA868277LA01HU	Human	ELISA, IHC
<u>NRN1L Antibody</u>	CSB-PA670967LA01HU	Human	ELISA, IHC, IF
<u>NRP2 Antibody</u>	CSB-PA016092LA01HU	Human	ELISA, IHC, IF
<u>NRSN1 Antibody</u>	CSB-PA815587LA01HU	Human, Mouse	ELISA, WB, IF
<u>NRSN2 Antibody</u>	CSB-PA884422LA01HU	Human, Mouse	ELISA, WB, IF
<u>NRXN1 Antibody</u>	CSB-PA347832LA01HU	Human, Mouse	ELISA, WB, IF
<u>NRXN3 Antibody</u>	CSB-PA875730LA01HU	Human, Mouse	ELISA, WB
<u>NS Antibody</u>	CSB-PA356065LA01IE	Influenza A virus	ELISA
<u>NS Antibody</u>	CSB-PA604972LA01IER	Influenza A virus	ELISA, WB
<u>NS1 Antibody</u>	CSB-PA357083LA01HPM	Human parvovirus B19	ELISA, WB
<u>NSD2 Antibody</u>	CSB-PA026107LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>NSDHL Antibody</u>	CSB-PA618994LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>Nsf Antibody</u>	CSB-PA15457A0Rb	Rat	ELISA
<u>NSL1 Antibody</u>	CSB-PA839345LA01HU	Human	ELISA, IHC
<u>NSMAF Antibody</u>	CSB-PA856416LA01HU	Human	ELISA, IHC, IF
<u>NSMCE2 Antibody</u>	CSB-PA836260LA01HU	Human	ELISA, IHC
<u>NSUN2 Antibody</u>	CSB-PA626626LA01HU	Human, Rat, Mouse	ELISA, WB, IHC, IF
<u>NSUN3 Antibody</u>	CSB-PA867156LA01HU	Human	ELISA
<u>NSUN4 Antibody</u>	CSB-PA850270LA01HU	Human	ELISA, WB, IHC
<u>NSUN5 Antibody</u>	CSB-PA839386LA01HU	Human	ELISA, WB
<u>NSUN7 Antibody</u>	CSB-PA822769LA01HU	Human	ELISA, IHC
<u>NT5C1A Antibody</u>	CSB-PA866304YA01HU	Human	ELISA, IHC, IF
<u>NT5C2 Antibody</u>	CSB-PA23789A0Rb	Human, Mouse	ELISA, WB, IF
<u>NT5E Antibody</u>	CSB-PA10529A0Rb	Human	ELISA, WB, IHC
<u>NTAN1 Antibody</u>	CSB-PA822176LA01HU	Human	ELISA, IHC, IF
<u>NTF4 Antibody</u>	CSB-PA09929A0Rb	Human	ELISA, IHC, IF
<u>NTH1 Antibody</u>	CSB-PA18159A0Rb	Saccharomyces cerevisiae	ELISA
<u>NTHL1 Antibody</u>	CSB-PA016125LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>NTMT1 Antibody</u>	CSB-PA880111HA01HU	Human	ELISA
<u>NTN1 Antibody</u>	CSB-PA016127LA01HU	Human	ELISA, WB, IHC
<u>ntpcr Antibody</u>	CSB-PA24179A0Rb	Zebrafish	ELISA, WB
<u>NTRK1 Antibody</u>	CSB-PA016133LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>NTS Antibody</u>	CSB-PA016136HA01HU	Human, Mouse	ELISA, WB, IHC
<u>NTSR1 Antibody</u>	CSB-PA016137LA01HU	Human	ELISA, IF
<u>NTSR2 Antibody</u>	CSB-PA016138LA01HU	Human	ELISA, WB, IHC, IF
<u>NUAK1 Antibody</u>	CSB-PA016139LA01HU	Human	ELISA
<u>NUBP2 Antibody</u>	CSB-PA896918LA01HU	Human	ELISA, IHC, IF
<u>NUBPL Antibody</u>	CSB-PA819450LA01HU	Human	ELISA, WB, IHC
<u>NUCB2 Antibody</u>	CSB-PA016146HA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>Nuclease P1 Antibody</u>	CSB-PA328459LA01PEU	Penicillium citrinum	ELISA
<u>NUCD2 Antibody</u>	CSB-PA840986LA01HU	Human	ELISA
<u>NUCD3 Antibody</u>	CSB-PA811605LA01HU	Human	ELISA, WB
<u>NUDT1 Antibody</u>	CSB-PA016154LA01HU	Human	ELISA, IHC, IF
<u>NUDT11 Antibody</u>	CSB-PA836216LA01HU	Human	ELISA, WB, IHC
<u>NUDT13 Antibody</u>	CSB-PA769800LA01HU	Human	ELISA
<u>NUDT14 Antibody</u>	CSB-PA016159LA01HU	Human	ELISA
<u>NUDT15 Antibody</u>	CSB-PA885733LA01HU	Human	ELISA, IHC, IF
<u>NUDT21 Antibody</u>	CSB-PA016170LA01HU	Human	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>NUDT5 Antibody</u>	CSB-PA891535LA01HU	Human	ELISA, IHC, IF
<u>NUDT9 Antibody</u>	CSB-PA880117LA01HU	Human, Mouse	ELISA, WB
<u>NUGGC Antibody</u>	CSB-PA718072LA01HU	Human	ELISA, WB
<u>NUMA1 Antibody</u>	CSB-PA016185YA01HU	Human	ELISA, WB, IHC, IF
<u>NUMBL Antibody</u>	CSB-PA897595LA01HU	Human	ELISA, IHC, IF
<u>NUP107 Antibody</u>	CSB-PA016188LA01HU	Human, Mouse	ELISA, WB, IHC
<u>NUP153 Antibody</u>	CSB-PA016190LA01HU	Human	ELISA, IHC
<u>NUP210 Antibody</u>	CSB-PA016195LA01HU	Human	ELISA, IHC, IF
<u>NUP35 Antibody</u>	CSB-PA836735LA01HU	Human	ELISA, WB
<u>NUP37 Antibody</u>	CSB-PA847704LA01HU	Human	ELISA
<u>NUP43 Antibody</u>	CSB-PA844045LA01HU	Human	ELISA, WB, IHC
<u>NUP58 Antibody</u>	CSB-PA887135LA01HU	Human, Mouse	ELISA, WB, IHC
<u>NUP62 Antibody</u>	CSB-PA016204LA01HU	Human, Mouse	ELISA, WB, IHC, IP
<u>NUP85 Antibody</u>	CSB-PA874821LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>NUP93 Antibody</u>	CSB-PA822683LA01HU	Human	ELISA, IHC
<u>NUP98 Antibody</u>	CSB-PA016209LA01HU	Human	ELISA, WB, IF, IP
<u>NUPL2 Antibody</u>	CSB-PA016211LA01HU	Human	ELISA, IHC
<u>NUPR1 Antibody</u>	CSB-PA527662LA01HU	Human	ELISA, IHC, IF
<u>NUTM2F Antibody</u>	CSB-PA008184LA01HU	Human	ELISA
<u>NXF1 Antibody</u>	CSB-PA016218LA01HU	Human	ELISA, IF
<u>NXF3 Antibody</u>	CSB-PA016221LA01HU	Human	ELISA, IHC, IF
<u>NXF5 Antibody</u>	CSB-PA016223LA01HU	Human	ELISA, WB, IHC, IF
<u>NXN Antibody</u>	CSB-PA016224HA01HU	Human	ELISA, WB
<u>NXNL2 Antibody</u>	CSB-PA726124LA01HU	Human	ELISA, IHC
<u>NXPE1 Antibody</u>	CSB-PA847620LA01HU	Human	ELISA, WB, IHC
<u>NXPE3 Antibody</u>	CSB-PA846579LA01HU	Human	ELISA, IHC
<u>NXPH1 Antibody</u>	CSB-PA016227LA01HU	Human	ELISA, IHC
<u>NXPH2 Antibody</u>	CSB-PA016228LA01HU	Human	ELISA, IHC
<u>OAS1 Antibody</u>	CSB-PA640913LA01HU	Human	ELISA, IHC
<u>OAS3 Antibody</u>	CSB-PA897591LA01HU	Human	ELISA, IHC, IF
<u>OAZ3 Antibody</u>	CSB-PA016247LA01HU	Human	ELISA
<u>OBP2A Antibody</u>	CSB-PA873688HA01HU	Human	ELISA, IHC
<u>OBSCN Antibody</u>	CSB-PA729422LA01HU	Human	ELISA, IHC, IF
<u>OBSL1 Antibody</u>	CSB-PA016254LA01HU	Human	ELISA, WB, IHC
<u>OCA2 Antibody</u>	CSB-PA22759A0Rb	Human, Rat	ELISA, WB
<u>OCC1 Antibody</u>	CSB-PA851524LA01HU	Human	ELISA, IHC
<u>OCIAD1 Antibody</u>	CSB-PA882129LA01HU	Human, Mouse	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>OCN Antibody</u>	CSB-PA016263LA01HU	Human	ELISA, WB, IHC, IF
<u>Ocm Antibody</u>	CSB-PA016264HA01MO	Mouse	ELISA, WB
<u>OCRL Antibody</u>	CSB-PA016267LA01HU	Human	ELISA, IHC
<u>OCT4 Monoclonal Antibody</u>	CSB-MA018403A0m	Human, Mouse, Rat	ELISA, WB, IHC, IF, FC
<u>ODAPH Antibody</u>	CSB-PA613780LA01HU	Human	ELISA, IHC
<u>ODC1 Antibody</u>	CSB-PA00535A0Rb	Human	ELISA, WB, IHC, IF
<u>ODF1 Antibody</u>	CSB-PA613531LA01HU	Human	ELISA
<u>ODF3L1 Antibody</u>	CSB-PA814225LA01HU	Human	ELISA, WB, IHC
<u>ODR4 Antibody</u>	CSB-PA725554LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>OFCC1 Antibody</u>	CSB-PA809031LA01HU	Human, Mouse	ELISA, WB
<u>OFD1 Antibody</u>	CSB-PA016289LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>OGDH Antibody</u>	CSB-PA016307LA01HU	Human	ELISA, IHC
<u>OGFOD2 Antibody</u>	CSB-PA754021LA01HU	Human, Mouse	ELISA, WB
<u>OGG1 Antibody</u>	CSB-PA016313LA01HU	Human	ELISA, WB, IHC
<u>OGN Antibody</u>	CSB-PA016314LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>OGT Antibody</u>	CSB-PA016315LA01HU	Human	ELISA, IHC, IF
<u>OIH Antibody</u>	CSB-PA318020LA01CH	Gallus gallus	ELISA
<u>OIP5 Antibody</u>	CSB-PA016316LA01HU	Human	ELISA
<u>OLA1 Antibody</u>	CSB-PA016318HA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>OLE9 Antibody</u>	CSB-PA19957A0Rb	Olea europaea	ELISA, WB
<u>OLFM1 Antibody</u>	CSB-PA858725LA01HU	Human, Mouse	ELISA, WB, IHC
<u>OLFM4 Antibody</u>	CSB-PA764941LA01HU	Human, Rat	ELISA, WB, IHC
<u>OLFML1 Antibody</u>	CSB-PA740917LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>OLFML2A Antibody</u>	CSB-PA731560LA01HU	Human	ELISA, IHC
<u>OLFML2B Antibody</u>	CSB-PA734650LA01HU	Human, Mouse	ELISA, WB, IHC
<u>OLFML3 Antibody</u>	CSB-PA885708LA01HU	Human, Rat	ELISA, WB, IHC
<u>OLIG1 Antibody</u>	CSB-PA016328LA01HU	Human, Rat	ELISA, WB
<u>OLIG1 Antibody</u>	CSB-PA016328YA01HU	Human, Mouse, Rat	ELISA, WB
<u>olig2 Antibody</u>	CSB-PA24549A0Rb	Zebrafish	ELISA, WB
<u>OLR1 Antibody</u>	CSB-PA016331LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>OMA1 Antibody</u>	CSB-PA856923LA01HU	Human	ELISA, WB, IHC, IF
<u>omcB Antibody</u>	CSB-PA315528HA01DSB	Chlamydia trachomatis	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>OMD Antibody</u>	CSB-PA016333LA01HU	Human	ELISA, WB, IHC
<u>Omega-gliadin Antibody</u>	CSB-PA360898DA01TKR	Triticum monococcum	ELISA
<u>OMP Antibody</u>	CSB-PA016335LA01HU	Human	ELISA, IF
<u>omp Antibody</u>	CSB-PA706583LA01RMU	Rickettsia japonica	ELISA, WB
<u>ompA Antibody</u>	CSB-PA028238LA01ZBR	Chlamydia pecorum	ELISA
<u>ompA Antibody</u>	CSB-PA16479A0Rb	Salmonella typhi	ELISA
<u>ompA Antibody</u>	CSB-PA328425LA01DSA	Chlamydia trachomatis	ELISA
<u>ompA Antibody</u>	CSB-PA333258LA01DRZ	Chlamydophila pneumoniae	ELISA
<u>ompA Antibody</u>	CSB-PA706579LA01RMS	Rickettsia conorii	ELISA, WB
<u>ompA Antibody</u>	CSB-PA332793LA01HSZ	Haemophilus influenzae	ELISA, WB
<u>ompA Antibody</u>	CSB-PA356881LA01DSG	Chlamydia trachomatis serovar L2	ELISA, WB
<u>ompC Antibody</u>	CSB-PA356994LA01ENV	Escherichia coli	ELISA
<u>ompC Antibody</u>	CSB-PA848941LA01EOD	Escherichia coli O157:H7	ELISA, WB
<u>ompK Antibody</u>	CSB-PA349512LA01VFE	Vibrio parahaemolyticus serotype O3:K6	ELISA, WB
<u>ompT Antibody</u>	CSB-PA348292LA01EOD	Escherichia coli O157:H7	ELISA
<u>ompT Antibody</u>	CSB-PA357589LA01EVN	Escherichia coli	ELISA, WB
<u>ONECUT1 Antibody</u>	CSB-PA880058LA01HU	Human, Rat	ELISA, WB, IHC
<u>ONECUT2 Antibody</u>	CSB-PA016337LA01HU	Human	ELISA, WB, IHC
<u>OOSP2 Antibody</u>	CSB-PA769798LA01HU	Human	ELISA, IHC
<u>OPA1 Antibody</u>	CSB-PA016340LA01HU	Human, Mouse	ELISA, WB
<u>OPA1 Antibody</u>	CSB-PA016340HA01HU	Human	ELISA, IHC
<u>OPALIN Antibody</u>	CSB-PA822279LA01HU	Human	ELISA, IF
<u>OPCML Antibody</u>	CSB-PA621799LA01HU	Human	ELISA, WB
<u>OPN1MW Antibody</u>	CSB-PA016352LA01HU	Human	ELISA, WB, IF
<u>OPN1SW Antibody</u>	CSB-PA016354LA01HU	Human	ELISA, WB, IF
<u>OPN3 Antibody</u>	CSB-PA872435LA01HU	Human	ELISA, IHC, IF
<u>OPN5 Antibody</u>	CSB-PA740895LA01HU	Human	ELISA, IHC, IF
<u>oprF Antibody</u>	CSB-PA318417LA01EZS	Pseudomonas aeruginosa	ELISA, WB
<u>OPRK1 Antibody</u>	CSB-PA016359OA01HU	Human	ELISA, IHC, IF
<u>OPRL1 Antibody</u>	CSB-PA016360LA01HU	Human, Mouse	ELISA, WB, IF
<u>OPRM1 Antibody</u>	CSB-PA016361LA01HU	Human	ELISA, WB, IHC, IF
<u>OPTC Antibody</u>	CSB-PA016362LA01HU	Human	ELISA, IHC
<u>OR10W1 Antibody</u>	CSB-PA818793LA01HU	Human	ELISA
<u>OR1D2 Antibody</u>	CSB-PA016490LA01HU	Human	ELISA
<u>OR2H1 Antibody</u>	CSB-PA875629LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>OR2L13 Antibody</u>	CSB-PA822694LA01HU	Human, Mouse	ELISA, WB, IF
<u>OR3A1 Antibody</u>	CSB-PA016636LA01HU	Human	ELISA



Product Name	Code	Species Reactivity	Applications
<u>OR5L1 Antibody</u>	CSB-PA839894LA01HU	Human	ELISA, IHC
<u>ORAI1 Antibody</u>	CSB-PA846601LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ORC4 Antibody</u>	CSB-PA017234LA01HU	Human	ELISA, IHC, IF
<u>ORF1 Antibody</u>	CSB-PA327230LA01HGE	Hepatitis E virus genotype 1	ELISA
<u>ORF2 Antibody</u>	CSB-PA727851LA01HVZ	Hepatitis E virus genotype 1	ELISA
<u>ORF3 Antibody</u>	CSB-PA527126LA01HVZ	Hepatitis E virus genotype 1	ELISA
<u>ORF7 Antibody</u>	CSB-PA21639A0Rb	Porcine reproductive and respiratory syndrome virus	ELISA, WB
<u>ORM1 Antibody</u>	CSB-PA05709A0Rb	Human	ELISA, IHC
<u>ORM1 Antibody</u>	CSB-PA017237LA01HU	Human	ELISA, IHC
<u>ORM2 Antibody</u>	CSB-PA14549A0Rb	Human, Mouse	ELISA, WB
<u>ORMDL1 Antibody</u>	CSB-PA868379LA01HU	Human	ELISA, IHC, IF
<u>ORMDL2 Antibody</u>	CSB-PA690476LA01HU	Human	ELISA, IHC, IF
<u>ORMDL3 Antibody</u>	CSB-PA854016LA01HU	Human	ELISA, IHC, IF
<u>OS9 Antibody</u>	CSB-PA613402LA01HU	Human	ELISA, WB, IHC, IF
<u>OSBP Antibody</u>	CSB-PA017254LA01HU	Human	ELISA, IF
<u>OSBPL2 Antibody</u>	CSB-PA884442LA01HU	Human, Mouse	ELISA, WB, IHC
<u>OSBPL6 Antibody</u>	CSB-PA883645LA01HU	Human	ELISA, IHC
<u>OSBPL7 Antibody</u>	CSB-PA861186LA01HU	Human	ELISA, IHC, IF
<u>OSER1 Antibody</u>	CSB-PA865164LA01HU	Human, Mouse	ELISA, WB, IHC
<u>OSGIN2 Antibody</u>	CSB-PA897460LA01HU	Human	ELISA, WB, IHC
<u>OSM Antibody</u>	CSB-PA09694A0Rb	Human	ELISA, IHC
<u>Osm Antibody</u>	CSB-PA723970LA01RA	Rat	ELISA, IF
<u>OSM34 Antibody</u>	CSB-PA342308LA01DOA	Arabidopsis thaliana	ELISA, WB
<u>OSMR Antibody</u>	CSB-PA857869LA01HU	Human	ELISA, IHC, IF
<u>OSR2 Antibody</u>	CSB-PA822691LA01HU	Human	ELISA, WB, IHC
<u>OTC Antibody</u>	CSB-PA017270LA01HU	Human	ELISA, IHC
<u>OTOG Antibody</u>	CSB-PA017273LA01HU	Human	ELISA, IHC, IF
<u>OTOP1 Antibody</u>	CSB-PA752530LA01HU	Human	ELISA, IHC, IF
<u>OTUB1 Antibody</u>	CSB-PA03595A0Rb	Human	ELISA, WB, IHC
<u>OTUD3 Antibody</u>	CSB-PA719399LA01HU	Human	ELISA, WB, IF
<u>OTUD5 Antibody</u>	CSB-PA850298LA01HU	Human	ELISA, WB, IHC, IF
<u>OTUD6B Antibody</u>	CSB-PA847646LA01HU	Human	ELISA, WB, IHC
<u>OTUD7B Antibody</u>	CSB-PA753588LA01HU	Human	ELISA, IHC, IF
<u>OTULIN Antibody</u>	CSB-PA853391LA01HU	Human	ELISA, IF
<u>OTX1 Antibody</u>	CSB-PA017298LA01HU	Human	ELISA, IHC, IF
<u>OTX2 Antibody</u>	CSB-PA017299LA01HU	Human, Mouse	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>Outer capsid glycoprotein VP7 Antibody</u>	CSB-PA318053LA01ROH	Rotavirus A	ELISA
<u>Outer capsid protein VP4 Antibody</u>	CSB-PA320838LA01ROH	Rotavirus A	ELISA
<u>OV16 Antibody</u>	CSB-PA341251HA01OEG	Onchocerca volvulus	ELISA
<u>OVGP1 Antibody</u>	CSB-PA614260LA01HU	Human, Mouse	ELISA, WB
<u>OVOL1 Antibody</u>	CSB-PA017303LA01HU	Human	ELISA, WB, IHC
<u>OVOL2 Antibody</u>	CSB-PA887110HA01HU	Human	ELISA, IHC
<u>Ovomucoid Antibody</u>	CSB-PA00400E0Rb	Chicken	ELISA
<u>Ovotransferrin Antibody</u>	CSB-PA00390E0Rb	Chicken	ELISA
<u>OXCT1 Antibody</u>	CSB-PA017306LA01HU	Human	ELISA, WB, IHC
<u>OXER1 Antibody</u>	CSB-PA837442LA01HU	Human, Rat	ELISA, WB, IHC
<u>OXNAD1 Antibody</u>	CSB-PA839337LA01HU	Human, Mouse	ELISA, WB
<u>OXR1 Antibody</u>	CSB-PA854043LA01HU	Human	ELISA, WB, IHC
<u>OXSM Antibody</u>	CSB-PA868340LA01HU	Human	ELISA, IHC
<u>OXT Antibody</u>	CSB-PA017315LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>OXTR Antibody</u>	CSB-PA017316LA01HU	Human	ELISA, IF
<u>P Antibody</u>	CSB-PA303154LA01RAH	Rabies virus	ELISA, WB
<u>P Antibody</u>	CSB-PA324598LA01HDO	Hepatitis B virus genotype A2 subtype adw2	ELISA
<u>P Antibody</u>	CSB-PA768100LA01VBH	Vesicular stomatitis Indiana virus	ELISA, WB
<u>P2RX1 Antibody</u>	CSB-PA017319LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>P2RX4 Antibody</u>	CSB-PA859513LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>P2RX7 Antibody</u>	CSB-PA859103LA01HU	Human	ELISA, IF
<u>P2RY1 Antibody</u>	CSB-PA017326LA01HU	Human	ELISA, IHC
<u>P2RY13 Antibody</u>	CSB-PA860651LA01HU	Human	ELISA, WB, IF
<u>p37 Antibody</u>	CSB-PA323355HA01MLR	Mycoplasma hyorhinis	ELISA
<u>P3H2 Antibody</u>	CSB-PA815555HA01HU	Human, Mouse	ELISA, WB, IHC
<u>p46 Antibody</u>	CSB-PA313811HA01MWK	Mycoplasma hyopneumoniae	ELISA
<u>p46 Antibody</u>	CSB-PA314621LA01MLQ	Mycoplasma hyopneumoniae	ELISA
<u>p46 Antibody</u>	CSB-PA314621YA01MLQ	Mycoplasma hyopneumoniae	ELISA
<u>P4HA1 Antibody</u>	CSB-PA017339LA01HU	Human, Rat	ELISA, WB, IHC
<u>P4HA3 Antibody</u>	CSB-PA772021LA01HU	Human	ELISA, IHC
<u>P4HB Antibody</u>	CSB-PA00254A0Rb	Human, Mouse	ELISA, WB, IHC
<u>P4hb Antibody</u>	CSB-PA017342YA01MO	Mouse	ELISA
<u>P4HB Antibody</u>	CSB-PA823188YA01DXU	Cricetulus griseus	ELISA



Product Name	Code	Species Reactivity	Applications
<u>P4HTM Antibody</u>	CSB-PA878913LA01HU	Human	ELISA
<u>PA Antibody</u>	CSB-PA395880EA01ILR	Influenza A virus	ELISA
<u>PA2G4 Antibody</u>	CSB-PA891987LA01HU	Human	ELISA, WB, IHC, IF
<u>PACRGL Antibody</u>	CSB-PA850819LA01HU	Human	ELISA
<u>PACS1 Antibody</u>	CSB-PA017371LA01HU	Human	ELISA, WB
<u>PACS2 Antibody</u>	CSB-PA017372LA01HU	Human	ELISA, IHC
<u>PACSIN2 Antibody</u>	CSB-PA890765EA01HU	Human	ELISA, WB, IHC, IF, IP
<u>PACSIN3 Antibody</u>	CSB-PA892456LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>PACSIN3 Antibody</u>	CSB-PA892456YA01HU	Human	ELISA, WB
<u>PADI1 Antibody</u>	CSB-PA891543LA01HU	Human	ELISA
<u>PADI2 Antibody</u>	CSB-PA896493LA01HU	Human, Mouse	ELISA, WB
<u>PADI2 Antibody</u>	CSB-PA896493YA01HU	Human	ELISA, WB, IHC, IF
<u>PADI3 Antibody</u>	CSB-PA891552LA01HU	Human	ELISA, IHC
<u>PAEP Antibody</u>	CSB-PA017381LA01HU	Human	ELISA, IHC, IP
<u>PAF1 Antibody</u>	CSB-PA843285LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PAFAH1B1 Antibody</u>	CSB-PA017383LA01HU	Human	ELISA, IF
<u>PAFAH1B2 Antibody</u>	CSB-PA017384LA01HU	Human	ELISA, IHC, IF
<u>PAFAH1B3 Antibody</u>	CSB-PA619963LA01HU	Human	ELISA, IHC, IF
<u>PAFAH2 Antibody</u>	CSB-PA857448LA01HU	Human	ELISA, IHC, IF
<u>PAG 1 Antibody</u>	CSB-PA644917LA01BO	Bovine, Human, Rat	ELISA, WB
<u>PAG1 Antibody</u>	CSB-PA885746LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PAG2 Antibody</u>	CSB-PA637747LA01BO	Bovine	ELISA
<u>PAGE1 Antibody</u>	CSB-PA13219A0Rb	Human	ELISA, IHC
<u>PAGE2B Antibody</u>	CSB-PA708481LA01HU	Human	ELISA, IHC
<u>PAH Antibody</u>	CSB-PA017396LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PAIP2 Antibody</u>	CSB-PA859027LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PAK3 Antibody</u>	CSB-PA017407LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PAK4 Antibody</u>	CSB-PA05437A0Rb	Human, Rat	ELISA, WB, IHC, IF
<u>pal Antibody</u>	CSB-PA340792HA01LNY	Legionella pneumophila	ELISA
<u>PALM Antibody</u>	CSB-PA017413HA01HU	Human, Mouse	ELISA, WB, IHC
<u>PAM Antibody</u>	CSB-PA017417LA01HU	Human	ELISA, WB, IHC, IF
<u>PAM16 Antibody</u>	CSB-PA896711DA01HU	Human	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>PAN2 Antibody</u>	CSB-PA681391LA01HU	Human	ELISA, IHC, IF
<u>PANK1 Antibody</u>	CSB-PA017421LA01HU	Human	ELISA, IHC, IF
<u>PANK2 Antibody</u>	CSB-PA874850LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PANX1 Antibody</u>	CSB-PA857017LA01HU	Human	ELISA, WB, IHC, IF
<u>PANX2 Antibody</u>	CSB-PA839405LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>PANX3 Antibody</u>	CSB-PA839402LA01HU	Human	ELISA, WB, IF
<u>Paox Antibody</u>	CSB-PA804896LA01MO	Mouse	ELISA
<u>PAP1 Antibody</u>	CSB-PA529312LA01DOA	Arabidopsis thaliana	ELISA, WB
<u>Papain Antibody</u>	CSB-PA365471LA01DQH	Carica papaya	ELISA, WB
<u>PAPLN Antibody</u>	CSB-PA017431LA01HU	Human	ELISA
<u>PAPPA Antibody</u>	CSB-PA06919A0Rb	Human	ELISA, IHC
<u>PAPSS2 Antibody</u>	CSB-PA527262LA01HU	Human	ELISA, IF
<u>PAQR5 Antibody</u>	CSB-PA889141LA01HU	Human	ELISA, WB, IHC, IF
<u>PAQR7 Antibody</u>	CSB-PA768224LA01HU	Human	ELISA, IF
<u>PARD3 Antibody</u>	CSB-PA819481LA01HU	Human	ELISA, WB, IHC, IF
<u>PARD6A Antibody</u>	CSB-PA878839LA01HU	Human, Mouse	ELISA, WB, IF
<u>PARD6G Antibody</u>	CSB-PA863155LA01HU	Human, Mouse	ELISA, WB, IF
<u>PARG Antibody</u>	CSB-PA773045LA01HU	Human	ELISA, IF
<u>PARK7 Antibody</u>	CSB-PA03724A0Rb	Human, Mouse, Rat	ELISA, WB, IHC
<u>PARK7 Antibody</u>	CSB-PA03729A0Rb	Human	ELISA, WB, IHC, IF
<u>PARN Antibody</u>	CSB-PA017456LA01HU	Human	ELISA, WB, IHC, IF
<u>PARP14 Antibody</u>	CSB-PA677190LA01HU	Human	ELISA, IF
<u>PARP2 Antibody</u>	CSB-PA866224LA01HU	Human	ELISA, WB
<u>PARP3 Antibody</u>	CSB-PA017467EA01HU	Human	ELISA, WB, IHC
<u>PARP8 Antibody</u>	CSB-PA847623LA01HU	Human	ELISA, IHC, IF
<u>PARS2 Antibody</u>	CSB-PA745334LA01HU	Human	ELISA, IHC
<u>PARVA Antibody</u>	CSB-PA868328LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>Parvalbumin beta Antibody</u>	CSB-PA365766HA01GAA	Gadus morhua subsp	ELISA
<u>PASD1 Antibody</u>	CSB-PA816882LA01HU	Human	ELISA, WB
<u>PASK Antibody</u>	CSB-PA822296LA01HU	Human	ELISA, IF
<u>pat Antibody</u>	CSB-PA701139HA01FOQ	Streptomyces viridochromogenes	ELISA
<u>pat Antibody</u>	CSB-PA701139YA01FOQ	Streptomyces viridochromogenes	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>PAX2 Antibody</u>	CSB-PA11564A0Rb	Human	ELISA
<u>PAX3 Antibody</u>	CSB-PA017489HA01HU	Human	ELISA, WB, IHC, IF
<u>PAX5 Antibody</u>	CSB-PA20659A0Rb	Human	ELISA, WB, IF
<u>PAX8 Antibody</u>	CSB-PA017494LA01HU	Human	ELISA, WB
<u>PAX9 Antibody</u>	CSB-PA017495LA01HU	Human	ELISA, IHC
<u>PAXIP1 Antibody</u>	CSB-PA761652LA01HU	Human	ELISA, IHC
<u>PBLD Antibody</u>	CSB-PA017501LA01HU	Human, Mouse	ELISA, WB, IHC
<u>pbpA Antibody</u>	CSB-PA401992HA01CWW	Clostridium botulinum	ELISA
<u>pbpD Antibody</u>	CSB-PA334764LA01BRJ	Bacillus subtilis	ELISA
<u>pbpE Antibody</u>	CSB-PA333755HA01BRJ	Bacillus subtilis	ELISA, WB
<u>PBX1 Antibody</u>	CSB-PA017505LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>PBXIP1 Antibody</u>	CSB-PA856896HA01HU	Human	ELISA, WB
<u>PC Antibody</u>	CSB-PA017511LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PCBD1 Antibody</u>	CSB-PA017514LA01HU	Human	ELISA, IHC, IF
<u>PCBP1 Antibody</u>	CSB-PA613590LA01HU	Human	ELISA
<u>PCBP2 Antibody</u>	CSB-PA622995LA01HU	Human	ELISA, WB, IHC
<u>PCBP4 Antibody</u>	CSB-PA017521LA01HU	Human	ELISA
<u>PCCA Antibody</u>	CSB-PA017522LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PCCB Antibody</u>	CSB-PA017523LA01HU	Human	ELISA, WB, IHC
<u>PCDH11X Antibody</u>	CSB-PA866326LA01HU	Human	ELISA, IHC, IF
<u>PCDH20 Antibody</u>	CSB-PA822719LA01HU	Human	ELISA, IF
<u>PCDH7 Antibody</u>	CSB-PA017536LA01HU	Human	ELISA, IHC
<u>PCDHB16 Antibody</u>	CSB-PA889094LA01HU	Human	ELISA, IF
<u>PCDHGC3 Antibody</u>	CSB-PA892155LA01HU	Human	ELISA, IF
<u>PCED1A Antibody</u>	CSB-PA867128LA01HU	Human	ELISA, IHC
<u>PCF11 Antibody</u>	CSB-PA017603LA01HU	Human	ELISA, IHC
<u>PCGF1 Antibody</u>	CSB-PA863114LA01HU	Human	ELISA, IHC, IF
<u>PCGF2 Antibody</u>	CSB-PA017606LA01HU	Human	ELISA, IF
<u>PCGF6 Antibody</u>	CSB-PA871612LA01HU	Human	ELISA, WB
<u>PCID2 Antibody</u>	CSB-PA705514LA01HU	Human	ELISA, IHC, IF
<u>PCIF1 Antibody</u>	CSB-PA884471LA01HU	Human	ELISA, IHC
<u>PCK1 Antibody</u>	CSB-PA017613LA01HU	Human	ELISA, IHC, IF
<u>PCK2 Antibody</u>	CSB-PA619092LA01HU	Human, Mouse	ELISA, WB
<u>PCLAF Antibody</u>	CSB-PA012164DA01HU	Human	ELISA
<u>PCMP-H63 Antibody</u>	CSB-PA890450YA01DOA	Arabidopsis thaliana	ELISA
<u>PCMT1 Antibody</u>	CSB-PA017618LA01HU	Human	ELISA, IHC, IF
<u>PCNA Antibody</u>	CSB-PA01567A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>PCNA Antibody</u>	CSB-PA878702LA01DOA	Arabidopsis thaliana	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>PCNA Monoclonal Antibody</u>	CSB-MA000081M0m		ELISA
<u>PCOLCE Antibody</u>	CSB-PA614881LA01HU	Human	ELISA, IHC
<u>PCOLCE2 Antibody</u>	CSB-PA017632LA01HU	Human	ELISA, WB, IHC, IF
<u>PCOTH Antibody</u>	CSB-PA685351LA01HU	Human	ELISA, IHC, IF
<u>PCP4 Antibody</u>	CSB-PA017636HA01HU	Human	ELISA
<u>Pcsk1 Antibody</u>	CSB-PA017640LA01RA	Rat	ELISA
<u>Pcsk2 Antibody</u>	CSB-PA017642LA01RA	Rat, Human	ELISA, WB
<u>PCSK2 Antibody</u>	CSB-PA017642LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PCSK5 Antibody</u>	CSB-PA017644LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>Pcsk5 Antibody</u>	CSB-PA017644LA01MO	Mouse	ELISA
<u>PCSK6 Antibody</u>	CSB-PA017645LA01HU	Human	ELISA, IHC
<u>PCSK7 Antibody</u>	CSB-PA017646LA01HU	Human	ELISA, IHC
<u>PCSK9 Antibody</u>	CSB-PA07039A0Rb	Human	ELISA, IHC, IF
<u>PCSK9 Antibody</u>	CSB-PA017647LA01HU	Human	ELISA, WB, IHC, IF
<u>PCTP Antibody</u>	CSB-PA017651LA01HU	Human	ELISA, WB, IHC, IF
<u>PDCD1 Antibody</u>	CSB-PA619964LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PDCD10 Antibody</u>	CSB-PA861141HA01HU	Human	ELISA, IHC, IF
<u>PDCD1LG2 Antibody</u>	CSB-PA017667EA01HU	Human	ELISA
<u>PDCD6 Antibody</u>	CSB-PA04874A0Rb	Human	ELISA, IHC, IF
<u>PDCL2 Antibody</u>	CSB-PA836663LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>PDCL3 Antibody</u>	CSB-PA872479LA01HU	Human	ELISA, WB, IF
<u>PDE11A Antibody</u>	CSB-PA875726LA01HU	Human	ELISA, IHC, IF
<u>PDE1A Antibody</u>	CSB-PA22739A0Rb	Human	ELISA, IF
<u>PDE1B Antibody</u>	CSB-PA27759A0Rb	Human	ELISA, WB, IF
<u>PDE1C Antibody</u>	CSB-PA623901LA01HU	Human	ELISA, IF
<u>PDE2A Antibody</u>	CSB-PA017687LA01HU	Human	ELISA, WB, IHC, IF
<u>PDE4B Antibody</u>	CSB-PA21029A0Rb	Human	ELISA, IHC
<u>PDE6C Antibody</u>	CSB-PA23139A0Rb	Human	ELISA, WB, IF
<u>PDE6D Antibody</u>	CSB-PA526126LA01HU	Human	ELISA, WB, IHC
<u>PDE6H Antibody</u>	CSB-PA617840LA01HU	Human	ELISA, IHC
<u>PDE8A Antibody</u>	CSB-PA525315LA01HU	Human	ELISA, IHC, IF
<u>PDE8B Antibody</u>	CSB-PA529743LA01HU	Human	ELISA, IF
<u>PDE9A Antibody</u>	CSB-PA527985LA01HU	Human	ELISA, WB, IHC, IF
<u>pdeH Antibody</u>	CSB-PA336327LA01ENV	Escherichia coli	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>PDF Antibody</u>	CSB-PA872538LA01HU	Human	ELISA, WB, IHC
<u>Pdgfd Antibody</u>	CSB-PA852889LA01MO	Mouse, Human	ELISA, WB, IF
<u>PDGFRA Antibody</u>	CSB-PA017712LA01HU	Human, Mouse	ELISA, WB, IF
<u>PDGFRB Antibody</u>	CSB-PA017713LA01HU	Human	ELISA, IHC, IF
<u>PDHA1 Antibody</u>	CSB-PA017715LA01HU	Human	ELISA, IHC, IF
<u>Pdia2 Antibody</u>	CSB-PA017719HA01MO	Mouse, Human	ELISA, WB
<u>PDIA3 Antibody</u>	CSB-PA06284A0Rb	Human	ELISA
<u>PDIA4 Antibody</u>	CSB-PA017722LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PDIA5 Antibody</u>	CSB-PA614533LA01HU	Human	ELISA, WB, IHC
<u>PDIK1L Antibody</u>	CSB-PA854018LA01HU	Human	ELISA, IHC
<u>PDILT Antibody</u>	CSB-PA818726LA01HU	Human	ELISA, WB
<u>PDK1 Antibody</u>	CSB-PA613584LA01HU	Human, Mouse	ELISA, WB, IHC
<u>Pdk3 Antibody</u>	CSB-PA852852LA01MO	Mouse, Human, Rat	ELISA, WB
<u>PDK4 Antibody</u>	CSB-PA614983LA01HU	Human	ELISA, IHC
<u>PD-L1 Monoclonal Antibody</u>	CSB-MA878942A0m	Human	ELISA, WB, IHC, IF, FC
<u>PD-L1 Monoclonal Antibody</u>	CSB-MA878942A1m	Human	ELISA, WB, IHC, IF, FC
<u>PDLIM2 Antibody</u>	CSB-PA856962LA01HU	Human	ELISA, IHC
<u>PDLIM4 Antibody</u>	CSB-PA017735LA01HU	Human	ELISA, WB, IP
<u>PDLIM5 Antibody</u>	CSB-PA822223LA01HU	Human	ELISA, IHC, IF
<u>PDP1 Antibody</u>	CSB-PA878946LA01HU	Human	ELISA, WB, IHC
<u>PDP2 Antibody</u>	CSB-PA868395LA01HU	Human	ELISA
<u>PDRG1 Antibody</u>	CSB-PA878890HA01HU	Human	ELISA, IHC, IF
<u>PDS5A Antibody</u>	CSB-PA643713LA01HU	Human	ELISA, IHC, IF
<u>PDS5B Antibody</u>	CSB-PA889108LA01HU	Human	ELISA
<u>PDSS2 Antibody</u>	CSB-PA773060LA01HU	Human	ELISA, IHC
<u>PDXDC1 Antibody</u>	CSB-PA761311LA01HU	Human	ELISA, IF
<u>Pdxdp Antibody</u>	CSB-PA017749LA01MO	Mouse, Rat	ELISA, WB
<u>PDZD11 Antibody</u>	CSB-PA707840HA01HU	Human	ELISA, IHC
<u>PDZD3 Antibody</u>	CSB-PA774807LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PDZK1 Antibody</u>	CSB-PA722565HA01HU	Human	ELISA
<u>PDZRN3 Antibody</u>	CSB-PA891797LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PEAK1 Antibody</u>	CSB-PA875688LA01HU	Human	ELISA, IHC
<u>PEAR1 Antibody</u>	CSB-PA017765LA01HU	Human, Rat, Mouse	ELISA, WB, IHC, IF
<u>PEBP1 Antibody</u>	CSB-PA017766HA01HU	Human	ELISA, IHC, IF
<u>PEBP4 Antibody</u>	CSB-PA857025LA01HU	Human	ELISA, IHC
<u>PECAM1 Antibody</u>	CSB-PA017767LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>PECR Antibody</u>	CSB-PA017769LA01HU	Human	ELISA



Product Name	Code	Species Reactivity	Applications
<u>Pectate lyase 1 Antibody</u>	CSB-PA305864HA01JAN	<i>Juniperus ashei</i>	ELISA
<u>Pectate lyase 1 Antibody</u>	CSB-PA856862LA01CHAX	<i>Chamaecyparis obtusa</i>	ELISA
<u>Pectate lyase 1 Antibody</u>	CSB-PA871029LA01COAD	<i>Hesperocyparis arizonica</i>	ELISA, WB
<u>Pectate lyase 5 Antibody</u>	CSB-PA329754HA01BYC	<i>Ambrosia artemisiifolia</i>	ELISA
<u>PEG10 Antibody</u>	CSB-PA768203LA01HU	Human	ELISA, IHC
<u>PELO Antibody</u>	CSB-PA861124LA01HU	Human	ELISA, IHC
<u>PEMT Antibody</u>	CSB-PA017780LA01HU	Human	ELISA, WB, IF
<u>PEPD Antibody</u>	CSB-PA017784LA01HU	Human	ELISA, WB, IHC
<u>pepF Antibody</u>	CSB-PA346917LA01MLW	<i>Mycoplasma pneumoniae</i>	ELISA, WB
<u>pepP Antibody</u>	CSB-PA303407LA01MLW	<i>Mycoplasma pneumoniae</i>	ELISA, WB
<u>PER1 Antibody</u>	CSB-PA017786LA01HU	Human	ELISA, IHC
<u>PER2 Antibody</u>	CSB-PA017787LA01HU	Human	ELISA, IHC
<u>PER3 Antibody</u>	CSB-PA017788LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>Perlwapin Antibody</u>	CSB-PA306283LA01HAZ	<i>Haliotis laevigata</i>	ELISA, WB
<u>PES1 Antibody</u>	CSB-PA017791LA01HU	Human, Mouse	ELISA, WB, IF
<u>PET100 Antibody</u>	CSB-PA319771LA01HU	Human	ELISA, IHC, IF
<u>petE Antibody</u>	CSB-PA344627HA01PSD	<i>Prochlorothrix hollandica</i>	ELISA, WB
<u>PETE Antibody</u>	CSB-PA365424HA01FKI	<i>Spinacia oleracea</i>	ELISA, WB
<u>PEX1 Antibody</u>	CSB-PA017793LA01HU	Human	ELISA, IHC, IF
<u>PEX10 Antibody</u>	CSB-PA017794LA01HU	Human, Rat, Mouse	ELISA, WB
<u>PEX11A Antibody</u>	CSB-PA017795LA01HU	Human, Rat	ELISA, WB, IHC
<u>PEX11B Antibody</u>	CSB-PA017796LA01HU	Human	ELISA, WB, IHC, IF
<u>PEX13 Antibody</u>	CSB-PA849799LA01HU	Human, Rat	ELISA, WB, IHC
<u>PEX14 Antibody</u>	CSB-PA867033LA01DOA	<i>Arabidopsis thaliana</i>	ELISA, WB
<u>PEX26 Antibody</u>	CSB-PA767196LA01HU	Human	ELISA, WB, IHC, IP
<u>PEX5 Antibody</u>	CSB-PA017805LA01HU	Human	ELISA, WB, IHC
<u>PF0142 Antibody</u>	CSB-PA823641LA01FHW	<i>Pyrococcus furiosus</i>	ELISA, WB
<u>PF14_0076 Antibody</u>	CSB-PA745322LA01EWP	<i>Plasmodium falciparum</i>	ELISA
<u>PF4 Antibody</u>	CSB-PA017809LA01HU	Human	ELISA, IHC
<u>PFAS Antibody</u>	CSB-PA017811LA01HU	Human	ELISA, IHC
<u>PFDN1 Antibody</u>	CSB-PA01485A0Rb	Human	ELISA, IHC, IF
<u>PFDN2 Antibody</u>	CSB-PA017813LA01HU	Human	ELISA, WB, IHC, IF
<u>PFDN4 Antibody</u>	CSB-PA017814LA01HU	Human	ELISA
<u>PFDN5 Antibody</u>	CSB-PA01414A0Rb	Human	ELISA, WB
<u>PFDN5 Antibody</u>	CSB-PA01415A0Rb	Human	ELISA, IHC, IF
<u>PFKFB2 Antibody</u>	CSB-PA017818LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>PFKFB4 Antibody</u>	CSB-PA613704LA01HU	Human	ELISA, IF



Product Name	Code	Species Reactivity	Applications
<u>PFKL Antibody</u>	CSB-PA017821LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PFKM Antibody</u>	CSB-PA017822LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PFN1 Antibody</u>	CSB-PA04154A0Rb	Human	ELISA, WB
<u>pfo Antibody</u>	CSB-PA314820HA01CMB	Clostridium perfringens	ELISA, WB
<u>PG_1424 Antibody</u>	CSB-PA886515LA01EYA	Porphyromonas gingivalis	ELISA, WB
<u>PGA3 Antibody</u>	CSB-PA11289A0Rb	Human	ELISA
<u>PGAM5 Antibody</u>	CSB-PA839338LA01HU	Human	ELISA, IHC, IF
<u>PGAP2 Antibody</u>	CSB-PA883397LA01HU	Human	ELISA, IHC, IF
<u>PGC Antibody</u>	CSB-PA11299A0Rb	Human, Mouse	ELISA, WB
<u>PGD Antibody</u>	CSB-PA01155A0Rb	Human	ELISA, WB
<u>PGF Antibody</u>	CSB-PA06959A0Rb	Human	ELISA
<u>PGK1 Antibody</u>	CSB-PA00035A0Rb	Human	ELISA, WB
<u>PGK2 Antibody</u>	CSB-PA017859LA01HU	Human	ELISA, WB
<u>PGLYRP1 Antibody</u>	CSB-PA017862LA01HU	Human	ELISA, IHC
<u>PGLYRP2 Antibody</u>	CSB-PA846684LA01HU	Human	ELISA, IHC
<u>PGLYRP3 Antibody</u>	CSB-PA017864LA01HU	Human	ELISA, WB
<u>PGLYRP4 Antibody</u>	CSB-PA017865LA01HU	Human	ELISA, IHC
<u>PGM1 Antibody</u>	CSB-PA017866HA01HU	Human	ELISA, WB, IHC, IF
<u>PGM2L1 Antibody</u>	CSB-PA764776LA01HU	Human	ELISA, WB, IHC
<u>PGP Antibody</u>	CSB-PA017873LA01HU	Human	ELISA, IHC, IF
<u>PGR Antibody</u>	CSB-PA11385A0Rb	Human	ELISA
<u>PGS1 Antibody</u>	CSB-PA653479LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PHACTR1 Antibody</u>	CSB-PA017880LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>PHAX Antibody</u>	CSB-PA863999LA01HU	Human	ELISA, WB, IF
<u>PHB Antibody</u>	CSB-PA017885LA01HU	Human	ELISA, WB, IHC, IP
<u>PHC1 Antibody</u>	CSB-PA017891HA01HU	Human	ELISA, IHC
<u>PHEX Antibody</u>	CSB-PA017896LA01HU	Human	ELISA, IF
<u>PHF10 Antibody</u>	CSB-PA017898LA01HU	Human	ELISA, WB, IHC, IF
<u>PHF13 Antibody</u>	CSB-PA017902LA01HU	Human	ELISA, IHC
<u>PHF21A Antibody</u>	CSB-PA017911LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PHF5A Antibody</u>	CSB-PA017915EA01HU	Human	ELISA, IHC, IF
<u>PHF6 Antibody</u>	CSB-PA017917LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>PHF7 Antibody</u>	CSB-PA017918LA01HU	Human	ELISA, IHC
<u>PHGDH Antibody</u>	CSB-PA03255A0Rb	Human	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>PHGDH Antibody</u>	CSB-PA03257A0Rb	Human	ELISA, WB, IHC
<u>PHIP Antibody</u>	CSB-PA837866LA01HU	Human	ELISA
<u>PHKA1 Antibody</u>	CSB-PA017922LA01HU	Human	ELISA, IHC
<u>PHLDA1 Antibody</u>	CSB-PA845148LA01HU	Human	ELISA, IHC, IF
<u>PHLDA3 Antibody</u>	CSB-PA017932LA01HU	Human	ELISA, IF
<u>PHLDB1 Antibody</u>	CSB-PA773033LA01HU	Human	ELISA, IHC
<u>PHLDB2 Antibody</u>	CSB-PA803110LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PHLPII Antibody</u>	CSB-PA337346HA01EUQ	Phleum pratense	ELISA, WB
<u>PHOSPHO1 Antibody</u>	CSB-PA819461LA01HU	Human	ELISA, IHC
<u>Phospholipase D LamSicTox-alphaIC1 Antibody</u>	CSB-PA495649LA01LQU	Loxosceles amazonica	ELISA
<u>PHOX2A Antibody</u>	CSB-PA017940LA01HU	Human	ELISA, IHC
<u>PHPT1 Antibody</u>	CSB-PA017942DA01HU	Human	ELISA, WB, IHC, IF
<u>phr Antibody</u>	CSB-PA361494LA01FPZ	Synechococcus sp.	ELISA, WB
<u>phyA Antibody</u>	CSB-PA333934LA01DOZ	Aspergillus niger	ELISA
<u>PI3 Antibody</u>	CSB-PA017952HA01HU	Human	ELISA, IHC, IF
<u>PI4KAP1 Antibody</u>	CSB-PA818731LA01HU	Human, Mouse	ELISA, WB
<u>PIAS2 Antibody</u>	CSB-PA13299A0Rb	Human	ELISA
<u>PIAS4 Antibody</u>	CSB-PA822692LA01HU	Human	ELISA, IHC
<u>PICALM Antibody</u>	CSB-PA622660LA01HU	Human	ELISA, IHC, IF
<u>PICK1 Antibody</u>	CSB-PA05445A0Rb	Human, Mouse	ELISA, WB, IHC
<u>PID1 Antibody</u>	CSB-PA800096LA01HU	Human	ELISA, IHC, IF
<u>PIF1 Antibody</u>	CSB-PA863987LA01HU	Human	ELISA, IHC
<u>PIGA Antibody</u>	CSB-PA017965LA01HU	Human	ELISA, IHC, IF
<u>PIGC Antibody</u>	CSB-PA856402LA01HU	Human	ELISA, IF
<u>PIGH Antibody</u>	CSB-PA617916LA01HU	Human, Mouse	ELISA, WB
<u>PIGK Antibody</u>	CSB-PA842143LA01HU	Human	ELISA, WB, IHC
<u>PIGL Antibody</u>	CSB-PA897465LA01HU	Human	ELISA, IHC
<u>PIGN Antibody</u>	CSB-PA017977LA01HU	Human	ELISA, IHC
<u>PIGO Antibody</u>	CSB-PA855068LA01HU	Human	ELISA, IHC, IF
<u>PIGP Antibody</u>	CSB-PA017979LA01HU	Human	ELISA, IHC, IF
<u>PIGQ Antibody</u>	CSB-PA883583LA01HU	Human	ELISA, IF
<u>PIGR Antibody</u>	CSB-PA017981LA01BO	Bovine	ELISA, WB
<u>PIGS Antibody</u>	CSB-PA847222LA01HU	Human	ELISA, IF
<u>PIGT Antibody</u>	CSB-PA836170LA01HU	Human	ELISA, IHC
<u>PIGU Antibody</u>	CSB-PA884465PA01HU	Human	ELISA, IHC
<u>PIGV Antibody</u>	CSB-PA868321LA01HU	Human	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>PIGW Antibody</u>	CSB-PA773605LA01HU	Human, Rat	ELISA, WB
<u>PIGX Antibody</u>	CSB-PA848395LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PIGZ Antibody</u>	CSB-PA773038LA01HU	Human	ELISA, IF
<u>PIH1D1 Antibody</u>	CSB-PA889130LA01HU	Human	ELISA, IHC, IF, IP
<u>PIH1D3 Antibody</u>	CSB-PA878855LA01HU	Human	ELISA, IHC, IF
<u>PIK3AP1 Antibody</u>	CSB-PA747806LA01HU	Human	ELISA, IF
<u>PIK3C3 Antibody</u>	CSB-PA854108LA01HU	Human	ELISA, WB, IHC
<u>PIK3CA Antibody</u>	CSB-PA017998LA01HU	Human	ELISA, IHC, IF
<u>PIK3CB Antibody</u>	CSB-PA017999LA01HU	Human	ELISA, IF
<u>PIK3CG Antibody</u>	CSB-PA13309A0Rb	Human	ELISA, IF
<u>PIK3IP1 Antibody</u>	CSB-PA850293LA01HU	Human, Mouse	ELISA, WB, IF
<u>PIK3R1 Antibody</u>	CSB-PA018003LA01HU	Human	ELISA, IHC
<u>PIK3R2 Antibody</u>	CSB-PA018005LA01HU	Human	ELISA, WB, IHC
<u>PIK3R4 Antibody</u>	CSB-PA858716LA01HU	Human	ELISA, IHC, IF
<u>PIK3R5 Antibody</u>	CSB-PA841225LA01HU	Human	ELISA, IHC, IF
<u>PIK3R6 Antibody</u>	CSB-PA719469LA01HU	Human	ELISA, WB, IHC, IF
<u>PILRA Antibody</u>	CSB-PA892331LA01HU	Human	ELISA, IHC
<u>PILRB Antibody</u>	CSB-PA890735LA01HU	Human	ELISA, IHC, IF
<u>PIM1 Antibody</u>	CSB-PA02959A0Rb	Human	ELISA, IHC, IF
<u>PIM2 Antibody</u>	CSB-PA873725LA01HU	Human	ELISA, IHC, IF
<u>PIM3 Antibody</u>	CSB-PA771444LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PIN1 Antibody</u>	CSB-PA04619A0Rb	Human	ELISA, IHC, IP
<u>PIN4 Antibody</u>	CSB-PA861478LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>PINK1 Antibody</u>	CSB-PA863144LA01HU	Human	ELISA, IHC
<u>PINX1 Antibody</u>	CSB-PA850264LA01HU	Human, Mouse	ELISA, WB
<u>PIP Antibody</u>	CSB-PA018020LA01HU	Human	ELISA, IHC
<u>PIP4K2C Antibody</u>	CSB-PA819455LA01HU	Human	ELISA, IHC, IF
<u>PIP4P2 Antibody</u>	CSB-PA023855LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PIPOX Antibody</u>	CSB-PA868381LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>PIR Antibody</u>	CSB-PA018032LA01HU	Human	ELISA, IHC, IF
<u>PIRT Antibody</u>	CSB-PA314551LA01HU	Human	ELISA, IHC, IF
<u>PISD Antibody</u>	CSB-PA880077LA01HU	Human	ELISA
<u>PITPNA Antibody</u>	CSB-PA018035LA01HU	Human	ELISA, IHC, IF
<u>PITPNC1 Antibody</u>	CSB-PA892330LA01HU	Human	ELISA, IHC
<u>PITPNM1 Antibody</u>	CSB-PA01995A0Rb	Human	ELISA
<u>PITX2 Antibody</u>	CSB-PA860349LA01HU	Human	ELISA, WB, IHC, IF
<u>PITX3 Antibody</u>	CSB-PA010844LA01HU	Human, Mouse	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>PIWIL4 Antibody</u>	CSB-PA767195LA01HU	Human	ELISA, WB, IF
<u>Pka-C1 Antibody</u>	CSB-PA318954LA01DLU	Drosophila melanogaster	ELISA
<u>PKD1 Antibody</u>	CSB-PA23149A0Rb	Human	ELISA, IHC, IF
<u>PKD1L1 Antibody</u>	CSB-PA823467LA01HU	Human	ELISA, IHC, IF
<u>PKD1L2 Antibody</u>	CSB-PA773589LA01HU	Human	ELISA, IHC, IF
<u>PKD1L3 Antibody</u>	CSB-PA767197LA01HU	Human	ELISA, IHC, IF
<u>PKD2 Antibody</u>	CSB-PA618787LA01HU	Human	ELISA, IHC, IF
<u>PKD2L2 Antibody</u>	CSB-PA882158LA01HU	Human	ELISA, IHC, IF
<u>PKHD1 Antibody</u>	CSB-PA357550LA01HU	Human	ELISA, WB, IHC, IF
<u>PKIB Antibody</u>	CSB-PA883651LA01HU	Human	ELISA, IHC, IF
<u>PKIG Antibody</u>	CSB-PA23619A0Rb	Human	ELISA, IF
<u>PKLR Antibody</u>	CSB-PA018071LA01HU	Human	ELISA
<u>Pklr Antibody</u>	CSB-PA018071LA01RA	Rat, Human	ELISA, WB
<u>PKM Antibody</u>	CSB-PA018072DA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PKM Antibody</u>	CSB-PA02949A0Rb	Human	ELISA, IHC, IF
<u>PKMYT1 Antibody</u>	CSB-PA16029A0Rb	Human	ELISA, IHC
<u>PKN2 Antibody</u>	CSB-PA620966LA01HU	Human	ELISA, IHC, IF
<u>PKN3 Antibody</u>	CSB-PA744059LA01HU	Human	ELISA, WB, IHC
<u>PKP2 Antibody</u>	CSB-PA859116LA01HU	Human	ELISA, IHC
<u>PLA1A Antibody</u>	CSB-PA018084LA01HU	Human	ELISA, WB, IHC
<u>PLA2G10 Antibody</u>	CSB-PA018085LA01HU	Human	ELISA, IHC
<u>PLA2G16 Antibody</u>	CSB-PA018089HA01HU	Human	ELISA
<u>PLA2G16 Antibody</u>	CSB-PA018089LA01HU	Human	ELISA, IHC
<u>PLA2G2A Antibody</u>	CSB-PA11209A0Rb	Human	ELISA, IHC
<u>Pla2g2a Antibody</u>	CSB-PA321246LA01RA	Rat	ELISA
<u>PLA2G2D Antibody</u>	CSB-PA891566LA01HU	Human	ELISA, IHC
<u>PLA2G4A Antibody</u>	CSB-PA018097LA01HU	Human	ELISA, IHC, IF
<u>PLA2G4B Antibody</u>	CSB-PA018098LA01HU	Human	ELISA, WB, IHC, IF
<u>PLA2G4E Antibody</u>	CSB-PA666060LA01HU	Human	ELISA, IHC, IF
<u>PLA2G4F Antibody</u>	CSB-PA721207LA01HU	Human	ELISA, IHC
<u>Pla2g5 Antibody</u>	CSB-PA018103LA01MO	Mouse	ELISA
<u>Pla2g5 Antibody</u>	CSB-PA018103LA01RA	Rat	ELISA
<u>PLA2G7 Antibody</u>	CSB-PA08917A0Rb	Human	ELISA, WB, IHC
<u>PLA2G7 Antibody</u>	CSB-PA08919A0Rb	Human, Mouse	ELISA, WB, IHC
<u>PLA2R1 Antibody</u>	CSB-PA623780LA01HU	Human	ELISA, IHC
<u>PLAC1 Antibody</u>	CSB-PA867193LA01HU	Human	ELISA
<u>PLAC8 Antibody</u>	CSB-PA873705LA01HU	Human	ELISA, IHC, IF
<u>plac8.1 Antibody</u>	CSB-PA29731A0Rb	Zebrafish	ELISA
<u>plac8.2 Antibody</u>	CSB-PA29741A0Rb	Zebrafish	ELISA



Product Name	Code	Species Reactivity	Applications
<u>PLAC9 Antibody</u>	CSB-PA689337LA01HU	Human	ELISA, IF
<u>PLAGL1 Antibody</u>	CSB-PA892351LA01HU	Human	ELISA, WB
<u>PLAGL2 Antibody</u>	CSB-PA891980LA01HU	Human	ELISA, WB
<u>PLAT Antibody</u>	CSB-PA11819A0Rb	Human, Mouse	ELISA, WB, IHC
<u>PLAU Antibody</u>	CSB-PA14319A0Rb	Human, Mouse	ELISA, WB, IHC
<u>Plaur Antibody</u>	CSB-PA018122LA01RA	Rat, Human	ELISA, WB
<u>PLAUR Antibody</u>	CSB-PA07505A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>Plbd2 Antibody</u>	CSB-PA018125LA01MO	Mouse	ELISA, WB
<u>Plbd2 Antibody</u>	CSB-PA683927LA01RA	Rat, Mouse	ELISA, WB
<u>PLCB1 Antibody</u>	CSB-PA865102LA01HU	Human	ELISA, IHC, IF
<u>PLCD1 Antibody</u>	CSB-PA018130LA01HU	Human	ELISA, IHC
<u>PLCD4 Antibody</u>	CSB-PA866264LA01HU	Human	ELISA, IHC
<u>PLCG1 Antibody</u>	CSB-PA018134LA01HU	Human	ELISA, IF
<u>PLCG2 Antibody</u>	CSB-PA018135LA01HU	Human	ELISA, IHC
<u>PLCL1 Antibody</u>	CSB-PA624015LA01HU	Human	ELISA, IHC, IF
<u>PLCL2 Antibody</u>	CSB-PA892366LA01HU	Human	ELISA, IHC
<u>PLCXD1 Antibody</u>	CSB-PA889114LA01HU	Human	ELISA, WB, IHC
<u>PLCZ1 Antibody</u>	CSB-PA768240LA01HU	Human, Mouse	ELISA, WB
<u>PLD1 Antibody</u>	CSB-PA614406LA01HU	Human	ELISA, IF
<u>PLD2 Antibody</u>	CSB-PA018145LA01HU	Human	ELISA, IHC
<u>PLD3 Antibody</u>	CSB-PA018146LA01HU	Human	ELISA, IHC, IF
<u>PLD4 Antibody</u>	CSB-PA839298LA01HU	Human	ELISA, WB, IHC, IF
<u>PLD5 Antibody</u>	CSB-PA836686LA01HU	Human	ELISA, IHC
<u>PLD6 Antibody</u>	CSB-PA836649LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PLEC Antibody</u>	CSB-PA619965LA01HU	Human	ELISA, IHC, IF
<u>PLEK Antibody</u>	CSB-PA018152LA01HU	Human	ELISA
<u>PLEKHA7 Antibody</u>	CSB-PA750339LA01HU	Human	ELISA, IF
<u>PLEKHA8 Antibody</u>	CSB-PA856957LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PLEKHB2 Antibody</u>	CSB-PA846599LA01HU	Human	ELISA, IHC
<u>PLEKHG5 Antibody</u>	CSB-PA018174LA01HU	Human	ELISA, IF
<u>PLEKHJ1 Antibody</u>	CSB-PA868334LA01HU	Human	ELISA, IHC
<u>PLEKHO1 Antibody</u>	CSB-PA706633LA01HU	Human	ELISA, IF
<u>PLEKHS1 Antibody</u>	CSB-PA716341LA01HU	Human	ELISA, WB, IHC
<u>Plet1 Antibody</u>	CSB-PA851843HA01MO	Mouse, Human	ELISA, WB, IHC
<u>PLG Antibody</u>	CSB-PA018188LA01HU	Human	ELISA, IHC
<u>PLGRKT Antibody</u>	CSB-PA875713LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PLIN1 Antibody</u>	CSB-PA526447YA01HU	Human	ELISA, IHC
<u>PLIN3 Antibody</u>	CSB-PA013294LA01HU	Human	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>PLIN4 Antibody</u>	CSB-PA822286HA01HU	Human	ELISA, WB, IHC, IF
<u>PLIN5 Antibody</u>	CSB-PA804867LA01HU	Human, Mouse	ELISA, WB, IF
<u>PLK1 Antibody</u>	CSB-PA018193LA01HU	Human	ELISA, IHC, IF
<u>PLLP Antibody</u>	CSB-PA896501LA01HU	Human, Mouse	ELISA, WB, IF
<u>PLN Antibody</u>	CSB-PA018198LA01HU	Human	ELISA, IHC
<u>PLOD2 Antibody</u>	CSB-PA018200LA01HU	Human	ELISA, WB
<u>plp Antibody</u>	CSB-PA018202LA01OEI	Oncorhynchus mykiss	ELISA
<u>PLP1 Antibody</u>	CSB-PA018202LA01HU	Human	ELISA, WB, IHC, IF
<u>PLPBP Antibody</u>	CSB-PA018755LA01HU	Human	ELISA, IHC
<u>PLPP2 Antibody</u>	CSB-PA018416LA01HU	Human	ELISA, WB, IHC, IF
<u>PLPP3 Antibody</u>	CSB-PA018415LA01HU	Human	ELISA, IHC, IF
<u>PLPPR2 Antibody</u>	CSB-PA836219HA01HU	Human	ELISA, IHC, IF
<u>PLS1 Antibody</u>	CSB-PA623928LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>PLS3 Antibody</u>	CSB-PA018206LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PLTP Antibody</u>	CSB-PA018212EA01HU	Human	ELISA, IHC
<u>PLXNA1 Antibody</u>	CSB-PA14309A0Rb	Human	ELISA, IHC
<u>PLXNA4 Antibody</u>	CSB-PA862071LA01HU	Human	ELISA, WB
<u>PLXNB1 Antibody</u>	CSB-PA018222LA01HU	Human	ELISA, IHC, IF
<u>PLXNB2 Antibody</u>	CSB-PA018223LA01HU	Human	ELISA, WB
<u>PLXNC1 Antibody</u>	CSB-PA018225LA01HU	Human	ELISA, IF
<u>PLXND1 Antibody</u>	CSB-PA659725LA01HU	Human	ELISA, IHC, IF
<u>ply Antibody</u>	CSB-PA314690HA01FMW	Streptococcus pneumoniae serotype 4	ELISA
<u>PMAIP1 Antibody</u>	CSB-PA018229LA01HU	Human	ELISA, IHC, IF
<u>PMAIP3 Antibody</u>	CSB-PA344230HA01PI	Pig	ELISA
<u>PMEI Antibody</u>	CSB-PA021324LA01HU	Human	ELISA, WB, IHC
<u>PMM1 Antibody</u>	CSB-PA018237LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PMP2 Antibody</u>	CSB-PA07854A0Rb	Human	ELISA, WB
<u>PMP2 Antibody</u>	CSB-PA315130LA01HO	Horse	ELISA
<u>PMP22 Antibody</u>	CSB-PA018241LA01HU	Human	ELISA, IF
<u>pmpD Antibody</u>	CSB-PA530746LA01DSB	Chlamydia trachomatis	ELISA, WB
<u>PMRT15 Antibody</u>	CSB-PA812747LA01DOA	Arabidopsis thaliana	ELISA, WB
<u>PNCK Antibody</u>	CSB-PA754178LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>PNKD Antibody</u>	CSB-PA843154LA01HU	Human, Mouse	ELISA, WB
<u>PNLIP Antibody</u>	CSB-PA018262LA01HU	Human	ELISA, IHC
<u>PNLIPRP1 Antibody</u>	CSB-PA018263LA01DO	Dog	ELISA



Product Name	Code	Species Reactivity	Applications
<u>PNLIPRP2 Antibody</u>	CSB-PA018264LA01HU	Human	ELISA, WB, IHC
<u>PNLIPRP3 Antibody</u>	CSB-PA018265LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PNMA1 Antibody</u>	CSB-PA018266LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PNMA2 Antibody</u>	CSB-PA018267LA01HU	Human	ELISA, IHC
<u>PNMA5 Antibody</u>	CSB-PA853483HA01HU	Human	ELISA, IHC
<u>PNMT Antibody</u>	CSB-PA02939A0Rb	Human	ELISA, IHC, IF
<u>PNO1 Antibody</u>	CSB-PA885712LA01HU	Human	ELISA
<u>PNP Antibody</u>	CSB-PA355527HA01HU	Human	ELISA, WB
<u>PNPLA1 Antibody</u>	CSB-PA843295LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PNPLA2 Antibody</u>	CSB-PA836180LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PNPLA5 Antibody</u>	CSB-PA773803LA01HU	Human	ELISA, WB, IHC, IF
<u>PNPLA6 Antibody</u>	CSB-PA809011LA01HU	Human	ELISA, IHC, IF
<u>PNPLA8 Antibody</u>	CSB-PA868272LA01HU	Human	ELISA, IHC, IF
<u>PNPO Antibody</u>	CSB-PA885740LA01HU	Human	ELISA, IHC, IF
<u>PNRC2 Antibody</u>	CSB-PA018290LA01HU	Human	ELISA, WB, IHC
<u>pntB Antibody</u>	CSB-PA364617LA01EVN	Escherichia coli	ELISA
<u>POC5 Antibody</u>	CSB-PA843309LA01HU	Human	ELISA, IHC
<u>PODXL Antibody</u>	CSB-PA13447A0Rb	Human	ELISA, WB, IHC, IF
<u>PODXL Monoclonal Antibody</u>	CSB-MA134471A0m	Human, Mouse, Rat	ELISA, WB, IHC, IP, IF, FC
<u>POF1B Antibody</u>	CSB-PA840988LA01HU	Human, Rat	ELISA, WB
<u>POFUT2 Antibody</u>	CSB-PA837737LA01HU	Human	ELISA, WB, IHC
<u>POGK Antibody</u>	CSB-PA018298LA01HU	Human	ELISA, IHC, IF
<u>POGLUT1 Antibody</u>	CSB-PA818754LA01HU	Human	ELISA, WB, IHC
<u>POGZ Antibody</u>	CSB-PA018299LA01HU	Human	ELISA, IHC
<u>polA Antibody</u>	CSB-PA16667A0Rb	Thermus aquaticus	ELISA, WB
<u>POLA2 Antibody</u>	CSB-PA614517LA01HU	Human	ELISA
<u>POLB Antibody</u>	CSB-PA018302LA01HU	Human	ELISA, WB, IHC
<u>POLD1 Antibody</u>	CSB-PA018303LA01HU	Human	ELISA, WB, IHC, IP
<u>POLD2 Antibody</u>	CSB-PA018304LA01HU	Human	ELISA, IHC
<u>POLD3 Antibody</u>	CSB-PA617997LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>POLDIP2 Antibody</u>	CSB-PA896496LA01HU	Human	ELISA, IF
<u>POLDIP3 Antibody</u>	CSB-PA871610LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>POLE Antibody</u>	CSB-PA894382LA01HU	Human	ELISA, IHC, IF
<u>POLE4 Antibody</u>	CSB-PA889087LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>POLG Antibody</u>	CSB-PA018313LA01HU	Human	ELISA, WB
<u>POLG2 Antibody</u>	CSB-PA883399LA01HU	Human	ELISA
<u>POLH Antibody</u>	CSB-PA018321LA01HU	Human	ELISA, WB, IHC, IF
<u>POLI Antibody</u>	CSB-PA892157LA01HU	Human	ELISA, IHC, IF
<u>POLK Antibody</u>	CSB-PA887009LA01HU	Human, Mouse	ELISA, WB
<u>POLL Antibody</u>	CSB-PA871513LA01HU	Human, Rat	ELISA, WB, IF
<u>Pollen allergen Dac g 3 Antibody</u>	CSB-PA309006HA01DAC	Dactylis glomeRata	ELISA
<u>Pollen allergen Phl p 5a Antibody</u>	CSB-PA669216HA01EUQ	Phleum pratense	ELISA
<u>Pollen allergen Phl p 5b Antibody</u>	CSB-PA671435HA01EUQ	Phleum pRatense	ELISA, WB
<u>POLR1C Antibody</u>	CSB-PA018324LA01HU	Human, Mouse, Rat	ELISA, WB
<u>POLR2C Antibody</u>	CSB-PA018329LA01HU	Human, Mouse	ELISA, WB, IHC
<u>POLR2D Antibody</u>	CSB-PA018331LA01HU	Human	ELISA
<u>POLR2E Antibody</u>	CSB-PA018332LA01HU	Human	ELISA
<u>POLR2G Antibody</u>	CSB-PA018334LA01HU	Human	ELISA, WB
<u>POLR2I Antibody</u>	CSB-PA018336LA01HU	Human	ELISA, IHC
<u>POLR2J Antibody</u>	CSB-PA347112LA01HU	Human	ELISA
<u>POLR2K Antibody</u>	CSB-PA018340LA01HU	Human	ELISA
<u>POLR2M Antibody</u>	CSB-PA316685LA01HU	Human	ELISA, IF
<u>POLR3A Antibody</u>	CSB-PA018343LA01HU	Human	ELISA, IHC, ChIP
<u>POLR3B Antibody</u>	CSB-PA885742LA01HU	Human	ELISA, IHC
<u>POLR3C Antibody</u>	CSB-PA871578LA01HU	Human, Mouse	ELISA, WB
<u>POLR3E Antibody</u>	CSB-PA868333LA01HU	Human	ELISA, IHC
<u>POLR3F Antibody</u>	CSB-PA884440LA01HU	Human	ELISA, WB
<u>POLR3K Antibody</u>	CSB-PA018352HA01HU	Human	ELISA
<u>POLRMT Antibody</u>	CSB-PA018353LA01HU	Human	ELISA, IHC
<u>Polygalacturonase Antibody</u>	CSB-PA337345LA01DYO	Cryptomeria japonica	ELISA
<u>Polygalacturonase Antibody</u>	CSB-PA752116LA01CHAX	Chamaecyparis obtusa	ELISA
<u>POMC Antibody</u>	CSB-PA018363NJ01HU	Human	ELISA, IHC
<u>POMC Antibody</u>	CSB-PA018363OJ01HU	Human, Mouse	ELISA, WB, IHC
<u>POMGNT1 Antibody</u>	CSB-PA855513LA01HU	Human	ELISA, IHC
<u>POMK Antibody</u>	CSB-PA863985LA01HU	Human	ELISA, IHC
<u>POMP Antibody</u>	CSB-PA018365LA01HU	Human	ELISA, IHC, IF
<u>PON1 Antibody</u>	CSB-PA018369EA01HU	Human, Mouse	ELISA, WB, IHC
<u>PON1 Antibody</u>	CSB-PA018369LA01HU	Human	ELISA, IHC
<u>PON2 Antibody</u>	CSB-PA613586LA01HU	Human	ELISA, IF
<u>PON3 Antibody</u>	CSB-PA614883LA01HU	Human	ELISA, IHC
<u>ponA Antibody</u>	CSB-PA345793LA01DKK	Slime mold	ELISA, WB
<u>POP4 Antibody</u>	CSB-PA018373LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>POPDC3 Antibody</u>	CSB-PA881020LA01HU	Human	ELISA, IF
<u>POR Antibody</u>	CSB-PA06275A0Rb	Human	ELISA, WB, IHC
<u>porB Antibody</u>	CSB-PA518750LA01NGH	Neisseria meningitidis serogroup B / serotype 15	ELISA
<u>PORCN Antibody</u>	CSB-PA887958LA01HU	Human, Mouse	ELISA, WB, IHC
<u>pos5 Antibody</u>	CSB-PA891853HA01SXV	Schizosaccharomyces pombe	ELISA
<u>POSTN Antibody</u>	CSB-PA09905A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>POT1 Antibody</u>	CSB-PA868326LA01HU	Human	ELISA, IHC, IF
<u>POU1F1 Antibody</u>	CSB-PA018391LA01HU	Human	ELISA, IHC
<u>POU2AF1 Antibody</u>	CSB-PA620973LA01HU	Human	ELISA
<u>POU2F1 Antibody</u>	CSB-PA018393LA01HU	Human	ELISA, IF
<u>POU2F2 Antibody</u>	CSB-PA018394LA01HU	Human	ELISA, WB, IHC, IF
<u>POU3F1 Antibody</u>	CSB-PA018396LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>POU3F2 Antibody</u>	CSB-PA018397HA01HU	Human	ELISA, IHC
<u>POU3F4 Antibody</u>	CSB-PA018399LA01HU	Human	ELISA, IHC
<u>POU4F1 Antibody</u>	CSB-PA22749A0Rb	Human	ELISA, WB, IF
<u>POU5F1 Antibody</u>	CSB-PA11579A0Rb	Human	ELISA, WB, IHC
<u>POU5F1B Antibody</u>	CSB-PA23129A0Rb	Human	ELISA, WB, IHC
<u>poxB Antibody</u>	CSB-PA17967A0Rb	Escherichia coli	ELISA
<u>PPA1 Antibody</u>	CSB-PA614884LA01HU	Human	ELISA, IHC, IF
<u>PPA2 Antibody</u>	CSB-PA867137LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>PPAN Antibody</u>	CSB-PA868284LA01HU	Human	ELISA
<u>PPARA Antibody</u>	CSB-PA018421LA01HU	Human	ELISA, WB
<u>PPARG Antibody</u>	CSB-PA018424LA01HU	Human	ELISA, IHC, IF
<u>PPBP Antibody</u>	CSB-PA09014A0Rb	Human	ELISA
<u>PPCDC Antibody</u>	CSB-PA018433LA01HU	Human	ELISA, IHC, IF
<u>PPDK Antibody</u>	CSB-PA336270LA01EKM	Entamoeba histolytica	ELISA
<u>PPDPFL Antibody</u>	CSB-PA840995LA01HU	Human	ELISA, IHC
<u>PPEF1 Antibody</u>	CSB-PA018435LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PPFIBP2 Antibody</u>	CSB-PA847688LA01HU	Human	ELISA, WB
<u>PPIA Antibody</u>	CSB-PA07814A0Rb	Human, Mouse	ELISA, WB
<u>PPIAL4G Antibody</u>	CSB-PA018453LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PPIC Antibody</u>	CSB-PA018473LA01HU	Human	ELISA, IF
<u>PPIE Antibody</u>	CSB-PA890944LA01HU	Human	ELISA, IHC
<u>PPIH Antibody</u>	CSB-PA018478LA01HU	Human	ELISA, WB, IHC
<u>PPIL3 Antibody</u>	CSB-PA863962LA01HU	Human	ELISA, WB, IHC
<u>PPIL4 Antibody</u>	CSB-PA851964LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PPM1F Antibody</u>	CSB-PA018493LA01HU	Human, Mouse	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>PPM1H Antibody</u>	CSB-PA892345LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PPM1J Antibody</u>	CSB-PA686382LA01HU	Human	ELISA, WB, IHC
<u>PPM1K Antibody</u>	CSB-PA822697LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PPM1L Antibody</u>	CSB-PA722509LA01HU	Human	ELISA, IHC
<u>PPM1M Antibody</u>	CSB-PA842717LA01HU	Human	ELISA, IHC
<u>PPME1 Antibody</u>	CSB-PA018501LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PPOX Antibody</u>	CSB-PA018502LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>PPP1CA Antibody</u>	CSB-PA018503LA01HU	Human	ELISA, IHC, IF
<u>PPP1CB Antibody</u>	CSB-PA018504LA01HU	Human	ELISA, WB, IHC
<u>Ppp1cb Antibody</u>	CSB-PA018504LA01MO	Mouse, Rat	ELISA, WB
<u>PPP1CC Antibody</u>	CSB-PA018505LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PPP1R11 Antibody</u>	CSB-PA03834A0Rb	Human	ELISA, IHC
<u>PPP1R12C Antibody</u>	CSB-PA871627LA01HU	Human	ELISA, IHC, IF
<u>PPP1R13B Antibody</u>	CSB-PA853449LA01HU	Human	ELISA, IHC
<u>PPP1R13L Antibody</u>	CSB-PA851965OA01HU	Human	ELISA, IHC, IF
<u>PPP1R14A Antibody</u>	CSB-PA850252LA01HU	Human	ELISA, IHC
<u>PPP1R16B Antibody</u>	CSB-PA842762LA01HU	Human	ELISA, IHC
<u>PPP1R3B Antibody</u>	CSB-PA769802LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>PPP1R7 Antibody</u>	CSB-PA622999LA01HU	Human	ELISA, IHC
<u>PPP2CA Antibody</u>	CSB-PA018559LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PPP2CB Antibody</u>	CSB-PA00985A0Rb	Human	ELISA, WB, IHC, IF
<u>PPP2R2A Antibody</u>	CSB-PA018564LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>PPP2R2B Antibody</u>	CSB-PA018565LA01HU	Human	ELISA, IHC, IF
<u>PPP2R2D Antibody</u>	CSB-PA727750LA01HU	Human	ELISA, WB, IHC, IF
<u>PPP2R3A Antibody</u>	CSB-PA21129A0Rb	Human	ELISA, IHC
<u>PPP2R3B Antibody</u>	CSB-PA896533HA01HU	Human	ELISA
<u>PPP2R5E Antibody</u>	CSB-PA619075LA01HU	Human	ELISA, IHC, IF
<u>PPP3R1 Antibody</u>	CSB-PA018575LA01HU	Human	ELISA, WB
<u>PPP5C Antibody</u>	CSB-PA018582LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PPP6C Antibody</u>	CSB-PA018583LA01HU	Human	ELISA, IHC, IF
<u>ppsA Antibody</u>	CSB-PA332921HA01ENV	Escherichia coli	ELISA, WB
<u>PPT1 Antibody</u>	CSB-PA018587LA01HU	Human	ELISA, WB, IHC
<u>PPTC7 Antibody</u>	CSB-PA836783LA01HU	Human	ELISA, IHC
<u>PPX1 Antibody</u>	CSB-PA334459HA01SVG	Saccharomyces cerevisiae	ELISA, WB
<u>PQBP1 Antibody</u>	CSB-PA018594LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>PQLC2 Antibody</u>	CSB-PA744399LA01HU	Human	ELISA, IF
<u>PQLC3 Antibody</u>	CSB-PA854059LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PR10A Antibody</u>	CSB-PA381182LA01CXF	Coptis japonica	ELISA, WB
<u>PRA1 Antibody</u>	CSB-PA308545YA01CZD	Candida albicans	ELISA, WB
<u>PRAC1 Antibody</u>	CSB-PA856965LA01HU	Human	ELISA
<u>PRAME Antibody</u>	CSB-PA018603LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PRAMEF10 Antibody</u>	CSB-PA018605LA01HU	Human	ELISA, IF
<u>PRAMEF17 Antibody</u>	CSB-PA018612LA01HU	Human	ELISA, WB
<u>PRAMEF18 Antibody</u>	CSB-PA018613LA01HU	Human	ELISA, WB, IHC
<u>PRAMEF5 Antibody</u>	CSB-PA018622LA01HU	Human	ELISA
<u>PRAMEF6 Antibody</u>	CSB-PA018623LA01HU	Human	ELISA, WB
<u>PRB1 Antibody</u>	CSB-PA018629LA01HU	Human	ELISA, WB
<u>PRCP Antibody</u>	CSB-PA018636LA01HU	Human	ELISA, IF
<u>PRDM10 Antibody</u>	CSB-PA878858LA01HU	Human	ELISA, IHC
<u>PRDM11 Antibody</u>	CSB-PA882077LA01HU	Human	ELISA, WB, IHC, IF
<u>PRDM13 Antibody</u>	CSB-PA887981LA01HU	Human, Mouse	ELISA, WB
<u>PRDM14 Antibody</u>	CSB-PA863934LA01HU	Human	ELISA, IHC
<u>PRDM16 Antibody</u>	CSB-PA872534LA01HU	Human	ELISA, IHC, IF
<u>PRDM2 Antibody</u>	CSB-PA613379LA01HU	Human	ELISA, IHC
<u>PRDM4 Antibody</u>	CSB-PA890920LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PRDM5 Antibody</u>	CSB-PA868292LA01HU	Human, Mouse	ELISA, WB
<u>PRDX1 Antibody</u>	CSB-PA01165A0Rb	Human	ELISA, WB, IHC, IF
<u>PRDX2 Antibody</u>	CSB-PA339235LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PRDX3 Antibody</u>	CSB-PA018656LA01HU	Human, Mouse	ELISA, WB, IHC, IF, IP
<u>PRDX4 Antibody</u>	CSB-PA621653LA01HU	Human	ELISA, WB
<u>PRDX6 Antibody</u>	CSB-PA018659LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>PREB Antibody</u>	CSB-PA867205LA01HU	Human, Rat	ELISA, WB, IHC
<u>PRELID2 Antibody</u>	CSB-PA850836LA01HU	Human	ELISA
<u>PRELID3A Antibody</u>	CSB-PA842722LA01HU	Human	ELISA
<u>PRELID3B Antibody</u>	CSB-PA896708LA01HU	Human	ELISA, IHC
<u>PREP Antibody</u>	CSB-PA018664LA01HU	Human	ELISA, WB, IHC, IF
<u>PREPL Antibody</u>	CSB-PA018665LA01HU	Human	ELISA, IHC
<u>PREX1 Antibody</u>	CSB-PA848403LA01HU	Human	ELISA, IHC, IF
<u>PREX2 Antibody</u>	CSB-PA741122LA01HU	Human, Mouse	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>PRF Antibody</u>	CSB-PA15107A0Rb	Toxoplasma gondii	ELISA
<u>PRF1 Antibody</u>	CSB-PA018668LA01HU	Human	ELISA
<u>Prf1 Antibody</u>	CSB-PA018668LA01MO	Mouse, Human	ELISA, WB, IF
<u>PRG2 Antibody</u>	CSB-PA018670LA01HU	Human	ELISA, WB, IHC, IF
<u>prgJ Antibody</u>	CSB-PA331347LA01SXB	Salmonella typhimurium	ELISA, WB
<u>PRICKLE3 Antibody</u>	CSB-PA018678PA01HU	Human	ELISA, WB, IHC
<u>PRICKLE4 Antibody</u>	CSB-PA652689LA01HU	Human	ELISA, IHC
<u>PRIM1 Antibody</u>	CSB-PA018680LA01HU	Human	ELISA, IF
<u>PRIMPOL Antibody</u>	CSB-PA822255LA01HU	Human	ELISA, IHC, IF
<u>PRKAA1 Antibody</u>	CSB-PA618759LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PRKAA2 Antibody</u>	CSB-PA805325LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PRKAB1 Antibody</u>	CSB-PA889918LA01HU	Human	ELISA, IF
<u>PRKACA Antibody</u>	CSB-PA018688LA01HU	Human	ELISA, WB, IHC
<u>PRKACB Antibody</u>	CSB-PA018689LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>PRKACG Antibody</u>	CSB-PA018690LA01HU	Human	ELISA, WB, IHC
<u>PRKAG1 Antibody</u>	CSB-PA01149A0Rb	Human, Rat	ELISA, WB, IHC, IF
<u>PRKAR2A Antibody</u>	CSB-PA018697LA01HU	Human	ELISA, WB, IHC, IF
<u>PRKCA Antibody</u>	CSB-PA15859A0Rb	Human	ELISA, IHC, IF
<u>PRKCB Antibody</u>	CSB-PA15849A0Rb	Human, Rat	ELISA, WB, IHC, IF
<u>PRKCD Antibody</u>	CSB-PA15899A0Rb	Human	ELISA, IF
<u>PRKCG Antibody</u>	CSB-PA15869A0Rb	Human	ELISA, IF
<u>PRKCH Antibody</u>	CSB-PA15879A0Rb	Human, Mouse, Rat	ELISA, WB, IHC
<u>PRKCSH Antibody</u>	CSB-PA02925A0Rb	Human	ELISA, IHC
<u>PRKCZ Antibody</u>	CSB-PA15889A0Rb	Human	ELISA, IF
<u>PRKD1 Antibody</u>	CSB-PA622989LA01HU	Human	ELISA, IF
<u>PRKD3 Antibody</u>	CSB-PA018713LA01HU	Human, Rat	ELISA, WB, IHC
<u>PRKDC Antibody</u>	CSB-PA018714LA01HU	Human	ELISA, WB
<u>PRKN Antibody</u>	CSB-PA017451LA01HU	Human	ELISA, WB, IHC, IF
<u>PRKRA Antibody</u>	CSB-PA13669A0Rb	Human	ELISA, IHC, IF
<u>PRL Antibody</u>	CSB-PA12344A0Rb	Human	ELISA
<u>PRL Antibody</u>	CSB-PA018724YA01HO	Horse	ELISA, WB
<u>Prl Antibody</u>	CSB-PA018724LA01RA	Rat	ELISA, IHC
<u>PRLR Antibody</u>	CSB-PA018727LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>PRMT1 Antibody</u>	CSB-PA03025A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>PRMT1 Antibody</u>	CSB-PA859113EA01HU	Human	ELISA, IHC, IF
<u>PRMT3 Antibody</u>	CSB-PA018733LA01HU	Human	ELISA, WB, IHC, IF
<u>PRMT5 Antibody</u>	CSB-PA018734LA01HU	Human	ELISA, IHC, IF
<u>PRMT6 Antibody</u>	CSB-PA856972LA01HU	Human	ELISA, WB, IHC, IF
<u>PRMT6 Antibody</u>	CSB-PA856972HA01HU	Human	ELISA, IHC, IF
<u>PRMT7 Antibody</u>	CSB-PA885738LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PRMT9 Antibody</u>	CSB-PA764718LA01HU	Human	ELISA, IHC
<u>prn Antibody</u>	CSB-PA321224LA01BUA	Bordetella pertussis	ELISA
<u>Prnp Antibody</u>	CSB-PA018739LA01RA	Rat	ELISA
<u>PRNP Antibody</u>	CSB-PA018739LA01HU	Human	ELISA, IHC, IF
<u>PRO1 Antibody</u>	CSB-PA517332HA01EQB	Cynodon dactylon	ELISA, WB
<u>Probable terminase, ATPase subunit Antibody</u>	CSB-PA343952LA01HAN	Haemophilus phage HP1	ELISA
<u>PROC Antibody</u>	CSB-PA018742LA01HU	Human	ELISA, IHC, IF
<u>PROCR Antibody</u>	CSB-PA06329A0Rb	Human, Mouse	ELISA, WB, IHC
<u>PRODH Antibody</u>	CSB-PA018744LA01HU	Human, Rat	ELISA, WB, IHC
<u>Profilin Antibody</u>	CSB-PA525079HA01MQZ	Mercurialis annua	ELISA, WB
<u>PROK1 Antibody</u>	CSB-PA018746LA01HU	Human	ELISA, IHC
<u>Prok2 Antibody</u>	CSB-PA879446LA01MO	Mouse	ELISA, WB
<u>PROKR1 Antibody</u>	CSB-PA018748LA01HU	Human	ELISA, IHC
<u>PROKR2 Antibody</u>	CSB-PA018749YA01HU	Human	ELISA, IHC
<u>PROM1 Antibody</u>	CSB-PA018751LA01HU	Human	ELISA, WB, IHC
<u>PROP1 Antibody</u>	CSB-PA018753LA01HU	Human	ELISA
<u>PROS1 Antibody</u>	CSB-PA018754LA01HU	Human	ELISA, IHC, IF
<u>Prothoracicotropic hormone Antibody</u>	CSB-PA324609LA01BTT	Bombyx mori	ELISA
<u>PROX1 Antibody</u>	CSB-PA852905LA01HU	Human	ELISA, WB
<u>PRPF3 Antibody</u>	CSB-PA018763LA01HU	Human	ELISA
<u>PRPF39 Antibody</u>	CSB-PA769774LA01HU	Human	ELISA, IHC, IF
<u>PRPF40A Antibody</u>	CSB-PA018769LA01HU	Human	ELISA, WB, IHC
<u>PRPF8 Antibody</u>	CSB-PA744036LA01HU	Human	ELISA, IF
<u>PRPS2 Antibody</u>	CSB-PA018779LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>PRPSAP1 Antibody</u>	CSB-PA614534LA01HU	Human	ELISA, WB
<u>PRR11 Antibody</u>	CSB-PA836225LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PRR3 Antibody</u>	CSB-PA018794LA01HU	Human	ELISA



Product Name	Code	Species Reactivity	Applications
<u>PRRC1 Antibody</u>	CSB-PA853460LA01HU	Human	ELISA, IHC
<u>PRRG3 Antibody</u>	CSB-PA866327LA01HU	Human	ELISA, IF
<u>Prss1 Antibody</u>	CSB-PA018811LA01RA	Rat	ELISA
<u>PRSS1 Antibody</u>	CSB-PA12309A0Rb	Human	ELISA, IHC
<u>PRSS12 Antibody</u>	CSB-PA018812LA01HU	Human	ELISA, IHC
<u>PRSS16 Antibody</u>	CSB-PA882073LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>PRSS2 Antibody</u>	CSB-PA018814LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PRSS23 Antibody</u>	CSB-PA018817LA01HU	Human	ELISA, IHC
<u>Prss29 Antibody</u>	CSB-PA859566LA01MO	Mouse	ELISA
<u>PRSS3 Antibody</u>	CSB-PA12319A0Rb	Human	ELISA, WB, IHC
<u>PRSS35 Antibody</u>	CSB-PA839792LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PRSS56 Antibody</u>	CSB-PA317372LA01HU	Human, Rat	ELISA, WB, IHC
<u>PRSS8 Antibody</u>	CSB-PA623091LA01HU	Human, Mouse	ELISA, WB, IF
<u>PRTFDC1 Antibody</u>	CSB-PA889092LA01HU	Human	ELISA, IHC
<u>PRTN3 Antibody</u>	CSB-PA018830LA01HU	Human	ELISA, IHC
<u>PRUNE1 Antibody</u>	CSB-PA771430LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>PRY Antibody</u>	CSB-PA018834LA01HU	Human	ELISA, IHC
<u>PSAP Antibody</u>	CSB-PA018836DA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>Psat1 Antibody</u>	CSB-PA860801HA01MO	Mouse	ELISA
<u>PSAT1 Antibody</u>	CSB-PA018838LA01HU	Human	ELISA, IHC, IF
<u>PSBO Antibody</u>	CSB-PA318952HA01FKI	Spinacia oleracea	ELISA
<u>PSD Antibody</u>	CSB-PA018842LA01HU	Human	ELISA, WB, IHC
<u>PSD4 Antibody</u>	CSB-PA843332LA01HU	Human	ELISA, IHC
<u>PSEN1 Antibody</u>	CSB-PA018846EA01HU	Human	ELISA, WB, IHC
<u>PSEN2 Antibody</u>	CSB-PA018847LA01HU	Human, Rat	ELISA, WB
<u>PSENNEN Antibody</u>	CSB-PA878932OA01HU	Human	ELISA, IHC
<u>PSG11 Antibody</u>	CSB-PA891577LA01HU	Human	ELISA, WB
<u>PSG3 Antibody</u>	CSB-PA618081LA01HU	Human	ELISA, WB, IHC
<u>PSG6 Antibody</u>	CSB-PA018856LA01HU	Human	ELISA, WB
<u>PSG9 Antibody</u>	CSB-PA27499A0Rb	Human	ELISA, IHC, IF
<u>PSKH2 Antibody</u>	CSB-PA822291LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>PSMA1 Antibody</u>	CSB-PA018865LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>PSMA2 Antibody</u>	CSB-PA018866LA01HU	Human	ELISA, IHC, IF
<u>PSMA3 Antibody</u>	CSB-PA018867LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>PSMA4 Antibody</u>	CSB-PA018869LA01HU	Human	ELISA, IHC, IF
<u>PSMA5 Antibody</u>	CSB-PA018870LA01HU	Human, Rat	ELISA, WB, IHC
<u>PSMA7 Antibody</u>	CSB-PA01335A0Rb	Human, Mouse	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>PSMB1 Antibody</u>	CSB-PA01395A0Rb	Human	ELISA, WB
<u>PSMB2 Antibody</u>	CSB-PA01385A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>PSMB4 Antibody</u>	CSB-PA01374A0Rb	Human	ELISA, WB, IHC, IF
<u>PSMB5 Antibody</u>	CSB-PA010083LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>PSMB7 Antibody</u>	CSB-PA01345A0Rb	Human	ELISA, WB, IHC, IF
<u>PSMB8 Antibody</u>	CSB-PA018886LA01HU	Human	ELISA, IHC
<u>PSMB9 Antibody</u>	CSB-PA018887LA01HU	Human	ELISA, WB, IHC
<u>PSMC3 Antibody</u>	CSB-PA018891LA01HU	Human	ELISA, IHC, IF
<u>PSMC4 Antibody</u>	CSB-PA018894LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PSMC5 Antibody</u>	CSB-PA018895LA01HU	Human	ELISA, WB, IF
<u>PSMC6 Antibody</u>	CSB-PA018896LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>PSMD10 Antibody</u>	CSB-PA018899LA01HU	Human	ELISA, IHC, IF
<u>PSMD11 Antibody</u>	CSB-PA00575A0Rb	Human	ELISA, IHC, IF
<u>PSMD14 Antibody</u>	CSB-PA018904OA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>PSMD2 Antibody</u>	CSB-PA619755LA01HU	Human	ELISA, WB, IHC, IF
<u>PSMD4 Antibody</u>	CSB-PA018908HA01HU	Human	ELISA, WB, IHC
<u>PSMD6 Antibody</u>	CSB-PA614875PA01HU	Human	ELISA, WB, IHC, IF
<u>PSMD7 Antibody</u>	CSB-PA018912LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>PSMD9 Antibody</u>	CSB-PA018914LA01HU	Human	ELISA, IHC, IF
<u>PSME2 Antibody</u>	CSB-PA01134A0Rb	Human	ELISA
<u>PSME3 Antibody</u>	CSB-PA13674A0Rb	Human	ELISA, IHC, IF
<u>PSME4 Antibody</u>	CSB-PA618917LA01HU	Human	ELISA, WB, IHC
<u>PSMF1 Antibody</u>	CSB-PA852882HA01HU	Human	ELISA, IHC
<u>PSMG1 Antibody</u>	CSB-PA018921LA01HU	Human	ELISA, WB, IHC, IF
<u>PSMG3 Antibody</u>	CSB-PA861132LA01HU	Human	ELISA
<u>PSORS1C1 Antibody</u>	CSB-PA883431LA01HU	Human	ELISA, IHC
<u>PSPC1 Antibody</u>	CSB-PA018937LA01HU	Human	ELISA
<u>PSTK Antibody</u>	CSB-PA816881LA01HU	Human	ELISA, WB
<u>PSTPIP2 Antibody</u>	CSB-PA884494LA01HU	Human	ELISA, WB, IHC
<u>pstS1 Antibody</u>	CSB-PA14919A0Rb	Mycobacterium tuberculosis	ELISA
<u>PTAFR Antibody</u>	CSB-PA018946LA01HU	Human, Mouse	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>PTBP1 Antibody</u>	CSB-PA00359A0Rb	Human	ELISA, IHC
<u>PTBP2 Antibody</u>	CSB-PA892133LA01HU	Human	ELISA, WB
<u>PTBP3 Antibody</u>	CSB-PA020060NA01HU	Human	ELISA, WB, IHC, IP
<u>PTCD2 Antibody</u>	CSB-PA823887LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PTCH1 Antibody</u>	CSB-PA623823LA01HU	Human	ELISA, IHC, IF
<u>PTCH2 Antibody</u>	CSB-PA896928LA01HU	Human	ELISA, IHC, IF
<u>PTDSS1 Antibody</u>	CSB-PA04305A0Rb	Human	ELISA, WB, IHC, IF
<u>PTF1A Antibody</u>	CSB-PA018967LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>PTGDS Antibody</u>	CSB-PA018969HA01HU	Human	ELISA, IHC, IF
<u>PTGER3 Antibody</u>	CSB-PA018972LA01HU	Human	ELISA, IHC, IF
<u>PTGER4 Antibody</u>	CSB-PA018973LA01HU	Human, Rat	ELISA, WB, IF
<u>PTGIS Antibody</u>	CSB-PA018982LA01HU	Human	ELISA, IHC
<u>PTGR1 Antibody</u>	CSB-PA018983LA01HU	Human	ELISA, IHC, IF
<u>PTGR2 Antibody</u>	CSB-PA018984LA01HU	Human	ELISA, WB, IHC, IF
<u>PTGS2 Antibody</u>	CSB-PA018986LA01HU	Human	ELISA, WB
<u>PTH Antibody</u>	CSB-PA09774A0Rb	Human	ELISA, IHC
<u>PTH Antibody</u>	CSB-PA09777A0Rb	Human	ELISA, IHC
<u>pth Antibody</u>	CSB-PA358958HA01ENV	Escherichia coli	ELISA
<u>PTH1R Antibody</u>	CSB-PA018988LA01HU	Human, Rat, Mouse	ELISA, WB
<u>PTH2R Antibody</u>	CSB-PA018990LA01HU	Human	ELISA, IF
<u>PTHLH Antibody</u>	CSB-PA018991LA01HU	Human	ELISA, IHC
<u>PTK2 Antibody</u>	CSB-PA14445A0Rb	Chicken	ELISA
<u>PTK2 Antibody</u>	CSB-PA018994LA01HU	Human	ELISA, IHC, IF
<u>PTK2B Antibody</u>	CSB-PA613494LA01HU	Human	ELISA, IHC, IF
<u>PTK6 Antibody</u>	CSB-PA614448LA01HU	Human	ELISA, WB, IF
<u>PTK7 Antibody</u>	CSB-PA018996LA01HU	Human	ELISA, IHC
<u>PTMA Antibody</u>	CSB-PA01524A0Rb	Human	ELISA, IHC
<u>PTN Antibody</u>	CSB-PA10239A0Rb	Human	ELISA, IHC
<u>PTOV1 Antibody</u>	CSB-PA769809LA01HU	Human	ELISA, IHC, IF
<u>PTP4A1 Antibody</u>	CSB-PA846294LA01HU	Human	ELISA
<u>PTP4A2 Antibody</u>	CSB-PA613277LA01HU	Human	ELISA, IHC, IF
<u>PTPA Antibody</u>	CSB-PA620876LA01HU	Human	ELISA, WB, IHC, IF
<u>PTPMT1 Antibody</u>	CSB-PA819889LA01HU	Human, Rat	ELISA, WB, IHC
<u>PTPN11 Antibody</u>	CSB-PA11779A0Rb	Human, Mouse	ELISA, WB
<u>PTPN12 Antibody</u>	CSB-PA019026LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>PTPN13 Antibody</u>	CSB-PA614263LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>PTPN14 Antibody</u>	CSB-PA623010LA01HU	Human	ELISA, IF
<u>PTPN22 Antibody</u>	CSB-PA019036LA01HU	Human	ELISA, IF
<u>PTPN3 Antibody</u>	CSB-PA019040LA01HU	Human	ELISA, IHC, IF
<u>PTPN5 Antibody</u>	CSB-PA019042LA01HU	Human	ELISA, WB, IHC
<u>PTPN6 Antibody</u>	CSB-PA019043HA01HU	Human	ELISA, WB, IHC, IF
<u>PTPN9 Antibody</u>	CSB-PA019045LA01HU	Human	ELISA, IHC
<u>PTPRF Antibody</u>	CSB-PA019053LA01HU	Human	ELISA, IHC, IF
<u>PTPRK Antibody</u>	CSB-PA613588LA01HU	Human	ELISA, IHC, IF
<u>PTPRN Antibody</u>	CSB-PA613703LA01HU	Human	ELISA
<u>PTPRR Antibody</u>	CSB-PA624019LA01HU	Human	ELISA, IHC
<u>PTPRS Antibody</u>	CSB-PA613393LA01HU	Human	ELISA, IHC
<u>ptprsa Antibody</u>	CSB-PA24489A0Rb	Zebrafish	ELISA, WB
<u>PTPRZ1 Antibody</u>	CSB-PA019068LA01HU	Human	ELISA, IHC
<u>PTPRZ1 Antibody</u>	CSB-PA019068YA01HU	Human	ELISA, IHC, IF
<u>ptsH Antibody</u>	CSB-PA362538HA01BRJ	Bacillus subtilis	ELISA, WB
<u>PTTG1 Antibody</u>	CSB-PA019074LA01HU	Human	ELISA, WB
<u>PTTG1IP Antibody</u>	CSB-PA00859A0Rb	Human	ELISA, IHC, IF
<u>PTTG2 Antibody</u>	CSB-PA878936LA01HU	Human	ELISA, IF
<u>PTX3 Antibody</u>	CSB-PA019079LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>ptxA Antibody</u>	CSB-PA356423LA01Bua	Bordetella pertussis	ELISA, WB
<u>PUF60 Antibody</u>	CSB-PA883409LA01HU	Human	ELISA, WB, IHC
<u>PURA Antibody</u>	CSB-PA019083LA01HU	Human	ELISA, WB
<u>PURG Antibody</u>	CSB-PA891529LA01HU	Human	ELISA, WB
<u>PVALB Antibody</u>	CSB-PA019092HA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>pvdA Antibody</u>	CSB-PA706553HA01EZX	Pseudomonas aeruginosa	ELISA, WB
<u>PVRIG Antibody</u>	CSB-PA721905LA01HU	Human, Mouse	ELISA, WB, IF
<u>PXK Antibody</u>	CSB-PA770363LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>PXN Antibody</u>	CSB-PA019112OA01HU	Human	ELISA, WB, IHC
<u>PYC1 Antibody</u>	CSB-PA319972LA01SVG	Baker's yeast	ELISA
<u>PYCARD Antibody</u>	CSB-PA890936LA01HU	Human, Mouse	ELISA, WB, IHC
<u>PYCR1 Antibody</u>	CSB-PA019115HA01HU	Human	ELISA, IHC
<u>PYCR2 Antibody</u>	CSB-PA822188LA01HU	Human	ELISA, IHC
<u>PYGB Antibody</u>	CSB-PA10065A0Rb	Human, Mouse	ELISA, WB, IHC
<u>PYGL Antibody</u>	CSB-PA019122EA01HU	Human	ELISA, WB, IHC
<u>PYGO2 Antibody</u>	CSB-PA883586LA01HU	Human	ELISA, IHC, IF
<u>PYHIN1 Antibody</u>	CSB-PA746957LA01HU	Human, Mouse	ELISA, WB
<u>pyrE Antibody</u>	CSB-PA506933HA01LNO	Laribacter hongkongensis	ELISA, WB
<u>pyrF Antibody</u>	CSB-PA14337A0Rb	Helicobacter pylori	ELISA
<u>PYY Antibody</u>	CSB-PA12214A0Rb	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>PZP Antibody</u>	CSB-PA019131HA01HU	Human	ELISA
<u>QDPR Antibody</u>	CSB-PA019133LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>QPCT Antibody</u>	CSB-PA019135LA01HU	Human, Rat	ELISA, WB, IHC
<u>QPCTL Antibody</u>	CSB-PA868351LA01HU	Human	ELISA, IF
<u>QPRT Antibody</u>	CSB-PA621868LA01HU	Human	ELISA, IHC, IF
<u>QSOX1 Antibody</u>	CSB-PA019144LA01HU	Human, Rat	ELISA, WB, IHC
<u>QTRT1 Antibody</u>	CSB-PA887150LA01HU	Human	ELISA, WB, IHC
<u>R107.2 Antibody</u>	CSB-PA333729LA01CXY	Caenorhabditis elegans	ELISA
<u>R3HDML Antibody</u>	CSB-PA872491LA01HU	Human	ELISA
<u>RAB10 Antibody</u>	CSB-PA019152LA01HU	Human	ELISA, IHC, IF
<u>RAB11A Antibody</u>	CSB-PA019153HA01HU	Human	ELISA, IHC
<u>RAB11B Antibody</u>	CSB-PA614903LA01HU	Human, Rat	ELISA, WB, IHC
<u>RAB11FIP5 Antibody</u>	CSB-PA863141LA01HU	Human	ELISA, IHC
<u>RAB12 Antibody</u>	CSB-PA753710LA01HU	Human	ELISA, WB, IF
<u>RAB1A Antibody</u>	CSB-PA00335A0Rb	Human	ELISA, WB
<u>RAB1B Antibody</u>	CSB-PA872427LA01HU	Human	ELISA, IF
<u>RAB20 Antibody</u>	CSB-PA878910LA01HU	Human	ELISA, WB, IHC
<u>RAB21 Antibody</u>	CSB-PA891733LA01HU	Human	ELISA, IHC, IF
<u>RAB24 Antibody</u>	CSB-PA842627LA01HU	Human	ELISA, IF
<u>RAB25 Antibody</u>	CSB-PA019175LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>RAB26 Antibody</u>	CSB-PA892347LA01HU	Human, Rat	ELISA, WB, IHC
<u>RAB27A Antibody</u>	CSB-PA019177LA01HU	Human	ELISA, IHC, IF
<u>RAB27B Antibody</u>	CSB-PA019178HA01HU	Human	ELISA
<u>RAB28 Antibody</u>	CSB-PA019179LA01HU	Human	ELISA, IHC, IF
<u>RAB2A Antibody</u>	CSB-PA019180LA01HU	Human	ELISA
<u>RAB30 Antibody</u>	CSB-PA613602LA01HU	Human	ELISA, WB, IHC
<u>RAB32 Antibody</u>	CSB-PA618795LA01HU	Human, Mouse	ELISA, WB
<u>RAB33A Antibody</u>	CSB-PA613485LA01HU	Human	ELISA, IHC, IF
<u>RAB35 Antibody</u>	CSB-PA614887LA01HU	Human	ELISA, IHC, IF
<u>RAB37 Antibody</u>	CSB-PA842637LA01HU	Human	ELISA, IHC
<u>RAB38 Antibody</u>	CSB-PA019191LA01HU	Human	ELISA
<u>RAB39A Antibody</u>	CSB-PA614818LA01HU	Human	ELISA, IF
<u>RAB39B Antibody</u>	CSB-PA822194LA01HU	Human	ELISA, IHC, IF
<u>RAB3A Antibody</u>	CSB-PA019194LA01HU	Human	ELISA, IHC, IF
<u>RAB3B Antibody</u>	CSB-PA019195LA01HU	Human, Mouse	ELISA, WB, IHC
<u>RAB3C Antibody</u>	CSB-PA846609HA01HU	Human, Mouse	ELISA, WB
<u>RAB3GAP1 Antibody</u>	CSB-PA621861LA01HU	Human	ELISA, WB, IHC, IF
<u>RAB3GAP2 Antibody</u>	CSB-PA872480LA01HU	Human	ELISA, WB, IHC
<u>RAB3IL1 Antibody</u>	CSB-PA837424LA01HU	Human	ELISA



Product Name	Code	Species Reactivity	Applications
<u>RAB3IP Antibody</u>	CSB-PA853488LA01HU	Human	ELISA, IHC, IF
<u>RAB40AL Antibody</u>	CSB-PA019203LA01HU	Human	ELISA, IF
<u>RAB40B Antibody</u>	CSB-PA615534LA01HU	Human	ELISA, IHC
<u>RAB42 Antibody</u>	CSB-PA843268LA01HU	Human	ELISA, IF
<u>RAB4B Antibody</u>	CSB-PA019212LA01HU	Human, Mouse	ELISA, WB, IHC
<u>RAB5A Antibody</u>	CSB-PA019213HA01HU	Human	ELISA, IHC
<u>RAB5A Antibody</u>	CSB-PA019213EA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>RAB6A Antibody</u>	CSB-PA019216HA01HU	Human	ELISA, IHC
<u>RAB6B Antibody</u>	CSB-PA873642LA01HU	Human, Mouse	ELISA, WB, IHC
<u>RAB6C Antibody</u>	CSB-PA875643LA01HU	Human	ELISA, IHC
<u>RAB7A Antibody</u>	CSB-PA019219LA01HU	Human	ELISA, IHC, IF
<u>RAB7B Antibody</u>	CSB-PA019220LA01HU	Human	ELISA, IF
<u>RAB8B Antibody</u>	CSB-PA838820LA01HU	Human	ELISA, IHC
<u>RABEP1 Antibody</u>	CSB-PA614886LA01HU	Human	ELISA, WB, IHC
<u>RABEP2 Antibody</u>	CSB-PA880980LA01HU	Human	ELISA, IHC
<u>RABGAP1L Antibody</u>	CSB-PA732584LA01HU	Human, Mouse	ELISA, WB, IHC
<u>RABGEF1 Antibody</u>	CSB-PA883453LA01HU	Human	ELISA, IHC
<u>RABGGTA Antibody</u>	CSB-PA852898LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>RABGGTB Antibody</u>	CSB-PA019235HA01HU	Human	ELISA, IHC, IF
<u>RABIF Antibody</u>	CSB-PA019236LA01HU	Human, Mouse	ELISA, WB
<u>RAC1 Antibody</u>	CSB-PA04329A0Rb	Human	ELISA, IHC, IF
<u>RAC2 Antibody</u>	CSB-PA019248LA01HU	Human	ELISA, IHC
<u>RACGAP1 Antibody</u>	CSB-PA884432LA01HU	Human	ELISA, IHC, IF
<u>RACK1 Antibody</u>	CSB-PA009605LA01HU	Human	ELISA, IHC, IF
<u>RAD17 Antibody</u>	CSB-PA019253EA01HU	Human, Mouse	ELISA, WB, IHC
<u>RAD18 Antibody</u>	CSB-PA889102LA01HU	Human	ELISA, WB, IHC
<u>RAD23A Antibody</u>	CSB-PA03774A0Rb	Human	ELISA, WB, IHC
<u>RAD23B Antibody</u>	CSB-PA019260LA01HU	Human	ELISA, WB, IHC
<u>RAD51 Antibody</u>	CSB-PA019268LA01HU	Human, Mouse	ELISA, WB, IHC, IF, IP
<u>RAD54B Antibody</u>	CSB-PA896748LA01HU	Human	ELISA, IHC, IF
<u>RAD54L2 Antibody</u>	CSB-PA896892LA01HU	Human, Rat	ELISA, WB, IF
<u>RAD9A Antibody</u>	CSB-PA857455LA01HU	Human	ELISA, IHC
<u>RAD9B Antibody</u>	CSB-PA740946LA01HU	Human, Mouse	ELISA, WB
<u>RAET1E Antibody</u>	CSB-PA855052LA01HU	Human	ELISA, IHC, IF
<u>RAF1 Antibody</u>	CSB-PA019284LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>RAG2 Antibody</u>	CSB-PA345128LA01HU	Human	ELISA, IF
<u>RAI1 Antibody</u>	CSB-PA768772LA01HU	Human	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>RAI14 Antibody</u>	CSB-PA868376LA01HU	Human	ELISA, IHC
<u>RALA Antibody</u>	CSB-PA019296LA01HU	Human	ELISA, IHC, IF
<u>RALBP1 Antibody</u>	CSB-PA618981LA01HU	Human	ELISA, WB, IHC, IF
<u>RALGPS1 Antibody</u>	CSB-PA683358HA01HU	Human	ELISA
<u>RALY Antibody</u>	CSB-PA892332LA01HU	Human	ELISA, WB, IHC, IF
<u>RALYL Antibody</u>	CSB-PA771422LA01HU	Human	ELISA, IHC
<u>RAMP2 Antibody</u>	CSB-PA019305LA01HU	Human	ELISA, IHC, IF
<u>RANBP1 Antibody</u>	CSB-PA019308HA01HU	Human	ELISA, IHC, IF
<u>RANBP10 Antibody</u>	CSB-PA764955LA01HU	Human	ELISA, WB, IHC, IP
<u>RANBP17 Antibody</u>	CSB-PA872481LA01HU	Human	ELISA, WB, IHC
<u>RANBP2 Antibody</u>	CSB-PA019311LA01HU	Human	ELISA, IHC
<u>RANBP3L Antibody</u>	CSB-PA803130LA01HU	Human	ELISA, IHC, IF
<u>RANBP9 Antibody</u>	CSB-PA853501OA01HU	Human	ELISA, WB, IHC
<u>RAP1A Antibody</u>	CSB-PA019321LA01HU	Human	ELISA, IHC, IF
<u>RAP1B Antibody</u>	CSB-PA03269A0Rb	Human, Mouse	ELISA, WB, IHC
<u>RAPGEF1 Antibody</u>	CSB-PA613430NA01HU	Human	ELISA, WB
<u>RAPGEF2 Antibody</u>	CSB-PA897493LA01HU	Human	ELISA, IF
<u>RAPGEF5 Antibody</u>	CSB-PA821676LA01HU	Human	ELISA, IHC
<u>RARB Antibody</u>	CSB-PA019339LA01HU	Human	ELISA, IF
<u>RARG Antibody</u>	CSB-PA019340LA01HU	Human, Mouse	ELISA, WB, IHC
<u>RARRES2 Antibody</u>	CSB-PA10759A0Rb	Human	ELISA, IHC, IF
<u>RARRES3 Antibody</u>	CSB-PA892340LA01HU	Human	ELISA, IHC
<u>RARS Antibody</u>	CSB-PA019344LA01HU	Human	ELISA, WB, IHC, IF
<u>RASA1 Antibody</u>	CSB-PA019346LA01HU	Human, Mouse	ELISA, WB, IF
<u>RASA2 Antibody</u>	CSB-PA622992LA01HU	Human	ELISA, IF
<u>RASAL1 Antibody</u>	CSB-PA019351LA01HU	Human	ELISA, IF
<u>RASAL2 Antibody</u>	CSB-PA890908LA01HU	Human	ELISA, IHC, IF
<u>RASGEF1A Antibody</u>	CSB-PA854076LA01HU	Human	ELISA, IF
<u>RASGEF1B Antibody</u>	CSB-PA610807LA01HU	Human	ELISA, IHC
<u>RASGRP1 Antibody</u>	CSB-PA019362LA01HU	Human	ELISA, IF
<u>RASL10B Antibody</u>	CSB-PA842758LA01HU	Human	ELISA, IHC
<u>RASL11B Antibody</u>	CSB-PA859410LA01HU	Human	ELISA, IHC, IF
<u>RASL12 Antibody</u>	CSB-PA878926LA01HU	Human	ELISA
<u>RASSE1 Antibody</u>	CSB-PA885715LA01HU	Human	ELISA, IHC
<u>RASSE2 Antibody</u>	CSB-PA019374LA01HU	Human	ELISA, WB, IP
<u>RASSE3 Antibody</u>	CSB-PA801247LA01HU	Human	ELISA, IHC, IF
<u>RASSE4 Antibody</u>	CSB-PA880956LA01HU	Human	ELISA
<u>RASSE7 Antibody</u>	CSB-PA863596LA01HU	Human, Mouse	ELISA, WB, IF



Product Name	Code	Species Reactivity	Applications
<u>RAVER2 Antibody</u>	CSB-PA872547LA01HU	Human	ELISA, WB, IHC, IF
<u>RB1CC1 Antibody</u>	CSB-PA855061LA01HU	Human	ELISA, IF
<u>RB38 Antibody</u>	CSB-PA24729A0Rb	Chlamydomonas reinhardtii	ELISA, WB
<u>RBBP4 Antibody</u>	CSB-PA600201LA01HU	Human, Mouse	ELISA, WB, IHC
<u>RBBP6 Antibody</u>	CSB-PA768776LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>RBBP7 Antibody</u>	CSB-PA621959LA01HU	Human	ELISA, WB
<u>RBBP8 Antibody</u>	CSB-PA857871LA01HU	Human	ELISA, IHC
<u>RBBP9 Antibody</u>	CSB-PA019394LA01HU	Human	ELISA, WB, IHC
<u>RBCK1 Antibody</u>	CSB-PA880138LA01HU	Human	ELISA, IHC
<u>rbcL Antibody</u>	CSB-PA328678LA01GGV	Soybean	ELISA, WB
<u>rbcL Antibody</u>	CSB-PA517305LA01DOA	Arabidopsis thaliana	ELISA, WB
<u>RBFA Antibody</u>	CSB-PA836636LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>RBFOX2 Antibody</u>	CSB-PA019439LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>RBFOX3 Antibody</u>	CSB-PA407385LA01HU	Human	ELISA, IHC
<u>RBKS Antibody</u>	CSB-PA019397LA01HU	Human, Rat	ELISA, WB, IHC
<u>RBL1 Antibody</u>	CSB-PA019398LA01HU	Human, Mouse	ELISA, WB
<u>RBL2 Antibody</u>	CSB-PA600921LA01HU	Human	ELISA, IHC
<u>RBM10 Antibody</u>	CSB-PA019400LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>RBM11 Antibody</u>	CSB-PA019401LA01HU	Human	ELISA, IHC
<u>RBM12 Antibody</u>	CSB-PA865137LA01HU	Human	ELISA, IHC, IF
<u>RBM22 Antibody</u>	CSB-PA889126LA01HU	Human, Mouse	ELISA, WB, IHC
<u>rbm24 Antibody</u>	CSB-PA745088LA01DIL	Zebrafish	ELISA, WB
<u>RBM3 Antibody</u>	CSB-PA019420HA01HU	Human	ELISA, IHC
<u>RBM38 Antibody</u>	CSB-PA880939LA01HU	Human	ELISA, IHC, IF
<u>RBM4 Antibody</u>	CSB-PA880118LA01HU	Human, Mouse	ELISA, WB, IHC
<u>RBM41 Antibody</u>	CSB-PA850315LA01HU	Human	ELISA, WB, IHC
<u>RBM45 Antibody</u>	CSB-PA814203LA01HU	Human, Rat	ELISA, WB, IF
<u>rbm47 Antibody</u>	CSB-PA23979A0Rb	Zebrafish	ELISA
<u>RBM6 Antibody</u>	CSB-PA019435LA01HU	Human	ELISA
<u>RBM7 Antibody</u>	CSB-PA896729LA01HU	Human	ELISA
<u>RBMS1 Antibody</u>	CSB-PA019440LA01HU	Human	ELISA, WB, IF
<u>RBMS2 Antibody</u>	CSB-PA613593LA01HU	Human	ELISA, WB, IHC
<u>RBMS3 Antibody</u>	CSB-PA761548LA01HU	Human	ELISA, IF
<u>RBMY1A1 Antibody</u>	CSB-PA319143LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>RBMY1F Antibody</u>	CSB-PA614891LA01HU	Human	ELISA, WB, IHC
<u>RBP1 Antibody</u>	CSB-PA019480YA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>RBP3 Antibody</u>	CSB-PA019482LA01HU	Human	ELISA
<u>RBP4 Antibody</u>	CSB-PA019483YA01HU	Human	ELISA
<u>RBP4 Antibody</u>	CSB-PA07057A0Rb	Human	ELISA, WB, IHC
<u>RBP4 Antibody</u>	CSB-PA07059A0Rb	Human	ELISA
<u>RBP5 Antibody</u>	CSB-PA019484LA01HU	Human, Mouse	ELISA, WB
<u>RBPJ Antibody</u>	CSB-PA019486LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>RBPJL Antibody</u>	CSB-PA883361LA01HU	Human	ELISA, WB, IHC
<u>RBPMS Antibody</u>	CSB-PA04569A0Rb	Human	ELISA, IHC, IF
<u>RBPMS Antibody</u>	CSB-PA821892LA01HU	Human	ELISA, IHC
<u>rbsK Antibody</u>	CSB-PA18707A0Rb	Escherichia coli	ELISA, WB
<u>RBSN Antibody</u>	CSB-PA884441LA01HU	Human	ELISA
<u>rbtD Antibody</u>	CSB-PA360400HA01EKN	Enterobacter aerogenes	ELISA
<u>RC3H2 Antibody</u>	CSB-PA884510LA01HU	Human	ELISA
<u>RCAN1 Antibody</u>	CSB-PA019500LA01HU	Human	ELISA, IHC
<u>RCAN2 Antibody</u>	CSB-PA618879LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>RCBTB2 Antibody</u>	CSB-PA019504LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>RCE1 Antibody</u>	CSB-PA896686LA01HU	Human	ELISA, WB
<u>RCN3 Antibody</u>	CSB-PA856914LA01HU	Human, Mouse	ELISA, WB, IHC
<u>RCOR1 Antibody</u>	CSB-PA019515HA01HU	Human	ELISA, IHC
<u>RCOR3 Antibody</u>	CSB-PA019517LA01HU	Human	ELISA, IHC, IF
<u>RCVRN Antibody</u>	CSB-PA019519HA01HU	Human, Mouse	ELISA, WB
<u>RD3 Antibody</u>	CSB-PA800101LA01HU	Human, Mouse	ELISA, WB, IHC
<u>RDH11 Antibody</u>	CSB-PA019523LA01HU	Human	ELISA, WB, IHC
<u>RDH12 Antibody</u>	CSB-PA842727LA01HU	Human	ELISA, IHC
<u>recF Antibody</u>	CSB-PA358973LA01ENV	Escherichia coli	ELISA
<u>RECK Antibody</u>	CSB-PA019536LA01HU	Human	ELISA, IHC, IF
<u>recO Antibody</u>	CSB-PA364006HA01ENV	Escherichia coli	ELISA
<u>recQ Antibody</u>	CSB-PA17979A0Rb	Escherichia coli	ELISA, WB
<u>RECQL4 Antibody</u>	CSB-PA019538LA01HU	Human	ELISA, WB, IHC, IF
<u>recR Antibody</u>	CSB-PA358976HA01ENV	Escherichia coli	ELISA
<u>REEP2 Antibody</u>	CSB-PA887109LA01HU	Human	ELISA, IHC, IF
<u>REEP3 Antibody</u>	CSB-PA754137LA01HU	Human	ELISA, IHC, IF
<u>REEP4 Antibody</u>	CSB-PA875682LA01HU	Human	ELISA, WB, IHC, IF
<u>REG1A Antibody</u>	CSB-PA019546LA01HU	Human	ELISA, IHC
<u>REG3A Antibody</u>	CSB-PA019548EA01HU	Human	ELISA, IHC, IF
<u>Reg3a Antibody</u>	CSB-PA019548LA01RA	Rat	ELISA
<u>Reg3b Antibody</u>	CSB-PA326510LA01RA	Rat	ELISA



Product Name	Code	Species Reactivity	Applications
<u>Reg3g Antibody</u>	CSB-PA019549LA01RA	Rat, Mouse	ELISA, WB
<u>REG3G Antibody</u>	CSB-PA751009EA01HU	Human	ELISA
<u>REG4 Antibody</u>	CSB-PA887169LA01HU	Human	ELISA, IHC
<u>RELA Antibody</u>	CSB-PA03987A0Rb	Human, Mouse	ELISA, WB, IHC, IF, ChIP
<u>RELB Antibody</u>	CSB-PA019554LA01HU	Human, Mouse	ELISA, WB, IF
<u>RELN Antibody</u>	CSB-PA019557LA01HU	Human	ELISA, IHC
<u>RELT Antibody</u>	CSB-PA08179A0Rb	Human	ELISA, WB, IHC, IF
<u>REM1 Antibody</u>	CSB-PA019559LA01HU	Human	ELISA, IF
<u>REN Antibody</u>	CSB-PA12069A0Rb	Human, Mouse	ELISA, WB, IF
<u>Ren2 Antibody</u>	CSB-PA17059A0Rb	Mouse, Human	ELISA, WB
<u>repD Antibody</u>	CSB-PA360946HA01FKZ	Staphylococcus aureus	ELISA
<u>REPS1 Antibody</u>	CSB-PA856915LA01HU	Human	ELISA, IHC, IF
<u>RER1 Antibody</u>	CSB-PA019567LA01HU	Human	ELISA, IHC, IF
<u>REG Antibody</u>	CSB-PA850254LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>REST Antibody</u>	CSB-PA019571LA01HU	Human	ELISA, IF
<u>RETN Antibody</u>	CSB-PA019573LA01BO	Bovine	ELISA
<u>RETN Antibody</u>	CSB-PA07079A0Rb	Human	ELISA, IHC
<u>REV1 Antibody</u>	CSB-PA883372LA01HU	Human	ELISA, IHC
<u>REXO4 Antibody</u>	CSB-PA863931EA01HU	Human	ELISA, IHC
<u>RFC4 Antibody</u>	CSB-PA019591LA01HU	Human, Mouse	ELISA, WB, IF
<u>RFESD Antibody</u>	CSB-PA840551LA01HU	Human	ELISA, IHC
<u>RFK Antibody</u>	CSB-PA836167LA01HU	Human	ELISA, IHC, IF
<u>RFNG Antibody</u>	CSB-PA896922LA01HU	Human, Rat	ELISA, WB, IHC
<u>RFPL3 Antibody</u>	CSB-PA019600LA01HU	Human	ELISA, WB
<u>RFPL4B Antibody</u>	CSB-PA019603LA01HU	Human	ELISA, IHC, IF
<u>RFTN2 Antibody</u>	CSB-PA687456LA01HU	Human	ELISA, IHC, IF
<u>RFX1 Antibody</u>	CSB-PA019610LA01HU	Human	ELISA, WB, IHC
<u>RFX2 Antibody</u>	CSB-PA019611LA01HU	Human	ELISA, IHC, IF
<u>RFX3 Antibody</u>	CSB-PA019612LA01HU	Human	ELISA
<u>RFX6 Antibody</u>	CSB-PA812809LA01HU	Human	ELISA, IHC
<u>RFXANK Antibody</u>	CSB-PA019617LA01HU	Human	ELISA, WB, IHC, IF
<u>RGL2 Antibody</u>	CSB-PA019625LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>RGMA Antibody</u>	CSB-PA836186LA01HU	Human, Mouse	ELISA, WB, IHC
<u>Rgn Antibody</u>	CSB-PA019630LA01RA	Rat, Human	ELISA, WB
<u>RGN Antibody</u>	CSB-PA623001LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>RGPI Antibody</u>	CSB-PA838781LA01HU	Human	ELISA
<u>rgpA Antibody</u>	CSB-PA338957LA01PQP	Porphyromonas gingivalis	ELISA



Product Name	Code	Species Reactivity	Applications
<u>rgpB Antibody</u>	CSB-PA310587LA01EYA	Porphyromonas gingivalis	ELISA
<u>RGS1 Antibody</u>	CSB-PA019641LA01HU	Human	ELISA, IHC, IF
<u>RGS18 Antibody</u>	CSB-PA019649LA01HU	Human	ELISA, IHC, IF
<u>RGS19 Antibody</u>	CSB-PA019650LA01HU	Human	ELISA, IHC, IF
<u>RGS2 Antibody</u>	CSB-PA019651LA01HU	Human	ELISA, WB, IHC, IF
<u>RGS22 Antibody</u>	CSB-PA019654LA01HU	Human	ELISA
<u>RGS4 Antibody</u>	CSB-PA019656LA01HU	Human	ELISA, IHC, IF
<u>RGS6 Antibody</u>	CSB-PA019658LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>RGS7 Antibody</u>	CSB-PA019659LA01HU	Human	ELISA, IHC
<u>RHBDD3 Antibody</u>	CSB-PA27742A0Rb	Human	ELISA, WB, IHC, IF
<u>RHBDF1 Antibody</u>	CSB-PA856909LA01HU	Human	ELISA, IHC
<u>RHBDL2 Antibody</u>	CSB-PA019672LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>RHCG Antibody</u>	CSB-PA880059LA01HU	Human, Rat	ELISA, WB, IHC
<u>RHD Antibody</u>	CSB-PA019677LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>RHEB Antibody</u>	CSB-PA624022LA01HU	Human	ELISA, IHC
<u>RHEBL1 Antibody</u>	CSB-PA855034LA01HU	Human	ELISA, IHC
<u>rho Antibody</u>	CSB-PA16627A0Rb	Escherichia coli	ELISA
<u>rho Antibody</u>	CSB-PA16629A0Rb	Escherichia coli	ELISA
<u>RHOA Antibody</u>	CSB-PA00875A0Rb	Human	ELISA, IHC
<u>rhoaa Antibody</u>	CSB-PA24499A0Rb	Zebrafish	ELISA, WB
<u>rhoac Antibody</u>	CSB-PA24509A0Rb	Zebrafish	ELISA, WB
<u>RHOB Antibody</u>	CSB-PA019683LA01HU	Human	ELISA, IF
<u>RHOBTB3 Antibody</u>	CSB-PA019686LA01HU	Human	ELISA, WB
<u>RHOC Antibody</u>	CSB-PA019687LA01HU	Human	ELISA, IHC
<u>RHOD Antibody</u>	CSB-PA019688LA01HU	Human	ELISA, IF
<u>RHOF Antibody</u>	CSB-PA867192LA01HU	Human, Mouse	ELISA, WB, IHC
<u>RHOH Antibody</u>	CSB-PA618992LA01HU	Human	ELISA, IHC, IF
<u>RHOJ Antibody</u>	CSB-PA872495LA01HU	Human	ELISA, IHC
<u>RHOT2 Antibody</u>	CSB-PA810289LA01HU	Human	ELISA, IHC
<u>RHOU Antibody</u>	CSB-PA745331LA01HU	Human	ELISA, IHC, IF
<u>RHOXF2 Antibody</u>	CSB-PA861117LA01HU	Human	ELISA, IHC
<u>RHPN1-AS1 Antibody</u>	CSB-PA871590LA01HU	Human	ELISA, IHC
<u>RIBC1 Antibody</u>	CSB-PA836661LA01HU	Human	ELISA, IHC
<u>RIBC2 Antibody</u>	CSB-PA884469LA01HU	Human	ELISA, WB, IHC
<u>Riboflavin-binding protein Antibody</u>	CSB-PA360879YA01CH	Chicken	ELISA
<u>Ribosome-inactivating protein</u>	CSB-PA357829HA01TIF	Trichosanthes kirilowii	ELISA



Product Name	Code	Species Reactivity	Applications
<u>alpha-trichosanthin Antibody</u>			
<u>Ribosome-inactivating protein bryodin I Antibody</u>	CSB-PA339327LA01BMR	Bryonia dioica	ELISA
<u>Ribosome-inactivating protein charybdin Antibody</u>	CSB-PA308125HA01DEZ	Drimia maritima	ELISA
<u>Ribosome-inactivating protein cucurmosin Antibody</u>	CSB-PA307485LA01CUI	Cucurbita moschata	ELISA, WB
<u>Ribosome-inactivating protein gynostemmin Antibody</u>	CSB-PA688696LA01GAAZ	Gynostemma pentaphyllum	ELISA
<u>Ribosome-inactivating protein karasurin-C Antibody</u>	CSB-PA326146HA01TIF	Trichosanthes kirilowii	ELISA
<u>Ribosome-inactivating protein lychnin Antibody</u>	CSB-PA307648HA01LKN	Silene chalcedonica	ELISA
<u>Ribosome-inactivating protein PD-L1/PD-L2 Antibody</u>	CSB-PA306294LA01PJJ	Phytolacca dioica	ELISA
<u>RIC1 Antibody</u>	CSB-PA671121LA01HU	Human	ELISA, IHC, IF
<u>RIDA Antibody</u>	CSB-PA010744LA01HU	Human	ELISA, IHC, IF
<u>RILP Antibody</u>	CSB-PA850399LA01HU	Human	ELISA, IHC
<u>RIMS2 Antibody</u>	CSB-PA892173LA01HU	Human	ELISA, IF
<u>RIN1 Antibody</u>	CSB-PA019726LA01HU	Human	ELISA, IF
<u>RIN3 Antibody</u>	CSB-PA019728LA01HU	Human	ELISA, IF
<u>RINT1 Antibody</u>	CSB-PA747351LA01HU	Human	ELISA, IHC
<u>RIOK2 Antibody</u>	CSB-PA887136LA01HU	Human	ELISA, IHC, IF
<u>RIOX2 Antibody</u>	CSB-PA810267LA01HU	Human	ELISA, IHC, IF
<u>Ripk1 Antibody</u>	CSB-PA720181LA01MO	Mouse, Human	ELISA, IP
<u>RIPK1 Antibody</u>	CSB-PA618785LA01HU	Human, Rat, Mouse	ELISA, WB, IF
<u>RIPK2 Antibody</u>	CSB-PA019736LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>RIPK3 Antibody</u>	CSB-PA897497LA01HU	Human	ELISA, WB, IHC, IF
<u>RIPOR2 Antibody</u>	CSB-PA008285LA01HU	Human	ELISA, IHC
<u>RIPPLY2 Antibody</u>	CSB-PA735979LA01HU	Human	ELISA, IHC
<u>RIT1 Antibody</u>	CSB-PA846116LA01HU	Human	ELISA, IHC
<u>RIT2 Antibody</u>	CSB-PA860769LA01HU	Human	ELISA, IHC, IF
<u>RLBP1 Antibody</u>	CSB-PA019743LA01HU	Human	ELISA, IHC, IF
<u>RMDN1 Antibody</u>	CSB-PA850277LA01HU	Human	ELISA, WB, IHC
<u>RMDN3 Antibody</u>	CSB-PA836307LA01HU	Human	ELISA, IHC
<u>RMI1 Antibody</u>	CSB-PA875699LA01HU	Human, Mouse	ELISA, WB
<u>RMND1 Antibody</u>	CSB-PA882123LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>RMND5A Antibody</u>	CSB-PA875693LA01HU	Human	ELISA, IHC
<u>RMND5B Antibody</u>	CSB-PA856939LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>RNASE2 Antibody</u>	CSB-PA019794LA01HU	Human	ELISA, IHC
<u>RNASE3 Antibody</u>	CSB-PA12789A0Rb	Human	ELISA
<u>RNASE9 Antibody</u>	CSB-PA019800LA01HU	Human	ELISA, IHC
<u>RNASEH1 Antibody</u>	CSB-PA019801LA01HU	Human	ELISA, WB, IHC
<u>RNASEH2A Antibody</u>	CSB-PA019805LA01HU	Human	ELISA, IHC
<u>RNASEH2C Antibody</u>	CSB-PA837441LA01HU	Human	ELISA, WB, IHC, IF
<u>RNASEL Antibody</u>	CSB-PA313223LA01HU	Human	ELISA, IHC
<u>RND1 Antibody</u>	CSB-PA835694LA01HU	Human	ELISA, IHC, IF
<u>RND2 Antibody</u>	CSB-PA019812LA01HU	Human, Mouse	ELISA, WB
<u>RND3 Antibody</u>	CSB-PA019813LA01HU	Human	ELISA, IHC
<u>RNF10 Antibody</u>	CSB-PA019814LA01HU	Human, Mouse	ELISA, WB, IHC
<u>RNF103 Antibody</u>	CSB-PA019815LA01HU	Human	ELISA, IHC, IF
<u>RNF125 Antibody</u>	CSB-PA836205LA01HU	Human	ELISA, IHC
<u>RNF126 Antibody</u>	CSB-PA019828LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>RNF128 Antibody</u>	CSB-PA837448LA01HU	Human, Mouse	ELISA, WB, IF
<u>RNF133 Antibody</u>	CSB-PA855485LA01HU	Human	ELISA, WB, IHC, IF
<u>RNF138 Antibody</u>	CSB-PA851971LA01HU	Human	ELISA, IHC
<u>RNF141 Antibody</u>	CSB-PA019839LA01HU	Human	ELISA
<u>RNF148 Antibody</u>	CSB-PA019844LA01HU	Human	ELISA, WB, IHC, IF
<u>RNF168 Antibody</u>	CSB-PA814236LA01HU	Human	ELISA, IHC
<u>RNF181 Antibody</u>	CSB-PA019860LA01HU	Human	ELISA, IHC
<u>RNF182 Antibody</u>	CSB-PA822714LA01HU	Human	ELISA, IHC, IF
<u>RNF183 Antibody</u>	CSB-PA019862LA01HU	Human	ELISA, IHC, IF
<u>RNF19A Antibody</u>	CSB-PA868327LA01HU	Human	ELISA, IHC
<u>RNF2 Antibody</u>	CSB-PA859101LA01HU	Human, Mouse	ELISA, WB
<u>RNF20 Antibody</u>	CSB-PA736198LA01HU	Human	ELISA, IHC
<u>RNF213 Antibody</u>	CSB-PA717453LA01HU	Human	ELISA, IHC, IF
<u>RNF219 Antibody</u>	CSB-PA019878LA01HU	Human	ELISA, WB, IHC
<u>RNF220 Antibody</u>	CSB-PA719619LA01HU	Human	ELISA, IHC
<u>RNF24 Antibody</u>	CSB-PA019881LA01HU	Human	ELISA, IHC, IF
<u>RNF26 Antibody</u>	CSB-PA019883LA01HU	Human	ELISA, WB, IHC, IF
<u>RNF31 Antibody</u>	CSB-PA856927LA01HU	Human	ELISA, IHC, IF
<u>RNF39 Antibody</u>	CSB-PA019888LA01HU	Human	ELISA, IHC, IF
<u>RNF4 Antibody</u>	CSB-PA019889LA01HU	Human	ELISA, IHC, IF
<u>RNF40 Antibody</u>	CSB-PA019890LA01HU	Human	ELISA, IHC
<u>RNF43 Antibody</u>	CSB-PA019892LA01HU	Human	ELISA, IHC
<u>RNF6 Antibody</u>	CSB-PA897080LA01HU	Human	ELISA, IF



Product Name	Code	Species Reactivity	Applications
<u>RNF7 Antibody</u>	CSB-PA019897DA01HU	Human	ELISA, IHC
<u>RNU2 Antibody</u>	CSB-PA360432LA01UBD	Ustilago sphaerogena	ELISA, WB
<u>ROBO1 Antibody</u>	CSB-PA896760LA01HU	Human	ELISA, WB, IF
<u>ROCK1 Antibody</u>	CSB-PA614415HA01HU	Human	ELISA, WB, IHC, IF
<u>ROCK2 Antibody</u>	CSB-PA020059DA01HU	Human, Mouse	ELISA, WB
<u>rodA Antibody</u>	CSB-PA340427LA01NGS	Neosartorya fumigata	ELISA
<u>ROGDI Antibody</u>	CSB-PA867108LA01HU	Human	ELISA, WB, IHC
<u>ROMO1 Antibody</u>	CSB-PA020063LA01HU	Human	ELISA, IF
<u>ROPN1 Antibody</u>	CSB-PA888012LA01HU	Human	ELISA, IHC
<u>ROPN1B Antibody</u>	CSB-PA020065EA01HU	Human	ELISA, IHC
<u>ROPN1L Antibody</u>	CSB-PA856908LA01HU	Human	ELISA, IHC
<u>ROR1 Antibody</u>	CSB-PA020067LA01HU	Human	ELISA, IHC
<u>ROR2 Antibody</u>	CSB-PA21119A0Rb	Human	ELISA, IHC, IF
<u>RORA Antibody</u>	CSB-PA020069LA01HU	Human, Mouse	ELISA, WB
<u>RORB Antibody</u>	CSB-PA821694LA01HU	Human	ELISA, WB
<u>RORC Antibody</u>	CSB-PA020071LA01HU	Human	ELISA, WB, IHC
<u>ROS1 Antibody</u>	CSB-PA874445LA01DOA	Arabidopsis thaliana	ELISA
<u>ROS1 Antibody</u>	CSB-PA020072LA01HU	Human	ELISA, IF
<u>RPA1 Antibody</u>	CSB-PA020088LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>RPA2 Antibody</u>	CSB-PA01555A0Rb	Human	ELISA, WB, IHC, IF, IP
<u>RPA3 Antibody</u>	CSB-PA16234A0Rb	Human	ELISA, WB
<u>RPA3 Antibody</u>	CSB-PA16239A0Rb	Human	ELISA, IHC
<u>RPA4 Antibody</u>	CSB-PA16249A0Rb	Human	ELISA
<u>RPE Antibody</u>	CSB-PA853385LA01HU	Human	ELISA, IHC, IF
<u>RPE65 Antibody</u>	CSB-PA624103LA01HU	Human, Rat	ELISA, WB, IF
<u>RPGR Antibody</u>	CSB-PA838814LA01HU	Human	ELISA, WB, IHC, IF
<u>RPH3AL Antibody</u>	CSB-PA891975LA01HU	Human, Rat	ELISA, WB
<u>rpiA Antibody</u>	CSB-PA540419DA01ENT	Escherichia coli	ELISA, WB
<u>RPIA Antibody</u>	CSB-PA020104LA01HU	Human	ELISA, WB
<u>RPL10A Antibody</u>	CSB-PA05249A0Rb	Human	ELISA, IF
<u>RPL10L Antibody</u>	CSB-PA850330LA01HU	Human, Mouse	ELISA, WB, IHC
<u>RPL11 Antibody</u>	CSB-PA00065A0Rb	Human	ELISA, IHC, IF
<u>RPL12 Antibody</u>	CSB-PA020116LA01HU	Human	ELISA, IHC
<u>RPL13A Antibody</u>	CSB-PA05289A0Rb	Human	ELISA, IF
<u>RPL14 Antibody</u>	CSB-PA00565A0Rb	Human	ELISA, IHC
<u>RPL15 Antibody</u>	CSB-PA02435A0Rb	Human	ELISA, WB
<u>RPL17 Antibody</u>	CSB-PA02504A0Rb	Human	ELISA, WB
<u>RPL17 Antibody</u>	CSB-PA02505A0Rb	Human	ELISA, WB
<u>RPL19 Antibody</u>	CSB-PA020160LA01HU	Human, Mouse	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>RPL22L1 Antibody</u>	CSB-PA740724LA01HU	Human	ELISA, IHC
<u>RPL23A Antibody</u>	CSB-PA020187LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>RPL24 Antibody</u>	CSB-PA020202LA01HU	Human	ELISA, WB
<u>RPL26 Antibody</u>	CSB-PA020207LA01HU	Human	ELISA, IHC, IF
<u>RPL26L1 Antibody</u>	CSB-PA892361LA01HU	Human	ELISA, IHC
<u>RPL27 Antibody</u>	CSB-PA02345A0Rb	Human	ELISA, WB, IHC, IF
<u>RPL28 Antibody</u>	CSB-PA020217EA01HU	Human	ELISA, IHC, IF
<u>RPL3 Antibody</u>	CSB-PA02385A0Rb	Human	ELISA, IHC
<u>RPL30 Antibody</u>	CSB-PA020225DA01HU	Human	ELISA, WB, IHC
<u>RPL34 Antibody</u>	CSB-PA020243LA01HU	Human	ELISA, IHC
<u>RPL35A Antibody</u>	CSB-PA020248LA01HU	Human	ELISA, IHC
<u>RPL36AL Antibody</u>	CSB-PA01804A0Rb	Human	ELISA, WB
<u>RPL36AL Antibody</u>	CSB-PA01805A0Rb	Human	ELISA, WB
<u>RPL38 Antibody</u>	CSB-PA02405A0Rb	Human	ELISA, WB, IF
<u>RPL39L Antibody</u>	CSB-PA850285LA01HU	Human	ELISA, IHC
<u>RPL4 Antibody</u>	CSB-PA020289LA01HU	Human	ELISA, WB, IHC, IF
<u>RPL5 Antibody</u>	CSB-PA02495A0Rb	Human	ELISA, WB, IHC
<u>RPL7A Antibody</u>	CSB-PA020303LA01HU	Human, Mouse	ELISA, WB, IHC
<u>RPL7L1 Antibody</u>	CSB-PA735330LA01HU	Human	ELISA, IHC
<u>RPL8 Antibody</u>	CSB-PA00595A0Rb	Human	ELISA, WB
<u>rplL Antibody</u>	CSB-PA354049HA01SKX	Staphylococcus aureus	ELISA, WB
<u>RPLP0 Antibody</u>	CSB-PA020336LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>RPLP1 Antibody</u>	CSB-PA02364A0Rb	Human	ELISA, WB, IHC, IF
<u>RPN1 Antibody</u>	CSB-PA05124A0Rb	Human	ELISA, WB, IHC
<u>RPN13 Antibody</u>	CSB-PA517092HA01SVG	Saccharomyces cerevisiae	ELISA
<u>RPP14 Antibody</u>	CSB-PA020346LA01HU	Human	ELISA, WB
<u>RPP21 Antibody</u>	CSB-PA887986EA01HU	Human	ELISA, WB
<u>RPP38 Antibody</u>	CSB-PA020350LA01HU	Human	ELISA, IHC
<u>RPRD1A Antibody</u>	CSB-PA850404LA01HU	Human	ELISA, WB
<u>RPRD1B Antibody</u>	CSB-PA878854HA01HU	Human	ELISA, WB, IHC
<u>RPS10 Antibody</u>	CSB-PA02565A0Rb	Human	ELISA, WB
<u>RPS10P5 Antibody</u>	CSB-PA878849LA01HU	Human	ELISA, WB, IP
<u>RPS11 Antibody</u>	CSB-PA01874A0Rb	Human	ELISA, WB, IHC
<u>RPS12 Antibody</u>	CSB-PA02525A0Rb	Human	ELISA, WB
<u>RPS13 Antibody</u>	CSB-PA02485A0Rb	Human	ELISA, WB, IHC
<u>RPS14 Antibody</u>	CSB-PA04255A0Rb	Human	ELISA



Product Name	Code	Species Reactivity	Applications
<u>RPS15 Antibody</u>	CSB-PA020372LA01HU	Human	ELISA, WB, IHC, IF
<u>RPS16 Antibody</u>	CSB-PA02465A0Rb	Human	ELISA, WB
<u>RPS18 Antibody</u>	CSB-PA01645A0Rb	Human	ELISA, WB
<u>RPS19 Antibody</u>	CSB-PA00655A0Rb	Human	ELISA, WB, IHC
<u>RPS2 Antibody</u>	CSB-PA020393LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>RPS21 Antibody</u>	CSB-PA020399LA01HU	Human	ELISA, IHC, IF
<u>RPS23 Antibody</u>	CSB-PA020400LA01HU	Human	ELISA, IF
<u>RPS24 Antibody</u>	CSB-PA020402LA01HU	Human	ELISA, IHC
<u>RPS25 Antibody</u>	CSB-PA00644A0Rb	Human	ELISA
<u>RPS25 Antibody</u>	CSB-PA00645A0Rb	Human	ELISA
<u>RPS26 Antibody</u>	CSB-PA020407LA01HU	Human	ELISA, IHC
<u>RPS27A Antibody</u>	CSB-PA01639A0Rb	Human	ELISA, IF
<u>RPS28 Antibody</u>	CSB-PA020430LA01HU	Human	ELISA, IHC
<u>RPS29 Antibody</u>	CSB-PA020432LA01HU	Human	ELISA, IF
<u>RPS4Y1 Antibody</u>	CSB-PA020457LA01HU	Human	ELISA, IHC, IF
<u>RPS5 Antibody</u>	CSB-PA00729A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>RPS6 Antibody</u>	CSB-PA020463LA01HU	Human	ELISA, IHC, IF
<u>RPS6KA1 Antibody</u>	CSB-PA618984LA01HU	Human	ELISA, IF
<u>RPS6KA2 Antibody</u>	CSB-PA618003LA01HU	Human	ELISA, IHC, IF
<u>RPS6KA3 Antibody</u>	CSB-PA020466LA01HU	Human, Rat	ELISA, WB, IHC
<u>RPS6KA4 Antibody</u>	CSB-PA020467LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>RPS6KA5 Antibody</u>	CSB-PA020468LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>RPS6KA6 Antibody</u>	CSB-PA891943LA01HU	Human	ELISA, IHC
<u>RPS6KB1 Antibody</u>	CSB-PA020470LA01HU	Human	ELISA, IF, IP
<u>RPS6KB2 Antibody</u>	CSB-PA02905A0Rb	Human	ELISA, WB, IHC
<u>RPS7 Antibody</u>	CSB-PA020476LA01HU	Human	ELISA, IHC
<u>RPS8 Antibody</u>	CSB-PA020478LA01HU	Human, Mouse	ELISA, WB
<u>RPS9 Antibody</u>	CSB-PA00695A0Rb	Human	ELISA, WB
<u>RPSA Antibody</u>	CSB-PA020485LA01HU	Human, Mouse	ELISA, WB
<u>rpsB Antibody</u>	CSB-PA08644A0Rb	Escherichia coli	ELISA, WB
<u>RPTOR Antibody</u>	CSB-PA839774LA01HU	Human	ELISA, IHC, IF
<u>RPUSD4 Antibody</u>	CSB-PA836192LA01HU	Human	ELISA, WB, IHC
<u>RRAGA Antibody</u>	CSB-PA03465A0Rb	Human	ELISA, IHC, IF
<u>RRAS Antibody</u>	CSB-PA020513PA01HU	Human	ELISA, IHC
<u>RRM1 Antibody</u>	CSB-PA020518LA01HU	Human	ELISA, WB
<u>RRM1 Antibody</u>	CSB-PA020518YA01HU	Human	ELISA, WB
<u>RRM2 Antibody</u>	CSB-PA020519HA01HU	Human	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>RS1 Antibody</u>	CSB-PA020534EA01HU	Human	ELISA
<u>RS1 Antibody</u>	CSB-PA020534LA01HU	Human	ELISA, IHC
<u>RSPH14 Antibody</u>	CSB-PA887042LA01HU	Human	ELISA, WB, IHC
<u>RSPH3 Antibody</u>	CSB-PA769775LA01HU	Human	ELISA
<u>RSPH9 Antibody</u>	CSB-PA867129LA01HU	Human, Mouse	ELISA, WB, IHC
<u>RSPO3 Antibody</u>	CSB-PA887154LA01HU	Human	ELISA, IHC
<u>rsrIM Antibody</u>	CSB-PA318583LA01RLE	Rhodobacter sphaeroides	ELISA, WB
<u>RSU1 Antibody</u>	CSB-PA020558LA01HU	Human	ELISA, WB
<u>RTBDN Antibody</u>	CSB-PA020560LA01HU	Human	ELISA
<u>RTCA Antibody</u>	CSB-PA520755LA01HU	Human	ELISA, IHC
<u>rtcB Antibody</u>	CSB-PA341647YA01EVN	Escherichia coli	ELISA
<u>RTCB Antibody</u>	CSB-PA897546LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>RTEL1 Antibody</u>	CSB-PA020563LA01HU	Human	ELISA, IHC, IF
<u>RTF1 Antibody</u>	CSB-PA821675HA01HU	Human	ELISA, WB, IHC, IF
<u>RTKN2 Antibody</u>	CSB-PA810311LA01HU	Human	ELISA, IHC
<u>RTL8A Antibody</u>	CSB-PA871589LA01HU	Human	ELISA
<u>RTL8B Antibody</u>	CSB-PA622040LA01HU	Human	ELISA, IHC
<u>RTL8C Antibody</u>	CSB-PA008044LA01HU	Human	ELISA, WB, IHC
<u>RTN4IP1 Antibody</u>	CSB-PA848832LA01HU	Human	ELISA, IHC
<u>RTN4RL1 Antibody</u>	CSB-PA803121LA01HU	Human	ELISA, IHC, IF
<u>RTN4RL2 Antibody</u>	CSB-PA774805NA01HU	Human	ELISA, WB, IHC, IF
<u>RTP2 Antibody</u>	CSB-PA728807LA01HU	Human, Mouse	ELISA, WB
<u>RTP3 Antibody</u>	CSB-PA887106LA01HU	Human, Mouse	ELISA, WB, IHC, IP
<u>RTP4 Antibody</u>	CSB-PA842659LA01HU	Human, Mouse	ELISA, WB, IHC
<u>RTTN Antibody</u>	CSB-PA020581LA01HU	Human	ELISA, WB
<u>RUBCN Antibody</u>	CSB-PA856415LA01HU	Human	ELISA, IHC
<u>RuBisCO activase Antibody</u>	CSB-PA00170E0Rb	Spinach	ELISA, WB
<u>Rubredoxin Antibody</u>	CSB-PA360387HA01CMA	Clostridium pasteurianum	ELISA, WB
<u>RUFY3 Antibody</u>	CSB-PA765974LA01HU	Human	ELISA, WB, IF
<u>RUFY4 Antibody</u>	CSB-PA751157LA01HU	Human	ELISA, IHC
<u>RUNDC3A Antibody</u>	CSB-PA685437LA01HU	Human	ELISA, WB, IHC
<u>RUNX1 Antibody</u>	CSB-PA020592LA01HU	Human	ELISA, WB, IF
<u>RUNX1T1 Antibody</u>	CSB-PA726765HA01HU	Human	ELISA, WB, IHC
<u>RUNX2 Antibody</u>	CSB-PA14439A0Rb	Human	ELISA, WB
<u>RUNX3 Antibody</u>	CSB-PA020595LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>RUSC1 Antibody</u>	CSB-PA871585LA01HU	Human	ELISA, IHC, IF
<u>ruvA Antibody</u>	CSB-PA359062HA01ENV	Escherichia coli	ELISA



Product Name	Code	Species Reactivity	Applications
<u>RUVBL1 Antibody</u>	CSB-PA896483LA01HU	Human	ELISA, WB, IHC, IF
<u>RUVBL2 Antibody</u>	CSB-PA897459LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>ruvC Antibody</u>	CSB-PA364095HA01ENV	Escherichia coli	ELISA
<u>RWDD2A Antibody</u>	CSB-PA883450LA01HU	Human, Mouse	ELISA, WB, IHC
<u>RXFP1 Antibody</u>	CSB-PA875715LA01HU	Human	ELISA, IF
<u>RXFP2 Antibody</u>	CSB-PA819908LA01HU	Human	ELISA, WB, IHC, IF
<u>RXRG Antibody</u>	CSB-PA020614LA01HU	Human	ELISA, IHC
<u>RYK Antibody</u>	CSB-PA020617LA01HU	Human	ELISA, WB, IHC
<u>RYR1 Antibody</u>	CSB-PA020619LA01HU	Human	ELISA, IF
<u>RYR2 Antibody</u>	CSB-PA15139A0Rb	Human	ELISA, IHC, IF
<u>RYR3 Antibody</u>	CSB-PA020621LA01HU	Human	ELISA, IHC, IF
<u>S Antibody</u>	CSB-PA333052LA01BJO	Bovine coronavirus	ELISA
<u>S Antibody</u>	CSB-PA336163EA01HIY	Human coronavirus OC43	ELISA
<u>S Antibody</u>	CSB-PA322803LA01BJK	Bovine coronavirus	ELISA, WB
<u>S100A1 Antibody</u>	CSB-PA020622HA01HU	Human	ELISA, IHC, IF
<u>S100A10 Antibody</u>	CSB-PA020623LA01HU	Human	ELISA, IHC, IF
<u>S100A11 Antibody</u>	CSB-PA020624HA01HU	Human	ELISA, IHC
<u>S100A12 Antibody</u>	CSB-PA020626HA01HU	Human	ELISA, IHC
<u>S100A16 Antibody</u>	CSB-PA853418LA01HU	Human	ELISA, IHC, IF
<u>S100A2 Antibody</u>	CSB-PA020630LA01HU	Human	ELISA
<u>S100A4 Antibody</u>	CSB-PA020632HA01HU	Human	ELISA, WB, IHC
<u>S100A6 Antibody</u>	CSB-PA020634HA01HU	Human	ELISA, IHC
<u>S100A7 Antibody</u>	CSB-PA020635HA01HU	Human	ELISA, IHC, IF
<u>S100a8 Antibody</u>	CSB-PA020641DA01MO	Mouse	ELISA, WB
<u>S100A8 Antibody</u>	CSB-PA020641HA01HU	Human	ELISA, IHC
<u>S100a8 Antibody</u>	CSB-PA020641HA01MO	Mouse	ELISA
<u>S100a8 Antibody</u>	CSB-PA020641LA01RA	Rat	ELISA, WB
<u>S100A9 Antibody</u>	CSB-PA020642HA01HU	Human	ELISA, IHC, IF
<u>S100a9 Antibody</u>	CSB-PA020642HA01MO	Mouse	ELISA, WB, IP
<u>S100A9 Antibody</u>	CSB-PA020642LA01HU	Human	ELISA, IF
<u>S100B Antibody</u>	CSB-PA10677A0Rb	Human	ELISA
<u>S100g Antibody</u>	CSB-PA020644LA01RA	Rat	ELISA
<u>S100G Antibody</u>	CSB-PA020644LA01HU	Human	ELISA, IHC, IF
<u>S100P Antibody</u>	CSB-PA020645LA01HU	Human	ELISA, IHC, IF
<u>S100PBP Antibody</u>	CSB-PA822185LA01HU	Human	ELISA, IHC
<u>S1PR1 Antibody</u>	CSB-PA020650LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>S1PR2 Antibody</u>	CSB-PA020651LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>S1PR3 Antibody</u>	CSB-PA859932LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>S1PR5 Antibody</u>	CSB-PA867132LA01HU	Human, Mouse	ELISA, WB, IF
<u>SAA1 Antibody</u>	CSB-PA020656LA01DO	Dog	ELISA
<u>SAA1 Antibody</u>	CSB-PA14047A0Rb	Human	ELISA, WB, IHC
<u>SAA2 Antibody</u>	CSB-PA020657LA01HU	Human	ELISA, IHC
<u>sacIM Antibody</u>	CSB-PA523234LA01FNS	Streptomyces achromogenes	ELISA
<u>SACM1L Antibody</u>	CSB-PA868317LA01HU	Human	ELISA, IHC
<u>SACS Antibody</u>	CSB-PA873707LA01HU	Human	ELISA, IHC, IF
<u>SAG Antibody</u>	CSB-PA020669YA01HU	Human	ELISA
<u>SAGE1 Antibody</u>	CSB-PA020670LA01HU	Human	ELISA, IHC
<u>SAMD12 Antibody</u>	CSB-PA847657LA01HU	Human	ELISA, IHC
<u>SAMD3 Antibody</u>	CSB-PA843281LA01HU	Human	ELISA, IHC
<u>SAMD7 Antibody</u>	CSB-PA801801LA01HU	Human	ELISA, IHC
<u>SAMD8 Antibody</u>	CSB-PA822254LA01HU	Human	ELISA, IHC
<u>Samhd1 Antibody</u>	CSB-PA020691LA01MO	Mouse, Human	ELISA, WB
<u>SAMSN1 Antibody</u>	CSB-PA885722LA01HU	Human	ELISA, IHC
<u>SAP9 Antibody</u>	CSB-PA522117LA01CZF	Yeast	ELISA, WB
<u>SARAF Antibody</u>	CSB-PA853392HA01HU	Human	ELISA
<u>SARS Antibody</u>	CSB-PA04145A0Rb	Human	ELISA, WB, IHC
<u>SART3 Antibody</u>	CSB-PA624010HA01HU	Human	ELISA, IHC
<u>SASH1 Antibody</u>	CSB-PA020714LA01HU	Human	ELISA, IHC
<u>SAT1 Antibody</u>	CSB-PA02754A0Rb	Human	ELISA, WB, IHC
<u>SAT1 Antibody</u>	CSB-PA02759A0Rb	Human	ELISA, IHC, IF
<u>SAT2 Antibody</u>	CSB-PA836208LA01HU	Human	ELISA, IHC, IF
<u>SATB2 Antibody</u>	CSB-PA892170LA01HU	Human	ELISA, IHC
<u>SBDS Antibody</u>	CSB-PA897481LA01HU	Human	ELISA, IHC, IF
<u>SBF2 Antibody</u>	CSB-PA773047LA01HU	Human	ELISA, IHC, IF
<u>SBK2 Antibody</u>	CSB-PA020731LA01HU	Human	ELISA, IHC, IF
<u>SBNO2 Antibody</u>	CSB-PA896873LA01HU	Human	ELISA, IF
<u>scaIM Antibody</u>	CSB-PA528741LA01FNZ	Streptomyces caespitosus	ELISA
<u>SCAMP1 Antibody</u>	CSB-PA020750LA01HU	Human, Rat	ELISA, WB, IHC
<u>SCAMP2 Antibody</u>	CSB-PA020751LA01HU	Human	ELISA, IHC, IF
<u>SCAMP3 Antibody</u>	CSB-PA05195A0Rb	Human	ELISA, WB, IHC, IF
<u>SCAMP4 Antibody</u>	CSB-PA856883LA01HU	Human	ELISA, IHC, IF
<u>SCAMP5 Antibody</u>	CSB-PA855032LA01HU	Human	ELISA, IHC, IF
<u>SCAND2P Antibody</u>	CSB-PA861977LA01HU	Human	ELISA
<u>SCAPER Antibody</u>	CSB-PA863150LA01HU	Human	ELISA, IHC
<u>SCARA3 Antibody</u>	CSB-PA721489LA01HU	Human, Mouse	ELISA, WB, IHC
<u>SCARA5 Antibody</u>	CSB-PA020766LA01HU	Human	ELISA, IHC, IF
<u>SCARB1 Antibody</u>	CSB-PA845139LA01HU	Human	ELISA, IHC
<u>SCARB2 Antibody</u>	CSB-PA619859LA01HU	Human, Mouse	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>SCD Antibody</u>	CSB-PA020802LA01HU	Human	ELISA, IHC
<u>SCFD2 Antibody</u>	CSB-PA840978LA01HU	Human	ELISA, IHC
<u>SCG2 Antibody</u>	CSB-PA020808LA01HU	Human	ELISA, IHC, IF
<u>SCG3 Antibody</u>	CSB-PA845166LA01HU	Human	ELISA, IF
<u>SCGB1A1 Antibody</u>	CSB-PA020811LA01HU	Human	ELISA, IF
<u>SCGB1C1 Antibody</u>	CSB-PA837436LA01HU	Human	ELISA
<u>SCGB2A1 Antibody</u>	CSB-PA020816LA01HU	Human	ELISA
<u>Scgb3a2 Antibody</u>	CSB-PA846028LA01MO	Mouse	ELISA
<u>SCGN Antibody</u>	CSB-PA020821LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SCIMP Antibody</u>	CSB-PA754405LA01HU	Human	ELISA, IHC, IF
<u>SCIN Antibody</u>	CSB-PA020823LA01HU	Human, Mouse	ELISA, WB, IF
<u>SCLT1 Antibody</u>	CSB-PA020826LA01HU	Human	ELISA, IHC
<u>SCLY Antibody</u>	CSB-PA020827LA01HU	Human	ELISA, IHC
<u>SCMH1 Antibody</u>	CSB-PA846625LA01HU	Human	ELISA
<u>SCML2 Antibody</u>	CSB-PA891470LA01HU	Human	ELISA, IHC
<u>scn Antibody</u>	CSB-PA858805HA01SKY	Staphylococcus aureus	ELISA
<u>SCN10A Antibody</u>	CSB-PA750910LA01HU	Human	ELISA, WB, IF
<u>SCN11A Antibody</u>	CSB-PA883417LA01HU	Human	ELISA, IHC, IF
<u>SCN1B Antibody</u>	CSB-PA020835LA01HU	Human	ELISA, IHC, IF
<u>SCN4B Antibody</u>	CSB-PA818246LA01HU	Human	ELISA, IF
<u>SCN9A Antibody</u>	CSB-PA613605LA01HU	Human	ELISA, IHC, IF
<u>SCNM1 Antibody</u>	CSB-PA020847LA01HU	Human	ELISA, IHC
<u>SCNN1B Antibody</u>	CSB-PA020849LA01HU	Human	ELISA, WB, IHC, IF
<u>SCNN1G Antibody</u>	CSB-PA020851LA01HU	Human	ELISA, IHC, IF
<u>SCO1 Antibody</u>	CSB-PA020852LA01HU	Human	ELISA, IHC
<u>SCOC Antibody</u>	CSB-PA020854LA01HU	Human	ELISA, IHC, IF
<u>SCP Antibody</u>	CSB-PA321237HA01EFA	Epstein-Barr virus	ELISA
<u>SCP2 Antibody</u>	CSB-PA020856DA01HU	Human	ELISA, WB, IHC
<u>SCPEP1 Antibody</u>	CSB-PA867189LA01HU	Human, Mouse	ELISA, WB, IHC
<u>Serg1 Antibody</u>	CSB-PA526023HA01MO	Mouse	ELISA
<u>SCRG1 Antibody</u>	CSB-PA530448LA01HU	Human	ELISA, IHC
<u>SCRN3 Antibody</u>	CSB-PA607167LA01HU	Human	ELISA
<u>SCTR Antibody</u>	CSB-PA020865LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>SCUBE1 Antibody</u>	CSB-PA811616LA01HU	Human	ELISA, IF
<u>SCUBE2 Antibody</u>	CSB-PA865101LA01HU	Human, Rat	ELISA, WB, IHC
<u>SCYL3 Antibody</u>	CSB-PA809026LA01HU	Human	ELISA, IHC
<u>Scytovirin Antibody</u>	CSB-PA308458LA01SEH	Scytonema varium	ELISA
<u>SDC1 Antibody</u>	CSB-PA020888HA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>SDC4 Antibody</u>	CSB-PA327046HA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SDCBP Antibody</u>	CSB-PA020893LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>SDCBP2 Antibody</u>	CSB-PA872430LA01HU	Human	ELISA, IHC
<u>sdeA Antibody</u>	CSB-PA28219A0Rb	Legionella pneumophila subsp	ELISA, WB
<u>SDF2 Antibody</u>	CSB-PA020900LA01HU	Human	ELISA, IHC
<u>SDF4 Antibody</u>	CSB-PA861120LA01HU	Human	ELISA, WB
<u>SDHA Antibody</u>	CSB-PA01985A0Rb	Human	ELISA, WB, IHC, IF
<u>SDHA Antibody</u>	CSB-PA01989A0Rb	Human, Mouse, Rat, Zebrafish	ELISA, WB, IHC
<u>SDHAF1 Antibody</u>	CSB-PA406800LA01HU	Human	ELISA, IHC
<u>SDHAF2 Antibody</u>	CSB-PA885749LA01HU	Human	ELISA, IHC, IF
<u>SDHC Antibody</u>	CSB-PA859106LA01HU	Human	ELISA, IHC, IF
<u>SDR Antibody</u>	CSB-PA24399A0Rb	Zebrafish	ELISA, WB
<u>SDSL Antibody</u>	CSB-PA020927LA01HU	Human, Mouse	ELISA, WB, IHC
<u>SEC11A Antibody</u>	CSB-PA020931LA01HU	Human	ELISA, IHC
<u>SEC11C Antibody</u>	CSB-PA871609LA01HU	Human	ELISA, IHC
<u>SEC14L2 Antibody</u>	CSB-PA020937LA01HU	Human	ELISA, WB, IHC, IF
<u>SEC14L3 Antibody</u>	CSB-PA887012LA01HU	Human, Mouse	ELISA, WB, IHC
<u>SEC14L4 Antibody</u>	CSB-PA883374LA01HU	Human	ELISA, IHC
<u>SEC16B Antibody</u>	CSB-PA853439LA01HU	Human	ELISA, IHC
<u>SEC22A Antibody</u>	CSB-PA856954LA01HU	Human	ELISA, IHC, IF
<u>SEC22B Antibody</u>	CSB-PA020945LA01HU	Human	ELISA
<u>SEC23A Antibody</u>	CSB-PA618985LA01HU	Human	ELISA, IF
<u>SEC23B Antibody</u>	CSB-PA621874LA01HU	Human	ELISA, IHC
<u>SEC23IP Antibody</u>	CSB-PA896767LA01HU	Human	ELISA, WB, IHC, IF
<u>SEC24C Antibody</u>	CSB-PA020952LA01HU	Human	ELISA
<u>SEC31A Antibody</u>	CSB-PA020954LA01HU	Human	ELISA, IHC
<u>SEC61G Antibody</u>	CSB-PA03509A0Rb	Human	ELISA, IHC
<u>SECTM1 Antibody</u>	CSB-PA819898LA01HU	Human	ELISA, IHC, IF
<u>SEL1L Antibody</u>	CSB-PA866208LA01HU	Human	ELISA, IHC, IF
<u>SELE Antibody</u>	CSB-PA020975EA01PI	Pig	ELISA
<u>SELENOM Antibody</u>	CSB-PA837868LA01HU	Human	ELISA
<u>SELENON Antibody</u>	CSB-PA878943LA01HU	Human, Rat	ELISA, WB
<u>SELENOP Antibody</u>	CSB-PA021018LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>Selenot Antibody</u>	CSB-PA350436LA01MO	Mouse	ELISA



Product Name	Code	Species Reactivity	Applications
<u>SELENOW Antibody</u>	CSB-PA021034LA01HU	Human	ELISA, IHC
<u>SELP Antibody</u>	CSB-PA020978LA01HU	Human	ELISA, IF
<u>SEM1 Antibody</u>	CSB-PA021263LA01HU	Human	ELISA, IHC
<u>SEMA3A Antibody</u>	CSB-PA613506LA01HU	Human, Rat, Mouse	ELISA, WB, IHC
<u>SEMA3C Antibody</u>	CSB-PA860786LA01HU	Human	ELISA, IHC
<u>SEMA3F Antibody</u>	CSB-PA622647LA01HU	Human	ELISA, IHC
<u>SEMA4B Antibody</u>	CSB-PA868280LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>SEMA4D Antibody</u>	CSB-PA835707LA01HU	Human	ELISA, IF
<u>SEMA5A Antibody</u>	CSB-PA613413LA01HU	Human	ELISA, IHC, IF
<u>SEMA5B Antibody</u>	CSB-PA868389LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SEMA6A Antibody</u>	CSB-PA875655YA01HU	Human	ELISA, IHC, IF
<u>SEMA6B Antibody</u>	CSB-PA862010LA01HU	Human	ELISA, IHC, IF
<u>SEMA6D Antibody</u>	CSB-PA847710LA01HU	Human	ELISA, WB, IF
<u>SENP1 Antibody</u>	CSB-PA885787LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SENP2 Antibody</u>	CSB-PA888026LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SENP6 Antibody</u>	CSB-PA867110LA01HU	Human	ELISA, IF
<u>SENP7 Antibody</u>	CSB-PA805340LA01HU	Human	ELISA, IHC
<u>SENP8 Antibody</u>	CSB-PA822251HA01HU	Human	ELISA
<u>SEPSECS Antibody</u>	CSB-PA14389A0Rb	Human	ELISA, WB, IHC, IF
<u>SEPT11 Antibody</u>	CSB-PA878895LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>SEPT12 Antibody</u>	CSB-PA809019LA01HU	Human, Mouse	ELISA, WB, IF
<u>SEPT2 Antibody</u>	CSB-PA617994LA01HU	Human	ELISA, IHC
<u>SEPT3 Antibody</u>	CSB-PA021026LA01HU	Human, Mouse	ELISA, WB, IHC
<u>serB Antibody</u>	CSB-PA018938HA01ENV	Escherichia coli	ELISA
<u>SERF1A Antibody</u>	CSB-PA021040LA01HU	Human	ELISA, IHC
<u>SERHL2 Antibody</u>	CSB-PA862016LA01HU	Human	ELISA
<u>SERINC1 Antibody</u>	CSB-PA021046LA01HU	Human, Mouse	ELISA, WB, IF
<u>Serine protease inhibitor Kazal-type 2-like protein Antibody</u>	CSB-PA19037A0Rb	Cyprinus carpio	ELISA
<u>Serpina1 Antibody</u>	CSB-PA021053HA01RA	Rat	ELISA
<u>SERPINA1 Antibody</u>	CSB-PA12565A0Rb	Human	ELISA
<u>SERPINA11 Antibody</u>	CSB-PA801228LA01HU	Human	ELISA, WB
<u>Serpina3n Antibody</u>	CSB-PA835578LA01MO	Mouse	ELISA
<u>SERPINA5 Antibody</u>	CSB-PA021061LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>SERPINA6 Antibody</u>	CSB-PA021062EA01HU	Human	ELISA, WB, IHC, IF
<u>SERPINB1 Antibody</u>	CSB-PA021065LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>SERPINB12 Antibody</u>	CSB-PA822276LA01HU	Human	ELISA, IHC, IF
<u>SERPINB14 Antibody</u>	CSB-PA004301E01CH	Chicken	ELISA, WB
<u>Serpib2 Antibody</u>	CSB-PA021070HA01RA	Rat, Human	ELISA, WB
<u>SERPINB2 Antibody</u>	CSB-PA021070LA01HU	Human	ELISA, WB, IHC
<u>SERPINB3 Antibody</u>	CSB-PA11275A0Rb	Human	ELISA, IHC, IF
<u>SERPINB4 Antibody</u>	CSB-PA021072LA01HU	Human	ELISA, IHC, IF
<u>SERPINB8 Antibody</u>	CSB-PA021076EA01HU	Human	ELISA, WB, IHC
<u>SERPINC1 Antibody</u>	CSB-PA021079LA01HU	Human, Rat	ELISA, WB, IHC
<u>SERPINE1 Antibody</u>	CSB-PA11239A0Rb	Human	ELISA, WB, IHC, IF
<u>Serpinf1 Antibody</u>	CSB-PA021084DA01RA	Rat	ELISA
<u>SERPINF1 Antibody</u>	CSB-PA021084LA01HU	Human	ELISA, IHC
<u>SERPING1 Antibody</u>	CSB-PA021086LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>Serpinh1 Antibody</u>	CSB-PA021087LA01MO	Mouse, Human	ELISA, WB
<u>SERPINI1 Antibody</u>	CSB-PA859934LA01HU	Human	ELISA, IHC
<u>SERTAD1 Antibody</u>	CSB-PA883406LA01HU	Human	ELISA, WB, IHC
<u>SES1 Antibody</u>	CSB-PA357079LA01SVG	Saccharomyces cerevisiae	ELISA
<u>SESN1 Antibody</u>	CSB-PA021095LA01HU	Human, Mouse	ELISA, WB, IHC
<u>SESN2 Antibody</u>	CSB-PA021096NA01HU	Human	ELISA, IF
<u>SET Antibody</u>	CSB-PA875345LA01HU	Human	ELISA, WB
<u>SETBP1 Antibody</u>	CSB-PA896935LA01HU	Human	ELISA
<u>SETD1A Antibody</u>	CSB-PA021101HA01HU	Human	ELISA, IHC
<u>SETD3 Antibody</u>	CSB-PA769772LA01HU	Human	ELISA, IHC
<u>SETD5 Antibody</u>	CSB-PA871633LA01HU	Human	ELISA, IHC, IF
<u>SETD6 Antibody</u>	CSB-PA855042LA01HU	Human	ELISA, IHC
<u>SETD7 Antibody</u>	CSB-PA840975LA01HU	Human	ELISA, WB, IHC
<u>SETD9 Antibody</u>	CSB-PA836727LA01HU	Human	ELISA, WB
<u>SETX Antibody</u>	CSB-PA800097LA01HU	Human	ELISA, IHC
<u>SEZ6L2 Antibody</u>	CSB-PA761499LA01HU	Human	ELISA, IF
<u>SF3B1 Antibody</u>	CSB-PA021125LA01HU	Human	ELISA, WB, IHC
<u>SF3B6 Antibody</u>	CSB-PA897543LA01HU	Human	ELISA, IHC
<u>SFN Antibody</u>	CSB-PA021135LA01HU	Human, Rat, Mouse	ELISA, WB, IHC, IF
<u>SFPQ Antibody</u>	CSB-PA021136LA01HU	Human, Mouse	ELISA, WB
<u>SFR1 Antibody</u>	CSB-PA801256LA01HU	Human	ELISA, IHC
<u>SFRP1 Antibody</u>	CSB-PA021138LA01HU	Human, Mouse	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>SFRP2 Antibody</u>	CSB-PA021139LA01HU	Human	ELISA, IF
<u>SFRP4 Antibody</u>	CSB-PA021140LA01HU	Human	ELISA, WB, IHC
<u>SFRP5 Antibody</u>	CSB-PA021141LA01HU	Human	ELISA, IHC
<u>SFT2D2 Antibody</u>	CSB-PA021164LA01HU	Human, Mouse	ELISA, WB, IF
<u>SFTA3 Antibody</u>	CSB-PA021167LA01HU	Human	ELISA
<u>SFTPB Antibody</u>	CSB-PA18055A0Rb	Pig	ELISA
<u>SFTPC Antibody</u>	CSB-PA021174LA01BO	Bovine	ELISA
<u>SFTPC Antibody</u>	CSB-PA021174LA01HU	Human	ELISA, IHC
<u>SFXN2 Antibody</u>	CSB-PA853470LA01HU	Human	ELISA, IHC, IF
<u>SFXN4 Antibody</u>	CSB-PA744046LA01HU	Human	ELISA, WB, IHC, IF
<u>SFXN5 Antibody</u>	CSB-PA819464LA01HU	Human	ELISA
<u>SGCA Antibody</u>	CSB-PA623086LA01HU	Human	ELISA, WB, IHC, IF
<u>SGCB Antibody</u>	CSB-PA618082LA01HU	Human	ELISA, IHC, IF
<u>SGCE Antibody</u>	CSB-PA021185LA01HU	Human	ELISA, IHC
<u>SGCG Antibody</u>	CSB-PA618771LA01HU	Human	ELISA, IHC, IF
<u>SGF29 Antibody</u>	CSB-PA856928LA01HU	Human	ELISA, IHC
<u>SGK1 Antibody</u>	CSB-PA021189LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SGMS1 Antibody</u>	CSB-PA801243LA01HU	Human	ELISA, WB, IHC
<u>SGMS2 Antibody</u>	CSB-PA854159LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>SGO1 Antibody</u>	CSB-PA688764LA01HU	Human, Mouse	ELISA, WB, IF
<u>SGTB Antibody</u>	CSB-PA846614LA01HU	Human	ELISA, IHC
<u>SH2D1A Antibody</u>	CSB-PA021209LA01HU	Human	ELISA, IHC, IF
<u>SH2D1B Antibody</u>	CSB-PA021210LA01HU	Human	ELISA, IHC, IF
<u>SH2D3A Antibody</u>	CSB-PA880088OA01HU	Human	ELISA, IHC
<u>SH3BGR Antibody</u>	CSB-PA021219LA01HU	Human	ELISA
<u>SH3BGRL Antibody</u>	CSB-PA021220LA01HU	Human	ELISA, IHC
<u>SH3BGRL3 Antibody</u>	CSB-PA880952LA01HU	Human	ELISA, WB, IHC
<u>SH3BP1 Antibody</u>	CSB-PA23069A0Rb	Human	ELISA, IHC
<u>SH3BP2 Antibody</u>	CSB-PA021224LA01HU	Human	ELISA, IHC
<u>SH3BP5 Antibody</u>	CSB-PA021226LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SH3D19 Antibody</u>	CSB-PA711356LA01HU	Human	ELISA, IHC
<u>SH3GLB1 Antibody</u>	CSB-PA878084LA01HU	Human	ELISA, IF
<u>SH3GLB2 Antibody</u>	CSB-PA868296LA01HU	Human, Mouse	ELISA, WB, IHC
<u>SH3RF1 Antibody</u>	CSB-PA773800LA01HU	Human	ELISA, IF
<u>SH3TC2 Antibody</u>	CSB-PA819482LA01HU	Human	ELISA, IHC, IF
<u>SHANK2 Antibody</u>	CSB-PA892492LA01HU	Human	ELISA, IHC
<u>SHARPIN Antibody</u>	CSB-PA867119LA01HU	Human	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>SHBG Antibody</u>	CSB-PA12289A0Rb	Human, Mouse	ELISA, WB
<u>SHC1 Antibody</u>	CSB-PA021253LA01HU	Human	ELISA, IHC, IF
<u>SHC3 Antibody</u>	CSB-PA838779LA01HU	Human, Mouse	ELISA, WB
<u>SHCBP1 Antibody</u>	CSB-PA847699LA01HU	Human	ELISA
<u>SHD Antibody</u>	CSB-PA850314LA01HU	Human	ELISA, IHC
<u>SHE Antibody</u>	CSB-PA730142LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>SHF Antibody</u>	CSB-PA755509LA01HU	Human	ELISA, IF
<u>SHH Antibody</u>	CSB-PA623000LA01HU	Human	ELISA, IF
<u>shha Antibody</u>	CSB-PA835612LA01DIL	Zebrafish	ELISA
<u>SHISA5 Antibody</u>	CSB-PA818677LA01HU	Human	ELISA, IHC, IF, IP
<u>SHISA9 Antibody</u>	CSB-PA457545LA01HU	Human	ELISA, IF
<u>SHOC2 Antibody</u>	CSB-PA891468LA01HU	Human	ELISA, WB, IHC, IF
<u>SHPRH Antibody</u>	CSB-PA619905LA01HU	Human	ELISA, IHC, IF
<u>SHQ1 Antibody</u>	CSB-PA744134LA01HU	Human	ELISA, IHC
<u>SHROOM2 Antibody</u>	CSB-PA622684LA01HU	Human	ELISA, IHC, IF
<u>SHROOM3 Antibody</u>	CSB-PA844740LA01HU	Human	ELISA, IHC, IF
<u>SIAH2 Antibody</u>	CSB-PA021288LA01HU	Human	ELISA, IHC, IF
<u>Sialic acid-binding lectin Antibody</u>	CSB-PA321644LA01RJL	Rana japonica	ELISA
<u>SIGLEC10 Antibody</u>	CSB-PA842706LA01HU	Human	ELISA, IHC, IF
<u>SIGLEC15 Antibody</u>	CSB-PA761623LA01HU	Human, Mouse	ELISA, WB
<u>SIGLEC6 Antibody</u>	CSB-PA021300LA01HU	Human	ELISA, WB, IHC
<u>SIGLEC7 Antibody</u>	CSB-PA896485LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SIGMAR1 Antibody</u>	CSB-PA021320LA01HU	Human	ELISA, IHC
<u>SIK1 Antibody</u>	CSB-PA021321LA01HU	Human	ELISA, IHC, IF
<u>SIL1 Antibody</u>	CSB-PA875645LA01HU	Human, Mouse	ELISA, WB
<u>SIM2 Antibody</u>	CSB-PA621763LA01HU	Human	ELISA, IHC
<u>SIMC1 Antibody</u>	CSB-PA822768LA01HU	Human	ELISA
<u>SIN3A Antibody</u>	CSB-PA836303LA01HU	Human	ELISA, WB, IHC
<u>SIN3B Antibody</u>	CSB-PA021328LA01HU	Human	ELISA, IHC, IF
<u>SIPA1 Antibody</u>	CSB-PA021330LA01HU	Human	ELISA, WB, IHC, IF
<u>SIPA1L2 Antibody</u>	CSB-PA882179LA01HU	Human	ELISA, IHC
<u>SIRPA Antibody</u>	CSB-PA021334LA01HU	Human	ELISA, IHC, IF
<u>SIRPB1 Antibody</u>	CSB-PA021335LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SIRPB1 Antibody</u>	CSB-PA732939LA01HU	Human	ELISA, IHC, IF
<u>SIRPD Antibody</u>	CSB-PA021337LA01HU	Human, Mouse	ELISA, WB
<u>SIRT1 Antibody</u>	CSB-PA822202HA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>SIRT2 Antibody</u>	CSB-PA812889LA01HU	Human	ELISA, IHC, IF
<u>SIRT3 Antibody</u>	CSB-PA882102LA01HU	Human, Mouse	ELISA, WB, IF
<u>SIRT5 Antibody</u>	CSB-PA873680LA01HU	Human	ELISA, IHC, IF
<u>SIRT7 Antibody</u>	CSB-PA885703LA01HU	Human	ELISA, IHC, IF
<u>SIVA1 Antibody</u>	CSB-PA021347LA01HU	Human	ELISA, IHC
<u>SIX1 Antibody</u>	CSB-PA618986LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>SKA1 Antibody</u>	CSB-PA846589LA01HU	Human	ELISA, IHC
<u>SKA2 Antibody</u>	CSB-PA823890LA01HU	Human	ELISA, IHC
<u>SKA3 Antibody</u>	CSB-PA816899LA01HU	Human	ELISA, IHC, IF
<u>SKAP1 Antibody</u>	CSB-PA769799LA01HU	Human	ELISA
<u>skib Antibody</u>	CSB-PA24359A0Rb	Zebrafish	ELISA, WB
<u>SKIL Antibody</u>	CSB-PA021357LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>SKOR1 Antibody</u>	CSB-PA308064HA01HU	Human	ELISA, IHC
<u>SKOR2 Antibody</u>	CSB-PA009018PA01HU	Human	ELISA, IHC, IF
<u>SKP1 Antibody</u>	CSB-PA021360LA01HU	Human, Mouse	ELISA, WB
<u>SKP2 Antibody</u>	CSB-PA613392LA01HU	Human	ELISA, WB, IHC
<u>SLA Antibody</u>	CSB-PA613387LA01HU	Human	ELISA, WB, IHC, IF
<u>SLAIN1 Antibody</u>	CSB-PA850864LA01HU	Human	ELISA, WB, IHC
<u>SLAMF1 Antibody</u>	CSB-PA614403LA01HU	Human	ELISA, WB, IHC, IF
<u>SLAMF6 Antibody</u>	CSB-PA021368LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SLAMF8 Antibody</u>	CSB-PA021370LA01HU	Human	ELISA, IHC
<u>SLC10A2 Antibody</u>	CSB-PA622519LA01HU	Human	ELISA, IHC, IF
<u>SLC10A3 Antibody</u>	CSB-PA021375LA01HU	Human, Rat	ELISA, WB, IHC
<u>SLC10A6 Antibody</u>	CSB-PA665974LA01HU	Human	ELISA
<u>SLC11A1 Antibody</u>	CSB-PA021380LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>SLC11A2 Antibody</u>	CSB-PA021381LA01HU	Human	ELISA, WB, IHC, IF
<u>SLC12A2 Antibody</u>	CSB-PA021383LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SLC12A3 Antibody</u>	CSB-PA021384LA01HU	Human	ELISA, IF
<u>SLC12A4 Antibody</u>	CSB-PA892363LA01HU	Human	ELISA, IHC
<u>SLC12A5 Antibody</u>	CSB-PA863965LA01HU	Human	ELISA, IHC, IF
<u>SLC12A7 Antibody</u>	CSB-PA893907LA01HU	Human	ELISA, IHC, IF
<u>SLC12A9 Antibody</u>	CSB-PA866307LA01HU	Human	ELISA, IF
<u>SLC13A4 Antibody</u>	CSB-PA891949LA01HU	Human	ELISA, WB, IHC, IF
<u>SLC13A5 Antibody</u>	CSB-PA768239LA01HU	Human, Mouse, Rat	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>SLC14A1 Antibody</u>	CSB-PA615678LA01HU	Human	ELISA, IHC, IF
<u>SLC15A3 Antibody</u>	CSB-PA815578LA01HU	Human, Mouse	ELISA, WB, IHC
<u>SLC15A4 Antibody</u>	CSB-PA843279LA01HU	Human	ELISA, IHC, IF
<u>SLC16A10 Antibody</u>	CSB-PA851563HA01HU	Human	ELISA, IHC, IF
<u>SLC16A11 Antibody</u>	CSB-PA843323LA01HU	Human	ELISA, WB, IF
<u>SLC16A12 Antibody</u>	CSB-PA757957LA01HU	Human	ELISA, WB, IHC, IF
<u>SLC16A3 Antibody</u>	CSB-PA021410LA01HU	Human, Mouse	ELISA, WB, IF
<u>SLC16A7 Antibody</u>	CSB-PA021414LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>SLC16A8 Antibody</u>	CSB-PA021415LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SLC16A9 Antibody</u>	CSB-PA742402LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SLC17A3 Antibody</u>	CSB-PA021419LA01HU	Human	ELISA, IHC, IF
<u>SLC17A4 Antibody</u>	CSB-PA897087LA01HU	Human	ELISA, IF
<u>SLC18A1 Antibody</u>	CSB-PA021426LA01HU	Human, Mouse	ELISA, WB, IHC
<u>SLC18A2 Antibody</u>	CSB-PA021427LA01HU	Human, Mouse	ELISA, WB, IF
<u>SLC19A2 Antibody</u>	CSB-PA021430LA01HU	Human	ELISA, IHC
<u>SLC19A3 Antibody</u>	CSB-PA866333LA01HU	Human	ELISA, WB, IF
<u>SLC1A2 Antibody</u>	CSB-PA021433LA01HU	Human, Rat, Mouse	ELISA, WB, IHC, IF
<u>SLC1A6 Antibody</u>	CSB-PA021437LA01HU	Human	ELISA, WB, IHC, IF
<u>SLC1A7 Antibody</u>	CSB-PA021438EA01HU	Human	ELISA, IHC, IF
<u>SLC20A1 Antibody</u>	CSB-PA823882LA01HU	Human	ELISA, WB, IHC, IF
<u>SLC22A1 Antibody</u>	CSB-PA021442LA01HU	Human	ELISA, IHC
<u>SLC22A17 Antibody</u>	CSB-PA823881LA01HU	Human	ELISA, WB, IHC, IF
<u>SLC22A18AS Antibody</u>	CSB-PA847610LA01HU	Human	ELISA, IHC
<u>SLC22A2 Antibody</u>	CSB-PA021453LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>SLC22A23 Antibody</u>	CSB-PA021455LA01HU	Human, Mouse	ELISA, WB, IHC
<u>SLC22A24 Antibody</u>	CSB-PA850799LA01HU	Human	ELISA, IHC, IF
<u>SLC22A25 Antibody</u>	CSB-PA740875LA01HU	Human	ELISA, WB, IHC, IF
<u>SLC22A5 Antibody</u>	CSB-PA021460LA01HU	Human	ELISA, WB, IF
<u>SLC22A9 Antibody</u>	CSB-PA811607LA01HU	Human	ELISA, IHC, IF
<u>SLC23A1 Antibody</u>	CSB-PA887038LA01HU	Human, Mouse	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>SLC23A3 Antibody</u>	CSB-PA761358LA01HU	Human	ELISA, IHC, IF
<u>SLC24A5 Antibody</u>	CSB-PA754654LA01HU	Human	ELISA, WB
<u>SLC25A1 Antibody</u>	CSB-PA021475LA01HU	Human	ELISA, IHC, IF
<u>SLC25A10 Antibody</u>	CSB-PA883371LA01HU	Human, Mouse	ELISA, WB, IHC
<u>SLC25A12 Antibody</u>	CSB-PA021478LA01HU	Human, Rat, Mouse	ELISA, WB, IHC
<u>SLC25A14 Antibody</u>	CSB-PA021480LA01HU	Human	ELISA, IHC, IF
<u>SLC25A18 Antibody</u>	CSB-PA880945LA01HU	Human	ELISA, WB, IHC, IF
<u>SLC25A23 Antibody</u>	CSB-PA871580LA01HU	Human, Mouse	ELISA, WB
<u>SLC25A24 Antibody</u>	CSB-PA743993LA01HU	Human, Rat, Mouse	ELISA, WB, IHC, IF
<u>SLC25A25 Antibody</u>	CSB-PA760850LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SLC25A26 Antibody</u>	CSB-PA761690LA01HU	Human	ELISA, IHC, IF
<u>SLC25A32 Antibody</u>	CSB-PA872477LA01HU	Human, Mouse	ELISA, WB
<u>SLC25A33 Antibody</u>	CSB-PA874800LA01HU	Human	ELISA, IHC, IF
<u>SLC25A35 Antibody</u>	CSB-PA668593LA01HU	Human	ELISA, IHC
<u>SLC25A36 Antibody</u>	CSB-PA853398LA01HU	Human	ELISA, WB, IHC, IF
<u>SLC25A37 Antibody</u>	CSB-PA873700LA01HU	Human	ELISA, IHC
<u>SLC25A38 Antibody</u>	CSB-PA846607LA01HU	Human	ELISA, IF
<u>SLC25A4 Antibody</u>	CSB-PA021510LA01HU	Human	ELISA, IHC, IF
<u>SLC25A40 Antibody</u>	CSB-PA840561LA01HU	Human	ELISA, IHC, IF
<u>SLC25A46 Antibody</u>	CSB-PA856895LA01HU	Human, Mouse	ELISA, WB, IF
<u>SLC25A5 Antibody</u>	CSB-PA021518LA01HU	Human	ELISA, WB, IHC
<u>SLC25A51 Antibody</u>	CSB-PA875649LA01HU	Human, Mouse	ELISA, WB, IF
<u>SLC25A52 Antibody</u>	CSB-PA662246LA01HU	Human	ELISA
<u>SLC25A6 Antibody</u>	CSB-PA021520LA01HU	Human	ELISA, WB, IHC
<u>SLC26A2 Antibody</u>	CSB-PA021525LA01HU	Human	ELISA, IHC, IF
<u>SLC26A4 Antibody</u>	CSB-PA021527LA01HU	Human	ELISA, WB, IHC, IF
<u>SLC26A6 Antibody</u>	CSB-PA861164LA01HU	Human	ELISA, IHC, IF
<u>SLC26A8 Antibody</u>	CSB-PA853496LA01HU	Human	ELISA, IHC, IF
<u>SLC26A9 Antibody</u>	CSB-PA741974LA01HU	Human	ELISA, IHC, IF
<u>SLC27A1 Antibody</u>	CSB-PA740768LA01HU	Human	ELISA, IF
<u>SLC27A2 Antibody</u>	CSB-PA021534LA01HU	Human	ELISA
<u>SLC28A1 Antibody</u>	CSB-PA021539LA01HU	Human	ELISA, IHC, IF
<u>SLC28A3 Antibody</u>	CSB-PA884504LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SLC29A1 Antibody</u>	CSB-PA857875LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SLC2A2 Antibody</u>	CSB-PA13329A0Rb	Human	ELISA, IF



Product Name	Code	Species Reactivity	Applications
<u>SLC2A3 Antibody</u>	CSB-PA021553LA01HU	Human	ELISA, IHC, IF
<u>SLC2A4 Antibody</u>	CSB-PA021557EA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SLC2A4RG Antibody</u>	CSB-PA021558LA01HU	Human	ELISA, IHC, IF
<u>SLC2A9 Antibody</u>	CSB-PA873638LA01HU	Human, Mouse	ELISA, WB, IHC
<u>SLC30A10 Antibody</u>	CSB-PA747695LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SLC30A3 Antibody</u>	CSB-PA859940LA01HU	Human	ELISA, IHC
<u>SLC30A6 Antibody</u>	CSB-PA744010LA01HU	Human	ELISA, WB, IHC, IF
<u>SLC30A7 Antibody</u>	CSB-PA836732LA01HU	Human, Mouse	ELISA, WB, IHC
<u>SLC30A8 Antibody</u>	CSB-PA818247LA01HU	Human	ELISA, IHC, IF
<u>SLC31A2 Antibody</u>	CSB-PA021577LA01HU	Human	ELISA, IHC, IF
<u>SLC32A1 Antibody</u>	CSB-PA863984LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SLC33A1 Antibody</u>	CSB-PA13339A0Rb	Human, Mouse	ELISA, WB, IHC
<u>SLC34A2 Antibody</u>	CSB-PA021581LA01HU	Human	ELISA, IHC
<u>SLC34A3 Antibody</u>	CSB-PA854015LA01HU	Human	ELISA, IHC, IF
<u>SLC35A2 Antibody</u>	CSB-PA021584LA01HU	Human, Mouse	ELISA, WB, IF
<u>SLC35B1 Antibody</u>	CSB-PA021588LA01HU	Human	ELISA, WB, IF
<u>SLC35B4 Antibody</u>	CSB-PA853373LA01HU	Human	ELISA, IF
<u>SLC35C1 Antibody</u>	CSB-PA839285LA01HU	Human	ELISA, IHC
<u>SLC35D2 Antibody</u>	CSB-PA765728LA01HU	Human	ELISA, IHC
<u>SLC35D3 Antibody</u>	CSB-PA021596LA01HU	Human	ELISA, IHC, IF
<u>SLC35F2 Antibody</u>	CSB-PA809010LA01HU	Human	ELISA, IHC, IF
<u>SLC35F5 Antibody</u>	CSB-PA840985LA01HU	Human	ELISA, WB, IHC, IF
<u>SLC35F6 Antibody</u>	CSB-PA843146LA01HU	Human	ELISA, IHC, IF
<u>SLC36A2 Antibody</u>	CSB-PA672077LA01HU	Human	ELISA, IHC, IF
<u>SLC36A3 Antibody</u>	CSB-PA674326LA01HU	Human	ELISA, IHC, IF
<u>SLC37A1 Antibody</u>	CSB-PA021610LA01HU	Human	ELISA, IHC, IF
<u>SLC37A2 Antibody</u>	CSB-PA819478LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SLC38A2 Antibody</u>	CSB-PA847210LA01HU	Human	ELISA, WB, IHC, IF
<u>SLC38A5 Antibody</u>	CSB-PA837853LA01HU	Human	ELISA, IF
<u>SLC38A6 Antibody</u>	CSB-PA818274LA01HU	Human, Mouse	ELISA, WB, IF
<u>SLC38A7 Antibody</u>	CSB-PA889118HA01HU	Human	ELISA, IHC, IF
<u>SLC38A9 Antibody</u>	CSB-PA836713LA01HU	Human	ELISA, WB, IHC, IF
<u>SLC39A1 Antibody</u>	CSB-PA17954A0Rb	Human	ELISA
<u>SLC39A10 Antibody</u>	CSB-PA891736LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>SLC39A11 Antibody</u>	CSB-PA854023LA01HU	Human	ELISA, IF
<u>Slc39a13 Antibody</u>	CSB-PA644804LA01RA	Rat	ELISA
<u>SLC39A8 Antibody</u>	CSB-PA021637LA01HU	Human	ELISA, IHC, IF
<u>SLC3A1 Antibody</u>	CSB-PA737352LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>SLC3A2 Antibody</u>	CSB-PA021640LA01HU	Human	ELISA, WB
<u>SLC40A1 Antibody</u>	CSB-PA865089LA01HU	Human	ELISA, IF
<u>SLC41A1 Antibody</u>	CSB-PA810274LA01HU	Human	ELISA, IHC, IF
<u>SLC41A2 Antibody</u>	CSB-PA021643LA01HU	Human	ELISA, IHC, IF
<u>SLC43A3 Antibody</u>	CSB-PA854092LA01HU	Human	ELISA, IHC, IF
<u>SLC44A4 Antibody</u>	CSB-PA681506LA01HU	Human	ELISA, IF
<u>SLC45A2 Antibody</u>	CSB-PA890941LA01HU	Human	ELISA, WB
<u>SLC45A3 Antibody</u>	CSB-PA836239LA01HU	Human	ELISA, IHC, IF
<u>SLC46A2 Antibody</u>	CSB-PA871607LA01HU	Human, Mouse	ELISA, WB, IF
<u>SLC46A3 Antibody</u>	CSB-PA773784LA01HU	Human	ELISA, WB, IHC, IF
<u>SLC4A5 Antibody</u>	CSB-PA874837LA01HU	Human	ELISA, WB, IF
<u>SLC4A8 Antibody</u>	CSB-PA649721LA01HU	Human	ELISA, WB, IHC, IF
<u>SLC4A9 Antibody</u>	CSB-PA842743LA01HU	Human	ELISA, IHC, IF
<u>SLC51A Antibody</u>	CSB-PA768212LA01HU	Human	ELISA, IHC, IF
<u>SLC51B Antibody</u>	CSB-PA771440LA01HU	Human	ELISA, IF
<u>SLC52A1 Antibody</u>	CSB-PA878904LA01HU	Human	ELISA, IHC, IF
<u>SLC52A2 Antibody</u>	CSB-PA881008LA01HU	Human	ELISA, IHC, IF
<u>SLC52A3 Antibody</u>	CSB-PA873618LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SLC5A1 Antibody</u>	CSB-PA021675LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SLC5A10 Antibody</u>	CSB-PA021676LA01HU	Human	ELISA, IHC, IF
<u>SLC5A2 Antibody</u>	CSB-PA021679LA01HU	Human	ELISA, IHC, IF
<u>SLC5A4 Antibody</u>	CSB-PA23589A0Rb	Human	ELISA, IF
<u>SLC5A5 Antibody</u>	CSB-PA835713LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>SLC5A8 Antibody</u>	CSB-PA836676LA01HU	Human	ELISA, IF
<u>SLC6A13 Antibody</u>	CSB-PA889104LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SLC6A15 Antibody</u>	CSB-PA875657LA01HU	Human	ELISA, IHC, IF
<u>SLC6A16 Antibody</u>	CSB-PA863929LA01HU	Human	ELISA, IHC, IF
<u>SLC6A19 Antibody</u>	CSB-PA731656LA01HU	Human	ELISA, IHC, IF
<u>SLC6A3 Antibody</u>	CSB-PA021701LA01HU	Human, Mouse	ELISA, WB, IF
<u>SLC6A6 Antibody</u>	CSB-PA021704LA01HU	Human	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>SLC6A7 Antibody</u>	CSB-PA04209A0Rb	Human	ELISA, IHC
<u>SLC6A8 Antibody</u>	CSB-PA04679A0Rb	Human, Mouse	ELISA, WB, IF
<u>SLC7A1 Antibody</u>	CSB-PA021709LA01HU	Human	ELISA, IHC, IF
<u>SLC7A14 Antibody</u>	CSB-PA844711LA01HU	Human	ELISA, IHC
<u>SLC7A2 Antibody</u>	CSB-PA021714LA01HU	Human	ELISA, IHC
<u>SLC7A3 Antibody</u>	CSB-PA855505LA01HU	Human	ELISA, IF
<u>SLC7A6 Antibody</u>	CSB-PA849761LA01HU	Human, Mouse	ELISA, WB, IHC
<u>SLC7A7 Antibody</u>	CSB-PA892349LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SLC9A1 Antibody</u>	CSB-PA021726LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SLC9A3R2 Antibody</u>	CSB-PA01759A0Rb	Human	ELISA, WB, IHC
<u>SLCO1A2 Antibody</u>	CSB-PA021741LA01HU	Human	ELISA, IF
<u>SLCO1B1 Antibody</u>	CSB-PA896932LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>SLCO1B3 Antibody</u>	CSB-PA882065LA01HU	Human	ELISA, IHC, IF
<u>SLCO1C1 Antibody</u>	CSB-PA868355LA01HU	Human	ELISA, IHC, IF
<u>SLCO2A1 Antibody</u>	CSB-PA856440LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>SLCO3A1 Antibody</u>	CSB-PA883432LA01HU	Human	ELISA, IHC, IF
<u>SLFN11 Antibody</u>	CSB-PA021760LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>Sln9 Antibody</u>	CSB-PA106533LA21MO	Mouse	ELISA, WB
<u>SLIT3 Antibody</u>	CSB-PA021769LA01HU	Human	ELISA, IHC
<u>SLITRK3 Antibody</u>	CSB-PA021772LA01HU	Human	ELISA
<u>SLITRK4 Antibody</u>	CSB-PA815560LA01HU	Human	ELISA, IHC
<u>SLITRK5 Antibody</u>	CSB-PA021774LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SLITRK6 Antibody</u>	CSB-PA867155LA01HU	Human	ELISA, IHC, IF
<u>SLPI Antibody</u>	CSB-PA07154A0Rb	Human	ELISA, IHC, IF
<u>SLU7 Antibody</u>	CSB-PA021783LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SLURP1 Antibody</u>	CSB-PA021784LA01HU	Human	ELISA, IHC, IF
<u>SMAD1 Antibody</u>	CSB-PA618998LA01HU	Human	ELISA, WB, IHC, IF
<u>SMAD1 Antibody</u>	CSB-PA618998HA01HU	Human	ELISA, IHC, IF
<u>SMAD2 Antibody</u>	CSB-PA15557A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>SMAD3 Antibody</u>	CSB-PA021788LA01HU	Human	ELISA, WB
<u>SMAD4 Antibody</u>	CSB-PA619768LA01HU	Human	ELISA, WB, IHC, IP
<u>SMAD5 Antibody</u>	CSB-PA859108LA01HU	Human, Mouse	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>SMAD7 Antibody</u>	CSB-PA021793LA01HU	Human	ELISA, WB, IF
<u>SMAP1 Antibody</u>	CSB-PA811625LA01HU	Human	ELISA, IHC
<u>SMARCA1 Antibody</u>	CSB-PA021798LA01HU	Human	ELISA, IHC
<u>SMARCA2 Antibody</u>	CSB-PA021799LA01HU	Human	ELISA, IHC, ChIP
<u>SMARCA4 Antibody</u>	CSB-PA021801LA01HU	Human	ELISA, IHC, IF
<u>SMARCAL1 Antibody</u>	CSB-PA885774LA01HU	Human, Mouse	ELISA, WB, IF
<u>SMARCB1 Antibody</u>	CSB-PA623654LA01HU	Human	ELISA, IHC
<u>SMARCD1 Antibody</u>	CSB-PA846627LA01HU	Human	ELISA, WB, IHC, IF
<u>SMARCD2 Antibody</u>	CSB-PA846114LA01HU	Human, Mouse	ELISA, WB, IHC
<u>SMARCD3 Antibody</u>	CSB-PA754358LA01HU	Human, Mouse	ELISA, WB
<u>SMARCE1 Antibody</u>	CSB-PA839275LA01HU	Human	ELISA, IHC
<u>SMC1B Antibody</u>	CSB-PA854104LA01HU	Human	ELISA, IHC
<u>SMC2 Antibody</u>	CSB-PA021814LA01HU	Human	ELISA, IF
<u>SMC4 Antibody</u>	CSB-PA885725YA01HU	Human	ELISA, IHC, IF
<u>SMC5 Antibody</u>	CSB-PA810294LA01HU	Human	ELISA, WB
<u>SMC6 Antibody</u>	CSB-PA836300LA01HU	Human	ELISA, IHC
<u>SMCO1 Antibody</u>	CSB-PA623935LA01HU	Human	ELISA, IHC
<u>SMCO3 Antibody</u>	CSB-PA003101LA01HU	Human	ELISA, IF
<u>SMCP Antibody</u>	CSB-PA021820DA01HU	Human	ELISA, IHC
<u>SMG6 Antibody</u>	CSB-PA773032LA01HU	Human	ELISA, IHC, IF
<u>SMG7 Antibody</u>	CSB-PA842128LA01HU	Human	ELISA, IHC, IF
<u>SMIM1 Antibody</u>	CSB-PA450720LA01HU	Human	ELISA, IHC, IF
<u>SMIM14 Antibody</u>	CSB-PA857012LA01HU	Human	ELISA, IHC, IF
<u>SMIM19 Antibody</u>	CSB-PA856922LA01HU	Human	ELISA, IHC, IF
<u>SMN1 Antibody</u>	CSB-PA614982LA01HU	Human	ELISA, WB, IHC, IF
<u>SMNDC1 Antibody</u>	CSB-PA021840LA01HU	Human	ELISA, IHC, IP
<u>SMO Antibody</u>	CSB-PA857876LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SMOC1 Antibody</u>	CSB-PA875673LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>Smox Antibody</u>	CSB-PA021844LA01MO	Mouse	ELISA
<u>Smox Antibody</u>	CSB-PA021844YA01MO	Mouse, Human	ELISA, WB
<u>SMOX Antibody</u>	CSB-PA021844LA01HU	Human	ELISA, IHC, IF
<u>SMPD2 Antibody</u>	CSB-PA021846LA01HU	Human	ELISA, IHC, IF
<u>SMPD3 Antibody</u>	CSB-PA878920LA01HU	Human, Mouse	ELISA, IHC, IP
<u>SMPDL3A Antibody</u>	CSB-PA842119LA01HU	Human	ELISA, IHC
<u>SMS Antibody</u>	CSB-PA021854LA01HU	Human	ELISA, IHC, IF
<u>SMT3 Antibody</u>	CSB-PA615489HA01SVG	Saccharomyces cerevisiae	ELISA
<u>SMT3 Antibody</u>	CSB-PA615489LA01SVG	Saccharomyces cerevisiae	ELISA
<u>SMUG1 Antibody</u>	CSB-PA706636LA01HU	Human	ELISA



Product Name	Code	Species Reactivity	Applications
<u>SMURF1 Antibody</u>	CSB-PA875720LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SMURF2 Antibody</u>	CSB-PA872532EA01HU	Human	ELISA, IHC
<u>SMYD2 Antibody</u>	CSB-PA885705LA01HU	Human	ELISA, IHC, IF
<u>Snaclec rhodocytin subunit alpha Antibody</u>	CSB-PA888242LA01CBG	Calloselasma rhodostoma	ELISA, WB
<u>SNAI1 Antibody</u>	CSB-PA021867LA01HU	Human	ELISA, IHC
<u>SNAI2 Antibody</u>	CSB-PA021869LA01HU	Human	ELISA, IHC, IF
<u>SNAP23 Antibody</u>	CSB-PA021871LA01HU	Human	ELISA, IHC, IF
<u>SNAPC1 Antibody</u>	CSB-PA623083LA01HU	Human	ELISA
<u>Snca Antibody</u>	CSB-PA021912DA01RA	Rat	ELISA
<u>SNCA Antibody</u>	CSB-PA06994A0Rb	Human	ELISA, WB, IHC
<u>SNCA Antibody</u>	CSB-PA06997A0Rb	Human	ELISA, IHC, IF
<u>SNCA Antibody</u>	CSB-PA021912OA01HU	Human	ELISA, IHC, IF
<u>SNCAIP Antibody</u>	CSB-PA896546LA01HU	Human	ELISA, IHC, IF
<u>SNCG Antibody</u>	CSB-PA021915HA01HU	Human, Mouse	ELISA, WB, IHC
<u>SND1 Antibody</u>	CSB-PA745330LA01HU	Human	ELISA, IHC, IF
<u>SNED1 Antibody</u>	CSB-PA840582LA01HU	Human	ELISA, IHC
<u>SNF8 Antibody</u>	CSB-PA856945LA01HU	Human	ELISA, IF
<u>SNPH Antibody</u>	CSB-PA022317LA01HU	Human	ELISA
<u>SNRK Antibody</u>	CSB-PA868301NA01HU	Human	ELISA, WB
<u>SNRNP25 Antibody</u>	CSB-PA863129LA01HU	Human	ELISA
<u>Snrnp35 Antibody</u>	CSB-PA883986LA01MO	Mouse, Human	ELISA, WB
<u>SNRNP40 Antibody</u>	CSB-PA822196LA01HU	Human	ELISA
<u>SNRNP48 Antibody</u>	CSB-PA750321LA01HU	Human	ELISA
<u>SNRPA Antibody</u>	CSB-PA00974A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>SNRPB2 Antibody</u>	CSB-PA362460LA01HU	Human	ELISA, IHC
<u>SNRPD1 Antibody</u>	CSB-PA022332LA01HU	Human	ELISA, WB, IHC, IF
<u>SNRPD2 Antibody</u>	CSB-PA02635A0Rb	Human	ELISA, WB
<u>SNRPE Antibody</u>	CSB-PA022336LA01HU	Human	ELISA, WB, IHC, IF
<u>SNRPG Antibody</u>	CSB-PA01535A0Rb	Human	ELISA, WB, IHC, IF
<u>SNTB1 Antibody</u>	CSB-PA621699LA01HU	Human	ELISA, IF
<u>SNX10 Antibody</u>	CSB-PA896538LA01HU	Human	ELISA, IHC
<u>SNX12 Antibody</u>	CSB-PA891787LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>SNX14 Antibody</u>	CSB-PA806844LA01HU	Human	ELISA, IHC, IF
<u>SNX15 Antibody</u>	CSB-PA882091LA01HU	Human	ELISA, IHC, IF
<u>SNX16 Antibody</u>	CSB-PA022360LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>SNX17 Antibody</u>	CSB-PA614877LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SNX20 Antibody</u>	CSB-PA022366EA01HU	Human, Mouse	ELISA, WB, IHC
<u>SNX22 Antibody</u>	CSB-PA04059A0Rb	Human	ELISA, IHC
<u>SNX24 Antibody</u>	CSB-PA022369HA01HU	Human	ELISA, IHC
<u>SNX27 Antibody</u>	CSB-PA842704LA01HU	Human	ELISA, WB, IHC
<u>SNX3 Antibody</u>	CSB-PA01324A0Rb	Human	ELISA
<u>SNX3 Antibody</u>	CSB-PA01325A0Rb	Human	ELISA
<u>SNX33 Antibody</u>	CSB-PA851969LA01HU	Human	ELISA, IHC
<u>SNX5 Antibody</u>	CSB-PA872015LA01HU	Human	ELISA, WB, IHC, IF
<u>SNX6 Antibody</u>	CSB-PA890766LA01HU	Human	ELISA, WB, IHC
<u>SNX9 Antibody</u>	CSB-PA896917LA01HU	Human	ELISA, WB, IHC, IF
<u>SOCS1 Antibody</u>	CSB-PA022388LA01HU	Human	ELISA, IHC
<u>SOCS2 Antibody</u>	CSB-PA022389EA01HU	Human	ELISA, IF
<u>SOCS3 Antibody</u>	CSB-PA022392EA01HU	Human	ELISA, IHC, IF
<u>Socs3 Antibody</u>	CSB-PA022392LA01RA	Rat, Human	ELISA, WB
<u>SOCS7 Antibody</u>	CSB-PA022396LA01HU	Human, Rat	ELISA, WB, IF
<u>Sod1 Antibody</u>	CSB-PA022397LA01MO	Mouse, Rat, Human	ELISA, WB, IHC, IP
<u>Sod1 Antibody</u>	CSB-PA022397YA01MO	Mouse	ELISA
<u>SOD1 Antibody</u>	CSB-PA02864A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>SOD2 Antibody</u>	CSB-PA022398LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF, IP
<u>SOD3 Antibody</u>	CSB-PA022399LA01HU	Human	ELISA, IHC
<u>sodB Antibody</u>	CSB-PA852873LA01NGS	Neosartorya fumigata	ELISA
<u>SOHLH2 Antibody</u>	CSB-PA878909LA01HU	Human	ELISA, IHC
<u>SORBS3 Antibody</u>	CSB-PA022406LA01HU	Human	ELISA, WB, IHC
<u>SORCS1 Antibody</u>	CSB-PA848843LA01HU	Human	ELISA, WB, IHC, IF
<u>SORCS2 Antibody</u>	CSB-PA822282LA01HU	Human	ELISA, IHC, IF
<u>SORD Antibody</u>	CSB-PA022410HA01HU	Human	ELISA, IHC
<u>SORL1 Antibody</u>	CSB-PA821688LA01HU	Human	ELISA, IHC
<u>Sort1 Antibody</u>	CSB-PA022412LA01RA	Rat	ELISA
<u>SORT1 Antibody</u>	CSB-PA859102LA01HU	Human	ELISA, IHC
<u>SOS1 Antibody</u>	CSB-PA726788LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>SOS2 Antibody</u>	CSB-PA22327A0Rb	Human	ELISA, WB, IHC, IF
<u>SOSTDC1 Antibody</u>	CSB-PA022416EA01HU	Human	ELISA



Product Name	Code	Species Reactivity	Applications
<u>sox10 Antibody</u>	CSB-PA20859A0Rb	Zebrafish	ELISA, WB
<u>SOX10 Antibody</u>	CSB-PA022418OA01HU	Human	ELISA, IHC
<u>SOX11 Antibody</u>	CSB-PA022419LA01HU	Human, Rat, Mouse	ELISA, WB, IHC
<u>SOX14 Antibody</u>	CSB-PA022422LA01HU	Human	ELISA, WB, IHC, IF
<u>SOX17 Antibody</u>	CSB-PA887989LA01HU	Human	ELISA, IHC
<u>SOX2 Antibody</u>	CSB-PA16539A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>SOX21 Antibody</u>	CSB-PA896749LA01HU	Human	ELISA, IF
<u>SOX3 Antibody</u>	CSB-PA022429LA01HU	Human	ELISA, WB, IHC
<u>SOX30 Antibody</u>	CSB-PA022430LA01HU	Human	ELISA, IHC, IF
<u>SOX4 Antibody</u>	CSB-PA022431LA01HU	Human	ELISA, WB, IHC
<u>SOX8 Antibody</u>	CSB-PA022436LA01HU	Human	ELISA, WB
<u>SOX9 Antibody</u>	CSB-PA022437LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SP110 Antibody</u>	CSB-PA022440LA01HU	Human	ELISA, WB, IHC, IF
<u>SP3 Antibody</u>	CSB-PA22589A0Rb	Human	ELISA, IHC
<u>SP6 Antibody</u>	CSB-PA668784LA01HU	Human	ELISA, IHC
<u>SPA17 Antibody</u>	CSB-PA12124A0Rb	Human	ELISA, IHC
<u>SPACA1 Antibody</u>	CSB-PA864023LA01HU	Human	ELISA, IHC
<u>SPACA9 Antibody</u>	CSB-PA839312LA01HU	Human	ELISA
<u>SPAG6 Antibody</u>	CSB-PA022467LA01HU	Human	ELISA, IF
<u>SPAG8 Antibody</u>	CSB-PA858730LA01HU	Human	ELISA, IF
<u>SPAG9 Antibody</u>	CSB-PA022470LA01HU	Human	ELISA, IHC, IF
<u>SPAM1 Antibody</u>	CSB-PA022471LA01HU	Human	ELISA
<u>SPANXC Antibody</u>	CSB-PA878921LA01HU	Human	ELISA, IHC
<u>SPANXD Antibody</u>	CSB-PA819444LA01HU	Human	ELISA, IHC
<u>SPATA16 Antibody</u>	CSB-PA871597LA01HU	Human	ELISA, IHC
<u>SPATA19 Antibody</u>	CSB-PA770355LA01HU	Human, Mouse	ELISA, WB
<u>SPATA46 Antibody</u>	CSB-PA716346LA01HU	Human	ELISA, IHC
<u>SPATA8 Antibody</u>	CSB-PA764872LA01HU	Human	ELISA
<u>SPATC1L Antibody</u>	CSB-PA867117LA01HU	Human	ELISA, IHC, IF
<u>SPATS1 Antibody</u>	CSB-PA676576LA01HU	Human	ELISA, IHC
<u>SPATS2L Antibody</u>	CSB-PA868324HA01HU	Human	ELISA, WB
<u>SPC24 Antibody</u>	CSB-PA836712LA01HU	Human	ELISA, IHC
<u>SPC25 Antibody</u>	CSB-PA881019LA01HU	Human	ELISA, IHC
<u>SPCS1 Antibody</u>	CSB-PA896543LA01HU	Human	ELISA, IF
<u>SPCS2 Antibody</u>	CSB-PA617993LA01HU	Human	ELISA, WB, IHC, IF
<u>SPCS3 Antibody</u>	CSB-PA022517LA01HU	Human, Mouse	ELISA, WB, IHC
<u>SPDEF Antibody</u>	CSB-PA022518LA01HU	Human, Mouse	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>SPDL1 Antibody</u>	CSB-PA842662LA01HU	Human	ELISA, WB, IHC
<u>speA Antibody</u>	CSB-PA14507A0Rb	Escherichia coli	ELISA
<u>speB Antibody</u>	CSB-PA676488LA01SBAF	Streptococcus pyogenes serotype M28	ELISA, WB
<u>speB Antibody</u>	CSB-PA314886LA01SMR	Streptococcus pyogenes	ELISA, WB
<u>SPEM1 Antibody</u>	CSB-PA822704LA01HU	Human	ELISA, IHC
<u>SPERT Antibody</u>	CSB-PA843308LA01HU	Human, Mouse, Rat	ELISA, WB
<u>SPG11 Antibody</u>	CSB-PA842692LA01HU	Human	ELISA, IF
<u>SPG21 Antibody</u>	CSB-PA889158LA01HU	Human	ELISA, IHC, IF
<u>SPHK1 Antibody</u>	CSB-PA022564LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SPI1 Antibody</u>	CSB-PA16619A0Rb	Human	ELISA
<u>SPIC Antibody</u>	CSB-PA818708LA01HU	Human, Mouse	ELISA, WB, IHC
<u>SPIN2A Antibody</u>	CSB-PA859525LA01HU	Human	ELISA
<u>SPIN3 Antibody</u>	CSB-PA708489LA01HU	Human	ELISA, WB
<u>SPIN4 Antibody</u>	CSB-PA704087LA01HU	Human	ELISA, IHC, IF
<u>SPINK1 Antibody</u>	CSB-PA022575LA01HU	Human	ELISA, IHC, IF
<u>SPINK5 Antibody</u>	CSB-PA889077LA01HU	Human	ELISA, IHC
<u>SPINK6 Antibody</u>	CSB-PA744263LA01HU	Human	ELISA, IHC
<u>SPINK9 Antibody</u>	CSB-PA707717LA01HU	Human	ELISA, IHC
<u>SPINT2 Antibody</u>	CSB-PA022585LA01HU	Human, Mouse	ELISA, WB, IF
<u>SPIRE2 Antibody</u>	CSB-PA855490LA01HU	Human	ELISA, IHC, IF
<u>SPNS2 Antibody</u>	CSB-PA815559LA01HU	Human	ELISA, IF
<u>SPO11 Antibody</u>	CSB-PA897506LA01HU	Human	ELISA, WB
<u>SPOCK1 Antibody</u>	CSB-PA730785LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SPOP Antibody</u>	CSB-PA022601LA01HU	Human	ELISA, IHC
<u>SPOPL Antibody</u>	CSB-PA022602LA01HU	Human	ELISA, WB
<u>Sporamin A Antibody</u>	CSB-PA320416HA01INB	Ipomoea batatas	ELISA
<u>Spp1 Antibody</u>	CSB-PA022603YA01MO	Mouse	ELISA, WB
<u>SPP1 Antibody</u>	CSB-PA06904A0Rb	Human	ELISA, IHC
<u>SPPL2A Antibody</u>	CSB-PA855050LA01HU	Human	ELISA, WB, IHC, IF
<u>SPPL2B Antibody</u>	CSB-PA823457LA01HU	Human	ELISA, IF
<u>SPPL3 Antibody</u>	CSB-PA848402LA01HU	Human	ELISA, WB, IF
<u>SPRED1 Antibody</u>	CSB-PA773600LA01HU	Human	ELISA, WB, IHC, IF
<u>SPRR1A Antibody</u>	CSB-PA022610LA01HU	Human	ELISA, IHC, IF
<u>SPRR1B Antibody</u>	CSB-PA022611LA01HU	Human	ELISA, IHC
<u>SPRR2A Antibody</u>	CSB-PA842658LA01HU	Human	ELISA
<u>SPRTN Antibody</u>	CSB-PA863937LA01HU	Human	ELISA, IHC
<u>SPRY1 Antibody</u>	CSB-PA022621LA01HU	Human, Mouse	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>SPRY2 Antibody</u>	CSB-PA022622LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SPRY4 Antibody</u>	CSB-PA866335LA01HU	Human	ELISA, IF
<u>SPT16 Antibody</u>	CSB-PA022955LA01DOA	Arabidopsis thaliana	ELISA
<u>SPTBN1 Antibody</u>	CSB-PA733819LA01HU	Human	ELISA, IHC, IF
<u>SPTLC2 Antibody</u>	CSB-PA022640LA01HU	Human	ELISA, WB, IHC, IF
<u>SPTLC3 Antibody</u>	CSB-PA865142LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SPZ1 Antibody</u>	CSB-PA887147LA01HU	Human, Mouse	ELISA, WB, IHC
<u>SQSTM1 Antibody</u>	CSB-PA615696LA01HU	Human	ELISA, IHC, IF, IP
<u>src Antibody</u>	CSB-PA24479A0Rb	Zebrafish	ELISA, WB
<u>SREBF1 Antibody</u>	CSB-PA022657LA01HU	Human	ELISA, WB, IHC, IF
<u>SREBF2 Antibody</u>	CSB-PA666241LA01HU	Human	ELISA, IHC, IF
<u>SRF Antibody</u>	CSB-PA02857A0Rb	Human, Mouse	ELISA, WB, IHC
<u>SRFBP1 Antibody</u>	CSB-PA839864LA01HU	Human	ELISA
<u>SRGAP1 Antibody</u>	CSB-PA770358LA01HU	Human	ELISA, IF
<u>SRGAP2B Antibody</u>	CSB-PA20639A0Rb	Human, Mouse	ELISA, WB, IF
<u>SRGAP2C Antibody</u>	CSB-PA20519A0Rb	Human	ELISA, IHC, IF
<u>SRGAP3 Antibody</u>	CSB-PA022663LA01HU	Human	ELISA, IF
<u>SRGN Antibody</u>	CSB-PA022664EA01HU	Human	ELISA, IHC
<u>SRM Antibody</u>	CSB-PA01115A0Rb	Human	ELISA, WB, IHC
<u>SRMS Antibody</u>	CSB-PA862011LA01HU	Human	ELISA, IHC, IF
<u>SRP14 Antibody</u>	CSB-PA022673LA01HU	Human	ELISA, IF
<u>SRP19 Antibody</u>	CSB-PA022674HA01HU	Human	ELISA, WB, IHC
<u>SRP19 Antibody</u>	CSB-PA022674LA01HU	Human	ELISA, IHC, IF
<u>SRP54 Antibody</u>	CSB-PA022675LA01HU	Human, Mouse	ELISA, WB, IF
<u>SRP72 Antibody</u>	CSB-PA334457HA01SVG	Saccharomyces cerevisiae	ELISA
<u>SRPK2 Antibody</u>	CSB-PA022683LA01HU	Human	ELISA, WB, IHC
<u>SRPK3 Antibody</u>	CSB-PA892486LA01HU	Human	ELISA, IHC
<u>SRPRB Antibody</u>	CSB-PA896911LA01HU	Human	ELISA, WB, IF
<u>SRPX Antibody</u>	CSB-PA022689LA01HU	Human	ELISA, WB, IHC
<u>SRR Antibody</u>	CSB-PA022691LA01HU	Human, Mouse	ELISA, WB, IHC
<u>SRRM2 Antibody</u>	CSB-PA892369LA01HU	Human	ELISA, IHC
<u>SRRM4 Antibody</u>	CSB-PA012292LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>SRSF10 Antibody</u>	CSB-PA525638LA01HU	Human	ELISA, IHC
<u>SRSF9 Antibody</u>	CSB-PA00214A0Rb	Human, Mouse	ELISA, WB
<u>Sry Antibody</u>	CSB-PA022699LA01MO	Mouse	ELISA
<u>SS18L1 Antibody</u>	CSB-PA022701LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ssaA1 Antibody</u>	CSB-PA683042LA01FLB	Staphylococcus aureus	ELISA



Product Name	Code	Species Reactivity	Applications
<u>SSB Antibody</u>	CSB-PA022705HA01HU	Human, Mouse	ELISA, WB
<u>ssb Antibody</u>	CSB-PA770616LA01BQE	Bacillus anthracis	ELISA
<u>SSBP3 Antibody</u>	CSB-PA880121LA01HU	Human	ELISA, IHC
<u>SSH1 Antibody</u>	CSB-PA837881LA01HU	Human	ELISA, IHC, IF
<u>SSH3 Antibody</u>	CSB-PA837447LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>SSMEM1 Antibody</u>	CSB-PA845158LA01HU	Human	ELISA, IHC, IF
<u>sspB Antibody</u>	CSB-PA859203EA01SKX	Staphylococcus aureus	ELISA
<u>SSR2 Antibody</u>	CSB-PA022718LA01HU	Human	ELISA, IHC
<u>SSR3 Antibody</u>	CSB-PA892358LA01HU	Human	ELISA, WB, IHC, IF
<u>SSRP1 Antibody</u>	CSB-PA022721LA01HU	Human	ELISA, WB, IF
<u>SSSCA1 Antibody</u>	CSB-PA022722HA01HU	Human	ELISA, WB, IHC
<u>SST Antibody</u>	CSB-PA022723NJ01HU	Human	ELISA, IHC
<u>SST Antibody</u>	CSB-PA022723OJ01HU	Human	ELISA, IHC
<u>SSTR2 Antibody</u>	CSB-PA022725LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SSU72 Antibody</u>	CSB-PA865091LA01HU	Human	ELISA, IF
<u>SSX1 Antibody</u>	CSB-PA614964LA01HU	Human	ELISA, IHC, IF
<u>SSX2 Antibody</u>	CSB-PA613674LA01HU	Human	ELISA
<u>SSX2IP Antibody</u>	CSB-PA897526LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>ST3GAL2 Antibody</u>	CSB-PA618099LA01HU	Human	ELISA
<u>ST3GAL3 Antibody</u>	CSB-PA606136LA01HU	Human, Mouse	ELISA, WB
<u>ST3GAL4 Antibody</u>	CSB-PA612121LA01HU	Human	ELISA, IHC
<u>ST5 Antibody</u>	CSB-PA022758LA01HU	Human	ELISA, WB, IHC
<u>ST5 Antibody</u>	CSB-PA022758NA01HU	Human	ELISA, IHC
<u>ST6GALNAC1 Antibody</u>	CSB-PA868312LA01HU	Human	ELISA, WB, IHC
<u>ST6GALNAC2 Antibody</u>	CSB-PA883452LA01HU	Human	ELISA, WB
<u>ST6GALNAC3 Antibody</u>	CSB-PA836724LA01HU	Human	ELISA
<u>ST6GALNAC6 Antibody</u>	CSB-PA850251LA01HU	Human	ELISA, IHC
<u>ST8SIA4 Antibody</u>	CSB-PA821630LA01HU	Human	ELISA
<u>ST8SIA5 Antibody</u>	CSB-PA022778LA01HU	Human	ELISA, WB, IHC
<u>STAC Antibody</u>	CSB-PA859100LA01HU	Human, Mouse	ELISA, WB, IF
<u>STAC3 Antibody</u>	CSB-PA839373LA01HU	Human, Rat	ELISA, WB
<u>STAG3 Antibody</u>	CSB-PA022787LA01HU	Human	ELISA, IHC
<u>STAMBP Antibody</u>	CSB-PA022794HA01HU	Human	ELISA, WB, IHC
<u>STAP1 Antibody</u>	CSB-PA890756LA01HU	Human	ELISA
<u>STAR Antibody</u>	CSB-PA022798LA01HU	Human	ELISA, IHC
<u>STARD10 Antibody</u>	CSB-PA896706LA01HU	Human, Mouse	ELISA, WB
<u>STARD13 Antibody</u>	CSB-PA896714LA01HU	Human	ELISA, IF
<u>STARD7 Antibody</u>	CSB-PA03375A0Rb	Human	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>STARD9 Antibody</u>	CSB-PA885805LA01HU	Human	ELISA, WB, IHC
<u>STAT1 Antibody</u>	CSB-PA022810LA01HU	Human	ELISA, IHC
<u>STAT1 Antibody</u>	CSB-PA022810HA01HU	Human	ELISA, IHC, IF
<u>STAT2 Antibody</u>	CSB-PA022811LA01HU	Human	ELISA, IHC, IF
<u>STAT3 Antibody</u>	CSB-PA022812LA01HU	Human	ELISA, IHC, IF, IP
<u>STAT4 Antibody</u>	CSB-PA022813EA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>STAT5B Antibody</u>	CSB-PA022815LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF, IP
<u>STAT6 Antibody</u>	CSB-PA022816LA01HU	Human	ELISA, IHC, IF
<u>STAU1 Antibody</u>	CSB-PA13759A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>STC1 Antibody</u>	CSB-PA09889A0Rb	Human	ELISA, IHC
<u>STE13 Antibody</u>	CSB-PA339395LA01SVG	Saccharomyces cerevisiae	ELISA, WB
<u>STEAP1 Antibody</u>	CSB-PA13354A0Rb	Human	ELISA, WB, IHC, IF
<u>STEAP3 Antibody</u>	CSB-PA717655LA01HU	Human	ELISA, WB, IHC, IF
<u>STH Antibody</u>	CSB-PA022827LA01HU	Human	ELISA, IHC
<u>STIL Antibody</u>	CSB-PA618008LA01HU	Human	ELISA, IHC, IF
<u>STIM1 Antibody</u>	CSB-PA022829LA01HU	Human	ELISA, IHC
<u>STIM2 Antibody</u>	CSB-PA022830LA01HU	Human, Mouse	ELISA, WB, IHC
<u>STIP1 Antibody</u>	CSB-PA022831HA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>STK17B Antibody</u>	CSB-PA022838LA01HU	Human, Mouse	ELISA, WB, IHC
<u>STK24 Antibody</u>	CSB-PA00249A0Rb	Human	ELISA, WB, IHC, IF
<u>STK25 Antibody</u>	CSB-PA022842LA01HU	Human	ELISA, IHC
<u>STK33 Antibody</u>	CSB-PA861179LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>STK35 Antibody</u>	CSB-PA851552LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>STK38L Antibody</u>	CSB-PA896491LA01HU	Human	ELISA, WB, IHC
<u>STK4 Antibody</u>	CSB-PA621644NA01HU	Human	ELISA, IHC, IF
<u>STMN1 Antibody</u>	CSB-PA022857LA01HU	Human	ELISA, IHC, IF
<u>STMN2 Antibody</u>	CSB-PA853100LA01HU	Human	ELISA, IHC, IF
<u>STOML1 Antibody</u>	CSB-PA887004LA01HU	Human	ELISA, IF
<u>STOML3 Antibody</u>	CSB-PA844707LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>STOX1 Antibody</u>	CSB-PA765089LA01HU	Human	ELISA, WB, IHC
<u>STPG2 Antibody</u>	CSB-PA847628LA01HU	Human	ELISA, WB, IHC
<u>STRA6 Antibody</u>	CSB-PA866301LA01HU	Human	ELISA, IHC, IF
<u>STRA8 Antibody</u>	CSB-PA767212LA01HU	Human	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>STRADB Antibody</u>	CSB-PA874868EA01HU	Human	ELISA, WB, IHC
<u>STRBP Antibody</u>	CSB-PA847225LA01HU	Human	ELISA, WB, IF
<u>Streptavidin Antibody</u>	CSB-PA18194A0Rb	Streptomyces avidinii	ELISA
<u>Streptavidin Antibody</u>	CSB-PA332849YA01SNO	Streptomyces avidinii	ELISA
<u>STT3A Antibody</u>	CSB-PA022885LA01HU	Human	ELISA, IHC
<u>STT3B Antibody</u>	CSB-PA823456LA01HU	Human	ELISA, IHC
<u>STUB1 Antibody</u>	CSB-PA892480LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>STX10 Antibody</u>	CSB-PA022888HA01HU	Human	ELISA, IHC
<u>STX10 Antibody</u>	CSB-PA022888LA01HU	Human	ELISA, IHC
<u>STX19 Antibody</u>	CSB-PA822703LA01HU	Human	ELISA, WB
<u>STX1A Antibody</u>	CSB-PA022895LA01HU	Human	ELISA, WB
<u>STX4 Antibody</u>	CSB-PA623656LA01HU	Human	ELISA, IF
<u>STX6 Antibody</u>	CSB-PA022902DA01HU	Human	ELISA, IHC
<u>STX6 Antibody</u>	CSB-PA022902LA01CH	Chicken, Mouse	ELISA, WB
<u>STX8 Antibody</u>	CSB-PA892159LA01HU	Human	ELISA, IHC
<u>stxB Antibody</u>	CSB-PA300809HA01EKZ	Enterobacteria H19B	phage ELISA
<u>stxB2 Antibody</u>	CSB-PA357663HA01ECC	Enterobacteria 933W	phage ELISA
<u>STXBP1 Antibody</u>	CSB-PA022905LA01HU	Human, Rat, Mouse	ELISA, WB, IHC, IF
<u>STXBP3 Antibody</u>	CSB-PA022907LA01HU	Human, Mouse	ELISA, WB, IHC
<u>STXBP5 Antibody</u>	CSB-PA722568LA01HU	Human	ELISA, IHC, IF
<u>STXBP6 Antibody</u>	CSB-PA850881LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>STYX Antibody</u>	CSB-PA851966LA01HU	Human	ELISA, WB
<u>SUB1 Antibody</u>	CSB-PA022915HA01HU	Human	ELISA, IHC
<u>SUCLG1 Antibody</u>	CSB-PA022919LA01HU	Human	ELISA, IHC
<u>SUCLG2 Antibody</u>	CSB-PA846636LA01HU	Human	ELISA, WB, IHC
<u>SUCNR1 Antibody</u>	CSB-PA022921LA01HU	Human	ELISA, WB, IHC, IF
<u>sulA Antibody</u>	CSB-PA365383LA01ENV	Escherichia coli	ELISA
<u>SULF1 Antibody</u>	CSB-PA814217LA01HU	Human	ELISA, WB, IHC, IF
<u>SULF2 Antibody</u>	CSB-PA810283LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>SULT1A2 Antibody</u>	CSB-PA022934LA01HU	Human	ELISA, IHC
<u>SULT1A3 Antibody</u>	CSB-PA00887A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>SULT1B1 Antibody</u>	CSB-PA022937EA01HU	Human, Mouse	ELISA, WB, IHC
<u>Sumo tag Monoclonal</u>	CSB-MA000132M0m	All	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>Antibody</u>			
<u>Sumo tag Monoclonal Antibody</u>	CSB-MA000131M0m	specific for sumo tag fusion protein	ELISA, WB
<u>SUMO2 Antibody</u>	CSB-PA05009A0Rb	Human	ELISA, IHC
<u>SUMO3 Antibody</u>	CSB-PA05509A0Rb	Human	ELISA, IHC, IF
<u>SUN1 Antibody</u>	CSB-PA24689A0Rb	Arabidopsis thaliana	ELISA, WB
<u>SUN2 Antibody</u>	CSB-PA883386LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>SUOX Antibody</u>	CSB-PA022954HA01HU	Human, Rat	ELISA, WB, IHC
<u>SUPT16H Antibody</u>	CSB-PA835618LA01HU	Human	ELISA, IHC
<u>SUPT3H Antibody</u>	CSB-PA022956LA01HU	Human, Mouse	ELISA, WB, IHC
<u>SUPV3L1 Antibody</u>	CSB-PA814232LA01HU	Human	ELISA, WB, IHC
<u>surA Antibody</u>	CSB-PA359693HA01ENV	Escherichia coli	ELISA, WB
<u>SUS1 Antibody</u>	CSB-PA335680LA01SVV	Saccharum officinarum	ELISA
<u>SUSD4 Antibody</u>	CSB-PA733164LA01HU	Human	ELISA, IHC
<u>SUSD5 Antibody</u>	CSB-PA022971LA01HU	Human	ELISA, WB, IHC, IF
<u>SUSD6 Antibody</u>	CSB-PA835675LA01HU	Human	ELISA, IHC, IF
<u>SUV39H1 Antibody</u>	CSB-PA022972LA01HU	Human	ELISA, IHC
<u>SUZ12 Antibody</u>	CSB-PA621859LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>SV2A Antibody</u>	CSB-PA022978LA01HU	Human	ELISA, IHC
<u>SV2B Antibody</u>	CSB-PA022979LA01HU	Human, Mouse	ELISA, WB, IF
<u>SV2C Antibody</u>	CSB-PA022980LA01HU	Human	ELISA, IHC
<u>SVBP Antibody</u>	CSB-PA836653LA01HU	Human	ELISA
<u>SVOPL Antibody</u>	CSB-PA818698LA01HU	Human, Mouse, Rat	ELISA, WB, IF
<u>SWAP70 Antibody</u>	CSB-PA887026LA01HU	Human, Mouse	ELISA, WB
<u>SWI5 Antibody</u>	CSB-PA636328LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>SYCE3 Antibody</u>	CSB-PA003820LA01HU	Human	ELISA, IHC
<u>SYCN Antibody</u>	CSB-PA022989LA01HU	Human	ELISA, IHC
<u>SYF2 Antibody</u>	CSB-PA022996LA01HU	Human	ELISA
<u>SYK Antibody</u>	CSB-PA022998LA01HU	Human	ELISA, IHC, IF
<u>SYN1 Antibody</u>	CSB-PA023002DA01HU	Human, Mouse	ELISA, WB, IHC
<u>SYNC Antibody</u>	CSB-PA862034LA01HU	Human	ELISA
<u>SYNCRIP Antibody</u>	CSB-PA03334A0Rb	Human	ELISA, WB, IHC
<u>SYNE2 Antibody</u>	CSB-PA841001LA01HU	Human	ELISA, IHC, IF
<u>SYNE4 Antibody</u>	CSB-PA847615LA01HU	Human	ELISA
<u>SYNGAP1 Antibody</u>	CSB-PA857006LA01HU	Human	ELISA, IF
<u>SYNJ1 Antibody</u>	CSB-PA023016LA01HU	Human	ELISA, WB, IHC, IF
<u>SYNJ2BP Antibody</u>	CSB-PA023018LA01HU	Human, Rat	ELISA, WB, IHC, IP



Product Name	Code	Species Reactivity	Applications
<u>SYNM Antibody</u>	CSB-PA023019LA01HU	Human	ELISA, IHC, IF
<u>SYNPR Antibody</u>	CSB-PA023023LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>SYT1 Antibody</u>	CSB-PA023029LA01HU	Human	ELISA
<u>SYT10 Antibody</u>	CSB-PA754474LA01HU	Human	ELISA, IHC
<u>SYT6 Antibody</u>	CSB-PA716356LA01HU	Human	ELISA, IF
<u>SYTL3 Antibody</u>	CSB-PA687190LA01HU	Human	ELISA, WB, IHC
<u>SYVNI Antibody</u>	CSB-PA803115LA01HU	Human	ELISA, IHC, IF
<u>SZT2 Antibody</u>	CSB-PA735962LA01HU	Human	ELISA, IHC
<u>TAAR5 Antibody</u>	CSB-PA023057LA01HU	Human, Rat, Mouse	ELISA, WB, IHC, IF
<u>TAB3 Antibody</u>	CSB-PA850806LA01HU	Human	ELISA, IHC, IF
<u>TAC1 Antibody</u>	CSB-PA023062EA01HU	Human	ELISA
<u>TAC1 Antibody</u>	CSB-PA023062LA01HU	Human	ELISA, WB, IHC
<u>TACC3 Antibody</u>	CSB-PA896927LA01HU	Human	ELISA, IHC, IF
<u>TACR1 Antibody</u>	CSB-PA09674A0Rb	Human	ELISA, IHC, IF
<u>TADA2B Antibody</u>	CSB-PA768204LA01HU	Human	ELISA, IHC
<u>TADA3 Antibody</u>	CSB-PA03494A0Rb	Human	ELISA
<u>TAF12 Antibody</u>	CSB-PA624102LA01HU	Human	ELISA, IHC
<u>TAF13 Antibody</u>	CSB-PA614893LA01HU	Human	ELISA, IHC
<u>TAF15 Antibody</u>	CSB-PA856431LA01HU	Human, Rat, Mouse	ELISA, WB, IHC, IF
<u>TAF1A Antibody</u>	CSB-PA614895LA01HU	Human	ELISA, WB, IHC
<u>TAF2 Antibody</u>	CSB-PA747395LA01HU	Human	ELISA, IHC
<u>TAF4B Antibody</u>	CSB-PA849782LA01HU	Human	ELISA, IHC, IF
<u>TAF6 Antibody</u>	CSB-PA023095LA01HU	Human	ELISA, IHC
<u>TAF6L Antibody</u>	CSB-PA896758LA01HU	Human	ELISA, IHC
<u>TAF7L Antibody</u>	CSB-PA688972LA01HU	Human	ELISA, IHC, IF
<u>TAF8 Antibody</u>	CSB-PA773805LA01HU	Human, Mouse, Rat	ELISA, WB
<u>TAGAP Antibody</u>	CSB-PA818676LA01HU	Human	ELISA, IHC
<u>TAGLN Antibody</u>	CSB-PA00804A0Rb	Human	ELISA, WB, IHC
<u>tal1 Antibody</u>	CSB-PA20969A0Rb	Zebrafish	ELISA, WB
<u>TAL1 Antibody</u>	CSB-PA023110LA01HU	Human, Rat	ELISA, WB, IF
<u>TALDO1 Antibody</u>	CSB-PA04275A0Rb	Human	ELISA, WB, IHC
<u>TANGO2 Antibody</u>	CSB-PA743560LA01HU	Human	ELISA
<u>TANK Antibody</u>	CSB-PA023116LA01HU	Human	ELISA, IHC, IF
<u>TAOK1 Antibody</u>	CSB-PA752101LA01HU	Human	ELISA, WB, IHC, IF
<u>TAOK2 Antibody</u>	CSB-PA892341LA01HU	Human	ELISA, IHC
<u>TAP2 Antibody</u>	CSB-PA22539A0Rb	Human	ELISA, IHC
<u>TAPBPL Antibody</u>	CSB-PA887144LA01HU	Human	ELISA, IHC
<u>tar Antibody</u>	CSB-PA356999HA01ENV	Escherichia coli	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>TARDBP Antibody</u>	CSB-PA023129EA01HU	Human	ELISA, WB, IHC, IF
<u>TARS2 Antibody</u>	CSB-PA023132LA01HU	Human	ELISA, IHC
<u>TAS1R2 Antibody</u>	CSB-PA840577LA01HU	Human	ELISA, IF
<u>TAS2R13 Antibody</u>	CSB-PA889153LA01HU	Human	ELISA, IHC, IF
<u>TAS2R14 Antibody</u>	CSB-PA868362LA01HU	Human	ELISA, IHC
<u>TAS2R38 Antibody</u>	CSB-PA023155LA01HU	Human	ELISA, IHC, IF
<u>TAS2R41 Antibody</u>	CSB-PA023159LA01HU	Human	ELISA, IHC, IF
<u>TAS2R42 Antibody</u>	CSB-PA766403LA01HU	Human	ELISA, IHC
<u>TAT Antibody</u>	CSB-PA023175LA01HU	Human	ELISA, WB, IF
<u>TAX1BP1 Antibody</u>	CSB-PA771448LA01HU	Human	ELISA, WB
<u>TAX1BP3 Antibody</u>	CSB-PA04445A0Rb	Human	ELISA, WB
<u>TAZ Antibody</u>	CSB-PA023181EA01HU	Human	ELISA, IHC, IF
<u>TBC1D1 Antibody</u>	CSB-PA803114LA01HU	Human	ELISA, IHC
<u>TBC1D14 Antibody</u>	CSB-PA889187LA01HU	Human	ELISA, WB, IHC
<u>TBC1D15 Antibody</u>	CSB-PA823452LA01HU	Human	ELISA, WB, IHC, IF
<u>TBC1D19 Antibody</u>	CSB-PA843274LA01HU	Human	ELISA, WB, IHC
<u>TBC1D20 Antibody</u>	CSB-PA846594LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TBC1D22A Antibody</u>	CSB-PA848818LA01HU	Human, Rat	ELISA, WB
<u>TBC1D23 Antibody</u>	CSB-PA865143LA01HU	Human	ELISA, IHC
<u>TBC1D24 Antibody</u>	CSB-PA890752LA01HU	Human	ELISA, IHC, IF
<u>TBC1D25 Antibody</u>	CSB-PA662139LA01HU	Human	ELISA, IHC
<u>TBC1D26 Antibody</u>	CSB-PA769776LA01HU	Human, Mouse	ELISA, WB
<u>TBC1D3 Antibody</u>	CSB-PA814244LA01HU	Human	ELISA
<u>TBC1D31 Antibody</u>	CSB-PA822197LA01HU	Human	ELISA, IHC
<u>TBC1D3C Antibody</u>	CSB-PA740248LA01HU	Human	ELISA
<u>TBC1D4 Antibody</u>	CSB-PA023218LA01HU	Human	ELISA, IHC, IF
<u>TBC1D7 Antibody</u>	CSB-PA873720LA01HU	Human	ELISA, IHC, IF
<u>TBC1D9B Antibody</u>	CSB-PA724138LA01HU	Human	ELISA, IHC, IF
<u>TBCA Antibody</u>	CSB-PA023225HA01HU	Human	ELISA, IHC
<u>TBCB Antibody</u>	CSB-PA03169A0Rb	Human	ELISA, WB, IHC, IF, IP
<u>TBCC Antibody</u>	CSB-PA023227LA01HU	Human	ELISA, WB, IHC
<u>TBCD Antibody</u>	CSB-PA883599LA01HU	Human	ELISA, IHC, IF
<u>TBCEL Antibody</u>	CSB-PA722183LA01HU	Human, Rat	ELISA, WB
<u>TBK1 Antibody</u>	CSB-PA023232LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TBK1 Antibody</u>	CSB-PA890690NA01HU	Human	ELISA, WB, IF
<u>TBKBP1 Antibody</u>	CSB-PA023233LA01HU	Human	ELISA, IF
<u>TBLIX Antibody</u>	CSB-PA023234LA01HU	Human	ELISA, WB, IF



Product Name	Code	Species Reactivity	Applications
<u>TBP Antibody</u>	CSB-PA023239HA01HU	Human	ELISA, IHC, IF
<u>TBPL1 Antibody</u>	CSB-PA023240LA01HU	Human	ELISA, IHC
<u>TBX1 Antibody</u>	CSB-PA023245YA01HU	Human	ELISA, IF
<u>Tbx15 Antibody</u>	CSB-PA023247LA01MO	Mouse	ELISA
<u>TBX15 Antibody</u>	CSB-PA839414HA01HU	Human	ELISA, WB
<u>fbx16 Antibody</u>	CSB-PA20979A0Rb	Zebrafish	ELISA
<u>TBX18 Antibody</u>	CSB-PA023248LA01HU	Human	ELISA
<u>TBX2 Antibody</u>	CSB-PA613386LA01HU	Human, Rat	ELISA, WB
<u>TBX21 Antibody</u>	CSB-PA890744LA01HU	Human, Mouse	ELISA, WB
<u>TBX3 Antibody</u>	CSB-PA023255LA01HU	Human	ELISA, IHC
<u>TBX4 Antibody</u>	CSB-PA023256LA01HU	Human	ELISA, WB
<u>fbx5a Antibody</u>	CSB-PA881254LA01DIL	Zebrafish	ELISA
<u>TBX6 Antibody</u>	CSB-PA023258LA01HU	Human	ELISA, IHC
<u>TCAF2 Antibody</u>	CSB-PA008015LA01HU	Human	ELISA, IHC, IF
<u>TCEA3 Antibody</u>	CSB-PA023271LA01HU	Human	ELISA, IHC, IF
<u>TCEAL7 Antibody</u>	CSB-PA863109LA01HU	Human	ELISA, WB, IHC
<u>TCEAL8 Antibody</u>	CSB-PA818261LA01HU	Human	ELISA, IHC
<u>TCEAL9 Antibody</u>	CSB-PA887044LA01HU	Human	ELISA
<u>TCEANC2 Antibody</u>	CSB-PA842718LA01HU	Human	ELISA, IHC
<u>TCERG1 Antibody</u>	CSB-PA023291LA01HU	Human, Rat	ELISA, WB, IHC
<u>TCF15 Antibody</u>	CSB-PA619636LA01HU	Human	ELISA, IHC
<u>TCF20 Antibody</u>	CSB-PA887024LA01HU	Human	ELISA, IHC
<u>TCF3 Antibody</u>	CSB-PA023301HA01HU	Human	ELISA, IHC
<u>TCF4 Antibody</u>	CSB-PA023302HA01HU	Human, Mouse	ELISA, WB, IHC
<u>TCF7L1 Antibody</u>	CSB-PA884626OA01HU	Human	ELISA, WB
<u>TCHH Antibody</u>	CSB-PA023309OA01HU	Human	ELISA, IHC, IF
<u>TCHP Antibody</u>	CSB-PA874804LA01HU	Human	ELISA, IHC, IF
<u>TCIRG1 Antibody</u>	CSB-PA615690LA01HU	Human	ELISA, IHC, IF
<u>TCL1A Antibody</u>	CSB-PA023313LA01HU	Human	ELISA, IHC, IF
<u>TCL1B Antibody</u>	CSB-PA023314LA01HU	Human	ELISA, IHC
<u>TCN1 Antibody</u>	CSB-PA023317LA01HU	Human	ELISA, IHC
<u>TCN2 Antibody</u>	CSB-PA023318LA01HU	Human	ELISA, WB, IF
<u>TCP1 Antibody</u>	CSB-PA023320LA01HU	Human	ELISA, IHC, IF
<u>TCP10 Antibody</u>	CSB-PA614256HA01HU	Human	ELISA
<u>TCP10L Antibody</u>	CSB-PA823465LA01HU	Human	ELISA, IHC
<u>TCP11L2 Antibody</u>	CSB-PA847634LA01HU	Human	ELISA, WB, IHC
<u>TCTA Antibody</u>	CSB-PA023331LA01HU	Human, Mouse	ELISA, WB, IHC
<u>TCTE1 Antibody</u>	CSB-PA689339LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TCTEX1D2 Antibody</u>	CSB-PA845157LA01HU	Human	ELISA, IHC
<u>TCTN3 Antibody</u>	CSB-PA747352LA01HU	Human, Rat	ELISA, WB, IF
<u>teyA Antibody</u>	CSB-PA331446HA01BRJ	Bacillus subtilis	ELISA



Product Name	Code	Species Reactivity	Applications
<u>TDGF1 Antibody</u>	CSB-PA10764A0Rb	Human	ELISA
<u>TDGF1 Antibody</u>	CSB-PA10769A0Rb	Human	ELISA
<u>TDRD1 Antibody</u>	CSB-PA874834LA01HU	Human	ELISA, IHC, IF
<u>TDRD12 Antibody</u>	CSB-PA707440LA01HU	Human	ELISA, IHC
<u>TDRD3 Antibody</u>	CSB-PA863996LA01HU	Human	ELISA, WB, IHC
<u>TDRD9 Antibody</u>	CSB-PA850866LA01HU	Human	ELISA, IHC
<u>TEAD3 Antibody</u>	CSB-PA860344LA01HU	Human	ELISA, IHC
<u>TEAD4 Antibody</u>	CSB-PA618010LA01HU	Human	ELISA, WB, IHC
<u>TECPR1 Antibody</u>	CSB-PA773603LA01HU	Human	ELISA, IHC, IF
<u>TECPR2 Antibody</u>	CSB-PA023369LA01HU	Human	ELISA, IHC, IF
<u>TECTB Antibody</u>	CSB-PA857003LA01HU	Human	ELISA, IF
<u>TEF Antibody</u>	CSB-PA608815LA01HU	Human	ELISA, WB
<u>TEFM Antibody</u>	CSB-PA857011LA01HU	Human	ELISA, IHC
<u>TEKT3 Antibody</u>	CSB-PA871598LA01HU	Human	ELISA, IHC
<u>TEKT4 Antibody</u>	CSB-PA855486LA01HU	Human	ELISA
<u>TEKT5 Antibody</u>	CSB-PA856979LA01HU	Human, Mouse	ELISA, WB
<u>TENM1 Antibody</u>	CSB-PA892339LA01HU	Human	ELISA, IHC
<u>TENT5B Antibody</u>	CSB-PA856891LA01HU	Human	ELISA, WB
<u>TENT5C Antibody</u>	CSB-PA736203LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TEP1 Antibody</u>	CSB-PA18049A0Rb	Human	ELISA, IHC
<u>TERB1 Antibody</u>	CSB-PA822742LA01HU	Human	ELISA, IHC
<u>TERF1 Antibody</u>	CSB-PA023386EA01HU	Human	ELISA, IHC
<u>TERF2 Antibody</u>	CSB-PA620883LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TERT Antibody</u>	CSB-PA10669A0Rb	Human	ELISA, IHC
<u>TERT Antibody</u>	CSB-PA023391LA01DO	Dog	ELISA
<u>TES Antibody</u>	CSB-PA890679LA01HU	Human	ELISA, WB, IHC, IF
<u>TES-26 Antibody</u>	CSB-PA347478HA01THA	Toxocara canis	ELISA
<u>TESC Antibody</u>	CSB-PA023393HA01HU	Human	ELISA, IF
<u>TESPA1 Antibody</u>	CSB-PA012205LA01HU	Human	ELISA, IHC, IF
<u>TET1 Antibody</u>	CSB-PA847709HA01HU	Human	ELISA, IHC
<u>TET3 Antibody</u>	CSB-PA023398LA01HU	Human	ELISA, IHC
<u>tetR Antibody</u>	CSB-PA366163HA01ENL	Escherichia coli	ELISA
<u>tetX Antibody</u>	CSB-PA356416LA01CMN	Clostridium tetani	ELISA
<u>TEX12 Antibody</u>	CSB-PA023402EA01HU	Human	ELISA
<u>TEX13A Antibody</u>	CSB-PA883626LA01HU	Human	ELISA, IHC
<u>TEX35 Antibody</u>	CSB-PA732917LA01HU	Human	ELISA, WB, IHC
<u>TEX43 Antibody</u>	CSB-PA765062LA01HU	Human	ELISA, IHC
<u>TEX47 Antibody</u>	CSB-PA851537LA01HU	Human	ELISA, WB, IHC
<u>TF Antibody</u>	CSB-PA00250E0Rb	Human	ELISA, WB



Product Name	Code	Species Reactivity	Applications
<u>Tf Antibody</u>	CSB-PA00470E0Rb	Rat	ELISA, WB
<u>Tf Antibody</u>	CSB-PA856363LA01MO	Mouse	ELISA, IHC
<u>TFAM Antibody</u>	CSB-PA023413LA01HU	Human	ELISA, WB, IHC
<u>TFAP2B Antibody</u>	CSB-PA765731LA01HU	Human	ELISA, WB
<u>TFAP2C Antibody</u>	CSB-PA842150LA01HU	Human	ELISA, IHC, IF
<u>TFAP2E Antibody</u>	CSB-PA761516LA01HU	Human	ELISA, WB
<u>TFB1M Antibody</u>	CSB-PA023420LA01HU	Human	ELISA, WB, IHC, IF
<u>TFB2M Antibody</u>	CSB-PA867154LA01HU	Human	ELISA, IHC
<u>TFDP2 Antibody</u>	CSB-PA623906LA01HU	Human, Mouse	ELISA, WB, IHC
<u>TFDP3 Antibody</u>	CSB-PA702134LA01HU	Human	ELISA, IHC, IF
<u>TFE3 Antibody</u>	CSB-PA023428LA01HU	Human	ELISA, WB, IHC
<u>TFEB Antibody</u>	CSB-PA023429EA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>TFEC Antibody</u>	CSB-PA023430LA01HU	Human	ELISA, WB
<u>TFF3 Antibody</u>	CSB-PA023433YA01HU	Human	ELISA
<u>TFF3 Antibody</u>	CSB-PA11935A0Rb	Human	ELISA
<u>TFF3 Antibody</u>	CSB-PA023433LA01HU	Human	ELISA, IHC
<u>TFPI2 Antibody</u>	CSB-PA023438LA01HU	Human	ELISA, IHC
<u>TFR2 Antibody</u>	CSB-PA14729A0Rb	Human	ELISA, IHC, IF
<u>TFRC Antibody</u>	CSB-PA07219A0Rb	Human	ELISA, WB, IHC, IF
<u>TG Antibody</u>	CSB-PA12325A0Rb	Human	ELISA, IHC
<u>TGDS Antibody</u>	CSB-PA023444LA01HU	Human, Rat	ELISA, WB, IHC
<u>TGFA Antibody</u>	CSB-PA07294A0Rb	Human	ELISA, WB, IHC, IF
<u>TGFB1 Antibody</u>	CSB-PA023446LA01BO	Bovine	ELISA
<u>TGFB1 Antibody</u>	CSB-PA023446LA01CH	Chicken	ELISA
<u>tgfb1 Antibody</u>	CSB-PA023446LA01OEI	Oncorhynchus mykiss, Human	ELISA, WB
<u>TGFB1 Antibody</u>	CSB-PA07309A0Rb	Human	ELISA, IHC
<u>Tgfb1 Antibody</u>	CSB-PA023446LA01MO	Mouse, Human	ELISA, WB, IP
<u>TGFB1I1 Antibody</u>	CSB-PA13049A0Rb	Human	ELISA, WB, IHC
<u>TGFB2 Antibody</u>	CSB-PA07315A0Rb	Human	ELISA, WB, IHC
<u>TGFB2 Antibody</u>	CSB-PA07319A0Rb	Human, Mouse	ELISA, WB, IHC
<u>TGFB3 Antibody</u>	CSB-PA023449LA01HU	Human	ELISA, IHC
<u>TGFB3 Antibody</u>	CSB-PA17749A0Rb	Human	ELISA, IHC
<u>TGFBR1 Antibody</u>	CSB-PA023451LA01HU	Human	ELISA, WB, IHC, IF
<u>TGFBR2 Antibody</u>	CSB-PA023452LA01HU	Human	ELISA, IHC, IF
<u>TGFBR3L Antibody</u>	CSB-PA27472A0Rb	Human	ELISA, IF
<u>TGFBRAP1 Antibody</u>	CSB-PA819888LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>TGIF1 Antibody</u>	CSB-PA623005LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TGIF2 Antibody</u>	CSB-PA884421LA01HU	Human	ELISA, IHC, IF
<u>TGM2 Antibody</u>	CSB-PA11809A0Rb	Human	ELISA, WB, IHC, IF
<u>TGM4 Antibody</u>	CSB-PA023464LA01HU	Human	ELISA, IHC, IF
<u>TGS1 Antibody</u>	CSB-PA847219LA01HU	Human, Mouse	ELISA, WB, IHC
<u>TH Antibody</u>	CSB-PA023470LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>th Antibody</u>	CSB-PA24419A0Rb	Zebrafish	ELISA, WB
<u>THADA Antibody</u>	CSB-PA751092LA01HU	Human	ELISA, IHC, IF
<u>THAP11 Antibody</u>	CSB-PA846612LA01HU	Human	ELISA
<u>THAP12 Antibody</u>	CSB-PA018719LA01HU	Human	ELISA, WB, IHC
<u>Thbs1 Antibody</u>	CSB-PA023487LA01MO	Mouse	ELISA
<u>THBS1 Antibody</u>	CSB-PA07349A0Rb	Human	ELISA
<u>Thbs2 Antibody</u>	CSB-PA023488HA01MO	Mouse	ELISA
<u>Thbs3 Antibody</u>	CSB-PA023489DA01MO	Mouse	ELISA
<u>Thbs3 Antibody</u>	CSB-PA023489EA01MO	Mouse, Human	ELISA, WB
<u>THEMIS Antibody</u>	CSB-PA850780LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>THEMIS2 Antibody</u>	CSB-PA735984LA01HU	Human	ELISA, WB, IHC, IF
<u>THI2.1 Antibody</u>	CSB-PA365586LA01VDO	Viscum album	ELISA, WB
<u>Thiaminase-1 Antibody</u>	CSB-PA332698HA01PAI	Paenibacillus thiaminolyticus	ELISA, WB
<u>thiO Antibody</u>	CSB-PA521738LA01BRJ	Bacillus subtilis	ELISA
<u>THNSL2 Antibody</u>	CSB-PA803155LA01HU	Human	ELISA, WB, IHC
<u>THOC3 Antibody</u>	CSB-PA822232LA01HU	Human	ELISA, IHC
<u>THOC5 Antibody</u>	CSB-PA613425LA01HU	Human	ELISA, IHC, IF
<u>THOC6 Antibody</u>	CSB-PA769793LA01HU	Human	ELISA, WB
<u>THOC7 Antibody</u>	CSB-PA746911LA01HU	Human	ELISA
<u>Thpo Antibody</u>	CSB-PA07649A0Rb	Rat, Mouse, Human	ELISA, WB
<u>THRA Antibody</u>	CSB-PA023510LA01HU	Human	ELISA, WB, IHC
<u>THRB Antibody</u>	CSB-PA023514LA01HU	Human, Mouse	ELISA, WB, IHC
<u>Thrombin-like enzyme crotalase Antibody</u>	CSB-PA520546LA01DYB	Crotalus adamanteus	ELISA
<u>Thrsp Antibody</u>	CSB-PA726789HA01MO	Mouse	ELISA, IF
<u>THRSP Antibody</u>	CSB-PA846095LA01HU	Human	ELISA, IHC
<u>THTPA Antibody</u>	CSB-PA883600LA01HU	Human	ELISA, IHC
<u>THY1 Antibody</u>	CSB-PA023525LA01HU	Human	ELISA
<u>TIAL1 Antibody</u>	CSB-PA023529HA01HU	Human	ELISA, IF
<u>TIAM1 Antibody</u>	CSB-PA614387LA01HU	Human	ELISA, IHC, IF
<u>TIAM2 Antibody</u>	CSB-PA818235LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>TIE1 Antibody</u>	CSB-PA07369A0Rb	Human	ELISA, IHC, IF
<u>TIFA Antibody</u>	CSB-PA839301LA01HU	Human	ELISA
<u>fig Antibody</u>	CSB-PA484773LA01EOB	Escherichia coli O127:H6	ELISA
<u>TIGAR Antibody</u>	CSB-PA889078LA01HU	Human	ELISA, IHC, IF
<u>TIGD4 Antibody</u>	CSB-PA818256LA01HU	Human	ELISA, IHC
<u>TIGD7 Antibody</u>	CSB-PA757502LA01HU	Human	ELISA, IHC
<u>TIGIT Antibody</u>	CSB-PA675446LA01HU	Human	ELISA, IF
<u>TIMD4 Antibody</u>	CSB-PA850304LA01HU	Human, Mouse	ELISA, WB, IF
<u>TIMELESS Antibody</u>	CSB-PA892360LA01HU	Human	ELISA, IHC
<u>TIMM10 Antibody</u>	CSB-PA023548LA01HU	Human	ELISA, IHC
<u>TIMM10B Antibody</u>	CSB-PA896532LA01HU	Human, Mouse	ELISA, WB, IHC
<u>TIMM8B Antibody</u>	CSB-PA897567LA01HU	Human	ELISA, IHC
<u>TIMM9 Antibody</u>	CSB-PA897296LA01HU	Human	ELISA, IHC
<u>TIMMDC1 Antibody</u>	CSB-PA885685LA01HU	Human, Mouse	ELISA, WB, IHC
<u>Timp1 Antibody</u>	CSB-PA023560LA01MO	Mouse	ELISA
<u>TIMP1 Antibody</u>	CSB-PA023560YA01HU	Human, Mouse	ELISA, WB, IHC
<u>Timp1 Antibody</u>	CSB-PA16719A0Rb	Rat	ELISA
<u>TIMP2 Antibody</u>	CSB-PA07405A0Rb	Human	ELISA, WB
<u>TIMP3 Antibody</u>	CSB-PA11749A0Rb	Human	ELISA, WB
<u>TIMP3 Antibody</u>	CSB-PA023562LA01HU	Human	ELISA, IHC
<u>TIMP4 Antibody</u>	CSB-PA023563LA01BO	Bovine	ELISA
<u>TIMP4 Antibody</u>	CSB-PA07419A0Rb	Human	ELISA, IHC, IF
<u>TINAGL1 Antibody</u>	CSB-PA875630LA01HU	Human	ELISA, IHC, IF
<u>TINF2 Antibody</u>	CSB-PA874799LA01HU	Human	ELISA, IHC
<u>TIPARP Antibody</u>	CSB-PA801800LA01HU	Human	ELISA, WB, IHC
<u>TIPRL Antibody</u>	CSB-PA023569LA01HU	Human	ELISA
<u>tir Antibody</u>	CSB-PA514896HA01EOF	Escherichia coli O157:H7	ELISA
<u>TJP1 Antibody</u>	CSB-PA023573OA01HU	Human	ELISA, IF
<u>TJP1 Antibody</u>	CSB-PA023573HA01HU	Human	ELISA, IF
<u>TJP2 Antibody</u>	CSB-PA880071NA01HU	Human	ELISA, WB, IHC, IF
<u>TK Antibody</u>	CSB-PA023577HA01HWY	Human herpesvirus 1	ELISA
<u>TK1 Antibody</u>	CSB-PA023577HA01HU	Human	ELISA, WB, IHC
<u>TKFC Antibody</u>	CSB-PA659468LA01HU	Human	ELISA, IHC, IF
<u>TKTL2 Antibody</u>	CSB-PA872424LA01HU	Human	ELISA, IHC
<u>TLE4 Antibody</u>	CSB-PA523988LA01HU	Human	ELISA, IF
<u>TLK2 Antibody</u>	CSB-PA801231LA01HU	Human	ELISA, IHC, IF
<u>TLL1 Antibody</u>	CSB-PA023595LA01HU	Human	ELISA, IHC, IF
<u>TLN1 Antibody</u>	CSB-PA05189A0Rb	Human	ELISA, IHC, IF
<u>TLN2 Antibody</u>	CSB-PA897282HA01HU	Human	ELISA, WB, IHC
<u>TLR1 Antibody</u>	CSB-PA023599LA01HU	Human	ELISA, IHC, IF
<u>TLR2 Antibody</u>	CSB-PA023601LA01HU	Human	ELISA



Product Name	Code	Species Reactivity	Applications
<u>TLR3 Antibody</u>	CSB-PA023602LA01HU	Human	ELISA, IHC, IF
<u>TLR5 Antibody</u>	CSB-PA023604LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TLR6 Antibody</u>	CSB-PA023605LA01HU	Human	ELISA, WB, IHC, IF
<u>Tlr7 Antibody</u>	CSB-PA023606EA01MO	Mouse, Human	ELISA, WB
<u>TLR7 Antibody</u>	CSB-PA16789A0Rb	Human	ELISA, IF
<u>TLR9 Antibody</u>	CSB-PA023608LA01HU	Human	ELISA, IHC, IF
<u>TM2D3 Antibody</u>	CSB-PA863106LA01HU	Human	ELISA, IHC, IF
<u>TM4SF1 Antibody</u>	CSB-PA023615LA01HU	Human	ELISA, IHC, IF
<u>TM4SF20 Antibody</u>	CSB-PA709608LA01HU	Human	ELISA, IHC, IF
<u>TM4SF5 Antibody</u>	CSB-PA023620LA01HU	Human	ELISA, IHC
<u>TM6SF2 Antibody</u>	CSB-PA023622LA01HU	Human	ELISA, IHC, IF
<u>TM7SF3 Antibody</u>	CSB-PA023624LA01HU	Human	ELISA, IHC, IF
<u>TMBIM1 Antibody</u>	CSB-PA853376LA01HU	Human, Mouse	ELISA, WB, IF
<u>TMC2 Antibody</u>	CSB-PA855056LA01HU	Human	ELISA, IF
<u>TMCC1 Antibody</u>	CSB-PA023641LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TMCO1 Antibody</u>	CSB-PA891784LA01HU	Human	ELISA, IHC, IF
<u>TMCO5A Antibody</u>	CSB-PA843283LA01HU	Human	ELISA, WB, IHC, IF
<u>TMED2 Antibody</u>	CSB-PA05049A0Rb	Human, Mouse	ELISA, WB, IHC
<u>TMED3 Antibody</u>	CSB-PA897107LA01HU	Human	ELISA, WB, IHC
<u>TMED6 Antibody</u>	CSB-PA840991LA01HU	Human	ELISA, IHC
<u>TMED9 Antibody</u>	CSB-PA04804A0Rb	Human, Mouse	ELISA, WB
<u>TMEM100 Antibody</u>	CSB-PA023668LA01HU	Human	ELISA, IF
<u>TMEM106A Antibody</u>	CSB-PA023673LA01HU	Human	ELISA, IF
<u>TMEM109 Antibody</u>	CSB-PA03899A0Rb	Human	ELISA, IHC, IF
<u>TMEM11 Antibody</u>	CSB-PA023679LA01HU	Human	ELISA, IHC
<u>TMEM115 Antibody</u>	CSB-PA023683LA01HU	Human	ELISA, IHC, IF
<u>TMEM120B Antibody</u>	CSB-PA023688LA01HU	Human	ELISA, IHC, IF
<u>TMEM127 Antibody</u>	CSB-PA023694LA01HU	Human	ELISA, IHC, IF
<u>TMEM131L Antibody</u>	CSB-PA012218NA01HU	Human	ELISA, IF
<u>TMEM132D Antibody</u>	CSB-PA023702LA01HU	Human	ELISA, IF
<u>TMEM139 Antibody</u>	CSB-PA023709LA01HU	Human	ELISA, IHC, IF
<u>TMEM140 Antibody</u>	CSB-PA023710LA01HU	Human	ELISA, IHC, IF
<u>TMEM141 Antibody</u>	CSB-PA023711EA01HU	Human	ELISA, WB, IHC
<u>TMEM143 Antibody</u>	CSB-PA023713LA01HU	Human	ELISA, IHC, IF
<u>TMEM14A Antibody</u>	CSB-PA023720LA01HU	Human	ELISA, WB, IHC, IF
<u>TMEM14B Antibody</u>	CSB-PA023721EA01HU	Human	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>TMEM150A Antibody</u>	CSB-PA774796LA01HU	Human	ELISA, IF
<u>TMEM154 Antibody</u>	CSB-PA023729LA01HU	Human	ELISA, WB, IHC
<u>TMEM159 Antibody</u>	CSB-PA850262LA01HU	Human	ELISA, IF
<u>TMEM161A Antibody</u>	CSB-PA023735LA01HU	Human	ELISA, WB, IHC, IF
<u>TMEM161B Antibody</u>	CSB-PA023736LA01HU	Human	ELISA, IHC, IF
<u>TMEM165 Antibody</u>	CSB-PA023739LA01HU	Human	ELISA, WB, IHC, IF
<u>TMEM173 Antibody</u>	CSB-PA023754LA01HU	Human	ELISA, WB, IHC, IF
<u>TMEM174 Antibody</u>	CSB-PA023755LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TMEM176B Antibody</u>	CSB-PA023758LA01HU	Human	ELISA, WB, IHC, IF
<u>TMEM177 Antibody</u>	CSB-PA023759LA01HU	Human	ELISA, WB, IHC, IF
<u>TMEM184B Antibody</u>	CSB-PA023770LA01HU	Human	ELISA, IHC, IF
<u>TMEM184C Antibody</u>	CSB-PA023771LA01HU	Human	ELISA, IHC, IF
<u>TMEM186 Antibody</u>	CSB-PA023774LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TMEM19 Antibody</u>	CSB-PA023779LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TMEM192 Antibody</u>	CSB-PA023784LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>TMEM196 Antibody</u>	CSB-PA023788LA01HU	Human	ELISA
<u>TMEM2 Antibody</u>	CSB-PA23279A0Rb	Human	ELISA, IHC, IF
<u>TMEM205 Antibody</u>	CSB-PA023799LA01HU	Human	ELISA, IHC
<u>TMEM207 Antibody</u>	CSB-PA023801LA01HU	Human	ELISA, IHC, IF
<u>TMEM209 Antibody</u>	CSB-PA023803LA01HU	Human	ELISA, IHC, IF
<u>TMEM213 Antibody</u>	CSB-PA023807LA01HU	Human	ELISA, IHC, IF
<u>TMEM214 Antibody</u>	CSB-PA023808LA01HU	Human	ELISA, WB, IP
<u>TMEM217 Antibody</u>	CSB-PA023811LA01HU	Human	ELISA, WB, IF
<u>TMEM222 Antibody</u>	CSB-PA023817LA01HU	Human	ELISA, WB, IHC, IF
<u>TMEM231 Antibody</u>	CSB-PA862028LA01HU	Human	ELISA, IF
<u>TMEM237 Antibody</u>	CSB-PA822287OA01HU	Human	ELISA, WB
<u>TMEM240 Antibody</u>	CSB-PA719389LA01HU	Human	ELISA, IF
<u>TMEM243 Antibody</u>	CSB-PA874810LA01HU	Human	ELISA, IHC, IF
<u>TMEM25 Antibody</u>	CSB-PA023821LA01HU	Human	ELISA
<u>TMEM255A Antibody</u>	CSB-PA708482LA01HU	Human, Mouse	ELISA, WB, IF
<u>TMEM255B Antibody</u>	CSB-PA819893LA01HU	Human	ELISA, IHC, IF
<u>TMEM256 Antibody</u>	CSB-PA847619LA01HU	Human	ELISA, IF



Product Name	Code	Species Reactivity	Applications
<u>TMEM26 Antibody</u>	CSB-PA023822LA01HU	Human	ELISA, IF
<u>TMEM267 Antibody</u>	CSB-PA610832LA01HU	Human	ELISA, IF
<u>TMEM27 Antibody</u>	CSB-PA881018EA01HU	Human	ELISA, IHC, IF
<u>TMEM270 Antibody</u>	CSB-PA744250LA01HU	Human	ELISA, IHC, IF
<u>TMEM31 Antibody</u>	CSB-PA023828LA01HU	Human	ELISA, IHC, IF
<u>TMEM33 Antibody</u>	CSB-PA023829LA01HU	Human	ELISA, IHC, IF
<u>TMEM38B Antibody</u>	CSB-PA878900LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TMEM40 Antibody</u>	CSB-PA023836LA01HU	Human	ELISA, IHC, IF
<u>TMEM41A Antibody</u>	CSB-PA023837LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TMEM43 Antibody</u>	CSB-PA023840LA01HU	Human, Mouse	ELISA, WB, IHC
<u>TMEM44 Antibody</u>	CSB-PA023841HA01HU	Human	ELISA, IHC, IF
<u>TMEM45A Antibody</u>	CSB-PA023842LA01HU	Human	ELISA, IHC, IF
<u>TMEM47 Antibody</u>	CSB-PA023845LA01HU	Human	ELISA, IHC, IF
<u>TMEM52B Antibody</u>	CSB-PA680179LA01HU	Human	ELISA, IHC
<u>TMEM53 Antibody</u>	CSB-PA023853LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TMEM56 Antibody</u>	CSB-PA023857LA01HU	Human	ELISA, IHC, IF
<u>TMEM59 Antibody</u>	CSB-PA023860LA01HU	Human	ELISA, WB, IF
<u>TMEM59L Antibody</u>	CSB-PA023861LA01HU	Human	ELISA, WB, IHC, IF
<u>TMEM61 Antibody</u>	CSB-PA023863LA01HU	Human	ELISA, IHC, IF
<u>TMEM63A Antibody</u>	CSB-PA023865LA01HU	Human	ELISA, IF
<u>TMEM64 Antibody</u>	CSB-PA023868LA01HU	Human, Mouse	ELISA, WB, IHC
<u>TMEM65 Antibody</u>	CSB-PA023869LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TMEM67 Antibody</u>	CSB-PA689243LA01HU	Human	ELISA, IHC, IF
<u>TMEM68 Antibody</u>	CSB-PA023872LA01HU	Human, Mouse	ELISA, WB, IF
<u>TMEM71 Antibody</u>	CSB-PA023875LA01HU	Human	ELISA, WB, IHC
<u>TMEM79 Antibody</u>	CSB-PA023881LA01HU	Human	ELISA, IHC
<u>TMEM82 Antibody</u>	CSB-PA023885LA01HU	Human	ELISA, IHC, IF
<u>TMEM87B Antibody</u>	CSB-PA023891LA01HU	Human	ELISA, IHC
<u>TMEM88 Antibody</u>	CSB-PA023892LA01HU	Human	ELISA, IF
<u>TMEM9 Antibody</u>	CSB-PA023894LA01HU	Human	ELISA, IHC
<u>TMEM91 Antibody</u>	CSB-PA023896LA01HU	Human	ELISA, WB, IHC, IF
<u>TMEM95 Antibody</u>	CSB-PA023899LA01HU	Human	ELISA, IHC, IF
<u>TMEM98 Antibody</u>	CSB-PA023901LA01HU	Human	ELISA
<u>TMF1 Antibody</u>	CSB-PA023904LA01HU	Human	ELISA, WB, IHC
<u>TMIGD2 Antibody</u>	CSB-PA850263LA01HU	Human	ELISA, IHC, IF
<u>tmk Antibody</u>	CSB-PA007226HA01BQG	Bacillus anthracis	ELISA



Product Name	Code	Species Reactivity	Applications
<u>TMLHE Antibody</u>	CSB-PA023908LA01HU	Human	ELISA, WB, IHC
<u>TMPRSS11B Antibody</u>	CSB-PA771428LA01HU	Human	ELISA, IHC
<u>TMPRSS11D Antibody</u>	CSB-PA023918LA01HU	Human	ELISA, IHC
<u>Tmprss15 Antibody</u>	CSB-PA018824LA01MO	Mouse, Human, Rat	ELISA, WB
<u>Tmprss4 Antibody</u>	CSB-PA844988LA01MO	Mouse	ELISA
<u>TMPRSS6 Antibody</u>	CSB-PA808531LA01HU	Human	ELISA, IHC
<u>TMPRSS7 Antibody</u>	CSB-PA745764LA01HU	Human	ELISA, IHC
<u>TMSB4X Antibody</u>	CSB-PA023934LA01HU	Human	ELISA, IHC
<u>TMSB4Y Antibody</u>	CSB-PA023935LA01HU	Human	ELISA
<u>TMTC1 Antibody</u>	CSB-PA808537LA01HU	Human, Mouse	ELISA, WB, IF
<u>TMTC4 Antibody</u>	CSB-PA713124LA01HU	Human	ELISA, IHC
<u>TMX2 Antibody</u>	CSB-PA023951LA01HU	Human	ELISA, IHC, IF
<u>TMX3 Antibody</u>	CSB-PA846645LA01HU	Human	ELISA, IHC
<u>TNC Antibody</u>	CSB-PA11504A0Rb	Human	ELISA, IHC
<u>TNC Antibody</u>	CSB-PA11509A0Rb	Human	ELISA, IHC
<u>TNF Antibody</u>	CSB-PA023955DA01SH	Sheep	ELISA, WB
<u>TNF Antibody</u>	CSB-PA07427A0Rb	Human	ELISA, WB, IHC, IF
<u>TNF Antibody</u>	CSB-PA08474A0Rb	Bovine	ELISA, WB
<u>TNFAIP1 Antibody</u>	CSB-PA622686LA01HU	Human	ELISA, IHC, IF
<u>TNFAIP2 Antibody</u>	CSB-PA023957LA01HU	Human, Mouse	ELISA, WB
<u>TNFAIP3 Antibody</u>	CSB-PA023958LA01HU	Human, Mouse	ELISA, WB, IF
<u>TNFAIP6 Antibody</u>	CSB-PA023959LA01HU	Human	ELISA, IHC
<u>TNFAIP8L1 Antibody</u>	CSB-PA837859LA01HU	Human	ELISA, IF
<u>TNFAIP8L2 Antibody</u>	CSB-PA754193LA01HU	Human	ELISA, IF
<u>TNFRSF10A Antibody</u>	CSB-PA023964LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>TNFRSF10B Antibody</u>	CSB-PA08165A0Rb	Human	ELISA, IF
<u>TNFRSF10C Antibody</u>	CSB-PA08189A0Rb	Human	ELISA, IF
<u>TNFRSF10D Antibody</u>	CSB-PA08209A0Rb	Human	ELISA, IHC
<u>TNFRSF11A Antibody</u>	CSB-PA08467A0Rb	Human	ELISA, WB, IF
<u>TNFRSF11B Antibody</u>	CSB-PA08229A0Rb	Human	ELISA, IHC
<u>TNFRSF13C Antibody</u>	CSB-PA853495LA01HU	Human	ELISA
<u>TNFRSF18 Antibody</u>	CSB-PA023975LA01HU	Human	ELISA, IHC, IF
<u>TNFRSF19 Antibody</u>	CSB-PA878878LA01HU	Human	ELISA, WB, IF
<u>TNFRSF1A Antibody</u>	CSB-PA023977LA01HU	Human	ELISA, WB, IHC, IF
<u>TNFRSF25 Antibody</u>	CSB-PA839000LA01HU	Human	ELISA, WB
<u>Tnfrsf4 Antibody</u>	CSB-PA023981HA01MO	Mouse	ELISA
<u>TNFRSF9 Antibody</u>	CSB-PA08109A0Rb	Human	ELISA, WB, IF
<u>TNFSF10 Antibody</u>	CSB-PA07475A0Rb	Human	ELISA, IHC, IF
<u>Tnfsf11 Antibody</u>	CSB-PA023986LA01MO	Mouse	ELISA, IF



Product Name	Code	Species Reactivity	Applications
<u>Tnfsf11 Antibody</u>	CSB-PA07964A0Rb	Mouse	ELISA
<u>TNFSF11 Antibody</u>	CSB-PA023986LA01HU	Human	ELISA, IHC, IF
<u>TNFSF12 Antibody</u>	CSB-PA023987LA01HU	Human	ELISA, IHC, IF
<u>TNFSF13B Antibody</u>	CSB-PA05767A0Rb	Human	ELISA
<u>TNFSF14 Antibody</u>	CSB-PA06825A0Rb	Human	ELISA, WB, IF
<u>TNFSF9 Antibody</u>	CSB-PA08047A0Rb	Human	ELISA, IHC, IF
<u>TNIK Antibody</u>	CSB-PA023998LA01HU	Human	ELISA, IHC, IF
<u>TNIP1 Antibody</u>	CSB-PA619959LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>TNKS Antibody</u>	CSB-PA024004LA01HU	Human	ELISA, WB, IHC, IF
<u>TNKS1BP1 Antibody</u>	CSB-PA871634LA01HU	Human	ELISA, IHC, IF
<u>TNKS2 Antibody</u>	CSB-PA867136LA01HU	Human	ELISA, IHC, IF
<u>TNMD Antibody</u>	CSB-PA024007LA01HU	Human, Mouse	ELISA, WB, IF
<u>TNN Antibody</u>	CSB-PA892498LA01HU	Human	ELISA, IHC
<u>TNNC1 Antibody</u>	CSB-PA024009LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>TNNI1 Antibody</u>	CSB-PA024011LA01HU	Human	ELISA, IF
<u>TNNI3 Antibody</u>	CSB-PA07936A0Rb	Human	ELISA
<u>TNNI3 Antibody</u>	CSB-PA07939A0Rb	Human, Rat	ELISA, IHC, IF, IP
<u>Tnni3 Antibody</u>	CSB-PA341785LA01MO	Mouse, Rat	ELISA, WB
<u>Tnnt2 Antibody</u>	CSB-PA024016HA01MO	Mouse	ELISA
<u>TNNT2 Antibody</u>	CSB-PA10047A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>TNP1 Antibody</u>	CSB-PA024018LA01HU	Human	ELISA, IF
<u>TNP2 Antibody</u>	CSB-PA024019LA01HU	Human	ELISA, IHC, IF
<u>TNPO3 Antibody</u>	CSB-PA757558LA01HU	Human	ELISA, IF
<u>TNR Antibody</u>	CSB-PA340307LA01HU	Human	ELISA, IHC
<u>TNRC6A Antibody</u>	CSB-PA839859LA01HU	Human	ELISA, IHC, IF
<u>TNRC6B Antibody</u>	CSB-PA892489LA01HU	Human	ELISA, IF
<u>TNRC6C Antibody</u>	CSB-PA862070LA01HU	Human	ELISA, WB, IHC
<u>TNS4 Antibody</u>	CSB-PA810319LA01HU	Human	ELISA, WB, IHC
<u>TOB1 Antibody</u>	CSB-PA024037LA01HU	Human	ELISA, WB, IHC, IF
<u>TOLLIP Antibody</u>	CSB-PA875640LA01HU	Human	ELISA, IHC
<u>TOM1 Antibody</u>	CSB-PA024043LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TOM1L1 Antibody</u>	CSB-PA024044LA01HU	Human	ELISA, IHC, IF
<u>TONSL Antibody</u>	CSB-PA836224LA01HU	Human	ELISA, IF
<u>TOP1 Antibody</u>	CSB-PA024056LA01HU	Human	ELISA, IF
<u>TOP1MT Antibody</u>	CSB-PA853372LA01HU	Human	ELISA, IHC
<u>TOP2B Antibody</u>	CSB-PA717547LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>TOR1A Antibody</u>	CSB-PA024067LA01HU	Human	ELISA, IHC, IF
<u>TOX Antibody</u>	CSB-PA024073LA01HU	Human	ELISA, WB, IHC
<u>TOX2 Antibody</u>	CSB-PA024074LA01HU	Human	ELISA, WB
<u>TOX3 Antibody</u>	CSB-PA024075LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>toxA Antibody</u>	CSB-PA324524HA01DUN	Peptoclostridium difficile	ELISA
<u>toxB Antibody</u>	CSB-PA321589HA01DUN	Peptoclostridium difficile	ELISA
<u>Toxin zeta Antibody</u>	CSB-PA363578LA01SMF	Streptococcus agalactiae	ELISA
<u>Tp53 Antibody</u>	CSB-PA024077LA01MO	Mouse	ELISA
<u>TP53 Antibody</u>	CSB-PA07889A0Rb	Human, Mouse, Rat	ELISA, WB, IHC, ChIP
<u>Tp53 Antibody</u>	CSB-PA15509A0Rb	Rat, Mouse	ELISA, WB
<u>TP53I11 Antibody</u>	CSB-PA024082LA01HU	Human	ELISA, IF
<u>TP53INP1 Antibody</u>	CSB-PA822174LA01HU	Human	ELISA, IHC, IF
<u>TP53INP2 Antibody</u>	CSB-PA812888LA01HU	Human	ELISA, IF
<u>TP53TG5 Antibody</u>	CSB-PA897086LA01HU	Human	ELISA, IHC
<u>tp63 Antibody</u>	CSB-PA23919A0Rb	Zebrafish	ELISA, WB
<u>TP73 Antibody</u>	CSB-PA024092OA01HU	Human	ELISA, WB, IHC
<u>Tpbg Antibody</u>	CSB-PA024093LA01RA	Rat	ELISA
<u>TPD52L2 Antibody</u>	CSB-PA024098LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>Tph1 Antibody</u>	CSB-PA024100HA01MO	Mouse	ELISA, IF
<u>TPH1 Antibody</u>	CSB-PA024100LA01HU	Human, Mouse	ELISA, WB, IF
<u>TPH2 Antibody</u>	CSB-PA809000LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>TPI1 Antibody</u>	CSB-PA13499A0Rb	Human, Mouse	ELISA, WB, IHC
<u>TPK1 Antibody</u>	CSB-PA024103EA01HU	Human	ELISA, IF
<u>TPM2 Antibody</u>	CSB-PA13387A0Rb	Human	ELISA
<u>TPM3 Antibody</u>	CSB-PA024106LA01HU	Human	ELISA, IHC, IF
<u>TPMT Antibody</u>	CSB-PA00585A0Rb	Human	ELISA, WB, IHC, IP
<u>TPO Antibody</u>	CSB-PA12279A0Rb	Human	ELISA, IHC
<u>TPP1 Antibody</u>	CSB-PA024113LA01HU	Human	ELISA, IHC
<u>TPPP2 Antibody</u>	CSB-PA024116LA01HU	Human	ELISA, WB, IHC, IF
<u>TPPP3 Antibody</u>	CSB-PA024117LA01HU	Human	ELISA, IHC
<u>TPRKB Antibody</u>	CSB-PA896710LA01HU	Human	ELISA, WB, IHC
<u>Tpsab1 Antibody</u>	CSB-PA024127LA01MO	Mouse	ELISA
<u>TPSB2 Antibody</u>	CSB-PA024128LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>Tpsb2 Antibody</u>	CSB-PA024128LA01MO	Mouse, Rat	ELISA, WB
<u>Tpsb2 Antibody</u>	CSB-PA024128LA01RA	Rat	ELISA
<u>TPST1 Antibody</u>	CSB-PA024132LA01HU	Human, Mouse, Rat	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>TPST2 Antibody</u>	CSB-PA024133LA01HU	Human	ELISA, IHC
<u>Tpt1 Antibody</u>	CSB-PA024134HA01MO	Mouse	ELISA, IHC
<u>TPT1 Antibody</u>	CSB-PA024134LA01HU	Human	ELISA, IHC, IF
<u>TPT1P8 Antibody</u>	CSB-PA884505LA01HU	Human	ELISA, IF
<u>TPTE Antibody</u>	CSB-PA024136LA01HU	Human	ELISA, WB, IF
<u>TPTE2 Antibody</u>	CSB-PA757850LA01HU	Human	ELISA, IHC
<u>TPX2 Antibody</u>	CSB-PA892151LA01HU	Human	ELISA, WB, IHC
<u>TRABD2B Antibody</u>	CSB-PA406235LA01HU	Human	ELISA, IHC, IF
<u>TRAC Antibody</u>	CSB-PA024144LA01HU	Human	ELISA
<u>TRADD Antibody</u>	CSB-PA621879EA01HU	Human, Rat	ELISA, WB, IHC, IF, IP
<u>TRAF2 Antibody</u>	CSB-PA024147LA01HU	Human	ELISA, IHC, IF
<u>TRAF3 Antibody</u>	CSB-PA024148LA01HU	Human	ELISA, WB, IF
<u>TRAF3IP3 Antibody</u>	CSB-PA896865LA01HU	Human, Mouse, Rat	ELISA, WB, IF
<u>TRAF5 Antibody</u>	CSB-PA024153LA01HU	Human	ELISA, IHC
<u>TRAF6 Antibody</u>	CSB-PA024154LA01HU	Human	ELISA, IHC, IF
<u>TRAF7 Antibody</u>	CSB-PA757671LA01HU	Human	ELISA, IHC, IF
<u>TRAFD1 Antibody</u>	CSB-PA024156LA01HU	Human	ELISA
<u>TRAK1 Antibody</u>	CSB-PA891465LA01HU	Human	ELISA, IHC
<u>TRAM2 Antibody</u>	CSB-PA624011LA01HU	Human	ELISA, IHC, IF
<u>Transposase for transposon Tn552 Antibody</u>	CSB-PA323588LA01FKZ	Staphylococcus aureus	ELISA
<u>TRAP1 Antibody</u>	CSB-PA03604A0Rb	Human, Mouse	ELISA, WB, IF, IP
<u>TRAPPC11 Antibody</u>	CSB-PA801799LA01HU	Human	ELISA, IF
<u>TRAPPC2L Antibody</u>	CSB-PA892461LA01HU	Human, Mouse	ELISA, WB, IHC
<u>TRAPPC3 Antibody</u>	CSB-PA024229LA01HU	Human	ELISA, IHC
<u>TRAPPC6A Antibody</u>	CSB-PA024232LA01HU	Human	ELISA, IHC
<u>TRAPPC6B Antibody</u>	CSB-PA803112LA01HU	Human, Mouse	ELISA, WB, IHC
<u>TRAPPC8 Antibody</u>	CSB-PA27999A0Rb	Human	ELISA, IHC
<u>TRAPPC9 Antibody</u>	CSB-PA842741LA01HU	Human	ELISA, IHC, IF
<u>TRAT1 Antibody</u>	CSB-PA764808LA01HU	Human	ELISA, IHC
<u>TRBC1 Antibody</u>	CSB-PA024291LA01HU	Human	ELISA, IHC, IF
<u>TRBV7-9 Antibody</u>	CSB-PA361265LA01HU	Human	ELISA, WB, IHC
<u>TRDMT1 Antibody</u>	CSB-PA675694LA01HU	Human, Rat	ELISA, WB, IHC
<u>TRDN Antibody</u>	CSB-PA024398LA01HU	Human	ELISA, IF
<u>TREH Antibody</u>	CSB-PA024402LA01HU	Human	ELISA, IF
<u>Trem2 Antibody</u>	CSB-PA024405LA01MO	Mouse, Human	ELISA, WB
<u>TREM2 Antibody</u>	CSB-PA024405LA01HU	Human	ELISA, IHC, IF
<u>TREML1 Antibody</u>	CSB-PA771477LA01HU	Human	ELISA
<u>TREML2 Antibody</u>	CSB-PA729201LA01HU	Human	ELISA, IHC
<u>TRHR Antibody</u>	CSB-PA024441LA01HU	Human	ELISA, IF
<u>TRIB1 Antibody</u>	CSB-PA836296LA01HU	Human	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>TRIB2 Antibody</u>	CSB-PA838778LA01HU	Human	ELISA, IHC, IF
<u>TRIB3 Antibody</u>	CSB-PA847220LA01HU	Human	ELISA, IHC
<u>TRIM11 Antibody</u>	CSB-PA856931LA01HU	Human	ELISA, IHC
<u>TRIM15 Antibody</u>	CSB-PA887181LA01HU	Human	ELISA, WB, IHC
<u>TRIM16 Antibody</u>	CSB-PA024453LA01HU	Human	ELISA
<u>TRIM17 Antibody</u>	CSB-PA897559LA01HU	Human	ELISA, WB, IHC
<u>Trim21 Antibody</u>	CSB-PA733816LA01MO	Mouse	ELISA
<u>TRIM21 Antibody</u>	CSB-PA024457LA01HU	Human	ELISA, WB, IHC
<u>TRIM24 Antibody</u>	CSB-PA024460LA01HU	Human	ELISA
<u>TRIM25 Antibody</u>	CSB-PA617909LA01HU	Human	ELISA, WB, IHC, IF
<u>TRIM26 Antibody</u>	CSB-PA619638LA01HU	Human	ELISA, IHC
<u>TRIM28 Antibody</u>	CSB-PA05254A0Rb	Human	ELISA, WB, IHC, IP
<u>TRIM3 Antibody</u>	CSB-PA024467LA01HU	Human	ELISA, IHC
<u>TRIM33 Antibody</u>	CSB-PA891796LA01HU	Human	ELISA, IHC, IF
<u>TRIM35 Antibody</u>	CSB-PA892166LA01HU	Human	ELISA, IHC, IF
<u>TRIM36 Antibody</u>	CSB-PA868285LA01HU	Human	ELISA, IHC
<u>TRIM37 Antibody</u>	CSB-PA024474LA01HU	Human	ELISA, IHC, IF
<u>TRIM38 Antibody</u>	CSB-PA024475LA01HU	Human	ELISA, WB
<u>TRIM39 Antibody</u>	CSB-PA881029LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TRIM4 Antibody</u>	CSB-PA866336LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>TRIM43 Antibody</u>	CSB-PA836189LA01HU	Human, Mouse	ELISA, WB
<u>TRIM45 Antibody</u>	CSB-PA862044LA01HU	Human	ELISA, IHC
<u>TRIM48 Antibody</u>	CSB-PA818248LA01HU	Human	ELISA
<u>TRIM5 Antibody</u>	CSB-PA874861LA01HU	Human, Mouse	ELISA, WB, IHC
<u>TRIM51 Antibody</u>	CSB-PA871567LA01HU	Human	ELISA, WB, IHC
<u>TRIM54 Antibody</u>	CSB-PA861180LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>TRIM55 Antibody</u>	CSB-PA880142LA01HU	Human, Mouse	ELISA, WB, IHC
<u>TRIM59 Antibody</u>	CSB-PA808999LA01HU	Human	ELISA, IHC
<u>TRIM65 Antibody</u>	CSB-PA750884LA01HU	Human	ELISA, WB, IHC, IF
<u>TRIM71 Antibody</u>	CSB-PA652403LA01HU	Human	ELISA, IHC
<u>TRIM74 Antibody</u>	CSB-PA801235LA01HU	Human	ELISA, IHC
<u>TRIM8 Antibody</u>	CSB-PA871629LA01HU	Human	ELISA, IHC, IF
<u>TRIM9 Antibody</u>	CSB-PA880153LA01HU	Human	ELISA, IHC
<u>TRIML1 Antibody</u>	CSB-PA822741LA01HU	Human, Mouse	ELISA, WB
<u>TRIP10 Antibody</u>	CSB-PA624029LA01HU	Human	ELISA, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>TRIP13 Antibody</u>	CSB-PA623008LA01HU	Human	ELISA, WB, IHC
<u>TRIP4 Antibody</u>	CSB-PA618991LA01HU	Human	ELISA, WB, IHC, IF
<u>TRIT1 Antibody</u>	CSB-PA867142LA01HU	Human	ELISA, WB, IHC
<u>TRMT112 Antibody</u>	CSB-PA883416LA01HU	Human	ELISA, IHC
<u>TRMT6 Antibody</u>	CSB-PA891716LA01HU	Human	ELISA, WB, IHC
<u>TRNAU1AP Antibody</u>	CSB-PA882127LA01HU	Human	ELISA, WB, IHC
<u>TRNP1 Antibody</u>	CSB-PA761223LA01HU	Human	ELISA, IHC
<u>TRNT1 Antibody</u>	CSB-PA836282LA01HU	Human	ELISA, WB
<u>TRO Antibody</u>	CSB-PA025069LA01HU	Human	ELISA, IF
<u>Tropomyosin Antibody</u>	CSB-PA889027HA01DRR	Charybdis feriata	ELISA
<u>TROVE2 Antibody</u>	CSB-PA13749A0Rb	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>TRPC1 Antibody</u>	CSB-PA025075LA01HU	Human	ELISA, IHC
<u>TRPM1 Antibody</u>	CSB-PA773791LA01HU	Human	ELISA, IF
<u>TRPM2 Antibody</u>	CSB-PA007779LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TRPM3 Antibody</u>	CSB-PA884520LA01HU	Human	ELISA, IF
<u>TRPM8 Antibody</u>	CSB-PA768757LA01HU	Human	ELISA, IHC, IF
<u>TRPV1 Antibody</u>	CSB-PA822774LA01HU	Human	ELISA, IHC, IF
<u>TRPV2 Antibody</u>	CSB-PA897571LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TRPV3 Antibody</u>	CSB-PA854113LA01HU	Human	ELISA, IF
<u>TRPV6 Antibody</u>	CSB-PA861992LA01HU	Human	ELISA, IHC
<u>TRUB1 Antibody</u>	CSB-PA851978LA01HU	Human	ELISA
<u>trxA Antibody</u>	CSB-PA18667A0Rb	Escherichia coli	ELISA, WB
<u>trxC Antibody</u>	CSB-PA314858HA01ENV	Escherichia coli	ELISA
<u>Trypsin Antibody</u>	CSB-PA018814YA01PI	Pig	ELISA, WB
<u>Trypsin Antibody</u>	CSB-PA18879A0Rb	Pig	ELISA
<u>Trypsin inhibitor DE-3 Antibody</u>	CSB-PA362840HA01EGO	Erythrina caffra	ELISA
<u>Trypsin inhibitor DE-3 Antibody</u>	CSB-PA362840YA01EGO	Erythrina caffra	ELISA
<u>TSACC Antibody</u>	CSB-PA842631LA01HU	Human	ELISA, IHC
<u>TSBP1 Antibody</u>	CSB-PA713105HA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TSC1 Antibody</u>	CSB-PA849764LA01HU	Human	ELISA, IHC, IF
<u>TSC2 Antibody</u>	CSB-PA343033LA01HU	Human, Mouse	ELISA, WB, IF
<u>TSC22D1 Antibody</u>	CSB-PA623011LA01HU	Human	ELISA
<u>TSC22D3 Antibody</u>	CSB-PA860343LA01HU	Human	ELISA, WB, IHC
<u>TSEN15 Antibody</u>	CSB-PA840989LA01HU	Human	ELISA, IHC
<u>TSGA10 Antibody</u>	CSB-PA863172LA01HU	Human	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>TSGA10IP Antibody</u>	CSB-PA660950LA01HU	Human	ELISA
<u>TSGA13 Antibody</u>	CSB-PA842737LA01HU	Human	ELISA, WB, IHC
<u>TSHB Antibody</u>	CSB-PA11769A0Rb	Human	ELISA
<u>TSHR Antibody</u>	CSB-PA11189A0Rb	Human	ELISA, IF
<u>TSKU Antibody</u>	CSB-PA025140LA01HU	Human	ELISA, WB
<u>TSLP Antibody</u>	CSB-PA025141LA01HU	Human, Mouse	ELISA, WB, IHC
<u>TSNARE1 Antibody</u>	CSB-PA822269LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TSNAX Antibody</u>	CSB-PA858718LA01HU	Human	ELISA, WB, IHC, IF
<u>TSPAN1 Antibody</u>	CSB-PA025146LA01HU	Human	ELISA, IHC, IF
<u>TSPAN12 Antibody</u>	CSB-PA025149LA01HU	Human	ELISA, IHC, IF
<u>TSPAN14 Antibody</u>	CSB-PA847712LA01HU	Human	ELISA, IHC, IF
<u>TSPAN15 Antibody</u>	CSB-PA025152LA01HU	Human	ELISA, IF
<u>TSPAN16 Antibody</u>	CSB-PA892335LA01HU	Human	ELISA, IF
<u>TSPAN18 Antibody</u>	CSB-PA853503LA01HU	Human	ELISA, IF
<u>TSPAN31 Antibody</u>	CSB-PA622522LA01HU	Human	ELISA, IHC, IF
<u>TSPAN32 Antibody</u>	CSB-PA842746LA01HU	Human	ELISA, IHC, IF
<u>TSPAN33 Antibody</u>	CSB-PA768209LA01HU	Human	ELISA, IHC, IF
<u>TSPAN5 Antibody</u>	CSB-PA025163LA01HU	Human	ELISA, WB, IF
<u>TSPAN6 Antibody</u>	CSB-PA025164LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TSPAN7 Antibody</u>	CSB-PA025165LA01HU	Human	ELISA, IHC, IF
<u>TSPAN8 Antibody</u>	CSB-PA025166LA01HU	Human	ELISA, IHC, IF
<u>TSPAN9 Antibody</u>	CSB-PA025167LA01HU	Human	ELISA, IHC, IF
<u>TSPEAR Antibody</u>	CSB-PA848817LA01HU	Human	ELISA, WB
<u>TSPO Antibody</u>	CSB-PA536576LA01HU	Human	ELISA, IHC
<u>TSR1 Antibody</u>	CSB-PA651478LA01HU	Human	ELISA, IF
<u>TSR2 Antibody</u>	CSB-PA836166HA01HU	Human	ELISA, WB, IHC, IF
<u>TST Antibody</u>	CSB-PA620980LA01HU	Human, Rat	ELISA, WB, IHC
<u>TSTA3 Antibody</u>	CSB-PA00935A0Rb	Human	ELISA, WB, IHC, IF
<u>TTC1 Antibody</u>	CSB-PA01104A0Rb	Human	ELISA, WB
<u>TTC12 Antibody</u>	CSB-PA862040LA01HU	Human, Mouse	ELISA, WB, IHC
<u>TTC16 Antibody</u>	CSB-PA850870LA01HU	Human	ELISA, IHC
<u>TTC19 Antibody</u>	CSB-PA724947LA01HU	Human	ELISA, IHC, IF
<u>TTC26 Antibody</u>	CSB-PA025219LA01HU	Human, Mouse	ELISA, WB
<u>TTC27 Antibody</u>	CSB-PA757569LA01HU	Human	ELISA, WB, IHC
<u>TTC30A Antibody</u>	CSB-PA803138LA01HU	Human, Mouse	ELISA, WB, IHC
<u>TTC30B Antibody</u>	CSB-PA818703LA01HU	Human	ELISA, WB, IHC
<u>TTC32 Antibody</u>	CSB-PA686303LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>TTC33 Antibody</u>	CSB-PA747492LA01HU	Human	ELISA, IHC
<u>TTC38 Antibody</u>	CSB-PA719094LA01HU	Human, Mouse	ELISA, WB, IHC
<u>TTC39B Antibody</u>	CSB-PA713348LA01HU	Human	ELISA, WB
<u>TTC5 Antibody</u>	CSB-PA822676LA01HU	Human, Mouse	ELISA, WB
<u>TTC9C Antibody</u>	CSB-PA836671LA01HU	Human	ELISA, WB, IHC
<u>TTF1 Antibody</u>	CSB-PA11677A0Rb	Human	ELISA
<u>TFI1 Antibody</u>	CSB-PA012188LA01HU	Human	ELISA, IHC, IF
<u>TFI2 Antibody</u>	CSB-PA764687LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>TTK Antibody</u>	CSB-PA025251LA01HU	Human	ELISA, IHC, IF
<u>TTL Antibody</u>	CSB-PA025252LA01HU	Human	ELISA
<u>TTLL2 Antibody</u>	CSB-PA863137LA01HU	Human	ELISA, IHC
<u>TTLL3 Antibody</u>	CSB-PA23329A0Rb	Human	ELISA, IHC
<u>TTLL5 Antibody</u>	CSB-PA715737LA01HU	Human	ELISA
<u>TTLL6 Antibody</u>	CSB-PA818727LA01HU	Human	ELISA, WB, IHC
<u>TTLL8 Antibody</u>	CSB-PA025265LA01HU	Human	ELISA, IHC, IF
<u>TTLL9 Antibody</u>	CSB-PA667484LA01HU	Human	ELISA, IHC
<u>TTN Antibody</u>	CSB-PA15079A0Rb	Human	ELISA, IHC, IF
<u>TTR Antibody</u>	CSB-PA00340E0Gp	Human	ELISA
<u>Ttr Antibody</u>	CSB-PA025270DA01MO	Mouse	ELISA, IHC
<u>Ttr Antibody</u>	CSB-PA025270HA01MO	Mouse	ELISA
<u>TTR Antibody</u>	CSB-PA025270HA01MOV	Macaca fascicularis	ELISA
<u>TTR Antibody</u>	CSB-PA12077A0Rb	Human	ELISA, IHC
<u>TTYH1 Antibody</u>	CSB-PA867139LA01HU	Human	ELISA, WB, IF
<u>TTYH2 Antibody</u>	CSB-PA871565LA01HU	Human	ELISA, WB
<u>TUBA1A</u>	CSB-MA754656A0m	Human, Rabbit, Rat, Mouse	ELISA, WB, IHC, IF, FC, IP
<u>TUBA1A Antibody</u>	CSB-PA754656LA01HU	Human	ELISA, IHC
<u>TUBA1B Antibody</u>	CSB-PA04189A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>TUBA1C Antibody</u>	CSB-PA025308HA01HU	Human, Mouse	ELISA, WB, IHC, IP
<u>TUBA3C Antibody</u>	CSB-PA621695LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>TUBA3E Antibody</u>	CSB-PA754255LA01HU	Human	ELISA, WB, IHC
<u>TUBA4A Antibody</u>	CSB-PA09499A0Rb	Human, Mouse	ELISA, WB, IHC, IF, IP
<u>tuba8l2 Antibody</u>	CSB-PA20829A0Rb	Zebrafish	ELISA
<u>TUBAL3 Antibody</u>	CSB-PA025315LA01HU	Human	ELISA, IHC
<u>TUBB Antibody</u>	CSB-PA025318NJ01HU	Human	ELISA, IHC
<u>TUBB Antibody</u>	CSB-PA025318OJ01HU	Human, Zebrafish	ELISA, WB
<u>TUBB1 Antibody</u>	CSB-PA867148LA01HU	Human, Mouse	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>TUBB2A Antibody</u>	CSB-PA025320LA01HU	Human	ELISA, WB
<u>TUBB3 Antibody</u>	CSB-PA03874A0Rb	Human, Mouse	ELISA, WB, IHC
<u>TUBB3 Antibody</u>	CSB-PA03879A0Rb	Human	ELISA, WB, IHC, IF
<u>TUBB4A Antibody</u>	CSB-PA025324LA01HU	Human	ELISA, WB, IHC, IF
<u>TUBB6 Antibody</u>	CSB-PA880105LA01HU	Human	ELISA, WB, IHC
<u>TUBG1 Antibody</u>	CSB-PA025334LA01HU	Human	ELISA, WB, IHC
<u>TUBG2 Antibody</u>	CSB-PA873636LA01HU	Human, Mouse	ELISA, WB, IHC
<u>TUBGCP4 Antibody</u>	CSB-PA871402LA01HU	Human	ELISA, WB, IHC, IF, IP
<u>Tubulin beta-1 chain Antibody</u>	CSB-PA18289A0Rb	Gallus gallus	ELISA
<u>tuf1 Antibody</u>	CSB-PA025342LA01AUD	Acinetobacter baumannii	ELISA
<u>TUFM Antibody</u>	CSB-PA02199A0Rb	Human	ELISA, WB
<u>TULP3 Antibody</u>	CSB-PA025348LA01HU	Human	ELISA, IHC
<u>TUSC3 Antibody</u>	CSB-PA025352LA01HU	Human, Rat	ELISA, WB, IHC
<u>TWF1 Antibody</u>	CSB-PA622511LA01HU	Human	ELISA, IHC, IF
<u>TWIST1 Antibody</u>	CSB-PA025358LA01HU	Human	ELISA, WB
<u>TWNK Antibody</u>	CSB-PA842754LA01HU	Human, Mouse	ELISA, WB, IHC
<u>TXK Antibody</u>	CSB-PA025362LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>TXN Antibody</u>	CSB-PA02814A0Rb	Human	ELISA, WB, IHC, IF
<u>TXNDC11 Antibody</u>	CSB-PA761361LA01HU	Human	ELISA, IHC, IF
<u>TXNDC12 Antibody</u>	CSB-PA025371HA01HU	Human	ELISA
<u>TXNDC8 Antibody</u>	CSB-PA718200LA01HU	Human	ELISA, IHC
<u>TXNDC9 Antibody</u>	CSB-PA025382LA01HU	Human	ELISA, IHC
<u>Txnip Antibody</u>	CSB-PA803849LA01MO	Mouse	ELISA
<u>TXNIP Antibody</u>	CSB-PA880966LA01HU	Human	ELISA, IF
<u>TXNL1 Antibody</u>	CSB-PA025384HA01HU	Human	ELISA, IHC
<u>TXNL4A Antibody</u>	CSB-PA025386LA01HU	Human	ELISA, IF
<u>TXNRD1 Antibody</u>	CSB-PA613705LA01HU	Human	ELISA, WB, IHC
<u>TYK2 Antibody</u>	CSB-PA025391LA01HU	Human	ELISA, IHC, IF
<u>TYMP Antibody</u>	CSB-PA025392LA01HU	Human	ELISA, IHC
<u>TYMS Antibody</u>	CSB-PA025393LA01HU	Human	ELISA, WB, IF, IP
<u>TYR Antibody</u>	CSB-PA025394LA01HU	Human	ELISA, IHC, IF
<u>TYRO3 Antibody</u>	CSB-PA025396LA01HU	Human	ELISA, IHC, IF
<u>TYRP1 Antibody</u>	CSB-PA025399LA01HU	Human	ELISA, IF
<u>TYW5 Antibody</u>	CSB-PA003860LA01HU	Human	ELISA, IF
<u>U2AF2 Antibody</u>	CSB-PA025408LA01HU	Human	ELISA, WB
<u>UAP1 Antibody</u>	CSB-PA618067LA01HU	Human	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>UBA1 Antibody</u>	CSB-PA00149A0Rb	Human	ELISA, WB, IF
<u>UBA5 Antibody</u>	CSB-PA875638LA01HU	Human	ELISA, WB
<u>UBA52 Antibody</u>	CSB-PA025416LA01HU	Human	ELISA, IHC, IF
<u>UBA6 Antibody</u>	CSB-PA025420LA01HU	Human	ELISA, IHC, IF
<u>UBA7 Antibody</u>	CSB-PA025421LA01HU	Human	ELISA
<u>UBAC2 Antibody</u>	CSB-PA854093LA01HU	Human	ELISA, IHC
<u>UBAP1 Antibody</u>	CSB-PA873701LA01HU	Human	ELISA, IHC
<u>UBAP2L Antibody</u>	CSB-PA622754LA01HU	Human	ELISA, IHC, IF
<u>UBASH3A Antibody</u>	CSB-PA025427LA01HU	Human	ELISA, WB, IF
<u>UBB Antibody</u>	CSB-PA04837A0Rb	Human	ELISA, IHC
<u>UBC Antibody</u>	CSB-PA03742A0Rb	Human	ELISA, WB
<u>UBD Antibody</u>	CSB-PA025435LA01HU	Human	ELISA, IHC, IF
<u>UBE2C Antibody</u>	CSB-PA025440LA01HU	Human	ELISA, WB, IHC
<u>UBE2D1 Antibody</u>	CSB-PA025443DA01HU	Human, Mouse	ELISA, WB, IHC
<u>UBE2D2 Antibody</u>	CSB-PA025444LA01HU	Human	ELISA, IHC, IF
<u>UBE2D3 Antibody</u>	CSB-PA025445HA01HU	Human, Rat	ELISA, WB
<u>UBE2DNL Antibody</u>	CSB-PA810279LA01HU	Human	ELISA
<u>UBE2E3 Antibody</u>	CSB-PA839280LA01HU	Human	ELISA, IF
<u>UBE2G1 Antibody</u>	CSB-PA025453LA01HU	Human	ELISA, IHC, IF
<u>UBE2G2 Antibody</u>	CSB-PA025454LA01HU	Human	ELISA, WB, IHC
<u>UBE2I Antibody</u>	CSB-PA00474A0Rb	Human	ELISA, WB, IHC
<u>UBE2K Antibody</u>	CSB-PA09667A0Rb	Human	ELISA, IHC, IF
<u>UBE2K Antibody</u>	CSB-PA09669A0Rb	Human	ELISA
<u>UBE2L3 Antibody</u>	CSB-PA025463LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>UBE2L6 Antibody</u>	CSB-PA025466LA01HU	Human	ELISA, IF
<u>UBE2M Antibody</u>	CSB-PA025468OA01HU	Human	ELISA, IHC
<u>UBE2N Antibody</u>	CSB-PA025470LA01HU	Human	ELISA, IHC, IF
<u>UBE2O Antibody</u>	CSB-PA866339LA01HU	Human	ELISA, WB
<u>UBE2Q2 Antibody</u>	CSB-PA845153HA01HU	Human	ELISA, WB, IHC
<u>UBE2T Antibody</u>	CSB-PA878841LA01HU	Human	ELISA, WB, IHC
<u>UBE2V1 Antibody</u>	CSB-PA025482LA01HU	Human	ELISA, IHC
<u>UBE2V2 Antibody</u>	CSB-PA025484HA01HU	Human	ELISA, IHC
<u>UBE2V2 Antibody</u>	CSB-PA025484DA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>UBE2W Antibody</u>	CSB-PA842638LA01HU	Human	ELISA
<u>UBE3A Antibody</u>	CSB-PA025488LA01HU	Human	ELISA, IHC, IF
<u>UBE3C Antibody</u>	CSB-PA621871LA01HU	Human	ELISA
<u>UBIAD1 Antibody</u>	CSB-PA23029A0Rb	Human	ELISA, IF
<u>UBL3 Antibody</u>	CSB-PA025498LA01HU	Human	ELISA, IHC
<u>UBL4A Antibody</u>	CSB-PA025499LA01HU	Human	ELISA
<u>UBL5 Antibody</u>	CSB-PA880150LA01HU	Human	ELISA



Product Name	Code	Species Reactivity	Applications
<u>UBL7 Antibody</u>	CSB-PA850429LA01HU	Human	ELISA, IHC
<u>UBLCP1 Antibody</u>	CSB-PA851974LA01HU	Human	ELISA, WB, IHC
<u>UBQLN1 Antibody</u>	CSB-PA892477LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>UBQLN2 Antibody</u>	CSB-PA883393LA01HU	Human	ELISA, IHC, IF
<u>UBQLNL Antibody</u>	CSB-PA811630LA01HU	Human	ELISA, IHC
<u>UBR1 Antibody</u>	CSB-PA812883LA01HU	Human	ELISA, IHC, IF
<u>UBR2 Antibody</u>	CSB-PA809001LA01HU	Human, Mouse	ELISA, WB
<u>UBR3 Antibody</u>	CSB-PA757959LA01HU	Human	ELISA, IHC, IF
<u>UBTD1 Antibody</u>	CSB-PA884502LA01HU	Human	ELISA
<u>UBTD2 Antibody</u>	CSB-PA845144LA01HU	Human	ELISA
<u>UBXN10 Antibody</u>	CSB-PA836251LA01HU	Human	ELISA, IHC
<u>UBXN2A Antibody</u>	CSB-PA025534LA01HU	Human	ELISA, IF
<u>UBXN2B Antibody</u>	CSB-PA614827LA01HU	Human	ELISA, IHC, IF
<u>UBXN6 Antibody</u>	CSB-PA874859LA01HU	Human, Mouse	ELISA, WB, IHC
<u>UBXN8 Antibody</u>	CSB-PA025539LA01HU	Human	ELISA, IHC
<u>UCHL1 Antibody</u>	CSB-PA11637A0Rb	Human	ELISA
<u>UCK1 Antibody</u>	CSB-PA025546LA01HU	Human	ELISA, IHC
<u>UCN2 Antibody</u>	CSB-PA822298LA01HU	Human	ELISA, IHC
<u>UCP1 Antibody</u>	CSB-PA025554LA01HU	Human	ELISA, IHC
<u>UCP2 Antibody</u>	CSB-PA025555EA01HU	Human	ELISA, IHC
<u>UCP3 Antibody</u>	CSB-PA025556LA01HU	Human	ELISA, WB, IHC, IF
<u>UFC1 Antibody</u>	CSB-PA897483LA01HU	Human	ELISA
<u>UFM1 Antibody</u>	CSB-PA025560LA01HU	Human	ELISA, IHC
<u>UFSP2 Antibody</u>	CSB-PA885731LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>UGDH Antibody</u>	CSB-PA025567HA01HU	Human	ELISA
<u>ugl Antibody</u>	CSB-PA865671HA01BGAD	Bacillus sp.	ELISA, WB
<u>UGP2 Antibody</u>	CSB-PA620986HA01HU	Human, Mouse	ELISA, WB, IHC
<u>UGT1A1 Antibody</u>	CSB-PA025570LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>UGT1A6 Antibody</u>	CSB-PA025579LA01HU	Human	ELISA, IHC
<u>UGT1A9 Antibody</u>	CSB-PA025582LA01HU	Human	ELISA, IHC, IF
<u>UGT2A3 Antibody</u>	CSB-PA761491LA01HU	Human	ELISA, IHC
<u>UGT2B4 Antibody</u>	CSB-PA025597LA01HU	Human	ELISA, IHC
<u>UGT2B7 Antibody</u>	CSB-PA025598LA01HU	Human	ELISA, IHC, IF
<u>UGT8 Antibody</u>	CSB-PA618101LA01HU	Human	ELISA
<u>UHMK1 Antibody</u>	CSB-PA840554LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>UHRF1BP1L Antibody</u>	CSB-PA025606LA01HU	Human	ELISA
<u>UHRF2 Antibody</u>	CSB-PA839395LA01HU	Human	ELISA, IF



Product Name	Code	Species Reactivity	Applications
<u>UIMC1 Antibody</u>	CSB-PA857019LA01HU	Human, Rat	ELISA, WB, IHC
<u>UL48 Antibody</u>	CSB-PA356268HA01HWZ	Human herpesvirus 1	ELISA
<u>ULBP1 Antibody</u>	CSB-PA887177LA01HU	Human	ELISA, IHC, IF
<u>ULK1 Antibody</u>	CSB-PA025612LA01HU	Human	ELISA, IF
<u>ULK2 Antibody</u>	CSB-PA811629LA01HU	Human	ELISA, WB, IHC
<u>UMODL1 Antibody</u>	CSB-PA025617LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>UMPS Antibody</u>	CSB-PA728898LA01HU	Human	ELISA, IHC, IF
<u>UNC13B Antibody</u>	CSB-PA025622LA01HU	Human	ELISA, WB, IHC, IF
<u>UNC13C Antibody</u>	CSB-PA843314LA01HU	Human	ELISA, IHC, IF
<u>UNC13D Antibody</u>	CSB-PA758222LA01HU	Human	ELISA, IHC, IF
<u>UNC45B Antibody</u>	CSB-PA812884LA01HU	Human	ELISA, IHC
<u>UNC50 Antibody</u>	CSB-PA687496LA01HU	Human	ELISA, WB, IF
<u>UNC5B Antibody</u>	CSB-PA025629LA01HU	Human	ELISA, IHC, IF
<u>UNC5D Antibody</u>	CSB-PA740925LA01HU	Human	ELISA, IF
<u>UNC79 Antibody</u>	CSB-PA603875LA01HU	Human	ELISA, IHC, IF
<u>UNC80 Antibody</u>	CSB-PA854025LA01HU	Human	ELISA, WB, IF
<u>UNC93A Antibody</u>	CSB-PA773046LA01HU	Human	ELISA, WB, IHC
<u>Uncharacterized protein Antibody</u>	CSB-PA27779A0Rb	Human	ELISA, WB
<u>UNG Antibody</u>	CSB-PA025641LA01HU	Human	ELISA, IHC, IF
<u>UPF1 Antibody</u>	CSB-PA838819LA01HU	Human	ELISA, IF
<u>UPF3B Antibody</u>	CSB-PA883646LA01HU	Human	ELISA, IHC, IF
<u>UPP2 Antibody</u>	CSB-PA025660LA01HU	Human, Mouse	ELISA, WB
<u>UPRT Antibody</u>	CSB-PA839297LA01HU	Human	ELISA, IHC
<u>UQCR11 Antibody</u>	CSB-PA03785A0Rb	Human	ELISA, IF
<u>UQCFS1 Antibody</u>	CSB-PA04645A0Rb	Human	ELISA, WB, IHC
<u>UQCRO Antibody</u>	CSB-PA03855A0Rb	Human	ELISA, WB
<u>ureA Antibody</u>	CSB-PA14347A0Rb	Helicobacter pylori	ELISA, WB
<u>ureB Antibody</u>	CSB-PA14359A0Rb	Helicobacter pylori	ELISA, WB
<u>URI1 Antibody</u>	CSB-PA529740LA01HU	Human	ELISA, WB, IHC, IF
<u>UROS Antibody</u>	CSB-PA025679LA01HU	Human	ELISA
<u>USB1 Antibody</u>	CSB-PA22889A0Rb	Human	ELISA, IHC
<u>USE1 Antibody</u>	CSB-PA889157LA01HU	Human, Mouse	ELISA, WB, IF, IP
<u>USF1 Antibody</u>	CSB-PA025681HA01HU	Human	ELISA, WB, IHC
<u>USF2 Antibody</u>	CSB-PA621888LA01HU	Human, Rat	ELISA, WB, IF
<u>USH1C Antibody</u>	CSB-PA896761HA01HU	Human	ELISA
<u>USH1G Antibody</u>	CSB-PA676572LA01HU	Human	ELISA
<u>USHBP1 Antibody</u>	CSB-PA839815LA01HU	Human	ELISA, IHC, IF
<u>USMG5 Antibody</u>	CSB-PA836231LA01HU	Human	ELISA, IF



Product Name	Code	Species Reactivity	Applications
<u>USO1 Antibody</u>	CSB-PA025694LA01HU	Human	ELISA, IHC, IF
<u>USP10 Antibody</u>	CSB-PA613516LA01HU	Human, Rat	ELISA, WB, IHC, IF, IP
<u>USP12 Antibody</u>	CSB-PA025699LA01HU	Human, Mouse	ELISA, WB
<u>USP14 Antibody</u>	CSB-PA025704LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>USP17L1 Antibody</u>	CSB-PA755915LA01HU	Human	ELISA, WB
<u>USP17L3 Antibody</u>	CSB-PA406793LA01HU	Human	ELISA, WB
<u>USP17L5 Antibody</u>	CSB-PA427555LA01HU	Human	ELISA, WB, IHC
<u>USP17L6P Antibody</u>	CSB-PA757690LA01HU	Human	ELISA, WB
<u>USP18 Antibody</u>	CSB-PA891559LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>USP2 Antibody</u>	CSB-PA025710LA01HU	Human, Mouse	ELISA, WB, IF
<u>USP20 Antibody</u>	CSB-PA896876LA01HU	Human	ELISA, IHC
<u>USP22 Antibody</u>	CSB-PA891982LA01HU	Human	ELISA, IHC, IF
<u>USP24 Antibody</u>	CSB-PA890777LA01HU	Human	ELISA, IHC
<u>USP25 Antibody</u>	CSB-PA883400LA01HU	Human	ELISA, IHC, IF
<u>USP3 Antibody</u>	CSB-PA896756LA01HU	Human	ELISA, IHC
<u>USP30 Antibody</u>	CSB-PA765106LA01HU	Human, Mouse	ELISA, WB, IHC
<u>USP31 Antibody</u>	CSB-PA747834LA01HU	Human	ELISA, IHC
<u>USP33 Antibody</u>	CSB-PA840583LA01HU	Human	ELISA, WB, IHC, IF
<u>USP36 Antibody</u>	CSB-PA878959LA01HU	Human	ELISA, IHC
<u>USP37 Antibody</u>	CSB-PA801223LA01HU	Human	ELISA, IHC
<u>USP4 Antibody</u>	CSB-PA614393LA01HU	Human	ELISA, WB, IHC, IF
<u>USP45 Antibody</u>	CSB-PA741097LA01HU	Human	ELISA, IHC
<u>USP46 Antibody</u>	CSB-PA025738LA01HU	Human	ELISA, IF
<u>USP47 Antibody</u>	CSB-PA856964LA01HU	Human	ELISA, IHC
<u>USP48 Antibody</u>	CSB-PA774808LA01HU	Human	ELISA, WB, IHC, IF
<u>USP5 Antibody</u>	CSB-PA025742LA01HU	Human, Mouse	ELISA, WB, IHC
<u>USP6 Antibody</u>	CSB-PA025747LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>USP6NL Antibody</u>	CSB-PA842149LA01HU	Human	ELISA, IHC
<u>USP7 Antibody</u>	CSB-PA849973LA01HU	Human	ELISA, IHC, IF
<u>USP9X Antibody</u>	CSB-PA856610LA01HU	Human	ELISA, IHC, IF
<u>uspF Antibody</u>	CSB-PA339887HA01ENV	Escherichia coli	ELISA
<u>UST Antibody</u>	CSB-PA897525LA01HU	Human	ELISA, IHC
<u>UTP11 Antibody</u>	CSB-PA897542LA01HU	Human	ELISA, IHC
<u>UTP15 Antibody</u>	CSB-PA840580LA01HU	Human	ELISA, WB
<u>UTP18 Antibody</u>	CSB-PA897566LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>UTP23 Antibody</u>	CSB-PA874796LA01HU	Human, Mouse	ELISA, WB, IHC
<u>UTP6 Antibody</u>	CSB-PA868357EA01HU	Human	ELISA
<u>UTY Antibody</u>	CSB-PA025775LA01HU	Human	ELISA, IF
<u>UVRAG Antibody</u>	CSB-PA882187HA01HU	Human	ELISA, WB, IHC, IF
<u>UVSSA Antibody</u>	CSB-PA012267LA01HU	Human	ELISA, IHC
<u>uvsY Antibody</u>	CSB-PA366179LA01EDZ	Enterobacteria phage T4	ELISA
<u>UXS1 Antibody</u>	CSB-PA822754LA01HU	Human, Mouse	ELISA, WB
<u>VAC14 Antibody</u>	CSB-PA600170LA01HU	Human	ELISA, IF
<u>vacA Antibody</u>	CSB-PA14379A0Rb	Helicobacter pylori	ELISA
<u>vacA Antibody</u>	CSB-PA345927LA01B2M	Helicobacter pylori	ELISA
<u>vacA Antibody</u>	CSB-PA345927LA01HUV	Helicobacter pylori	ELISA
<u>vacA Antibody</u>	CSB-PA669657LA01HUS	Helicobacter pylori	ELISA
<u>VACWR062 Antibody</u>	CSB-PA325116LA01VAI	Vaccinia virus	ELISA
<u>VAMP1 Antibody</u>	CSB-PA025780LA01HU	Human	ELISA, IHC
<u>VAMP4 Antibody</u>	CSB-PA025783LA01HU	Human, Mouse	ELISA, WB
<u>VAMP5 Antibody</u>	CSB-PA025784LA01HU	Human	ELISA, IHC
<u>VAMP7 Antibody</u>	CSB-PA025785LA01CH	Chicken	ELISA, IHC
<u>VAMP7 Antibody</u>	CSB-PA025785LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>VARS2 Antibody</u>	CSB-PA735951LA01HU	Human	ELISA, IHC
<u>VASH1 Antibody</u>	CSB-PA745337HA01HU	Human	ELISA, IHC
<u>VASH2 Antibody</u>	CSB-PA769783EA01HU	Human	ELISA, IHC
<u>VASH2 Antibody</u>	CSB-PA769783LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>VASN Antibody</u>	CSB-PA025796OA01HU	Human	ELISA, WB, IHC, IF
<u>VAT1L Antibody</u>	CSB-PA875723LA01HU	Human	ELISA, WB, IHC
<u>VAV3 Antibody</u>	CSB-PA892140LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>VAX1 Antibody</u>	CSB-PA025806LA01HU	Human	ELISA, WB
<u>VAX2 Antibody</u>	CSB-PA025807LA01HU	Human	ELISA, WB
<u>VCAM1 Antibody</u>	CSB-PA025809LA01HU	Human	ELISA
<u>VCL Antibody</u>	CSB-PA13779A0Rb	Human	ELISA, WB
<u>VCP Antibody</u>	CSB-PA025813LA01HU	Human	ELISA, WB, IHC
<u>vcp Antibody</u>	CSB-PA772050LA01DIL	Zebrafish	ELISA, WB
<u>VDAC1 Antibody</u>	CSB-PA865954LA01DOA	Arabidopsis thaliana	ELISA, WB
<u>VDAC1 Antibody</u>	CSB-PA025821PA01HU	Human	ELISA, IHC
<u>VDAC2 Antibody</u>	CSB-PA025824LA01HU	Human	ELISA, IHC
<u>VEGFA Antibody</u>	CSB-PA025833LA01PI	Pig	ELISA
<u>VEGFA Antibody</u>	CSB-PA07529A0Rb	Human	ELISA, WB
<u>Vegfa Antibody</u>	CSB-PA08249A0Rb	Rat	ELISA, WB, IHC
<u>vegfaa Antibody</u>	CSB-PA524826LA01DIL	Zebrafish	ELISA



Product Name	Code	Species Reactivity	Applications
<u>VEGFB Antibody</u>	CSB-PA025834LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>Vegfc Antibody</u>	CSB-PA025835LA01MO	Mouse	ELISA
<u>VEGFC Antibody</u>	CSB-PA07545A0Rb	Human	ELISA, WB, IHC
<u>Vegfc Antibody</u>	CSB-PA07725A0Rb	Mouse	ELISA
<u>VEGFD Antibody</u>	CSB-PA07554A0Rb	Human	ELISA, IHC
<u>Vertebrate ancient opsin Antibody</u>	CSB-PA517517LA01SWI	Salmo salar	ELISA, WB
<u>VGF Antibody</u>	CSB-PA025847LA01HU	Human	ELISA, IHC
<u>VGLL1 Antibody</u>	CSB-PA858733LA01HU	Human	ELISA, IHC
<u>VGLL2 Antibody</u>	CSB-PA836691LA01HU	Human	ELISA, IHC
<u>VGLL4 Antibody</u>	CSB-PA613487LA01HU	Human	ELISA, IHC
<u>VHLL Antibody</u>	CSB-PA761420LA01HU	Human	ELISA, IHC
<u>viaA Antibody</u>	CSB-PA16519A0Rb	Salmonella typhi	ELISA
<u>VII Antibody</u>	CSB-PA302477LA01ECY	Enterobacteria phage M13	ELISA
<u>VIM Antibody</u>	CSB-PA025857LA01HU	Human	ELISA, WB, IHC, IF
<u>VIPAS39 Antibody</u>	CSB-PA888005LA01HU	Human, Mouse	ELISA, WB, IF
<u>VIPR2 Antibody</u>	CSB-PA025860LA01HU	Human	ELISA, WB
<u>virD2 Antibody</u>	CSB-PA03944A0Rb	Agrobacterium tumefaciens	ELISA
<u>virE2 Antibody</u>	CSB-PA357288LA01AYS	Agrobacterium fabrum	ELISA
<u>vjbR Antibody</u>	CSB-PA651339YA01BMN	Brucella abortus	ELISA
<u>VKORC1 Antibody</u>	CSB-PA22899A0Rb	Human	ELISA, IHC
<u>VKORC1L1 Antibody</u>	CSB-PA818673LA01HU	Human	ELISA, IHC
<u>VLDLR Antibody</u>	CSB-PA025865LA01HU	Human	ELISA, IHC, IF
<u>VMO1 Antibody</u>	CSB-PA800218LA01HU	Human	ELISA, IHC, IF
<u>VMP1 Antibody</u>	CSB-PA856940LA01HU	Human, Mouse	ELISA, WB, IHC
<u>VN1R1 Antibody</u>	CSB-PA861973LA01HU	Human	ELISA, WB, IF
<u>VN1R2 Antibody</u>	CSB-PA847711LA01HU	Human	ELISA, WB, IF
<u>VN1R5 Antibody</u>	CSB-PA801813LA01HU	Human	ELISA, IHC, IF
<u>VNN1 Antibody</u>	CSB-PA025883LA01HU	Human	ELISA, WB, IHC, IF
<u>VOPPI Antibody</u>	CSB-PA03544A0Rb	Human	ELISA, WB
<u>VOPPI Antibody</u>	CSB-PA03545A0Rb	Human	ELISA, WB, IF
<u>VP1 Antibody</u>	CSB-PA859915LA01CID	Chicken anemia virus	ELISA, WB
<u>VP2 Antibody</u>	CSB-PA302471LA01CID	Chicken anemia virus	ELISA, WB
<u>VP3 Antibody</u>	CSB-PA860325DA01CID	Chicken anemia virus	ELISA, WB
<u>VP3 Antibody</u>	CSB-PA346512LA01CIF	Chicken anemia virus	ELISA, WB
<u>VP35 Antibody</u>	CSB-PA856262LA01RCT	Reston ebolavirus	ELISA
<u>VP40 Antibody</u>	CSB-PA810348HA01RCJ	Reston ebolavirus	ELISA
<u>VP40 Antibody</u>	CSB-PA762349HA01ZAT	Zaire ebolavirus	ELISA, WB
<u>VPREB3 Antibody</u>	CSB-PA891724LA01HU	Human	ELISA



Product Name	Code	Species Reactivity	Applications
<u>VPS16 Antibody</u>	CSB-PA863958LA01HU	Human, Rat	ELISA, WB, IHC
<u>VPS18 Antibody</u>	CSB-PA882173LA01HU	Human	ELISA, IHC
<u>VPS25 Antibody</u>	CSB-PA887108LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>VPS26B Antibody</u>	CSB-PA676824LA01HU	Human, Rat, Mouse	ELISA, WB, IHC, IF
<u>VPS29 Antibody</u>	CSB-PA890661LA01HU	Human	ELISA, WB, IHC, IF
<u>VPS33A Antibody</u>	CSB-PA850260LA01HU	Human	ELISA, IHC, IF
<u>VPS36 Antibody</u>	CSB-PA801240OA01HU	Human	ELISA, IHC
<u>VPS39 Antibody</u>	CSB-PA836235LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>VPS45 Antibody</u>	CSB-PA882092LA01HU	Human, Mouse	ELISA, WB, IHC
<u>VPS4B Antibody</u>	CSB-PA025915LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>VPS50 Antibody</u>	CSB-PA856958LA01HU	Human	ELISA, IHC
<u>VPS51 Antibody</u>	CSB-PA883424LA01HU	Human	ELISA, WB
<u>VPS52 Antibody</u>	CSB-PA850779LA01HU	Human	ELISA, IHC
<u>vpx Antibody</u>	CSB-PA320847HA01SHA	Simian immunodeficiency virus	ELISA
<u>vpx Antibody</u>	CSB-PA323941HA01SHD	Simian immunodeficiency virus	ELISA
<u>vpx Antibody</u>	CSB-PA325553LA01HLU	Human immunodeficiency virus type 2 subtype A	ELISA
<u>vpx Antibody</u>	CSB-PA356976LA01HLZ	Human immunodeficiency virus type 2 subtype A	ELISA, WB
<u>VRK2 Antibody</u>	CSB-PA768235LA01HU	Human	ELISA, IHC, IF
<u>VSIG2 Antibody</u>	CSB-PA839343EA01HU	Human	ELISA, IHC
<u>VSIG4 Antibody</u>	CSB-PA896869LA01HU	Human	ELISA, IHC, IF
<u>VSIG8 Antibody</u>	CSB-PA722789LA01HU	Human	ELISA, IHC, IF
<u>VSIR Antibody</u>	CSB-PA884484LA01HU	Human, Mouse	ELISA, WB, IHC
<u>VSTM2L Antibody</u>	CSB-PA842721LA01HU	Human	ELISA, IHC
<u>VTA1 Antibody</u>	CSB-PA882062LA01HU	Human	ELISA, IHC
<u>VTCN1 Antibody</u>	CSB-PA13005A0Rb	Human	ELISA, IHC, IF
<u>VTG2 Antibody</u>	CSB-PA16369A0Rb	Gallus gallus	ELISA
<u>VTI1A Antibody</u>	CSB-PA842635LA01HU	Human	ELISA, IHC, IF
<u>VTI1B Antibody</u>	CSB-PA887014LA01HU	Human	ELISA
<u>VTN Antibody</u>	CSB-PA025944EA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>VTN Antibody</u>	CSB-PA029544EA01HU	Human	ELISA, WB, IHC
<u>vvhA Antibody</u>	CSB-PA323085LA01VFI	Vibrio vulnificus	ELISA, WB
<u>VWA3B Antibody</u>	CSB-PA681377LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>VWA5A Antibody</u>	CSB-PA025954LA01HU	Human	ELISA, WB
<u>VWA5B2 Antibody</u>	CSB-PA818694HA01HU	Human	ELISA, IHC
<u>VWF Antibody</u>	CSB-PA12779A0Rb	Human	ELISA, IHC
<u>WAC Antibody</u>	CSB-PA861133HA01HU	Human	ELISA, IHC, IF
<u>WAPL Antibody</u>	CSB-PA772025LA01HU	Human	ELISA, WB, IHC
<u>WARS Antibody</u>	CSB-PA025965LA01HU	Human	ELISA, IHC
<u>WAS Antibody</u>	CSB-PA025967LA01HU	Human	ELISA, WB, IHC, IF
<u>WASF1 Antibody</u>	CSB-PA852883LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>WASF2 Antibody</u>	CSB-PA896553LA01HU	Human	ELISA, IHC, IF
<u>WASF3 Antibody</u>	CSB-PA891467LA01HU	Human	ELISA, WB, IHC, IF
<u>WASHC3 Antibody</u>	CSB-PA896504LA01HU	Human	ELISA
<u>WBP2 Antibody</u>	CSB-PA842628LA01HU	Human	ELISA, IHC
<u>WDFY1 Antibody</u>	CSB-PA816890LA01HU	Human, Mouse	ELISA, WB
<u>WDFY2 Antibody</u>	CSB-PA850405LA01HU	Human	ELISA, WB, IF
<u>WDFY3 Antibody</u>	CSB-PA809030LA01HU	Human	ELISA, IHC, IF
<u>WDFY4 Antibody</u>	CSB-PA747798LA01HU	Human	ELISA, IHC
<u>Wdpcp Antibody</u>	CSB-PA804911EA01MO	Mouse	ELISA, WB
<u>WDPCP Antibody</u>	CSB-PA528516LA01HU	Human	ELISA, IHC, IF
<u>WDR13 Antibody</u>	CSB-PA861995LA01HU	Human	ELISA, IHC
<u>WDR18 Antibody</u>	CSB-PA880110LA01HU	Human	ELISA
<u>WDR19 Antibody</u>	CSB-PA836733LA01HU	Human	ELISA, IHC, IF
<u>WDR24 Antibody</u>	CSB-PA836297LA01HU	Human	ELISA, IHC
<u>WDR27 Antibody</u>	CSB-PA026011LA01HU	Human	ELISA, IHC
<u>WDR34 Antibody</u>	CSB-PA842668LA01HU	Human	ELISA, WB, IHC
<u>WDR35 Antibody</u>	CSB-PA879071LA01HU	Human	ELISA, IHC
<u>WDR37 Antibody</u>	CSB-PA897530LA01HU	Human	ELISA, WB, IHC
<u>Wdr37 Antibody</u>	CSB-PA805557LA01MO	Mouse, Human	ELISA, WB, IHC
<u>WDR41 Antibody</u>	CSB-PA864011LA01HU	Human	ELISA, IF
<u>WDR49 Antibody</u>	CSB-PA810271LA01HU	Human	ELISA, WB, IHC
<u>WDR5 Antibody</u>	CSB-PA026038HA01HU	Human	ELISA
<u>WDR53 Antibody</u>	CSB-PA772027LA01HU	Human	ELISA, IHC
<u>WDR54 Antibody</u>	CSB-PA881001LA01HU	Human, Mouse	ELISA, WB, IHC
<u>WDR59 Antibody</u>	CSB-PA764810LA01HU	Human	ELISA, WB, IF
<u>WDR60 Antibody</u>	CSB-PA845154HA01HU	Human	ELISA, IHC, IF
<u>WDR61 Antibody</u>	CSB-PA872416LA01HU	Human, Mouse	ELISA, WB, IHC
<u>WDR62 Antibody</u>	CSB-PA026050LA01HU	Human	ELISA, IHC
<u>WDR66 Antibody</u>	CSB-PA848398LA01HU	Human	ELISA, IHC
<u>WDR73 Antibody</u>	CSB-PA744047LA01HU	Human	ELISA, WB, IHC, IF
<u>WDR74 Antibody</u>	CSB-PA747557LA01HU	Human	ELISA, WB, IHC



Product Name	Code	Species Reactivity	Applications
<u>WDR77 Antibody</u>	CSB-PA859030LA01HU	Human	ELISA, IHC
<u>WDR83 Antibody</u>	CSB-PA863110LA01HU	Human	ELISA, IHC, IF
<u>WDR88 Antibody</u>	CSB-PA744395LA01HU	Human	ELISA, WB, IHC
<u>WDR89 Antibody</u>	CSB-PA842671LA01HU	Human	ELISA, IHC
<u>WDR91 Antibody</u>	CSB-PA026077LA01HU	Human	ELISA, IHC
<u>WDR92 Antibody</u>	CSB-PA822265LA01HU	Human	ELISA, IHC
<u>WDTC1 Antibody</u>	CSB-PA822706LA01HU	Human, Mouse	ELISA, WB
<u>WDYHV1 Antibody</u>	CSB-PA850306LA01HU	Human, Mouse	ELISA, WB
<u>WEE1 Antibody</u>	CSB-PA16059A0Rb	Human	ELISA
<u>WFDC11 Antibody</u>	CSB-PA818778LA01HU	Human	ELISA, IHC
<u>WFDC2 Antibody</u>	CSB-PA09794A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>WFDC8 Antibody</u>	CSB-PA818226LA01HU	Human	ELISA, IHC, IF
<u>WFDC9 Antibody</u>	CSB-PA822776LA01HU	Human	ELISA
<u>WFIKKN1 Antibody</u>	CSB-PA842729OA01HU	Human	ELISA, IHC
<u>Wheat Gliadin Antibody</u>	CSB-PA00450E0Rb	Plant	ELISA, WB
<u>WHSC1L1 Antibody</u>	CSB-PA880145LA01HU	Human	ELISA, IHC
<u>WIF1 Antibody</u>	CSB-PA026113LA01HU	Human	ELISA, IHC
<u>WIPE2 Antibody</u>	CSB-PA823478LA01HU	Human	ELISA, IHC
<u>WISP1 Antibody</u>	CSB-PA026119LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>WISP2 Antibody</u>	CSB-PA026120LA01HU	Human	ELISA, IHC
<u>WISP3 Antibody</u>	CSB-PA026121LA01HU	Human, Mouse, Rat	ELISA, WB, IHC
<u>WIZ Antibody</u>	CSB-PA026123LA01HU	Human	ELISA, IHC
<u>WLS Antibody</u>	CSB-PA716359LA01HU	Human	ELISA, WB, IHC
<u>WNK1 Antibody</u>	CSB-PA887978LA01HU	Human	ELISA, IHC, IF
<u>WNK2 Antibody</u>	CSB-PA896509LA01HU	Human	ELISA, IHC, IF
<u>WNK3 Antibody</u>	CSB-PA861176LA01HU	Human	ELISA, IF
<u>WNK4 Antibody</u>	CSB-PA842690LA01HU	Human	ELISA, IHC, IF
<u>WNT1 Antibody</u>	CSB-PA026128LA01HU	Human	ELISA, IHC
<u>WNT10A Antibody</u>	CSB-PA884424LA01HU	Human	ELISA, IHC
<u>Wnt3 Antibody</u>	CSB-PA026135LA01MO	Mouse	ELISA, WB
<u>WNT3A Antibody</u>	CSB-PA15549A0Rb	Human, Mouse	ELISA, WB, IHC
<u>WNT4 Antibody</u>	CSB-PA026137LA01HU	Human	ELISA, WB, IHC
<u>WNT5A Antibody</u>	CSB-PA026138LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>WNT7A Antibody</u>	CSB-PA026141LA01HU	Human	ELISA, IHC, IF
<u>wnt8 Antibody</u>	CSB-PA340936LA01XBE	Xenopus laevis	ELISA
<u>WNT8A Antibody</u>	CSB-PA863951LA01HU	Human	ELISA, IHC
<u>wnt8b Antibody</u>	CSB-PA026144LA01DIL	Zebrafish	ELISA, WB
<u>WNT9B Antibody</u>	CSB-PA026146LA01HU	Human	ELISA, IHC
<u>WRB Antibody</u>	CSB-PA026148LA01HU	Human, Mouse	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>WRKY40 Antibody</u>	CSB-PA874370LA01DOA	Arabidopsis thaliana	ELISA
<u>WSB1 Antibody</u>	CSB-PA897589LA01HU	Human, Rat	ELISA, WB, IHC
<u>WSS1 Antibody</u>	CSB-PA334485LA01SVG	Saccharomyces cerevisiae	ELISA
<u>WT1 Antibody</u>	CSB-PA11609A0Rb	Human	ELISA, IHC
<u>WTAP Antibody</u>	CSB-PA618970LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>WWC2 Antibody</u>	CSB-PA734879LA01HU	Human	ELISA, IHC, IF
<u>WWC3 Antibody</u>	CSB-PA890929NA01HU	Human	ELISA, WB, IHC, IF
<u>WVOX Antibody</u>	CSB-PA873704LA01HU	Human, Rat	ELISA, WB, IHC, IF, IP
<u>WWP2 Antibody</u>	CSB-PA026168LA01HU	Human	ELISA, IHC
<u>X Antibody</u>	CSB-PA323774LA01WAW	Woodchuck hepatitis B virus	ELISA
<u>X Antibody</u>	CSB-PA526046LA01HVP	Hepatitis B virus genotype D	ELISA
<u>X Antibody</u>	CSB-PA865559LA01HEO	Hepatitis B virus genotype D subtype ayw	ELISA
<u>X Antibody</u>	CSB-PA355969LA01HEM	Hepatitis B virus genotype D subtype ayw	ELISA
<u>XAB2 Antibody</u>	CSB-PA862073LA01HU	Human	ELISA, WB, IHC, IF
<u>XAF1 Antibody</u>	CSB-PA764111LA01HU	Human	ELISA, IHC, IF
<u>XAGE1A Antibody</u>	CSB-PA872552LA01HU	Human	ELISA, IHC
<u>XAGE2 Antibody</u>	CSB-PA839331LA01HU	Human	ELISA
<u>XPB1 Antibody</u>	CSB-PA05319A0Rb	Human	ELISA, IHC, IF
<u>XG Antibody</u>	CSB-PA026190LA01HU	Human, Mouse	ELISA, WB
<u>XIAP Antibody</u>	CSB-PA026193LA01HU	Human	ELISA, IHC
<u>XIRP2 Antibody</u>	CSB-PA026196LA01HU	Human	ELISA, IHC
<u>XKR8 Antibody</u>	CSB-PA875681LA01HU	Human	ELISA, WB, IHC, IF
<u>XPC Antibody</u>	CSB-PA026217LA01HU	Human	ELISA
<u>XPNPEP1 Antibody</u>	CSB-PA878859HA01HU	Human, Mouse	ELISA, WB, IHC
<u>XPO1 Antibody</u>	CSB-PA026221LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>XPR1 Antibody</u>	CSB-PA871381LA01HU	Human	ELISA, IHC, IF
<u>XRCC3 Antibody</u>	CSB-PA026231LA01HU	Human	ELISA
<u>XRCC5 Antibody</u>	CSB-PA026233LA01HU	Human	ELISA, WB, IHC
<u>XRCC5 Antibody</u>	CSB-PA026233HA01HU	Human	ELISA, WB, IHC, IP
<u>XRCC6 Antibody</u>	CSB-PA01617A0Rb	Human	ELISA, WB, IHC,



Product Name	Code	Species Reactivity	Applications
			IF
<u>XRRA1 Antibody</u>	CSB-PA747397LA01HU	Human, Mouse	ELISA, WB, IHC
<u>YAF2 Antibody</u>	CSB-PA810297LA01HU	Human	ELISA, IHC
<u>YARS Antibody</u>	CSB-PA00465A0Rb	Human	ELISA, WB, IHC, IF
<u>YBEY Antibody</u>	CSB-PA003773LA01HU	Human	ELISA, IHC
<u>ydhF Antibody</u>	CSB-PA300544HA01ENV	Escherichia coli	ELISA
<u>yebF Antibody</u>	CSB-PA329240HA01ENV	Escherichia coli	ELISA, WB
<u>YES1 Antibody</u>	CSB-PA026253LA01HU	Human	ELISA, WB, IHC, IF
<u>YIF1B Antibody</u>	CSB-PA688542LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>YIPF3 Antibody</u>	CSB-PA863928LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>YIPF4 Antibody</u>	CSB-PA861129LA01HU	Human	ELISA, WB, IF
<u>YIPF5 Antibody</u>	CSB-PA842626LA01HU	Human, Mouse	ELISA, WB, IF
<u>YIPF6 Antibody</u>	CSB-PA836203LA01HU	Human	ELISA, IHC, IF
<u>YKT6 Antibody</u>	CSB-PA026265LA01HU	Human	ELISA, WB, IHC
<u>yoaJ Antibody</u>	CSB-PA522026HA01BRJ	Bacillus subtilis	ELISA
<u>yoaJ Antibody</u>	CSB-PA522026YA01BRJ	Bacillus subtilis	ELISA
<u>YOD1 Antibody</u>	CSB-PA713352LA01HU	Human	ELISA, IF
<u>yopM Antibody</u>	CSB-PA324657HA01YAS	Yersinia pestis	ELISA, WB
<u>YPEL1 Antibody</u>	CSB-PA026270LA01HU	Human	ELISA, IHC
<u>YPEL3 Antibody</u>	CSB-PA026272LA01HU	Human	ELISA, IHC
<u>YPEL5 Antibody</u>	CSB-PA026274LA01HU	Human	ELISA, IHC
<u>YTHDC1 Antibody</u>	CSB-PA853467LA01HU	Human, Mouse	ELISA, WB
<u>YTHDF2 Antibody</u>	CSB-PA896902LA01HU	Human	ELISA, WB, IHC, IF
<u>YTHDF3 Antibody</u>	CSB-PA800225LA01HU	Human	ELISA
<u>yvgO Antibody</u>	CSB-PA520104LA01BRJ	Bacillus subtilis	ELISA
<u>YWHAB Antibody</u>	CSB-PA01544A0Rb	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>YWHAE Antibody</u>	CSB-PA026287DA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>YWHAG Antibody</u>	CSB-PA026288LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IP
<u>YWHAH Antibody</u>	CSB-PA02174A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>YWHAZ Antibody</u>	CSB-PA01459A0Rb	Human, Mouse	ELISA, WB, IHC, IF
<u>YY1API Antibody</u>	CSB-PA884488LA01HU	Human	ELISA, WB, IHC, IF



Product Name	Code	Species Reactivity	Applications
<u>Zap70 Antibody</u>	CSB-PA026303LA01MO	Mouse, Human	ELISA, WB, IF
<u>zapA Antibody</u>	CSB-PA423441HA01BOL	Bacillus pumilus	ELISA, WB
<u>ZAR1 Antibody</u>	CSB-PA653624LA01DOA	Arabidopsis thaliana	ELISA, WB
<u>ZBED1 Antibody</u>	CSB-PA02715A0Rb	Human	ELISA, WB, IHC
<u>ZBED3 Antibody</u>	CSB-PA839344LA01HU	Human	ELISA, IHC
<u>ZBED6CL Antibody</u>	CSB-PA842669LA01HU	Human	ELISA, IHC
<u>ZBTB16 Antibody</u>	CSB-PA21019A0Rb	Human, Mouse	ELISA, WB, IHC
<u>ZBTB18 Antibody</u>	CSB-PA026604LA01HU	Human	ELISA, WB, IHC
<u>ZBTB32 Antibody</u>	CSB-PA897538LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ZBTB33 Antibody</u>	CSB-PA769768LA01HU	Human	ELISA, IF
<u>ZBTB43 Antibody</u>	CSB-PA026335LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ZBTB46 Antibody</u>	CSB-PA803124HA01HU	Human, Mouse	ELISA, WB, IHC
<u>ZBTB6 Antibody</u>	CSB-PA621889LA01HU	Human	ELISA, WB, IHC
<u>ZBTB8OS Antibody</u>	CSB-PA816894LA01HU	Human	ELISA, WB, IHC, IF
<u>ZC3H10 Antibody</u>	CSB-PA846649LA01HU	Human	ELISA
<u>ZC3H12A Antibody</u>	CSB-PA701732LA01HU	Human	ELISA, IHC
<u>ZC3H14 Antibody</u>	CSB-PA740793LA01HU	Human	ELISA, WB, IHC
<u>ZC3H7A Antibody</u>	CSB-PA818245LA01HU	Human	ELISA, IHC
<u>ZC4H2 Antibody</u>	CSB-PA873628LA01HU	Human, Mouse	ELISA, WB, IF
<u>ZCCHC13 Antibody</u>	CSB-PA837862LA01HU	Human	ELISA, IHC
<u>ZCCHC17 Antibody</u>	CSB-PA878834LA01HU	Human	ELISA
<u>ZCCHC4 Antibody</u>	CSB-PA887984LA01HU	Human	ELISA, WB, IHC
<u>ZDHHC1 Antibody</u>	CSB-PA837848LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ZDHHC12 Antibody</u>	CSB-PA822219LA01HU	Human	ELISA, WB, IHC, IF
<u>ZDHHC14 Antibody</u>	CSB-PA811637LA01HU	Human	ELISA, WB, IHC, IF
<u>ZDHHC15 Antibody</u>	CSB-PA842720LA01HU	Human	ELISA, IHC, IF
<u>ZDHHC2 Antibody</u>	CSB-PA883437LA01HU	Human	ELISA, IHC, IF
<u>ZDHHC20 Antibody</u>	CSB-PA722800LA01HU	Human	ELISA, IF
<u>ZDHHC23 Antibody</u>	CSB-PA818262LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ZDHHC3 Antibody</u>	CSB-PA882141LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ZDHHC4 Antibody</u>	CSB-PA873612LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ZDHHC5 Antibody</u>	CSB-PA861196LA01HU	Human	ELISA, IHC, IF
<u>ZDHHC6 Antibody</u>	CSB-PA862029LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ZEB1 Antibody</u>	CSB-PA026424LA01HU	Human	ELISA
<u>ZFAND2A Antibody</u>	CSB-PA847647LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>ZFAND6 Antibody</u>	CSB-PA739009HA01HU	Human	ELISA, IHC, IF
<u>ZFP36 Antibody</u>	CSB-PA026448LA01HU	Human	ELISA, IF
<u>ZFP36L2 Antibody</u>	CSB-PA026450LA01HU	Human, Rat	ELISA, WB
<u>ZFP42 Antibody</u>	CSB-PA856985LA01HU	Human, Mouse	ELISA, WB, IF
<u>ZFP57 Antibody</u>	CSB-PA865138LA01HU	Human, Mouse	ELISA, WB
<u>ZFP64 Antibody</u>	CSB-PA865094LA01HU	Human, Mouse	ELISA, WB
<u>ZFP69B Antibody</u>	CSB-PA026913LA01HU	Human	ELISA, IHC
<u>ZFP82 Antibody</u>	CSB-PA843135LA01HU	Human	ELISA
<u>ZFR Antibody</u>	CSB-PA842700LA01HU	Human	ELISA, IF
<u>ZFX Antibody</u>	CSB-PA026467LA01HU	Human	ELISA, IHC
<u>ZFYVE26 Antibody</u>	CSB-PA734654LA01HU	Human	ELISA, IHC, IF
<u>ZFYVE9 Antibody</u>	CSB-PA026477LA01HU	Human	ELISA, WB, IF
<u>ZG16 Antibody</u>	CSB-PA528909LA01HU	Human	ELISA, IHC
<u>ZGRF1 Antibody</u>	CSB-PA803153LA01HU	Human	ELISA, IF
<u>ZHX1-C8orf76 Antibody</u>	CSB-PA856963LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ZIC4 Antibody</u>	CSB-PA850841LA01HU	Human	ELISA, IHC, IF
<u>ZIC5 Antibody</u>	CSB-PA850433LA01HU	Human	ELISA, IHC
<u>ZIM2 Antibody</u>	CSB-PA885781LA01HU	Human	ELISA, IHC
<u>ZIP1 Antibody</u>	CSB-PA339148LA01SVG	Saccharomyces cerevisiae	ELISA, WB
<u>ZKSCAN2 Antibody</u>	CSB-PA730959LA01HU	Human	ELISA, IHC
<u>ZKSCAN3 Antibody</u>	CSB-PA861122LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ZKSCAN4 Antibody</u>	CSB-PA836168EA01HU	Human	ELISA
<u>ZKSCAN5 Antibody</u>	CSB-PA897093LA01HU	Human	ELISA, IHC
<u>ZKSCAN7 Antibody</u>	CSB-PA026552LA01HU	Human	ELISA
<u>ZKSCAN8 Antibody</u>	CSB-PA026568LA01HU	Human	ELISA
<u>ZMAT2 Antibody</u>	CSB-PA13517A0Rb	Human	ELISA
<u>ZMAT4 Antibody</u>	CSB-PA867172LA01HU	Human	ELISA, IHC
<u>ZMAT5 Antibody</u>	CSB-PA866211LA01HU	Human	ELISA, WB, IHC
<u>ZMIZ1 Antibody</u>	CSB-PA890749LA01HU	Human	ELISA, IHC, IF
<u>ZMPSTE24 Antibody</u>	CSB-PA026504LA01HU	Human	ELISA, IHC, IF
<u>ZMYM3 Antibody</u>	CSB-PA614520LA01HU	Human	ELISA, IF
<u>ZMYND11 Antibody</u>	CSB-PA621870LA01HU	Human, Mouse	ELISA, WB, IF
<u>ZMYND19 Antibody</u>	CSB-PA836202LA01HU	Human	ELISA, IHC
<u>ZNF114 Antibody</u>	CSB-PA026523LA01HU	Human	ELISA, IHC
<u>ZNF131 Antibody</u>	CSB-PA026530LA01HU	Human	ELISA, IHC
<u>ZNF133 Antibody</u>	CSB-PA026532LA01HU	Human	ELISA, WB, IHC, IF
<u>ZNF155 Antibody</u>	CSB-PA026546LA01HU	Human	ELISA, IHC
<u>ZNF165 Antibody</u>	CSB-PA026551LA01HU	Human	ELISA, IHC
<u>ZNF175 Antibody</u>	CSB-PA026556LA01HU	Human	ELISA, IHC
<u>ZNF177 Antibody</u>	CSB-PA026558LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>ZNF184 Antibody</u>	CSB-PA026563LA01HU	Human	ELISA, IHC
<u>ZNF19 Antibody</u>	CSB-PA026567LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ZNF2 Antibody</u>	CSB-PA026572LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ZNF207 Antibody</u>	CSB-PA026578LA01HU	Human	ELISA, IF
<u>ZNF212 Antibody</u>	CSB-PA026582LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ZNF215 Antibody</u>	CSB-PA026585LA01HU	Human	ELISA, WB, IHC
<u>ZNF221 Antibody</u>	CSB-PA026589LA01HU	Human	ELISA, IHC
<u>ZNF223 Antibody</u>	CSB-PA026591LA01HU	Human	ELISA
<u>ZNF224 Antibody</u>	CSB-PA026592LA01HU	Human	ELISA, IHC, IF
<u>ZNF226 Antibody</u>	CSB-PA026594LA01HU	Human	ELISA, IHC
<u>ZNF227 Antibody</u>	CSB-PA026595LA01HU	Human	ELISA
<u>ZNF230 Antibody</u>	CSB-PA026598LA01HU	Human	ELISA, IHC
<u>ZNF232 Antibody</u>	CSB-PA026599LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ZNF252P-AS1 Antibody</u>	CSB-PA612563LA01HU	Human	ELISA, IHC
<u>ZNF26 Antibody</u>	CSB-PA026619LA01HU	Human	ELISA, IHC
<u>ZNF260 Antibody</u>	CSB-PA026620LA01HU	Human	ELISA, IHC
<u>ZNF266 Antibody</u>	CSB-PA026623LA01HU	Human	ELISA, IHC
<u>ZNF268 Antibody</u>	CSB-PA026625LA01HU	Human	ELISA, WB, IF
<u>ZNF274 Antibody</u>	CSB-PA850299LA01HU	Human	ELISA, IHC
<u>ZNF281 Antibody</u>	CSB-PA026637LA01HU	Human	ELISA, IHC
<u>ZNF282 Antibody</u>	CSB-PA026638LA01HU	Human	ELISA, IHC
<u>ZNF285 Antibody</u>	CSB-PA822271LA01HU	Human	ELISA
<u>ZNF287 Antibody</u>	CSB-PA026644LA01HU	Human	ELISA, IHC
<u>ZNF3 Antibody</u>	CSB-PA026651LA01HU	Human	ELISA, IHC
<u>ZNF300 Antibody</u>	CSB-PA026653LA01HU	Human	ELISA
<u>ZNF32 Antibody</u>	CSB-PA026667LA01HU	Human	ELISA
<u>ZNF329 Antibody</u>	CSB-PA026676LA01HU	Human	ELISA
<u>ZNF331 Antibody</u>	CSB-PA026680LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ZNF333 Antibody</u>	CSB-PA026681LA01HU	Human	ELISA
<u>ZNF335 Antibody</u>	CSB-PA026683LA01HU	Human	ELISA, WB
<u>ZNF346 Antibody</u>	CSB-PA026691EA01HU	Human	ELISA, WB
<u>ZNF365 Antibody</u>	CSB-PA754628LA01HU	Human	ELISA, IHC, IF
<u>ZNF366 Antibody</u>	CSB-PA026702LA01HU	Human	ELISA, IHC, IF
<u>ZNF367 Antibody</u>	CSB-PA026703LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ZNF37A Antibody</u>	CSB-PA026704LA01HU	Human	ELISA, WB
<u>ZNF384 Antibody</u>	CSB-PA03455A0Rb	Human	ELISA, WB, IHC
<u>ZNF385C Antibody</u>	CSB-PA026712LA01HU	Human	ELISA
<u>ZNF396 Antibody</u>	CSB-PA026719LA01HU	Human	ELISA, IHC
<u>ZNF397 Antibody</u>	CSB-PA026720LA01HU	Human	ELISA, IHC
<u>ZNF398 Antibody</u>	CSB-PA026722LA01HU	Human	ELISA, WB
<u>ZNF404 Antibody</u>	CSB-PA026724LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>ZNF408 Antibody</u>	CSB-PA026728LA01HU	Human	ELISA, IHC
<u>ZNF410 Antibody</u>	CSB-PA026730LA01HU	Human, Mouse	ELISA, WB
<u>ZNF420 Antibody</u>	CSB-PA026737LA01HU	Human	ELISA, IHC
<u>ZNF423 Antibody</u>	CSB-PA026738OA01HU	Human	ELISA, IHC, IF
<u>ZNF432 Antibody</u>	CSB-PA026746LA01HU	Human	ELISA, IHC
<u>ZNF433 Antibody</u>	CSB-PA026747LA01HU	Human	ELISA
<u>ZNF440 Antibody</u>	CSB-PA026754LA01HU	Human	ELISA
<u>ZNF454 Antibody</u>	CSB-PA026765LA01HU	Human	ELISA, IHC
<u>ZNF462 Antibody</u>	CSB-PA026768LA01HU	Human	ELISA, IHC
<u>ZNF484 Antibody</u>	CSB-PA026781LA01HU	Human	ELISA, IHC
<u>ZNF486 Antibody</u>	CSB-PA026783LA01HU	Human	ELISA
<u>ZNF488 Antibody</u>	CSB-PA026785LA01HU	Human	ELISA, WB, IHC, IF
<u>ZNF490 Antibody</u>	CSB-PA026786LA01HU	Human	ELISA, IHC
<u>ZNF502 Antibody</u>	CSB-PA026795LA01HU	Human	ELISA, IHC
<u>ZNF513 Antibody</u>	CSB-PA026804LA01HU	Human	ELISA, IHC
<u>ZNF529 Antibody</u>	CSB-PA026817LA01HU	Human	ELISA, WB, IHC, IF
<u>ZNF530 Antibody</u>	CSB-PA026818LA01HU	Human	ELISA
<u>ZNF546 Antibody</u>	CSB-PA026827LA01HU	Human	ELISA, WB
<u>ZNF554 Antibody</u>	CSB-PA026834LA01HU	Human	ELISA, WB
<u>ZNF555 Antibody</u>	CSB-PA026835LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ZNF559 Antibody</u>	CSB-PA026839LA01HU	Human	ELISA, IHC
<u>ZNF560 Antibody</u>	CSB-PA026841LA01HU	Human	ELISA, IHC
<u>ZNF562 Antibody</u>	CSB-PA026843LA01HU	Human	ELISA, WB, IHC
<u>ZNF563 Antibody</u>	CSB-PA026844LA01HU	Human	ELISA, WB, IHC
<u>ZNF564 Antibody</u>	CSB-PA026845LA01HU	Human	ELISA, IHC
<u>ZNF582 Antibody</u>	CSB-PA026864LA01HU	Human	ELISA, IHC
<u>ZNF583 Antibody</u>	CSB-PA026865LA01HU	Human	ELISA, IHC
<u>ZNF586 Antibody</u>	CSB-PA026869LA01HU	Human	ELISA, IHC
<u>ZNF592 Antibody</u>	CSB-PA03634A0Rb	Human	ELISA
<u>ZNF595 Antibody</u>	CSB-PA026875LA01HU	Human	ELISA
<u>ZNF596 Antibody</u>	CSB-PA026876LA01HU	Human	ELISA, IHC, IF
<u>ZNF606 Antibody</u>	CSB-PA026886LA01HU	Human	ELISA
<u>ZNF610 Antibody</u>	CSB-PA026890LA01HU	Human	ELISA, WB
<u>ZNF614 Antibody</u>	CSB-PA026893LA01HU	Human	ELISA, WB, IHC
<u>ZNF621 Antibody</u>	CSB-PA026899LA01HU	Human	ELISA
<u>ZNF622 Antibody</u>	CSB-PA026900LA01HU	Human	ELISA, IHC
<u>ZNF624 Antibody</u>	CSB-PA026902LA01HU	Human	ELISA, IHC, IF
<u>ZNF625 Antibody</u>	CSB-PA026903LA01HU	Human, Mouse	ELISA, WB
<u>ZNF626 Antibody</u>	CSB-PA026904LA01HU	Human	ELISA, WB, IHC
<u>ZNF627 Antibody</u>	CSB-PA026905LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>ZNF644 Antibody</u>	CSB-PA026914LA01HU	Human	ELISA, IHC
<u>ZNF662 Antibody</u>	CSB-PA026928LA01HU	Human	ELISA, IHC
<u>ZNF670 Antibody</u>	CSB-PA026936LA01HU	Human	ELISA, WB, IHC
<u>ZNF679 Antibody</u>	CSB-PA026945LA01HU	Human	ELISA, IHC
<u>ZNF680 Antibody</u>	CSB-PA026946LA01HU	Human	ELISA
<u>ZNF688 Antibody</u>	CSB-PA026952LA01HU	Human	ELISA, IHC
<u>ZNF689 Antibody</u>	CSB-PA026953LA01HU	Human	ELISA, WB
<u>ZNF695 Antibody</u>	CSB-PA026957LA01HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
<u>ZNF699 Antibody</u>	CSB-PA026960LA01HU	Human	ELISA, WB, IHC
<u>ZNF703 Antibody</u>	CSB-PA026966LA01HU	Human, Rat	ELISA, WB, IHC, IF
<u>ZNF707 Antibody</u>	CSB-PA026975LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ZNF71 Antibody</u>	CSB-PA885697LA01HU	Human	ELISA, IHC
<u>ZNF738 Antibody</u>	CSB-PA843334LA01HU	Human	ELISA, IHC
<u>ZNF75A Antibody</u>	CSB-PA027013LA01HU	Human, Mouse	ELISA, WB
<u>ZNF765 Antibody</u>	CSB-PA027021LA01HU	Human	ELISA
<u>ZNF767P Antibody</u>	CSB-PA741699LA01HU	Human	ELISA, IHC
<u>ZNF77 Antibody</u>	CSB-PA027025LA01HU	Human	ELISA, IHC
<u>ZNF774 Antibody</u>	CSB-PA027030LA01HU	Human	ELISA, IHC
<u>ZNF791 Antibody</u>	CSB-PA027047LA01HU	Human	ELISA, IHC
<u>ZNF800 Antibody</u>	CSB-PA027053LA01HU	Human	ELISA, IHC
<u>ZNF816 Antibody</u>	CSB-PA610840LA01HU	Human	ELISA, IHC, IF
<u>ZNF821 Antibody</u>	CSB-PA027066LA01HU	Human	ELISA, WB, IHC
<u>ZNF829 Antibody</u>	CSB-PA027071LA01HU	Human	ELISA, WB
<u>ZNF844 Antibody</u>	CSB-PA027086LA01HU	Human	ELISA, IHC
<u>ZNF91 Antibody</u>	CSB-PA02005A0Rb	Human	ELISA, IHC
<u>ZNF98 Antibody</u>	CSB-PA027106LA01HU	Human	ELISA, IF
<u>ZNRD1 Antibody</u>	CSB-PA885791HA01HU	Human	ELISA, WB
<u>ZNRF3 Antibody</u>	CSB-PA890933HA01HU	Human	ELISA, IF
<u>zot Antibody</u>	CSB-PA334425LA01VEX	Vibrio cholerae serotype O1	ELISA, WB
<u>ZP4 Antibody</u>	CSB-PA621518LA01HU	Human	ELISA, IHC
<u>ZPBP Antibody</u>	CSB-PA866271LA01HU	Human	ELISA
<u>ZRANB1 Antibody</u>	CSB-PA883378LA01HU	Human	ELISA, WB, IHC
<u>ZRANB2 Antibody</u>	CSB-PA027127HA01HU	Human	ELISA, IHC
<u>ZRANB3 Antibody</u>	CSB-PA705058LA01HU	Human	ELISA, IHC
<u>ZRSR2 Antibody</u>	CSB-PA618014LA01HU	Human, Mouse	ELISA, WB, IHC, IF
<u>ZRSR2P1 Antibody</u>	CSB-PA614898LA01HU	Human	ELISA, IHC
<u>ZSCAN16 Antibody</u>	CSB-PA862018LA01HU	Human	ELISA
<u>ZSCAN20 Antibody</u>	CSB-PA027136LA01HU	Human	ELISA, IHC



Product Name	Code	Species Reactivity	Applications
<u>ZSCAN21 Antibody</u>	CSB-PA13789A0Rb	Human	ELISA
<u>ZSCAN29 Antibody</u>	CSB-PA809002LA01HU	Human	ELISA, IHC
<u>ZSCAN32 Antibody</u>	CSB-PA026748LA01HU	Human	ELISA, WB
<u>ZSCAN9 Antibody</u>	CSB-PA026569LA01HU	Human, Mouse	ELISA, WB, IHC
<u>ZSWIM6 Antibody</u>	CSB-PA884622LA01HU	Human	ELISA, IHC, IF
<u>ZXDC Antibody</u>	CSB-PA648494LA01HU	Human	ELISA
<u>α1-acidglycoprotein Antibody</u>	CSB-PA00330E0Rb	Human	ELISA, WB