



The catalogue of CUSABIO Histone Antibodies

Target	Product Name	Code	Tested Species	Tested Applications
CENPA	CENPA (Ab-7) Antibody	CSB-PA005205OA07nphHU	Human	ELISA, IHC
H1FO	H1FO (Ab-101) Antibody	CSB-PA010087OA101nme1HU	Human, Rat	ELISA, WB, IHC, IF, ChIP
H1FO	H1FO (Ab-107) Antibody	CSB-PA010087PA107nme1HU	Human, Mouse, Rat	ELISA, WB, IF
H1FO	H1FO (Ab-11) Antibody	CSB-PA010087PA11nme1HU	Human, Rat	ELISA, WB
H1FO	H1FO (Ab-81) Antibody	CSB-PA010087OA81nme1HU	Human	ELISA, IHC, IF, ChIP
H1FO	Mono-methyl-H1FO (K101) Antibody	CSB-PA010087OA101me1HU	Human, Rat	ELISA, WB, ICC, IF, ChIP
H1FO	Mono-methyl-H1FO (K11) Antibody	CSB-PA010087PA11me1HU	Human	ELISA, IF
H1FO	Mono-methyl-H1FO (K81) Antibody	CSB-PA010087OA81me1HU	Human	ELISA, ICC, IF, ChIP
H2AFX	H2AFX (Ab-139) Antibody	CSB-PA010097OA139nphHU	Human	ELISA, WB, IHC, IF, ChIP
H2AFX	H2AFX (Ab-142) Antibody	CSB-PA010097OA142nphHU	Human	ELISA, WB, IHC, IF
H2AFX	Phospho-H2AFX (S139) Antibody	CSB-PA010097OA139phHU	Human	ELISA, WB, IF, ChIP
H2AFX	Phospho-H2AFX (Y142) Antibody	CSB-PA010097OA142phHU	Human	ELISA, IF
H2AFZ	H2AFZ (Ab-11) Antibody	CSB-PA010100PA11nacHU	Human	ELISA, IF
H2AFZ	H2AFZ (Ab-4) Antibody	CSB-PA010100OA04nacHU	Human, Mouse	ELISA, WB, IHC, IF, IP, ChIP
H2AFZ	H2AFZ (Ab-4) Antibody	CSB-PA010100OA04nme1HU	Human	ELISA, WB, IHC, IF, ChIP
H2AFZ	H2AFZ (Ab-7) Antibody	CSB-PA010100PA07nme1HU	Human	ELISA, IF



Target	Product Name	Code	Tested Species	Tested Applications
H2AFZ	H2AFZ (Ab-7) Antibody	CSB-PA010100OA07nacHU	Human	ELISA, WB, IHC, IF
H2AFZ	Acetyl-H2AFZ (K11) Antibody	CSB-PA010100PA11acHU	Human	ELISA, WB, IF, CHIP
H2AFZ	Acetyl-H2AFZ (K4) Antibody	CSB-PA010100OA04acHU	Human	ELISA, ICC, IF, IP, CHIP
H2AFZ	Mono-methyl-H2AFZ (K4) Antibody	CSB-PA010100OA04me1HU	Human	ELISA, ICC, IF, CHIP
H2AFZ	Mono-methyl-H2AFZ (K7) Antibody	CSB-PA010100PA07me1HU	Human	ELISA, ICC, IF
H2AFZ	Acetyl-H2AFZ (K7) Antibody	CSB-PA010100OA07acHU	Human	ELISA, ICC, IF, IP, CHIP
H3F3A	H3F3A (Ab-79) Antibody	CSB-PA010109OA79nme2HU	Human, Rat	ELISA, WB, CHIP
H3F3A	Di-methyl-H3F3A (K79) Antibody	CSB-PA010109OA79me2HU	Human	ELISA, WB, CHIP
HIST1H1A	HIST1H1A (Ab-21) Antibody	CSB-PA010376PA21nme1HU	Human	ELISA, IF
HIST1H1A	HIST1H1A (Ab-21) Antibody	CSB-PA010376PA21nacHU	Human	ELISA, IF, CHIP
HIST1H1A	Acetyl-HIST1H1A (K21) Antibody	CSB-PA010376PA21acHU	Human	ELISA, IF
HIST1H1B	HIST1H1B (Ab-10) Antibody	CSB-PA010377PA10nphHU	Human	ELISA, IHC
HIST1H1B	HIST1H1B (Ab-137) Antibody	CSB-PA010377PA137nphHU	Human	ELISA, WB, IHC, IF
HIST1H1B	HIST1H1B (Ab-154) Antibody	CSB-PA010377PA154nphHU	Human	ELISA, IHC
HIST1H1B	HIST1H1B (Ab-16) Antibody	CSB-PA010377PA16nacHU	Human	ELISA, WB, IHC, IF
HIST1H1B	HIST1H1B (Ab-17) Antibody	CSB-PA010377PA17nphHU	Human, Mouse, Rat	ELISA, WB, IHC, IF
HIST1H1B	HIST1H1B (Ab-188) Antibody	CSB-PA010377PA188nphHU	Human	ELISA, WB, IHC, IF
HIST1H1B	HIST1H1B (Ab-48) Antibody	CSB-PA010377PA48nacHU	Human	ELISA, IHC



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H1B	HIST1H1B (Ab-66) Antibody	CSB-PA010377PA66nforHU	Human	ELISA, IHC
HIST1H1B	Acetyl-HIST1H1B (K16) Antibody	CSB-PA010377PA16achHU	Human	ELISA, IF, ChIP
HIST1H1B	Phospho-HIST1H1B (S17) Antibody	CSB-PA010377PA17phHU	Human	ELISA, WB, ICC, IF
HIST1H1B	Phospho-HIST1H1B (S188) Antibody	CSB-PA010377PA188phHU	Human, Rat	ELISA, WB, IF
HIST1H1B	Phospho-HIST1H1B (T154) Antibody	CSB-PA010377PA154phHU	Human	ELISA, WB, IF, ChIP
HIST1H1C	HIST1H1C (Ab-105) Antibody	CSB-PA010378OA105nme1HU	Human	ELISA, IF
HIST1H1C	HIST1H1C (Ab-109) Antibody	CSB-PA010378OA109nforHU	Human, Rat	ELISA, WB, IHC, IF
HIST1H1C	HIST1H1C (Ab-116) Antibody	CSB-PA010378OA116nhibHU	Human	ELISA, IHC
HIST1H1C	HIST1H1C (Ab-118) Antibody	CSB-PA010378PA118nme1HU	Human	ELISA, IHC, IF
HIST1H1C	HIST1H1C (Ab-158) Antibody	CSB-PA010378PA158ncrHU	Human	ELISA, IF
HIST1H1C	HIST1H1C (Ab-16) Antibody	CSB-PA010378PA16nachHU	Human, Mouse, Rat	ELISA, WB, IF, ChIP
HIST1H1C	HIST1H1C (Ab-164) Antibody	CSB-PA010378OA164nphHU	Human, Rat	ELISA, WB, IHC, ChIP
HIST1H1C	HIST1H1C (Ab-167) Antibody	CSB-PA010378OA167ncrHU	Human, Mouse, Rat	ELISA, WB, IHC, IF
HIST1H1C	HIST1H1C (Ab-168) Antibody	CSB-PA010378OA168nhibHU	Human	ELISA, IHC, IF
HIST1H1C	HIST1H1C (Ab-186) Antibody	CSB-PA010378PA186nme1HU	Human, Mouse, Rat	ELISA, WB, IP
HIST1H1C	HIST1H1C (Ab-210) Antibody	CSB-PA010378OA210nhibHU	Human	ELISA, IHC, IF
HIST1H1C	HIST1H1C (Ab-22) Antibody	CSB-PA010378OA22nhibHU	Human	ELISA, ChIP
HIST1H1C	HIST1H1C (Ab-25) Antibody	CSB-PA010378OA25nhibHU	Human	ELISA, IHC, IF



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H1C	HIST1H1C (Ab-26) Antibody	CSB-PA010378OA26nhibHU	Human	ELISA, IHC
HIST1H1C	HIST1H1C (Ab-33) Antibody	CSB-PA010378PA33ncrHU	Human, Mouse	ELISA, WB, IHC, IF
HIST1H1C	HIST1H1C (Ab-45) Antibody	CSB-PA010378OA45nme1HU	Human	ELISA, IHC, IF, ChIP
HIST1H1C	HIST1H1C (Ab-45) Antibody	CSB-PA010378OA45nme2HU	Human, Mouse, Rat	ELISA, WB, IHC, ChIP
HIST1H1C	HIST1H1C (Ab-51) Antibody	CSB-PA010378PA51nhibHU	Human	ELISA, IF
HIST1H1C	HIST1H1C (Ab-62) Antibody	CSB-PA010378PA62nacHU	Human	ELISA, IF
HIST1H1C	HIST1H1C (Ab-62) Antibody	CSB-PA010378OA62nsucHU	Human	ELISA, IHC
HIST1H1C	HIST1H1C (Ab-62) Antibody	CSB-PA010378OA62ncrHU	Human	ELISA, IHC
HIST1H1C	HIST1H1C (Ab-74) Antibody	CSB-PA010378OA74nhibHU	Human	ELISA, IHC
HIST1H1C	HIST1H1C (Ab-80) Antibody	CSB-PA010378OA80nhibHU	Human	ELISA, IHC, IF
HIST1H1C	HIST1H1C (Ab-89) Antibody	CSB-PA010378OA89nacHU	Human	ELISA, IHC
HIST1H1C	HIST1H1C (Ab-96) Antibody	CSB-PA010378PA96nacHU	Human, Mouse, Rat	ELISA, WB, IF
HIST1H1C	HIST1H1C (Ab-96) Antibody	CSB-PA010378OA96ncrHU	Human	ELISA, IHC, IF
HIST1H1C	HIST1H1C (Ab-96) Antibody	CSB-PA010378OA96nme1HU	Human, Mouse, Rat	ELISA, WB, IHC, IF, ChIP
HIST1H1C	HIST1H1C (Ab-96) Antibody	CSB-PA010378OA96nhibHU	Human	ELISA, WB, IHC, IF
HIST1H1C	HIST1H1C (Ab-96) Antibody	CSB-PA010378OA96nforHU	Human	ELISA, WB, IHC, IF
HIST1H1C	Mono-methyl-HIST1H1C (K105) Antibody	CSB-PA010378OA105me1HU	Human	ELISA, IF
HIST1H1C	Formyl-HIST1H1C (K109) Antibody	CSB-PA010378OA109forHU	Human, Rat	ELISA, WB, ICC, IF



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H1C	2-hydroxyisobutyryl-HIST1H1C (K109) Antibody	CSB-PA010378OA109hibHU	Human	ELISA, WB, ICC, IF, ChIP
HIST1H1C	2-hydroxyisobutyryl-HIST1H1C (K116) Antibody	CSB-PA010378OA116hibHU	Human	ELISA, ICC, IF, ChIP
HIST1H1C	Mono-methyl-HIST1H1C (K118) Antibody	CSB-PA010378PA118me1HU	Human	ELISA, IF, ChIP
HIST1H1C	2-hydroxyisobutyryl-HIST1H1C (K128) Antibody	CSB-PA010378OA128hibHU	Human	ELISA, ICC
HIST1H1C	2-hydroxyisobutyryl-HIST1H1C (K135) Antibody	CSB-PA010378OA135hibHU	Human	ELISA, IF
HIST1H1C	Crotonyl-HIST1H1C (K158) Antibody	CSB-PA010378PA158crHU	Human	ELISA, ICC, IF, ChIP
HIST1H1C	2-hydroxyisobutyryl-HIST1H1C (K158) Antibody	CSB-PA010378OA158hibHU	Human	ELISA, WB, ICC, ChIP
HIST1H1C	Formyl-HIST1H1C (K159) Antibody	CSB-PA010378QA159forHU	Human	ELISA, IF
HIST1H1C	Acetyl-HIST1H1C (K16) Antibody	CSB-PA010378PA16acHU	Human	ELISA, ICC, IF, ChIP
HIST1H1C	2-hydroxyisobutyryl-HIST1H1C (K168) Antibody	CSB-PA010378OA168hibHU	Human	ELISA, ICC, IF
HIST1H1C	Mono-methyl-HIST1H1C (K186) Antibody	CSB-PA010378PA186me1HU	Human, Mouse, Rat	ELISA, WB, IF, ChIP
HIST1H1C	2-hydroxyisobutyryl-HIST1H1C (K210) Antibody	CSB-PA010378OA210hibHU	Human	ELISA, ICC
HIST1H1C	2-hydroxyisobutyryl-HIST1H1C (K22) Antibody	CSB-PA010378OA22hibHU	Human	ELISA, ICC, IF
HIST1H1C	2-hydroxyisobutyryl-HIST1H1C (K25) Antibody	CSB-PA010378OA25hibHU	Human	ELISA, ICC, IF
HIST1H1C	2-hydroxyisobutyryl-HIST1H1C (K26) Antibody	CSB-PA010378OA26hibHU	Human, Rat	ELISA, WB, IF
HIST1H1C	Mono-methyl-HIST1H1C (K45) Antibody	CSB-PA010378OA45me1HU	Human	ELISA, WB, ICC, IF
HIST1H1C	Di-methyl-HIST1H1C (K45) Antibody	CSB-PA010378OA45me2HU	Human	ELISA, WB, ICC, IP, ChIP
HIST1H1C	Formyl-HIST1H1C (K45) Antibody	CSB-PA010378OA45forHU	Human	ELISA, ICC, IF



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H1C	2-hydroxyisobutyryl-HIST1H1C (K51) Antibody	CSB-PA010378PA51hibHU	Human	ELISA, IF
HIST1H1C	Acetyl-HIST1H1C (K62) Antibody	CSB-PA010378PA62achU	Human	ELISA, ICC, IF, CHIP
HIST1H1C	2-hydroxyisobutyryl-HIST1H1C (K63) Antibody	CSB-PA010378PA63hibHU	Human	ELISA, IF
HIST1H1C	Formyl-HIST1H1C (K74) Antibody	CSB-PA010378OA74forHU	Human	ELISA, CHIP
HIST1H1C	Acetyl-HIST1H1C (K74) Antibody	CSB-PA010378OA74achU	Human	ELISA, IF, CHIP
HIST1H1C	2-hydroxyisobutyryl-HIST1H1C (K80) Antibody	CSB-PA010378OA80hibHU	Human	ELISA, CHIP
HIST1H1C	Formyl-HIST1H1C (K84) Antibody	CSB-PA010378OA84forHU	Human	ELISA, IF
HIST1H1C	Acetyl-HIST1H1C (K84) Antibody	CSB-PA010378OA84achU	Human	ELISA, WB, ICC, CHIP
HIST1H1C	Acetyl-HIST1H1C (K96) Antibody	CSB-PA010378PA96achU	Human	ELISA, WB, ICC, IF, IP, CHIP
HIST1H1C	Mono-methyl-HIST1H1C (K96) Antibody	CSB-PA010378OA96me1HU	Human, Mouse, Rat	ELISA, WB, IF, CHIP
HIST1H1C	Formyl-HIST1H1C (K96) Antibody	CSB-PA010378OA96forHU	Human, Rat	ELISA, WB, ICC, IF
HIST1H1C	Phospho-HIST1H1C (T164) Antibody	CSB-PA010378OA164phHU	Human	ELISA, IF
HIST1H1D	HIST1H1D (Ab-106) Antibody	CSB-PA010379OA106nbhbHU	Human	ELISA, WB, IHC
HIST1H1D	HIST1H1D (Ab-140) Antibody	CSB-PA010379PA140nforHU	Human, Mouse, Rat	ELISA, WB, IHC, IF
HIST1H1D	HIST1H1D (Ab-146) Antibody	CSB-PA010379PA146nphHU	Human, Mouse, Rat	ELISA, WB, IHC, IF
HIST1H1D	HIST1H1D (Ab-16) Antibody	CSB-PA010379PA16nme1HU	Human	ELISA, IF
HIST1H1D	HIST1H1D (Ab-179) Antibody	CSB-PA010379OA179nphHU	Human	ELISA, IF, CHIP
HIST1H1D	HIST1H1D (Ab-34) Antibody	CSB-PA010379OA34nbhbHU	Human	ELISA, IHC, IF



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H1D	HIST1H1D (Ab-64) Antibody	CSB-PA010379OA64nbhbHU	Human	ELISA, IHC
HIST1H1D	HIST1H1D (Ab-85) Antibody	CSB-PA010379OA85nbhbHU	Human	ELISA, IHC, CHIP
HIST1H1D	HIST1H1D (Ab-90) Antibody	CSB-PA010379OA90nbhbHU	Human	ELISA, IHC, IF
HIST1H1D	β-hydroxybutyryl-HIST1H1D (K106) Antibody	CSB-PA010379OA106bhbHU	Human	ELISA, WB, ICC
HIST1H1D	Mono-methyl-HIST1H1D (K16) Antibody	CSB-PA010379PA16me1HU	Human	ELISA, IF
HIST1H1D	β-hydroxybutyryl-HIST1H1D (K75) Antibody	CSB-PA010379OA75bhbHU	Human, Rat	ELISA, WB, ICC
HIST1H1D	Phospho-HIST1H1D (T146) Antibody	CSB-PA010379PA146phHU	Human	ELISA, WB, IF
HIST1H1D	Phospho-HIST1H1D (T179) Antibody	CSB-PA010379QA179phHU	Human	ELISA, IF, ChIP
HIST1H1E	HIST1H1E (Ab-145) Antibody	CSB-PA010380OA145nphHU	Human	ELISA, ChIP
HIST1H1E	HIST1H1E (Ab-147) Antibody	CSB-PA010380PA147nme1HU	Human	ELISA, IF
HIST1H1E	HIST1H1E (Ab-16) Antibody	CSB-PA010380PA16nme2HU	Human	ELISA, IF
HIST1H1E	HIST1H1E (Ab-16) Antibody	CSB-PA010380PA16nachHU	Human	ELISA, WB, IHC, IF
HIST1H1E	HIST1H1E (Ab-17) Antibody	CSB-PA010380OA17nphHU	Human, Rat	ELISA, WB, IHC, IF
HIST1H1E	HIST1H1E (Ab-25) Antibody	CSB-PA010380PA25nme1HU	Human, Mouse	ELISA, WB, IHC, IF, ChIP
HIST1H1E	HIST1H1E (Ab-26) Antibody	CSB-PA010380PA26nphHU	Human	ELISA, WB, IHC, IF
HIST1H1E	HIST1H1E (Ab-3) Antibody	CSB-PA010380PA03nphHU	Human	ELISA, WB, IHC, IF
HIST1H1E	HIST1H1E (Ab-33) Antibody	CSB-PA010380PA33nachHU	Human	ELISA, WB, IHC, IF
HIST1H1E	HIST1H1E (Ab-35) Antibody	CSB-PA010380PA35nphHU	Human	ELISA, IF



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H1E	HIST1H1E (Ab-45) Antibody	CSB-PA010380PA45nacHU	Human	ELISA, IF
HIST1H1E	HIST1H1E (Ab-45) Antibody	CSB-PA010380OA45nacHU	Human	ELISA, WB, IHC, IF
HIST1H1E	HIST1H1E (Ab-51) Antibody	CSB-PA010380PA51nacHU	Human	ELISA, IF
HIST1H1E	HIST1H1E (Ab-62) Antibody	CSB-PA010380PA62nforHU	Human	ELISA, WB, IHC, IF
HIST1H1E	HIST1H1E (Ab-63) Antibody	CSB-PA010380PA63nacHU	Human	ELISA, IHC
HIST1H1E	HIST1H1E (Ab-63) Antibody	CSB-PA010380PA63nme1HU	Human	ELISA, IHC
HIST1H1E	HIST1H1E (Ab-63) Antibody	CSB-PA010380PA63ncrHU	Human	ELISA, IHC
HIST1H1E	HIST1H1E (Ab-63) Antibody	CSB-PA010380PA63nforHU	Human	ELISA, IHC
HIST1H1E	Mono-methyl-HIST1H1E (K147) Antibody	CSB-PA010380PA147me1HU	Human	ELISA, ICC, IF
HIST1H1E	Di-methyl-HIST1H1E (K16) Antibody	CSB-PA010380PA16me2HU	Human	ELISA, ICC, IF, CHIP
HIST1H1E	Acetyl-HIST1H1E (K16) Antibody	CSB-PA010380PA16acHU	Human	ELISA, ICC, IF, CHIP
HIST1H1E	Mono-methyl-HIST1H1E (K25) Antibody	CSB-PA010380PA25me1HU	Human	ELISA, ICC, IF
HIST1H1E	Acetyl-HIST1H1E (K33) Antibody	CSB-PA010380PA33acHU	Human	ELISA, ICC, IF, CHIP
HIST1H1E	Acetyl-HIST1H1E (K51) Antibody	CSB-PA010380PA51acHU	Human	ELISA, ICC, IF, CHIP
HIST1H1E	Acetyl-HIST1H1E (K63) Antibody	CSB-PA010380PA63acHU	Human	ELISA, ICC, IF, CHIP
HIST1H1E	Phospho-HIST1H1E (S35) Antibody	CSB-PA010380PA35phHU	Human	ELISA, ICC, IF
HIST1H1E	Phospho-HIST1H1E (T145) Antibody	CSB-PA010380OA145phHU	Human	ELISA, IF
HIST1H1E	Phospho-HIST1H1E (T17) Antibody	CSB-PA010380OA17phHU	Human, Mouse	ELISA, WB, IF



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H1E	Phospho-HIST1H1E (T3) Antibody	CSB-PA010380PA03phHU	Human	ELISA, IF
HIST1H2AG	HIST1H2AG (Ab-101) Antibody	CSB-PA010389PA101nogHU	Human	ELISA, IHC, IF
HIST1H2AG	HIST1H2AG (Ab-11) Antibody	CSB-PA010389OA11nme2HU	Human	ELISA, WB, IF
HIST1H2AG	HIST1H2AG (Ab-118) Antibody	CSB-PA010389PA118ncrHU	Human	ELISA, WB, IHC, IF
HIST1H2AG	HIST1H2AG (Ab-118) Antibody	CSB-PA010389OA118nme1HU	Human, Mouse	ELISA, WB, IP, ChIP
HIST1H2AG	HIST1H2AG (Ab-118) Antibody	CSB-PA010389OA118nforHU	Human, Mouse	ELISA, WB, IHC, IF, IP
HIST1H2AG	HIST1H2AG (Ab-118) Antibody	CSB-PA010389OA118nhibHU	Human	ELISA, IHC, IF, ChIP
HIST1H2AG	HIST1H2AG (Ab-118) Antibody	CSB-PA010389OA118nbhbHU	Human	ELISA, IHC
HIST1H2AG	HIST1H2AG (Ab-119) Antibody	CSB-PA010389PA119ncrHU	Human	ELISA, IHC, IF, ChIP
HIST1H2AG	HIST1H2AG (Ab-119) Antibody	CSB-PA010389OA119nbhbHU	Human	ELISA, IHC, IF
HIST1H2AG	HIST1H2AG (Ab-125) Antibody	CSB-PA010389PA125ncrHU	Human	ELISA, IHC, IF
HIST1H2AG	HIST1H2AG (Ab-13) Antibody	CSB-PA010389OA13nacHU	Human, Mouse, Rat	ELISA, WB, IHC, ChIP
HIST1H2AG	HIST1H2AG (Ab-29) Antibody	CSB-PA010389OA29nme2HU	Human	ELISA, WB, IHC, IF
HIST1H2AG	HIST1H2AG (Ab-36) Antibody	CSB-PA010389PA36ncrHU	Human	ELISA, IHC, IF
HIST1H2AG	HIST1H2AG (Ab-36) Antibody	CSB-PA010389PA36nacHU	Human	ELISA, WB, IHC, IF, ChIP
HIST1H2AG	HIST1H2AG (Ab-36) Antibody	CSB-PA010389OA36nhibHU	Human, Mouse	ELISA, WB, IHC, IF
HIST1H2AG	HIST1H2AG (Ab-36) Antibody	CSB-PA010389OA36nsuchHU	Human	ELISA, WB, IHC
HIST1H2AG	HIST1H2AG (Ab-36) Antibody	CSB-PA010389OA36nbhbHU	Human	ELISA, WB



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H2AG	HIST1H2AG (Ab-5) Antibody	CSB-PA010389PA05nacHU	Human	ELISA, IF
HIST1H2AG	HIST1H2AG (Ab-5) Antibody	CSB-PA010389OA05nacHU	Human	ELISA, WB, IF
HIST1H2AG	HIST1H2AG (Ab-5) Antibody	CSB-PA010389OA05nhibHU	Human, Mouse, Rat	ELISA, WB, IHC, IF
HIST1H2AG	HIST1H2AG (Ab-5) Antibody	CSB-PA010389OA05nbhbHU	Human	ELISA, IHC
HIST1H2AG	Histone-H2A(Ab-5) Antibody	CSB-PA010389OA05nbutHU	Human	ELISA, WB, IHC, IF
HIST1H2AG	HIST1H2AG (Ab-74) Antibody	CSB-PA010389OA74nacHU	Human, Rat	ELISA, WB, IP
HIST1H2AG	HIST1H2AG (Ab-9) Antibody	CSB-PA010389OA09nhibHU	Human	ELISA, WB, IHC, IF
HIST1H2AG	HIST1H2AG (Ab-9) Antibody	CSB-PA010389OA09nsucHU	Human	ELISA, WB, IHC, IF
HIST1H2AG	HIST1H2AG (Ab-9) Antibody	CSB-PA010389OA09nme1HU	Human	ELISA, CHIP
HIST1H2AG	HIST1H2AG (Ab-9) Antibody	CSB-PA010389OA09nacHU	Human	ELISA, IHC, IF
HIST1H2AG	HIST1H2AG (Ab-9) Antibody	CSB-PA010389OA09nbhbHU	Human	ELISA, IHC
HIST1H2AG	HIST1H2AG (Ab-95) Antibody	CSB-PA010389PA95nforHU	Human	ELISA, IHC, IF
HIST1H2AG	HIST1H2AG (Ab-95) Antibody	CSB-PA010389PA95nhibHU	Human	ELISA, IHC
HIST1H2AG	HIST1H2AG (Ab-95) Antibody	CSB-PA010389OA95nsucHU	Human	ELISA, IHC
HIST1H2AG	HIST1H2AG (Ab-95) Antibody	CSB-PA010389OA95ncrHU	Human	ELISA, IHC, IF, CHIP
HIST1H2AG	Crotonyl-HIST1H2AG (K118) Antibody	CSB-PA010389PA118crHU	Human	ELISA, WB, ICC, IF, CHIP
HIST1H2AG	Formyl-HIST1H2AG (K118) Antibody	CSB-PA010389OA118forHU	Human, Rat	ELISA, WB, IF
HIST1H2AG	Crotonyl-HIST1H2AG (K125) Antibody	CSB-PA010389PA125crHU	Human	ELISA, IF, CHIP



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H2AG	Acetyl-HIST1H2AG (K13) Antibody	CSB-PA010389OA13acHU	Human, Mouse, Rat	ELISA, WB, IHC, IF, IP
HIST1H2AG	Acetyl-HIST1H2AG (K15) Antibody	CSB-PA010389OA15acHU	Human	ELISA, ICC, IF
HIST1H2AG	Acetyl-HIST1H2AG (K36) Antibody	CSB-PA010389PA36acHU	Human	ELISA, WB, ICC, IF, ChIP
HIST1H2AG	2-hydroxyisobutyryl-HIST1H2AG (K36) Antibody	CSB-PA010389OA36hibHU	Human	ELISA, ICC
HIST1H2AG	β-hydroxybutyryl-HIST1H2AG (K36) Antibody	CSB-PA010389OA36bhbHU	Human, Rat	ELISA, WB
HIST1H2AG	Acetyl-HIST1H2AG (K5) Antibody	CSB-PA010389OA05acHU	Human	ELISA, WB, ICC, IF
HIST1H2AG	2-hydroxyisobutyryl-HIST1H2AG (K74) Antibody	CSB-PA010389OA74hibHU	Human, Rat	ELISA, WB, ICC
HIST1H2AG	Acetyl-HIST1H2AG (K74) Antibody	CSB-PA010389OA74acHU	Human, Rat	ELISA, WB, ICC, IF
HIST1H2AG	Mono-methyl-HIST1H2AG (K9) Antibody	CSB-PA010389OA09me1HU	Human	ELISA, ICC, IF, ChIP
HIST1H2AG	2-hydroxyisobutyryl-HIST1H2AG (K9) Antibody	CSB-PA010389OA09hibHU	Human, Rat	ELISA, WB
HIST1H2AG	Acetyl-HIST1H2AG (K9) Antibody	CSB-PA010389OA09acHU	Human, Rat	ELISA, WB, IF
HIST1H2AG	Formyl-HIST1H2AG (K95) Antibody	CSB-PA010389PA95forHU	Human	ELISA, ICC, IF
HIST1H2AG	2-hydroxyisobutyryl-HIST1H2AG (K95) Antibody	CSB-PA010389PA95hibHU	Human	ELISA, IF
HIST1H2AG	Di-methyl-HIST1H2AG (R11) Antibody	CSB-PA010389OA11me2HU	Human	ELISA, IF
HIST1H2AG	Mono-methyl-HIST1H2AG (R29) Antibody	CSB-PA010389OA29me1HU	Human	ELISA, IF
HIST1H2AG	Di-methyl-HIST1H2AG (R29) Antibody	CSB-PA010389OA29me2HU	Human	ELISA, ChIP
HIST1H2BB	HIST1H2BB (Ab-16) Antibody	CSB-PA010402PA16nacHU	Human	ELISA, IF
HIST1H2BB	HIST1H2BB (Ab-16) Antibody	CSB-PA010402OA16nacHU	Human, Mouse, Rat	ELISA, WB, IF, ChIP



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H2BB	HIST1H2BB (Ab-5) Antibody	CSB-PA010402PA05nacHU	Human	ELISA, IHC, IF
HIST1H2BB	HIST1H2BB (Ab-5) Antibody	CSB-PA010402NA05nacHU	Human	ELISA, IF
HIST1H2BB	Acetyl-HIST1H2BB (K16) Antibody	CSB-PA010402PA16acHU	Human	ELISA, WB, ICC
HIST1H2BB	Acetyl-HIST1H2BB (K16) Antibody	CSB-PA010402OA16acHU	Human	ELISA, ICC, IF, ChIP
HIST1H2BB	Acetyl-HIST1H2BB (K5) Antibody	CSB-PA010402PA05acHU	Human	ELISA, ICC, IF
HIST1H2BB	Acetyl-HIST1H2BB (K5) Antibody	CSB-PA010402NA05acHU	Human	ELISA, ICC, IF, ChIP
HIST1H2BC	HIST1H2BC (Ab-108) Antibody	CSB-PA010403OA108nforHU	Human	ELISA, WB, IHC, IP
HIST1H2BC	HIST1H2BC (Ab-108) Antibody	CSB-PA010403OA108nhibHU	Human, Rat	ELISA, WB, ChIP
HIST1H2BC	HIST1H2BC (Ab-108) Antibody	CSB-PA010403OA108nacHU	Human	ELISA, IF
HIST1H2BC	HIST1H2BC (Ab-11) Antibody	CSB-PA010403OA11ncrHU	Human	ELISA, IHC, IF
HIST1H2BC	HIST1H2BC (Ab-11) Antibody	CSB-PA010403OA11nacHU	Human	ELISA, IHC, IF
HIST1H2BC	HIST1H2BC (Ab-11) Antibody	CSB-PA010403OA11nbhbHU	Human	ELISA, IHC, IF
HIST1H2BC	HIST1H2BC (Ab-116) Antibody	CSB-PA010403OA116nforHU	Human, Rat	ELISA, WB, IHC, IF
HIST1H2BC	HIST1H2BC (Ab-116) Antibody	CSB-PA010403OA116nhibHU	Human, Rat	ELISA, WB
HIST1H2BC	HIST1H2BC (Ab-116) Antibody	CSB-PA010403OA116nme1HU	Human, Mouse	ELISA, WB, IP, ChIP
HIST1H2BC	HIST1H2BC (Ab-116) Antibody	CSB-PA010403OA116nbhbHU	Human	ELISA, WB
HIST1H2BC	HIST1H2BC (Ab-116) Antibody	CSB-PA010403OA116nacHU	Human	ELISA, WB, IHC, IF
HIST1H2BC	HIST1H2BC (Ab-116) Antibody	CSB-PA010403OA116ngluHU	Human, Rat	ELISA, WB, IHC



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H2BC	HIST1H2BC (Ab-12) Antibody	CSB-PA010403OA12ncrHU	Human	ELISA, WB, IF
HIST1H2BC	HIST1H2BC (Ab-12) Antibody	CSB-PA010403OA12nme1HU	Human, Mouse, Rat	ELISA, WB, CHIP
HIST1H2BC	HIST1H2BC (Ab-12) Antibody	CSB-PA010403OA12nacHU	Human, Mouse, Rat	ELISA, WB, IHC, IF
HIST1H2BC	HIST1H2BC (Ab-12) Antibody	CSB-PA010403OA12nhibHU	Human	ELISA, WB, IHC
HIST1H2BC	Histone H2B (Ab-12) Antibody	CSB-PA010403OA12nbutHU	Human	ELISA, IHC, IF
HIST1H2BC	HIST1H2BC (Ab-120) Antibody	CSB-PA010403OA120nforHU	Human, Rat	ELISA, WB, CHIP
HIST1H2BC	HIST1H2BC (Ab-120) Antibody	CSB-PA010403OA120nhibHU	Human, Rat	ELISA, WB, IHC
HIST1H2BC	HIST1H2BC (Ab-120) Antibody	CSB-PA010403OA120ngluHU	Human	ELISA, IHC
HIST1H2BC	HIST1H2BC (Ab-120) Antibody	CSB-PA010403OA120nbhbHU	Human, Rat	ELISA, WB, IHC
HIST1H2BC	Histone H2B (Ab-120) Antibody	CSB-PA010403OA120nacHU	Human	ELISA, IHC
HIST1H2BC	HIST1H2BC (Ab-14) Antibody	CSB-PA010403OA14nphHU	Human, Mouse, Rat	ELISA, WB, CHIP
HIST1H2BC	HIST1H2BC (Ab-15) Antibody	CSB-PA010403OA15nme1HU	Human	ELISA, IHC
HIST1H2BC	HIST1H2BC (Ab-15) Antibody	CSB-PA010403OA15nacHU	Human, Rat	ELISA, WB, IHC
HIST1H2BC	HIST1H2BC (Ab-16) Antibody	CSB-PA010403OA16ncrHU	Human, Rat	ELISA, WB, IHC, IF
HIST1H2BC	HIST1H2BC (Ab-20) Antibody	CSB-PA010403OA20nme1HU	Human, Mouse	ELISA, WB
HIST1H2BC	HIST1H2BC (Ab-20) Antibody	CSB-PA010403OA20ncrHU	Human	ELISA, IHC, IF
HIST1H2BC	HIST1H2BC (Ab-20) Antibody	CSB-PA010403OA20nbutHU	Human, Rat	ELISA, WB, IHC, IF
HIST1H2BC	HIST1H2BC (Ab-20) Antibody	CSB-PA010403OA20nhibHU	Human, Mouse	ELISA, WB, IHC, IF



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H2BC	HIST1H2BC (Ab-20) Antibody	CSB-PA010403OA20nbhbHU	Human	ELISA, IHC, IF
HIST1H2BC	HIST1H2BC (Ab-23) Antibody	CSB-PA010403OA23ncrHU	Human	ELISA, WB, IHC, IF
HIST1H2BC	HIST1H2BC (Ab-23) Antibody	CSB-PA010403OA23nhibHU	Human	ELISA, IHC
HIST1H2BC	HIST1H2BC (Ab-34) Antibody	CSB-PA010403OA34ncrHU	Human	ELISA, IHC
HIST1H2BC	HIST1H2BC (Ab-34) Antibody	CSB-PA010403OA34nsuchU	Human	ELISA, IHC, IF
HIST1H2BC	HIST1H2BC (Ab-34) Antibody	CSB-PA010403OA34nbhbHU	Human	ELISA, IF
HIST1H2BC	HIST1H2BC (Ab-34) Antibody	CSB-PA010403OA34nhibHU	Human	ELISA, IF
HIST1H2BC	HIST1H2BC (Ab-36) Antibody	CSB-PA010403QA36nogHU	Human	ELISA, IHC
HIST1H2BC	HIST1H2BC (Ab-5) Antibody	CSB-PA010403OA05ncrHU	Human	ELISA, IHC, IF
HIST1H2BC	HIST1H2BC (Ab-5) Antibody	CSB-PA010403OA05nforHU	Human, Rat	ELISA, WB, IHC, IF
HIST1H2BC	HIST1H2BC (Ab-5) Antibody	CSB-PA010403OA05nhibHU	Human	ELISA, IHC
HIST1H2BC	HIST1H2BC (Ab-5) Antibody	CSB-PA010403OA05nbhbHU	Human	ELISA, IHC, IF
HIST1H2BC	HIST1H2BC (Ab-5) Antibody	CSB-PA010403OA05ngluHU	Human	ELISA, IHC, IF
HIST1H2BC	HIST1H2BC (Ab-5) Antibody	CSB-PA010403OA05nbutHU	Human	ELISA, WB, IHC
HIST1H2BC	HIST1H2BC (Ab-79) Antibody	CSB-PA010403OA79nme1HU	Human	ELISA, IHC, IF
HIST1H2BC	HIST1H2BC (Ab-85) Antibody	CSB-PA010403OA85nme1HU	Human	ELISA, WB
HIST1H2BC	HIST1H2BC (Ab-85) Antibody	CSB-PA010403OA85nsuchU	Human	ELISA, IHC
HIST1H2BC	HIST1H2BC (Ab-85) Antibody	CSB-PA010403OA85nbhbHU	Human	ELISA, WB, IF



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H2BC	HIST1H2BC (Ab-85) Antibody	CSB-PA010403OA85nacHU	Human	ELISA, IHC
HIST1H2BC	HIST1H2BC (Ab-99) Antibody	CSB-PA010403OA99nme1HU	Human	ELISA, IHC
HIST1H2BC	2-hydroxyisobutyryl-HIST1H2BC (K108) Antibody	CSB-PA010403OA108hibHU	Human, Rat	ELISA, WB, IF
HIST1H2BC	Acetyl-HIST1H2BC (K108) Antibody	CSB-PA010403OA108acHU	Human	ELISA, ICC, IF, ChIP
HIST1H2BC	Crotonyl-HIST1H2BC (K11) Antibody	CSB-PA010403OA11crHU	Human	ELISA, WB, IF
HIST1H2BC	Acetyl-HIST1H2BC (K11) Antibody	CSB-PA010403OA11acHU	Human, Rat	ELISA, WB, IHC, ChIP
HIST1H2BC	β-hydroxybutyryl-HIST1H2BC (K11) Antibody	CSB-PA010403OA11bhbHU	Human, Rat	ELISA, WB
HIST1H2BC	Formyl-HIST1H2BC (K116) Antibody	CSB-PA010403OA116forHU	Human	ELISA, WB, IF
HIST1H2BC	Mono-methyl-HIST1H2BC (K116) Antibody	CSB-PA010403OA116me1HU	Human	ELISA, WB
HIST1H2BC	Acetyl-HIST1H2BC (K116) Antibody	CSB-PA010403OA116acHU	Human, Rat	ELISA, WB, IF, ChIP
HIST1H2BC	Glutaryl-HIST1H2BC (K116) Antibody	CSB-PA010403OA116gluHU	Human, Rat	ELISA, WB, IF
HIST1H2BC	Crotonyl-HIST1H2BC (K12) Antibody	CSB-PA010403OA12crHU	Human	ELISA, WB, IF, IP, ChIP
HIST1H2BC	Mono-methyl-HIST1H2BC (K12) Antibody	CSB-PA010403OA12me1HU	Human, Mouse, Rat	ELISA, WB, ICC, ChIP
HIST1H2BC	Acetyl-HIST1H2BC (K12) Antibody	CSB-PA010403OA12acHU	Human, Rat	ELISA, WB, ICC, IF
HIST1H2BC	2-hydroxyisobutyryl-HIST1H2BC (K12) Antibody	CSB-PA010403OA12hibHU	Human, Rat	ELISA, WB, ICC, IF
HIST1H2BC	2-hydroxyisobutyryl-HIST1H2BC (K120) Antibody	CSB-PA010403OA120hibHU	Human, Rat	ELISA, WB
HIST1H2BC	Glutaryl-HIST1H2BC (K120) Antibody	CSB-PA010403OA120gluHU	Human, Rat	ELISA, WB



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H2BC	β-hydroxybutyryl-HIST1H2BC (K120) Antibody	CSB-PA010403OA120bhbHU	Human	ELISA, IF
HIST1H2BC	acetyl-Histone H2B (K120) Antibody	CSB-PA010403OA120acHU	Human	ELISA, ICC, IF, ChIP
HIST1H2BC	Mono-methyl-HIST1H2BC (K15) Antibody	CSB-PA010403OA15me1HU	Human	ELISA, WB, ICC
HIST1H2BC	Acetyl-HIST1H2BC (K15) Antibody	CSB-PA010403OA15acHU	Human, Rat	ELISA, WB, IF, IP
HIST1H2BC	Crotonyl-HIST1H2BC (K16) Antibody	CSB-PA010403OA16crHU	Human	ELISA, WB, ICC, IF, ChIP
HIST1H2BC	Mono-methyl-HIST1H2BC (K20) Antibody	CSB-PA010403OA20me1HU	Human	ELISA, WB, ChIP
HIST1H2BC	Crotonyl-HIST1H2BC (K20) Antibody	CSB-PA010403OA20crHU	Human, Mouse, Rat	ELISA, WB, IF, ChIP
HIST1H2BC	Butyryl-HIST1H2BC (K20) Antibody	CSB-PA010403OA20butHU	Human, Rat	ELISA, WB, ICC, IF
HIST1H2BC	2-hydroxyisobutyryl-HIST1H2BC (K20) Antibody	CSB-PA010403OA20hibHU	Human, Rat	ELISA, WB, IF
HIST1H2BC	Acetyl-HIST1H2BC (K20) Antibody	CSB-PA010403OA20acHU	Human, Rat	ELISA, WB, ICC, IF, IP, ChIP
HIST1H2BC	β-hydroxybutyryl-HIST1H2BC (K20) Antibody	CSB-PA010403OA20bhbHU	Human, Rat	ELISA, WB, ICC
HIST1H2BC	Crotonyl-HIST1H2BC (K23) Antibody	CSB-PA010403OA23crHU	Human	ELISA, WB, IF, IP, ChIP
HIST1H2BC	2-hydroxyisobutyryl-HIST1H2BC (K24) Antibody	CSB-PA010403OA24hibHU	Human, Mouse, Rat	ELISA, WB, ICC, IF
HIST1H2BC	Acetyl-HIST1H2BC (K24) Antibody	CSB-PA010403OA24acHU	Human, Mouse, Rat	ELISA, WB, IHC, IP
HIST1H2BC	2-hydroxyisobutyryl-HIST1H2BC (K34) Antibody	CSB-PA010403OA34hibHU	Human, Rat	ELISA, WB, ICC, IF
HIST1H2BC	Glutaryl-HIST1H2BC (K5) Antibody	CSB-PA010403OA05gluHU	Human, Rat	ELISA, WB
HIST1H2BC	Butyryl-HIST1H2BC (K5) Antibody	CSB-PA010403OA05butHU	Human, Rat	ELISA, WB, IF, IP, ChIP
HIST1H2BC	Acetyl-HIST1H2BC (K85) Antibody	CSB-PA010403OA85acHU	Human	ELISA, IHC, ChIP



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H2BC	Phospho-HIST1H2BC (S14) Antibody	CSB-PA010403OA14phHU	Human	ELISA, WB
HIST1H2BC	Hydroxyl-HIST1H2BC (Y37) Antibody	CSB-PA010403OA37ohHU	Human, Rat	ELISA, WB
HIST1H3A	HIST1H3A (Ab-10) Antibody	CSB-PA010418PA10nphHU	Human, Mouse, Rat	ELISA, WB, IHC, IF, CHIP
HIST1H3A	HIST1H3A (Ab-10) Antibody	CSB-PA010418OA10nphHU	Human, Mouse, Rat	ELISA, WB, IHC, IF, CHIP
HIST1H3A	HIST1H3A (Ab-10) Antibody	CSB-PA010418OA10nogHU	Human, Rat	ELISA, WB, IHC, IF, CHIP
HIST1H3A	HIST1H3A (Ab-11) Antibody	CSB-PA010418OA11nphHU	Human	ELISA, WB, IHC, IF
HIST1H3A	HIST1H3A (Ab-115) Antibody	CSB-PA010418OA115nachU	Human, Mouse, Rat	ELISA, WB
HIST1H3A	HIST1H3A (Ab-122) Antibody	CSB-PA010418OA122nhibHU	Human, Mouse, Rat	ELISA, WB, IHC, IP
HIST1H3A	HIST1H3A (Ab-122) Antibody	CSB-PA010418OA122nforHU	Human, Mouse, Rat	ELISA, WB, IP
HIST1H3A	HIST1H3A (Ab-122) Antibody	CSB-PA010418OA122nbutHU	Human, Mouse	ELISA, WB, IP
HIST1H3A	HIST1H3A (Ab-122) Antibody	CSB-PA010418OA122nachU	Human, Mouse, Rat	ELISA, WB, IHC, IF, IP
HIST1H3A	HIST1H3A (Ab-122) Antibody	CSB-PA010418OA122nbhbHU	Human, Mouse, Rat	ELISA, WB, IF
HIST1H3A	Histone H3 (Ab-128) Antibody	CSB-PA010418OA128nme1HU	Human	ELISA, WB, IHC, IP
HIST1H3A	HIST1H3A (Ab-14) Antibody	CSB-PA010418PA14nachU	Human, Mouse, Rat	ELISA, WB, IHC, IF, CHIP
HIST1H3A	HIST1H3A (Ab-14) Antibody	CSB-PA010418OA14nachU	Human, Mouse, Rat	ELISA, WB, IF, CHIP
HIST1H3A	HIST1H3A (Ab-14) Antibody	CSB-PA010418OA14nsuchU	Human	ELISA, IHC, IF
HIST1H3A	HIST1H3A (Ab-14) Antibody	CSB-PA010418OA14nhibHU	Human	ELISA, IF
HIST1H3A	HIST1H3A (Ab-17) Antibody	CSB-PA010418PA17nme1HU	Human	ELISA, WB, IHC, IF, CHIP



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H3A	HIST1H3A (Ab-17) Antibody	CSB-PA010418OA17nme1HU	Human	ELISA, WB, IHC, IF, ChIP
HIST1H3A	HIST1H3A (Ab-17) Antibody	CSB-PA010418PA17nme2HU	Human, Mouse, Rat	ELISA, WB, IHC, IF, ChIP
HIST1H3A	HIST1H3A (Ab-17) Antibody	CSB-PA010418OA17nme2HU	Human	ELISA, IHC, IF
HIST1H3A	HIST1H3A (Ab-18) Antibody	CSB-PA010418PA18nacHU	Human	ELISA, IHC, IF
HIST1H3A	HIST1H3A (Ab-18) Antibody	CSB-PA010418OA18nacHU	Human	ELISA, IHC, IF
HIST1H3A	HIST1H3A (Ab-18) Antibody	CSB-PA010418OA18nforHU	Human	ELISA, IHC, IF
HIST1H3A	HIST1H3A (Ab-18) Antibody	CSB-PA010418OA18ncrHU	Human	ELISA, WB
HIST1H3A	HIST1H3A (Ab-18) Antibody	CSB-PA010418OA18nbutHU	Human, Mouse, Rat	ELISA, WB, IHC, IF, IP
HIST1H3A	HIST1H3A (Ab-18) Antibody	CSB-PA010418OA18nbhbHU	Human, Mouse, Rat	ELISA, WB, IHC, IF
HIST1H3A	HIST1H3A (Ab-23) Antibody	CSB-PA010418PA23nacHU	Human, Rat	ELISA, WB, IF
HIST1H3A	HIST1H3A (Ab-23) Antibody	CSB-PA010418OA23nacHU	Human, Mouse	ELISA, WB, IHC, IF
HIST1H3A	HIST1H3A (Ab-23) Antibody	CSB-PA010418OA23ncrHU	Human	ELISA, IHC, IF
HIST1H3A	HIST1H3A (Ab-23) Antibody	CSB-PA010418OA23nforHU	Human, Mouse, Rat	ELISA, WB, IHC, IF, IP
HIST1H3A	HIST1H3A (Ab-23) Antibody	CSB-PA010418OA23nhibHU	Human, Mouse, Rat	ELISA, WB, IHC, ChIP
HIST1H3A	HIST1H3A (Ab-23) Antibody	CSB-PA010418OA23nbhbHU	Human	ELISA, IHC, ChIP
HIST1H3A	HIST1H3A (Ab-23) Antibody	CSB-PA010418OA23nsuchHU	Human	ELISA, IHC, IF
HIST1H3A	HIST1H3A (Ab-23) Antibody	CSB-PA010418OA23nbutHU	Human, Mouse, Rat	ELISA, WB, IHC, IF, ChIP
HIST1H3A	HIST1H3A (Ab-23) Antibody	CSB-PA010418OA23npropHU	Human	ELISA, WB, IHC, IF



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H3A	HIST1H3A (Ab-27) Antibody	CSB-PA010418PA27nme2HU	Human	ELISA, IHC, IF
HIST1H3A	HIST1H3A (Ab-27) Antibody	CSB-PA010418OA27nme2HU	Human	ELISA, IHC, IF
HIST1H3A	HIST1H3A (Ab-27) Antibody	CSB-PA010418PA27nme1HU	Human, Mouse, Rat	ELISA, WB, IHC, IF, CHIP
HIST1H3A	HIST1H3A (Ab-27) Antibody	CSB-PA010418OA27nme1HU	Human	ELISA, IHC
HIST1H3A	HIST1H3A (Ab-27) Antibody	CSB-PA010418PA27nme3HU	Human	ELISA, WB, IHC, IF
HIST1H3A	HIST1H3A (Ab-27) Antibody	CSB-PA010418OA27nme3HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
HIST1H3A	HIST1H3A (Ab-27) Antibody	CSB-PA010418OA27nhibHU	Human	ELISA, IHC
HIST1H3A	HIST1H3A (Ab-27) Antibody	CSB-PA010418OA27nacHU	Human, Mouse, Rat	ELISA, WB, IF
HIST1H3A	HIST1H3A (Ab-27) Antibody	CSB-PA010418OA27nbutHU	Human, Mouse, Rat	ELISA, WB, IHC, IF, CHIP
HIST1H3A	HIST1H3A (Ab-27) Antibody	CSB-PA010418OA27nbhbHU	Human, Mouse	ELISA, WB, IHC
HIST1H3A	HIST1H3A (Ab-28) Antibody	CSB-PA010418PA28nacHU	Human	ELISA, IHC, IF
HIST1H3A	HIST1H3A (Ab-28) Antibody	CSB-PA010418OA28nacHU	Human, Mouse, Rat	ELISA, WB, IHC, IF
HIST1H3A	HIST1H3A (Ab-28) Antibody	CSB-PA010418OA28nphHU	Human	ELISA, IHC, IF
HIST1H3A	HIST1H3A (Ab-3) Antibody	CSB-PA010418OA03nphHU	Human, Rat	ELISA, WB, IHC, IF, IP
HIST1H3A	HIST1H3A (Ab-36) Antibody	CSB-PA010418PA36nme2HU	Human	ELISA, IHC, IF
HIST1H3A	HIST1H3A (Ab-36) Antibody	CSB-PA010418OA36nme2HU	Human	ELISA, IHC, IF
HIST1H3A	HIST1H3A (Ab-36) Antibody	CSB-PA010418PA36nme1HU	Human, Mouse	ELISA, WB, IHC, IF, CHIP
HIST1H3A	HIST1H3A (Ab-36) Antibody	CSB-PA010418OA36nme1HU	Human, Mouse, Rat	ELISA, WB, IHC, IF, CHIP



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H3A	HIST1H3A (Ab-36) Antibody	CSB-PA010418PA36nme3HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
HIST1H3A	HIST1H3A (Ab-36) Antibody	CSB-PA010418OA36nme3HU	Human, Mouse, Rat	ELISA, WB, IHC, IF, ChIP
HIST1H3A	HIST1H3A (Ab-36) Antibody	CSB-PA010418OA36nacHU	Human	ELISA, WB, IHC
HIST1H3A	HIST1H3A (Ab-37) Antibody	CSB-PA010418OA37nacHU	Human	ELISA, WB, IF, ChIP
HIST1H3A	HIST1H3A (Ab-4) Antibody	CSB-PA010418PA04nme1HU	Human	ELISA, WB, IHC, IF
HIST1H3A	HIST1H3A (Ab-4) Antibody	CSB-PA010418OA04nme1HU	Human, Mouse	ELISA, WB, IHC, IF, ChIP
HIST1H3A	HIST1H3A (Ab-4) Antibody	CSB-PA010418PA04nme2HU	Human, Mouse	ELISA, WB, IHC, IF, ChIP
HIST1H3A	HIST1H3A (Ab-4) Antibody	CSB-PA010418OA04nme2HU	Human, Mouse	ELISA, WB, IHC, IF, ChIP
HIST1H3A	HIST1H3A (Ab-4) Antibody	CSB-PA010418PA04nme3HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
HIST1H3A	HIST1H3A (Ab-4) Antibody	CSB-PA010418NA04nme3HU	Human, Mouse, Rat	ELISA, WB, IHC, IF, IP
HIST1H3A	HIST1H3A (Ab-4) Antibody	CSB-PA010418OA04nacHU	Human, Rat, Mouse	ELISA, WB, IHC, IF, ChIP
HIST1H3A	HIST1H3A (Ab-4) Antibody	CSB-PA010418OA04ncrHU	Human	ELISA, IHC, IF
HIST1H3A	HIST1H3A (Ab-4) Antibody	CSB-PA010418OA04nhibHU	Human, Mouse, Rat	ELISA, WB, IHC, IF
HIST1H3A	HIST1H3A (Ab-4) Antibody	CSB-PA010418OA04nbhbHU	Human	ELISA, IHC
HIST1H3A	Histone H3 (Ab-42) Antibody	CSB-PA010418OA42nme2HU	Human	ELISA, IHC
HIST1H3A	HIST1H3A (Ab-45) Antibody	CSB-PA010418OA45nphHU	Human, Rat, Mouse	ELISA, WB, IHC, IF
HIST1H3A	HIST1H3A (Ab-56) Antibody	CSB-PA010418OA56nme1HU	Human	ELISA, WB
HIST1H3A	HIST1H3A (Ab-56) Antibody	CSB-PA010418OA56nhibHU	Human	ELISA, IHC, IF



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H3A	HIST1H3A (Ab-56) Antibody	CSB-PA010418OA56nbhbHU	Human	ELISA, IHC
HIST1H3A	HIST1H3A (Ab-6) Antibody	CSB-PA010418OA06nphHU	Human, Rat	ELISA, WB, IHC, IP, ChIP
HIST1H3A	Histone H3 (Ab-63) Antibody	CSB-PA010418OA63nme1HU	Human	ELISA, IHC
HIST1H3A	HIST1H3A (Ab-64) Antibody	CSB-PA010418OA64nacHU	Human	ELISA, IHC
HIST1H3A	HIST1H3A (Ab-64) Antibody	CSB-PA010418OA64nme1HU	Human	ELISA, IF
HIST1H3A	HIST1H3A (Ab-64) Antibody	CSB-PA010418OA64nforHU	Human	ELISA, IHC
HIST1H3A	HIST1H3A (Ab-79) Antibody	CSB-PA010418OA79nacHU	Human	ELISA, WB, IHC
HIST1H3A	HIST1H3A (Ab-79) Antibody	CSB-PA010418OA79nme1HU	Human	ELISA, WB, IHC
HIST1H3A	HIST1H3A (Ab-79) Antibody	CSB-PA010418OA79nbhbHU	Human, Mouse	ELISA, WB, IHC
HIST1H3A	HIST1H3A (Ab-9) Antibody	CSB-PA010418PA09nacHU	Human, Mouse	ELISA, WB, IHC, IF, ChIP
HIST1H3A	HIST1H3A (Ab-9) Antibody	CSB-PA010418OA09nacHU	Human	ELISA, WB, IHC, IF
HIST1H3A	HIST1H3A (Ab-9) Antibody	CSB-PA010418PA09nme1HU	Human	ELISA, WB, IHC, IF
HIST1H3A	HIST1H3A (Ab-9) Antibody	CSB-PA010418OA09nme1HU	Human	ELISA, WB, IHC, IF, ChIP
HIST1H3A	HIST1H3A (Ab-9) Antibody	CSB-PA010418PA09nme2HU	Human, Mouse	ELISA, WB, IF
HIST1H3A	HIST1H3A (Ab-9) Antibody	CSB-PA010418OA09nme2HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
HIST1H3A	HIST1H3A (Ab-9) Antibody	CSB-PA010418PA09nme3HU	Human, Mouse, Rat	ELISA, WB, IF, ChIP
HIST1H3A	HIST1H3A (Ab-9) Antibody	CSB-PA010418OA09nme3HU	Human, Mouse, Rat	ELISA, WB, IHC, IF
HIST1H3A	HIST1H3A (Ab-9) Antibody	CSB-PA010418OA09nhhbHU	Human, Mouse, Rat	ELISA, WB, IHC, IF



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H3A	HIST1H3A (Ab-9) Antibody	CSB-PA010418OA09nbhbHU	Human	ELISA, WB
HIST1H3A	HIST1H3A (Ab-9) Antibody	CSB-PA010418OA09ncrHU	Human	ELISA, WB, IHC, IF
HIST1H3A	HIST1H3A (Ab-9) Antibody	CSB-PA010418OA09nbuthHU	Human	ELISA, WB, IHC
HIST1H3A	Acetyl-HIST1H3A (K115) Antibody	CSB-PA010418OA115achHU	Human, Rat	ELISA, WB, ICC
HIST1H3A	2-hydroxyisobutyryl-HIST1H3A (K122) Antibody	CSB-PA010418OA122hibHU	Human, Mouse, Rat	ELISA, WB
HIST1H3A	Formyl-HIST1H3A (K122) Antibody	CSB-PA010418OA122forHU	Human, Mouse, Rat	ELISA, WB
HIST1H3A	Butyryl-HIST1H3A (K122) Antibody	CSB-PA010418OA122butHU	Human, Rat	ELISA, WB, ChIP
HIST1H3A	Acetyl-HIST1H3A (K122) Antibody	CSB-PA010418OA122achHU	Human, Mouse, Rat	ELISA, WB, ICC
HIST1H3A	β-hydroxybutyryl-HIST1H3A (K122) Antibody	CSB-PA010418OA122bhbHU	Human, Rat	ELISA, WB, ICC
HIST1H3A	Acetyl-HIST1H3A (K14) Antibody	CSB-PA010418PA14achHU	Human	ELISA, IF
HIST1H3A	Acetyl-HIST1H3A (K14) Antibody	CSB-PA010418OA14achHU	Human	ELISA, ICC, IF
HIST1H3A	2-hydroxyisobutyryl-HIST1H3A (K14) Antibody	CSB-PA010418OA14hibHU	Human	ELISA, IF, ChIP
HIST1H3A	Acetyl-HIST1H3A (K18) Antibody	CSB-PA010418OA18achHU	Human	ELISA, ICC
HIST1H3A	Crotonyl-HIST1H3A (K18) Antibody	CSB-PA010418OA18crHU	Human	ELISA, IF, ChIP
HIST1H3A	Formyl-HIST1H3A (K18) Antibody	CSB-PA010418OA18forHU	Human	ELISA, IF
HIST1H3A	2-hydroxyisobutyryl-HIST1H3A (K18) Antibody	CSB-PA010418OA18hibHU	Human, Rat	ELISA, WB, ICC, IF
HIST1H3A	Butyryl-HIST1H3A (K18) Antibody	CSB-PA010418OA18butHU	Human, Mouse, Rat	ELISA, WB, ICC, IF
HIST1H3A	β-hydroxybutyryl-HIST1H3A (K18) Antibody	CSB-PA010418OA18bhbHU	Human, Rat	ELISA, WB, ICC, IF, ChIP



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H3A	2-hydroxyisobutyryl-HIST1H3A (K18) Antibody	CSB-PA010418PA18hibHU	Human	ELISA, WB, ICC, IF, ChIP
HIST1H3A	Acetyl-HIST1H3A (K23) Antibody	CSB-PA010418OA23achU	Human	ELISA, WB, ICC, IF
HIST1H3A	Formyl-HIST1H3A (K23) Antibody	CSB-PA010418OA23forHU	Human, Rat	ELISA, WB, IF
HIST1H3A	2-hydroxyisobutyryl-HIST1H3A (K23) Antibody	CSB-PA010418OA23hibHU	Human, Rat	ELISA, WB, IF, ChIP
HIST1H3A	β-hydroxybutyryl-HIST1H3A (K23) Antibody	CSB-PA010418OA23bhbHU	Human, Rat	ELISA, WB, ChIP
HIST1H3A	Butyryl-HIST1H3A (K23) Antibody	CSB-PA010418OA23butHU	Human, Rat	ELISA, WB, IHC, IF, IP, ChIP
HIST1H3A	Propionyl-HIST1H3A (K23) Antibody	CSB-PA010418OA23propHU	Human, Mouse, Rat	ELISA, WB, ICC, IF, ChIP
HIST1H3A	Di-methyl-HIST1H3A (K27) Antibody	CSB-PA010418OA27me2HU	Human	ELISA, IF
HIST1H3A	Mono-methyl-HIST1H3A (K27) Antibody	CSB-PA010418PA27me1HU	Human	ELISA, ICC, IF
HIST1H3A	Tri-methyl-HIST1H3A (K27) Antibody	CSB-PA010418PA27me3HU	Human	ELISA, IF
HIST1H3A	Tri-methyl-HIST1H3A (K27) Antibody	CSB-PA010418OA27me3HU	Human, Mouse, Rat	ELISA, WB, ICC, IF
HIST1H3A	2-hydroxyisobutyryl-HIST1H3A (K27) Antibody	CSB-PA010418OA27hibHU	Human, Rat	ELISA, WB, ICC, IF
HIST1H3A	Acetyl-HIST1H3A (K27) Antibody	CSB-PA010418OA27achU	Human, Mouse, Rat	ELISA, WB, IHC, IF, IP, ChIP
HIST1H3A	Butyryl-HIST1H3A (K27) Antibody	CSB-PA010418OA27butHU	Human, Mouse, Rat	ELISA, WB, IHC, IF, ChIP
HIST1H3A	β-hydroxybutyryl-HIST1H3A (K27) Antibody	CSB-PA010418OA27bhbHU	Human, Rat	ELISA, WB, ICC, ChIP
HIST1H3A	Di-methyl-HIST1H3A (K36) Antibody	CSB-PA010418PA36me2HU	Human, Mouse, Rat	ELISA, WB, IF
HIST1H3A	Di-methyl-HIST1H3A (K36) Antibody	CSB-PA010418OA36me2HU	Human	ELISA, IF
HIST1H3A	Mono-methyl-HIST1H3A (K36) Antibody	CSB-PA010418PA36me1HU	Human, Mouse, Rat	ELISA, WB, ICC, IF



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H3A	Mono-methyl-HIST1H3A (K36) Antibody	CSB-PA010418OA36me1HU	Human, Mouse, Rat	ELISA, WB, IF
HIST1H3A	Tri-methyl-HIST1H3A (K36) Antibody	CSB-PA010418PA36me3HU	Human, Mouse, Rat	ELISA, WB, ICC, IF
HIST1H3A	Tri-methyl-HIST1H3A (K36) Antibody	CSB-PA010418OA36me3HU	Human, Mouse, Rat	ELISA, WB, IF
HIST1H3A	Acetyl-HIST1H3A (K36) Antibody	CSB-PA010418OA36acHU	Human, Rat	ELISA, WB, IF, ChIP
HIST1H3A	2-hydroxyisobutyryl-HIST1H3A (K36) Antibody	CSB-PA010418OA36hibHU	Human	ELISA, IF
HIST1H3A	Acetyl-HIST1H3A (K37) Antibody	CSB-PA010418OA37acHU	Human, Rat	ELISA, WB, IF, ChIP
HIST1H3A	Mono-methyl-HIST1H3A (K4) Antibody	CSB-PA010418PA04me1HU	Human, Mouse, Rat	ELISA, WB, IF
HIST1H3A	Mono-methyl-HIST1H3A (K4) Antibody	CSB-PA010418OA04me1HU	Human	ELISA, IF
HIST1H3A	Di-methyl-HIST1H3A (K4) Antibody	CSB-PA010418PA04me2HU	Human	ELISA, ICC, IF
HIST1H3A	Di-methyl-HIST1H3A (K4) Antibody	CSB-PA010418OA04me2HU	Human	ELISA, IF
HIST1H3A	Tri-methyl-HIST1H3A (K4) Antibody	CSB-PA010418PA04me3HU	Human, Mouse, Rat	ELISA, WB, ICC, IF
HIST1H3A	Tri-methyl-HIST1H3A (K4) Antibody	CSB-PA010418NA04me3HU	Human, Mouse, Rat	ELISA, WB, ICC, IF
HIST1H3A	Acetyl-HIST1H3A (K4) Antibody	CSB-PA010418OA04acHU	Human	ELISA, WB, IF, IP, ChIP
HIST1H3A	Crotonyl-HIST1H3A (K4) Antibody	CSB-PA010418OA04crHU	Human, Mouse, Rat	ELISA, WB, IP, ChIP
HIST1H3A	2-hydroxyisobutyryl-HIST1H3A (K4) Antibody	CSB-PA010418OA04hibHU	Human, Rat	ELISA, WB, ICC, IF
HIST1H3A	β-hydroxybutyryl-HIST1H3A (K4) Antibody	CSB-PA010418OA04hbHU	Human, Rat	ELISA, WB, ICC, ChIP
HIST1H3A	2-hydroxyisobutyryl-HIST1H3A (K56) Antibody	CSB-PA010418OA56hibHU	Human	ELISA, WB, IF
HIST1H3A	Crotonyl-HIST1H3A (K56) Antibody	CSB-PA010418OA56crHU	Human, Mouse, Rat	ELISA, WB



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H3A	Acetyl-HIST1H3A (K56) Antibody	CSB-PA010418OA56acHU	Human, Rat	ELISA, WB
HIST1H3A	β-hydroxybutyryl-HIST1H3A (K56) Antibody	CSB-PA010418OA56bhbHU	Human	ELISA, ICC, IF
HIST1H3A	Malonyl-HIST1H3A (K56) Antibody	CSB-PA010418OA56malHU	Human, Rat	ELISA, WB
HIST1H3A	Acetyl-HIST1H3A (K64) Antibody	CSB-PA010418OA64acHU	Human	ELISA, IF
HIST1H3A	Mono-methyl-HIST1H3A (K64) Antibody	CSB-PA010418OA64me1HU	Human	ELISA, WB
HIST1H3A	Acetyl-HIST1H3A (K79) Antibody	CSB-PA010418OA79acHU	Human	ELISA, IF, ChIP
HIST1H3A	β-hydroxybutyryl-HIST1H3A (K79) Antibody	CSB-PA010418OA79bhbHU	Human, Rat	ELISA, WB, ICC, IF
HIST1H3A	Butyryl-HIST1H3A (K79) Antibody	CSB-PA010418OA79butHU	Human, Rat	ELISA, WB, ICC
HIST1H3A	Acetyl-HIST1H3A (K9) Antibody	CSB-PA010418PA09acHU	Human	ELISA, WB, ICC, IF, ChIP
HIST1H3A	Acetyl-HIST1H3A (K9) Antibody	CSB-PA010418OA09acHU	Human	ELISA, ICC
HIST1H3A	Mono-methyl-HIST1H3A (K9) Antibody	CSB-PA010418PA09me1HU	Human, Mouse, Rat	ELISA, WB, ICC, IF
HIST1H3A	Mono-methyl-HIST1H3A (K9) Antibody	CSB-PA010418OA09me1HU	Human, Mouse, Rat	ELISA, WB, IF
HIST1H3A	Di-methyl-HIST1H3A (K9) Antibody	CSB-PA010418PA09me2HU	Human, Mouse, Rat	ELISA, WB, IF
HIST1H3A	Di-methyl-HIST1H3A (K9) Antibody	CSB-PA010418OA09me2HU	Human, Mouse, Rat	ELISA, WB, IF
HIST1H3A	Tri-methyl-HIST1H3A (K9) Antibody	CSB-PA010418PA09me3HU	Human	ELISA, IF
HIST1H3A	Tri-methyl-HIST1H3A (K9) Antibody	CSB-PA010418OA09me3HU	Human	ELISA, IF
HIST1H3A	β-hydroxybutyryl-HIST1H3A (K9) Antibody	CSB-PA010418OA09bhbHU	Human, Rat	ELISA, WB, ICC, ChIP
HIST1H3A	Crotonyl-HIST1H3A (K9) Antibody	CSB-PA010418OA09crHU	Human, Mouse, Rat	ELISA, WB, ICC, IP, ChIP



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H3A	Butyryl-HIST1H3A (K9) Antibody	CSB-PA010418OA09butHU	Human, Rat	ELISA, WB, IHC, IF, ChIP
HIST1H3A	Mono-methyl-HIST1H3A (R17) Antibody	CSB-PA010418OA17me1HU	Human	ELISA, IF
HIST1H3A	Di-methyl-HIST1H3A (R17) Antibody	CSB-PA010418PA17me2HU	Human	ELISA, IF
HIST1H3A	Di-methyl-HIST1H3A (R17) Antibody	CSB-PA010418OA17me2HU	Human	ELISA, IF
HIST1H3A	phospho-HIST1H3A (S10) Antibody	CSB-PA010418PA10phHU	Human	ELISA, WB, IHC, IF, IP
HIST1H3A	Acetyl-HIST1H3A (S28) Antibody	CSB-PA010418OA28acHU	Human	ELISA, IF
HIST1H3A	Phospho-HIST1H3A (S28) Antibody	CSB-PA010418OA28phHU	Human	ELISA, IF, ChIP
HIST1H3A	Acetyl-HIST1H3A (T22) Antibody	CSB-PA010418OA22acHU	Human	ELISA, ChIP
HIST1H3A	Phospho-HIST1H3A (T3) Antibody	CSB-PA010418OA03phHU	Human, Mouse, Rat	ELISA, WB, IP, ChIP
HIST1H3A	Phospho-HIST1H3A (T6) Antibody	CSB-PA010418OA06phHU	Human	ELISA, IF, ChIP
HIST1H3A	Phospho-HIST1H3A (T80) Antibody	CSB-PA010418OA80phHU	Human	ELISA, IF, ChIP
HIST1H3A	Phospho-HIST1H3A (Y41) Antibody	CSB-PA010418OA41phHU	Human	ELISA, IF
HIST1H4A	HIST1H4A (Ab-1) Antibody	CSB-PA010429OA01nphHU	Human, Mouse, Rat	ELISA, WB, IHC, IF, ChIP
HIST1H4A	HIST1H4A (Ab-12) Antibody	CSB-PA010429PA12nacHU	Human, Mouse	ELISA, WB, IHC, IF, IP, ChIP
HIST1H4A	HIST1H4A (Ab-12) Antibody	CSB-PA010429OA12nacHU	Human	ELISA, WB, IHC, IF, ChIP
HIST1H4A	HIST1H4A (Ab-12) Antibody	CSB-PA010429OA12nbioHU	Human	ELISA, IHC
HIST1H4A	HIST1H4A (Ab-12) Antibody	CSB-PA010429OA12ncrHU	Human	ELISA, WB, IF, IP
HIST1H4A	HIST1H4A (Ab-12) Antibody	CSB-PA010429OA12nsucHU	Human	ELISA, IF



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H4A	HIST1H4A (Ab-12) Antibody	CSB-PA010429OA12npropHU	Human	ELISA, WB, IHC, IP, ChIP
HIST1H4A	HIST1H4A (Ab-12) Antibody	CSB-PA010429OA12nbutHU	Human	ELISA, WB, IHC, IP, ChIP
HIST1H4A	HIST1H4A (Ab-12) Antibody	CSB-PA010429OA12nhibHU	Human	ELISA, WB, IHC, IF, IP
HIST1H4A	HIST1H4A (Ab-12) Antibody	CSB-PA010429OA12nbhbHU	Human	ELISA, WB, IHC, IP
HIST1H4A	HIST1H4A (Ab-16) Antibody	CSB-PA010429PA16nacHU	Human	ELISA, IHC, IF
HIST1H4A	HIST1H4A (Ab-16) Antibody	CSB-PA010429NA16nacHU	Human	ELISA, WB, IHC, IF
HIST1H4A	HIST1H4A (Ab-16) Antibody	CSB-PA010429OA16ncrHU	Human	ELISA, IHC, IF
HIST1H4A	HIST1H4A (Ab-16) Antibody	CSB-PA010429OA16nbutHU	Human	ELISA, IF, ChIP
HIST1H4A	HIST1H4A (Ab-16) Antibody	CSB-PA010429OA16npropHU	Human	ELISA, WB, IHC
HIST1H4A	HIST1H4A (Ab-20) Antibody	CSB-PA010429OA20nme3HU	Human	ELISA, IHC, IF, ChIP
HIST1H4A	HIST1H4A (Ab-3) Antibody	CSB-PA010429PA03nme1HU	Human	ELISA, WB, IF, ChIP
HIST1H4A	HIST1H4A (Ab-3) Antibody	CSB-PA010429OA03nme1HU	Human	ELISA, WB, IF, ChIP
HIST1H4A	HIST1H4A (Ab-31) Antibody	CSB-PA010429OA31npropHU	Human, Mouse	ELISA, WB, IHC
HIST1H4A	HIST1H4A (Ab-31) Antibody	CSB-PA010429OA31nhibHU	Human, Rat	ELISA, WB, IHC
HIST1H4A	HIST1H4A (Ab-31) Antibody	CSB-PA010429OA31nsucHU	Human	ELISA, WB, IHC, ChIP
HIST1H4A	HIST1H4A (Ab-31) Antibody	CSB-PA010429OA31nacHU	Human	ELISA, IHC
HIST1H4A	HIST1H4A (Ab-31) Antibody	CSB-PA010429OA31nbhbHU	Human	ELISA, IF
HIST1H4A	HIST1H4A (Ab-35) Antibody	CSB-PA010429OA35nme1HU	Human	ELISA, IHC, IF, ChIP



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H4A	HIST1H4A (Ab-44) Antibody	CSB-PA010429OA44npropHU	Human	ELISA, WB
HIST1H4A	HIST1H4A (Ab-44) Antibody	CSB-PA010429OA44nbutHU	Human	ELISA, IHC
HIST1H4A	HIST1H4A (Ab-44) Antibody	CSB-PA010429OA44nhibHU	Human	ELISA, IHC
HIST1H4A	HIST1H4A (Ab-5) Antibody	CSB-PA010429OA05nacHU	Human, Mouse	ELISA, WB, IHC, IF
HIST1H4A	HIST1H4A (Ab-5) Antibody	CSB-PA010429OA05ncrHU	Human, Rat	ELISA, WB, IHC, IF, IP, ChIP
HIST1H4A	HIST1H4A (Ab-5) Antibody	CSB-PA010429OA05nme1HU	Human	ELISA, IF, ChIP
HIST1H4A	HIST1H4A (Ab-5) Antibody	CSB-PA010429OA05nbutHU	Human, Mouse	ELISA, WB, IHC, IF, IP
HIST1H4A	HIST1H4A (Ab-5) Antibody	CSB-PA010429OA05npropHU	Human, Mouse	ELISA, WB, IHC, IF, IP, ChIP
HIST1H4A	HIST1H4A (Ab-5) Antibody	CSB-PA010429OA05nhibHU	Human	ELISA, WB, IHC, IF, IP
HIST1H4A	HIST1H4A (Ab-5) Antibody	CSB-PA010429OA05nbhbHU	Human	ELISA, WB, IHC, IF, IP
HIST1H4A	HIST1H4A (Ab-55) Antibody	CSB-PA010429OA55nme1HU	Human	ELISA, WB, IHC
HIST1H4A	HIST1H4A (Ab-59) Antibody	CSB-PA010429OA59nme1HU	Human	ELISA, WB, IHC, IF, ChIP
HIST1H4A	HIST1H4A (Ab-59) Antibody	CSB-PA010429OA59nforHU	Human	ELISA, IHC
HIST1H4A	HIST1H4A (Ab-59) Antibody	CSB-PA010429OA59nhibHU	Human	ELISA, IHC
HIST1H4A	HIST1H4A (Ab-77) Antibody	CSB-PA010429OA77nme1HU	Human	ELISA, IHC, IF
HIST1H4A	HIST1H4A (Ab-77) Antibody	CSB-PA010429OA77nsuchHU	Human	ELISA, IHC
HIST1H4A	HIST1H4A (Ab-77) Antibody	CSB-PA010429OA77nhibHU	Human	ELISA, WB, IHC
HIST1H4A	HIST1H4A (Ab-77) Antibody	CSB-PA010429OA77nbutHU	Human, Mouse	ELISA, WB, IHC



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H4A	HIST1H4A (Ab-77) Antibody	CSB-PA010429OA77nbhbHU	Human	ELISA, IHC
HIST1H4A	HIST1H4A (Ab-79) Antibody	CSB-PA010429OA79nforHU	Human	ELISA, WB, IHC
HIST1H4A	HIST1H4A (Ab-79) Antibody	CSB-PA010429OA79nhibHU	Human	ELISA, IHC, ChIP
HIST1H4A	HIST1H4A (Ab-79) Antibody	CSB-PA010429OA79npropHU	Human	ELISA, WB, IHC, IF
HIST1H4A	HIST1H4A (Ab-79) Antibody	CSB-PA010429OA79nbutHU	Human	ELISA, IHC
HIST1H4A	HIST1H4A (Ab-8) Antibody	CSB-PA010429PA08nacHU	Human, Rat	ELISA, WB, IHC, IF, IP, ChIP
HIST1H4A	HIST1H4A (Ab-8) Antibody	CSB-PA010429OA08nacHU	Human, Mouse	ELISA, WB, IHC, IF
HIST1H4A	HIST1H4A (Ab-8) Antibody	CSB-PA010429OA08ncrHU	Human	ELISA, WB, IHC, IF
HIST1H4A	HIST1H4A (Ab-8) Antibody	CSB-PA010429OA08nbutHU	Human, Mouse, Rat	ELISA, WB, IHC, IF
HIST1H4A	HIST1H4A (Ab-8) Antibody	CSB-PA010429OA08npropHU	Human	ELISA, WB, IHC, IF, IP, ChIP
HIST1H4A	HIST1H4A (Ab-8) Antibody	CSB-PA010429OA08nhibHU	Human	ELISA, WB, IHC, IF
HIST1H4A	HIST1H4A (Ab-8) Antibody	CSB-PA010429OA08nbhbHU	Human	ELISA, WB, IHC, IF, ChIP
HIST1H4A	HIST1H4A (Ab-88) Antibody	CSB-PA010429OA88nohHU	Human	ELISA, IHC, IF
HIST1H4A	HIST1H4A (Ab-91) Antibody	CSB-PA010429OA91nacHU	Human	ELISA, IHC, IF
HIST1H4A	HIST1H4A (Ab-91) Antibody	CSB-PA010429OA91npropHU	Human	ELISA, IHC, IF
HIST1H4A	HIST1H4A (Ab-91) Antibody	CSB-PA010429OA91nbutHU	Human, Mouse	ELISA, WB, IHC, ChIP
HIST1H4A	HIST1H4A (Ab-91) Antibody	CSB-PA010429OA91nforHU	Human	ELISA, IHC, IF
HIST1H4A	HIST1H4A (Ab-91) Antibody	CSB-PA010429OA91nbhbHU	Human	ELISA, IHC



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H4A	Acetyl-HIST1H4A (K12) Antibody	CSB-PA010429PA12acHU	Human	ELISA, ICC, IF, ChIP
HIST1H4A	Acetyl-HIST1H4A (K12) Antibody	CSB-PA010429OA12acHU	Human	ELISA, ICC, IF, ChIP
HIST1H4A	Crotonyl-HIST1H4A (K12) Antibody	CSB-PA010429OA12crHU	Human, Mouse, Rat	ELISA, WB, IF
HIST1H4A	Succinyl-HIST1H4A (K12) Antibody	CSB-PA010429OA12sucHU	Human	ELISA, WB, IF
HIST1H4A	Formyl-HIST1H4A (K12) Antibody	CSB-PA010429OA12forHU	Human, Rat	ELISA, WB, ICC, IF
HIST1H4A	Butyryl-HIST1H4A (K12) Antibody	CSB-PA010429OA12butHU	Human	ELISA, WB, ICC, IF, IP, ChIP
HIST1H4A	2-hydroxyisobutyryl-HIST1H4A (K12) Antibody	CSB-PA010429PA12hibHU	Human	ELISA, WB, ICC, IF, IP, ChIP
HIST1H4A	Acetyl-HIST1H4A (K16) Antibody	CSB-PA010429NA16acHU	Human	ELISA, WB, ICC, IF, ChIP
HIST1H4A	Crotonyl-HIST1H4A (K16) Antibody	CSB-PA010429OA16crHU	Human, Mouse, Rat	ELISA, WB, ICC, IF
HIST1H4A	Butyryl-HIST1H4A (K16) Antibody	CSB-PA010429OA16butHU	Human	ELISA, WB, IF, ChIP
HIST1H4A	Propionyl-HIST1H4A (K16) Antibody	CSB-PA010429OA16propHU	Human	ELISA, WB, ChIP
HIST1H4A	2-hydroxyisobutyryl-HIST1H4A (K16) Antibody	CSB-PA010429OA16hibHU	Human, Rat	ELISA, WB, IF, IP, ChIP
HIST1H4A	Acetyl-HIST1H4A (K16) Antibody	CSB-PA010429PA16acHU	Human	ELISA, ChIP
HIST1H4A	Tri-methyl-HIST1H4A (K20) Antibody	CSB-PA010429OA20me3HU	Human	ELISA, WB
HIST1H4A	Acetyl-HIST1H4A (K20) Antibody	CSB-PA010429OA20acHU	Human	ELISA, WB, IF
HIST1H4A	Propionyl-HIST1H4A (K31) Antibody	CSB-PA010429OA31propHU	Human, Mouse, Rat	ELISA, WB, ChIP
HIST1H4A	Formyl-HIST1H4A (K31) Antibody	CSB-PA010429OA31forHU	Human, Rat	ELISA, WB
HIST1H4A	2-hydroxyisobutyryl-HIST1H4A (K31) Antibody	CSB-PA010429OA31hibHU	Human, Rat	ELISA, WB, IF



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H4A	Butyryl-HIST1H4A (K31) Antibody	CSB-PA010429OA31butHU	Human	ELISA, ICC, IF
HIST1H4A	Acetyl-HIST1H4A (K31) Antibody	CSB-PA010429OA31acHU	Human, Rat	ELISA, WB, ICC, IF, ChIP
HIST1H4A	Acetyl-HIST1H4A (K5) Antibody	CSB-PA010429PA05acHU	Human	ELISA, ChIP
HIST1H4A	Acetyl-HIST1H4A (K5) Antibody	CSB-PA010429OA05acHU	Human	ELISA, IF
HIST1H4A	Crotonyl-HIST1H4A (K5) Antibody	CSB-PA010429OA05crHU	Human, Mouse, Rat	ELISA, WB, IF, IP, ChIP
HIST1H4A	Mono-methyl-HIST1H4A (K5) Antibody	CSB-PA010429OA05me1HU	Human	ELISA, IF, ChIP
HIST1H4A	Butyryl-HIST1H4A (K5) Antibody	CSB-PA010429OA05butHU	Human, Mouse, Rat	ELISA, WB, IF, IP, ChIP
HIST1H4A	Propionyl-HIST1H4A (K5) Antibody	CSB-PA010429OA05propHU	Human, Mouse, Rat	ELISA, WB, ICC, IP, ChIP
HIST1H4A	2-hydroxyisobutyryl-HIST1H4A (K5) Antibody	CSB-PA010429OA05hibHU	Human	ELISA, WB, ICC, IF
HIST1H4A	β hydroxybutyryl-HIST1H4A (K5) Antibody	CSB-PA010429OA05bhbHU	Human	ELISA, WB
HIST1H4A	Mono-methyl-HIST1H4A (K59) Antibody	CSB-PA010429OA59me1HU	Human	ELISA, WB, IF
HIST1H4A	Formyl-HIST1H4A (K77) Antibody	CSB-PA010429OA77forHU	Human, Rat	ELISA, WB, IHC, ChIP
HIST1H4A	2-hydroxyisobutyryl-HIST1H4A (K77) Antibody	CSB-PA010429OA77hibHU	Human, Rat	ELISA, WB
HIST1H4A	Acetyl-HIST1H4A (K77) Antibody	CSB-PA010429OA77acHU	Human, Mouse, Rat	ELISA, WB
HIST1H4A	Formyl-HIST1H4A (K79) Antibody	CSB-PA010429OA79forHU	Human, Rat	ELISA, WB, ICC
HIST1H4A	Acetyl-HIST1H4A (K79) Antibody	CSB-PA010429OA79acHU	Human	ELISA, ICC, IF, ChIP
HIST1H4A	Acetyl-HIST1H4A (K8) Antibody	CSB-PA010429PA08acHU	Human	ELISA, ICC, IF, ChIP
HIST1H4A	Acetyl-HIST1H4A (K8) Antibody	CSB-PA010429OA08acHU	Human	ELISA, ICC, IF



Target	Product Name	Code	Tested Species	Tested Applications
HIST1H4A	Crotonyl-HIST1H4A (K8) Antibody	CSB-PA010429OA08crHU	Human, Mouse, Rat	ELISA, WB, ICC, IF, ChIP
HIST1H4A	Butyryl-HIST1H4A (K8) Antibody	CSB-PA010429OA08butHU	Human, Mouse, Rat	ELISA, WB, ICC, IF, ChIP
HIST1H4A	Propionyl-HIST1H4A (K8) Antibody	CSB-PA010429OA08prophHU	Human	ELISA, WB, IF, IP, ChIP
HIST1H4A	2-hydroxyisobutyryl-HIST1H4A (K8) Antibody	CSB-PA010429OA08hibHU	Human	ELISA, WB, ICC, IF
HIST1H4A	Acetyl-HIST1H4A (K91) Antibody	CSB-PA010429OA91acHU	Human	ELISA, WB, IF
HIST1H4A	2-hydroxyisobutyryl-HIST1H4A (K91) Antibody	CSB-PA010429OA91hibHU	Human	ELISA, IF
HIST1H4A	Formyl-HIST1H4A (K91) Antibody	CSB-PA010429OA91forHU	Human	ELISA, ICC, IF
HIST1H4A	β-hydroxybutyryl-HIST1H4A (K91) Antibody	CSB-PA010429OA91bhbHU	Human, Rat	ELISA, WB
HIST1H4A	Mono-methyl-HIST1H4A (R3) Antibody	CSB-PA010429PA03me1HU	Human	ELISA, IF
HIST1H4A	Mono-methyl-HIST1H4A (R3) Antibody	CSB-PA010429OA03me1HU	Human	ELISA, WB, IF
HIST1H4A	Mono-methyl-HIST1H4A (R35) Antibody	CSB-PA010429OA35me1HU	Human	ELISA, IF
HIST1H4A	Mono-methyl-HIST1H4A (R67) Antibody	CSB-PA010429OA67me1HU	Human	ELISA, ICC
HIST1H4A	O-GlcNAcyl-HIST1H4A (S47) Antibody	CSB-PA010429OA47ogHU	Human, Rat	ELISA, WB
HIST1H4A	Hydroxyl-HIST1H4A (Y88) Antibody	CSB-PA010429OA88ohHU	Human, Rat	ELISA, WB, IF