

Recombinant Hepatitis C virus NS3 protein



Catalog Number: CSB-RP184674Ba

Product Name:	Recombinant Hepatitis C virus NS3 protein
Catalog Number:	CSB-RP184674Ba
Mol. Weight:	72kD
Source:	E.coli expressed
Purity:	>95%(SDS-PAGE)
Storage Buffer:	20mM Tris-HCl ,pH8.0,50% Glycerol
Storage :	Store at -20℃, for extended storage, conserve at -20℃ or -80℃.
Notes :	Repeated freezing and thawing is not recommended. Store working aliquots at 4℃ for up to one week.
AA sequence:	APITAYAQQTRGLLSTIVTSLTGRDTNENCGEVQVLSTATQSFLGTAVNGVMW TVYHGAGSKTISGPQGPVNQMYTNVDQDLVGWPAPPGVRS LAPCTCGASD LYLVTRHADVVPVRRRGDTRGALLSPRPIS TLKGSSGGPLLCPMGHAAGIFR AAVCTRGVAKAVDFVPVESLETTMRSPVFTDNSTPPAVPQTYQVAHLHAPTG SGKSTKVPAAYAAQGYKVLVLNPSVAATLGFGAYMSKAYGIDPNIRTGVRTITT GAPITYSTY GKFLADGGCSGGAYDIIICDECHSTDATTILGIGTALDQAETAGM RLVVLATATPPGSVTTPHSNIEEVALPTTGEVPFYGRAIPLGLIKGGRHLIFCH SKKKCDELAGQLTSLGLNAVAYYRGLDVSVIPTSGDVVVCATDALMTGFTGD FDSVIDCNTSVIQTIDFSLDPTFSIETTTVPQDAVSR SQRRGRTRGRGLGIYRY VTPGERPSGIFDTAVLCECYDAGCAWYELTPAETTVRLRAYFNTPGLPVCQD HLEFWESVFTGLTNIDGHFLSQTKQSGENFPYLVAYQATVCARALAPPPSWD TIWKCLTRLKPTLHGPTPLLYRLGSVQNEVT LTHPITKYIMACMSADLEVVT
References:	"Mechanisms of HCV Genotype 3-induced hepatocellular steatosis." Piodi A., Chouteau P., Barbotte L., Brillet R., Hezode C., Dhumeaux D., Pawlotsky J.M.