



Human HMGB1 cDNA Clone

Catalog Number: CSB-CL010553HU1

Synonyms	HMG3, SBP-1
Accession Number	BC0033781
cDNA Size	648
Definition	Homo sapiens high-mobility group box 1, mRNA (cDNA clone MGC:5223 IMAGE:2901382), complete cds.
recombination site	attL1, attL2
GenelD	3146
Vector	pENTR223.1 (spectinomycin)
Sequence	ATGGGCAAAGGAGATCCTAAGAAGCCGAGAGGCAAATGTCATCATATGCATTTTTTGTG CAAACCTTGTCCGGGAGGAGCATAAGAAGAAGCACCCAGATGCTTCAGTCAACTTCTCAG AGTTTTCTAAGAAGTGCTCAGAGAGGTGGAAGACCATGTCTGCTAAAGAGAAAGGAAAA TTTGAAGATATGGCAAAGGCCGACAAGGCCCGTTATGAAAGAGAAATGAAAACCTATATC CCTCCCAAAGGGGAGACAAAAAGAAGTTCAAGGATCCCAATGCACCCAAGAGGCCTC CTTCGGCCTTCTTCTTCTGCTCTGAGTATCGCCCAAAATCAAAGGAGAACATCCT GGCCTGTCCATTGGTGATGTTGCGAAGAACTGGGAGAGATGTGGAATAACACTGCTG CAGATGACAAGCAGCCTTATGAAAAGAAGGCTGCGAAGCTGAAGGAAAAATACGAAAA GGATATTGCTGCATATCGAGCTAAAGGAAAGCCTGATGCAGCAAAAAGGGAGTTGTCA AGGCTGAAAAAAGCAAGAAAAAGAAGGAAGAGGAGGAAGATGAGGAAGATGAAGAGG ATGAGGAGGAGGAGGAAGATGAAGAAGATGAAGATGAAGAAGAAGATGATGATGATGAA TAA