



Human TRPM1 cDNA Clone

Catalog Number: CSB-CL773791HU

Synonyms	LTRPC1
Accession Number	BC058286
cDNA Size	411
Definition	Homo sapiens transient receptor potential cation channel, subfamily M, member 1, mRNA (cDNA clone IMAGE:6673722), complete cds.
Insert site	rrnB_T1, rrnB_T2
GeneID	4308
Vector	pENTR223.1 (spectinomycin)
Sequence	ATGAAAGACTCTAACAGGTGTTGCTGTGGCCAGTTCACCAACCAGCATATCCCCCCTC TGCCAAGTGCAACACCCAGCAAAAATGAAGAGGAAAACAAACAGGTGGAGACTCAGC CTGAGAAATGGTCTGTTGCCAAGCACACCCAGAGCTACCCAACAGATTCCTATGGAGT TCTTGAATTCCAGGGTGGCGGATATTCCAATAAAGCCATGGTGAGAAAGGCATTCAGA CATGGTGCCACTAGGATCACAGCTTTCATTGGCGGCCAGTCTCCCAGCCCCAAACTG CAGATACCTGGTCTTCTTCATGGCTGTGGCTCAATCTTCTAGATATTTTCATTGAAAAAC CAAGAGATATATCTGTGCACATGGCTTTTAGCCATGAGGCTTGAAACTGGACACCAC TGTA