



## Human P2RY12 cDNA Clone

Catalog Number: CSB-CL861997HU

<b>Synonyms</b>	ADPG-R, HORK3, P2T(AC), P2Y12, P2Y(AC), P2Y(ADP), P2Y(cyc), SP1999
<b>Accession Number</b>	BC017898
<b>cDNA Size</b>	1029
<b>Definition</b>	Homo sapiens purinergic receptor P2Y, G-protein coupled, 12, mRNA (cDNA clone MGC:23802 IMAGE:4251263), complete cds.
<b>Insert site</b>	rrnB_T1, rrnB_T2
<b>GeneID</b>	64805
<b>Vector</b>	pENTR223.1 (spectinomycin)
<b>Sequence</b>	ATGCAAGCCGTCGACAATCTCACCTCTGCGCCTGGGAACACCAGTCTGTGCACCAGA GACTACAAAATCACCCAGGTCCTCTTCCCCTGCTCTACACTGCTCCTGTTTTTTGTTGG ACTTATCACAAATGGCCTGGCGATGAGGATTTTCTTTCAAATCCGGAGTAAATCAA TTATTATTTTTCTTAAGAACACAGTCATTTCTGATCTTCTCATGATTCTGACTTTCCATT CAAATTCTTAGTGATGCCAAACTGGGAACAGGACCACTGAGAACTTTTGTGTGTCAA GTTACCTCCGTCATATTTATTTTACAATGTATATCAGTATTTCACTCCTGGGACTGATAA CTATCGATCGCTACCAGAAGACCACCAGGCCATTTAAACATCCAACCCCAAAAATCTC TTGGGGGCTAAGATTCTCTCTGTTGTCATCTGGGCATTCATGTTCTTACTCTCTTTGCC TAACATGATTCTGACCAACAGGCAGCCGAGAGACAAGAATGTGAAGAAATGCTCTTTC CTAAATCAGAGTTCGGTCTAGTCTGGCATGAAATAGTAAATTACATCTGTCAAGTCATT TTCTGGATTAATTTCTTAATTGTTATTGTATGTTATACACTCATTACAAAAGAACTGTACC GGTCATACGTAAGAACGAGGGGTGTAGGTAAAGTCCCAGGAAAAGGTGAACGTCA AAGTTTTTCATTATCATTGCTGTATTCTTTATTTGTTTTGTTCCCTTTCCATTTTGCCCGAATT CCTTACACCCTGAGCCAAACCCGGGATGTCTTTGACTGCACTGCTGAAAATACTCTGT TCTATGTGAAAGAGAGCACTCTGTGGTTAACTTCCTTAAATGCATGCCTGGATCCGTTTC ATCTATTTTTTCTTTGCAAGTCCTTCAGAAATTCCTTGATAAGTATGCTGAAGTGCCCC AATTCTGCAACATCTCTGTCCCAGGACAATAGGAAAAAAGAACAGGATGGTGGTGACC CAAATGAAGAGACTCCAATGTAA