




Recombinant Lymphocyte antigen 6E(Ly6e) for Mus musculus (Mouse)

Catalog Number: CSB-EP717533MO

Product Name:	Recombinant Lymphocyte antigen 6E(Ly6e) for Mus musculus (Mouse)
Alternative names:	Stem cell antigen 2
Catalog Number:	CSB-EP717533MO
Relevance :	Involved in T-cell development.
Mol. Weight:	26kD
Product Info :	His-tagged
Source:	E.coli derived
Images	
Purity:	>90%(SDS-PAGE)
Storage Buffer:	20mM Tris-HCl, 0.5M Arg, PH 8.0,50% glycerol
Storage :	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Notes :	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
AA sequence:	LMCFSDTDQKNNINCLWVSCQEKDHYCITLSAAAGFGNVNLYTLNKGCSPICPSENVNLN LGVASVNSYCCQSSFCNFSA
References:	"Isolation of a cDNA encoding thymic shared antigen-1. A new member of the Ly6 family with a possible role in T cell development." Macneil I., Kennedy J., Godfrey D.I., Jenkins N.A., Masciantonio M., Mineo C., Gilbert D.J., Copeland N.G., Boyd R.L., Zlotnik A. J. Immunol. 151:6913-6923(1993)