

Mouse anti- Cortisol monoclonal antibody

Catalog Number: CSB-MA003711I0m



Synonym Names	Hydrocortisone
Product type	Primary antibodies
Description	Mouse monoclonal to Cortisol
Clonality	monoclonal
Isotype	IgG1
Purity	Caprylic Acid Ammonium Sulfate Precipitation purified
Conjugate	Non-conjugated
Storage buffer	Preservative: 0.03% Proclin 300 Constituents: 0.01M PBS, PH 7.4
Storage	Shipped at -20°C or -80°C. Avoid repeated freeze.
Form	Freeze-dried powder
Raised in	Mouse
Clone Number	003711
Tested applications	ELISA, Immune colloidal gold, Immunochromatography IC50: 39.8ppb Detect Range: 31.25-2000ppb